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Abbott

Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Sligo, Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Lisa Murphy, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Finisklin Business Park, Sligo, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Authorised Official

24/06/2019

Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Sligo, Ireland

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Signed: Lisa Murphy

Date: 24 Jun 2018

Lisa Murphy,
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division
Finisklin Business Park, Sligo, Ireland

Alinity i

CEA Controls

CEA

REF 07P6210

G71251R03

C7P620

Test Kit Package Changes: Revised September 2019

NAME

Alinity i CEA Controls (also referred to as CEA Ctrls)

INTENDED USE

The Alinity i CEA Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer (reagents, calibrators and instrument), when used for the quantitative measurement of Carcinoembryonic Antigen (CEA) in human serum or plasma.

For additional information, refer to the Alinity i CEA reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

The **CONTROL L**, **CONTROL M**, and **CONTROL H** contain CEA (human) prepared in phosphate buffer with protein (bovine) stabilizer. Preservatives: Antimicrobial Agents.

The following concentration ranges may be used for individual replicate control specifications on the Alinity i analyzer:

Control	Quantity	CEA (CONG) (ng/mL)	RANGE (ng/mL)
CONTROL L	1 x 8.0 mL	5	3.3 - 6.8
CONTROL M	1 x 8.0 mL	20	13.0 - 27.0
CONTROL H	1 x 8.0 mL	100	65.0 - 135.0

NOTE: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Control lot
- Reagent lot
- Calibrator lot
- Instrument

PRECAUTIONS

- **IPD**
- For *In Vitro* Diagnostic Use
- **RE ONLY**

Safety Precautions

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. The controls contain antigen derived from a human cell line. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that this product and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴

The following warnings and precautions apply to Disease Components: CONTROL L , CONTROL M , and CONTROL H		
WARNING		
H316	Causes mild irritation	
H332	Irritation on inhalation	
P303+P361+P531	In case of contact with skin, wash immediately with plenty of soap and water	

For additional information, refer to Section 12.2 of the OSHA Hazard Communication 29 CFR 1910.1200 (HCS) 2012 have been implemented.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically and dispense 4 drops of the low control, 4 drops of the medium control, and 4 drops of the high control into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, if calibration does not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria, or if controls do not meet the appropriate criteria.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization, *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI), *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
REF	List Number

Other Symbols	
CN	Control Number
CONC	Concentration
CONTROL L	Control Low, Medium, High (L,M,H)
DISTRIBUTED IN THE USA BY	Distributed in the USA by
INFORMATION FOR USA ONLY	Information provided for United States classification only
PRODUCT OF IRELAND	Product of Ireland
RANGE	Range
Rx ONLY	For use by or on the order of a physician only (applicable to USA classification only).

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott Ireland
Diagnostics Division
Finiskin Business Park
Sligo
Ireland
+353-71-9171712

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, Abbott Park, IL 60064, USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised September 2018.

©2016, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity I CEA Controls (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 26 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transfer and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

90/269/EEC, EN ISO 9001, ISO 13485, EN ISO 18113-1, EN ISO 18113-2

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbott diagnostics.com or contact your local representative.

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1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization, *Laboratory Biosafety Manual*, 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline—Fourth Edition*. CLSI Document M29-A4, Wayne, PA: CLSI; 2014.

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DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Diagnostic Division

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