Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Sligo, Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Lisa Murphy, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Finisklin Business Park, Sligo, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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SLIGO CHAMBER OF COMMERCE (INC)

CERTIFIED CORRECT IN ALL GOOD FAITH

Abbott

Abbott Ireland Diagnostics Division,
Finisklin Business Park,
Sligo, Ireland

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Read Highlig

Instructions
cannot be g
instructions
WARNING:
determined
due to diff
results rep
the identit
assay me
monitorin
levels se
carried o
baseline

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Signed: Lisa Murphy

Date: 24 Jun 2018

Lisa Murphy,
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division
Finisklin Business Park, Sligo, Ireland

■ **NAI**
Alinity i
■ **INT**
The Al
immun
Carcin
Alinity
The A
mana
CEA
■ **S**
Car
Fre
as
a f
rac
de
ar
e
T
c

Read Highlighted Changes: Revised September 2018.

REF 07P6220

REF 07P6230

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

WARNING: The concentration of CEA in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the CEA assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining CEA levels serially is changed, additional sequential testing should be carried out. Before changing assays, the laboratory MUST confirm baseline values for patients being serially monitored.

NAME

Alinity i CEA Reagent Kit

INTENDED USE

The Alinity i CEA assay is a chemiluminescent microparticle immunoassay (CMIA) used for the quantitative determination of Carcinoembryonic Antigen (CEA) in human serum and plasma on the Alinity i analyzer.

The Alinity i CEA assay is to be used as an aid in the prognosis and management of cancer patients in whom changing concentrations of CEA are observed.

SUMMARY AND EXPLANATION OF THE TEST

Carcinoembryonic antigen (CEA), first described in 1965 by Gold and Freedman,¹ is a tumor associated antigen. CEA was characterized as a glycoprotein of approximately 200 000 molecular weight with a β -electrophoretic mobility.^{2,3} Subsequent development of a radioimmunoassay (RIA) by Thomson, *et al*⁴ made it possible to detect the very low concentrations of CEA in blood, other body fluids, and also in normal and diseased tissues.⁵⁻⁷ Two years later, Hansen, *et al*⁸ developed a modified RIA for CEA.

The result of clinical studies to date indicate that CEA, although originally thought to be specific for digestive tract cancers, may also be elevated in other malignancies and in some nonmalignant disorders.⁹⁻¹⁵

CEA testing can have significant value in the monitoring of patients with diagnosed malignancies in whom changing concentrations of CEA are observed. A persistent elevation in circulating CEA following treatment is strongly indicative of occult metastatic and/or residual disease.¹⁶⁻²⁰

A persistently rising CEA value may be associated with progressive malignant disease and a poor therapeutic response.²¹⁻²³ A declining CEA value is generally indicative of a favorable prognosis and a good response to treatment.^{21, 23, 24} Patients who have low pretherapy CEA levels may later show elevations in the CEA level as an indication of progressive disease.²⁵

Clinical relevance of the CEA assay has been shown in the follow-up management of patients with colorectal, gastric, breast, lung, prostatic, pancreatic, and ovarian carcinoma.^{18, 24, 26-31}

Follow-up studies of patients with colorectal, breast, and lung carcinoma suggest that the preoperative CEA level has prognostic significance.³²⁻³⁵

CEA testing is not recommended as a screening procedure to detect cancer in the general population; however, use of the CEA test as an adjunctive test in predicting prognosis and as an aid in the management of cancer patients has been widely accepted.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay is a two-step immunoassay for the quantitative determination of CEA in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology. Sample and anti-CEA coated paramagnetic microparticles are combined and incubated. The CEA present in the sample binds to the anti-CEA coated microparticles. The mixture is washed. Anti-CEA acridinium-labeled conjugate is added to create a reaction mixture and incubated. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of CEA in the sample and the RLUs detected by the system optics.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity i CEA Reagent Kit 07P62

NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P6220	07P6230
Tests per cartridge	100	600
Number of cartridges per kit	2	2
Tests per kit	200	1200
MICROPARTICLES	6.6 mL	32.1 mL
CONJUGATE	6.1 mL	31.6 mL

MICROPARTICLES Anti-CEA (mouse, monoclonal) coated microparticles in TRIS buffer with protein (bovine) stabilizer. Minimum concentration: 0.1% solids. Preservative: antimicrobial agents.

CONJUGATE Anti-CEA (mouse, monoclonal) acridinium-labeled conjugate in phosphate buffer with protein (bovine) stabilizer. Minimum concentration: 0.8 μ g/mL. Preservative: antimicrobial agents.

Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.³⁶⁻³⁹

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
 - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.
- After mixing, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 1 hour before use.
Onboard	System Temperature	30 days	
Opened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity i CEA assay file must be installed on the Alinity i analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum
	Serum separator
Plasma	Sodium heparin
	Lithium heparin
	Potassium EDTA

- Plasma specimens collected in lithium or sodium heparin have been shown to exhibit an average of 7% to 8% higher results compared to corresponding serum results.
- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter or which are hazy or cloudy in appearance.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Centrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	2 to 8°C	7 days	If testing will be delayed more than 24 hours, serum or plasma should be removed from the clot, serum separator, or red blood cells.

If testing will be delayed more than 7 days, specimens should be stored frozen at -20°C or colder.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P62 Alinity i CEA Reagent Kit

Materials Required but not Provided

- Alinity i CEA assay file
- 07P6201 Alinity i CEA Calibrators
- 07P6210 Alinity i CEA Controls or other commercially available controls
- 09P1540 Alinity i Multi-Assay Manual Diluent
- Alinity Trigger Solution
- Alinity Pre-Trigger Solution
- Alinity i-series Concentrated Wash Buffer

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 10
 - Priority:
 - Sample volume for first test: 60 µL
 - Sample volume for each additional test from same sample cup: 10 µL
 - ≤ 3 hours on the reagent and sample manager:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 10 µL
 - > 3 hours on the reagent and sample manager:
 - Replace with a fresh aliquot of sample.
- Refer to the Alinity i CEA calibrator package insert and/or Alinity i CEA control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.

- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with a CEA value exceeding 1500.00 ng/mL are flagged with the code ("> 1500.00 ng/mL") and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

The system performs a 1:10 dilution of the sample and automatically calculates the concentration by multiplying the result by the dilution factor.

Manual Dilution Procedure

Suggested dilution: 1:100.

An additional 1:10 dilution may be made if needed. It is recommended that dilutions not exceed 1:1000.

Add 20 µL of the sample to 1980 µL of Alinity i Multi-Assay Manual Diluent.

The operator must enter the dilution factor in the Specimen or Control tab of the Create Order screen. The system will use this dilution factor to automatically calculate the concentration of the sample and report the result. The result should be > 4.00 ng/mL before the dilution factor is applied.

If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result. If a diluted sample result is less than 4.00 ng/mL, do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Each assay control must be tested to evaluate the assay calibration. Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of this package insert.
 - If statistically-based quality control limits are not available, then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

The recommended control requirement for the Alinity i CEA assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each clinically relevant control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (range) for the laboratory. Sources of variation that should be included in this study in order to be representative of future system performance include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules (if applicable)

- Data points collected at different times of the day

Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines, for general quality control recommendations.⁴⁰

- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity i-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.⁴¹

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity i-series Operations Manual.

RESULTS

Calculation

The Alinity i CEA assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration and results.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity i-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in ng/mL which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity i CEA assay is 1.73 to 1500.00 ng/mL.

LIMITATIONS OF THE PROCEDURE

- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits which employ mouse monoclonal antibodies.^{42, 43} Alinity i CEA reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.⁴⁴
- The Alinity i CEA assay should not be used as a cancer screening test.

- Patients with confirmed carcinoma frequently have a pretreatment CEA level in the same range as healthy individuals. Elevations in circulating CEA levels may be observed in smokers as well as patients with nonmalignant disease. For these reasons, a serum or plasma CEA level, regardless of value, should not be interpreted as absolute evidence for the presence or absence of malignant disease. The CEA level should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

EXPECTED VALUES

This study was performed on the ARCHITECT i System.

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

The distribution of ARCHITECT CEA values determined in 1141 specimens is shown in the following table.

Distribution of ARCHITECT CEA Values					
	Number of Subjects	Percent (%)			
		0 - 3 (ng/mL)	>3 - 5 (ng/mL)	>5 - 10 (ng/mL)	>10 (ng/mL)
<u>Healthy Subjects</u>					
Smokers	159	74.2	18.2	6.9	0.6
Non-smokers	149	83.2	11.4	5.4	0.0
Total	308	78.6	14.9	6.2	0.3
<u>Nonmalignant Disease</u>					
Ulcerative Colitis	50	72.0	20.0	4.0	4.0
Rectal Polyps	78	83.3	10.3	5.1	1.3
Pulmonary	60	61.7	20.0	13.3	5.0
Cirrhosis	110	47.3	30.0	15.5	7.3
Hepatitis	60	70.0	16.7	11.7	1.7
Renal	20	60.0	15.0	15.0	10.0
<u>Malignant Disease</u>					
Colorectal	150	24.0	10.7	10.7	54.7
Gastric	37	62.2	5.4	10.8	21.6
Pulmonary	110	47.3	19.1	9.1	24.5
Mammary	117	62.4	11.1	10.3	16.2
Ovarian	41	78.0	7.3	2.4	12.2

In this study, 93.5% of healthy subjects (n=308) had CEA values of 5.00 ng/mL or less.

The distribution table above for malignant disease is derived primarily from patients representing both active (clinical evidence of disease progression) and inactive (no clinical evidence of disease progression) disease states. When changing CEA assay methods in the course of monitoring a patient, additional sequential testing should be carried out to confirm baseline values.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity i analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2. Testing was conducted using 1 lot of the Alinity i CEA Reagent Kit, 1 lot of the Alinity i CEA Calibrators, 1 lot of the Alinity i CEA Controls and 1 instrument. Three buffer-based controls and 2 human serum panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.⁴⁵

by individuals.
ed in smokers
these
f value,
the presence
uld be used
evaluation

Sample	n	Mean (ng/mL)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Low Control	120	4.75	0.121	2.5	0.140	2.9
Medium Control	120	18.99	0.414	2.2	0.503	2.6
High Control	120	97.64	2.182	2.2	2.748	2.8
Panel 1	120	462.62	11.379	2.5	12.937	2.8
Panel 2	120	1266.38	31.595	2.5	37.701	3.0

^a Includes within-run, between-run, and between-day variability.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2. Testing was conducted using 3 lots of the Alinity i CEA Reagent Kit on each of 2 instruments over a minimum of 3 days. The maximum observed Limit of Blank (LoB), Limit of Detection (LoD), and Limit of Quantitation (LoQ) values are summarized below.⁴⁶

	ng/mL
LoB ^a	0.30
LoD ^b	0.42
LoQ ^c	1.73

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a total allowable error of 25% was met.

Linearity

A study was performed based on guidance from CLSI EP06-A.⁴⁷ This assay is linear across the measuring interval of 1.73 to 1500.00 ng/mL.

Interference

This study was performed on the ARCHITECT i System.

The specificity of the ARCHITECT CEA assay was determined by testing sera containing the compounds listed below. These compounds showed less than 10% interference in the CEA assay at the levels indicated.

Potentially Interfering Substance	Interferent Level
Bilirubin	≤ 20 mg/dL
Hemoglobin	≤ 500 mg/dL
Total Protein	≤ 12 g/dL
Triglycerides	≤ 3000 mg/dL

Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.⁴⁸

	Units	n	Correlation Coefficient	Intercept	Slope	Concentration Range
Alinity i CEA vs ARCHITECT CEA	Serum ng/mL	187	1.00	0.14	0.97	1.87 - 1372.22

Carryover

No detectable carryover (-0.02 ng/mL) was observed when a sample containing 42 372 ng/mL of CEA was assayed.

High Dose Hook

This study was performed on the ARCHITECT i System.

High dose hook is a phenomenon whereby very high level specimens may read within the dynamic range of the assay. For the CEA assay, no high dose hook effect was observed when samples containing up to approximately 60 000 ng/mL of CEA were assayed.

BIBLIOGRAPHY

- Gold P, Freedman SO. Demonstration of Tumor-Specific Antigens in Human Colonic Carcinomata by Immunological Tolerance and Absorption Techniques. *J Exp Med* 1964;121:439.
- Krupey J, Gold P, Freedman SO. Physicochemical Studies of the Carcinoembryonic Antigens of the Human Digestive System. *J Exp Med* 1968;183:387.
- Krupey J, Wilson T, Feedman SO, et al. The Preparation of Purified Carcinoembryonic Antigen of the Human Digestive System from Large Quantities of Tumor Tissue. *Immunochem* 1972;9(6):617-622.
- Thomson DMP, Krupey J, Freedman SO, et al. The Radioimmunoassay of Circulating Carcinoembryonic Antigen of the Human Digestive System. *Proc Natl Acad Sci USA* 1969;64:161-167.
- Zamcheck N. Carcinoembryonic Antigen: Quantitative Variations in Circulating Levels in Benign and Malignant Digestive Tract Diseases. *Adv Intern Med* 1974;19:413.
- Go VLM, Ammon HV, Holtermuller KH, et al. Quantification of Carcinoembryonic Antigen-Like Activities in Normal Human Gastrointestinal Secretions. *Cancer* 1975;36:2346-2350.
- Khoo SK, Warner NL, Lie JT, et al. Carcinoembryonic Antigenic Activity of Tissue Extracts: A Quantitative Study of Malignant and Benign Neoplasms. Cirrhotic Liver, Normal Adult and Fetal Organs. *Int J Cancer* 1973;11(3):681-687.
- Hansen HJ, Lance KP, Krupey J. Demonstration of an Ion Sensitive Antigenic Site on Carcinoembryonic Antigen Using Zirconyl Phosphate Gel. *Clin Res* 1971;19:143.
- Oehr P, Schlosser T, Adolphs HD. Applicability of an Enzymatic Test for the Determination of CEA in Serum and CEA-Like Products in Urine of Patients with Bladder Cancer. *Tumor Diagn* 1980;1:40.
- Ng WW, Tong KJ, Tam TN, et al. Clinical Values of CA 19-9, CA125, and CEA in Malignant Obstructive Jaundice. *Chung Hua I Hsueh Tsa Chih (Taipei)* 1995;55(6):438-446.
- Jerzersek B, Cervek J, Rudolf Z, et al. Clinical Evaluation of Potential Usefulness of CEA, CA 15-3, and MCA in Follow-up Breast Cancer Patients. *Cancer Letters* 1996;110(1-2):137-144.
- Wollenberg B, Jan V, Schmit UM, et al. CYFRA 21-1 is not Superior to SCC Antigen and CEA in Head and Neck Squamous Cell Cancer. *Anticancer Res* 1996;16(5b):3117-3124.
- Cerwenka H, Aigner R, Quehenberger F, et al. Preoperative Differential Diagnosis of Benign and Malignant Pancreatic Lesions - The Value of Pancreatic Secretory Trypsin Inhibitor, Procarboxypeptidase B, CA 19-9 and CEA. *Hepato-Gastroenterology* 1997;44(16):1117-1121.
- Pastor A, Menendez R, Cremades JM, et al. Diagnostic Value of SCC, CEA and CYFRA 21-1 in Lung Cancer: A Bayesian Analysis. *Eur Respir J* 1997;10(3):603-609.
- Maestranzi S, Przemioslo R, Mitchell H, et al. The Effect of Benign and Malignant Liver Disease on the Tumour Markers CA 19-9 and CEA. *Ann Clin Biochem* 1998;35:99-103.
- Duk JM, Aalders JG, Fluere GJ, et al. Tumor Markers CA 125, Squamous Cell Carcinoma Antigen, and Carcinoembryonic Antigen in Patients with Adenocarcinoma of the Uterine Cervix. *Obstet Gynecol* 1989;73(4):661-668.
- Munck-Wikland E, Kuylenstierna R, Lindholm T, et al. Carcinoembryonic Antigen, CA 19-9, and CA 50 in Monitoring Human Squamous Cell Carcinoma of the Esophagus. *Anticancer Res* 1990;10(3):703-708.
- Jager W, Kramer S, Palapelas V, et al. Breast Cancer and Clinical Utility of CA 15-3 and CEA. *Scand J Clin Lab Invest Suppl* 1995;221:87-92.
- Ebert W, Hoppe M, Muley T, et al. Monitoring of Therapy in Inoperable Lung Cancer Patients by Measurement of CYFRA 21-1, TPA-M, TPS, CEA, and NSE. *Anticancer Res* 1997;17(4B): 2875-2878.
- King J, Caplehorn JRM, Ross WB, et al. High Serum Carcinoembryonic Antigen Concentrations in Patients with Colorectal Liver Metastases is Associated with Poor Cell-Mediated Immunity, Which is Predictive of Survival. *Br J Surg* 1997;84:1382-1385.
- Noda M, Kusunoki M, Yanagi H, et al. Serum Carcinoembryonic Antigen (CEA) Correlates with the Survival Time During 5-FU Hepatic Arterial Infusion Chemotherapy for Unresectable Colorectal Hepatic Metastase - Clinical Study. *International Journal of Oncology* 1996;9(4):741-746.

22. Korenaga D, Saeki H, Mawatari K, et al. Serum Carcinoembryonic Antigen Concentration Doubling Time Correlates with Tumor Biology and Life Expectancy in Patients with Recurrent Gastrointestinal Carcinoma. *Archives of Surgery* 1997;132(2):188-194.
23. Nakayama T, Watanabe M, Teramoto T, et al. Slope Analysis of CA 19-9 and CEA for Predicting Recurrence in Colorectal Cancer Patients. *Anticancer Res* 1997;17(2b):1379-1382.
24. Lokich JJ, Zamcheck N, Lowenstein M. Sequential Carcinoembryonic Antigen Levels in the Therapy of Metastatic Breast Cancer. *Ann Intern Med* 1978;39:902.
25. Yasue M, Sakamoto J, Teramukai S, et al. Prognostic Values of Preoperative and Postoperative CEA and CA 19-9 Levels in Pancreatic Cancer. *Pancreas* 1994;9(6):735-740.
26. Zamcheck N, Martin EW. Factors Controlling the Circulating CEA Levels in Pancreatic Cancer. Some Clinical Correlations. *Cancer* 1981;47:1620.
27. Alsabti EAK, Kamel A. Carcinoembryonic Antigen (CEA) in Patients with Malignant and Non-Malignant Disease. *Neoplasma* 1979;26:603.
28. Khoo SK, Whitaker S, Jones I, et al. Predictive Value of Serial Carcinoembryonic Antigen Levels in Long-Term Follow-Up of Ovarian Cancer. *Cancer* 1979;43:2471.
29. Chevinsky AH. CEA in Tumors of Other than Colorectal Origin. *Semin Surg Oncol* 1991;7(3):162-166.
30. Diez M, Gomez A, Hernando F, et al. Serum CEA, CA125, and SCC Antigens and Tumor Recurrence in Resectable Non-Small Cell Lung Cancer. *Int J Biol Markers* 1995;10(1):5-10.
31. Bast RC, Bates S, Bredt AB, et al. Clinical Practice Guidelines for the Use of Tumor Markers in Breast and Colorectal Cancer. *J Clin Oncol* 1996;14(10):2843-2877.
32. Wanebo HJ, Rao B, Pinsky CM, et al. Preoperative Carcinoembryonic Antigen Level as a Prognostic Indicator in Colorectal Cancer. *N Engl J Med* 1978;299(9):448-451.
33. Concannon JP, Dalbow MH, Hodgson SE, et al. Prognostic Value of Preoperative Carcinoembryonic Antigen (CEA) Plasma Levels in Patients with Bronchogenic Carcinoma. *Cancer* 1978;42:1477-1483.
34. Ikeda Y, Mori M, Kajiyama K, et al. Indicative Value of Carcinoembryonic Antigen (CEA) for Liver Recurrence Following Curative Resection of Stage II and III Gastric Cancer. *Hepatogastroenterol* 1996;43(11):1281-1287.
35. Gebauer G, Muller-Ruchholtz W. Tumor Marker Concentrations in Normal and Malignant Tissues of Colorectal Cancer Patients and Their Prognostic Relevance. *Anticancer Res* 1997;17(4a):2731-2734.
36. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
37. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
38. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
39. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
40. Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures; Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
41. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
42. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
43. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
44. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
45. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
46. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

47. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
48. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.





Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

■ Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
	Serial number

Other Symbols

	Conjugate
	Distributed in the USA by
	Information needed for United States of America only
	Inversions Performed
	Microparticles
	Product of Ireland
	For use by or on the order of a physician only (applicable to USA classification only).

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Finisklin Business Park
Sligo
Ireland
+353-71-9171712

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised September 2018.

©2017, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i CEA Reagent Kit (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8°C. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 18 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, EN ISO 9001, ISO 13485, EN ISO 18113-1, EN ISO 18113-2

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

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- Gold P, Freedman SO. Demonstration of Tumor-Specific Antigens in Human Colonic Carcinomata by Immunological Tolerance and Absorption Techniques. *J Exp Med* 1964;121:439.
- Krupey J, Gold P, Freedman SO. Physicochemical Studies of the Carcinoembryonic Antigens of the Human Digestive System. *J Exp Med* 1968;183:387.
- Krupey J, Wilson T, Freedman SO, et al. The Preparation of Purified Carcinoembryonic Antigen of the Human Digestive System from Large Quantities of Tumor Tissue. *Immunochem* 1972;9(6):617-622.
- Thomson DMP, Krupey J, Freedman SO, et al. The Radioimmunoassay of Circulating Carcinoembryonic Antigen of the Human Digestive System. *Proc Natl Acad Sci USA* 1969;64:161-167.
- Zamcheck N. Carcinoembryonic Antigen: Quantitative Variations in Circulating Levels in Benign and Malignant Digestive Tract Diseases. *Adv Intern Med* 1974;19:413.
- Go VLM, Ammon HV, Holtermuller KH, et al. Quantification of Carcinoembryonic Antigen-Like Activities in Normal Human Gastrointestinal Secretions. *Cancer* 1975;36:2346-2350.
- Khoo SK, Warner NL, Lie JT, et al. Carcinoembryonic Antigenic Activity of Tissue Extracts: A Quantitative Study of Malignant and Benign Neoplasms. Cirrhotic Liver, Normal Adult and Fetal Organs. *Int J Cancer* 1973;11(3):681-687.
- Hansen HJ, Lance KP, Krupey J. Demonstration of an Ion Sensitive Antigenic Site on Carcinoembryonic Antigen Using Zirconyl Phosphate Gel. *Clin Res* 1971;19:143.
- Oehr P, Schlosser T, Adolphs HD. Applicability of an Enzymatic Test for the Determination of CEA in Serum and CEA-Like Products in Urine of Patients with Bladder Cancer. *Tumor Diagn* 1980;1:40.
- Ng WW, Tong KJ, Tam TN, et al. Clinical Values of CA 19-9, CA125, and CEA in Malignant Obstructive Jaundice. *Chung Hua I Hsueh Tsa Chih (Taipei)* 1995;55(6):438-446.

11. Jerzersek B, Cervek J, Rudolf Z, et al. Clinical Evaluation of Potential Usefulness of CEA, CA 15-3, and MCA in Follow-up Breast Cancer Patients. *Cancer Letters* 1996;110(1-2):137-144.
12. Wollenberg B, Jan V, Schmit UM, et al. CYFRA 21-1 is not Superior to SCC Antigen and CEA in Head and Neck Squamous Cell Cancer. *Anticancer Res* 1996;16(5b):3117-3124.
13. Cerwenka H, Aigner R, Quehenberger F, et al. Preoperative Differential Diagnosis of Benign and Malignant Pancreatic Lesions - The Value of Pancreatic Secretory Trypsin Inhibitor, Procarboxypeptidase B, CA 19-9 and CEA. *Hepato-Gastroenterology* 1997;44(16):1117-1121.
14. Pastor A, Menendez R, Cremades JM, et al. Diagnostic Value of SCC, CEA and CYFRA 21.1 in Lung Cancer: A Bayesian Analysis. *Eur Respir J* 1997;10(3):603-609.
15. Maestranzi S, Przemioslo R, Mitchell H, et al. The Effect of Benign and Malignant Liver Disease on the Tumour Markers CA 19-9 and CEA. *Ann Clin Biochem* 1998;35:99-103.
16. Duk JM, Aalders JG, Flueren GJ, et al. Tumor Markers CA 125, Squamous Cell Carcinoma Antigen, and Carcinoembryonic Antigen in Patients with Adenocarcinoma of the Uterine Cervix. *Obstet Gynecol* 1989;73(4):661-668.
17. Munck-Wikland E, Kuylensstierna R, Lindholm T, et al. Carcinoembryonic Antigen, CA 19-9, and CA 50 in Monitoring Human Squamous Cell Carcinoma of the Esophagus. *Anticancer Res* 1990;10(3):703-708.
18. Jager W, Kramer S, Palapelas V, et al. Breast Cancer and Clinical Utility of CA 15-3 and CEA. *Scand J Clin Lab Invest Suppl* 1995;221:87-92.
19. Ebert W, Hoppe M, Muley T, et al. Monitoring of Therapy in Inoperable Lung Cancer Patients by Measurement of CYFRA 21-1, TPA-M, TPS, CEA, and NSE. *Anticancer Res* 1997;17(4B): 2875-2878.
20. King J, Caplehorn JRM, Ross WB, et al. High Serum Carcinoembryonic Antigen Concentrations in Patients with Colorectal Liver Metastases is Associated with Poor Cell-Mediated Immunity, Which is Predictive of Survival. *Br J Surg* 1997;84:1382-1385.
21. Noda M, Kusunoki M, Yanagi H, et al. Serum Carcinoembryonic Antigen (CEA) Correlates with the Survival Time During 5-FU Hepatic Arterial Infusion Chemotherapy for Unresectable Colorectal Hepatic Metastase - Clinical Study. *International Journal of Oncology* 1996;9(4):741-746.
22. Korenaga D, Saeki H, Mawatari K, et al. Serum Carcinoembryonic Antigen Concentration Doubling Time Correlates with Tumor Biology and Life Expectancy in Patients with Recurrent Gastrointestinal Carcinoma. *Archives of Surgery* 1997;132(2):188-194.
23. Nakayama T, Watanabe M, Teramoto T, et al. Slope Analysis of CA 19-9 and CEA for Predicting Recurrence in Colorectal Cancer Patients. *Anticancer Res* 1997;17(2b):1379-1382.
24. Lokich JJ, Zamcheck N, Lowenstein M. Sequential Carcinoembryonic Antigen Levels in the Therapy of Metastatic Breast Cancer. *Ann Intern Med* 1978;39:902.
25. Yasue M, Sakamoto J, Teramukai S, et al. Prognostic Values of Preoperative and Postoperative CEA and CA 19-9 Levels in Pancreatic Cancer. *Pancreas* 1994;9(6):735-740.
26. Zamcheck N, Martin EW. Factors Controlling the Circulating CEA Levels in Pancreatic Cancer. Some Clinical Correlations. *Cancer* 1981;47:1620.
27. Alsabti EAK, Kamel A. Carcinoembryonic Antigen (CEA) in Patients with Malignant and Non-Malignant Disease. *Neoplasma* 1979;26:603.
28. Khoo SK, Whitaker S, Jones I, et al. Predictive Value of Serial Carcinoembryonic Antigen Levels in Long-Term Follow-Up of Ovarian Cancer. *Cancer* 1979;43:2471.
29. Chevinsky AH. CEA in Tumors of Other than Colorectal Origin. *Semin Surg Oncol* 1991;7(3):162-166.
30. Diez M, Gomez A, Hernandez F, et al. Serum CEA, CA125, and SCC Antigens and Tumor Recurrence in Resectable Non-Small Cell Lung Cancer. *Int J Biol Markers* 1995;10(1):5-10.
31. Bast RC, Bates S, Bredt AB, et al. Clinical Practice Guidelines for the Use of Tumor Markers in Breast and Colorectal Cancer. *J Clin Oncol* 1996;14(10):2843-2877.
32. Wanebo HJ, Rao B, Pinsky CM, et al. Preoperative Carcinoembryonic Antigen Level as a Prognostic Indicator in Colorectal Cancer. *N Engl J Med* 1978;299(9):448-451.
33. Concannon JP, Dalbow MH, Hodgson SE, et al. Prognostic Value of Preoperative Carcinoembryonic Antigen (CEA) Plasma Levels in Patients with Bronchogenic Carcinoma. *Cancer* 1978;42:1477-1483.
34. Ikeda Y, Mori M, Kajiyama K, et al. Indicative Value of Carcinoembryonic Antigen (CEA) for Liver Recurrence Following Curative Resection of Stage II and III Gastric Cancer. *Hepatogastroenterol* 1996;43(11):1281-1287.
35. Gebauer G, Muller-Ruchholtz W. Tumor Marker Concentrations in Normal and Malignant Tissues of Colorectal Cancer Patients and Their Prognostic Relevance. *Anticancer Res* 1997;17(4a):2731-2734.
36. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
37. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
38. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
39. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
40. Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
41. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
42. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
43. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
44. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
45. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
46. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
47. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
48. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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