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SPECIFIC PERFORMANCE CHARACTERISTICS

Each cell sample is shown to have a negative direct agglutinin test.

When properly stored and used according to standard procedures, these reagents will demonstrate the appropriate agglutinins / antibodies specified in the reagent description.

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This reagent has been validated for use on the TANGO[®] optimo and TANGO infinity[®] automated blood grouping test system, for use with Erytype and Erytype S as well as for the manual tube technique. Its suitability for use by other technologies cannot be guaranteed. Users must confirm reagent suitability before using alternative techniques.

BIOGRAPHY

Technical Manual, 18th ed. Bethesda, MD: American Association of Blood Banks, 2014.

Standards for Blood Banks and Transfusion Services, 29th ed. Bethesda, MD: American Association of Blood Banks, 2014.

DATE OF ISSUE

February 2015

For additional information please contact your local distributor.

Distributor

Alba Bioscience Limited
Ellen's Glen Road
Edinburgh
Scotland, UK
EH17 7QT

Distributor

Bio-Rad Medical Diagnostics GmbH
Industriestr. 1
D-63303 Dreieich, Germany

Tel:

+49 6103 3130 0

Fax:

+49 6103 3130 724

E-Mail:

techsupport.bmd@bio-rad.com

www.bio-rad.com

Alba Bioscience Limited
doing business as

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BIBLIOGRAFIA

¹ Technical Manual. 18^a ed. Bethesda, MD. American Association of Blood Banks, 2014.

Standards for Blood Banks and Transfusion Services. 25^a ed. Bethesda, MD. American Association of Blood Banks, 2014.

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Control Set QC

Simulated Whole Blood Controls

Use as controls of blood grouping reagents

IN VITRO DIAGNOSTIC USE

PACKAGE SIZE

816034 VOL 4 x 5 ml Control Set QC

SUMMARY AND EXPLANATION

The purpose of daily quality assurance in the blood bank is to confirm the reliability of the test system. The test system includes reagents, test procedures and equipment. Testing known samples is an accepted method of quality control. If expected test results are observed, procedures are being performed accurately and reagents and equipment are performing properly. If unexpected results are observed, the problem may be due to improper test performance, faulty equipment, contamination or deterioration of reagents. The source of the problem should be determined and resolved before patient test results are reported.

Control Set QC provides a means of confirming the reactivity of routinely used reagents and are to be tested on each day of use. Observation of expected test results with Control Set QC will confirm the reactivity of anti-A, anti-B, anti-AB, anti-D, anti-C, anti-e, anti-K and anti-K.

INTERPRETATION OF LABELING SYMBOLS

Batch code

Use by (YYYY-MM-DD)

Storage temperature limitation (2...8°C)

In vitro diagnostic medical device

Consult instructions for use

Manufacturer

INTENDED USE

Control Set QC is intended for use as a control for blood grouping reagents used in the TANGO® optimo and TANGO infinity® blood grouping system, for Erytype and Erytype S as well as for manual tube tests.

PRINCIPLE OF PROCEDURE

The procedures used with these reagents are based on the principle of agglutination. Normal human red blood cells will agglutinate in the presence of antibody directed against antigens on those red blood cells. No agglutination indicates the absence of the demonstrable antigen or antibody.

The simulated whole blood controls provided in Control Set QC confirm the reactivity of the reagents used for ABO, Rh and K determinations, as well as reverse grouping cells.

REAGENT DESCRIPTION

Control Set QC is prepared from red blood cells collected from blood donors. Each individual donation contains the appropriate ABO, Rh and K blood group antigens and also the appropriate ABO blood group antibodies. ABO antibodies are of monoclonal origin.

The concentration of red blood cells in each of the controls has been adjusted to 20-25%. The red blood cells are suspended in a preservative solution to retard haemolysis and bacterial contamination.

Vial	Blood group								Antibody
	ABO	RH	RH1 D	RH2 C	RH3 E	RH4 c	RH5 e	KEL1 K	
1	A1	R1R1	+	+	-	-	+	+	Anti-B
2	B	R1R2	+	+	+	+	+	-	Anti-A
3	O	R2R2	+	-	+	+	-	-	Anti-A, Anti-B
4	AB	rr	-	-	-	+	+	-	-

STORAGE CONDITIONS

Do not transfer these reagent red cells to another container as this could result in spillage or contamination. Store at 2...8°C. Do not use if obviously discolored or haemolysed. Do not use beyond the notified expiry date. Slight discoloration in the supernatant is normal. After opening the vial the product can be stored under proper storage conditions (2...8°C) until the expiry date.

PRECAUTIONS FOR USE AND DISPOSAL

The preservative solution has been specially formulated to preserve red cell integrity and antigenicity and contains the following components - trisodium citrate, citric acid, dextrose, inosine, neomycin sulphate (0.103g/L) and chloramphenicol (0.349g/L).

Chloramphenicol is classified as a carcinogen and neomycin sulphate is classified as an irritant.

CAUTION: Each donor unit was tested in accordance with current European and FDA required tests and found to be reactive for HBsAg, Anti-HIV1/2, Anti-HCV and syphilis.

As this reagent is of human origin care must be taken during use and disposal as there is a potential infection. Known test methods can offer assurance that products derived from human blood will not transmit infectious agents. This reagent is for in vitro diagnostic use only.

TEST PROCEDURES

General information

When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer.

Control Set QC is intended to simulate normal blood samples. The samples contained in Control Set QC must be kept at room temperature (18...25°C) and should be tested by following standard procedures in accordance with the Directia Use accompanying each reagent used routinely.

Vials require to be centrifuged prior to first use. Centrifugation must be performed according to the users standard laboratory practice for the centrifugation of patient samples. Subsequent upright, refrigerated storage eliminates the requirement for further centrifugation unless mixing has occurred.

Additional Material and Reagents Required

For automated, semi-automated or manual testing:

- TANGO® optimo or TANGO infinity®
- Erytypecell-A1, -A2, -B, -0 [REF] 816006
- Bromelain Erytype™ ready for use [REF] 806220
- Isotonic saline
- Pipettes (drop volume 40-50 µl)
- Shaker e.g. [REF] 845 225
- Photometer for colorimetric and agglutination reading e.g. [REF] 841 010
- Erytype™ or Erytype™ S
- Erytypecell A, & B [REF] 810056100
- Centrifuge
- Glass tubes
- Reading mirror e.g. REF 845 310

For manual tube tests:

- Glass test tubes
- Timer
- Pipettes
- Centrifuge

QUALITY CONTROL

This is a quality control reagent and its satisfactory performance when used by the recommended techniques requires an adequate level of control.

Quality control of testing must be performed according to the local state, federal and accreditation requirements.

INTERPRETATION OF RESULTS

The following table illustrates the expected results in tests with Control Set QC and routine blood bank reagents

Component of Biotech Control Set QC	Reagent Under Test	Expected Test Results*	Reagent Under Test	Expected Test Results*
Vial 1	Anti-A	+	Anti-D	+
	Anti-B	0	Anti-C	+
	Anti-AB	+	Anti-E	0
	A ₁ cells	0	Anti-c	0
	A ₂ cells	0	Anti-e	+
	B cells	+	Anti-K	+
	O cells	0		
Vial 2	Anti-A	0	Anti-D	+
	Anti-B	+	Anti-C	+
	Anti-AB	+	Anti-E	+
	A ₁ cells	+	Anti-c	+
	A ₂ cells	+	Anti-e	+
	B cells	0	Anti-K	0
	O cells	0		
Vial 3	Anti-A	0	Anti-D	+
	Anti-B	0	Anti-C	0
	Anti-AB	0	Anti-E	+
	A ₁ cells	+	Anti-c	0
	A ₂ cells	+	Anti-e	+
	B cells	+	Anti-K	0
	O cells	0		
Vial 4	Anti-A	+	Anti-D	0
	Anti-B	+	Anti-C	0
	Anti-AB	+	Anti-E	0
	A ₁ cells	0	Anti-c	+
	A ₂ cells	0	Anti-e	+
	B cells	0	Anti-K	- 0
	O cells	0		

*Discrepant results must be investigated further.

PERFORMANCE LIMITATIONS

Improper techniques may invalidate the results obtained with this product.

False positive or false negative results can occur due to contamination of test materials, improper reaction temperature, improper storage of materials and omission of test reagents.

Individual laboratory procedures may affect the final reaction strength observed in tests performed with Control Set QC.

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Appendix to Instructions for Use of in vitro diagnostic medical device (for use only in Russian Federation)

Medical Device Name

Control Set QC for immunohematology tests, 4 x 5 ml/1 pkg.

Note:

The following name of the product may be used in the text - Control Set QC.

Devices used together with the medical device	For testing by automatic method • TANGO infinity automatic immunohematology analyzer, with accessories ([REF] 850 000 010) • Erytypecell-A₁, -A₂, -B, -0 reagent red blood cells for reverse typing in a set, 4 x 10 ml/1pkg. ([REF] 816 006) - Erytypecell-A₁ reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 041) - Erytypecell-B reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 043) - Erytypecell-0 reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 044) • Bromelin Erytype ready for use solution for enzyme treatment of red blood cells, 10 x 50 ml/1 pkg. ([REF] 806 220) - Erytype S AB0+D microtestplates (single strips) for AB0- and Rh-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806105) - Erytype S ABD Confirmation microtestplates (single strips) for AB0- and Rh-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806110) - Erytype S ABD+Rev. A1, B microtestplates (single strips) for AB0- and D-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806127) - Erytype S Rh+K single microtestplates (single strips) for Rh- and Kell-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806170) - Erytype S Rh+K microtestplates (single strips) for Rh- and Kell-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806181) For testing by semi-automatic / manual method: - Erytypecell-A₁, -A₂, -B, -0 reagent red blood cells for reverse typing in a set, 4 x 10 ml/1pkg. ([REF] 816 006) - Erytypecell-A₁ reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 041) - Erytypecell-B reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 043) - Erytypecell-0 reagent red blood cells for reverse typing, 10 ml/1 vial. ([REF] 816 044) - Bromelin Erytype ready for use solution for enzyme treatment of red blood cells, 10 x 50 ml/1 pkg. ([REF] 806 220) - Erytype S AB0+D microtestplates (single strips) for AB0- and Rh-typing
---	---

	of red blood cells, 10 pcs./1 pkg. ([REF] 806105) - Erytype S ABD Confirmation microtestplates (single strips) for AB0- and Rh-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806110) - Erytype S ABD+Rev. A1, B microtestplates (single strips) for AB0- and D-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806127) - Erytype S Rh+K single microtestplates (single strips) for Rh- and Kell-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806170) - Erytype S Rh+K microtestplates (single strips) for Rh- and Kell-typing of red blood cells, 10 pcs./1 pkg. ([REF] 806181) • Isotonic saline solution • Pipettes (drop volume 40-50µl) • Shaker • Photometer for colorimetric and agglutination reading • Centrifuge • Glass tubes • Reading mirror
--	--

Intended Use

Control Set QC for immunohematology tests, 4 x 5 ml/1 pkg. is intended for use as a control for blood grouping reagents used in the TANGO infinity blood grouping system, for Erytype S.

Intended User

Trained qualified laboratory professionals.

Physical Properties

Appearance - reddish, turbid liquid with sedimentation

Device Chemical Composition

The product contains material of human origin.





Component	Amount
Dextrose	2 % v/v
Trisodium citrate	0.8 % v/v
Sodium chloride	0.42 % v/v
Citric acid	0.05 % v/v
Inosine	0.09 % v/v
Chloramphenicol	0.349g/L
Neomycin sulphate	0.103g/L
Methylated spirits (95%)	0.2 % v/v

Human red blood cells	20- 25 % v/v
Water	71-76 % v/v
Monoclonal antibody	Variable*

*Variable depending antibodies specificity and potency

Each donor unit underwent approved method testing to demonstrate non reactivity for HBsAg, Anti-HIV1/2, Anti-HCV and syphilis.

Description of symbols

Symbol	Description
	Manufacturer
	Distributor logo
	The device complies with the Directive 98/79/EC on in-vitro diagnostic medical devices
	<i>In vitro</i> diagnostic medical device
	Refer to instructions for use
	Storage temperature range
	Batch code
	Use by
	reference number
	product packaging is able to be recycled
	2D barcode Data Matrix
	UDI: (01) - Device Identifier (17) –Shelf life (10) – batch code

Shipping Content

Shipping package includes 4 tubes with control material and instructions for use.

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Control Set QC is intended to simulate normal blood samples. The samples contained in Control Set QC must be used at room temperature (18...25°C) and should be tested by following standard procedures in accordance with the Directions for Use accompanying each reagent used routinely.

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Analytical Characteristics

The product allows performing reproducible control of reagents.

Inter- and intra-lot reproducibility of the product has been evaluated during the internal study. All test results met the acceptance criteria with all positive reactions achieving the required grades and all negative tests giving a clear unequivocal result.

Specific Performance Characteristics

Each cell sample is shown to have a negative direct antiglobulin test.

When properly stored and used according to standard procedures, these reagents will demonstrate the appropriate antigens / antibodies specified in the reagent description.

The Procedure and Interpretation of Results must be followed closely to ensure the accuracy of the test results. Each laboratory should have a program that will train personnel on the proper use and handling of the product.

This reagent has been validated for use on the TANGO infinity automated blood grouping test system, for use with Erytype S. Its suitability for use by other technologies cannot be guaranteed. Users must confirm reagent suitability before using alternative techniques.

Metrological Traceability

Not applicable.

Shipping

The product must be transported in the original package that protects products from external impacts in all kinds of vehicles in accordance with the shipping rules applicable to the transport of this type. Shipping temperature: 2-8 °C.

Loading and unloading and transportation should strictly comply with requirements of handling marks printed on the packaging labels and on the transport packaging labels.

Do not use this medical device if transport and storage conditions are not met.

Storage

Unopened vial should be stored at 2-8 °C in refrigerated chambers or in refrigerators till expiry date printed on the label.

Opened vial should be stored at 2-8 °C in refrigerated chambers or in refrigerators till expiry date printed on the label. (For manual/semiautomatic method only).

Do not use this medical device if storage conditions are not met.

Disposal

This product must be handled and discarded in accordance with all local, regional and national regulations. All parts of this medical device, any samples and materials used and produced in the performance of tests, including expired reagents should be disposed as potentially epidemiologically dangerous waste (class B).

Warranty

The manufacturer warrants constant product characteristics stated in the documentation when the product is used in accordance with the IFU and for its intended use specified by the manufacturer.

For all questions and support please contact manufacturer's authorized representative on the territory of the Russian Federation:

ООО "Bio-Rad Laboratorii"

Address: 105064, Moscow, Nizhnij Susalnyj pereulok, 5, building 5A.

Phone: +7 (495) 721-14-04.

Lt Reinise
28 June 2018

Lindsey Helen Robinson
Solicitor and Notary Public

DWF LLP
No 2 Lochrin Square
96 Fountainbridge
Edinburgh
EH3 9QA
OX ED15

Alba Bioscience Limited
doing business as
Quotient

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