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Bimeo PRO

User manual

Approved

Aleš Hribar
CEO

bimeo

Kinestica d.o.o. has patents or pending patent applications covering the subject matter contained in this document. The furnishing of this document does not give you any license to such patents. This device is protected under the patents.

No part of this manual or Bimeo PRO software may be reproduced in any form or by any means, including, without limitation, electronic or mechanical, such as photocopying or recording, or by any information storage and retrieval systems, without the express written consent of Kinestica d.o.o. Specifications are subject to change without notice.

Copyright © 2019 Kinestica d.o.o. All rights reserved.

Version: 003

Version date: 21.11.2019

Contact information

Kinestica d.o.o.
Podmilščakova 46,
1000 Ljubljana
Slovenia
info@kinestica.com
www.kinestica.com

Important safety instructions

Ignoring these warnings may cause severe injury or death as a result of incorrect operation.

Ignoring these cautions may cause moderate injury or property damage as a result of incorrect operation.

	WARNING
	Before using this product, please read the warnings and safety instructions manual. For safe and correct use of Bimeo PRO device, read this manual before use. Keep the manual and safety information in a convenient place for future reference.

	ATTENTION
	- Never leave the system unattended while charging. - There are hazardous voltages inside the Bimeo charger. The casing provides for safety and must not be opened or removed. - Use only original accessories from the manufacturer.

	WARNING
	- In case of heating up during charging, immediately disconnect the charger from the device. - Only use provided charger to charge the Bimeo PRO. - Do not disassemble the Bimeo PRO system. No serviceable parts are to be found inside the device. Maintenance must be carried out by a professionally skilled person only. - Do not dispose of system or its batteries in a fire or trash incinerator, or leave it in hot places. - Do not immerse the battery in water or seawater and keep the system in a cool dry environment during stands by period. - Do not use it in a location where static electricity is rich. - In case the electrolytes get into the eyes due to the leakage of battery, do not rub the eyes! Rinse the eyes with clean running water and seek medical attention immediately. Otherwise, it may injure eyes or cause loss of sight.

	WARNING
	The surroundings of the patient workspace should NOT CONTAIN FERROMAGNETIC components in order not to deteriorate the performance of the sensing devices.

	WARNING
	Do not place the sensing devices (Bimeo master, Bimeo sense) close to PERMANENT MAGNETS or devices that could produce HIGH MAGNETIC FIELD in order not to magnetize the sensing devices. This also includes laptops, loudspeakers and fMRI.

	ATTENTION
	In the event of exposure to HIGH MAGNETIC FIELD, the device must be calibrated. Contact your dealer or the manufacturer about calibrating the device.

	ATTENTION
	The patient should not be touching the bottom surface of the Bimeo PAD with the arm. The patient should not wear long sleeve shirt. The sleeves could also disturb the position sensing of the Bimeo PAD. When turning on the Bimeo PAD do not place any objects onto the bottom surface of the Bimeo PAD sensor. The Bimeo disc or Bimeo double disc can be placed onto the sensor when the status light starts blinking.

Safety and maintenance recommendations

Important labels

	CE MARKING
	TYPE OF PROTECTION AGAINST ELECTRIC SHOCK ACCORDING TO EN 60601-1: CLASS II
	TYPE BF APPLIED PART
	NON-IONIZING RADIATION
	WEEE MARKING
	SEE THE INSTRUCTIONS FOR USE
S/N:	SERIAL NUMBER

EMC statements

Bimeo PRO needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in accompanying EMC declaration. Portable and mobile RF communications equipment can affect the Bimeo PRO device.

Therapist responsibility

Therapist are responsible to warn technicians or others if the device does not work correctly or they suspect the device is not working as intended. The therapist must assess whether or how long the patient should perform therapy. Make sure that the straps are not too tight or they cause skin irritation. Straps are consumables and can be replaced if they are damaged or worn out.

Therapy duration

When setting the duration of the task or setting the number of tasks in therapy program consider the patient's ability to avoid excessive exhaustion of the patient.

	ATTENTION
	Set the therapy duration according to the patient's ability to avoid excessive exhaustion of the patient. When using the Bimeo PRO system, it is the operator's responsibility to:

	- take the time to properly instruct all patients especially those who have difficulty following and carrying out instructions. It is highly recommended that such patients be continually coached and verbally guided throughout the session, - set the appropriate range-of-motion according to the patient's ability, - closely monitor all patients using the system.
--	---

Sound level setting

Bimeo PRO software includes audio sounds that are presented to the patient via system's audio device. Set the volume control at a low setting. Slowly increase the sound until you can hear it comfortably and clearly. Once you have established a comfortable sound level, set it and leave it there. Refer to the user manual of the sound device for the detailed instructions on setting the sound level.

	WARNING
	- Set the volume level of the audio device to a minimum, and then adjust the volume gradually. - Sudden exposure to loud noises could cause hearing damage.

Cleaning the Bimeo PRO system

It is recommended to clean and disinfect the Bimeo PRO device regularly with 70% or less Ethyl Alcohol solution. Do not use aggressive cleaning agents, volatile liquids or excessive force when cleaning the device. Do not wash the device or immerse it in water. Also, do not use petrol, thinners or similar solvents to clean the device. Regular cleaning of the device prevents cross-contamination between patients or cross-contamination in a single patient.

	ATTENTION
	Clean the device on regular basis with 70% or less solution Ethyl alcohol.

	WARNING
	- Do not use aggressive cleaning agents, volatile liquids or excessive force when cleaning the device. - Do not wash the device or immerse it in water.

Service and maintenance intervals

It is recommended to calibrate arm sensors of the Bimeo PRO device once per year. Contact your dealer or the manufacturer about calibrating the device. Bimeo PAD device does not need to be calibrated during the product life time. The device can be serviced or calibrated only by manufacturer or service personnel certified by manufacturer.

	WARNING
	- Do not open the device. The device does not contain any parts which can be replaced by a user. Do not alter or adapt the device.

Product life and storage

If correctly used and maintained, the device will have a minimum service life of 5 years. You can store the device in suitable conditions for a maximum of 2 years. The battery should be half-empty during the storage. When using the device after storage, we recommend you to subject the device to a thorough maintenance check and calibrate it.

	ATTENTION
	Do not dispose of the Bimeo PRO device in domestic waste. Always dispose of electronic appliances at the designated collection facilities.

Warranty information

A warranty period of one year applies to the device, starting from the date of putting into operation at the end user site. Products under warranty claim should not show signs of deliberate damage. Products should be operated within normal use conditions specified in the User's manual.

Guidance and manufacturer's declaration

Special precautions concerning electromagnetic compatibility (EMC) must be taken for all medical electrical equipment.

- All medical electrical equipment must be installed and put into service in accordance with the EMC information provided in this document.
- Portable and mobile RF communications equipment can affect the behavior of medical electrical equipment.

The Bimeo PRO system complies with all applicable and required standards for electromagnetic interference.

- They do not normally affect nearby equipment and devices.
- They are not normally affected by nearby equipment and devices.

The Bimeo PRO device is intended for use in the electromagnetic environment specified below. The customer or the user of the Bimeo PRO device should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 2	The Bimeo PRO device must emit electromagnetic energy to perform its intended function. Nearby electronic equipment can be affected.
RF emissions CISPR 11	Class B	The Bimeo PRO device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

The Bimeo PRO device is intended for use in the electromagnetic environment specified below. The customer or the user of the Bimeo PRO device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Electromagnetic compatibility (EMC) of the Bimeo PRO device

The Bimeo PRO device is intended for use in the electromagnetic environment specified below. The customer or the user of the Bimeo PRO device should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Bimeo PRO device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Conducted RF IEC 61000-4-6	Not applicable	Not applicable	Not applicable
Radiated RF IEC 61000-4-3	3 Vrms 80 MHz to 2.5 GHz	3 Vrms	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts according to the transmitter manufacturer and d is the recommended separation distance in meters. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Bimeo PRO device is used exceeds the applicable RF compliance level above, the Bimeo PRO device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Bimeo PRO device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

The Bimeo PRO device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Bimeo PRO device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Bimeo PRO device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power (meters) of transmitter	Separation distance (meters) according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (w) according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Medical device classification MDD	Class I (Annex IX, Rule 12)
Battery	Battery type: rechargeable lithium polymer Total capacity, arm sensors: 880 mAh Total capacity, PAD sensor: 2000 mAh Battery autonomy: 6 hours
Charger	AC/DC power supply: • DUE 15WCP1-5.0, UE15WCP1-050 Input: 90-240 V AC/50-60 Hz/500 mA Output: 5 V DC/2.0 A max • or MW GEM12I05 Input: 100-240 V AC/50-60 Hz/400 mA Output: 5 V DC/2.4 A
Protection type	Compliant with standards: • EN 60601-1:2006/A12:2014 General requirements for basic safety and essential performance, • EN 60601-1-2:2007 Electromagnetic compatibility - Requirements and tests,
Temperature and humidity range	Working environment: 10 °C to 35 °C, 30 % to 85 % relative air humidity, IPX0 protection, Transport: -20 °C to 60 °C, up to 85 % relative air humidity Storage: 0 °C to 30 °C, up to 85 % relative air humidity
Communication	Wireless 2.400 GHz – 2.500 GHz, max. 5 mW
Force sensor	Three-axes, Range: 100 N, non-measuring function
Inertial measurement sensors	Accelerometer, gyroscope, magnetometer, non-measuring function
Bimeo PAD measurement principle	IR position detection, non-measuring function

Processor	Dual Core 2.0 GHz or equivalent processor
Memory	2 GB RAM
Disk space	1 GB free hard disk space
Graphic	DirectX 9.0c compliant video card with 128 MB of RAM
Screen resolution	1280x720
Operating system	Microsoft Windows 7 SP1/8/8.1/10

Processor	Quad-core Intel i5 or equivalent AMD CPU
Memory	4 GB RAM
Disk space	1 GB free hard disk space
Graphic	DirectX 9.0c compatible NVIDIA or AMD ATI video card with 512 MB of RAM
Screen resolution	1920x1080
Operating system	Microsoft Windows 8/8.1/10 Professional

The System for recovery of motion activity, motion coordination of upper limbs using interactive programs Bimeo PRO, also referred as Bimeo PRO System, Bimeo PRO, addresses a number of important issues in rehabilitation: it can be used in clinical environments enabling intensive therapy of prolonged duration, it provides objective assessment of motor functions and it includes a number of features, which are common for robot supported rehabilitation. Bimeo PRO is a sensor based device designed for unimanual and bimanual therapy. The Bimeo PRO therapy approach exploits the preserved motor functions of the affected or unaffected limb of the patient, thus it does not require any external actuators. The less affected limb can directly support the movement of the affected limb, thus taking over the function of a robot manipulator. Therefore, the system addresses also the problem of intention detection. The manipulator that moves the affected arm is actually the patient's limb (affected or unaffected), which is under the patient's voluntary control. Thus, the intended movement is always known and correctly supported, while the support is evaluated and fed back to the patient through various modalities in order to optimize the therapy. The Bimeo PRO concept can be combined with most other therapeutic approaches (EMG based biofeedback, electrical stimulation, constraint-induced movement therapy, repetitive task training, virtual reality, robot and other). Sensors embedded into the system enable objective functions assessment (movement trajectory, active or passive range of motion, maximal forces, support of the unaffected limb during rehabilitation and all other derived quantities) and provide simultaneous feedback to the patient and the therapist.

Bimanual therapy differs from unimanual and this can be efficiently exploited for implementation of new methods and devices for therapy. During the bimanual therapy interaction occurs between the unaffected limb, which functions are usually well preserved, and the affected limb, which functions are (partially) lost and need to be regained. In this regard bimanual therapy can be considered as unimanual therapy with a robot where the role of the robot is replaced by the unaffected limb function. The unaffected limb thus assumes the role of the manipulator. An important feature of robots is their ability to objectively measure functions of the affected limb, which the unaffected limb in the role of the manipulator does not allow. Therefore, bimanual therapy needs to be augmented with measurement, analysis and decision-making capabilities.

The Bimeo PAD device is an accessory to Bimeo PRO device. Bimeo PAD is a sensor based device designed for unimanual and bimanual therapy on surface and offers better position indication than Bimeo PRO system itself.

Important labels.....	5
EMC statements.....	5
Therapist responsibility.....	5
Therapy duration.....	5
Sound level setting.....	6
Cleaning the Bimeo PRO system.....	6
Service and maintenance intervals.....	6
Product life and storage.....	7
Warranty information.....	7
Guidance and manufacturer's declaration.....	7
Guidance and manufacturer's declaration – Electromagnetic emissions.....	8
Guidance and manufacturer's declaration – Electromagnetic immunity.....	8
Guidance and manufacturer's declaration – Electromagnetic immunity.....	9
Recommended separation distances between portable and mobile RF communications equipment and the Bimeo PRO device.....	10
Minimum PC system requirements.....	11
Recommended PC system requirements.....	11
Bimeo PRO parts.....	16
Descriptions of parts.....	18
Bimeo PRO sensing units.....	20
Bimeo PRO therapy Modes.....	20
Bimeo PAD.....	22
Indications.....	22
Contraindications for the Bimeo PRO use.....	22
Risk factors.....	23
Expected and predictable side effects.....	23

Information on circumstances when user should consult with medical worker.....	23
Chapter 2. Bimeo PRO software installation and configuration	
Installing Bimeo PRO software.....	24
Uninstalling Bimeo PRO software.....	27
Bimeo PRO license.....	27
Troubleshooting.....	28
The Bimeo PRO drivers are not installed.....	28
Chapter 3. Bimeo PRO sensors and attachments	
Bimeo PRO sensors.....	28
Troubleshooting.....	29
Charging the Bimeo PRO system.....	29
Bimeo PRO setup.....	30
Installing Bimeo sense units.....	30
Mounting Bimeo attachments.....	32
Chapter 4. Bimeo PAD	
Troubleshooting.....	33
Charging the Bimeo PAD.....	33
Bimeo PAD setup.....	34
Positioning Bimeo PAD.....	34
Troubleshooting.....	34
Chapter 5. Bimeo software	
General overview.....	36
Main user interface window.....	37
Language setting area.....	37
Bimeo sensors status indicators.....	37
Patient related information area.....	38
Add new patient or select an existing patient from the database.....	38
Edit patient information and configure therapy program for the patient.....	40
Analysis of patient's performance.....	42
Export of patient's performance.....	43
Therapy assessment parameters.....	45
Therapy related settings.....	47
Therapy modality.....	48
Movement modality selection.....	49

Therapy and assessment scenarios.....	53
Therapy program selection.....	55
Visual representation of therapy settings.....	55
Therapy control functions.....	56
Therapy settings status bar.....	56
Initialization procedure.....	58
Range-of-motion settings.....	60
General scenario user interface.....	62
Assessment scenarios.....	66
Reaching.....	66
Tracking.....	67
Labyrinth.....	68
Rotation.....	69
Therapy scenarios.....	70
Kitchen (ADL task).....	70
Fruit-picker (ADL task).....	71
Puzzle (cognitive task).....	72
Memory (cognitive task).....	72
Slide (cognitive task).....	73
Ghost scenario (game).....	74
Space (game).....	74
Tic Toc (game).....	75
Rocket (game).....	76
Racer (game).....	77
Snowball (game).....	77

2.1.1. System for recovery of motion activity

Bimeo PRO is a system for recovery of motion activity, motion coordination of upper limbs using interactive programs.

Bimeo PRO is a passive device consisting only of wireless inertial and force sensors. In order to promote affected arm movement the system exploits the preserved patient's motor functions. The less affected limb can assist the movements of the more affected limb, while its contribution is continuously measured and evaluated. Data collected from sensors attached to the user is wirelessly transmitted to a computer. The user interacts with a virtual environment adapted for motor rehabilitation therapy. Voluntary activation of the affected limb is maximized by adaptation of the therapy scenario. An expert system acquires and evaluates information on the status and progress of the user. Information is immediately available to clinicians in order to optimize therapy parameters. The Bimeo PRO system was designed to combine intensive and safe unimanual and bimanual therapy with the capability of objective assessment and cognitive visual and auditory feedback that characterize rehabilitation robots.

2.1.2. System components

System for recovery of motion activity, motion coordination of upper limbs using interactive programs Bimeo PRO, consists of:

1. Manipulator Bimeo Master no more than 3 pcs. (option).
2. Additional manipulator Bimeo Servant no more than 3 pcs. (option).
3. Upper arm sensor no more than 3 pcs. (option).
4. Forearm sensor no more than 3 pcs. (option).
5. Wireless receiver no more than 3 pcs. (option).
6. Support disc for training of one arm no more than 3 pcs. (option).
7. Double support disc for training of two arms no more than 3 pcs. (option).
8. Pad for training on table Bimeo PAD no more than 3 pcs. (option).
9. Anti-slip support no more than 3 pcs. (option).
10. USB flash card with software no more than 3 pcs. (option).
11. Charger no more than 3 pcs. (option).
12. Additional charger for Bimeo PAD no more than 3 pcs. (option).
13. Triple micro USB charging cable no more than 3 pcs. (option).
14. Micro USB cable for additional charger no more than 3 pcs. (option).
15. Storage case (option).
16. Set of straps - S, M, or L size – (includes 4 straps per set) no more than 2 set (option).
17. Table for training – 1 pcs. (option).
18. User manual – 1 pcs.

Figure 1 Bimeo PRO parts


Figure 2 Bimeo PAD and therapy table

- (1) Manipulator Bimeo Master (hereinafter – Bimeo master),
- (2) Additional manipulator Bimeo Servant (hereinafter – Bimeo servant),
- (3) Bimeo sense upper arm sensor,
- (4) Bimeo sense forearm sensor,
- (5) Wireless receiver (hereinafter – Bimeo receiver),
- (6) Support disc for training of one arm (hereinafter – Bimeo disc),
- (7) Double support disc for training of two arms (hereinafter – Bimeo double disc),
- (8) Pad for training on table Bimeo PAD (hereinafter – Bimeo PAD),
- (9) Anti-slip support,
- (10) USB flash card with software,
- (11) Charger,
- (12) Additional charger for Bimeo PAD,
- (13) Triple micro USB charging cable,
- (14) Micro USB cable for additional charger,
- (15) Storage case,
- (16) Set of straps - S, M, or L size,
- (17) Table for training.

Part Nr.	Name Description	Dimensions	Mass
(1)	Manipulator Bimeo Master Main active sensing unit held in patient's affected hand, used in every therapy mode.	13 cm x 12 cm x 7 cm $\pm 10\%$	192 g $\pm 10\%$
(2)	Additional manipulator Bimeo Servant Therapy attachment held in patient's unaffected hand, used in bimanual therapy modes.	12 cm x 12 cm x 6 cm $\pm 10\%$	75 g $\pm 10\%$
(3)	Bimeo sense upper arm sensor Active sensing unit attached to patient's affected upper arm.	6 cm x 5 cm x 1.3 cm $\pm 20\%$, strap length 38 cm $\pm 20\%$	38 g $\pm 20\%$
(4)	Bimeo sense forearm sensor Active sensing unit attached to patient's affected	6 cm x 5 cm x 1.3 cm, $\pm 20\%$	35 g $\pm 20\%$

	forearm.	strap length 30 cm $\pm$ 20%	
(5)	Wireless receiver USB key enables communication between Bimeo sensing units and software package.	7.5 cm x 2.5 cm x 1.4 cm $\pm$ 20%	16 g $\pm$ 20%
(6)	Support disc for training of one arm Coupling element for joining different therapy attachments, also used for unimanual training on support surface.	12.5 cm x 12.5 cm x 1.5 cm $\pm$ 10%	150 g $\pm$ 10%
(7)	Double support disc for training of two arms Therapy attachment used for bimanual training on support surface and for training of wrist.	30 cm x 12.5 cm x 3 cm $\pm$ 10%	255 g $\pm$ 10%
(8)	Pad for training on table Bimeo PAD An extension sensor for surface therapy mode.	60 cm x 77 cm x 1.7 cm $\pm$ 10%	4100 g $\pm$ 10%
(9)	Anti-slip support Therapy attachment that prevents bimeo double disc to slip during training of wrist.	30 cm x 12.5 cm x 0.3 cm $\pm$ 20%	9 g $\pm$ 20%
(10)	USB flash card with software Bimeo PRO software for training on PC.	4 cm x 1.3 cm x 0.5 cm $\pm$ 20%	6 g $\pm$ 20%
(11)	Charger Charger for charging bimeo master and bimeo sense units.	4.5 cm x 7.5 cm x 8 cm $\pm$ 20%, cable length 2 m $\pm$ 20%	140 g $\pm$ 20%
(12)	Additional charger for Bimeo PAD Additional charger for charging bimeo PAD sensor.	4 cm x 9 cm x 9 cm $\pm$ 20%, cable length 1 m $\pm$ 20%	120 g $\pm$ 20%
(13)	Triple micro USB charging cable Cable for charging with three uUSB connectors to charge Bimeo sensors	cable length 40 cm $\pm$ 20%	27 g $\pm$ 20%
(14)	Micro USB cable for additional charger Cable for charging with one uUSB connector to charge Bimeo PAD.	cable length 11 cm $\pm$ 20%	9 g $\pm$ 20%
(15)	Storage case Aluminum case for transportation and storage of bimeo.	47 cm x 37 cm x 11 cm $\pm$ 10%	2530 g $\pm$ 10%
(16)	Set of straps - S, M, or L size Additional straps for Bimeo master, Bimeo sense upper arm and Bimeo sense forearm sensor.	S size lengths: 32 cm, 30 cm, 35 cm 27 cm $\pm$ 10% M size lengths: 32 cm, 30 cm, 38 cm 30 cm $\pm$ 10% L size lengths: 32 cm, 30 cm, 42 cm 34 cm $\pm$ 10%	60 g $\pm$ 10% 61 g $\pm$ 10% 62 g $\pm$ 10%
(17)	Table for training Height-adjustable therapy table to support the patient during therapy. The table size and design are specially adjusted for the use of Bimeo PRO device and the attachment Bimeo PAD. Tabletop is laminated MDF and the frame material is aluminium. User manual Bimeo PRO instructions for use.	Height: 640 mm to 960 mm $\pm$ 10% Width: 1200 mm $\pm$ 10% Depth: 900 mm $\pm$ 10%	30 kg $\pm$ 10%

3.2.2. Bimeo PRO system

The main unit of the Bimeo PRO system is the ergonomically shaped and soft rubber covered handle, Bimeo master, that the user comfortably holds in his/her affected arm (Figure 3). The handle houses a custom built three-axes force sensor (range 100 N in all axes) and a three-dimensional inertial measurement unit (IMU - consisting of a gyroscope, an accelerometer and a magnetometer). The Bimeo master can be combined with various adapters that enable different therapy modes (Figure 3).

In addition to the main unit, the Bimeo PRO system includes two IMUs, Bimeo senses, to be attached on the upper arm and the forearm of the exercising patient. The two units measure orientation of the respective arm segment and combined with the information from the handle provide complete kinematic data for the arm.

The Bimeo PAD sensor is an additional sensor that indicates the position of the arm in the support surface therapy mode more accurately. The therapy in the support surface mode can be performed with or without the Bimeo PAD sensor.

3.2.3. Application modes

The sensor-based system for unimanual and bimanual therapy enables diverse motor function and assessment tasks with only minor system adaptations (for example, by changing the Bimeo master attachment). The application modes are shown in Figure 4.

During bimanual therapy the patient holds the device with both hands (Figure 4a). The following tasks can be performed: 1) passive scratching of the more affected limb with the support of the less affected limb where all interaction forces and postures are monitored, 2) bimanual reaching and tracking in three-dimensional space where the less affected limb provides support and assistance (or resistance) to the more affected limb, and 3) isometric therapy by generating opposing forces on the device (pushing or pulling) with both arms.

During unimanual therapy the patient holds the device in the affected arm (Figure 4b). The following tasks can be performed: 1) measurement of the limb active range of motion by voluntarily positioning the arm (or individual joint) to its minimal and maximal positions, and 2) reaching and tracking in three-dimensional space augmented with virtual reality scenarios (simulation of activities of daily living) or real objects.

Therapy on support surface (Figure 4c) enables unimanual or bimanual partially restrained movement, where the surface provides weight support for the arm and resistive force to the movement through friction between the surface and the disc adapter.

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Bimeo PAD

Therapy on a support surface including Bimeo PAD enables unimanual or bimanual partially restrained movement, where the surface provides weight support for the arm and resistive force to the movement through friction between the surface and the disc adapter. In a bimanual therapy the patient holds the device with both hands (Figure 5a). In unimanual therapy the patient holds the device in the affected arm (Figure 5b).


Figure 5 The Bimeo PAD used in support surface: a) bimanual therapy; b) unimanual therapy.

Clinical application

Indications

The Bimeo PRO is intended for rehabilitation of patients with disability in upper extremity, for example patients with hemiplegia, hemiparesis or tetraparesis due to stroke, spinal cord injury, traumatic brain injury, cerebral palsy, multiple sclerosis, amyotrophic lateral sclerosis or alike. Patients who suffered upper extremity injury or surgical procedure of an upper limb also profit from rehabilitation enhanced by the Bimeo PRO device. Eligible patients should have preserved at least minimal motor function in one of the upper limbs (both limbs can be impaired but at least in one limb minimal motor function is required). Device is intended for use in adults and children age 7 and above, if they understand instructions and their hands fit into the main Bimeo PRO sensor.

Contraindications for the Bimeo PRO use

The decision as to whether a patient is suitable for therapy with the Bimeo PRO device comes under the remit of the prescribing physician, who has sole medical responsibility for the therapy. As part of this, the prescribing physician must evaluate in each individual case, possible risks and side-effects of the therapy against the benefit gained from it.

The lists with contraindications is as follows, but does not claim to be exhaustive:

- Certain implanted ferromagnetic materials or active devices within the body (such as pacemakers or pumps)

- Bone instability (fractures, severe osteoporosis)
- Non-stable vital signs: pulmonary or cardio-circulatory contraindications, etc.
- Inability to sustain posture
- Inability to comply with instructions or procedure
- Severe ataxia
- Oculomotor Deficits (e.g. inability of vertical and horizontal eye tracking)
- Sensors cannot be fitted to the relevant arm

Contraindications

- Cardiovascular conditions (e.g. heart failure, subacute phase after myocardial infarction, autonomic failure),
- Severe balance deficits with increased risk of falling
- Cognitive deficits
- Epileptic seizure or photosensitivity
- Retinal disease which may be exacerbated by watching video games.

Additional information

None known.

Additional information for the user

The user is a medical worker.

Version

Installation

The purpose of this section is to provide users of the Bimeo PRO device with a simple procedure for installation of Bimeo PRO software under Windows 7, Windows 8 or Windows 10 operating system. The Bimeo PRO setup wizard will install Bimeo PRO software and Bimeo PRO drivers.

Installing Bimeo PRO software

The installation package comes as an all-in-one executable installer wizard. Launch the executable `setup.exe` by double-clicking the setup file. The **Welcome to the setup wizard window** will be displayed, click **NEXT** to continue (Figure 6).


Figure 6 Welcome screen of the setup wizard

Select the installation folder where the setup will install Bimeo PRO software. Click **NEXT** to continue (Figure 7). Ensure that the disk that you want to select has enough space to store all data.


Figure 7 Select installation folder

Select the folder where the setup will install Bimeo PRO Start Menu shortcuts. Click **NEXT** to continue (Figure 8).


Figure 8 Select Start menu folder

By selecting checkbox **Additional language** and by choosing a language, you add the language which will be available. Click **NEXT** to continue (Figure 9).


Figure 9 Language selection

Click **INSTALL** to start Bimeo PRO software installation (Figure 10). Installation may take a few minutes.


Figure 10 Review of install settings

The **Welcome to the Device Driver Installation Wizard** will be shown during installation (Figure 11). Click **NEXT** to continue.

Figure 11 Welcome screen of the driver installation wizard

Select **I accept this agreement** after you have read the License Agreement. Click **NEXT** to install the Bimeo PRO drivers.




Figure 12 Driver installation: License Agreement.

After drivers are successfully installed the screen in Figure 13 is shown. Click **FINISH**.

When the installation is complete the screen in Figure 14 is shown. To end the installation click **Finish**. You can now insert the Bimeo receiver into the USB port.

You can uninstall Bimeo PRO software by clicking Start button -> Programs -> Bimeo PRO -> Uninstall Bimeo PRO. Confirm uninstall by clicking **Yes**. When the setup finishes, click **OK**. The Bimeo PRO software is uninstalled.

The patient database is not deleted when uninstalling Bimeo PRO software and remains in current user's directory.

You must enter the software license key after the installation, the first time when the Bimeo software is started. The software license key can be found in the storage case under Bimeo sense upper-arm sensor. The software license key is linked to each Wireless receiver and can be only used with Wireless receiver provided within the storage case.

In case when the Bimeo PRO drivers are not installed automatically, you can run driver setup manually by double-clicking dpinst-amd64.exe (64-bit operating system) or dpinst-x86.exe (32-bit operating system) in »BimeoDriver\CDM« folder of Bimeo PRO installation folder (default: C:\Program Files\Bimeo PRO\BimeoDriver\CDM). Follow the setup wizard to install drivers.

Item	Description
1	Charging port
2	Power switch with status light
3	The arrow pointing toward the patient

To turn on the Bimeo sense unit turn on the Power switch marked as number 2 in Figure 15. The status light will turn blue for five seconds. After five seconds the status of the Bimeo sense unit is indicated as a short blink of blue light every two seconds. To turn off the Bimeo sense unit turn off the Power switch. The blue blinking light will turn off afterwards.

Item	Description
1	Charging port
2	Power switch with status light

To turn on the Bimeo master turn on the Power switch marked as number 1 in Figure 16. The button status light will turn blue for five seconds. After five seconds the status of the Bimeo master is indicated as a short blink of blue light every two seconds. To turn off the Bimeo master turn off the Power switch. The blue light will remain off.

THE BIMEO MASTER OR IMU TURNS OFF BY ITSELF

If the battery voltage is too low for normal operation, the unit will turn off automatically. Please charge the system immediately to prevent self-discharging of the battery causing the battery to lose its performances.

THE STATUS LIGHT DOES NOT TURN ON AFTER SWITCHING ON THE POWER BUTTON – THE SYSTEM DOES NOT POWER UP

The battery voltage is too low. Please charge the system immediately to prevent self-discharging of the battery causing the battery to lose its performances.

Plug the Bimeo charger in a wall power outlet. To charge the Bimeo master, connect one micro USB cable of Bimeo charger to Bimeo master charging port. To charge the Bimeo sense unit, connect one micro USB cable of Bimeo charger to the unit's charging port. Make sure that the system is turned off prior charging.

When the system is charging, the status light of each unit will turn white. The charging cycle lasts approximately 2 hours. The charging of each unit is complete when the white status light turns off.

If the unit is charging for more than four hours and the status light is still white, please disconnect the USB cable. In such a case the battery of the system has lost some of its

performances and cannot be completely charged. The battery should only be replaced by authorized personnel.

Important information

The Bimeo PRO system consists of three units. Since all units measure orientation of the segment that is attached to, it is important that each unit is correctly oriented and attached on the corresponding segment of the arm of the exercising patient. The incorrect placement or orientation of the units prevents the system to provide a complete kinematic data for the arm.



Attaching the Bimeo PRO system

Two Bimeo sense units must be attached on the upper arm and the forearm. To distinguish where to attach each unit, observe the mark on the unit and the attachment belt.

- Insert the Bimeo sense unit with marked upper arm into attachment strap with larger diameter. Insert the Bimeo sense unit with marked forearm into attachment strap with smaller diameter.
- The unit that has attachment strap with larger diameter must be attached on the upper arm and the unit with the smaller diameter of the attachment strap must be attached on the forearm.
- It is important to correctly orient the Bimeo sense unit. Figure 17 shows the correct placement and orientation of the IMU units.

The Bimeo sense units must be oriented in a way that the exercising patient is able to see the "figure" on the Bimeo sense unit from his/her view perspective and the arrow is pointing toward the patient.

The same rules apply if the system is attached on the left hand (see Figure 18).

- When attaching the Bimeo sense units, make sure that the straps are tightened to a comfortable level. Do not tighten too much.
- After attachment is done, make sure that all the units are powered on.

4.2.2. Attaching the Bimeo PRO system

The Bimeo PRO system can be used to perform bimanual and unimanual tasks. Some of the tasks require exchanging and combining of provided adapters – attachments. The software will provide visual information of which adapter should be used for a selected task.

- When attaching the Bimeo master or Bimeo servant, note the correct orientation of the key (see Figure 19).


Figure 19. Bimeo master/servant unit

- To fixate the Bimeo master (or Bimeo servant) to Bimeo disc or Bimeo double disc, rotate the Bimeo master (Bimeo servant) 90 degrees clockwise until the Bimeo master (Bimeo servant) is locked (see Figure 20 a).
- To release the Bimeo master (or Bimeo servant) from Bimeo disc or Bimeo double disc rotate the Bimeo master (Bimeo servant) 90 degrees counter clockwise until the Bimeo master (Bimeo servant) is released (see Figure 20 b).



Figure 20. Attaching and releasing the Bimeo master/servant unit from the Bimeo disc or Bimeo double disc

Figure 20

Item	Description
1	Charging port
2	Power switch with status light

To turn on the Bimeo PAD turn on the Power switch marked as number 2 in Figure 21. The status light will turn blue for five seconds. After five seconds the status of the Bimeo PAD is indicated as a short blink of blue light every two seconds.

To turn off the Bimeo PAD turn off the Power switch. The blue blinking light will turn off afterwards.

When turning on the Bimeo PAD do not place any objects onto the bottom surface of the Bimeo PAD sensor. The Bimeo disc or Bimeo double disc can be placed onto the sensor when the status light starts blinking.

THE BIMEO PAD TURNS OFF BY ITSELF

If the battery voltage is too low for normal operation, the unit will turn off automatically. Please charge the system immediately to prevent self-discharging of the battery causing the battery to lose its performances.

THE STATUS LIGHT DOES NOT TURN ON AFTER SWITCHING ON THE POWER BUTTON – THE SYSTEM DOES NOT POWER UP

The battery voltage is too low. Please charge the system immediately to prevent self-discharging of the battery causing the battery to lose its performances.

Plug the Bimeo charger in a wall power outlet. To charge the Bimeo PAD, connect micro USB cable of Bimeo charger to Bimeo PAD charging port.


When the system is charging, the status light of will turn white. The charging cycle lasts approximately 3 to 4 hours. The charging of each unit is complete when the white status light turns off.

If the unit is charging for more than six hours and the status light is still white, please disconnect the USB cable. In such a case the battery of the system has lost some of its performances and cannot be completely charged. The battery should only be replaced by authorized personnel.

Bimeo PAD can only be used together with Bimeo arm sensors. Attach the Bimeo PRO sensors as described before.

The Bimeo PAD must be positioned in front of the patient on the working surface - table.

- It is important to correctly orient the Bimeo PAD on the table. Figure 22 shows the correct placement and orientation.

THE CURSOR ON THE SCREEN JUMPS AROUND DURING THERAPY WITH BIMEO PAD

If the cursor on the screen jumps around there could be two main problems.

- The Bimeo receiver is losing communication signal to Bimeo PAD. Try changing the Bimeo receiver (PC) position and orientation to see if the cursor on the screen moves normally. You could connect the Bimeo receiver to the USB extension cable to find more suitable place for Bimeo receiver.

- The patient is touching the Bimeo PAD sensor bottom surface with the elbow or hand (Figure 23). The patient should not be touching the bottom surface of the Bimeo PAD. The patient should not wear long sleeve shirt. The sleeves could also disturb the position sensing.


Figure 23. The patient is touching the Bimeo PAD sensor bottom surface with the elbow or hand.

1. Introduction

Bimeo PRO software is related to the Bimeo PRO hardware. Its main functionalities are:

1. User interface for all tasks related to therapy and patient evaluation,
2. Database of patients,
3. Therapy tasks,
4. Evaluation tasks,
5. Definition of therapy programs for individual patients,
6. Definition of patient's range of motion,
7. Evaluation of patient's performance,
8. Export of patient's results.

Bimeo PRO software is organized as follows:

- a. Patient related information
- b. Therapy related settings
- c. Therapy setting status bar

Main user interface window

Main user interface window (Figure 24) is divided into 7 functional display areas:

- A. Language setting area (A) for setting the user interface language
- B. Bimeo sensors status indicators (B) indicate status of Bimeo sensors
- C. Patient related information (C) area is used to configure and display all patient related data
- D. Therapy related settings are used in three modes
 - a. Configuration of therapy task
 - b. Configuration of therapy task to be added to the Therapy program
 - c. Selection of therapy task for which patient's data is analyzed
- E. Visual representation of therapy settings (E) is composed of three images for better presentation
- F. Therapy control functions (F) are important for conducting the therapy
- G. Therapy settings status bar (G) shows instantaneous status of all parameters relevant for therapy


Figure 24 Main user interface window: A) language settings, B) Bimeo sensors status indicators, C) patient related information, D) therapy related settings, E) visual representation of therapy settings, F) therapy control functions, and G) therapy settings status bar.

LANGUAGE SETTING AREA

Pressing the button with the country flag in the language setting area (Figure 24 A) enables language selection for the Bimeo PRO. The default language is defined in relation to installation setting. By sequentially clicking on the country we can switch between English and additional language.

BIMEO SENSORS STATUS INDICATORS

Bimeo sensor status indicators (Figure 24 B) indicate activity of the three Bimeo sensors

1. Bimeo PAD sensor – the indicator on top,
2. Bimeo upper arm sensor – the second indicator,

3. Bimeo forearm sensor – the third indicator,
4. Bimeo master – the indicator below.

When signal strength and battery status are displayed next to the sensor image (as in the case of forearm and upper arm sensors in Figure 24 B), the sensor operates normally. If a red cross is displayed next to the sensor (as in the case of Bimeo master in Figure 24 B), the sensor is either turned off or malfunctioning. If all sensors operate normally, the Bimeo sensors status indicators are not displayed on the screen during therapy while they are always displayed in the main user interface.

The Bimeo PAD sensor status is shown only if the Bimeo PAD is turned on. Then the therapy on the support surface can be performed on the Bimeo PAD. If there is no Bimeo PAD sensor indicator, the arm position is estimated from the Bimeo upper arm, Bimeo forearm and Bimeo master sensors.

3.2.3.2. Patient menu and patient data

All aspects of patient settings and data analysis can be accessed from the patient menu positioned at the left part of the main user interface window (Figure 24 and Figure 25). Three functions can be accessed from the patient menu

- A. Add new patient or select an existing patient from the database (A),
- B. Edit patient information and configure therapy program for the specific patient (B),
- C. Visualize patient performance from all previous therapy sessions and export patient results (C).


Figure 25. Patient menu. 1 – Add new patient or select an existing patient from the database (A), 2 – Edit patient information and configure therapy program for the specific patient (B), 3 – Visualize patient performance from all previous therapy sessions and export patient results (C).

3.2.3.2.1. Add new patient

New patient can be added to the database or an existing patient can be selected from the database by clicking on the **ADD NEW PATIENT** button (Figure 25 A).

By clicking on the **ADD NEW PATIENT** button a submenu appears with the search field for finding patients in the database (Figure 26). The found patient's names are shown right to the search field. A patient can be selected by clicking on the name of the patient.

Below the search field is the **ADD NEW PATIENT** button that opens a new window for entering patient's data (Figure 27).

When a new patient is added, the following data can be entered:

1. Name
2. Surname
3. Date of birth
4. Gender
5. Height [cm] (enter only the number in the entry field)
6. Weight [kg] (enter only the number in the entry field)
7. Injury type (e.g. stroke, SCI, ...)
8. Year of injury occurrence
9. Affected side (left side, right side or both sides)
10. Clinical score; up to four different clinical scores can be selected from the clinical score menu (Figure 27 right side) that opens by clicking on the clinical score button; clinical scores need to be entered as numerical values to be plotted in a time series plot (tracking of changes of clinical scores). Initial clinical score values can be inserted when a new patient is added to the database. Clinical scores can be later updated from the clinical score menu. Clinical scores do not affect the Bimeo therapy and are not evaluated during therapy. Clinical scores should separately be done by therapists.

Name, surname, height, weight and affected side are mandatory entries. A new patient is entered into the database by clicking the **NEW PATIENT** button below. By clicking the **EXIT** button the process of inserting a new patient is aborted and the patient is not inserted into the database.


Figure 27: BimeoPRO software interface (Screenshot of the BimeoPRO software interface)

BimeoPRO EXIT

BIMANUAL THERAPY

VERTICAL PLANE

ASSESSMENT

THERAPY

THERAPY PROGRAM

PATIENT

EDIT

DATA ANALYSIS

NEW PATIENT

Enter personal details

Name:

Surname:

Date of birth:

Gender:

Enter physical details

Height [cm]:

Weight [kg]:

Injury type:

Year of injury:

Affected side: ☐ Left ☒ Right

Left arm clinical evaluation

Right arm clinical evaluation

Clinical evaluation

✓ SYSTEM ON

✗ SYSTEM INITIALIZED

✗ PATIENT

✗ THERAPY MODALITY

✗ MOVEMENT MODALITY

✗ ASSESSMENT

✗ THERAPY

By clicking on the **EDIT** button a window for editing patient's data opens (Figure 28).


The following parameters can be edited

1. Name
2. Surname
3. Date of birth
4. Gender
5. Height [cm] (enter only the number in the entry field)
6. Weight [kg] (enter only the number in the entry field)
7. Injury type (e.g. stroke, SCI, ...)
8. Year of injury occurrence
9. Affected side (left side, right side or both sides)
10. Clinical score; up to four different clinical scores can be selected from the **CLINICAL EVALUATION** menu (Figure 27 right side) that opens by clicking on the **CLINICAL EVALUATION** button; clinical scores need to be entered as numerical values to be plotted in a time series plot (tracking of changes of clinical scores). If clinical scores are updated from the **CLINICAL EVALUATION** menu, the previous values are stored in the database. Thus, the entire history of clinically evaluated performance is available for analysis of patient's progress.

Edited patient's data is entered into the database by clicking the **ENTER** button below. By clicking the **CANCEL** button the process of editing patient is aborted and the updated information is NOT inserted into the database.

Patient can be deleted from the database by clicking the **DELETE PATIENT** button. This option is enabled when editing patient's data and if at least two patients are currently in


database. When deleting a patient a pop-up window appears to confirm the deletion of the current patient.

By clicking on the **EDIT** button in the **PATIENT** section a **EDIT PATIENT** window is opened (Figure 29). Up to six different therapy programs can be configured for each patient. Each therapy program consists of a sequence of up to six therapy tasks. Each task can be configured independently. The **Therapy program** consists of the following areas:

1. **Therapy program selection**; six buttons enable selection of one of the six therapy programs: green icon behind the program number indicates the current active therapy program that is being configured, yellow icon behind the program number indicates a therapy program that is configured and active for patient's therapy, red icon behind the program number indicates a therapy program that is not active for patient's therapy.
2. **Therapy program name**; enables setting of the therapy program name by which the it is identified during therapy.
3. **Therapy program side**; the therapy program can be specified either for the left or the right side; only those therapy programs that correspond to the current setup (left or right arm initialized) are available during therapy.
4. **Therapy program sequence**; tasks in the therapy program can be presented as in the list of tasks pertaining to the therapy program or in a random order (when random sequence radio button is active).
5. **Therapy program auto-switch**; the therapy program automatically switches to the next task in the therapy program if this option is enabled. Otherwise the main menu is shown and the user has to manually click **ENTER** button.

6. **Configuration** button is used for configuring the therapy program; therapy tasks are defined by using the **Configuration** button (Figure 33); see chapter **Configuration** below!
7. **Configuration** button; radio buttons are used to activate or deactivate a specific task within the therapy program.
8. **Configuration** button is used for configuring the specific task within the therapy program; parameters are taken from the current settings defined by **Configuration** button (Figure 33); see chapter **Configuration** below!
9. **Configuration** button of each task in the therapy program can be changed with **Configuration** buttons. The **Configuration** button is show at the bottom-right corner of the window.
10. **Configuration** button stores the therapy program settings into the database.
11. **Configuration** button cancels the configuration of the therapy program.

An important aspect of intensive therapy is the analysis of patient's performance through consecutive therapy sessions. Data from each therapy session is stored into the database and available for later analysis. By clicking on **Configuration** button a new window for visualization of therapy results opens (Figure 30). The data analysis window is split into the following areas

1. **Configuration** button is the main window that presents data in a graphical form.
2. **Configuration** button on top of time series plot selects analysis of data for the left arm.
3. **Configuration** button on top of time series plot selects analysis of data for the right arm.
4. **Configuration** button on the right side enables selection of a parameter to be plotted (all parameters are explained in the **Configuration** section below) on the time series plot (by clicking on a button with the parameter name). At the same time numerical values for all parameters for a specific therapy

session are shown. The specific therapy session is indicated by date and time of the therapy session at the upper-right corner of the **Therapy assessment window**. The same therapy session is also indicated with a red square marker in the time series plot.

5. **Slider** below the time series plot enables selection of therapy session for which numerical values of individual parameters are presented in the **Therapy assessment window** on the right. Use the slider by pulling it left and right to see the red marker move across the blue markers indicating individual therapy session result.
6. **FINISH ANALYSIS** button terminates the data analysis.

Data analysis is performed for each therapy task separately. The therapy task for which patient's performance is analyzed is specified with the four buttons in the Therapy related settings area of the Main user interface window (see chapter Therapy related settings below).


Figure 30 Analysis of patient's performance window

EXPORT OF PATIENT'S PERFORMANCE

The therapy data of each patient can be exported to PDF file. First, the patient has to be selected. When the patient is selected, the export button is accessible in analysis of the patient's performance window (Figure 30).

The destination folder, where the exported PDF file is going to be saved, can be selected in the lower right corner of the analysis window (Figure 30). The available options are desktop and documents folder. The therapy data can be exported by clicking **EXPORT PATIENT'S RESULTS** button. A new PDF is created in the selected folder. The filename consists of patient's name, surname and date of exportation («name_surname_date.pdf»). If the PDF file with the same filename already exists in the same directory, it is overwritten.

Figure 31 shows an example of an exported file. The patient trained four different therapy tasks. The date of document's creation and the page number is in the upper left corner. Patient's information is located in the first table. It contains patient's name, surname, date of

birth, gender, height, weight, injury type, year of injury, affected side and the time period of therapy.

The next tables contain therapy data and analysis obtained from different therapy tasks. The therapy tasks are defined in the second row. There are as many tables as there are different therapy tasks played. To include the particular therapy data in the exported file, the therapy task must last at least 20 seconds. To include a particular therapy task in the exported file, the same therapy task has to be performed at least 2 times.

bimeo

Date: 19.04.2015
Page: 1

Name	Surname	Date of birth	Gender	Height [cm]	Weight [kg]	Injury type	Year of injury	Affected side
John	Doe	23.09.1990	m	175	190	sh	2015	RIGHT ARM
Exercises			From:	To:				
			09.02.2015	20.02.2015				

SELECTED ARM	MANUAL MODE			VERTICAL PLANE		JOINT MODE		GAME MODE		NUM. OF EXERCISES	
RIGHT ARM	BIMANUAL TRAINING			VERTICAL PLANE		2		REACHING		8	
	QI	PER	OPT	VEL	DEV	SMT		SUP	MVT	THT	SC
MIN/MAX	55/77	64/100	37/77	80/02	57/77	74/01		35/66	14/19	184/0	40/0/90/0
MEAN/STD	71/0.0	98/0.5	60/1.7	86/0.1	69/0.6	84/0.5		50/1.1	15/0.1	280/0/397	81/0/15/6
IV/TREND	162/0.23	9/1/0.13	57/0.40	87/0.02	64/0.16	80/0.13		16/0.39	18/-0.03	254/0	66/3/4.21
										10.00	

SELECTED ARM	MANUAL MODE		VERTICAL PLANE		JOINT MODE		GAME MODE		NUM. OF EXERCISES	
RIGHT ARM	BIMANUAL TRAINING		VERTICAL PLANE		2		TRACKING		3	
	QI	PER	OPT	VEL	DEV	SMT	SUP	MVT	THT	SC
MIN/MAX	30/49	45/78	31/57	60/65	70/84	04/14	13/27	63/63	184/0	125/0
MEAN/STD	41/0.6	58/1.4	42/1.1	63/0.7	78/0.6	08/0.4	22/0.6	63/0.0	264/0/56.6	190/7/49.4
IV/TREND	16/0.57	44/1.41	30/1.26	65/-0.25	71/0.69	04/0.39	23/-0.04	63/0.00	204/0	131/2
									60.00	59.50

Page 2

SELECTED ARM	MANUAL MODE		VERTICAL PLANE		JOINT MODE		GAME MODE		NUM. OF EXERCISES	
RIGHT ARM	BIMANUAL TRAINING VERTICAL PLANE						TRACKING		5	
	QI	PER	OPT	VEL	DEV	SMT	SUP	MVT	THT	SC
MIN/MAX	55/69	67/82	57/75	56/71	83/85	38/51	13/58	52/63	304/0	215/0
MEAN/STD	60/0.5	69/0.7	65/0.7	65/0.7	84/0.1	57/0.8	25/1.7	50/0.5	304/0/20	239/6/17.5
IV/TREND	56/0.17	77/-0.29	68/-0.13	57/0.37	84/-0.01	44/0.42	06/0.95	61/-0.31	304/0/0.00	254/0/-7.20

SELECTED ARM	MANUAL MODE			VERTICAL PLANE		JOINT MODE		GAME MODE		NUM. OF EXERCISES	
RIGHT ARM	BIMANUAL TRAINING			VERTICAL PLANE				LABYRINTH		9	
	QI	PER	OPT	VEL	DEV	SMT	SUP	MVT	THT	SC	
MIN/MAX	00/81	00/80	00/95	00/100	00/86	00/100	00/70	00/123	43/0/304/0	00/45/0	
MEAN/STD	64/2.5	63/2.3	74/2.8	88/3.1	58/3.1	76/2.9	55/2.0	91/3.4	261/7/89.9	281/13.2	
IV/TREND	10/0.63	4/0.57	45/0.72	68/0.50	10/0.98	52/0.60	42/0.34	81/0.24	177/5	125/4.15	

Figure 31: An example of exported patient's results file. The therapy data is shown for each therapy task separately. The abbreviations are defined below.

Figure 32 shows an example of a Movement QI plot in the exported file. The plot is shown for four assessment scenarios and for left and right arm, respectively. Therefore, there could be maximum 8 plots. To include the particular assessment scenario in the exported file, the assessment scenario must be performed at least 2 times.

MOVEMENT QI (LABYRINTH. RIGHT ARM)


Figure 2. Movement quality index (QI) for the right arm in a labyrinth task.

Table 2. Therapy assessment parameters. The parameters are computed from the real-time measurements available from the Bimeo sensors. These parameters provide objective information about the quality of movement of the patient

Abbreviation	Parameter
QI	Movement quality index
PER	Performance
OPT	Optimality
VEL	Velocity
DEV	Deviation
SMT	Smoothness
SUP	Support
MVT	Movement time
THT	Therapy time
SC	Score
MIN/MAX	Minimum value / Maximum value
MEAN/STD	Mean value / Standard deviation
IV/TREND	Initial value / Trend of the parameter

Therapy assessment parameters are computed from the real-time measurements available from the Bimeo sensors. These parameters provide objective information about the quality of movement of the patient (Table 2).

Result of the therapy session in terms of the number of achieved scores	
(Assessment & Therapy)	
[score/s]	
(Assessment & Therapy)	Result of the therapy session in terms of the number of achieved scores normalized to the therapy time in order to simplify comparison between sessions and between patients
[0-10]	
(Assessment only)	Non-weighted average value of the following parameters (general movement quality index for the entire therapy session): Performance, Optimality, Velocity, Deviation, Smoothness and Support
[0-10]	
(Assessment only)	"Fitts law" information index that defines quality of movement in terms of accuracy of hitting the target and the speed of executing

[0-10]
(Assessment only)

the task; for the labyrinth this index includes information on the number of hits with the labyrinth wall; index value of 10 indicates an approximate performance value for a healthy person during a normal execution of the task.

Ratio between the length of an ideal path and the actually completed path while executing the task; the parameter indicates which trajectory the patient chose for completing the move (for example, during the reaching task the ideal path is a straight line between the initial and the final point on the trajectory, during the tracking task the ideal path is the circle, in the labyrinth the ideal trajectory is defined with the midline through the labyrinth path); value of 10 indicates the tracking of the ideal path that would be expected from an average healthy person.

[0-10]
(Assessment only)

Movement quality expressed as the average movement speed (actual travelled path normalized by the movement time); value of 10 indicates speed that would be expected from an average healthy person performing the task with a normal speed.

[0-10]
(Assessment only)

Error is defined in terms of deviations from the ideal trajectory of the move; in case of labyrinth it also includes information about how much path was completed inside the labyrinth wall; value of 10 indicates the error from the ideal trajectory that would be expected from an average healthy person.

[0-10]
(Assessment only)

Smoothness defines how smoothly or continuously was the movement completed; movements of healthy subjects are smooth; patient's movements are composed of accelerations and decelerations; value of 10 indicates an approximate smoothness value for an average healthy person.

[0-10]
(Assessment only)

Index of support provided by the less affected arm to the more affected arm; support value of 10 means that the affected arm was able to complete the task without the support of the less affected arm; value close to 0 indicates that the affected arm behaved passively and the less affected arm contributed to motion.

[s]
(Assessment only)

Average movement time achieved within the therapy tasks (for example, time required to complete a single move from the center to the outward position in a reaching task or a time to complete one trajectory in a tracking task).

[s]
(Assessment & Therapy)

Duration of an individual therapy task (used for computation of normalized score).

[m]
(Assessment & Therapy)

Range of motion in the selected plane of interest for therapy.

(Therapy only)

Highest achieved level in multi-level therapy tasks.

(Assessment & Therapy)

Clinically evaluated score value 1.

(Assessment & Therapy)

Clinically evaluated score value 2.

(Assessment & Therapy)

Clinically evaluated score value 3.


Chapter 4 Clinically evaluated score value 4. (Assessment & Therapy)

All therapy related settings are defined using the five buttons on the left side of the **Therapy** window. These five buttons allow the setting of:

- A. Therapy modality (unimanual or bimanual therapy) (A),
- B. Movement modality (B),
- C. Assessment scenarios (C),
- D. Therapy scenarios (D),
- E. Select a therapy program from a list of predefined therapy programs (E).

The same five buttons are used for configuring a therapy task to be added to a therapy program (see chapter [Therapy task configuration](#)).

The same five buttons are used for selecting a therapy task for which data is to be analyzed (see chapter [Therapy task selection](#)).

The same five buttons are used for selecting a therapy task for which data is to be analyzed (see chapter [Therapy task selection](#)).

THERAPY MODALITY

The Bimeo PRO system enables unimanual and bimanual therapy. The **Therapy modality selection** (Figure 34) enables switching between the two therapy approaches. This selection influences combination of attachments to be used with the Bimeo master, preferred movement modality selection, as well as the set of scenarios available for therapy.


Figure 34 Manual selection mode

The bimanual mode enables therapy with both arms in free space or on a support surface. Attachments to be used for bimanual therapy are shown in Figure 35 and Figure 36. The two configurations for bimanual therapy are shown in Figure 37 and Figure 38.


Figure 39 Components of the unimanual therapy attachment


Figure 40 Components of the unimanual therapy attachment

Figure 41 Components of the unimanual therapy attachment

Figure 42 Components of the unimanual therapy attachment


Figure 43 Components of the unimanual therapy attachment

The unimanual mode enables therapy with a single arm in free space, either on a support surface or on a support sphere (wrist therapy). Attachments to be used for unimanual therapy are shown in Figure 39 and Figure 40. The three configurations for unimanual therapy are shown in Figure 41, Figure 42 and Figure 43.

Figure 41 Components of the unimanual therapy attachment


Figure 40 Attachments for unimanual therapy on a sphere

Figure 42 Components of the unimanual therapy attachment

Figure 43 Components of the unimanual therapy attachment

Figure 41 Components of the unimanual therapy attachment

The Bimeo PRO system enables therapy in a three-dimensional space. However, to simplify representation of information on a two-dimensional screen, preferred movement modalities were defined and can be selected from the **Movement modality menu**.


Figure 44 Preferred movement modality selection

In **Bimanual therapy** mode three preferred movement modalities for therapy are available

1. **Vertical plane** mode (Figure 45) – preferred movements are parallel to the subject's coronal plane at any distance from the body; the more affected arm is supported by the activity of the less affected arm,
2. **Horizontal plane** mode (Figure 46) – preferred movements are parallel to the subject's transverse plane at any height; the more affected arm is supported by the activity of the less affected arm,
3. **Support surface** (Figure 47) – preferred movements are constrained to the support surface parallel to the subject's transverse plane at any height; the more affected arm is supported by the activity of the less affected arm.

Figure 45 Bimanual vertical plane mode

Figure 46 Bimanual horizontal plane mode


Figure 47 Bimanual support surface mode

In **Unimanual therapy** mode five preferred movement modalities for therapy are available

1. **Vertical plane** mode (Figure 48) – preferred movements are parallel to the subject's coronal plane at any distance from the body,
2. **Horizontal plane** mode (Figure 49) – preferred movements are parallel to the subject's transverse plane at any height,
3. **Support surface** mode (Figure 50) – preferred movements are constrained to the support surface parallel to the subject's transverse plane at any height,
4. **Support sphere** mode (Figure 51) – preferred movements are constrained to the support sphere attached to a fixed surface; the therapy focuses on activity of wrist joint,
5. **Joint mode** (Figure 52) – therapy focuses on activation of individual arm joints (Figure 53)
 - a. Shoulder flexion/extension
 - b. Shoulder abduction/adduction
 - c. Shoulder internal/external rotation
 - d. Elbow flexion/extension
 - e. Wrist pronation/supination
 - f. Wrist ulnar/radial deviation
 - g. Wrist flexion/extension


Figure 48 Unimanual vertical plane mode


Figure 49 Unimanual horizontal plane mode


Figure 50 Unimanual support surface mode


Figure 51 Unimanual support sphere mode


Figure 52 Unimanual individual joint mode


Figure 53 Selection of different joint modes

The measured arm postures and forces are inputs to a virtual environment. The subject interacts with the virtual environment by manipulating or moving the Bimeo master. Virtual reality scenarios are suitable for unimanual as well as bimanual modes. Though the Bimeo PRO device enables movement in full three-dimensional space, the therapy tasks are limited to horizontal and vertical planes in order to simplify presentation of information on a two-dimensional screen.

Depending on the selected scenario and the device configuration, not all therapy and assessment scenarios are available. The availability of therapy task depends on the possibility of measuring arm postures and interaction forces.

Therapy and assessment scenarios can be selected from the main menu (Figure 54 and Figure 55).



Figure 55 Therapy scenarios

If therapy programs are predefined for the specific patient and the specific arm (see paragraph above), they can be selected from the **Therapy Program** submenu. In the submenu all available therapy programs for the selected patient are listed. If therapy programs are defined for left and right arm and right Bimeo PRO system is initialized, only therapy programs for the right arm will be displayed (and vice versa if the left Bimeo PRO system is initialized). Program names appear as they were defined during the configuration of the therapy program.

Once the therapy program is selected therapy tasks are presented to the patient in the order specified in the therapy program (if random sequence was selected then the order of therapy tasks is randomly determined). Once the therapy task is completed, the configuration for the next therapy task will be automatically configured. If autoplay is disabled each therapy task needs to be started manually (see section **Therapy Task Configuration** below). The therapy program is terminated once all preconfigured tasks are completed or the **EXIT** button in the **Therapy Program** menu is pressed.

The current therapy configuration settings are represented by three images displayed in the middle of the **Therapy Program** submenu (Figure 57). The three images visually demonstrate what was selected from the **Therapy Program** submenu (see above)

- A. The first image (A) displays the required configuration of the Bimeo master; the required attachment needs to be fitted onto the Bimeo master before the therapy starts,
- B. The second image (B) illustrates the selected preferred movement modality; this gives the patient a visual representation of required movements so that the patient can move the arm in the required posture before the therapy starts,

- C. The third image (C) illustrates the selected therapy or assessment scenario; the patient can mentally prepare for the therapy task before the therapy starts.

Images are automatically replaced when therapy configuration settings are changed manually or from within the therapy program.


Figure 57. The three images (A, B, C) illustrate the selected therapy or assessment scenario; the patient can mentally prepare for the therapy task before the therapy starts.

Figure 58. The control interface.

The four buttons in the middle of the control interface are the control interface for the therapy (Figure 58).

- The button (A) is the interface to the initialization window (see the initialization window section below); Bimeo PRO system initialization is required upon the start of the Bimeo application before the first therapy session can be conducted; the system can be reinitialized several times during the running of the Bimeo application, however only one initialization is strictly required,
- The button (B) is the interface to the range of motion window (see the range of motion window section below); the button is activated after the system has been initialized,
- The button (C). The therapy time for each task can be set with the button. The task time is shown below the button.
- The button (D) starts the selected therapy task; the button is NOT active until the therapy parameters are properly defined; clicking on the button opens the task instructions (see the task instructions section and individual scenario instructions sections below).


Figure 58. The four buttons in the middle of the control interface are the control interface for the therapy. The button (A) is the interface to the initialization window; the button (B) is the interface to the range of motion window; the button (C) is the interface to the duration window; the button (D) is the interface to the task instructions window.

Therapy settings status bar indicates the correctness of therapy settings (Figure 59). Seven configuration states are indicated either with a green checkmark or with a red cross:

- A. **Right arm status** (A) – if no Bimeo master is connected, a red cross is displayed; if system is connected, **Right arm** or **Left arm** is displayed indicating the status of the initialized arm,
- B. **System initialized** (B) indicates whether the system is initialized (green checkmark) or not (red cross),
- C. **Patient** (C) indicates whether a new patient or an existing patient is selected (green checkmark) for therapy or not (red cross),
- D. **Therapy modality** (D) indicates whether bimanual or unimanual therapy is selected (green checkmark); if neither is selected a red cross is displayed ,
- E. **Movement modality** (E) indicates whether a preferred movement direction is selected (green checkmark) or not (red cross),
- F. **Assessment** (F) indicates that an assessment scenario is selected (green checkmark) or not (red cross),
- G. **Therapy** (G) indicates that a therapy scenario is selected (green checkmark) or not (red cross),

In general all indicators need to be marked with a green checkmark before the therapy session can start. The only exceptions are **Therapy modality** and **Movement modality** that indicate two alternative options; only one at a time can be marked with a green checkmark, the other is marked with a red cross.

The Bimeo PRO system consists of three sensors to be attached to the patient's arm. These sensors measure relative arm orientation to the patient's trunk. In order to keep the system simple there is no sensor attached to the torso. Therefore an initialization procedure is required to determine the approximate torso orientation in space before the therapy can start. This is done by the initialization procedure. **Proper system initialization before the start of the therapy is essential for correct performance of the Bimeo PRO system.**

1. The computer screen (laptop) should be placed on a flat surface (table) approximately 1 m away from the patient.
2. The flat surface should **NOT contain ferromagnetic** components in order not to deteriorate the performance of the sensing devices.
3. The patient should be seated behind the flat surface with his/her frontal plane aligned with the plane of the computer screen. This will simplify the interpretation of information presented on the computer screen.
4. The forearm and upper arm Bimeo sensors should be attached to the patient's arm. The attachment of the two sensors on the subject's arm is not strictly required for initialization procedure, nevertheless it is recommended.
5. All Bimeo sensors need to be turned on (Figure 60).
6. **IMPORTANT: Bimeo master must be put on the flat surface in horizontal position and the finger grasping part of the Bimeo master facing toward the computer screen (Figure 61). The Bimeo master should be aligned with the frontal plane of the patient. Any attachment can be attached to the Bimeo master, not necessarily the one shown in Figure 61.**
7. On the **Initialization** window click the **Turn on sensors** button (Figure 58A).
8. If not all Bimeo sensors are turned on, **Initialization** window will appear instructing the user to turn the Bimeo system on (Figure 60). In the Initialization window sensors status is always shown. All sensors need to be marked with a green checkmark. If a red cross is displayed next to the sensor image, a sensor is either turned off or malfunctioning. The click on the **Turn on sensors** button leads to the **Initialization** window without completing the initialization.
9. When all Bimeo sensors are turned on the system will automatically proceed to **Initialization** window as shown in Figure 61. Follow the instructions for positioning the Bimeo master for initialization. Select the desired arm for therapy – left arm or right arm. The system will be initialized for selected arm. If you want to change the arm, repeat the initialization procedure.
10. The initialization starts and takes a few seconds to complete.
11. When the initialization procedure is completed the system returns automatically to the **Initialization** window.
12. The system can be reinitialized several times during the therapy session. However, only one initialization procedure is strictly required.
13. The initialization procedure is necessary every time when Bimeo software is started or when the training arm is changed.


Figure 60 Initialization window: instructions for turning the Bimeo PRO system on


Figure 61 Initialization window: instructions for positioning the Bimeo master as shown in the figure; the Bimeo master needs to be oriented in a way that the finger grasping part of the Bimeo master is facing the direction of the computer screen; the Bimeo master needs to be parallel with the frontal plane of the patient.

	WARNING The surroundings of the patient workspace should NOT CONTAIN FERROMAGNETIC components in order not to deteriorate the performance of the sensing devices.
	WARNING Do not place the sensing devices (Bimeo master, Bimeo sense) close to PERMANENT MAGNETS or devices that could produce HIGH MAGNETIC FIELD in order not to magnetize the sensing devices. This also includes laptops, loudspeakers and fMRI.
	ATTENTION In the event of exposure to HIGH MAGNETIC FIELD, the device must be calibrated. Contact your dealer or the manufacturer about calibrating the device.

Range-of-motion settings


Range-of-motion is a critical therapy related parameter that needs to be defined for each patient, specifically for bimanual and unimanual therapy and specifically for each preferred movement modality. Therefore, settings for each therapy modality are stored into the database. Before the new therapy session range-of-motion values are retrieved from the database and do not need to be redefined.

Range-of-motion is defined as a rectangle in the preferred plane of motion (Figure 62) or as the lower and upper limit of the joint motion (Figure 63). **All therapy and assessment scenarios are confined to the specified range-of-motion and no movements are required beyond the specified range-of-motion.**

The range-of-motion setup window can be accessed directly from the **Main user interface window** (Figure 58) or from the **Scenario user interface window** (Figure 64).

The range-of-motion for the preferred movement plane can be defined using the **Range-of-motion setting window** as shown in Figure 62. There are several options to simplify the definition of patient-specific range-of-motion

- A. Before setting the range-of motion the preferred central position of the arm needs to be defined. Range-of-motion is then defined relative to this central point. The central point is defined by placing the patient's arm in the preferred position and clicking on the white circular button marked with a blue cross (A),
- B. After the central point has been defined the appropriate range-of-motion can simply be defined by dragging the mouse while holding the left mouse button over the white circular button with blue arrows at the top-right corner of the rectangle indicating the range-of-motion area (B) (the button is blinking when this mode is active),
- C. The current extent of the range-of-motion is displayed as an approximate value in centimeters in the **Range of motion** window at the top-right corner of the screen (C),
- D. The white rectangle with currently selected scenario in background indicates the approximate range of motion,

[illegible]

- A. The lower value for the range-of-motion can be defined by clicking (holding the left mouse button) and pulling the arrow (A); the lower value for the range-of-motion cannot exceed the upper value for the range-of-motion,
- B. The upper value for the range-of-motion can be defined by clicking (holding the left mouse button) and pulling the arrow (B); the lower value for the range-of-motion cannot exceed the upper value for the range-of-motion,
- C. The range-of-motion can be defined manually by rotating the joint to the lower and upper value after clicking the hand button (C); the automatic setting of the range-of-

- motion can be combined with manual moving of arrows (A) and (B) (for example if the range-of-motion needs to be reduced),
- D. The current orientation of the joint is represented by a red rotating bar (D),
 - E. The setting circle provides an approximate value of range-of-motion in degrees (E),
 - F. Range-of-motion limits are displayed in the upper-left part of the screen (F),
 - G. By clicking the **CLOSE** button (H) the setting of range-of-motion is terminated and the application returns either to the **Main user interface window** or to **Scenario user interface window** depending on where the range-of-motion setup was triggered.


Figure 63 Joint range-of-motion settings window: A) setting for lower joint limit, B) setting for upper joint limit, C) button (hand) for setting the joint range-of-motion by moving the patient's arm, D) display of joint range-of-motion settings, E) circle with degrees indicating range-of-motion, F) upper and lower joint limit, and G) button for closing the joint range-of-motion setup window.

General scenario user interface

All Bimeo PRO therapy scenarios follow the same concept in terms of user interface. The general user interface is therefore described in this section, while individual therapy scenarios are described in their respective sections below. The general scenario user interface is shown in Figure 64 and consists of the following elements:

- A. The central part of the scenario user interface (A) is reserved for the presentation of the scenario. Different scenarios are available for therapy and assessment (see details about therapy and assessment scenario in the sections below),
- B. **Left force bar** (B) indicates activity of the left arm in relation to activity of the right arm; if the left arm is more active in terms of force (power) than the right arm, the left force bar shows the level of activity; if the left arm is the affected arm holding the Bimeo master, the left bar is colored green when active (Figure 65, green color

indicates correct activation when affected arm generates more power than the unaffected arm); if the left arm is the unaffected arm, the left bar is colored red when active (red color indicates prevailing activation of the unaffected arm); force bars can be turned on or off within the scenario settings parameters and are just available in Bimanual mode in vertical or horizontal plane (Figure 67).

- C. **Force bar (C)** indicates activity of the right arm in relation to activity of the left arm; if the right arm is more active in terms of force (power) than the left arm, the right force bar shows the level of activity; if the right arm is the affected arm holding the Bimeo master, the right bar is colored green when active (green color indicates correct activation when affected arm generates more power than the unaffected arm); if the right arm is the unaffected arm, the right bar is colored red when active (Figure 66, red color indicates prevailing activation of the unaffected arm), force bars can be turned on or off within the scenario settings parameters and are just available in Bimanual mode in vertical or horizontal plane (Figure 67),
- D. **Top bar (D)** indicates elapsed time since the start of the therapy session (minutes : seconds), current scenario level and achieved scores during the therapy session,
- E. **Reset button (E)** is used for resetting the center position of the arm workspace during the therapy; the resetting is done by positioning the arm in the preferred central position and clicking the **Reset** button,
- F. **Range-of-motion button (F)** provides the interface to the range-of-motion setting window that enables redefinition of the preferred arm range-of-motion (see section **Range-of-motion** above),
- G. **Settings button (G)** opens the popup window for setting the scenario-specific parameters (Figure 67); these include, for example, therapy duration, force-bars visibility, scenario difficulty, etc.,
- H. **Terminate button (H)** terminates the scenario and returns to the **Therapy** interface
- I. **Assessment display bar (I)** display bar provides information about the quality of movement during the therapy in terms of **Assessment scores** (Table 2); the display is only available for assessment scenarios (see **Therapy scenarios** section above); the display can be switched between the assessment scores for the last move, the median assessment scores for the entire therapy session and no display (default),
- J. **Assessment display button (J)** for switching between different real-time assessment display options.

Scenario parameters, including range-of-motion, therapy duration and difficulty, are stored for each individual patient and are retrieved from the database before the next therapy session. The range-of-motion parameters are stored specifically for each therapy modality (unimanual or bimanual) and preferred movement modality.

The two vertical bars in Figure 64 indicate relative forces applied by the unaffected and/or affected arms. When the subject applies force in the direction of movement or against gravity with the affected arm (thus supporting the unaffected arm), the bar on the affected side turns green with the bar height proportional to the applied force. On the other hand, when the unaffected arm supports the movement and compensates the weight of the affected arm, the bar on the unaffected side turns red. The two bars thus provide cognitive feedback to the patient about the activity of the affected arm.


Figure 64 Scenario user interface window: A) central part with scenario, B) left force bar, C) right force bar, D) elapsed time since the start of the therapy session, level and achieved scores, E) set center position button, F) range-of-motion-settings button, G) scenario settings button, H) exit button, I) performance indices real-time display (only assessment scenarios) and J) toggle between different performance indices real-time displays (only assessment scenarios).


Figure 65 Green bar on the affected side indicates higher activity of the affected limb


Figure 66 Red bar on the non-affected side indicates higher activity of the non-affected limb

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The main value of assessment scenarios is the objective assessment of patient's motor capabilities. Nevertheless, these scenarios can also be used for therapy.

The scenario promotes straight line movements in the selected plane. Movements begin in a starting position in the center of the screen, which corresponds to a predefined starting posture of the arm. Starting position is marked with a green circular object. Eight virtual objects – smileys – are appearing along a circular pattern with a predefined radius. Only one smiley is displayed for each movement. Smileys are displayed in a random sequence and are evenly distributed along the circular pattern. When the patient completes the movement from the center to the target smiley, the smiley disappears and the patient needs to return to the center position. The movement is then repeated with a new smiley. Center position and circular pattern radius can be customized to patient's needs.

- Patient gets 10 points for every successful move.

- Smiley appearance position is selected by pressing **Smiley** button. Following options can be chosen
 - All directions - smileys appearing randomly around the entire circle,
 - Upper half - smileys appearing randomly on the upper half of the circle,
 - Lower half - smileys appearing randomly on the lower half of the circle,
 - Left half - smileys appearing randomly on the left half of the circle,
 - Right half - smileys appearing randomly on the right half of the circle,
 - Upper right quarter – smileys appearing randomly on the upper right quarter of the circle.
- **Waiting time** is defined as the maximum time in which the patient's movement has to be completed. If the patient does not complete movement inside **Waiting time** interval, a new smiley appears and the patients need to return back to center position and start new movement. Waiting time can be adjusted from 1 s to 20 s.
- The vertical bars, which indicate relative forces applied by the unaffected and/or affected arms, can be enabled or disabled by pressing **Vertical bars** button only in bimanual therapy in vertical plane or horizontal plane.
- Therapy duration for selected scenario can be defined by changing **Therapy duration**. **Therapy duration** can be adjusted from 30 s to 20 min.


Figure 68 Reaching


Figure 69 Reaching settings

TRACKING

The scenario promotes circular motion in the selected plane. The scenario requires tracking of a reference point (the target) that moves on a round trajectory with a predefined speed. Speed and direction of tracking can be adjusted according to the patient's needs. The task is scored according to the tracking accuracy in relation to the center of the target.

Scoring

Patient gets points every second:

- 6 points for position in the center of the target,
- 4 points for position to the middle of the target,
- 2 points for position to the edge of the target and
- No points for position outside of the target.

Settings

- **Speed** of the target can be adjusted from 0.1 to 2 where 2 corresponds to the maximum speed of the target and 0.1 corresponds to the minimum speed of the target.
- Direction of the target motion can be changed by pressing **CHANGE DIRECTION** button.
- The two vertical bars which indicate relative forces applied by the unaffected and/or affected arms can be enabled or disabled by pressing **SHOW FORCE BARS** button.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 70 Tracking


Figure 71 Tracking settings

LABYRINTH

The scenario focuses on complex movements in the selected plane. The goal of the scenario is to complete the labyrinth path as quickly as possible with the least possible number of contacts with the labyrinth wall. When the path is completed, a new labyrinth is generated and the scenario is repeated. Path width and complexity are customizable scenario parameters.

Scoring

- Patient gets 10 points for every labyrinth path completed.
- Extra 20 points are earned completing without the contact with the labyrinth wall.

Settings

- Labyrinth difficulty level can be set by changing **Difficulty** parameter. Three options are available: easy, medium and difficult. We recommend starting with an easy level and then adjusting it according to the patient's ability.
- The two vertical bars which indicate relative forces applied by the unaffected and/or affected arms can be enabled or disabled by pressing **SHOW FORCE BARS** button.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 72 Labyrinth


Figure 73 Labyrinth settings

ROTATION

The scenario promotes movements in individual joints. In contrast to other scenarios, this task does not require multi-joint coordination. It rather focuses on therapy and assessment of isolated movements in individual joints. The green line corresponds to the actual joint angle, while blue lines indicate customizable joint rotation limits. The subject must rotate the limb segment alternately in the two target orientations (blue lines), which are pre-set according to abilities of the exercising person.

Scoring

- Patient gets 5 points for every successful move.

Settings

- The horizontal bar which indicates relative joint torques applied by the affected arm can be enabled or disabled by pressing **SHOW FORCE BARS** button.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 74 Rotation


Figure 75 Rotation settings

The therapy scenarios are divided into:

- activities-of-daily-living tasks,
- cognitive tasks and
- games.

The main focus of therapy scenarios is performing the activities-of-daily-living, cognitive tasks or just playing fun games. Due to high variability of required movements, these scenarios are not suitable for motor performance assessment. Nevertheless, basic assessment parameters are available also for these scenarios.

Kitchen scenario mimics an activity of daily living. By following recipes displayed on the screen the task is to use multiple ingredients in a proper order to prepare different dishes. First, the patient needs to place hands to the position where ingredients are stored. Grasping in bimanual mode in the vertical or horizontal plane is performed by applying force on the Bimeo master. Automatic grasping is enabled in unimanual therapy and bimanual therapy on support surface. Then the ingredient needs to be brought onto the kitchen table and poured into the pot. The task consists of several difficulty levels. While preparing dishes the insects need to be avoided from level two.

- Patient gets 10 points for every successful performed instruction.
- If patient completes the recipe, he/she gets 100 points.
- If patient touches the insect, he/she loses 20 points.
- If patient hits the insect with the flyswatter, he/she gets 50 points.

- Difficulty level of the scenario can be set manually by adjusting **difficulty level** parameter. The level is automatically increased when the certain score is reached.
- Amount of grasping force can be adjusted by changing **grasping force** parameter. **grasping force** parameter can be set from 3 N to 20 N.
- Therapy duration for selected scenario can be defined by changing **therapy duration** parameter. **therapy duration** can be adjusted from 30 s to 20 min. When the therapy duration is over, the patient has still 1 minute to finish the last recipe. After 1 minute, the scenario is finished.


Figure 76 Kitchen

Figure 77 Kitchen settings

FRUIT-PICKER (ADL TASK)

The scenario represents an activity of daily living. Apples need to be collected from the tree and put into the basket. Hands first need to be positioned over an apple. Grasping in bimanual mode in the vertical or horizontal plane is performed by applying force on the Bimeo master. Then the apple is brought into the basket by moving the arms and simultaneously applying force on the Bimeo master. At the same time birds and worms need to be avoided. Automatic grasping is enabled in unimanual therapy and bimanual therapy on support surface. The scenario consists of several difficulty levels.

Scoring

- Patient gets 10 points for every successful placed apple in the basket.
- If the animal touches the apple, the patient loses 20 points.

Settings

- Difficulty level of the scenario can be set manually by adjusting **Level** parameter. The level is automatically increased when the certain score is reached.
- Amount of grasping force can be adjusted by changing **Squeezing force** parameter. **Squeezing force** parameter can be set from 3 N to 20 N.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min. When the therapy duration is over, the patient has still 1 minute to finish the last tree. After 1 minute, the scenario is finished.


Figure 78 Fruit-picker



Figure 79 Fruit-picker settings

The scenario represents a cognitive task. A picture with the missing part is presented to the patient. There are four different parts of the same picture on the left and right side. The patient needs to find the part that is missing in the picture. The patient selects the part by moving the cursor to the missing part and holding the cursor on the missing part for two seconds. If the correct part is selected, the whole picture is shown. The scenario consists of two difficulty levels.

- Patient gets 10 points for every correct answer.

- Puzzle scenario has two difficulties levels. Difficulty level can be selected by changing **Difficulty** parameter. Options are easy and hard.
- Therapy duration for selected scenario can be defined by changing **Duration** parameter. **Duration** can be adjusted from 30 s to 20 min.

The scenario represents a cognitive task. Several cards that are faced down are presented to the patient. The patient needs to find the pair of cards with the same figure. The patient turns the card by moving the cursor to the selected card and holding the cursor on the card for two seconds. The same action is required to turn the second card. If the pair is found, the cards are removed. Otherwise, the cards are turned back. The goal of the scenario is to clear all the cards. The scenario consists of three difficulty levels: four pairs of cards, six pairs of cards and eight pairs of cards.

- Patient gets 10 points for every found pair.
- Patient gets 50 points when all pairs are found.

- Memory scenario has three difficulties levels. Difficulty level can be selected by changing **Difficulty** parameter. Options are easy, medium and hard.

- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 82 Memory

Figure 83 Memory settings

SLIDE (COGNITIVE TASK)

The scenario represents a cognitive task. Several parts of the picture are presented to the patient. There is also an empty space with grey color. The patient needs to move the parts to compose the picture. Only the parts next to the grey space can be moved. The patient moves the selected part by moving the cursor to the part and holding the cursor on the part for one second. The selected part is then moved to the empty space. When all the parts are in their correct positions, the whole picture is shown.

Scoring

- Patient gets 5 points for every moved part.
- Patient gets 50 points when the picture is composed (easy difficulty).
- Patient gets 300 points when the picture is composed (medium difficulty).
- Patient gets 800 points when the picture is composed (hard difficulty).

Settings

- Slide scenario has three difficulties levels. Difficulty level can be selected by changing **Difficulty** parameter. Options are easy, medium and hard.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.

Figure 84 Slide


Figure 85 Slide settings

GHOST SCENARIO (GAME)

The scenario represents an abstract game with several difficulty levels. The cat runs into the tree stump and the ghost has to save the cat from dogs that are tracking it. The user moves the red cursor and the ghost follows the cursor. To scare the dog, move the ghost in front of the dog. Those dogs who wear sun glasses or are hiding behind the bush cannot see the ghost, so the ghost must perform a scary sound. The user must squeeze the Bimeo device to apply the sound in bimanual mode in vertical or horizontal plane. The game consists of several difficulty levels. There is just one dog at a time in first two levels while there are up to three dogs at a time in higher levels.

Scoring

- Patient gets 10 points for every scared dog.
- Patient loses 5 points if the dog barks at the cat.

Settings

- Difficulty level of the scenario can be set manually by adjusting **Level** parameter. Levels range from 1 to 10, level 1 being the easiest and level 10 the most difficult. The level is automatically increased when the certain score is reached.
- Amount of squeezing force can be adjusted by changing **Squeezing force** parameter. **Squeezing force** parameter can be set from 3 N to 20 N.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 86 Ghost


Figure 87 Ghost settings

SPACE (GAME)

Space is an abstract two-dimensional arcade game. The aim is to shoot or crash meteorites with an in-game spaceship (Figure 88). The spaceship position and orientation is controlled by the patient's hands movement. Shooting is performed by applying force on the Bimeo master with both of the patient's hands. Shooting is disabled in unimanual therapy and in bimanual therapy on support surface. The patient can destroy the meteorites by hitting them with the spaceship.

Scoring

- 10 points are earned for every meteorite shot down.

- 5 points are earned if the space ship collides with a meteorite.

Settings

- Difficulty level of the game can be set manually by adjusting **Level** parameter. Levels range from 1 to 11, level 1 being the easiest and level 11 the most difficult. The level is automatically increased when the certain score is reached.
- Amount of shooting force can be adjusted by changing **Squeezing force** parameter. **Squeezing force** parameter can be set from 3 N to 20 N.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 88 Space


Figure 89 Space settings

TIC TOC (GAME)

Tic toc is a classic two-dimensional sport game. It simulates table tennis in which a patient controls an in-game paddle by moving it horizontally across the bottom of the screen. Patients use the paddles to hit a ball back and forth and compete against a computer controlled opponent. The aim is for a patient to earn more points than the opponent.

Scoring

- 50 points are earned when the opponent fails to return the ball.
- 10 points are lost when the patient fails to return the ball.

Settings

- Tic toc game has three difficulties levels. Difficulty level can be selected by changing **Difficulty** parameter. Options are easy, medium and hard.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 90 Tic toc


Figure 91 Tic toc settings

ROCKET (GAME)

The game represents an abstract game with several difficulty levels. The user moves the rocket between planets. The user must deliver the packages from the central planet to the selected outer planet. To move the package with the rocket, the patient has to apply force on the Bimeo while moving in bimanual mode in vertical or horizontal plane. In unimanual and bimanual on support surface mode auto grasping is enabled. The game consists of several difficulty levels. The user has to avoid the comets that comes in the workspace. The enemy alien shall not steal with the package. In higher levels, the user can change between red and blue rocket. The red rocket can deliver red and brown packages while the blue rocket can deliver the blue and brown packages. The rocket can be switched by moving to the change rocket sign at the top and applying force on the Bimeo master. Automatic squeezing is enabled in unimanual therapy modality.

Scoring

- Patient gets 20 points for every delivered package.
- Patient loses 10 points if the enemy takes the package.
- Patient loses 10 points if the comet hit the rocket.

Settings

- Difficulty level of the game can be set manually by adjusting **Level** parameter. Levels range from 1 to 9, level 1 being the easiest and level 9 the most difficult. The level is automatically increased when the certain score is reached.
- Amount of squeezing force can be adjusted by changing **Squeezing force** parameter. **Squeezing force** parameter can be set form 3 N to 20 N.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 92 Rocket


Figure 93 Rocket settings

RACER (GAME)

Racer is a classic action racing game. The patient controls the car by moving it horizontally across the bottom of the screen. Patient needs to avoid the obstacles and collect the coins. The game consists of several difficulty levels with different speed of the car.

Scoring

- 20 points are earned when the coin is collected.
- 10 points are lost when the patient hit the obstacle.

Settings

- Racer game has three difficulties levels. Difficulty level can be selected by changing **Difficulty** parameter. Options are easy, medium and hard.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 94 Racer scenario


Figure 95 Racer scenario settings

SNOWBALL (GAME)

The game represents an abstract game with several difficulty levels. The user must hit the dwarfs and animals with the snowball. To throw the snowball, the patient has to apply force on the Bimeo while moving in bimanual therapy in vertical and horizontal plane. Automatic throwing is enabled in unimanual therapy and bimanual therapy on support surface.

Scoring

- Patient gets 5 points for every dwarf's hat hit.
- Patient gets 10 points for every dwarf hit.
- Patient gets 20 points for every snowballs pile hit.
- Patient gets 30 points for every rabbit hit.
- Patient gets 40 points for every deer hit.
- Patient can get 50 points bonus for house hit.

Settings

- Difficulty level of the game can be set manually by adjusting **Level** parameter. Levels range from 1 to 8, level 1 being the easiest and level 8 the most difficult.
- Amount of squeezing force can be adjusted by changing **Squeezing force** parameter. **Squeezing force** parameter can be set from 3 N to 20 N. The level is automatically increased when the certain score is reached.
- Therapy duration for selected scenario can be defined by changing **Duration**. **Duration** can be adjusted from 30 s to 20 min.


Figure 96 Snowball


Figure 97 Snowball settings

Delrus Ltd.
Kutuzovskaya 43
Moscow, Russia
125080

info@delrus.ru
www.delrus.ru

Russian authorized representative

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