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LH Robinson

11 March 2020

Lindsay Helen Robinson
Solicitor and Notary Public

DWF LLP
No 2 Lochrin Square
96 Fountainbridge
Edinburgh
EH3 9QA
DX ED15

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«APPROVED»

Alba Bioscience Ltd.
James Hamilton Way,
Penicuik, EH26 0BF,
United Kingdom

Date: 11 March 2020

Name: Janine Robb

Janine Robb
(signature)

Position: Head of Technical Support

Stamp:

Alba Bioscience Limited
doing business as
Quotient

OPERATIONAL DOCUMENTATION
on medical device

**Quality Control Kit «Ortho Confidence WB» - erythrocytes and
antibodies to erythrocyte antigens**

Manufacturer:

Alba Bioscience Ltd. James Hamilton Way, Penicuik, EH26 0BF, United Kingdom

2020

I. GENERAL INFORMATION

1. Name of medical device

Quality Control Kit «Ortho Confidence WB» - erythrocytes and antibodies to erythrocyte antigens

Composition:

- «Ortho Confidence WB 1» Quality control - 6,5 ml in vial with white cover – 1 unit.
- «Ortho Confidence WB 2» Quality control - 6,5 ml in vial with red cover – 1 unit.
- «Ortho Confidence WB 3» Quality control - 6,5 ml in vial with blue cover – 1 unit.
- «Ortho Confidence WB 4» Quality control - 6,5 ml in vial with black cover – 1 unit.
- Information insert.

Short names that will be used thereafter in the text:

Quality Control Kit ORTHO CONFIDENCE WB

The manufacturer's abbreviation for the test: ORTHO CONFIDENCE WB

2. Name and address of manufacturer

(Alba Bioscience Ltd. James Hamilton Way, Penicuik, EH26 0BF, United Kingdom)

II. INSTRUCTION FOR USE.

See Attachment.

(Attachment 1
to Operational Documentation)

Instruction
On
Quality Control Kit «Ortho Confidence WB»

TRANSFUSION MEDICINE

INSTRUCTION

Quality Control Kit «Ortho Confidence WB» - erythrocytes and antibodies to erythrocyte antigens

REF

6842785

Intended Use

For *in vitro* diagnostic use only

Quality Control Kit «Ortho Confidence WB» - erythrocytes and antibodies to erythrocyte antigens

Intended for use as qualitative controls to assure the reliability of basic blood bank serological reagents and automated/manual blood grouping systems using column agglutination techniques.

Summary and Explanation

The purpose of daily quality assurance in the blood bank is to confirm the reliability of the test system. The test system includes reagents, test procedures and equipment. Testing known samples is an accepted method of quality control. If expected test results are observed, procedures are being performed accurately and reagents and equipment are performing properly. If unexpected results are observed, the problem may be due to improper test performance, faulty equipment or contamination or deterioration of reagents. The source of the problem should be determined and resolved before patient/donor test results are reported.

The ORTHO CONFIDENCE™ WB Quality Control Kit provides a means of confirming the reactivity of routinely used reagents and is to be tested on each day of use. Observation of expected test results with the ORTHO CONFIDENCE™ WB Quality Control Kit will confirm the reactivity of ABO, Rh, and K antibodies, as well as reverse grouping cells and reagent red blood cells used for antibody detection.

Principle of Procedure

The procedures used with these reagents are based on the principle of agglutination. Normal human red blood cells will agglutinate in the presence of the appropriate antibody directed against antigens on those red blood cells. No agglutination indicates the absence of the demonstrable antigen or antibody.

Reagents

The ORTHO CONFIDENCE™ WB Quality Control Kit is prepared from red blood cells collected from blood donors. ABO and anti-Fy^a antibodies are of monoclonal origin and anti-D and anti-c antibodies are of polyclonal origin. Polyclonal anti-D source material is prepared from an initial concentration of <10 IU/mL to a concentration of approximately 0.1 IU/mL.

The concentration of red blood cells in each of the samples is 15% ± 2%. The red cells are suspended in a preservative solution to retard haemolysis and bacterial contamination.

- Vial 1 - Group A₁B D+C+E+c+e+ (R₁R₂)
- Vial 2 - Group O D+C+E-c-e+ (R₁R₁) K+, anti-A, anti-B, anti-c
- Vial 3 - Group A D+C-E+c+e- (R₂R₂) Fy(a-), anti-B, anti-Fy^a
- Vial 4 - Group B D-C-E-c+e+ (rr) K-, anti-A, anti-D (concentration of approximately 0.1 IU/mL)

Storage Requirements

In refrigerators at a temperature of 2 to 8 ° C with a relative humidity of 20 to 80%.

Shelf life

Vials 1-4	Conditions	Stability
Unopened	Refrigerated, at 2-8 ° C	To expiration date. Shelf life is indicated on the package.
Opened	Refrigerated, at 2-8 ° C	14 days

Transportation

In refrigerators at a temperature of Shipping temperature: 2 to 8 ° C with a relative humidity of 20 to 80%.

Pub. No. J39796_RU_RU

Precautions

- Do not transfer any of the vial contents to another container as this could result in spillage or contamination.
- Do not use beyond the product expiry date.
- Once the product vials are opened, the reagent's performance characteristics are maintained for 14 days.
- Do not freeze.
- Do not use if marked haemolysis is observed. Slight discoloration in the supernatant is normal over the shelf life of the product. Extreme turbidity, precipitation or haemolysis of the red blood cells may indicate product alteration.
- Do not use if evidence of contamination is present.
- Replace vial caps when not in use.
- Allow the vials to warm to room temperature (18-25 °C) prior to use.

Caution: Source material from which this product is derived was found non-reactive for HBsAg, Anti-HIV 1/2 and Anti-HCV. No known test methods can offer assurance that products derived from human blood will not transmit infectious disease. Appropriate care should be taken in the use and disposal of this product.

Procedure

This reagent has been validated for use manually on the ORTHO BioVue® System, and for use on ORTHO AutoVue® System (Innova or Ultra), ORTHO VISION® Analyzer, and ORTHO VISION® Max Analyzer automated blood grouping test systems and therefore its suitability for use by other technologies cannot be guaranteed. Customers who choose to use commercial antisera and red cell reagents in an off-label manner must ensure the test method is appropriate by validating its intended use. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer.

The ORTHO CONFIDENCE™ WB Quality Control Kit is intended to simulate normal blood samples. The samples contained in the ORTHO CONFIDENCE™ WB Quality Control Kit must be used at room temperature (18-25 °C) and should be tested by following standard procedures in accordance with the applicable Instructions for Use for routine reagent use.

Vials are required to be centrifuged prior to first use. Centrifugation must be performed according to the user's standard laboratory practice for the centrifugation of patient samples. Subsequent upright, refrigerated storage eliminates the requirement for further centrifugation unless mixing has occurred.

Materials Provided

- «Ortho Confidence WB 1» Quality control - Quality control material containing erythrocytes of the A1B group (according to the ABO system); Rh phenotype erythrocytes D+C+E+c+e+ (R1R2); 6,5 ml in vial – 1 unit.
- «Ortho Confidence WB 2» Quality control - Quality control material containing erythrocytes of the group O, Rh phenotype erythrocytes D+C+E-c+e+ (R1R1), Kell positive erythrocytes - Kell-K+; antibodies (according to the ABO system): anti-A, anti-B; alloimmune anti-erythrocyte antibodies - anti-c; 6,5 ml in vial – 1 unit.
- «Ortho Confidence WB 3» Quality control - Quality control material containing erythrocytes of group A, Rh phenotype D+C-E+c+e- (R2R2) Fy(a-) negative erythrocytes; antibodies (according to the ABO system): anti-B, alloimmune anti-erythrocyte antibodies - anti-Fya; 6,5 ml in vial – 1 unit.
- «Ortho Confidence WB 4» Quality control Kell (K-) - Quality control material containing erythrocytes group-B; Rh phenotype D-C-E-c+e+ (rr), Kell-negative erythrocytes – Kell-K-; antibodies (according to the ABO system): anti-A, alloimmune anti-erythrocyte antibodies - anti-D (concentration of approximately 0.1 IU/mL); 6,5 ml in vial – 1 unit.
- Information insert sheet.

Additional Materials Required but not Provided

For Automated Testing:

- ORTHO BioVue® System Cassettes
- Reagent Red Blood Cells
- ORTHO AutoVue® System (Innova or Ultra)/ORTHO VISION® Analyzer/ORTHO VISION® Max Analyzer

For Manual Testing:

- ORTHO BioVue® System Cassettes
- Reagent Red Blood Cells
- Manual or Electronic Pipettes
- Pipette Tips
- Timer
- Test Tubes
- ORTHO™ Workstation/ORTHO BioVue® System Centrifuge and ORTHO BioVue® Incubator

TRANSFUSION MEDICINE

INSTRUCTION

Interpretation of Results

Interpretation of Results

The following table illustrates the expected results in tests with ORTHO CONFIDENCE™ WB Quality Control Kit and routine blood bank reagents.

Positive Result (in the table designated +): Agglutination of the red blood cells is a positive test result and indicates the presence of the corresponding antigen or antibody.

Negative Result (in the table designated 0): No agglutination of the red blood cells is a negative test result and indicates the corresponding antigen or antibody is not demonstrable.

Component of ORTHO CONFIDENCE™ WB Quality Control Kit	Reagent Under Test	Expected Test Results*
Vial 1	Anti-A	+
	Anti-B	+
	Anti-A,B	+
	Anti-D	+
	A ₁ cells	0
	A ₂ cells	0
	B cells	0
	O cells	0
	Screening cell 1 (H ₁ H ₁)	0
	Screening cell 2 (R ₂ R ₂)	0
	Screening cell 3 (rr)	0
	Anti-C	+
	Anti-E	+
	Anti-c	+
	Anti-e	+
Vial 2	Anti-A	0
	Anti-B	0
	Anti-A,B	0
	Anti-D	+
	A ₁ cells	+
	A ₂ cells	+
	B cells	+
	O cells	0
	Screening cell 1 (H ₁ H ₁)	0
	Screening cell 2 (R ₂ R ₂)	+
	Screening cell 3 (rr)	+
	Anti-C	+
	Anti-E	0
	Anti-c	0
	Anti-e	+
Vial 3	Anti-A	+
	Anti-B	0
	Anti-A,B	+
	Anti-D	+
	A ₁ cells	0
	A ₂ cells	0
	B cells	+
	O cells	0
	Screening cell 1 (H ₁ H ₁)	**
	Screening cell 2 (R ₂ R ₂)	**
	Screening cell 3 (rr)	**
	Anti-C	0
	Anti-E	+
	Anti-c	+
	Anti-e	0

* Discrepant results must be investigated further.

** Results depend upon antigen configuration.

Pub. No. J39796_RU_RU

TRANSFUSION MEDICINE INSTRUCTION

Quality Control

Component of ORTHO CONFIDENCE™ WB Quality Control Kit	Reagent Under Test	Expected Test Results*
Vial 4	Anti-A	0
	Anti-B	+
	Anti-A,B	+
	Anti-D	0
	A ₁ cells	+
	A ₂ cells	+
	B cells	0
	O cells	0
	Screening cell 1 (R ₁ R ₁)	+
	Screening cell 2 (R ₂ R ₂)	+
	Screening cell 3 (rr)	0
	Anti-C	0
	Anti-E	0
	Anti-c	+
	Anti-e	+
	Anti-K	0

* Discrepant results must be investigated further.

** Results depend upon antigen configuration.

Quality Control

This product is a quality control reagent and its satisfactory performance when used by the recommended techniques represents an adequate level of control.

Limitations of the Procedure

1. ORTHO CONFIDENCE™ WB Quality Control Kit is designed to be tested with reagents that have not been diluted.
2. Improper techniques may invalidate the results obtained with this product.
3. False positive or false negative results can occur due to contamination of test materials, improper reaction temperature, improper storage of materials, and omission of test reagents.
4. Individual laboratory procedures may affect the final reaction strength observed in tests performed with ORTHO CONFIDENCE™ WB Quality Control Kit.
5. Under certain circumstances, the anti-c specificity present in the plasma of vial 2 may react at room temperature with Affirmagen O cells when tested in conjunction with ORTHO BioVue® System Reverse Diluent cassettes, generating an unexpected weak positive reaction.

Performance Characteristics

- Each cell sample is shown to have a negative direct antiglobulin test.
- When properly stored and used according to standard procedures, these reagents will demonstrate the appropriate antigens / antibodies specified in the reagent description.
- The Procedure and Interpretation of Results must be followed closely to ensure the accuracy of the test results. Each laboratory should have a program that will train personnel on the proper use and handling of the product.
- Technical questions concerning this reagent should be directed to your local Ortho Care™ Technical Solutions Center.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

 Use by or Expiration Date (Year-Month-Day)	 Manufacturer	 Temperature Limitation
 LOT Lot Number	 IVD <i>In vitro</i> Diagnostic Medical Device	 Consult Instructions for Use
 REF Catalog Number or Product Code	 EC REP Authorised Representative in the European Community	

Pub. No. J39796_RU_RU

TRANSFUSION MEDICINE

INSTRUCTIONS

General Information

Guaranties

Ortho Clinical Diagnostics Inc. guarantees that the product is applicable for the intended use. The product is ready for use before the indicated shelf life and in the event of compliance with all storage, transport, package integrity and use procedures specified in the instructions for use.

Standards

The production is certified for compliance with the requirements of the international standard ISO 13485: quality management system. Medical devices comply with the National Standards of the Russian Federation and regulatory documents in terms of safety and efficiency requirements.

Potential users

Users - specialists in the field of laboratory diagnostics:

- a doctor of clinical laboratory diagnostics,
- medical laboratory technician (paramedic laboratory assistant),
- laboratory specialist (including with biological education)

Scope of application

Laboratory of IN VITRO diagnostics.

The medical device is intended for professional use. Training on work with the product is provided by Ortho. The medical product is used manually with the ORTHO BioVue system, as well as with automated systems for determining the blood group ORTHO AutoVue (Innova or Ultra), and with automated blood group testing systems - ORTHO VISION and ORTHO VISION Max analyzers.

Requirements for technical maintenance.

Maintenance and repair is not applicable.

Methods of sterilization.

It is not required sterilization

Software of medical device

It is absent.

Utilization

Care must be taken when handling materials obtained from human blood. All hematological materials should be considered as potentially infected. Disposal of waste must be in accordance with international requirements, federal and local legislation.

According to SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the management of medical waste" the product can be classified in class B.

Disposal or destruction, disinfection of reagent kits should be carried out in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for handling medical waste" and MU 287-113 "Methodological instructions for disinfection, pre-sterilization cleaning and sterilization of medical devices".

Measures in the event of a spillage of the reagent: wipe the contaminated surface with absorbent material (for example: fabric, fleece, etc.). Thoroughly clean the contaminated surface.

Pub. No. J39796_RU_RU

TRANSFUSION MEDICINE

INSTRUCTIONS

Summary of Revisions


Summary of Revisions

Version	Date of Revision	Description of Technical Changes
7	12-12-2019	Glossary of Symbols - Added EC Representative symbol Back Page - Added EC Representative address, Changed Notified Body number Registration certificate number added Changed Publication Number of IFU

Authorized representative of the manufacturer in the territory of the Russian Federation:

Limited Liability Company Ortho-Clinical Diagnostics
100 bld, Litera A, A., Novaya str, Skolkovo, Odintsovo district, Moscow region,
143026, Russia

Registration certificate number: № RZN 2019/9224 from 13.11.2019

Quotient Suisse S.A.
Unit B1 Terre Bonne Business Park
Route de Crassier 13
Eysins 1262, Switzerland

 **Alba Bioscience Limited**
James Hamilton Way, Penicuik, EH26 0BF, UK

Distributed by:

Ortho-Clinical Diagnostics, Inc.
1001 US Highway 202
Raritan, NJ 08869 USA

ORTHO CONFIDENCE™ WB Quality Control Kit is a
trademark of Ortho Clinical Diagnostics,
© Ortho Clinical Diagnostics, 2010-2019

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