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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Chengdu Besmile Biotechnology Co., Ltd on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhou Liping

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(Date: Feb. 10, 2020)

证书查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Instruction for use

Approved by : Moushan Liu Date : 2019.12.04
Signature : Moushan Liu Position : Manager Representative
Director : Xinzheng Yan Date : 2019.12.04

Chengdu Besmile Biotechnology Co., Ltd.

Instruction for use

Product name: Ceramic Zirconium Discs for CAD / CAM systems with accessories

Model specification:

Ceramic Zirconium Discs:

HT; ST; ST-preshade; SHT; SHTML; SHT-preshade; TT; TTML

98×10 mm, 98×12 mm, 98×14 mm, 98×16 mm, 98×18 mm, 98×20 mm, 98×22 mm, 98×25 mm, 95×10 mm, 95×12 mm, 95×14 mm, 95×16 mm, 95×18 mm, 95×20 mm, 95×22 mm, 95×25 mm, 89×71×12 mm, 89×71×14 mm, 89×71×16 mm, 89×71×18 mm, 89×71×20 mm, 89×71×22 mm, 89×71×25 mm)

Accessories (Coloring liquids) :

- Suitable for Ceramic Zirconium Discs HT, ST, SHT, TT:

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4; 1M1, 1M2, 2L1.5, 2L2.5, 2M1, 2M2, 2M3, 2R1.5, 2R2.5, 3L1.5, 3L2.5, 3M1, 3M2, 3M3, 3R1.5, 3R2.5, 4L1.5, 4L2.5, 4M1, 4M2, 4M3, 4R1.5, 4R2.5, 5M1, 5M2, 5M3: (50ml, 100ml)

- Suitable for all Ceramic Zirconium Discs:

O1, O2, O3, G1, G2, G3, P1, P2, P3, V1, V2, V3, TO1, TO2, TO3, TO4: (20ml)

Purpose: Ceramic Zirconium Discs are used for the manufacture of zirconium crowns, bridges and crowns for implants. Coloring liquids are used to color crowns and bridges before sintering.

Structure composition: The product contains Ceramic Zirconium Discs and Coloring liquids, in which Ceramic Zirconium Discs for all-ceramic dentures are mainly composed of zirconia (including hafnium oxide and yttrium oxide) ($ZrO_2 + HfO_2 + Y_2O_3$): $\geq 99\%$; yttrium oxide (Y_2O_3): 4.5% - 10.0%; hafnium oxide (HfO_2): $\leq 5\%$; alumina (Al_2O_3): $\leq 0.5\%$; other oxides: $\leq 0.5\%$. The Coloring liquids are mainly composed of purified water, ferric chloride, Erbium Chloride and polyethylene glycol.

Expected use: Ceramic Zirconium Discs is used for clinical patients with tooth defects and tooth loss. Coloring liquid is used for coloring Ceramic Zirconium Discs.

Contraindication: Contraindicated Use of Zirconia Material Allergy and Heavy Metal Allergy.

Product performance

Ceramic Zirconium Discs:

Physical and chemical properties	Density	$> 2.8 \text{ g/cm}^3$
	Sintered Density	$> 6.02 \text{ g/cm}^3$
	Flexural strength	$\geq 500 \text{ MPa}$
	Chemical solubility	$< 100 \text{ } \mu\text{g/cm}^2$
	Radioactivity	$\leq 1.0 \text{ Bq/g}$
	CTE	$(10.5 \pm 0.5) \times 10^{-6} \text{ K}^{-1}$
	Veneering porcelain bond strength	$> 25 \text{ MPa}$

Sintering temperature program:

Initial temperature (°C)	Heating time (min)	Temperature (°C)	Total (h)
Natural temperature	90	300	1.5
300	240	1480-1550	5.5
1480-1550	90	1480-1550	7.0
1480-1550	160	800	9.7
800	Natural cooling	Natural temperature	-

Coloring liquids:

Physical and chemical properties	Specification requirements	Specification deviation of coloring liquid should be (± 2 ml).
	Appearance	The coloring liquid should be homogeneous and impurity-free liquid.
	Post-dyeing comparison requirements	After dyeing, there was no significant difference in color between dental specimens and VITA veneer colorimeter or standard dental specimens.
	Heavy metal content	Heavy metal content (in pb) less than 10 mg/L.

Operation instruction

Ceramic Zirconium Discs usually used CAD/CAM (computer aided design / computer aided manufacturing) method of making porcelain dentures, production process is as follows:

- 1) Scanning the teeth or dental digital, access to dental three-dimensional data.
- 2) In accordance with the dental data CAD design, design of porcelain block processing model.
- 3) The use of CNC machine tools in accordance with the porcelain block processing model for CAM manufacturing, all-ceramic crowns made of denture.
- 4) Porcelain denture crown by sintering in air or vacuum sintering, to meet the use requirements of strength and appearance.
- 5) With porcelain powder in all-ceramic crowns of denture surface glazing, forming porcelain porcelain bonding, to create a full porcelain denture.

Coloring liquids:**Immersion method:**

- 1) Shake the coloring liquid well before dyeing. Choose the appropriate container and pour in the appropriate amount of coloring liquid.
- 2) Zirconia was lightly put into coloring liquid with plastic tweezers and soaked for 10 seconds.
- 3) Take out zirconia and dry it for 30-60 minutes (it is recommended to use infrared lamp or oven and place it in the position of 8-100mm under infrared lamp or at 80-90 C of oven).
- 4) The dyed zirconia was compared with VITA*16 color system, VITA*26 color system and Special colors (O1, O2, O3, G1, G2, G3, P1, P2, P3, V1, V2, V3, TO1, TO2, TO3, TO4).

Brushing method:

- 1) Shake the coloring liquid well before dyeing. Choose the appropriate container and pour in the appropriate amount of coloring liquid.
- 2) Zirconia was dipped in coloring liquid with a ceramic pen and brushed 1-2 times.
- 3) Zirconia is dried for 30-60 minutes (it is recommended to use infrared lamp or oven and place it in the position of 8-100mm under infrared lamp or at the temperature of 80-90 degrees C in the oven).
- 4) The dyed zirconia was compared with VITA*16 color system, VITA*26 color system and Special colors (O1, O2, O3, G1, G2, G3, P1, P2, P3, V1, V2, V3, TO1, TO2, TO3, TO4).

Storage requirements

The zirconia blanks must be stored in its original packaging in range of 10-50°C.

The coloring liquids must be stored in its original packaging, the best temperature for its storage is 10-30°C.

Transportation requirements

Transportation of zirconia blanks and coloring liquids should avoid damage from collision, extrusion and vibration.

The transport environment temperature should be maintained in range of 10-40 °C.

Temperature operation: 32-42 °C.

Sterility, method of sterilization

No sterility.

Not sterilable.

The medical device is disposable, not reused.

Warranty obligations

Ceramic Zirconium Discs shelf life: 5 years from date of manufacture.

Coloring liquids Shelf life: 1 year from date of manufacture.

Our products' warranty covers all product quality issues caused by material failure or production processes.

During the warranty period, all products with unqualified product quality can be replaced free of charge.

For specific replacement condition please contact with local dealers.

Chengdu Besmile Biotechnology Co., Ltd. is not responsible for any direct, indirect or eventual damage or delay caused by:

- Already used products;
- Damage caused by abnormal use beyond the specified conditions of use;
- Replace or remove original product label or manufacturing mark;
- Improper use.

Disposal considerations

Ceramic Zirconium Discs is a non-hazardous material. Disposal should be in accordance with applicable federal, state and local regulations. State or local regulations or restrictions may differ from federal regulations.

Remained coloring liquids should be diluted (10ml coloring liquids with 1L water) and disposal in accordance with applicable federal, state and local regulations. State or local regulations or restrictions may differ from federal regulations.

Packing instructions:

Each Zirconia blank with instruction for use is packed independently in a paper box, one blank in each box. Inside the box around the blank is PE foam, outside the box is bubble bag.

The coloring liquid is packed separately in plastic bottles.

Caution:

- 1) Storage, transportation and processing process to ensure good health environment and do not squeeze, bump.
- 2) The coloring liquid should be stored separately to avoid direct exposure to sunlight.
- 3) The dye is weak acidic and the skin is susceptible to allergies. If it splashes into the eyes, please rinse it with clean water immediately and see a doctor in time.
- 4) Don't dilute the coloring liquid by adding water to avoid affecting the color stability.
- 5) When soaking, make sure that zirconia is completely soaked in the Coloring Liquid.
- 6) Pneumoconiosis caused by inhalation of ceramic dust.
- 7) During the sintering process, the operator must wear appropriate protective apparatus.
- 8) Without sintering ceramic can not be applied directly to the body or oral denture repair; this product is required by professional artisans, by professional doctors with inlay, adjustment and wear.
- 9) Equipment or materials of the same linear expansion coefficient as our department when used with other materials or instruments.
- 10) Disposal of waste by the user in accordance with the requirements of the regulations.
- 11) This product needs to be operated by a professional dentist.

Index of CE symbol:

	Date of Manufacturer		Do not reuse	Lot	Batch number
	Warning, please refer to the instructions in the annex		Manufacturer		Keep dry
	European union authorization representative		0197	The product meets the basic requirements of European medical devices directive 93/42/EEC	

Chengdu Besmile Biotechnology Co. Ltd.

Address: No.1666, Chenglong Road Section Two,
Economic-technical Development Zone,
Chengdu, Sichuan 610100, China

Tel: +86-28-689078

TH Stomatology, LLC

Address: 125040, Moscow, st. Pravdy, 3/1,
premises 9, room 1-7

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