APPROVED BY

Sophie REHAULT
Quality and Regulatory Affairs Director

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INSTRUCTIONS FOR USE

IVA Introducer Sheaths, Sterile, Single-Use, Design Options

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These Instructions for Use have been prepared to provide information for registering a medical device in the Russian Federation. The information is provided as required in accordance with the Rules for State Registration of Medical Devices (Decree of the Government of the Russian Federation No. 1416 dated December 27, 2012). The contents of the Instructions for Use comply with the requirements of Order of the Ministry of Healthcare of the Russian Federation No. 11н "On Approval of Requirements to the Contents of Technical and Operational Documentation of the Manufacturer (Producer) of the Medical Device" dated January 19, 2017.

The device is developed, manufactured and sold in accordance with the requirements of EN ISO 13485:2016 + AC : 2017-07; EN ISO 13485:2016 + AC:2016; ISO:2016 (Certificate No. 513975 MP2016) and is CE marked (EC certificate number: 513975 MR2 annex II excluding section 4).

1. Product name:

IVA introducer sheaths, Sterile, Single-Use, Design Options:

I. IVA introducer sheaths, Short, Sterile, Single-Use, Design Options:

- IVA3F;

- Instructions for use (supplied with each device) – 1 pc.

II. IVA introducer sheaths, Long, Sterile, Single-Use, Design Options:

- IVA3F50;

- IVA5F45;

- IVA5F41.45;

- IVA5F96.45;

- IVA6F80;

- IVA6F80ST;

- IVA6F80ST_MP;

- IVA7V41.45;

- IVA7F96.45;

- IVA7F96.90;

- Instructions for use (supplied with each device) – 1 pc.

Device classification in accordance with the Nomenclature Classifier of Medical Devices approved by Order of the Ministry of Healthcare of the Russian Federation No. 4н dated June 6, 2012 – 192010.

OKPD2 – 32.50.13.190 – Instruments and devices used for medical purposes, not included in other groups.

Class depending on the potential risk degree of medical application in accordance with the Nomenclature Classifier of Medical Devices approved by Order of the Ministry of Healthcare of the Russian Federation No. 4н dated June 6, 2012 – 3.

Type of contact with the human body – short-term (less than 24 hours) contact with circulating blood.

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3. Intended use of the device and principles of operation

3.1. Intended use of the device: The device is intended for delivery of therapeutic or diagnostic devices to the vascular system. They enable the insertion and manipulation of such devices as guidewires, guide catheters, hyper-selective neuroradiological catheters, angioplasty catheters, or interventional radiology catheters.

3.2. Principle of operation: IVA introducer sheaths, Sterile, Single-Use, Design Options, are inserted into an artery upon completing a Seldinger procedure to provide safe blood vessel access.

The sheath is used to insert catheters or other devices for endoluminal procedures, such as angioplasty. Fluoroscopy can be applied to monitor the catheter position and advance it. Contrast media can be injected to visualize the organs. The device can be used for such interventional procedures as thermal ablation, angioplasty, embolization, or biopsy.

When the procedure is completed, the sheath is removed. In some cases, sealing devices can be used for post procedure puncture closure.

a) Insert the needle percutaneously into the femoral artery. The needle is carefully inserted until blood appears.


Figure 1. Needle insertion.

b) Insert the angiography guidewire and then advance it up to the artery. After positioning the guidewire, remove the needle.


Figure 2. Guidewire insertion.

c) Then insert the sheath with the dilator.

Figure 3. Sheath insertion.

d) Remove the guidewire and the dilator. The sheath remains in the blood vessel to ensure patency and prevent injuries.


Figure 4. Dilator removal.

e) Then you can proceed with insertion of intravascular devices.

4. Description of basic functional characteristics:

The device is intended for delivery of therapeutic or diagnostic devices to the vascular system. They enable insertion and manipulation of such devices as guidewires, guide catheters, hyper-selective neuroradiological catheters, angioplasty catheters, or interventional radiology catheters.

As for the raw materials and the manufacturing process, the sheath and the dilator of IVA introducer sheaths, Sterile, Single-Use, Design Option, Short, and IVA introducer sheaths, Sterile, Single-Use, Design Option, Long, are the same. However, there are differences between the two device types, as provided below:

- Dilator:

- Different dimensions and length.

- Sheath:

- Different dimensions and length
- Availability of a 6F (800 mm) soft tip (+/- angular) on the sheath end in the design option: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

- Sheath and dilator protector:

- **Sheath and dilator protector** – a PVC protective tube covering the dilator and the sheath to prevent the device from kinking during its transportation and storage up to its use by a physician.
- Not available in IVA Short IVA3F and IVA Long IVA7F96.90 design options.
- The protector does not contact the patient's blood, mucosal membranes or skin.


- Haemostatic valve:

- Non-removable (IVA Short) vs. removable (IVA Long).

- Color code:

- The anti-kinking sleeve is not exposed to blood (IVA Short) unlike the dilator with indirect blood contact (IVA Long).

- Formed (IVA Short) vs. pressed (IVA Long) size symbol.
- **Guidewire:**
 - Not available in the design option: IVA Long.
- **Needle** - designed for puncture during the introduction of the Introducer. Insert the needle percutaneously into the femoral artery. The needle is carefully inserted until blood appears. After installing the Guidewire, the needle is removed.
- **The needle protector** is a protective cap made of polyethylene put on the needle to prevent it from being damaged during transportation and storage up to its use by a physician.
- Not available in the design option: IVA Long.

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Guidewire protector
Straight guidewire
Needle

Figure 5. MD appearance: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short, unpacked.

Figure 6. MD appearance: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long, unpacked.

5. Description of main assemblies and parts:

5.1. IVA Introducer Sheaths, Sterile, Single-Use, Design Option, Short, consist of the following main parts: 1. Sheath; 2. Dilator; 3. Sheath hub; 4. Dilator hub; 5. Haemostatic valve; 6. Anti-kinking sleeve; 7. Transverse branch; 8. Straight guidewire; 9. Guidewire protector; 10. Needle.

Note: The sheath hub consists of a) a body; b) a ring; c) a connector.

The transverse branch consists of a) a side arm; b) a stopcock consisting of a body and a 3-way key.

This device consists of a dilator and a sheath. On the proximal end of the sheath, there is a non-removable silicone haemostatic valve connected to an extension infusion line with a transverse branch which prevents blood outflow and air aspiration.

The dilator sheath is used to dilate the vessel and ensure insertion of the dilator followed by the applied medical device. The dilator sheath is made of high density polyethylene (HDPE), bismuth carbonate and blue pigment. The dilator has a long, rounded tip to ensure non-traumatic insertion. The dilator is secured in the valve to prevent leakages during insertion.

5.2. IVA Introducer Sheaths, Sterile, Single-Use, Design Option, Long, consist of the following main parts: 1. Dilator sheath; 2. +/- soft tip; 3. Dilator; 4. Sheath; 5. Haemostatic valve; 6. Transverse branch.

Note: As with IVA Introducer sheaths, Sterile, Single-Use, Design Options, Short, on the proximal end of the dilator sheath, there is a non-removable silicone haemostatic valve encased in a transparent hub and connected to an extension infusion line with a stopcock (side arm) which prevents blood outflow and air aspiration. Unlike short models, long IVA introducer sheaths have removable haemostatic valves. The guidewire is not supplied with long IVA devices.

IVA Introducer sheaths, Sterile, Single-Use, Design Option, Long, are available in different diameters (French scale) and lengths. For convenience, the proximal ends of the sheath and the dilator have color marking and coding specifying the size. IVA introducer sheaths, Sterile, Single-Use, Design Options, Long, are available with diameters from 3F to 8F and lengths from 500 mm to 960 mm. IVA introducer sheaths, Sterile, Single-Use, Design Options, Long, are not available in 4F size. Two 6F articles have a soft tip (ST) on the sheath end. The tip can be either angled (MP) or straight.

6. Indications and contraindications, side effects, precautions:

a. Indications

IVA introducer sheaths, Sterile, Single-Use, Design Options, are designed to provide vascular access for insertion of catheters for treatment or diagnostic purposes.

b. Contraindications

The use of IVA introducer sheaths, Sterile, Single-Use, Design Options, is contraindicated for neonates, premature infants, infants, children, pregnant and breast-feeding women.

c. Potential complications and side effects

Potential complications include, but are not limited to:

- Haematoma at the insertion site;
- Embolism;
- Infection;
- Vascular wall perforation;
- Intima dissection;
- Thrombus formation.

d. Precautions

- Do not use if package is opened or damaged. The devices are sterile if package is undamaged.
- Percutaneous introducer sheaths are intended for single use only. Do not reuse. Any reuse of the medical device may result in a high risk of contamination to the patient and performance degradation.
- Do not resterilize these devices.
- Store in a dry and dark place at room temperature. Do not use past the expiration date.
- These devices should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology and/or cardiology.
- Never move an intravascular device if resistance is met until the causes of resistance have been determined by means of angiographic control. Movement against resistance may result in damage to the catheter or injury to the patient.
- Follow the instructions for use of the catheters inserted through the introducer for diagnostic or therapeutic purposes.
- Follow the instructions for use of the drugs injected into the percutaneous introducer sheath.
- Never use an injector on the side injection device for contrast injection.

e. Information on potential consumers

The device should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology. *Conditions of use:* Medical treatment facilities and healthcare centers, only when prescribed by a medical specialist.

7. Technical specifications of the medical device:

Table 1. MD technical specifications: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short.

Part No.	Dimensions	Color marking of the anti-kinking sleeve	Effective length of the sheath (mm)	Inner diameter of the sheath (mm)	Guidewire length (mm)	Guidewire diameter (mm)	NeedleØ, (mm)	Weight (g)	Minimum tear strength between the connector hub and the sheath, N	Minimum tear strength between the connector hub and the dilator, N
IVA3F	3F	White	70	1.13 ± 0.02	300	0.46	21G	9.3	10	5

Unspecified acceptable values: $\pm 1\%$

Note: IVA percutaneous introducer sheath, short, is supplied with a stainless steel puncture needle, a straight stainless steel guidewire, a dilator and a sheath made of high density polyethylene.

The short introducer family includes only one design option: with 3F diameter and effective length of 70 mm. The anti-kinking sleeve located between the sheath and haemostatic valve is color coded and formed with the model size to facilitate its identification.

Table 2. MD Technical Specifications: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

Part No.	Size (French scale)	Color code of the dilator hub and the sheath	Effective length of the sheath (mm)	Sheath ID/OD (mm)	Dilator ID/OD (mm)	Dilator effective length (mm)	Weight (g)	Min. tear force of the sheath, N	Min. tear force of the dilator, N
IVA3F50	3F	White	500	1.13/1.45	0.56/1.05	575	6.9	10	5
IVA5F41.45	5F	Grey	410	1.85/2.40	1.1/1.75	485	6.9	15	10
IVA5F45			450			525	6.9		
IVA5F96.45			960			1035	6.9		
IVA6F80			800			525	14.6		
IVA6F80ST	6F	Green	800	2.20/2.80	1.1/2.10	725	14.6	15	15
IVA6F80ST MP			800			875	14.6		
IVA7F41.45			410			485	14.6		
IVA7F96.45	7F	Orange	960	2.55/3.15	1.1/2.45	1035	14.6	15	15
IVA7F96.90			960			1035	6.9		

Unspecified acceptable values: $\pm 1\%$

Note: The dilator and sheath hubs have specific color marking and size stamps for identification purposes.

7.1. Device appearance: The device surface has no foreign particles, bubbles or stains. The surface is smooth and free from technological and external defects. The metal parts of the device should not show signs of corrosion affecting the functioning of the products and biocompatibility.

These devices are used for arterial and venous procedures that require percutaneous insertion of intravascular devices. They enable insertion and manipulation of such devices as guidewires, guide catheters, hyper-selective catheters for neuroradiology, angioplasty, or interventional radiology. The hub is tightly fitted around the tube ensuring mechanical integrity.

Requirements for Radiopaque: due to radiopaque materials, bismuth trioxide and bismuth carbonate included in the components of the product "Shell and expander, soft tip" of the product are radiopaque, which facilitates the introduction of the product and allows the use of fluoroscopic control.

Table 3. Requirements for Needle.

Parameter	Index
OD, mm	from 0.800 to 0.830
ID, mm	0.54 -0.02 / + 0.01
Effective length mm	42 $\pm$ 1.5
Sharpening angle	18° $\pm$ 2
Roughness parameter	no more than Ra 2.5 μ m
Rockwell hardness:	52-54 HRC

Parameter	Index
Radius of blunting of the working part, no more	0.03 mm
Connection type	Taper not less than 6% according to ISO 594-1
Connection strength, not less	20 N

Note: The device is intended for delivery of therapeutic or diagnostic devices to the vascular system. They enable insertion and manipulation of such devices as guidewires, guide catheters, hyper-selective neuroradiological catheters, angioplasty catheters, or interventional radiology catheters.

As for the raw materials and the manufacturing process, the sheath and the dilator of IVA introducer sheaths, Sterile, Single-Use, Design Option, Short, and IVA introducer sheaths, Sterile, Single-Use, Design Option, Long, are the same. However, there are differences between the two device types, as provided below:

- Dilator:

- Different dimensions and length.

- Sheath and dilator protector:

- **Sheath and dilator protector** – a PVC protective tube covering the dilator and the dilator sheath to prevent the device from kinking during its transportation and storage up to its use by a physician.

- Not available in IVA Short IVA3F and IVA Long IVA7F96.90 design options.

The protector does not contact the patient's blood, mucosal membranes or skin.

- Sheath:

- Different dimensions and length.

- Availability of a 6F (800 mm) soft tip (+/- angular) on the sheath end in the design option: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

- Haemostatic valve:

- Non-removable (IVA Short) vs. removable (IVA Long).

- Color code:

- The anti-kinking sleeve is not exposed to blood (IVA Short) unlike the dilator with indirect blood contact (IVA Long).

- Formed (IVA Short) vs. pressed (IVA Long) size symbol.

- Guidewire:

- Not available in the design option: IVA Long.

- **Needle** - designed for puncture during the introduction of the Introducer. Insert the needle percutaneously into the femoral artery. The needle is carefully inserted until blood appears. After installing the Guidewire, the needle is removed.

- **The needle protector** is a protective cap made of polyethylene put on the needle to prevent it from being damaged during transportation and storage up to its use by a physician.

- Not available in the design option: IVA Long.

7.2. MD schematic view: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short

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Side arm

Straight guidewire

Figure 7. Device structure: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short:

1. Sheath; 2. Dilator; 3. Sheath hub; 4. Dilator hub; 5. Haemostatic valve; 6. Anti-kinking sleeve; 7. Transverse branch; 8. Straight guidewire; 9. Guidewire protector (see Fig. 5); 10. Needle (see Fig. 5).



Figure 8. Schematic drawing of the 3-way key and the device body.

7.3. MD schematic view: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long

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Transverse branch

Side arm

Figure 9. Device structure: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long:

1. Sheath and dilator; 2. +/- soft tip; 3. Dilator hub; 4. Sheath hub; 5. Haemostatic valve; 6. Transverse branch.

7.4. MD drawing: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short, with technical specifications


Figure 10. MD drawing: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short, with technical specifications.

Table 4. MD Technical Specifications: IVA introducer sheaths, Sterile, Single-Use, Design Option, Short.

Dilator length, mm	Sheath length, mm	Straight guidewire length, mm	Sheath retractable dilator length, mm
A	B	C	D
129±3	70±3	300±3	15±0.5

IVA Introducer sheaths, Sterile, Single-Use, Design Option, Short	Dilator		Sheath	
	ID, mm	OD, mm	ID, mm	OD, mm
IVA3F	0.56	1.05	1.13	1.40

7.5. MD drawing: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long, with technical specifications


Table 6. MD Technical Specifications: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

Page 12/21

IVA Introducer sheaths, Sterile, Single-Use, Design Option, Long	Dilator length, mm	Sheath length, mm	Sheath retractable dilator length, mm	+/- soft tip, mm (Only applicable for models: IVA6F80ST, IVA6F80ST_MP)
IVA6F80ST	875+5/-0	800+5/-0		50±0.5
IVA6F80ST_MP	875+5/-0	800+5/-0		

Table 7. MD Technical Specifications: IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

IVA Introducer sheaths, Sterile, Single-Use, Design Option, Long	Dilator		Sheath	
Design option	OD, mm	ID, mm	ID, mm	OD, mm
IVA3F50	1.05	0.56	1.13	1.45
IVA5F45	1.75	1.10	1.85	2.40
IVA6F80	2.10	1.10	2.20	2.80
IVA5F41.45	1.75	1.10	1.85	2.40
IVA5F96.45				
IVA7F41.45	2.45		2.55	3.15
IVA7F96.45				
IVA7F96.90				
IVA6F80ST	2.10		2.20	2.80
IVA6F80ST MP				

Unspecified acceptable values: ±1%

8. List of the medical device materials:

Table 8. The materials used in IVA introducer sheaths, Sterile, Single-Use, Design Option, Short.

Name of the device or component and accessory (under study)			Materials used in manufacture	Manufacturer
Dilator and Sheath			High density polyethylene (low-pressure) (Brand: Alcudia 4810 S25K)	GAZECHIM
			Blue pigment (Brand: PIGBLEU30P410)	PLASTI COULEUR
			Bismuth carbonate (CAS # 5892-10-4)	5N PLUS FRANCE
Dilator hub			Polypropylene (Brand: PPH5060) Blue pigment (Brand: PIGBLEU30P410)	SPEEDPLAST PLASTICOCULEUR
Sheath hub	Body	Natural polycarbonate (Brand: LEXAN 124R112)		GE PLASTICS FRANCE
	Ring	ABS (acrylonitrile butadiene styrene) (Brand: TERLUX 2802TR)		SPEEDPLAST
	Connector	Polyamide (Brand: Grilamid TR90)		ESSE INTERNATIONAL
Haemostatic valve			Natural silicone (CAS# 10097-28-6)	BALTON
Anti-kinking sleeve			Thermoplastic (Brand: F4STE (TPE))	KRAIBURG TPE
			White pigment (Brand: PIGBLANC10020)	PLASTICOCULEUR
Transverse branch	Side arm		PVC (Brand: gfm/4)	BALTON
	Stopcock	Body	Natural polycarbonate (Brand: LEXAN 124R112)	
		3-way key	Polyethylene (CAS# 9002-88-4)	
			White pigment (Brand: POLYWHITE 8000)	PLASTICOCULEUR

Name of the device or component and accessory (under study)	Materials used in manufacture	Manufacturer
Straight guidewire	Stainless steel (Brand: 304V)	GUIDEWIRE TECHNOLOGIES
Guidewire protector	PVC (Brand: VH-611/3)	INEOS Compounds Aycliffe Ltd, United Kingdom
Needle	Stainless steel (Brand: SUS 304)	UNISIS
Needle protector	Polyethylene (CAS# 9002-88-4)	UNISIS, Japan
Blister*	PVC (Brand: VH-621/9)	MKF-ERGIS Sp. z o.o., Poland
Individual package	Tyvek® 1059B and ABS plastic (CYCOLAC EX58F)	SUDPACK MEDICA, France

Note: Pigment ratio is typically more than 1 g per 5 kg of polymer.

*The device is packed in a blister and then in an individual package.

Table 9. The materials used in IVA introducer sheaths, Sterile, Single-Use, Design Option, Long.

Name of the device or component and accessory (under study)			Materials used in manufacture	Manufacturer
Dilator and Sheath			High density polyethylene (low-pressure) (Brand: Alcudia 4810 S25K)	GAZECHIM
			Blue pigment (Brand: PIGBLEU30P410)	PLASTI COULEUR
			Bismuth carbonate (CAS # 5892-10-4)	5N PLUS FRANCE
+/- soft tip			Polypropylene (Brand: PPH5060)	Total Petrochemicals (Manufacturer) Speedplast (Distributor)
			Thermoplastic (Brand: TF4STE)	TPE Kraiburg GmbH
			Bismuth trioxide (CAS # 1304-76-3)	5N PLUS FRANCE
			Red pigment (Brand: PIGROUGE60P260)	PLASTICOULEUR
Dilator hub			Polypropylene (Brand: PPH5060)	TOTAL Research and Technology Feluy, Belgium
			White pigment (Brand: PIGBLANC10020) Black pigment (Brand: PIGNOIR) Dark-Green pigment (Brand: PIGVERTFONCE70P) Pale-Green pigment (Brand: PIGVERTCLAIR70P) Yellow pigment (Brand: PIGJAUNE) Red pigment (Brand: PIGROUGE60P260)	PLASTICOULEUR
Sheath and dilator protector			Polypropylene (Brand: PPH5060)	TOTAL Research and Technology Feluy, Belgium
Sheath hub	Body		Polyamide (Brand: Grilamid TR90)	ESSE INTERNATIONAL
	Ring		Natural polycarbonate (Brand: LEXAN 124R112)	GE PLASTICS FRANCE
	Connector		Natural polycarbonate (Brand: LEXAN 124R112)	GE PLASTICS FRANCE
Haemostatic valve			Natural silicone (CAS# 10097-28-6)	BALTON
			PVC (Brand: gfm/4)	BALTON
Transverse branch	Side arm		Natural polycarbonate (Brand: LEXAN 124R112)	GE PLASTICS FRANCE
	Stopcock	Body		
		3-way key	Polyethylene (CAS# 9002-88-4)	LM Plastiques, France
			White pigment (Brand: POLYWHITE 8000)	PLASTICOULEUR
+/- soft tip casing*			Polyethylene (CAS# 9002-88-4)	LM Plastiques, France
Individual package			Tyvek® 1059B and ABS plastic (CYCOLAC EX58F)	SUDPACK MEDICA, France

Note: Pigment ratio is typically over 1 g per 5 kg of polymer. The haemostatic valve and the sheath proximal end with color marking does not come into contact with circulating blood, but can have indirect contact with blood.

*The +/- soft tip casing is used in the packaging only for IVA introducer sheaths, Sterile, Single-Use, Design Option, Long of IVA7F96.90 to hold in place the 90° curve.

9. Description of accessories and medical devices used with the reported device:

IVA Introducer Sheaths, Sterile, Single-Use, Design Options, are invasive devices with a 6 % taper Luer connector used with any devices for insertion of endovascular surgery devices. These devices are used for arterial and venous procedures that require percutaneous insertion of intravascular devices. They enable insertion and manipulation of such devices as guidewires, guide catheters, hyper-selective catheters for neuroradiology, angioplasty, or interventional radiology. The hub is tightly fitted around the tube ensuring mechanical integrity.

10. Information on the procedure of installation and other actions required to commission the medical device:

10.1. IVA introducer sheath preparation:

Check the package integrity: it should not be opened or damaged. Open the pouch and remove various components using an aseptic technique. Inspect them to ensure they are not damaged. Use sterile saline solution to check the system patency and remove air from it. Insert the dilator into the sheath.

10.2. IVA introducer sheath insertion:

After performing local anesthesia (if applicable):

- Puncture the selected vessel wall with an appropriately sized needle.
- Insert the flexible part of the guidewire into the needle lumen. When using a J-wire, use a straightener.
- Remove the needle when the guidewire is inserted into the vessel lumen. If necessary, apply pressure to the puncture site.
- Connect the flushing line with the 3-way key of the stopcock on the introducer sheath. Fill the introducer sheath with heparinized saline and remove air.
- Prepare the dilator using a syringe filled with heparinized saline.
- Insert the dilator into the introducer sheath completely.

CAUTION! Insert the dilator in the center of the sheath valve. Forced insertion outside the valve center may cause damage to the valve and blood leakage.

- Place the dilator with the sheath on the guidewire. When advancing the sheath with the dilator percutaneously, apply slight rotation.
- Verify the introducer position angiographically.
- Remove the guidewire and the dilator while holding the introducer sheath in place.
- When introducing, manipulating or removing the catheter from the sheath, always hold the sheath in place.

10.3. IVA introducer sheath removal:

After removing the catheter for treatment or diagnostic purposes:

- Slowly retract the introducer sheath moving in parallel to the skin surface.
- After removing the sheath, apply pressure to the puncture site.

11. Information on the conditions of transportation, storage and operation of the medical device:

Transportation conditions: Ambient temperature from -10 °C to +35 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight.

Storage conditions: Store the device in its original package in a dry and clean place, away from heat sources and direct sunlight, at controlled temperature. The shelf life is indicated on each device label. Ambient temperature from +15 °C to +27 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa.

Operation conditions: Ambient temperature from +20 °C to +25 °C, relative humidity from 35 % to 65 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight. The device is temporarily operated inside the human body at a temperature ranging from +32 °C to +42 °C. The device is operated inside the human body for less than 24 hours.

12. Information on the use of a human medicinal product, animal or human origin materials in the medical device:

Not applicable. IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not contain any medicinal products, pharmaceutical substances, animal or human origin materials.

13. Information on the procedure for installation, assembly, setting and calibration required to commission the medical device:

Do not use if package is opened or damaged. The devices are sterile if package is undamaged.

- IVA Introducer Sheaths, Sterile, Single-Use, Design Options, are intended for single use only. Do not reuse. Any reuse of the medical device may result in a high risk of contamination to the patient and performance degradation.
- Do not resterilize these devices.
- Store in a dry and dark place at room temperature.
- Do not use past the expiration date.
- These devices should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology and/or cardiology.
- Follow the instructions for use provided with compatible devices to be used together with this device in the field of interventional neuroradiology and/or interventional radiology and/or cardiology.
- Follow the instructions for use of devices used in combination with the medical device: IVA Introducer Sheaths, Sterile, Single-Use, Design Options.

14. Requirements to the premises intended for installation of the medical device:

IVA Introducer Sheaths, Sterile, Single-Use, Design Options, shall be used in the premises with ambient temperature from +20 °C to +25 °C, relative humidity from 35 % to 65 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight. The device is temporarily operated inside the human body at a temperature ranging from +32 °C to +42 °C. The device is operated inside the human body for less than 24 hours.

15. Information for checking correctness of installation of the medical device and possibility of safe operation/use, including:

a) Scope and frequency of maintenance, including cleaning and disinfection of the medical device:

IVA Introducer Sheaths, Sterile, Single-Use, Design Options, shall not be maintained, cleaned or disinfected.

b) The list of information, keys and access passwords required for assembly and maintenance of the medical device provided by the manufacturer:

IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not contain any software, information, keys or access passwords required for assembly and are not subject to maintenance.

c) List of consumables, as well as procedure for use and replacement thereof:

IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not contain any consumables.

d) Need to calibrate the medical device to ensure appropriate and safe operation throughout the service life:

Not applicable. MD "IVA Introducer Sheaths, Sterile, Single-Use, Design Options" does not require calibration to ensure appropriate and safe operation throughout the service life.

e) Methods of mitigating risks related to the installation, calibration or maintenance of the medical device:

To mitigate the risks related to installation of MD "IVA Introducer Sheaths, Sterile, Single-Use, Design Options", read carefully the Instructions for Use provided with this device and concurrently used devices. For additional information please refer to Clause 9 "Information on the Procedure for Installation, Assembly, Setting, and Calibration Required to Commission the Medical Device".

The device is intended for single use and does not require maintenance. Calibration is also not required, all the necessary settings are preset by the manufacturer.

f) Information on assembly, adjustment, setting, calibration and other actions required to commission the medical device and for correct use (operation) thereof:

Not applicable. MD "IVA Introducer Sheaths, Sterile, Single-Use, Design Options" does not require assembly, adjustment, setting, calibration and other actions to commission the medical device and for correct use (operation) thereof.

g) Information on the list of main specifications for the use of the medical device, transportation and storage conditions:

Transportation conditions: Ambient temperature from -10 °C to +35 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight.

Storage conditions: Store the device in its original package in a dry and clean place, away from heat sources and direct sunlight, at controlled temperature. The shelf life is indicated on each device label. Ambient temperature from +15 °C to +27 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa.

16. Information on the sterile condition of the medical device, sterilization method and procedure in case of breach of sterile packaging:

Each MD "IVA Introducer Sheaths, Sterile, Single-Use, Design Options" is supplied as a single-use, sterile device in a terminally sterilized package and sterilized with ethylene oxide gas. The device should not be used and should be disposed of immediately if the package sterility is compromised. Resterilization of the device may result in a high risk of contamination to the patient and performance degradation.

17. Information on processing of the medical device for reuse:

The device is supplied sterile and should not be subject to processing, disinfection or resterilization. After use, the device should not be reused and should be disposed of immediately.

Reuse of the device may result in a high risk of contamination to the patient and performance degradation.

18. Information required to identify medical devices to obtain a safe combination and information on the known restrictions of concurrent use of medical devices:

Cl. 9 "Description of Accessories and Medical Devices Used with the Reported Device" specifies devices IVA Introducer Sheaths, Sterile, Single-Use, Design Options, can be used with. Please refer to the instructions for use for each BALT medical device for all the necessary information on the use and known limitations.

19. Information on the nature, type, intensity and distribution of radiation during the use of the medical device:

Not applicable. IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not require RF energy and therefore do not create electromagnetic fields and interference induced thereby.

20. Information on the precautions taken if

a) The medical device is defective, malfunctions or features deviations of operation when these issues may affect safety of the medical device, including visual defects:

Please refer to Clause 10 for the information on identifying potential failures, malfunctions or deviations in operation that may affect safety of the medical device, including those identified by means of visual inspection for the possibility of their elimination.

b) The medical device is exposed to external factors related to the use of the medical device with other medical devices:

Please refer to the instructions for use of the devices used with IVA Introducer Sheaths, Sterile, Single-Use, Design Options (Clause 9. Instructions for Use);

c) A risk of electromagnetic interference created by the medical device affecting other medical equipment:

Not applicable. IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not require RF energy and therefore do not create electromagnetic fields and interference induced thereby.

21. Warnings and precautions taken by the consumers when using the medical device containing a human medicinal product, carcinogenic, mutagenic or toxic materials, possible extraction or washout whereof results in sensitization, allergic reaction or negatively affects fertility:

Not applicable. IVA Introducer Sheaths, Sterile, Single-Use, Design Options, do not contain any medicinal products, pharmaceutical substances, animal or human origin materials, carcinogenic, mutagenic or toxic materials possible extraction or leaching whereof results in sensitization, allergic reaction or negatively affects fertility.

22. Warnings and precautions taken by the consumer when disposing of the medical device, accessories and consumables used therewith:

Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). Package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

23. Information on the circumstances requiring consultation of the consumer with a health care professional:

The device should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology. Conditions of use: Medical treatment facilities and healthcare centers, only when prescribed by a medical specialist. In view of the foregoing, it follows that in any case, the consumer should consult with a health care professional to determine the necessity to perform the procedure, and as a consequence, the use of MD "IVA Introducer Sheaths, Sterile, Single-Use, Design Options" prior to the cerebral embolization.

24. Information on the first issue or the latest revision of the operating instructions:

Not applicable. The information on the first issue and the latest revision of the operating instructions is missing.

25. Procedure and conditions of disposal or destruction of the medical device:

Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). Package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

26. Information printed on the product label:

The marking on the packaging of the final sterilization is identical to the marking on the individual cardboard packaging.

The label contains the following information:

- Device name and model;
- Company logo;
- Information on the manufacturer;
- Schematic drawing of the device;
- Quantity;
- Dilator ID/OD;
- Sheath OD;
- Dilator effective length;
- Lot number;
- Catalogue number;
- Date of manufacture;
- Expiration date;
- "Do not use if package is damaged" symbol;
- "Do not reuse" symbol;
- "Keep away from sunlight" symbol;
- "Consult instructions for use" symbol;
- "Caution! Consult accompanying documents" symbol;
- "Do not resterilize" symbol;
- "EO sterilized" symbol;
- "Keep dry" symbol;
- "CE mark" symbol;
- "QR code" symbol;
- "Non-pyrogenic" symbol;
- "European article number" symbol;
- "Caution! Fragile!" symbol;
- "Non-toxic" sign;

- Device name in different languages;
- Number of units in the package.
- Information about the authorized representative*;
- Registration Certificate No.*

*After the registration procedure.

Table 10. List of conventional symbols

	Catalogue number
	Lot number
	Use by
	Date of manufacture
	CE marking
	EO sterilized
	Caution! Consult instructions for use
	Do not reuse
	Do not use if package is damaged
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	European article number
	Manufacturer
	Do not re-sterilize
	QR code
	Caution! Fragile
	Non-pyrogenic

27. Guarantees

Shelf life – 5 years. Sterilization validity – 5 years.

Address for complaints:

Center for Perinatal Medicine, Limited Liability Company
Legal / Actual address: 20/4 Tovarishchesky per., 109004, Moscow
INN 7709063317, KPP 770901001, OKPO 40286236.

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