STATE OF NORTH CAROLINA

Department of The
Secretary of State


APOSTILLE

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1. COUNTRY: UNITED STATES OF AMERICA
2. This Public Document
has been signed by TIFFANNY A THOMAS
3. acting in the capacity of NOTARY PUBLIC
COMMISSION EXPIRES 9/5/2022
4. bears the seal/stamp of COUNTY OF GUILFORD, NORTH CAROLINA

CERTIFIED

5. at Raleigh, North Carolina
6. the 12TH DAY OF JUNE, 2019
7. by Secretary of State or Deputy Secretary of State, State of North Carolina
8. No. 9981
9. Seal/Stamp
10. Signature

Scan to verify online.

Claine J. Marshall
Claine J. Marshall

Certification# 8101615926

Verify the issuance of this Apostille at <http://www.sosnc.gov/aut>

Secretary of State

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and where applicable, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille is not valid for use anywhere within the United States of America, its territories or possessions. (Verify the issuance of this Apostille at <http://www.sosnc.gov/apostille>)

Endoscopic Wire Guides with Accessories

INSTRUCTION FOR USE

Endoscopic Wire Guides (ACRO-25-260; -25-450; -35-205; -35-260; -35-450)

Endoscopic Wire Guides (AWG2-25-260; AWG2-25-260-A; AWG2-25-450; AWG2-25-450-A; AWG2-35-205; AWG2-35-205-A; AWG2-35-260; AWG2-35-260-A; AWG2-35-450; AWG2-35-450-A)

Endoscopic Wire Guides (FS-LT-35-205; -35-260; -35-480)

Endoscopic Wire Guides (FSW-35)

Endoscopic Wire Guides (HYB-48015; -48025)

Endoscopic Wire Guides with adapter Pin-Vise (HYB-48015-A; -48025-A)

Endoscopic Wire Guides with adapter Olcott Torque Device (MET-25-480; -25-480-A; -35-480; -35-480-A)

Endoscopic Wire Guides with adapter Olcott Torque Device (METII -21-260; -21-480; -25-260; -25-480; -35-260; -35-480; -35-480-A; -35-600E)

INTENDED USE

This device is used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: perforation, allergic reaction to contrast or medication.

PRECAUTIONS

Use of this wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide.

This wire guide may be left in place when used under normal conditions in conjunction with, and according to the instruction for, a compatible Cook electrosurgical (ES) device.

TECHNICAL SPECIFICATIONS

Page 1 of 11

ACRO

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
ACRO-25-260	0.064/ 0.025	260	Straight	27
ACRO-25-450	0.064/ 0.025	450	Straight	27
ACRO-35-205	0.089/ 0.035	205	Straight	27
ACRO-35-260	0.089/ 0.035	260	Straight	27
ACRO-35-450	0.089/ 0.035	450	Straight	27

AWG2

Roughness, Ra - not more than 0.64 microns

All AWG2 wire guides evenly and consistently marked with 25 marks with an interval of 1 cm.

The distance from the tip to the first mark is 5 cm.

At a distance of 5 cm from the first mark there are 2 marks.

At a distance of 10 cm from the first mark there are 3 marks.

At a distance of 15 cm from the first mark there are 4 marks.

At a distance of 20 cm from the first mark there are five marks

Product	Diameter inch/ mm	Length cm	Tip Type	radius of curvature inch/ cm	Hydrophilic Tip Coating Length cm
tolerances	$\pm 0.04\text{mm}$ / ± 0.002 inch	+4cm / -5cm	-	-	$\pm 5\text{mm}$
AWG2-25-260	0.025/0.064	260	Straight	-	5.8
AWG2-25-450	0.025/0.064	450	Straight	-	5.8
AWG2-25-260-A	0.025/0.064	260	Angled	Between 0.432/1.1 and 1.000/2.54cm	5.8
AWG2-25-450-A	0.025/0.064	450	Angled	Between 0.432/1.1 and 1.000/ 2.54cm	5.8

Product	Diameter inch/ mm	Length cm	Tip Type	radius of curvature inch/ cm	Hydrophilic Tip Coating Length cm
tolerances	±0.10mm / ±0.004 inch	+4cm /-5cm	-	-	± 5mm
AWG2-35-205	0.035/0.089	205	Straight	-	5.8
AWG2-35-260	0.035/0.089	260	Straight	-	5.8
AWG2-35-450	0.035/0.089	450	Straight	-	5.8
AWG2-35-205-A	0.035/0.089	205	Angled	Between 0.432/1.1 and 1.000/2.54cm	5.8
AWG2-35-260-A	0.035/0.089	260	Angled	Between 0.432/1.1 and 1.000/2.54cm	5.8
AWG2-35-450-A	0.035/0.089	450	Angled	Between 0.432/1.1 and 1.000/2.54	5.8



Accessories

1. Wire guide locking FS-WL-O-S
2. Wire guide locking FS-WL-P-S

Wire guide locking accessories (FS-WL) are designed to provide a wire (ACRO, AWG2, DWG, FS-LT, FSW, HYB, RR, THSF, MET, METII) lock to secure wire placement.

The Wire Guide Locking Device (FS-WL) device is an accessory to be used with endoscopic biliary devices to lock the wire guide(s) in place during ERCP procedures. This device is used as described in the steps of the graphics shown in the Instruction for Use (IFU):

1. Press the upper portion of the Wire Guide Locking Device to allow bottom portion of device to open for easy placement onto the endoscope
2. Attach the Wire Guide Locking Device to endoscope accessory port.
3. Lock Wire Guide onto Wire Guide Locking Device to allow Wire Guide to remain in place for use of other devices if necessary.
4. Once procedure is complete, unlock Wire Guide from Wire Guide Locking Device and remove from endoscope accessory port.

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Stripcoat Length (CM)
HYB-48015-A	0.089/ 0.035	480	Straight	15
HYB-48015-A	0.089/ 0.035	480	Angled	15
HYB-48025	0.089/ 0.035	480	Straight	25
HYB-48025-A	0.089/ 0.035	480	Angled	25

HYB

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Stripcoat Length (CM)
FSW-35	0.089/ 0.035	205	Straight	5

FSW

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Stripcoat Length (CM)
FS-LT-35-205	0.089/ 0.035	205	Loop	N/A
FS-LT-35-260	0.089/ 0.035	260	Loop	N/A
FS-LT-35-480	0.089/ 0.035	480	Loop	N/A

FS-LT

MET

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
MET-25-480	0.064/ 0.025	480	Straight	5
MET-25-480-A	0.064/ 0.025	480	Angled	5
MET-35-480	0.089/ 0.035	480	Straight	5
MET-35-480-A	0.089/ 0.035	480	Angled	5

METII

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
METII-21-260	0.053/ 0.021	260	Straight	~5
METII-21-480	0.053/ 0.021	480	Straight	~5
METII-25-260	0.064/ 0.025	260	Straight	5
METII-25-480	0.064/ 0.025	480	Straight	5
METII-35-260	0.089/ 0.035	260	Straight	5
METII-35-480	0.089/ 0.035	480	Straight	5
METII-35-480-A	0.089/ 0.035	480	Angled	5
METII-35-600E	0.089/ 0.035	600	Straight	5


INSTRUCTIONS FOR USE

1. Prior to removing wire guide from holder, flush with 30 cc of sterile water. (See fig.)
2. Flush endoscope accessory channel and/or lumen of device with sterile water, **then insert wire guide floppy end first.**

Note: For best results, wire guide should be kept wet, if applicable.

3. Fluoroscopically monitor wire guide advancement in ductal system.
4. Proceed with exchange of compatible wire-guided accessories over wire guide.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

INSTRUCTION FOR USE

Endoscopic Wire Guides (DWG-25-260-S; -25-260-A-S; DWG-35-260-S; DWG-35-260-A-S)

Endoscopic Wire Guides (RR-18-260; -18-260-A)

Endoscopic Wire Guides with adapter Pin-Vise (RR-18-480; -18-480-A; -18-480-MUSC)

Endoscopic Wire Guides (THSF-35-400; -18-480; -21-480; -25-480; -35-480)

INTENDED USE

This device is used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Contraindications include those specific to ERCP.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with ERCP include, but are not limited to: perforation, allergic reaction to contrast or medication.

PRECAUTIONS

Use of this wire guide with metal tip ERCP devices may result in damage to the external coating and/or tip of the wire guide (except THSF wire guides).

WARNING

Remove this wire guide before applying electrical current.

TECHNICAL SPECIFICATIONS

DWG

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
DWG-25-260-A-S	0.064/ 0.025	260	Angled	260
DWG-25-260-S	0.064/ 0.025	260	Straight	260
DWG-35-260-A-S	0.089/ 0.035	260	Angled	260
DWG-35-260-S	0.089/ 0.035	260	Straight	260

RR

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
RR-18-260	0.046/ 0.018	260	Straight	N/A
RR-18-260-A	0.046/ 0.018	260	Angled	N/A
RR-18-480	0.046/ 0.018	480	Straight	N/A
RR-18-480-A	0.046/ 0.018	480	Angled	N/A
RR-18-480-MUSC	0.046/ 0.018	480	Straight	N/A

THSF

Product Code/Model	Wire Guide Diameter Cm/ Inch	Wire Guide Length (CM)	Wire Guide Tip Type	Slipcoat Length (CM)
THSF-18-480	0.046/ 0.018	480	Straight	N/A
THSF-21-480	0.053/ 0.021	480	Straight	N/A
THSF-25-480	0.064/ 0.025	480	Straight	N/A
THSF-35-400	0.089/ 0.035	400	Straight	N/A
THSF-35-480	0.089/ 0.035	480	Straight	N/A


INSTRUCTIONS FOR USE

1. Prior to removing wire guide from holder, flush with 30 cc of sterile water. (*See fig.*)
2. Flush endoscope accessory channel and/or lumen of device with sterile water, **then insert wire guide floppy end first.**

Note: For best results, wire guide should be kept wet.

3. Fluoroscopically monitor wire guide advancement in ductal system.
4. Proceed with exchange of compatible wire-guided accessories over wire guide.

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

INSTRUCTION FOR USE

Wire guide locking accessories FS-WL-O-S; FS-WL-P-S.

INTENDED USE

This device is an accessory to be used with endoscopic biliary devices to lock the wire guide(s) in place during ERCP procedures.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Notify Cook for return authorization.

Do not use this device for any purpose other than the stated intended use.

Store in a dry location, away from temperature extremes.

Do not reuse, reprocess or resterilize this device.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

None known.

POTENTIAL COMPLICATIONS

This device is used with Cook biliary devices. Read appropriate Instructions for Use for complications associated with the procedure and device being used.

TECHNICAL SPECIFICATIONS

Product Code/Model	Scope Compatibility
FS-WL-O-S	Olympus / Fujinon
FS-WL-P-S	Pentax

INSTRUCTIONS FOR USE

Apply water-soluble lubricant to top of cap as needed to ease passage of accessories.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Please notify Cook for return authorization.

Cook Medical Europe Ltd

O'HALLORAN ROAD

NATIONAL TECHNOLOGY PARK

LIMERICK, IRELAND

TEL: +353 61 334440 FAX: +353 61 334441

WWW.COOKMEDICAL.COM

SHELF LIFE

The "use by" date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

These devices have a three (3) year shelf life.

STERILIZATION DETAILS

The devices included in this technical data summary are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . Do not use the sterile product if there is doubt as to whether the product is sterile.

STORAGE RULES

These devices should be stored in a dry location, away from temperature extremes. ACRO wire guides should also be stored away from sunlight.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

Doris A. Hawks
Doris A. Hawks
Global Regulatory Affairs Specialist
Wilson-Cook Medical, Inc.

June 5, 2019
Date

State of North Carolina, Forsyth COUNTY

I, Tiffany A. Thomas a Notary Public certify that,

Doris A. Hawks personally came before me and acknowledged that she is Regulatory Affairs Specialist of Wilson-Cook Medical Inc. and acknowledged on behalf of the corporation, the due execution of the foregoing instrument.

Witness my hand and notarial seal this 5 day of JUNE 2019.

Tiffany A. Thomas
Notary Public

My Commission expires: September 5, 2022

Tiffany A Thomas
Notary Public, North Carolina
County of Guilford
My Commission Expires September 5, 2022

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