NIALL SHEEHY

NOTARY PUBLIC

6/7, Glentworth St.,
Limerick
Ireland

Telephone +353 61 414355


Facsimile + 353 61 414738

Email nsheehy@sellors.ie

TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glentworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with certificates as issued from time to time by Cook Medical Europe and that the original document hereunto annexed namely, **Instructions for Use-Endoscopic Needles** upon which I have signed my name and affixed my Seal of Office is a true and faithful original and that the original appears to me to be genuine.

AND I HAVE CERTIFIED as herein at the request of Cook Medical Europe having its place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.

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9. Seal / stamp:
Sceau / timbre
Sello / timbre 

10. Signature:
Signature:
Firma: *[Signature]*

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document. It does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.aifm.ie

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
Niall Sheehy
NIALL SHEEHY
Notary Public
6/7 Glentworth St., Limerick,
Commissioned For Life

04-12-2020

Williams Cystoscopic Injection Needle

(090001, 090001-S10, 090001-S26, 090001-S34, 090001-S9)


090001, 090001-S26, 090001-S34, 090001-S10, 090001-S9

The Williams Cystoscopic Injection Needle is used for cystoscopic injection into urethral strictures, bladder neck contracture, and the bladder wall. Williams Cystoscopic Injection Needle come in diameters 23 G /0.64 mm and 25 G /0.52 mm. The Williams Cystoscopic Injection Needle is a semi-rigid injection needle. French size, gage, length, and tip length will vary by configuration.


Parameters of FLLA (for 090001-S26 only):

Diameters: A – 0.050 inch/ 1.270 mm (± 0.002 inch/0.051mm); D – 0.210 inch/5.334mm (± 0.005 inch/0.127mm); E – 0.344 inch/ 8.738mm (± 0.005 inch/0.127mm);
Length: B – 0.792 inch/ 20.117 mm (± 0.01 inch/0.254mm); C – 0.060 inch/1.524mm (± 0.005 inch/0.127mm); F – 0.172 inch/ 4.369mm (± 0.005 inch/0.127mm); G – 0.281inch/7.137mm (± 0.005 inch/0.127mm); H – 0.312 inch/7.925mm (± 0.005 inch/0.127mm).

Parameters of grip wing fitting:

Lengths: A – 1.222 inch/ 31.038mm (± 0.01 inch/0.254mm); B – 0.60 inch / 15.24mm (± 0.005 inch/0.127mm).

Order number	Needle		Sheath	
	Diameter G/mm	Length mm	Diameter Fr/mm	Length cm
Tolerance	Indicated in the technical characteristics above	$\pm 10\%$	$\pm 10\%$	$\pm 10\%$
090001	23/0.64	8	5.0/1.67	35
090001-S26	23/0.64	8	5.0/1.67	45
090001-S34	23/0.64	8	3.7/1.23	35
090001-S10	23/0.64	8	3.7/1.23	35
090001-S9	25/0.52	8	5.0/1.67	35

MATERIALS

Component	Material	Patient Contacting
090001-S26		
Needle	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0.03; Cu-Max 1	Direct contact
Needle Catheter (sheath)	Polyurethane Radiopaque #20493	Direct contact
Leur lock adapter	Polyamide 6 # 40340 Blue colorant: CI Pigment Blue 29 CI Pigment Green 7	Indirect contact

Connector cap	Polyamide 6 # 40340 Blue colorant: CI Pigment Blue 29	Indirect contact
090001, 090001-S9		
Needle	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0.03; Cu-Max 1	Direct contact
Needle Catheter (sheath)	Polyurethane Radiopaque # 20113 Blue colorant: Titanium dioxide CI Pigment Blue 29 CI Pigment Blue 15 Hub: Polyurethane # 20060 White colorant: Titanium dioxide CI Pigment Blue 28	Direct contact
Coil	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0.03; Cu-Max 1	Direct contact
Grip wing fitting	Polyurethane # 20060 White colorant: Titanium dioxide CI Pigment Blue 28	Indirect contact
090001-S10, 090001-S34		
Needle	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0.03; Cu-Max 1	Direct contact
Needle Catheter (sheath)	Polyurethane Radiopaque # 20113 Blue colorant: Titanium dioxide CI Pigment Blue 29 CI Pigment Blue 15 Hub: Polyurethane # 20060 White colorant: Titanium dioxide CI Pigment Blue 28	Direct contact

Caution: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

INTENDED USE

Williams Cystoscopic Injection Needle used for cystoscopic injection into the urethra, bladder neck, and bladder wall.

CONTRAINDICATIONS

The device is contraindicated in the presence of conditions which create unacceptable risk to the patient during cystoscopy.

PRECAUTIONS

Do not force the device during the procedure; stop and determine the cause of resistance before proceeding.

This product is intended for use by physicians trained and experienced in cystoscopic techniques. Standard techniques should be employed.

POTENTIAL ADVERSE EVENTS

Complications associated with cystoscopy procedures may include, but are not limited to:

- Bleeding
- Urinary tract perforation
- Stricture
- Urinary retention

Additional complications associated with the urological injection needles may include, but are not limited to:

- Inadvertent puncture with tissue damage or bleeding
- Infection
- Inadvertent injection of the therapeutic agent.

INSTRUCTION FOR USE

1. Securely attach a pre-filled syringe to the Luer lock fitting on device, then flush air from the needle. Check for patency and leaks.
2. With the cystoscope already placed in the urethra, introduce the needle into the working channel of the cystoscope.
3. Advance the needle in short increments until it is visualized exiting the cystoscope.
NOTE: Needle should be advanced no further than 5 cm from the end of the scope.
4. Cystoscopically visualize the desired injection site and maintain position of cystoscope.
5. Advance the needle into the intended injection site.
6. Perform injection(s) to desired location(s).
7. Upon completion of injections, carefully remove the needle from the cystoscope.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PERFORMANCE CHARACTERISTICS

Endoscopic needles of various types, in packaging are stable at temperatures in the range of +10 to +40 °C.

Endoscopic needles of various types, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Endoscopic needles of various types are stable under sudden changes in temperature.

Endoscopic needles of various types are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Endoscopic needles of various types, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

BIOCOMPATIBILITY TESTING

Product Name/Description	Biocompatibility Test
Williams Cystoscopic Injection Needle	Cytotoxicity, Sensitization, Intracutaneous reactivity, Acute systemic toxicity, Hemocompatibility, Risk Assessment

FUNCTIONAL (TECHNICAL) TESTING

Test Description	Results
Burst test	Pass
Needle Elasticity	Pass
Needle bending strength	Pass
Tensile test	Pass

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SHELF LIFE

The “use by” date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Endoscopic Needles of Various Types - have a 3 year expiration date.

STERILIZATION DETAILS

Endoscopic Needles of Various Types are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

Endoscopic Needles of Various Types should be stored in a dry location, away from temperature extremes.

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USAGE CONDITIONS

Endoscopic Needles of Various Types are used in the patient care medical institutions.

SUPPLY DETAILS

Endoscopic Needles of Various Types are sold in peel-away packages together with the instructions for use. Peel-away packages are put in the corrugated board boxes. Boxes are being taped with the control tape to cover the sides or along the perimeter.

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
LABELING

Endoscopic Needles of Various Types are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

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- Company-manufacturer and manufacturing address;
- Device name (trade mark);
- Serial number;
- marking about CE certification;
- marking about sterilization;

- Address of the EU representative;
- Trade mark of the manufacturer;
- manufacturing date;
- Sizes;
- Schematic picture;
- Lot number;
- marking for "Single use only"
- marking for "Attention. See instruction for use"
- marking for "Keep dry"
- marking for "Keep away from sunlight"
- Shelf life

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For EU countries - in accordance with the guidance document by World Health Organization.

For US countries – in accordance with the US Department of Labor of Occupational Safety and Health Administration (OSHA)

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Facsimile + 353 61 414738

Email nshcehy@sellors.ie

*TO ALL TO WHOM THESE PRESENTS SHALL COME, I, Niall Sheehy, Notary Public, duly authorised, admitted, and sworn, practicing at 6/7 Glentworth St., Limerick, Republic of Ireland, DO HEREBY CERTIFY that I am familiar with certificates as issued from time to time by Cook Medical Europe and that the original document hereunto annexed namely, **Instructions for Use-Endoscopic Needles** upon which I have signed my name and affixed my Seal of Office is a true and faithful original and that the original appears to me to be genuine.*

*AND I HAVE CERTIFIED as herein at the request of Cook Medical Europe having its place of business at O'Halloran Road, National Technology Park, in the city of Limerick, Ireland, and for the benefit of the proper authorities in **Russian Federation** and for all to whom it shall come.*

IN FAITH AND TESTIMONY whereof I have hereunto set my hand and affixed my Seal of Office.

Niall Sheehy

*Niall Sheehy
Notary Public*

9th December 2020

Date

NIALL SHEEHY
Notary Public
6/7 Glentworth St., Limerick.
Commissioned For Life

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Pays/Páis:		IRELAND	
This public document Le présent acte public / El presente documento público			
2. has been signed by a été signé par ha sido firmado por		Niall Shreehy	
3. acting in the capacity of agissant en qualité de quien actúa en calidad de		Notary Public	
4. bears the seal / stamp of est revêtu du sceau / timbre de y está revestido del sello / timbre de			
Certified Attesté / Certificado			
5. at à / en	Cork	6. the le / el día	17-12-2020
7. by par / por	Dept. of Foreign Affairs & Trade		
8. No sous no bajo el número	4073142020		
9. Seal / stamp: Sceau / Sello	10. Signature: Signature: Firma:		

This Apostille certifies the authenticity of the signature and the capacity of the person who has signed the public document, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued. To verify the issuance of this Apostille, see www.austlii.edu.au/au/other/dfat/ie

INSTRUCTIONS FOR USE

Endoscopic Needles

PURPOSE

Intended for insertion through a compatible endoscope to perform various operations during endoscopy.

DEVICE DESCRIPTION


Table of conversions

Parameter	Conversion
1/3 Fr	1 mm
1 Inch	25.4 mm

Technical characteristics of the needles mentioned as components in this Technical Data Summary:

Connection of the needle tip and tube should not be broken under minimum force exercised on the tip and tube in the direction of their separation, which equals 40 N. Resilience: after force 10.0 N, maximum deviation – 0.45 mm. Straightness: deviation from tube and tip concentricity should not exceed more than 0.2 mm. Maximum allowed deviation from the tip axis 3°. When needle tip is deviated at the distance of 17.5 mm – needles should not break. Blunting radius of the working parts of the instruments should be no more than 0.03 mm. Roughness of the surface - no more than 0.16 mkm.


Tolerances for needles: 16 G /1.65 mm (± 0.4 mm); 19 G /1.07 mm (± 0.03 mm); 20 G /0.89 mm (± 0.03 mm); 21 G / 0.82 mm (± 0.01 mm); 22 G /0.71 mm (± 0.02 mm); 23 G /0.64 mm (± 0.03 mm); 25 G /0.52 mm (± 0.01 mm).

NIAL SHEEHY
Notary Public
6/7 Glentworth St., Limerick.
Commissioned For Life

09-12-2020

EchoTip Ultra Celiac Plexus Neurolysis Needle (ECHO-20-CPN)


ECHO-20-CPN

The EchoTip Celiac Plexus Neurolysis Needle (ECHO-20-CPN) is comprised of a 20G/0.89 mm needle enclosed in a sheath and with a handle to allow the user to adjust the extension of the needle and sheath. The distal end of the needle has a pencil point/cone tip with four side openings toward the tip to administer a uniform radial spray of neurolysing agent, diameter of the openings 0.038 inch/0.97 mm, distance between openings of 1.01 inch/25.65mm. The stylet has a diameter of 0.020 inch/0.508mm (+/- 0.001 inch/0.0254mm) length of stylet – 170.5cm (+3mm/-2mm). The device has a 6 Fr/2 mm outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an ultrasound endoscope and the handle allows the needle to be advanced through the sheath by a distance of 0- 8cm and sheath adjustment length of 5 cm. The tip of the needle is dimpled diameter 0.0045 inch/0.1143mm (± 0.0015 inch/0.0381mm) creating an echogenic patterning effect to enhance echogenicity of the tip. Laser cut, 660 dimples ($\pm 3\%$) allow for needle tip visualisation under ultrasonography.

Dimple parameters:

Diameter - 0.0045 inch/0.1143mm (± 0.0015 inch/0.0381mm)

A – 9.5 mm; B – 1.0 mm (± 0.75 mm)

The EchoTip Celiac Plexus Neurolysis Needle is used to inject neurolysing agents into the celiac plexus under the guidance of Endoscopic ultrasonography. The needle is advanced endoscopically through the stomach wall and into the Celiac Plexus, which is adjacent to the Aorta and inferior to the Super Mesenteric Artery, and Celiac Artery (trunk). The procedure will be performed under the guidance of Endoscopic ultrasonography. Injection allows neurolysis to take place, providing temporary pain relief for patients suffering from inoperable malignancies, for example, pancreatic cancer. Minimum endoscope accessory channel – 2.8 mm.

Product Code/Model	Sheath		Needle		Stylet	
	Diameter Fr/mm	Length cm	Diameter G/mm	Length cm	Diameter inch/mm	Length cm
Tolerance	+/- 10%	+/- 10%	Indicated in the technical characteristics above	+/- 10%	+/- 0.001 inch/ 0.0254mm	+3/-2 mm
ECHO-20-CPN	6/ 2	142.2	20/ 0.89	8	0.020 / 0.508	170.5


Parameters of handle:

A – 4.528 inch/ 115.011mm ($\pm 0.005/0.127$ mm); B – 0.591 inch/ 15.011mm ($\pm 0.003/0.076$ mm);
C – 0.504 inch/ 12.802 mm ($\pm 0.003/0.076$ mm); D – 0.384 inch/ 9.754mm ($\pm 0.015/0.381$ mm); E – 0.110 inch/ 2.794mm ($\pm 0.005/0.127$ mm).

Parameters of Luer-lock:

Length: A – 0.661 inch/ 16.789mm (± 0.005 inch/0.127mm); B – 0.285 inch/ 7.239 mm (± 0.005 inch/0.127mm); C – 0.171 inch / 4.343 mm (± 0.005 inch/0.127mm);
Diameters: D – 0.215 inch / 5.461 mm (± 0.002 inch/0.051mm); E – 0.200 inch/ 5.08 mm (± 0.000 inch/0.000mm; -0.005 inch/0.127mm); F – 0.0720 inch / 1.829 mm (± 0.005 inch/0.127mm).

MATERIALS

Component	Material	Patient Contacting
Needle assembly with stylet		
Needle	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0,045; S-Max 0,03; Cu-Max 1. Coating: S # 63148-62-9 Polydimethylsiloxane Dow Corning 360	Direct contact
Needle Catheter (sheath)	Polytetrafluoroethylene #110852 Sheath fixator: Polycarbonate Lexan 241R, Nickel CDA 360 Hub: Brass CDA 360	Direct contact
Stylet	AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0,03; Cu-Max 1.	Direct contact
Stylet cap	Polycarbonate: Lexan 121R-112 Colorants: Blue – Mevopur-Blue Grey - Mevopur-Grey	No contact
Handle assembly		
Handle	Polystyrene STYRON 484 Grey colorant: Q757-3-4, Blue colorant: DynaFlex D3202-1000-03 Luer lock adapter: Nickel plated brass CDA 360 Proximal stop: Polystyrene Dow 484W, Polycarbonate Lexan 121 Colorant: white WH1164 O-ring: BUNA-N	No contact
Locking ring	Polystyrene STYRON 484 Grey colorant: Q757-3-4 Hub: Brass CDA 360 Screw: AISI 304 Stainless Steel: C-Max 0.08; Cr-18 – 20; Fe-66.345-74; Mn-Max 2; Ni-8 – 10.5; P-Max 0.045; S-Max 0.03; Cu-Max 1.	No contact

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

INTENDED USE

EchoTip Ultra Celiac Plexus Neurolysis Needle (ECHO-20-CPN) is used to deliver neurolysing agents to the celiac plexus under guidance by endoscopic ultrasound.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired injection site.

Contraindications to celiac plexus neurolysis include, but are not limited to: coagulopathy, use of anticoagulants, aspirin or nonsteroidal anti-inflammatory drugs (NSAIDs).

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with gastrointestinal endoscopy and celiac plexus neurolysis include, but are not limited to: airway obstruction, allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, discomfort, fever, gastroparesis, hemorrhage, hypotension, infection, inflammation, lower extremity weakness, nausea, paraplegia, paresthesias, perforation, renal puncture, respiratory depression or arrest, sepsis, transient diarrhea, transient increasing pain and vomiting.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the endoscope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

A thorough understanding of the technical principles, clinical applications and risks associated with endoscopic ultrasound and celiac plexus neurolysis is necessary before using this device. The use of prophylactic antibiotics is recommended for endoscopic ultrasound- guided celiac plexus neurolysis. The use of general anesthesia may be needed in some cases.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

INSTRUCTION FOR USE

1. Identify the desired injection site by endoscopic ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting on the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 1). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring (Fig. 2). Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, to the desired injection site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional neurolysing injections are required later.
9. If required, inadvertent puncture of a blood vessel may be determined by attempting to aspirate blood. A Luer lock syringe may be attached to the Luer fitting on the needle handle for this purpose and then removed.
10. Attach the Luer lock fitting of a previously prepared syringe, containing local anaesthetic and neurolysing agents, securely to the Luer fitting on the needle handle. Instill volume required while observing the ultrasound view.
11. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
12. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
13. For additional neurolysing injections at alternative sites on the celiac plexus, disconnect the syringe from the Luer fitting on the needle handle and gently reinsert the stylet. Repeat Steps 2 through 12 of the "Instructions for Use". **Note:** Prior to reinserting stylet, wipe with saline or sterile

water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.


Fig. 1


Fig. 2

PERFORMANCE CHARACTERISTICS

Endoscopic needles, in packaging are stable at temperatures in the range of +10 to + 40 ° C.

Endoscopic needles, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Endoscopic needles are stable under sudden changes in temperature.

Endoscopic needles are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Endoscopic needles, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

BIOCOMPATIBILITY TESTING

Product Name/Description	Biocompatibility Test
EchoTip Ultra Celiac Plexus Neurolysis Needle (ECHO-20-CPN)	Cytotoxicity Study, Sensitization Study, Irritation / Intracutaneous Reactivity, Acute Systemic Toxicity, Intracutaneous Toxicity, Implantation, Evaluation of Equivalency and Biological Risk Assessment of PTFE stents

FUNCTIONAL (TECHNICAL) TESTING

Test Description	Results
Product Evaluation Testing	Pass
Compression Testing of Needle Tip	Pass
Tensile Testing of the Stylet Hub	Pass
Tensile Testing of Luer to Cannula Solder Joint	Pass
Tensile Testing of Strain Relief to Cannula Solder Joint	Pass
Tensile Testing of Sheath holder, Molded Hub to PTFE Tubing bond	Pass
Dimpling Validation	Pass
Product Dimension Checks	Pass

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated, Cook Ireland and Wilson-Cook Medical Inc. manufacture all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated and Wilson-Cook Medical Inc. manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2016) *Quality Management System*.

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SCHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9, floor 5, office 501.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Endoscopic Needles - have a 3 year expiration date.

STERILIZATION DETAILS

Endoscopic Needles are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

Endoscopic Needles should be stored in a dry location, away from temperature extremes.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

USAGE CONDITIONS

Endoscopic Needles are used in the patient care medical institutions.

SUPPLY DETAILS

Endoscopic Needles are sold in peel-away packages together with the instructions for use. Peel-away packages are put in the corrugated board boxes. Boxes are being taped with the control tape to cover the sides or along the perimeter.

PACKAGING PARAMETERS - of the product packaged in the peel-away pouch


RPN	Weight (gr)	Height (cm)	Width (cm)	Tolerance
ECHO-20-CPN	244.4	42.8	24.5	+/- 10%

LABELING

Endoscopic Needles are supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- Company-manufacturer and manufacturing address;
- Device name (trade mark);
- Serial number;
- marking about CE certification;
- marking about sterilization;
- Address of the EU representative;
- Trade mark of the manufacturer;
- manufacturing date;
- Sizes;
- Schematic picture;
- Lot number;
- marking for "Single use only"
- marking for "Attention. See instruction for use"
- Shelf life

REPAIR/MAINTENANCE

Endoscopic Needles packages are single use only, do not attempt to repair, carry out maintenance or reuse.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries - in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

Maryna Zhukava
Regulatory Affairs Specialist

03-Dec-2020
Date

COOK MEDICAL EUROPE LTD.
O'Halloran Road,
National Technology Park
Limerick, Ireland.
Company Registration No. 492272
VAT Registration No. 9776234K

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