19.02.19 131326

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APPROVED BY

Sophie REHAULT

Quality and Regulatory Affairs Director

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REHAULT

(seen exclusively to certify the above signature)

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USER MANUAL

CORAIL+ guiding sterile catheter for single use

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70196
28.10.2019

This user manual has been prepared to provide information for purposes of registration of a medical device in the territory of the Russian Federation. The information is provided in the amount required in accordance with the Rules for Registration of Medical Devices (Russian Government Decree No. 1416 dated 27.12.2012).

1. Product name:

CORAIL+ guiding sterile catheter for single use:

- CORAIL6F+;

- CORAIL8F+

- Instructions for use (supplied with each product).

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Development, production and sale of the device complies with the requirements of DIN EN ISO 13485:2016+AC:2017-07, EN ISO 13485:2016+AC:2016, ISO 13485:2016; the device bears CE-marking EC certificate number: 513975 MR2.

OKPD2 - 32.50.13.110 - Syringes, needles, catheters, cannulas and similar instruments;

Class depending on the degree of potential risk of medical use in accordance with the medical devices nomenclature classifier approved by Order of the Ministry of Health of the Russian Federation No. 4n - 3 dated June 6, 2012;

Device classification according to the Medical Device Nomenclature Classification approved by the Order of the Ministry of Health of the Russian Federation of June 6, 2012, No.4n – 131670.

Type of contact with the human body: short-term (less than 24 hours) contact with circulating blood.

This device is sterile. For single use.

Conditions of use: medical and medical-preventive health care institutions, only as prescribed by a medical specialist.

2. Purpose

The device is intended to facilitate the introduction of microcatheters for therapeutic and diagnostic use, as well as for percutaneous transluminal passage and placement of diagnostic / interventional devices (balloon dilating catheter, stent and delivery system, or embolization device) in the vascular system.

3. Device special properties

CORAIL+ balloon tipped guiding catheters are a coaxial catheter with a low pressure distal balloon, capable of reaching various diameters depending on the volume of injected liquid. The side entry is intended to inflate and deflate the dilation balloon, and the internal lumen provides a passage for various catheters. CORAIL+ balloon tipped guiding catheter is radiopaque and thermoformable and fulfils its role as a guiding catheter by guiding various catheters inside the vascular system. In addition, the balloon allows for temporary occlusion during various operations.

Technical characteristics of the devices are given in clause 4. Each device is packed in an individual package, each individual package is packed in a cardboard box, which also contains a user manual. For more information on packaging and labeling, see clause 5. Appearance of the device is shown in figures 1-4.


Figure 1: Device appearance without packaging CORAIL6F+


Figure 2: Appearance of the device without package CORAIL8F+

Figure 4: Appearance of the CORAIL+ device in the final package for sterilization


Figure 5: Device design

4. Technical specification

Table 1 – Materials of the devices

Component	Raw material CORAIL6F+	Raw material CORAIL8F+	Raw material CORAIL8F+LT3 0	Suppliers
Manipulator with Y- connector	Polypropylene PPH5060			SPEEDPLAST, Spain
Counter-current node ¥	Thermoplast TF6ST			KRAIBURG TPE, Germany
	HDPE Alcludia 4810			Gazechim Groupe, France
	Bismuth trioxide, CAS # 1304-76-3			5N PLUS FRANCE, France
External rigid tube	Polypropylene PPH5060			SPEEDPLAST, Spain

tube	Pigment Red 254 - CAS # 84632-65-5		Plasticouleurs, France
	Bismuth trioxide, CAS # 1304-76-3		5N PLUS FRANCE, France
External flexible tube	Polypropylene PPH5060		SPEEDPLAST, Spain
	Thermoplast TF4ST		KRAIBURG TPE, Germany
	Bismuth trioxide, CAS # 1304-76-3		5N PLUS FRANCE, France
External tip	-	Polypropylene PPH5060	SPEEDPLAST, Spain
	-	Thermoplast TF4ST	KRAIBURG TPE, Germany
	-	Pigment Red 254 - CAS # 84632-65-5	Plasticouleurs, France
	-	Yellow pigment (AZO Pigment Yellow) CAS # 66511-62-6	Plasticouleurs, France
Internal tube	Low density polyethylene (LDPE) Lacqtene1070		SPEEDPLAST, Spain
Internal flexible tube	Polypropylene PPH5060		SPEEDPLAST, Spain
	Thermoplast TF4STE	-	KRAIBURG TPE, Germany
	-	Chronoprene 25A	ADVANSOURCE Biomaterials, USA
	-	Melcometer orange Melcomaitres® 28413 PPHF	MELCOPLAST, France
	Yellow pigment (AZO Pigment Yellow) CAS # 66511-62-6		Plasticouleurs, France
	Pigment Red 254 - CAS # 84632-65-5		Plasticouleurs, France
Inner layer	Polytetrafluoroethylene (PTFE)	-	Accelent Inc. USA
	-	High density polyethylene HDPE Alcudia 4810	Gazechim Groupe, France
Braid ¥	Tungsten wire CAS # 7440-33-7	-	LUMA METALL, Sweden
	-	stainless steel AISI 304V	ALTHOFF-LOTTERS, Germany
Balloon	Thermoplast TF6ST		KRAIBURG TPE, Germany
ORX radiopaque marker ¥	Gold 99% (gold99,99), impurities 1% (not more) Lead - 0,003, Iron - 0,004, Antimony - 0,001, Bismuth - 0,002, Copper - 0,007, Silver - 0,008		Johnson Matthey, UK
Mandrine	stainless steel AISI 304V		ALTHOFF-LOTTERS, Germany
Protective rim	Polyethylene (Rigidex 5502 GA)		INEOS Olefins & Polymers, UK

(¥): components, not intended for direct or indirect contact with the body.

The differences between CORAIL6F + and CORAIL8F + are shown in Table 1. The inner tube and inner layer are made from different raw materials and pigments.

Table 3: Device technical characteristics

	CORAIL6F+	CORAIL8F+
Total length, mm	1035 ± 15	1035 ± 15
Effective length, mm	940 (+33/-20)	960 (+33/-20)
Inner diameter, mm	1.235	1.9
Outer diameter, mm	2.00	2.76
Wall thickness, mm	0.38	0.43
Overall dimensions of the balloon, (Max. Ø x length) mm	11x14	12x11
Maximum balloon volume, ml	1.20	
Maximum allowable pressure in the balloon, Bar	1	
Tensile strength, N	11,4 ± 0,8	
Overall dimensions of the distal tip (Length x Diameter)	2x1,52	2x2,25
Inner diameter (Main lumen), mm	1.23	1.90
Weight, g	5,65±0,29	
Connector type	Luer lock according to ISO 594	

Unspecified permissible values: ± 1%

Appearance of the device: The device surface has no foreign inclusions, bubbles or streaks. The surface is smooth, has no technological and surface defects and provides minimal trauma to vessels when being used. The catheter is biologically safe and compatible with biological human body fluids. Cannula and other parts of the catheter do not allow fluid leakage. The device is not implantable and is removed from the human body after being used.

5. Device range:

CORAIL+ guiding sterile catheter for single use:

- CORAIL6F+;
- CORAIL8F+;
- Instructions for use (supplied with each product).

6. Marking and packaging

1 CORAIL+ guiding sterile catheter for single use is supplied packaged in an individual Tyvek pack which is the final package for sterilization and in an individual cardboard package. Each device is supplied with a user manual. Each cardboard package contains 2 devices in the final package for sterilization.

6.1. Technical characteristics of the final package for sterilization:

One side of the packs is made of Tyvek® 2 FS UNCo 7502918M coating and the other side is made of PET/PE 65 in accordance with standard NF EN ISO11607-1 (2009). The final package for sterilization is glued together.

Table 4: Technical characteristics of the package

Basic properties	Results
Pack overall dimensions (LxWxT), mm	1260±5x90±5x0.3
Width of adhesive layer on the long side, mm	8
Width of adhesive layer on the wide side, mm	8
Thickness of adhesive layer, mm	0.3
Tyvek paper - Unit weight, g/m ²	59.5 (56.5- 62.5)
Tyvek paper - Thickness, µm	155 (75- 235)
Tyvek Paper - Porosity (Gurley Method)	22 s (9-35)
Tyvek Paper - Porosity (Bendtsen Method)	529 ml/min (155-885)
Tyvek paper - Mullenburst test	925 kPa (760-1,110)
Tyvek paper - MD breaking strength test	156 N/2.54 cm (135-177)
Tyvek paper - XD breaking strength test	157 N/2.54 cm (127 - 187)
Tyvek paper - Elmendorf MD breaking test	2.8 N (1.7-4.0)
Tyvek paper - Elmendorf XD breaking test	3.7 N (2.8-4.6)
Tyvek paper - Stratification	2.70 N/2.54 cm (2.05-3.25)
Tyvek paper - Hydrostatic flow pressure	145 cm H ₂ O (125- 165)
PET film - Thickness	65 µm ± 5%
PET film - Total weight	63 g/m ²

PET film - Breaking test	30 - 55 N/mm ²
PET film - Elongation at break	50- 100 %
PET film - Elongation at break	80- 160 %
PET film - O ₂ permeability	90 m ³ /m ² /bar
PET film - CO ₂ permeability	345 m ³ /m ² /bar
PET film - N ₂ permeability	50 m ³ /m ² /bar
PET film - Water vapor permeability	< 3.0 g/m ²
PET film - Sealing parameters	115-150°C
Sealing	> 1.5 N / 15 mm
Final packaging manufacturer	SUDPACK MEDICA, France

Unspecified permissible values: ±1%

Table 5: Technical characteristics of the individual cardboard package

Basic properties	Results
Box overall dimensions (LxWxT), mm	1260±4 x 90±3 x 40±3
Cardboard density	365 g/m ²
Material	Cardboard
Manufacturer	IMPRIMERIE LLP/BUIGNET Address: 29, avenue des Chataigniers – 95150 Taverny

Unspecified permissible values: ±1%

6.2. Appearance of the final package marking:

REF

CORAIL6F+

CONTENTS
1


CORAIL+ 6F

Cathéter Guide à Ballonnet/Guiding Catheter with Balloon
/Catéter Guia con Balon/Einführungs Ballonkatheter
/Catetere – guida a palloncino/Cateteres–guias com balão
/Cewniki prowadzące zbrojone z balonikiem
/Катетерные проводники с баллончиком
/Balonlu kilavuz kateter

100 cm

94 cm

2.03 mm/6F

Balloon inflation chart :

Vol (ml)	ø (mm)
0.1	5
0.3	8
0.5	9
0.7	10
1.0	11

Volume ballon max. :
Max. balloon volume :

1.20 ml

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1.20 mm (0.047")

ind.7

20170321

REF

CORAIL6F+

EAN 13:

3 700481 313381

CODE 128:

(01)03700481313381(17)220328(10)00283872

LOT

00283872

2017-03

2022-03

0459

STERILE

EO

BALT
EXTRUSION

Balt Extrusion
10, rue de la Croix-Vigneron
95160 Montmorency
France
Tél. +33 (0)1 39 89 46 41
Fax. +33 (0)1 34 17 03 46
www.balt.fr

REF

CORAIL6F+

(01)03700481313381(17)220328(10)00283872

LOT

00283872

REF

CORAIL6F+

(01)03700481313381(17)220328(10)00283872

LOT

00283872

REF

CORAIL6F+

(01)03700481313381(17)220328(10)00283872

LOT

00283872

REF

CORAIL6F+

(01)03700481313381(17)220328(10)00283872

LOT


00283872

Figure 5: CORAIL6F+ package and box marking

000 ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

70196
28.10.2019

Catheter Guide & Balloon/Guiding Catheter with Balloon
/Catheter Guia con Balon/Einführungs Ballonkatheter
/Catetere - guida a palloncino/Catetere-guia com balão
/Катетер-проводник с баллончиком
/Kateteri keuzut kateter

Max. balloon volume : **1.20 mL**
Volume ballon max. :

Inner diameter : **1.89 mm (.074")**
Diamètre interne : **1.89 mm (.074")**

Ball Extrusion
10, rue de la Croix-Vigneron
95160 Montmorency
France
Tél. +33 (0)1 39 89 46 41
Fax. +33 (0)1 34 17 03 40
www.ball.fr

Figure 6: CORAIL8F+ package and box marking;


Figure 7: Box layout

Marking on the final sterilization package is identical to marking on the individual cardboard package.

6.3. The marking contains:

- Quantity;
- Device name and model;
- Schematic representation of the device;
- The maximum volume of balloon;
- Reference number
- Serial number;
- Date of manufacture;
- Expiration date;
- Information on the manufacturer;

Symbol legends list

	Reference number
	Serial number
	Use before
	Date of manufacture
	CE label

STERILE EO	EO-sterilized
	Attention! Refer to user manual
	Prohibition on reuse
	Do not expose to sunlight
	Do not use if package is damaged
	Refer to user manual
	Keep dry!
EAN13:  3 700481 315538	European article number
	Manufacturer

Cardboard packages with the device are packed in a corrugated cardboard transport box. One transport box contains – no more than 10 cardboard packages.

6.4. Transport box marking contains:

- Device name and model;
- Information on the manufacturer;
- Quantity;
- Lot number;
- Serial number;
- Date of manufacture;
- Expiration date;
- "Do not expose to sunlight" symbol;
- "Keep dry!" symbol;

7. Transportation and storage

Transportation conditions: air temperature from -10°C to +35°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 - 1,060 hPa away from direct sunlight.

Storage conditions: store the device in the manufacturer's package in a dry, clean place, away from heat sources and direct sunlight, at controlled temperature. Shelf life is indicated

on the label of each device. Air temperature from +15°C to +27°C at relative humidity from 35% to 80% (non-condensing) and atmospheric pressure of 500 – 1,060 hPa.

Operating conditions: air temperature from +20°C to +25°C at relative humidity from 35% to 65% (non-condensing) and atmospheric pressure of 500 – 1,060 hPa away from direct sunlight. The device is temporary operated in the human body at the temperature from +36°C to +42°C. The device is operated in the human body for less than 24 hours.

8. Information on the device sterility

The device is supplied sterilized with Ethylene Oxide sterilization method.

9. Disposal information

The device shall be disposed in accordance with SanPiN 2.1.7.2790-10 for class B waste (epidemiologically hazardous waste). After the apparatus decontamination methods using physical methods and changing the waste appearance, excluding the possibility of their reuse, class B waste may be accumulated, temporarily stored, transported, destroyed and buried together with class A waste (epidemiologically safe waste approximated to solid household waste). A package of disinfected class B medical wastes shall have a mark indicating the completed waste disinfection.

10. Warranties

Shelf life - 5 years. Sterilization validity - 5 years.

Address for complaints submission:

Perinatal Medicine Center Limited Liability Company

Registered/Actual address: 109004, Moscow, Tovarishchesky per., 20, bldg. 4

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+33 (0)1 39 89 46 41

+33 (0)1 37 17 03 46

11. Data for development and production

Development, production and sale of the device complies with the requirements of DIN EN ISO 13485:2016+AC:2017-07, EN ISO 13485:2016+AC:2016, ISO 13485:2016; the device bears CE-marking EC certificate number: 513975 MR2.

12. General indications, contraindications and side effects

Indications:

Balloon tipped guiding catheter CORAIL+ fulfills its role as a guiding catheter by guiding various catheters inside the vascular system. In addition, the balloon allows for temporary occlusion during various operations.

Standard Luer connections provide a connection to the material.

Contraindications:

The use of CORAIL + catheters is contraindicated if the clinician concludes that a procedure using catheters can impair the patient's condition.

Not intended for embolectomy or angioplasty.

Not intended for application in coronary vessels.

Not intended for application in children or neonatals.

Potential complications:

Possible complications, connected with use of balloon catheters or with intravascular procedures include, without limitation, the following disorders:

- complications in the access site area;
- allergic reaction;
- puncture of the vessel or aneurysm;
- aneurysm rupture;
- death;
- embolism;
- hematoma at the access site;
- hemorrhage;
- infection;
- ischemia;
- neurological disorders;
- cerebral accident;
- transitory ischaemic attack;
- angiospasm;
- vessel dissection;
- vascular occlusion;
- vascular rupture;
- vascular thrombosis.


13. Application method.

Check the integrity of the package: it must not be opened or damaged. Open the package and remove the catheter from it, observing the aseptic rules. Check the condition of the catheter for damage.

13.1. Washing the balloon

- Fill the syringe with sufficient volume (10 ml), with a sterile radiopaque preparation (if necessary, diluted with sterile saline solution).
- Never use air or any gas to fill the balloon during washing and, especially, during the procedure.
- Fill the balloon with a radiopaque preparation after attaching the syringe to the side entry for inflation and deflation the balloon.

Place the catheter in such a way that the distal end of the catheter is pointing downward, as shown in the figure:

- Holding the syringe above the balloon, the end of which is directed downwards, inject the radiopaque preparation and draw air from the catheter deadspace.
- Repeat these steps several times until the air is completely removed.
- The balloon should be completely filled with radiopaque preparation, and the volume of the injected liquid should correspond to the diameter of the balloon.

13.2. Introduction of a balloon tipped guiding catheter CORAIL+

- Empty the balloon completely.
- Wash the lumen of the CORAIL+ catheter with sterile saline.
- Connect the valve Y-shaped connector to the CORAIL+ guiding catheter
- After making sure that the balloon is completely emptied, insert the washed CORAIL + catheter through the valve introducer using the removable protective cap of the balloon.
- Remove the protective cap of the balloon.
- Move the CORAIL+ guiding catheter up to the required area, performing the flushing.

13.3. Inflation the balloon

- After the procedure for washing the balloon, which is specified in paragraph 5.1.
- Connect the syringe, the volume of which corresponds to the required inflation diameter of the balloon.
- Slowly inject the radiopaque preparation (if necessary diluted with sterile saline solution) until the required diameter of the inflated balloon is reached (see device label).

13.4. Precautions for use.

- Do not use if the package is opened or damaged. If the package is intact, they are sterile.
- CORAIL+ catheters intended for single use only. Do not reuse. Any reuse of the device leads to a high risk of infection of the patient and a deterioration in the performance of the device.
- Do not resterilize these devices.
- Keep in a dry, dark place at room temperature.
- Observe the shelf life.
- These devices shall only be used by medical specialists in the sphere of interventional neuroradiology and/or interventional radiology.
- Observe the inflation volumes, never exceed the maximum inflation volume (see device label).

- Use only appropriate preparations for inflating the balloon (for example, 50:50 contrast-saline solution), it is recommended to use radiopaque preparations having the lowest viscosity to facilitate inflation and deflation.
- Balloon tipped guiding catheters CORAIL+ must be used under fluoroscopic control and with appropriate anticoagulants.
- If resistance occurs, never move the intravascular device without determining its causes by angiographic control. Movement in the event of resistance may result in damage to the catheter or injury to the patient.
- Observe user manual of preparations injected into a balloon tipped guiding catheters CORAIL+.
- Wash the guiding catheter throughout the procedure.
- Observe user manual of interventional catheters used with balloon tipped guiding catheters CORAIL+.

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