**Medtronic***Alleviating Pain · Restoring Health · Extending Life***CERTIFICATE OF AUTHENTICITY AND COPY**

State of Massachusetts
County of Essex

On this 22nd day of November, I certify that the attached document is a true, exact, complete, and unaltered copy made by me of the TRAcelet IFU M970663A001.

Signature of Document Custodian (Colleen Mullins)

(Official signature and seal of Notary)

Subscribed and sworn to before me this the
22nd day of November, 2016

Heather P Morose, Notary Public
Commonwealth of Massachusetts
My commission expires on June 22, 2023

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, USA
Tel: +1 763 514 4000

2016/NOV/18 at 7:51 a.m. Doc number: M970663A001 Rev. 1A [global]
Printspec x - 4.0 x 7.75 inches

Medtronic

TRAcelet™

Compression Device
Компрессионно устройство
Kompresní zařízení
Kompressionsenhet
Kompressionsgerät
Συσκευή συμπίεσης
Dispositivo de compresión
Kompressiooniseade
Puristuslaite
Dispositif de compression
Uređaj za komprimiranje
Kompressziós eszköz
Dispositivo di compressione
압축 장치
Užspaudimo įtaisas
Kompresijas ierīce
Druksysteem
Kompresjonsenhet
Urządzenie do kompresji
Dispositivo de compressão
Dispositivo de compressão
Dispozitiv de compresie
Компрессионное устройство
Zariadenie na kompresiu
Kompresijski pripomoček
Uređaj za kompresiju
Kompressionsanordning
Kompresyon Cihazı
Компресійний пристрій

Instructions for Use • Инструкции за употреба • Návod k použití •
Bruksanvisning • Gebrauchsanweisung • Οδηγίες χρήσης •
Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi •
Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • 사용 지침 •
Naudojimo instrukcijos • Lietošanas pamācība • Gebruiksaanwijzing •
Bruksanvisning • Instrukcja użytkowania • Instruções de utilização •
Instruções de utilização • Instrucțiuni de utilizare •
инструкция по применению • Pokyny na používanie • Navodila za uporabo •
Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları •
Інструкція з експлуатації • 使用說明

CE
0086

master:20160719, alldocMod:20161025, alldocMod2:20161024, alldocStyle_graph:20160504,
alldocStyle:20160504, master_pageset:20160504, languageLookup:20161110, table:20160113,
characterJoin:20110211, 3of9:20110211, 3of9_svg:20110211

2016/NOV/18 at 7:51 a.m. Doc number: M970663A001 Rev. 1A [global]
Printspec x - 4.0 x 7.75 inches

master:20160719, alldocMod:20161025, alldocMod2:20161024, alldocStyle_graph:20160504,
alldocStyle:20160504, master_pageset:20160504, languageLookup:20161110, table:20160113,
characterJoin:20110211, 3of9:20110211, 3of9_svg:20110211

Explanation of symbols on package labeling / Обяснение на символите от етикетите на опаковката / Vysvětlení symbolů uvedených na štítcích balení / Forklaring af symbolerne på emballagens mærkater / Erläuterung der Symbole auf dem Verpackungsetikett / Επεξήγηση των συμβόλων στη σήμανση της συσκευασίας / Explicación de los símbolos que aparecen en el etiquetado del envase / Pakendi siltidel olevate sümbolite tähendus / Pakkausetiketin symbolien selitykset / Explication des symboles des étiquettes sur l'emballage / Objašnjenje simbola na ambalažnoj etiketi / A csomagoláson látható szimbólumok jelentése / Spiegazione dei simboli presenti sulle etichette della confezione / 포장 레이블의 기호 설명 / Pakuotės etikečių simbolių reikšmės / Uz iepakojuma etiķetes attēloto simbolu skaidrojums / Verklaring van de symbolen op de verpakkingslabels / Forklaring av symboler på pakningsetiketten / Objaśnienie symboli znajdujących się na etykietach opakowania / Explicação dos símbolos nas etiquetas da embalagem / Explicação dos símbolos na documentação da embalagem / Legenda simbolurilor cuprinse în etichetele de pe ambalaj / Обяснение символов на этикетках упаковки / Vysvetlivky k symbolom na označení balenia / Razlaga simbolov na embalaži / Objašnjenje simbola na oznakama na pakovanju / Förklaring av symboler på förpackningsmärkningen / Ambalaj etiketi üzerindeki sembollerin açıklaması / Пояснения символов на маркуванні упаковки / 包装标签上的符号说明

Quantity / Количество / Množství / Antal / Menge / Ποσότητα / Cantidad / Kogus / Määrä / Quantité / Količina / Mennyiség / Quantità / 量 / Kiekis / Daudzums / Aantal / Antall / Ilość / Quantidade / Quantitate / Cantitate / Количество / Množstvo / Količina / Količina / Antal / Miktar / Кількість / 数量

Open here / Отворете тук / Zde otevřete / Åbnes her / Hier öffnen / Ανοίξτε εδώ / Abrir aquí / Avage siit / Avaa tästä / Ouvrir ici / Ovdje otvoriti / Itt nyílik / Aprire qui / 개봉 위치 / Atidaryti čia / Atvērt šeit / Hier openen / Åpnes her / Otwierać tutaj / Abrir aquí / Abrire aqui / Deschideți aici / Открывать здесь / Tu otvorit / Odrprie tukaj / Ovdje otvoriti / Öppna här / Buradan açın / Відкрити тут / 在此打开

For US audiences only / Само за САЩ / Pouze pro uživatele z USA / Gælder kun i USA / Gilt nur für Leser in den USA / Μόνο για πλάτες εντός των ΗΠΑ / Solo aplicable en EE. UU. / Ainult kasutajatele Ameerika Ühendriikides / Koskee vain Yhdysvaltoja / Ne s'applique qu'aux États-Unis / Samo za američko tržište / Csak egyesült államokbeli felhasználóknak / Esclusivamente per il mercato statunitense / 미국 내 사용자를 해당 / Tik JAV naudotojams / Tikai pircējiem ASV / Alleen van toepassing voor de VS / Gjelder kun USA / Dotyczy tylko odbiorców w USA / Aplicável apenas aos EUA / Aplicável apenas aos E.U.A. / Numai pentru clienții din SUA / Только для США / Len pre používateľov v USA / Samo za uporabnike v ZDA / Samo za korisnike u SAD / Gäller endast i USA / Yalnızca ABD'deki kullanıcılar için / Лише для користувачів із США / 仅供美国观众使用

Manufactured in / Произведено в / Vyrobeno v / Fremstillet i / Hergestellt in / Κατασκευάστηκε σε / Fabricado en / Tootmiskoht / Valmistuspaikka / Lieu de fabrication / Proizvedeno u / A gyártás helye / Fabbicato in / 製造 場所 / Pagaminta / Ražošanas vieta / Vervaardigd in / Produziert in / Miejsce produkcji / Fabricado em / Fabricado em / Fabricat în / Произведено в / Vyrobené v / Izdelano v / Proizvedeno u / Tillverkad i / Imalat yer i / Місце виготовлення / 製造地

IVM	Inflation volume maximum / Максимален обем на надуване / Maximální plnicí objem / Maks. oppustningskapacitet / Maximaales Befüllungsvolumen / Μέγιστος όγκος εμφύσησης / Volumen de inflado máximo / Maksimaalne täitmismah / Enimmäistäyttötilavuus / Volume de gonflage maximal / Maksimalni kapacitet napuhivanja / Maximális feltöltési térfogat / Volume massimo di gonfiaggio / 최대 팽창 용량 / Maksimalus pripildymo tūris / Maksimālais uzpildes tilpums / Maximaal vulvolume / Maksimālti fyllingsvolum / Maksymalna objętość napelniania / Volume máximo de insuflação / Volume máximo de insuflação / Volumul de umflare maxim / Максимальный объем раздувания / Maximálny objem nafuknutia / Največja prostornina napihovanja / Maksimalna zapremina za naduvavanje / Maximal fyllningsvolum / Maksimum şişirme hacmi / Максимальный объем роздувания / 最大膨張容積
L	Large / Голямо / Velký / Stor / Groß / Μεγάλη / Grande / Suur / Suuri / Grand / Veliki / Nagy / Large (Grande) / 대형 / Didelis / Liels / Groot / Stor / Duza / Grande / Grande / Mare / Большой / Velký / Velik / Veliki / Stor / Büyük / Великий / 大 / 大規模
R	Regular / Нормално / Normální / Normal / Normal / Κανονική / Regular / Standardsuurus / Tavallinen / Standard / Uobičajeni / Reguláris / Regular (Regolare) / 일반 / Normalus / Parasts / Normaal / Vanlig / Standardowa / Regular / Normal / Običnuit / Стандартный / Regulárny / Običajen / Regularni / Normal / Normal / Звичайний / 規則

TRAcelet™
Compression Device

1. Device Description

The compression device consists of a closure band (Figure 1, detail 1) and a syringe (Figure 1, detail 2).


Figure 1. Compression Device Components

The closure band consists of a compression balloon (3), an air injection port (4), a dial (5), a dial-lock (6), and dial markings (7).

2. Indications for Use

The compression device is used to assist patent hemostasis of the radial artery after a transradial procedure.

3. Contraindications

The use of the compression device is contraindicated in the following populations.

- patients with an abnormal Allen's test, abnormal Barbeau test, abnormal radial pulse, or insufficient dual arterial supply
- patients with infection or other serious skin diseases at the access site

4. Warnings

Caution: Read all instructions carefully. Failure to properly follow the instructions, warnings, and precautions may lead to serious consequences or injury to the patient.

- For single-patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a new risk of device contamination, which could result in patient injury, illness, or death. Cleaning, disinfection, and resterilization may compromise essential material and design characteristics of the device and lead to device failure.
- The syringe provided should be used only to inject air into the compression balloon of the closure band. If an introducer sheath is in place, do not inject air into the introducer sheath.
- When in use, do not put excessive load on the compression balloon, which could damage or alter the performance of the device.
- Do not inject more than a total of 18 mL of air. The closure band may be damaged if this volume is exceeded.
- Do not soak or wipe with agents containing organic solvents. Such agents may cause damage to the closure band.
- The patient should not be left unattended while the closure band is in use.

5. Precautions

- This device is intended for use by trained personnel.
- During the procedure, only use the syringe provided to adjust air pressure in the compression balloon.
- When connecting a syringe to the closure band, keep the syringe plunger in place to prevent loss of air compression in the compression balloon. Loss of air compression may cause bleeding.
- Ensure that no foreign particles enter the air valve when injecting air to avoid air leaking.
- If there is itching or redness of the skin during compression, stop using the closure band and treat appropriately.
- Dial markings do not indicate a specific pressure or compression balloon volume and should be used for relative reference only.

- Do not use this product if the unit package or the product has been damaged or soiled.
- Use this product immediately after opening the package.
- Avoid exposure to water, direct sunlight, extreme temperatures, or high humidity during storage.
- After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

6. Adverse Events

Depending on the patient's condition and the degree of compression balloon pressure, potential adverse events include, but are not limited to, the following:

- | | |
|-----------------------------------|---------------------|
| ■ Artery occlusion | ■ Numbness |
| ■ Arterial or venous embolization | ■ Infection |
| ■ Hemorrhage | ■ Allergic reaction |
| ■ Hypodermic hematoma | ■ Pain |

7. Instructions for Use

7.1. Use of the compression device during procedure

1. Ensure that the Δ dial marking lines up with the dial-lock.
2. After the transradial procedure is complete, withdraw the hub of the introducer sheath 2 to 3 cm from the access site.
3. Position the closure band over the wrist such that the access site is visible through the clear window and the center of the dial is aligned with the arteriotomy (Figure 2).

Note: Position the closure band so that the Medtronic logo is on the same side as the patient's little finger.


Figure 2. Positioning the Closure Band.

4. Guide the strap around the wrist and into the strap hook. Pull the strap so that it is stable, and secure the strap using the Velcro® fasteners.
Note: Ensure that there is at least 2 cm of Velcro® overlap when the strap is fastened.
5. Connect the air injection port to the syringe provided. Inflate the compression balloon by injecting approximately 13 mL of air.
Caution: Do not inject more than a total of 18 mL of air into the compression balloon.
6. Remove the introducer sheath according to the instructions for use provided with that device.
7. If bleeding is visible from the access site after removing the sheath, use the syringe provided to inject air into the compression balloon until there is no bleeding.
Caution: Do not inject more than a total of 18 mL of air into the compression balloon.
8. Using the syringe provided, carefully remove air from the compression balloon until bleeding is observed at the access site. Gradually reinject air into the compression balloon until bleeding at the access site stops.
9. Optional step: To confirm patent hemostasis, evaluate the patient's radial distal flow through appropriate means. Use the syringe provided to adjust air in the compression balloon until there is radial distal flow and bleeding from the access site has stopped.

7.2. Adjustment of the compression device during patient recovery

1. To decrease pressure on the access site, depress the dial-lock with one hand and simultaneously rotate the dial counterclockwise with the other hand (Figure 3, detail 1).

To increase pressure on the access site, depress the dial-lock with one hand and simultaneously rotate the dial clockwise with the other hand (Figure 3, detail 2).

Caution: Dial markings do not indicate a specific pressure or compression balloon volume and should be used for relative reference only.

Note: As the numerical value of the dial marking that is aligned with the dial lock decreases, pressure also decreases.

Figure 3. Using the Dial to Adjust Pressure.

1. Decreasing pressure
 2. Increasing pressure
2. If necessary, use the syringe provided or a sterile, standard Luer syringe to adjust air pressure in the compression balloon.
- Caution:** Do not inject more than a total of 18 mL of air into the compression balloon.
- Note:** Air injection volume and compression time may differ according to the patient's condition and medication regimen, volume of heparin administered, and the size of the access site. Adjust the pressure based on an evaluation of the access site.
3. Before removing the closure band, ensure that bleeding at the access site has stopped.

8. Disclaimer of Warranty

ALTHOUGH THE COMPRESSION DEVICE, HEREAFTER REFERRED TO AS "PRODUCT", HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC, INC., AND AFFILIATES (COLLECTIVELY, "MEDTRONIC") HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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
CERTIFICATE OF AUTHENTICITY AND SIGNATURE WITNESS

On this the 17th of July before me, the undersigned Notary Public, personally appeared Colleen Mullins provided to me through satisfactory evidence of identification, which was their employee identification number and/or drivers license number, to be the person whose name is signed on the preceding or attached document in my presence.

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Signature of Document Custodian

Subscribed and sworn to before me this the
17th day of July 2017

Heather P Morose

My commission expires on June 22, 2023

Colleen Mullins – Principal RA Specialist

Medtronic Inc.
710 Medtronic Parkway
Minneapolis, MN 55432, US.
Tel: +1 763 514 4000

Subscribed and sworn to me on this date – July 17th, 2017

Heather Morose – Notary Public

Appendix to the Instruction for use (Version # M970663A001)

TRAcelet Compression Device

1.1. With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity, location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;
Design Facility: Medtronic Vascular, 3576 Unocal Place, Santa Rosa, CA 95403, USA. Tel: 707-566-1143 Legal Manufacturer: Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432, USA www.medtronic.com, Tel: 763-514-4000 Manufacturing Facility: Medplast Medical Costa Rica, S.A., Parque Zona Franca Metropolitana, Edificio 2c Barreal de Heredia, 40101 Costa Rica
1.2. With respect to authorized representative of the manufacturer of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity

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Authorized Representative of Manufacturer:
“Medtronic” LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10
Tel: (495) 580-73-77, Fax: (495) 580-73-78
www.medtronic.ru

Holder of Registration Certificate in Russia
“Medtronic” LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10
Tel: (495) 580-73-77, Fax: (495) 580-73-78
www.medtronic.ru

2. Technical characteristics of the medical device

PRODUCT CHARACTERISTIC	SPECIFICATION
Effective Strap Length	Regular: 19.2 +1.0 / -0.5 cm Large: 25.2 +1.0 / -0.5 cm
Syringe to Anti-Lock Cap Tensile	<u>Tensile Force</u> 8.4 lbf min (37,37 H)
Syringe Luer 6% Conical Fitting	Conical Luer fitting on syringe shall meet requirements per Table A-1.
Check Valve Luer 6% Conical Fitting	Conical Luer fitting on check valve shall meet requirements per Table A-1.
Balloon Rupture (full assembly, post-flexural conditioning)	<u>Pressure</u> 21 psi minimum (144,8 kPa)
Side-Tube to Balloon Tensile	<u>Tensile Force</u> 5.1 lbf min (22,69 H)
Side-Tube to Check Valve Tensile	<u>Tensile Force</u> 5.1 lbf min (22,69 H)
Velcro to Layflat Tubing Weld Tensile	<u>Tensile Force</u> 4.1 lbf min (18,24 H)
Rivet Joint Tensile	<u>Tensile Force</u> 10.5 lbf min (46,71 H)
Dial Removal Torque	<u>Torque (to remove dial from frame with lock disengaged)</u> 7.5 in-lbs. min (0,85 H-m)
Dial to Threaded Window Shear Strength	<u>Torque</u> 40 in-lbs. min (4,52 H-m)
Balloon Contact Pressure Reduction with Dial	<u>Pressure reduction from nominal inflation</u> 60% min
Initial Balloon Inflation Pressure	<u>Contact Pressure at 13cc (Nominal Volume)</u> 3.15 ± 1.00 psi (21,72 kPa±6,89 kPa)
	<u>Contact Pressure at 18cc (Maximum Volume)</u> 4.50 ± 1.00 psi (31,03 kPa ±6,89 kPa)
Balloon Internal Pressure Over Time (full assembly)	<u>Pressure Decay over 4 Hours</u> 30% max
Dial Torque	<u>Torque (with lock disengaged)</u> 5 in-lbs. max (0,56 H-m)
Dial Lock Disengagement Force	<u>Force</u> 2.6 lbf max (11,57 H)


Available on request

SPECIFICATION
Device visibility shall aid placement and leakage assessment.
Device shall remain stable during inflation and dial rotation.
Pressure can be adjusted until leakage is stopped.
Pressure can be adjusted until leakage is stopped and distal flow Is observed.
The pulse in the ulnar artery can be assessed while the device is in place.
The dial lock should help to discourage patient tampering.

Syringe Characteristics:

Weight of syringe	13.6 GRAMS±0.68 GRAMS
Rated syringe capacity	20mL±1 mL
Dead space volume	Approximately .1mL
Syringe barrel diameter	Spec Sheet Below
Syringe barrel length	Spec Sheet Below
Total scale length	Spec Sheet Below
Division value	Spec Sheet Below
Force required to drive the plunger rod	No specification for this
Plunger rod length	Spec Sheet Below
Dimensions of plunger rod stop	Spec Sheet Below
Syringe tip location	Center, Luer Lock
Syringe tip port diameter.	Per ISO standard

MONOJECT 20mL Luer Lock Syringe, Nonsterile

1. Printed Barrel 2. Plunger 3. Plunger Tip 4. Lubrication - Silicone

Syringe Assembly:

- 20 mL Clear Molded Printed Barrel
- **Dimensions:** Overall Length: 3.88 Inches (9,855 sm±0.5 sm) Diameter: 0.850 Inches (2.159 sm±0.1 sm);
- **Barrel Printing:** Black Ink

Available on request

- Syringe Plunger, 20mL±1 mL, 4 Lobe
- **Material:** Radiation Grade Polypropylene PF535 for use in manufacture of medical devices

Dimensions: Overall Length: 3.948 Inches (10,03 sm±0.05 sm)

3.. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)

In accordance with IFU avoid exposure to water, direct sunlight, extreme temperatures, or high humidity during storage. Use at room temperature (18 ° C) not higher than 25 ° C, relative humidity of 75%, a pressure of from 86.0 to 106.0 kPa

Avoid exposure to water, direct sunlight, extreme temperatures, or high humidity during storage. Prior to use, the product should be stored in a cool, dry place:

- Temperature from + 20°C to +25°C;
- relative humidity of 20% - 80%;
- pressure from 100 kPa to 103 kPa.

The product should not be exposed to organic solvents, ionizing radiation, or ultraviolet light. Carefully inspect the sterile package for damage before opening. Do not use if the package has been damaged or opened.

The transport conditions are the same as storage condition.:

- Temperature from + 20°C to +25°C;
- relative humidity of 20% - 80%;

pressure from 100 kPa to 103 kPa.

4. List of national standards applied by the manufacturer

Document Number	Revision	Document Title
EN ISO 13485	2012/AC: 2012	Medical devices - Quality management systems - Requirements for regulatory purposes
ISO 594-1	1986(E)	Conical Fittings with a 6% (Luer) Taper for Syringes, Needles, and Certain other Medical Equipment – Part 1: General Requirements
ISO 14644-1	2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration

Available on request

ISO 14644-2	2015	Cleanrooms and associated controlled environments -- Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
EN ISO 14971	2012	Medical devices - Application of risk management to medical devices
EN 1041	2008/A1:2012	Information supplied by the manufacturer of medical devices
ISO 15223-1	2012	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 11607-1	2006/Amd 1: 2014	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
ISO 11607-2	2006/Amd 1: 2014	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
EN ISO 10993-1	2009/Cor 1:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5	2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
EN ISO 10993-7	2008/Cor 1: 2009	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals
ISO 10993-10	2013	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
EN ISO 10993-11	(BS EN) 2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
EN ISO 10993-12	2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials
ISO 11135	2014	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices
EN 556-1	2001/AC:2006	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices

Available on request

MEDDEV 2.7.1	Rev. 3 (2009)	Clinical Evaluation: A guide for manufacturers and notified bodies
MDD 93/42/EEC	Amended by 2007/47/EC	Medical Device Directive

5. Warnings and/or precautions taken by the consumer when disposing of the medical device, accessories and consumables used together with it (if any), including information about the infectious, microbial, environmental or physical hazard of the medical device, Procedure and conditions for the disposal or destruction of a medical device.

After use, this product may be a potential biohazard (class B). Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

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© 2016 Medtronic