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INT SYS 200
Installation Instructions
60722119 Rev. C

Introduction

These Instructions provide information for the installation of the MEDRAD® Intego PET Infusion System only and do not provide operation or safety related use information. Read all Warnings and Cautions thoroughly as provided in the MEDRAD® Intego Operation Manual before installing this device. Installation shall be performed by a qualified Bayer HealthCare Services representative or authorized dealer.

UNITED STATES Bayer Medical Care Inc. 1 Bayo Drive Millsboro, PA 19061-0780 U.S.A. Phone: +1 812 267 2400 +1 800 633 2231 Fax: +1 812 267 4220	EUROPE Bayer Medical Care B.V. Horsterweg 24 6159 AC Maastricht Airport The Netherlands Phone: +31 (0) 43 365 5101 Fax: +31 (0) 43 365 5590	JAPAN エーザイ 薬品株式会社 〒510-0001 大阪市北区梅田2-4-9 日本 電話: +81 (0) 6-6333-6250 FAX: +81 (0) 6-6344-2395
--	--	---

Definitions

The following are definitions of the terms WARNING, CAUTION and NOTE found throughout this document:

WARNING: Indicates that the information is a warning. Warnings advise of circumstances that could result in injury or death to the patient or operator. Read and understand the warnings before operating the system.

CAUTION: Indicates that the information is a caution. Cautions advise of circumstances that could result in damage to the device or improper functioning of the device. Read and understand the cautions before operating the system.

NOTE: Indicates that the information that follows is additional important information or a tip.

Warnings, cautions and notes are located throughout this document, where applicable, and in the MEDRAD® Intego Operation Manual.

Installing the System

The following is a list of tools required for the uncrating and set-up of the system:

- Utility Knife
- 10 mm & 19 mm sockets and ratchet
- 4 mm & 5 mm Allen wrench
- Phillips Screwdriver
- Pliers or Vice Grips
- Electrical tape
- Cesium calibration source to be supplied and used by hospital
- Cobalt calibration source to be supplied and used by hospital
- FDG Linearity Check source to be supplied and used by hospital

Uncrating the System

NOTE: Uncrate the system at the loading dock.

1. Cut the straps holding the cardboard sleeve with a utility knife.
2. Remove the ramp from the top of the cardboard sleeve and set it aside.
3. Remove the top foam and the boxes (accessories box and ionization chamber) that are underneath the foam on top of the system. Remove both boxes and set aside.
4. Remove the cardboard sleeve from around the system.

5. Place the ramp next to the bottom of crate as shown.

6. On the floor of the crate next to the system is a covered compartment holding the ramp screws. Remove this cover and remove the two screws from inside the compartment. Install these screws into each side of the ramp where it overlaps the bottom of the crate.

NOTE: The screws only need to be installed deep enough to prevent the ramp from moving.


7. Remove the rear cover to expose the brackets that secure the system to the shipping pallet.
 - a. Unscrew the manual override knob.
 - b. Remove the screws that secure the cable storage posts to the system using a 5 mm Allen wrench.

8. Remove the printer.
 - a. Disconnect the printer cables.
 - b. Remove the screws at the bottom of the mounting plate.
 - c. Lift the printer off of the top mounting flange.

9. Remove panel on the right side of the system.
 - a. Remove the two screws that secure the side panel to the rear of the system.

- b. Remove the two screws that secure the side panel to the front of the system.


10. Remove the mounting hardware from the front and rear of the system.
 - a. Using a 19 mm socket and ratchet, remove the two mounting bolts from the front of the system.

- b. Using a 19 mm socket and ratchet, remove the mounting bolts that secure the rear mounting brackets to the pallet at the rear of the system.
- c. Using a 10 mm socket and ratchet, remove the bolts that secure the mounting brackets to the rear of the system.

11. Open the Service Tray and setup the batteries.
 - a. Remove the 2 bottom screws from the bracket on the service tray using a 4 mm Allen wrench and pull it open.

- b. Install the fuses in the UPS Battery. The two fuses are in a plastic bag attached to the battery.

WARNING: Fire hazard: Patient injury could result from using incorrect fuses. To avoid an electrical fire, ensure that the correct type of fuse is used for replacement. The fuse must be replaced with a Type F 250 V, 5 A fuse by qualified personnel only.

- c. Connect the RED Motor Control Battery lead wire into the circuit breaker.

- d. Leaving the service tray open at this time will help in the installation of the ion chamber.

CAUTION: Component damage may occur if not installed properly. Ensure all connections are secure and do not get pinched when closing the service tray; DO NOT over tighten.

Quick SAS Installation Guide

Pressing the red stop button will stop all fluid movements during any portion of the procedure.
E.g. When infusing a patient, pressing the red stop button will pause the infusion. It will not discontinue the procedure.

WARNING: Air embolism hazard. Injury or death can result from infusing air.
Ensure air is expelled from the disposable sets before beginning infusion.

✓ Verify tubing is lying flat prior to closing shielded vault lid.


- 3. Insert needle cartridge clip properly.**
Tube side of clip should face toward black RP pump.


- 2. Install tubing through saline pump.**
Center tubing on rollers. Place luer against pump. Center tubing in door.

- 1. Fully spike saline bag.**

After installing RP vial:

- 10. Place collar properly.**
Ensure that the collar rests immediately outside the pump, opposite the hinge.


- 4. Place confluence into the holder.**

- 5. Verify that the dose calibrator is empty before installing the coil.**


- 6. Insert waste container into the shielded waste storage.**

- 7. Place tubing T-Connector into the holder.**
Note the orientation of the triangle in one corner.

- 8. Insert tubing into pinch valves.**
Depress pinch valve button, insert tubing fully, release button. Ensure tubing is straight. Install pinch valve covers onto pinch valves, if available.

- 9. Fully seat tubing into air detector.**

medrad® Intego
PET Infusion System

CE 0086 April 8, 2015 60765535 Rev. C

60755640 Rev G © 2011 - 2017 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.

Български
繁體中文
Hrvatski
Česky
Dansk
Nederlands
English
Eesti
Suomi
Français
Deutsch
Ελληνικά
Magyar
Italiano
한국어
Latvija
Lietuvių
Norsk
Polski
Portuguê (Brazil)
Portuguê (European)
Română
Русский
Srpski
Slovenský
Slovenščina
Español (EU)
Svenska
Türkçe
Türk

MEDRAD® Intego

If you require a printed copy, contact Customer Support.
For a full list of Bayer offices and authorized dealers,
go to <http://www.radiology.bayer.com/contact>

Bayer, the Bayer Cross, MEDRAD, Intego and MEDRAD
Intego are trademarks owned by and/or
registered to Bayer in the U.S. and/or other countries.
Acrobat Reader is needed to view the
information on this CD. www.adobe.com

Operation Manuals & SAS Quick Guides

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[illegible][illegible]

© 2014, 2015, 2017 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer. Bayer and the Bayer Cross are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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