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
CE 0483

B·R·A·H·M·S PIGF plus KRYPTOR Calibrator

IVD for professional use only	B·R·A·H·M·S is a registered trademark of B·R·A·H·M·S GmbH. Lumi4®-Tb is a registered trademark of Lumiphore, Inc., licensed for use by B·R·A·H·M·S.
TRACE™	Other product names in this document are used for identification purposes; they may be trademarks and/or registered trademarks of their respective companies.
	The manufacture and/or use of this product is covered by one or more of the following patents: AU2003265294, AU2009202176, CA2496253, EP1575416, EP2172220, EP2305301, EP2308507, HK1142535, HK1085368, HK1156521, HK1156233, IL166393, IN218/CHENP/2005, MXPA05000832
Date	
20.02.2020	Initial version

Intended Use

B·R·A·H·M·S PIGF plus KRYPTOR Calibrator to readjust the standard curve during the automated immunofluorescent assay for the quantitative determination of the concentration of human placental growth factor (PIGF) in human serum and EDTA plasma on the analyses of B·R·A·H·M·S KRYPTOR series.

Contraindications

none

Side effects

none

Instruments

B·R·A·H·M·S PIGF plus KRYPTOR may be used on:

- ☐ B·R·A·H·M·S KRYPTOR
- ☐ B·R·A·H·M·S KRYPTOR compact
- ☒ B·R·A·H·M·S KRYPTOR compact PLUS
- ☒ B·R·A·H·M·S KRYPTOR GOLD

Availability of instrument depends on registration in the country.

B·R·A·H·M·S PIGF plus KRYPTOR Calibrator CAL			
Is used together with B·R·A·H·M·S PIGF plus KRYPTOR kit (859.075)			
REF 85991		  see label for expiry date	
Name	Quantity	Quality	Description
B·R·A·H·M·S PIGF plus KRYPTOR CAL	VIALS 6	LYOPH	Defibrinated normal human plasma heat treated, recombinant human PIGF, kathon
Bar code card	1	ready for use	see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual. The bar code card contains information related to the calibrator lot including its concentration.

Other Requisites (Not supplied with the kit)	
Name	REF
B·R·A·H·M·S PIGF plus KRYPTOR	859.075
B·R·A·H·M·S KRYPTOR compact PLUS Consumables	
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89981
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89982
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89983
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89984
B·R·A·H·M·S KRYPTOR compact DILCUP	89985
B·R·A·H·M·S KRYPTOR compact REACT	89986
B·R·A·H·M·S KRYPTOR GOLD Consumables	
Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89991
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89992
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89993
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89994
B·R·A·H·M·S KRYPTOR compact DILCUP	89995

Precautions

	Follow the information given in the user manual of the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD and in this document.
	– Check the expiration dates specified on the label. – Reagents from different reagent kits should not be mixed. – Do not separate the elements of a unit. – Avoid any microbial contamination of the reagents.
	Raw materials of animal and/or human origin contained in the reagents have been tested and/or processed through virus inactivation procedures as far as possible. However, as it is impossible to strictly warrant that such products will not transmit any infectious agents, all raw materials of human or animal origin including the samples to be assayed must be treated as potentially infectious.
 	The generally acknowledged safety precautions and laboratory techniques must be observed when handling reagents and patient samples. – Do not pipet by mouth. – Wash hands after work.
 	– Wear protective clothing, protective gloves and safety glasses at work.
 	– Do not eat, drink or smoke in areas where samples or kit-reagents are handled.
	– Remove spills with absorbing paper. – All the material used for cleaning up must be disposed of as infectious laboratory waste. – Prevent from getting into sewage, water, ground.
	– Used reagent plates and reagent kits are disposed of as potential infectious laboratory waste according to local regulations. – Empty containers should be returned to local recyclers.

Calibration	CAL
• Reconstitute each vial with the volume distilled water (conductivity of less than 50 $\mu\text{S}/\text{cm}$ is recommended) indicated on the vial label. • Allow 15 min for the complete dissolution of the lyophilisate. • Mix gently after reconstitution. • Calibration must be carried out with every new reagent kit lot, it is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD in order to readjust the standard curve. • Calibration is automatically managed by instrument. • An aliquot of 70μL of calibrator is pipetted automatically by the instrument. • Use the calibrator/calibrators only once. • Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5 hours. • The calibrator bar code card must be read in for each new lot of calibrator. • For further information see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.	

Quality Control
Good laboratory practice requires that control samples are measured regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that the results be analysed using appropriate statistical methods. If desired, B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD can automatically check the quality of assays at intervals, by statistical analysis on the basis of Levey Jennings graphs. It is necessary to comply with national quality assurance guidelines for quantitative tests in the medical laboratory (current version). For instance, test accuracy and precision should be monitored by means of laboratory in house and/or commercially available control materials. If unacceptable control values are obtained, proceed as outlined in standard laboratory diagnostic procedures to determine the cause and implement corrective measures.

Storage conditions
<u>Storage (stability)</u> Product stability tests are carried out in accordance with the requirements of EN 23640: 2002 (Tests of the stability of reagents for in vitro diagnostics), as well as according to Guidelines of the Institute of Clinical and Laboratory Standards EP25-A (Evaluation of the stability of reagents for in vitro diagnostics).
<u>Storage (stability) after opening</u> The calibrator is stable for up to 5 hours on board B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD analyzer at room temperature.
Transportation conditions

The product should be transported in the original manufacturer's packaging, protecting products from external influences, by all kinds of transport in covered vehicles in accordance with the transportation rules applicable to any specific type of transport.

Heat impact test confirms stability at a temperature of 2 to 8°C for the entire shelf life and for 5 hours at a room temperature from 18 to 25°C on board B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD analyzer.

The calibrator in the packaging of the manufacturer should be stored at a temperature of 2 to 8°C for the entire shelf life in cold chambers or refrigerators, providing a regulated temperature regime with daily temperature recording).

During cargo handling and transportation, it is necessary to strictly comply with the requirements of handling signs applied on the package and the transport containers.

If the transportation and storage conditions are not followed, the medical device must not be used!

Shelf life

Maximum 24 months

Disposal

The medical device must be disposed according to local, regional and national laws. All samples and materials used and generated during research, including expired reagents, should be disposed according to SanPiN 2.1.7.2790-10 as potential infected medical waste (Class B).

Product	- Products must be disposed of as laboratory chemical waste in accordance with the local regulations. - Please contact the responsible authority. - Used reagent plates and reagent kits dispose of as potential infectious laboratory waste.
Contaminated packaging	- Dispose of according to local regulations - Recommendations: Remove the residues. Return the empty containers to local waste recycling stations.
Pollutes	- Remove pollutes with absorbing paper. - All materials used for cleaning up, must be disposed of as infectious laboratory waste.

When disposing and destroying medical devices that have expired, as well as reagent waste, it is necessary to act in accordance with hygienic standards and follow official instructions for the disposal and destruction of potentially infectious and chemical material adopted in the Russian Federation and in the given health facility.

The key issues in waste management strategies are optimization of tactics and approaches to final

products, as well as minimization of possible environmental pollution. Given the nature of the waste reagents that contain human biomaterial, thermal disposal and destruction is the most optimal at a temperature of 1000 C. The waste should be stored in an airtight container marked "Biological Hazard."

Employees should undergo professional training in occupational safety and hygiene.

Thermal disposal and destruction at 1,000 C is also allowed for expired medical device.

When using the kit, waste of classes A, B and G are generated, which are classified and destroyed (disposed of) in accordance with SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for medical waste management" and MU 287-113 "Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical products".

In order to avoid possible environmental pollution, expired reagents as well as waste from reagents used, must be disposed of in accordance with the local government and/or hospital regulations.

Proper destruction and disposal will support the conservation of natural resources and help protect human health and the environment.

Manufacturer's authorized representative

Thermo Fischer Scientific Joint Stock Company (Thermo Fischer Scientific JSC)

73, Kubinskaya St., Letter A, Bldg 1, St. Petersburg 196240

Tel.: (812) 703-42-15, facsimile: (812)703-42-16

Email: info.lcp.spb@thermofisher.com

Warranties

The Manufacturer guarantees the conformity of the characteristics of the product stated in the operating documentation, subject to use in accordance with the instructions and intended use, provided by the manufacturer.

For all questions, technical advice and support please contact the **manufacturer's authorized representative in the Russian Federation:**

Thermo Fischer Scientific Joint Stock Company (Thermo Fischer Scientific JSC), Russia

73, Kubinskaya St., Letter A, Bldg 1, St. Petersburg 196240

Tel.: (812) 703-42-15, facsimile: (812)703-42-16,

Email: info.lcp.spb@thermofisher.com

The list of applicable standards

1. EN ISO 15223-1:2016

Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied- Part 1: General requirements (ISO 15223-1:2016)

2. EN ISO 13485:2016

Medical devices - Quality management systems - Requirements for regulatory purposes

3. EN 13612:2002

Performance evaluation of in vitro diagnostic medical devices
EN 13612:2002/AC:2002

4. EN 13641:2002

Elimination or reduction of risk of infection related to in vitro diagnostic reagents

5. EN ISO 14971:2012

Medical devices - Application of risk management to medical devices
(ISO 14971:2007, Corrected version 2007-10-01)

6. EN ISO 17511:2003

In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials
(ISO 17511:2003)

7. EN ISO 18113-1:2011

In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements
(ISO 18113-2:2009)

8. EN ISO 18113-2:2011

In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use
(ISO 18113-2:2009)

9. EN ISO 23640:2015

Evaluation of stability of in vitro diagnostic reagents
(ISO 23640:2011)

10. IEC 62366:2015 / EN 62366-1:2015

Medical devices - Part 1: Application of usability engineering to medical devices (IEC 62366-1:2015 + COR1:2016);

■ Symbols

Symbols used in Instruction for Use and Product Labelling of B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
	Reference to the Intended use of the Medical Device		In Vitro Diagnostic Medical Device		Batch Code
	Contents		Calibrator		Control
	Buffer		B-R-A-H-M-S KRYPTOR SOLUTION 1/B-R-A-H-M-S KRYPTOR compact SOLUTION 1/ B-R-A-H-M-S Solution 1 KRYPTOR GOLD		B-R-A-H-M-S KRYPTOR SOLUTION 2/B-R-A-H-M-S KRYPTOR compact SOLUTION 2/ B-R-A-H-M-S Solution 2 KRYPTOR GOLD
	B-R-A-H-M-S KRYPTOR SOLUTION 3/B-R-A-H-M-S KRYPTOR compact SOLUTION 3/ B-R-A-H-M-S Solution 3 KRYPTOR GOLD		B-R-A-H-M-S KRYPTOR SOLUTION 4/B-R-A-H-M-S KRYPTOR compact SOLUTION 4/ B-R-A-H-M-S Solution 4 KRYPTOR GOLD		Bags contained
	Bags		Plates contained		Plates
	Vials contained		Vials		Val
	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL		Lyophilized, Freeze Dried		Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
	Registered Trade Mark		Trade Mark		Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		GHS05, Corrosion
	GHS07, Exclamation Mark		Trade Mark for TRACE™-technology		CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg. Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

The notary asked before if he was or is active in the matter which is object of this document outside of his notarial activity. The question was negated.

Hereby I certify that both members of the management team of the company B.R.A.H.M.S GmbH

1. Dr. Bernhard **C i o m m e r**
born January 28, 1955
domiciled in D-13593 Berlin
Folkungerstr. 21,
2. Dr. Jörg-Markus **P a u l**
born August 02, 1967,
domiciled in D-13503 Berlin
Reimerswalder Steig 22,

both commercially domiciled in 16761 Hennigsdorf, Neuendorfstrasse 25,
acknowledged their signatures before me.

The identity of the deponents is personally known to me.

At the same time I certify according to § 21 BNotO due to the presented certified extract from the statement of the Register of Companies of the District Court of Neuruppin dated March 09, 2020 that with HRB-No. 8872 NP B.R.A.H.M.S GmbH domiciled in Hennigsdorf is registered and that Dr. Bernhard Ciommer as well as Dr. Jörg-Markus Paul are registered as management team members and are both authorised to represent the company.

Berlin, in the business premises of the notary Lutz Baxmeier,
Berliner Strasse 32, D-13467 Berlin

March 10, 2020

Scroll no. 237/2020

Baxmeier
Notary in Berlin
Federal Republic of Germany

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