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
Instruction for Use

B·R·A·H·M·S sFlt-1 KRYPTOR

IVD for professional use only	B·R·A·H·M·S is a registered trademark of B·R·A·H·M·S GmbH. Lumi4®-Tb is a registered trademark of Lumiphore, Inc., licensed for use by B·R·A·H·M·S.
TRACE™	Other product names in this document are used for identification purposes; they may be trademarks and/or registered trademarks of their respective companies.
	The manufacture and/or use of this product is covered by one or more of the following patents: AU2003265294, AU2009202176, CA2496253, EP1575416, EP2172220, EP2305301, EP2308507, HK1142535, HK1085368, HK1156521, HK1156233, IL166393, IN218/CHENP/2005, MXPA05000832
Date	
20.02.2020	Initial release

Intended Use

B·R·A·H·M·S sFlt-1 KRYPTOR is an automated immunofluorescent assay for the quantitative determination of the concentration of sFlt-1 (soluble fms-like tyrosine kinase 1, also known as VEGF receptor-1) in human serum and EDTA plasma on the analyses of B·R·A·H·M·S KRYPTOR series.

Contraindications

none

Side effects

none

Instruments

B·R·A·H·M·S sFlt-1 KRYPTOR may be used on:

☐ B·R·A·H·M·S KRYPTOR

- ☐ B·R·A·H·M·S KRYPTOR compact
☒ B·R·A·H·M·S KRYPTOR compact PLUS
☒ B·R·A·H·M·S KRYPTOR GOLD

Availability of instrument depends on registration in the country.

Introduction

The B·R·A·H·M·S sFlt-1 KRYPTOR assay is intended to be used – together with the B·R·A·H·M·S PIGF plus KRYPTOR assay and additional clinical and diagnostic findings - as an aid in the diagnosis of pre-eclampsia and/or short-term prognosis of pre-eclampsia in pregnant women with suspicion of pre-eclampsia.

Pre-eclampsia, which affects about 2-8 % of pregnancies, is a major cause of perinatal and maternal morbidity and mortality. First clinical symptoms of pre-eclampsia such as new onset of hypertension and proteinuria are observed after 20 weeks of gestation. Clinically, pre-eclampsia may vary from mild to severe forms, requiring delivery before 34 weeks of gestation. The severe form of pre-eclampsia, HELLP syndrome (Hemolysis, elevated liver enzymes, low platelets), occurs in about 20 % of women affected by preeclampsia.[1] [2]

Although the cause of pre-eclampsia remains unclear, the syndrome may be initiated by an imbalance of placental factors that induce endothelial dysfunction. The antiangiogenic factor sFlt-1 (soluble fms-like tyrosine kinase) acts as potent antagonist of proangiogenic factors such as Placental Growth Factor (PlGF) by adhering to the receptor-binding domains, preventing the interaction with endothelial receptors and thereby inducing endothelial dysfunction.[3]

During pregnancy, PlGF levels increase progressively in first and second trimester and decrease towards term. In women with clinical pre-eclampsia, sFlt-1 levels are significantly increased while concentrations of circulating free PlGF are significantly decreased. Significantly reduced levels of free PlGF have been shown in pregnant women who subsequently develop preeclampsia.[4] [5]

Measurement of sFlt-1 levels together with PlGF levels in maternal serum starting in mid pregnancy can improve the current diagnostic workup of patients with suspected pre-eclampsia, which includes clinical symptoms, proteinuria and Uterine artery Doppler velocimetry. The sFlt-1/PlGF ratio is a tool for the diagnosis of pre-eclampsia by discriminating between preeclamptic pregnancies and other forms of pregnancy-related hypertensive disorders. In addition, the sFlt-1/PlGF ratio may be a valuable prognostic tool of subsequent maternal and fetal adverse outcome in women diagnosed with pre-eclampsia.[6][7][8][9]



Contents

Kit

B·R·A·H·M·S sFlt-1 KRYPTOR

CONT

REF 845.075

see label for expiry date

Name	Quantity	Quality	Description
K- conjugate	VIAL 1	ready for use	1 Anti-sFlt-1 monoclonal antibody conjugated with Lumi4®-Tb, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, dextran, kathon

XL-conjugate	VIAL 1	ready for use	1 Anti-sFit-1 monoclonal antibody conjugate with Cyanin5.5, buffer, bovine albumin, bovine immunoglobulins, mouse immunoglobulins, dextran, kathon
Bar code card	1	ready for use	see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD user manual. The bar code card contains all information needed for the registration of a new reagent lot.

Accessories	
Not supplied with the product	
Name	REF
B·R·A·H·M·S sFit-1 KRYPTOR CAL	84591
B·R·A·H·M·S sFit1- KRYPTOR QC	84592

Other Requisites	
Not supplied with the product	
B·R·A·H·M·S KRYPTOR compact PLUS Consumables	
Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89981
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89982
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89983
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89984
B·R·A·H·M·S KRYPTOR compact DILCUP	89985
B·R·A·H·M·S KRYPTOR compact REACT	89986
B·R·A·H·M·S KRYPTOR GOLD Consumables	
Name	REF
B·R·A·H·M·S KRYPTOR BUFFER	89970
B·R·A·H·M·S KRYPTOR compact SOLUTION 1	89991
B·R·A·H·M·S KRYPTOR compact SOLUTION 2	89992
B·R·A·H·M·S KRYPTOR compact SOLUTION 3	89993
B·R·A·H·M·S KRYPTOR compact SOLUTION 4	89994

B·R·A·H·M·S KRYPTOR compact DILCUP	89995
B·R·A·H·M·S KRYPTOR compact REACT	89996

Precautions	
	Follow the information given in the user manual of the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD/KRYPTOR GOLD and in this document.
	– Check the expiration dates specified on the label. – Reagents from different reagent kits should not be mixed. – Do not separate the elements of a unit. – Avoid any microbial contamination of the reagents.
	Raw materials of animal and/or human origin contained in the reagents have been tested and/or processed through virus inactivation procedures as far as possible. However, as it is impossible to strictly warrant that such products will not transmit any infectious agents, all raw materials of human or animal origin including the samples to be assayed must be treated as potentially infectious.
	The generally acknowledged safety precautions and laboratory techniques must be observed when handling reagents and patient samples. – Do not pipet by mouth. – Wash hands after work.
	– Wear protective clothing, protective gloves and safety glasses at work.
	– Do not eat, drink or smoke in areas where samples or kit-reagents are handled.
	– Remove spills with absorbing paper. – All the material used for cleaning up must be disposed of as infectious laboratory waste. – Prevent from getting into sewage, water, ground.
	– Used reagent plates and reagent kits are disposed of as potential infectious laboratory waste according to local regulations. – Empty containers should be returned to local recyclers.

Principle

The measurement principle of B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD is based on immunofluorescence and in particular on TRACE™ Technology (Time-Resolved Amplified Cryptate Emission), which measures the signal that is emitted from an immunocomplex with time delay. The basis of the TRACE™ Technology is non-radiative energy transfer from a donor (Cryptate) to an acceptor (XL).

The proximity of donor (cryptate) and acceptor (XL) when they are part of an immunocomplex and the spectral overlap between donor emission and acceptor absorption spectra on the one hand, intensify the fluorescent signal of the donor and on the other hand they extend the life span of the acceptor signal, permitting the measurement of temporally delayed fluorescence.

The assay is using two monoclonal antibodies which recognize sFlt1 antigen in the sample. One monoclonal antibody is labeled with a Cryptate and the second monoclonal antibody is labeled with XL. Cryptate conjugate, XL conjugate and patient sample are pipetted and dispensed in the reaction well. The two antibodies will bind the sFlt1 antigen present in the patient sample and the reaction well will be excited with the nitrogen laser at 337nm

Precise measuring of analyte concentration: When the sample is excited with a nitrogen laser at 337 nm, the donor (cryptate) emits a long-life fluorescent signal in the milli-second range at 620 nm, while the acceptor (XL) generates a short-life signal in the nanosecond-range at 707 nm. When the two components are bound in an immunocomplex, both the signal amplification and the prolongation of the life span of the acceptor signal occur at 707 nm, so that it can be measured over μ -seconds. This long-life signal is proportional to the concentration of the analyte to be measured.

Reliable prevention of interference: Non-specific signals, e.g. the signals of the short-life and unbound acceptor XL and the medium-specific interference signals conditional upon the natural fluorescence of the sample, are eliminated by temporal delay of the fluorescence measurement. The signal generated by the cryptate at 620 nm serves as an internal reference and is measured simultaneously with the long-life acceptor signal at 707 nm which is the specific signal. Interfering influences, e.g. from turbid sera, are automatically corrected by means of the internally calculated ratio of the intensities at these wavelengths.

Instructions

Sample volume:	8 μ L
Incubation time:	9 min
Results are given in:	pg/mL
Conversion factor:	not applicable
Direct measuring range	29...90 000 pg/mL
Measuring range with automatic dilution:	not applicable
Sample type	Serum, EDTA plasma (K2 or K3 EDTA)
Kit stability on board:	29 days
Calibrator:	2 points
Calibration stability:	15 days
Assay principle:	sandwich

- If the assay is carried out within 72 hours of blood collection, store samples at 2...8 °C. If not, samples shall be aliquoted and stored frozen (-20 °C).
- Citrated and heparin plasma should not be used.
- It is possible to use both K2 and K3 EDTA plasma tubes for plasma samples
- Place the sample in a tube suited for use on B·R·A·H·M·S KRYPTOR compact PLUS/ KRYPTOR GOLD (11–17 mm diameter). This might be the primary tube.
- The sample tube must contain a minimum volume, which will vary depending on the diameter of the sample tube. A 13 mm diameter tube will require an additional 150 µL of sample.
- sFlt-1 concentrations in patient samples are not expected to exceed the direct measuring range. Samples cannot be diluted.
- Icteric, hemolytic or hyperlipemic samples, or samples which are turbid or contain fibrin may yield imprecise results. Such samples are signaled by B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.

Conditions of preliminary sample processing

- Serum and EDTA (K2 or K3) plasma can be used
- When coagulation of whole blood is complete (after 30 min at room temperature), the sample should be centrifuged for at least 10 min at a minimum speed of 1500 g.
- The sample should be refrigerated (2°-8°C) not later than 2 hours after centrifugation (3 hours from blood collection) and up to 72 hours. If testing is not done within 72 hours, the serum has to be frozen at -20°C.
- Repeated freezing/thawing steps are not recommended


The operation and maintenance of B·R·A·H·M·S KRYPTOR compact PLUS/ KRYPTOR GOLD are described in the related user manual.

After it has been opened the reagent unit may be stored on the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD in the space provided. Each reagent unit is individually identified (bar code) and its maximum period of use after opening is controlled by the instrument.

To prepare a reagent unit, proceed as follows:

- Remove the security band from the reagent pack.
- Push in the lid by pressing it firmly (see diagram below).

Make sure to remove all the aluminium foil.

Opening the Kit

Note: Always handle the reagent unit carefully to avoid foaming or bubbles. Foam/bubbles may interfere with proper detection of reagent and/or correct dispensing of reagents. Allow the kit to stand 5...10 min within B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD prior to use.

A standard curve does not need to be established on B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.

The instrument stores the required information after reading the bar code from the reagent card. A recalibration must be carried out for every new reagent lot, and subsequently on a regular basis. B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD automatically indicates when a calibration is required. Both, the stored and recalibrated standard curve may be displayed on the screen.

The following steps are carried out:

- Conjugates and sample are dispensed into the reaction plate and the fluorescent signal emitted is measured periodically.
- After measurement of the fluorescent signal, the data obtained from the software are compared to the memorized standard curve.

Calibration	CAL
• Reconstitute each vial with the volume distilled water (conductivity of less than 50 $\mu\text{S}/\text{cm}$ is recommended) indicated on the vial label. • Allow 15 min for the complete dissolution of the lyophilisate. • Mix gently after reconstitution. • Calibration must be carried out with every new reagent kit lot, it is then repeated on a regular basis automatically managed by B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD in order to readjust the standard curve. • Calibration is automatically managed by instrument. • An aliquot of 8 μL of calibration is pipetted automatically pipetted by instrument. • Use the calibrator/calibrators only once. • Do not leave the calibrator/calibrators at room temperature or on the carousel for more than 5	

hours.

- The calibrator bar code card must be read in for each new lot of calibrator.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.

Control

CONTROL

- It is recommended that controls be run once a day but at least after each calibration.
- A control tube is directly processed like a sample tube.
- Reconstitute each vial with the volume of distilled water (conductivity of less than 50 $\mu\text{S}/\text{cm}$ is recommended) given on the vial label.
- Allow 15 min. for the complete dissolution of the lyophilisate.
- Mix gently after reconstitution.
- After reconstitution, do not keep a vial more than 5 hours at 18...25 °C or 24 hours at 2...8 °C.
- It is recommended that the contents of a reconstituted vial be divided into aliquots (volume of 175 μL is recommended to measure 2 replicates of control) which may then be stored frozen at -20 °C for a maximum period of 1 month.
- Use one of the tubes immediately for measurement.
- After thawing an aliquot, mix gently and use immediately for measurement.
- Once thawed, a control aliquot must not be refrozen.
- The bar code stick-on labels are used for identifying the controls when assayed on B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD.
- The control kit bar code card must be entered for each new lot of control.
- For further information see the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.

Quality Control

Good laboratory practice requires that control samples are measured regularly to ensure the quality of the results obtained. These samples must be processed exactly the same way as the assay samples, and it is recommended that the results be analysed using appropriate statistical methods. If desired, B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD can automatically check the quality of assays at intervals, by statistical analysis on the basis of Levey Jennings graphs.

It is necessary to comply with national quality assurance guidelines for quantitative tests in the medical laboratory (current version). For instance, test accuracy and precision should be monitored by means of laboratory in house and/or commercially available control materials. If unacceptable control values are obtained, proceed as outlined in standard laboratory diagnostic procedures to determine the cause and implement corrective measures.

Reference Range

A) Expected values

From EDTA plasma samples collected from 347 singleton pregnancies with a normal pregnancy outcome, levels of sFlt-1 and PIGF were determined in parallel and the sFlt-1/PIGF ratio was calculated.[10]

sFlt-1 KRYPTOR (pg/ml)

Gestational age (completed)	n	5 th percentile (pg/mL)	Median (pg/mL)	95 th percentile (pg/mL)
-----------------------------	---	------------------------------------	----------------	-------------------------------------

weeks)				
14-19	21	764	1 515	2 545
20-24	60	626	1 319	4 443
25-28	66	612	1 231	3 958
29-33	44	697	1 395	3 722
34-36	39	1023	2 971	8 492
37-delivery	117	961	3 079	7 524

PIGF plus KRYPTOR (pg/mL)

Gestational age (completed weeks)	n	5 th percentile (pg/mL)	Median (pg/mL)	95 th percentile (pg/mL)
14-19	21	43	101	348
20-24	60	183	391	868
25-28	66	145	459	1 000
29-33	44	71	507	1 157
34-36	39	41	170	752
37-delivery	117	37	131	750

sFlt-1/PIGF ratio

Gestational age (completed weeks)	n	5 th percentile	Median	95 th percentile
14-19	21	5	14	26
20-24	60	1	4	17
25-28	66	1	3	9
29-33	44	1	3	35
34-36	39	2	17	158
37-delivery	117	2	29	131

Note: It is recommended that each laboratory establish its own reference ranges based on representative patient collectives and/or test the validity of the manufacturer's commercial test kit data. The values above have been established on serum and are given only as an example

B) Clinical sensitivity and specificity

In an external study, using the B·R·A·H·M·S PIGF plus KRYPTOR and B·R·A·H·M·S sFlt-1 KRYPTOR assays in parallel on samples from 109 pregnant women with normal pregnancy outcome and 20 patients with preeclampsia, the published cut-off [7] for the sFlt-1/PIGF ratio of 85 was applied

to calculate sensitivity and specificity for this cohort. At a cut-off of 85 the sensitivity was calculated at 95 % and the specificity at 84 %.

All pregnancies were singleton pregnancies. Pre-eclampsia was defined as new onset of hypertension and proteinuria after 20 weeks of gestation. A pre-eclamptic pregnancy was defined as early-onset if clinical signs appeared before 34 weeks of gestation.

C) Method comparison

A comparison of the BRAHMS sFlt-1 KRYPTOR assay (y) with a commercially available automated sFlt-1 assay (x) using clinical samples gave the following correlation (pg/mL):

Number of samples measured: 150

Passing/Bablok

$$y = 0.93x - 267$$

$$r = 0.996$$

The sample concentrations were between 686 and 43598 pg/mL.

Assay Performance

Detection limit

The **Limit of Detection (LOD)**, and **Limit of Blank (LOB)** have been determined following CLSI EP17-A2 guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). The LOB has been assessed as being 16.7 pg/mL. The LOD has been assessed as being 22 pg/mL.

Sensitivity

The **Functional Assay Sensitivity (FAS)**, detected by inter-assay precision of 20 %, is estimated in accordance with CLSI EP05-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods). It has been assessed as being 29 pg/mL.

The **Limit of Quantitation (LOQ)** is the lowest reported concentration level with total error (imprecision and bias) at 40 % The LOQ is estimated in accordance with CLSI EP17-A guideline (Protocols for Determination of Limits of Detection and Limits of Quantitation). It has been assessed as being 34 pg/mL.

Specificity

The antibodies used in this assay are specific for the measurement of sFlt-1. The cross reactivity is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with sFlt-1 concentrations of 1 586 pg/mL or 10 248 pg/mL and found to cross-react at the concentrations indicated below:

Tested substance	Tested concentration	Observed cross reactivity
PIGF-1	10 000 pg/mL	No cross-reactivity
PIGF-2	10 000 pg/mL	No cross-reactivity
PIGF-3	10 000 pg/mL	No cross-reactivity
VEGF-A	10 000 pg/mL	No cross-reactivity
VEGF-B	10 000 pg/mL	No cross-reactivity
VEGF/PIGF heterodimer	10 000 pg/mL	No cross-reactivity

The **interference testing** is performed in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with sFlt-1 concentrations of 1 586 pg/mL or 10 248 pg/mL and found not to interfere at the concentrations indicated below:

Tested substance	Tested concentration
Acetaminophen	200 mg/L
Acetylsalicylic Acid	652 mg/L
L-Ascorbic acid	30 mg/L
Caffeine	59 mg/L
Calcium	200 mg/L
Ethanol	5% (v/v)
Folic Acid	18 µg/L
Gentamicin	10 mg/L
Heparin	3 000 U/L
Ibuprofen	500 mg/L
A-Methyldopa	15 mg/L
Metoprolol	5 mg/L
Iron	150 µg/dL
Dihydralazin	0.2 g/L
<u>Accuracy / Linearity</u>	
The linearity is determined in accordance with the CLSI EP6-A guideline (Evaluation of the Linearity of Quantitative Measurement Procedures); based on dilution with recombinant sFlt-1 the linearity of the assay is demonstrated from the LOQ up to 90 000 pg/mL. Patient samples cannot be diluted.	
<u>Accuracy / "High Dose Hook" Effect</u>	
No "High Dose Hook" Effect up to 400 000 pg/mL.	

Precision/Reproducibility for Intra Assay CV	
This was evaluated using samples of different concentrations in accordance with the CLSI EP5-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods).	
Concentration range	Within-assay CV%
599 pg/mL	1.0 %
2 696 pg/mL	0.8 %
31 173 pg/mL	0.8 %
50 452 pg/mL	0.9 %

Precision/ Reproducibility for Inter Assay CV	
This was evaluated using samples of different concentrations in accordance with the CLSI EP5-A2 guideline (Evaluation of Precision Performance of Quantitative Measurement Methods).	
Concentration range	Between-assay CV %
1 559 pg/mL	4.2 %
3 055 pg/mL	3.7 %

31 173 pg/mL	3.0 %
53 462 pg/mL	3.9 %

Disturbing Factors

This is determined in accordance with the CLSI EP7-A2 guideline (Interference Testing in Clinical Chemistry). The following substances were added to serum with PIGF concentrations of 1 586 pg/mL or 10 248 pg/mL and found not to interfere at the concentrations indicated below:

Factor	Description
Free unconjugated bilirubin	No significant effect up to 250 mg/L
Hemoglobin	No significant effect up to 5 g/L
Human serum albumin	No significant effect up to 60 g/L
Triglycerids	No significant effect up to 30 g/L

Traceability

An International sFlt-1- Reference Preparation is not available. Therefore the assay is calibrated using a reference preparation of recombinant human sFlt-1 as highest available standard. Dilutions of this preparation are checked by regression analysis and compared to theoretical expected results. Quality Control of new calibrator manufacturing ensures reproducibility and stability of calibration along time (according to internal procedure).

Storage conditions

Storage (stability)

Product stability tests are carried out according to EN 23640:2002 (In vitro diagnostic reagent stability tests) and also according to the CLSI Guidelines EP25-A (In vitro diagnostic reagent stability assessment). The kit should be stored in the manufacturer's package at the temperature range between 2 °C and 8 °C during the entire shelf life in cold chambers or refrigerators that secure the specified temperature mode with daily temperature recording.

Storage (stability) after opening

The reagent kit is stable after opening for 29 days in the B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD analyzer at a temperature of 2 to 8°C with a calibration period of 15 days.

Transportation conditions

The product should be transported in the manufacturer's original package that protect the product from external impact, transportation, means in roofed vehicles according to the Transportation Rules applicable.

Heat impact test confirms stability at the temperature between 2°C and 8°C during the entire shelf life and 15 days at a temperature between 2°C and 8°C on board B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD analyzer.

During cargo handling and transportation, it is necessary to strictly comply with the requirements of handling signs applied on the package and the transport containers.

If the transportation and storage conditions are not followed, the medical device must not be used!

Shelf life	
Maximum 18 months	
Disposal	
Disposal of the product must be carried out in accordance with the local, regional and national legislation.	
When disposing and destroying expired reagents, as well as waste and used reagents, it is necessary to act in accordance with the hygienic standards and follow official instructions for disposal and destruction of potentially infectious and chemical material applicable in the Russian Federation and in the specific health facility. The key issues in the waste management strategy are the optimization of tactics and approaches to final products, as well as minimization of possible environmental pollution. Given the nature of the used reagents waste material - containing human biomaterial, the most optimal is thermal disposal and destruction at a temperature of 1000°C. Waste should be stored in an airtight container marked "Biologically dangerous." Employees must undergo occupational safety and health training. For expired reagents we also allow the method of thermal utilization and destruction at a temperature of 1000°C. All samples and materials used and generated during research, including reagents that have expired, should be disposed of in accordance with SanPiN 2.1.7.2790-10 as potentially infected medical waste (class B). (When using the kit, waste of classes A, B and G are generated, which are classified and destroyed (disposed of) in accordance with SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for medical waste management" and MU 287-113 "Methodological guidelines for disinfection , pre-sterilization cleaning and sterilization of medical devices". In order to avoid possible environmental pollution, expired reagents as well as waste from reagents used must be disposed of in accordance with the local government and/or hospital regulations. Proper destruction and disposal will support the conservation of natural resources and help protect human health and the environment).	
Product	Products must be disposed of as laboratory chemical waste in accordance with the local regulations. - Please contact the responsible authority. - Used reagent plates and reagent kits should be disposed of as potentially infectious laboratory waste.
Contaminated packaging	- Dispose according to the local regulations. - Recommendations: Remove residues. Return empty containers to local waste recycling facilities.
Pollutes	- Remove pollutes with absorbing paper. - All materials used for cleaning up must be disposed of as infectious laboratory waste.

Manufacturer's authorized representative

Thermo Fischer Scientific Joint Stock Company (Thermo Fischer Scientific JSC)

73, Kubinskaya St., Letter A, Bldg 1, St. Petersburg 196240

Tel.: (812) 703-42-15, facsimile: (812)703-42-16

e-mail: info.lcp.spb@thermofisher.com

Warranties

The Manufacturer guarantees the conformity of the characteristics of the product stated in the operating documentation, subject to use in accordance with the instructions and intended use, provided by the manufacturer

For all questions, technical advice and support please contact the **manufacturer's authorized representative in the Russian Federation:**

Thermo Fischer Scientific Joint Stock Company (Thermo Fischer Scientific JSC), Russia

Address: 73, Kubinskaya St., Letter A, Bldg 1, St. Petersburg 196240

Tel.: (812) 703-42-15, facsimile: (812)703-42-16,

e-mail: info.lcp.spb@thermofisher.com.

Bibliography

Analyte

1. Ghulmiyyah L, et al (2012), Seminars in Perinatology; 36:56-59
2. Powe C, et al (2011), Circulation; 123: 2856-69
3. Levine RJ et al. (2004). N Engl J Med; 350:672-83.
4. Lam C, et al (2005). Hypertension; 46: 1077-85
5. De Vivo, et al. (2008). Acta Obstet et Gynecol; 87: 837-842
6. Thadhani R, et al. (2004). J Clin Endo; 89(2):770-775.
7. Verlohren S, et al (2012). Am J Obstet Gynecol. 206(1):58.e1-8.
8. Rana S, et al (2012). Circulation; 125 (7):911-919
9. Andersen L et al (2015). J Am Soc Hypertens; 9(2):86-96.
10. Salahuddin S et al (2016). Hypertens Pregnancy, 35(3):330-45

Technology

- a) Mathis G, Lehn JM.: Trace – Another Story of Time. Isotopics 1995; Vol. 9.
- b) Mathis G. ,Clin. Chem. 1993; 39:1953-9.
- c) B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD User Manual.
- d) B·R·A·H·M·S KRYPTOR compact PLUS/KRYPTOR GOLD Manual.

The list of applicable standards

1. EN ISO 15223-1:2016

Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied- Part 1: General requirements (ISO 15223-1:2016)

2. EN ISO 13485:2016

Medical devices - Quality management systems - Requirements for regulatory purposes

3. EN 13612:2002

Performance evaluation of in vitro diagnostic medical devices

EN 13612:2002/AC:2002

4. EN 13641:2002

Elimination or reduction of risk of infection related to in vitro diagnostic reagents

5. EN ISO 14971:2012

Medical devices - Application of risk management to medical devices

(ISO 14971:2007, Corrected version 2007-10-01)

6. EN ISO 17511:2003

In vitro diagnostic medical devices - Measurement of quantities in biological samples - Metrological traceability of values assigned to calibrators and control materials

(ISO 17511:2003)

7. EN ISO 18113-1:2011

In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 1: Terms, definitions and general requirements

(ISO 18113-2:2009)

8. EN ISO 18113-2:2011

In vitro diagnostic medical devices - Information supplied by the manufacturer (labelling) - Part 2: In vitro diagnostic reagents for professional use

(ISO 18113-2:2009)

9. EN ISO 23640:2015

Evaluation of stability of in vitro diagnostic reagents

(ISO 23640:2011)

10. IEC 62366:2015 / EN 62366-1:2015

Medical devices - Part 1: Application of usability engineering to medical devices (IEC 62366-1:2015 + COR1:2016);

■ Symbols

Symbols used in instruction for Use and Product Labelling of B-R-A-H-M-S KRYPTOR compact PLUS/KRYPTOR GOLD products.

Symbol	Usage	Symbol	Usage	Symbol	Usage
	Reference to the Intended use of the Medical Device		In Vitro Diagnostic Medical Device		Batch Code
	Contents		Calibrator		Control
	Buffer		B-R-A-H-M-S KRYPTOR SOLUTION 1/B-R-A-H-M-S KRYPTOR compact SOLUTION 1/ B-R-A-H-M-S Solution 1 KRYPTOR GOLD		B-R-A-H-M-S KRYPTOR SOLUTION 2/B-R-A-H-M-S KRYPTOR compact SOLUTION 2/ B-R-A-H-M-S Solution 2 KRYPTOR GOLD
	B-R-A-H-M-S KRYPTOR SOLUTION 3/B-R-A-H-M-S KRYPTOR compact SOLUTION 3/ B-R-A-H-M-S Solution 3 KRYPTOR GOLD		B-R-A-H-M-S KRYPTOR SOLUTION 4/B-R-A-H-M-S KRYPTOR compact SOLUTION 4/ B-R-A-H-M-S Solution 4 KRYPTOR GOLD		Bags contained
	Bags		Plates contained		Plates
	Vials contained		Vials		Val
	Use given Volume of distilled Water (conductivity of less than 50 µS/cm is recommended) for Reconstitution, e.g. 0.75 mL		Lyophilized, Freeze Dried		Reconstitute
	Manufacturer		Use by		Green Dot according to German Law
	Registered Trade Mark		Trade Mark		Article Number/Catalogue Number
	Contains sufficient for (Number of) tests, e.g. 50		Consult Instructions for Use		See Accompanying Compact Disk
	Biohazard		Wear Protective Gloves		Wear Safety Glasses
	Wash hands		General Regulatory Sign		General Prohibitive Sign
	Do not Smoke		Do not Eat and Drink		GHS05, Corrosion
	GHS07, Exclamation Mark		Trade Mark for TRACE™-technology		CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices
	CE Conformity Marking According to Directive 98/79/EC on In Vitro Diagnostic Medical Devices, Annex II with Reg. Number of Notified Body		Temperature Limitation		Do not Reuse
	Caution, consult accompanying documents		Accidental Release Measures		Waste

The notary asked before if he was or is active in the matter which is object of this document outside of his notarial activity. The question was negated.

Hereby I certify that both members of the management team of the company B.R.A.H.M.S GmbH

1. Dr. Bernhard **C i o m m e r**
born January 28, 1955
domiciled in D-13593 Berlin
Folkungerstr. 21,

2. Dr. Jörg-Markus **P a u l**
born August 02, 1967,
domiciled in D-13503 Berlin
Reimerswalder Steig 22,

both commercially domiciled in 16761 Hennigsdorf, Neuendorfstrasse 25,
acknowledged their signatures before me.

The identity of the deponents is personally known to me.

At the same time I certify according to § 21 BNotO due to the presented certified extract from the statement of the Register of Companies of the District Court of Neuruppin dated March 09, 2020 that with HRB-No. 8872 NP B.R.A.H.M.S GmbH domiciled in Hennigsdorf is registered and that Dr. Bernhard Ciommer as well as Dr. Jörg-Markus Paul are registered as management team members and are both authorised to represent the company.

Berlin, in the business premises of the notary Lutz Baxmeier,
Berliner Strasse 32, D-13467 Berlin

March 10, 2020

Scroll no. 233/2020

Baxmeier
Notary in Berlin
Federal Republic of Germany

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