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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HANGZHOU BIOTEST BIOTECH CO., LTD. on the annexed APPLICATION INSTRUCTION is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Chen Cheng

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(Date: Mar. 21, 2019)

证明书查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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B i o t e s T Hangzhou Biotest Biotech Co., Ltd.

APPLICATION INSTRUCTION

Set of reagents «Express test for troponin I assay in the human plasma, serum and whole blood by the immunochromatography method (CARDIAC TROPONIN I RAPID TEST)

• Designation

The set of reagents «Express test for troponin I assay in the human plasma, serum and whole blood by the immunochromatography method (CARDIAC TROPONIN I RAPID TEST) is designed for the qualitative determination of the cardiac troponin I (cTnI) with the purpose of diagnostics of the myocardial infarction (MI) in the whole blood, serum or plasma by the immunochromatography method.

For the single use of the appliance according to its purpose.

For the clinical laboratory diagnostics in vitro.

Potential customers of the appliance

Laboratory subdivisions of healthcare centers.

Professional level of potential users

Physician for the clinical laboratory diagnostics, medical laboratory technician.

Indications

Qualitative detection of cardiac troponin I (cTnI) to diagnose myocardial infarction (MI).

Contraindications

When used by specially trained personnel and taking into account the intended use of contraindications are not identified.

Possible side effects

When used by specially trained personnel and taking into account the intended use of side effects are not identified.

Target analyte description, data on its scientific validity

Cardiac Troponin I (cTnI) is a protein, contained in the cardiac muscle and having the molecular weight of 22.5 kDa. Troponin I is a part of the complex of three sub-units, including Troponin T and Troponin C. Together with tropomyosine this structural complex forms the basic component, regulating the calcium-sensitive activity of the AT-phase of the actomyosine in striated skeletal and cardiac muscles. After the myocardium damage Troponin I gets into blood 4-6 hours after the beginning of pain. The release of cTnI is equivalent to such of CK-MB, but while CK-MB levels begin normal 72 hours later, Troponin I levels remain increased during 6-10 days, providing the larger myocardium damage detection level. The high specificity of cTnI measurements in order to reveal myocardium damages can be demonstrated in such conditions as perioperative period, after marathons and the blunt injury of the chest. The cTnI release also has been registered at heart diseases, not referred to the acute myocardial infarction (AMI), such as unstable angina pectoris, congestive cardiac failure and ischemic diseases, referred to the coronary artery bypass grafting. Due to its high specificity and sensitivity in the myocardium tissue Troponin I lately became the most preferable biomarker of the myocardium infarction.

CARDIAC TROPONIN I RAPID TEST is a simple test, using the composition of particles, covered with antibodies to the anti -cTnI, and the capture reagent to reveal the cTnI in the whole blood, serum or plasma. The minimal detection level is 0.5 ng/ml.

Complete set

1.1. Test holder	Rectangular tablet 70x22x7 (length x width x thickness), made of white plastic with two windows: S-round – for the sample input, II- oblong (test area) with two markings: T (Test) – closer to the round window, C (Control) – closer to the upper part of the oblong window	1.1. Test holder
1.2. Buffer solution for the sample dilution	Transparent colorless liquid	1.2. Buffer solution for the sample dilution
1.3. Disposable droppers	Colorless transparent plastic dropper with balloon-like extension one the one end, volume 40 mcl	1.3. Disposable droppers

Note:

1. One package of laminated foil comprises 1 test holder and 1 disposable dropper with the desiccant.

2. The appliance does not contain any materials, directly or indirectly contacting the organism of the patient and personnel, using the appliance, while observing requirements of the maintenance documentation (application instruction).

3. The appliance does not contain pharmaceuticals or pharmaceutical substances.

Set components are packed in a box, the box contains the application instruction.

MAIN CONSUMER PROPERTIES

The basic set option allows to study 10 samples.

Appliance sterilization methods

The appliance does not require sterilization.

Appliance software

Absent.

Technical maintenance and repair of the appliance

The item does not require technical maintenance and cannot be repaired.

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The set of reagents "Express test for the troponin I assay in the human plasma, serum and whole blood by the immunochromatography method" (CARDIAC TROPONIN I RAPID TEST) is a chromatographic immunotest, including anti-cTnI-antibodies, conjugated with colloidal gold particles, antibodies to anti-cTnI on the test line and IgG antibodies on the control line.

During the testing the whole blood, serum or plasma sample reacts with particles, covered with antibodies to anti-cTnI, on the test strip. Then the mixture chromatographically migrates up on

the membrane, influenced by capillary forces, and reacts with anti -cTnI –antibodies on the membrane in the test line area.

If the sample contains the cardiac Troponin I (cTnI), the color line will show in the test line area, witnessing the positive result.

If the sample does not contain the cardiac Troponin I (cTnI), the color line will not be reflected in this area, what witnesses the negative result.

As procedure control the color line will always appear in the control line area, showing, that the due sample volume has been inserted and the membrane is capillary moisturized.

COLLECTION OF SAMPLES AND ITS PREPARATION

Serum or plasma(K2EDTA, 3,2% sodium citrate, sodium heparin, lithiumheparin) of the human blood, volume of 50 mcl, whole blood, volume of 75 mcl.

- To collect whole blood samples from the finger;
- Wash patient's hand with soap and warm water or clean it with the sponge with alcohol.

Let dry.

- Massage the hand without touching the place of puncture in the direction to the tip of the long or third finger.

- Puncture the skin with a sterile lancette. Wipe out first drops of blood appeared.

- Carefully rub the hand from the wrist to the finger in order the round blood drop appears.

- Place the whole finger blood sample in the test holder with the capillary pipe:

- Touch the blood drop with the end of the capillary pipe to fill the pipe for about 75 mcl.

Avoid air bubbles.

- Place the bolus in the upper part of the capillary pipe, then compress the bolus in order to force out the whole blood on the support for test holder samples.

- Immediately separate the serum or plasma from the blood to avoid the hemolysis. Use only clean non-hemolysed samples.

- The test should be performed immediately after the collection of samples. Do not leave samples at the room temperature for a long time. Serum and plasma samples can be stored at 2 - 8°C for 3. If the long time storage is necessary, samples are stored below -20°C. If the test is performed within 2 days since the blood collection, the whole blood from the vein should be stored at 2 - 8°C. Don't freeze whole blood samples. The test of the whole blood from the finger should be performed immediately.

- Before the test bring samples to the room temperature. Before the test frozen samples should thaw completely and be mixed. It is not allowed to freeze and defrost samples again.

- If samples are to be transported, it should be packed in cofnromity with local rules for the transportation of disease agents.

- Blood samples should be collected in conformity with the good practice by vein puncture method (GOST R 53079.4-2008 «Provision for the quality of clinical laboratory studies »).

• ANALYTIC AND DIAGNOSTIC FEATURES

- Analytic sensitivity (minimum detectable concentration) – 0,5ng/ml.
- The test sensitivity is determined by control biological indices. Containing the cardiac Troponin I (cTnI), as the per cent of samples, determined by the test as positive, and makes 100%.
- The set specificity is determined by control negative biological samples, not containing the cardiac Troponin I (cTnI), as the per cent of samples, determined by the set as negative, and makes 100%.
- The time for the test to get the determined color is NMT 20 min.

• SAFETY MEASURES

• Only for professional in vitro diagnostics. Don't use after the expiration of the shelf life.

- The test should be stored in the sealed package till the direct use.

It is prohibited to take food, drink or smoke in places for the work with samples.

All samples should be treated as samples, containing agents of infectious diseases. During all procedures observe provided safety measures for the work with microbiological sources of danger, as we;; as observe standard procedures for the due disposal of samples.

While working with samples wear protective clothes, including the laboratory coat, disposable gloves and goggles.

All used tests, samples and contaminated materials should be disposed in conformity with local norms.

Humidity and temperature can impact results.

Do not mix components from different batches.

Safety measures when operating the set – observance of rules "Instructions on measures for the prophylactics of spreading of infectious diseases while working and clinical & diagnostic laboratories of healthcare centers " (App. by the Ministry of Healthcare of USSR on January 17 1991), GOST R 52905–2007 (ISO 15190:2003).

The set is biologically safe, anyway, studied samples should be treated as potentially infect material, while strictly observing safety requirements at serologic studies of materials, potentially infected by agents of I and II pathogenicity groups (SR 1.32322-08).

The set of reagents in conformity with SRN 2.1.7.2790-10 «Sanitary & epidemiologic requirements for the treatment of medical waste » is referred to the class A – epidemiologically save waste, which composition is close to the one of domestic garbage.

Sets of reagents should be disposed of or destroyed, disinfected in conformity with SRN 2.1.7.2790-10 "Sanitary & epidemiologic requirements for the treatment of medical waste " and MI 287-113 "Methodic instructions for the disinfection, pre-sterilization cleaning and sterilization of medical appliances ".

The list of risks, qualified in the process of the risk analysis, as well as its management methods are represented by a separate document (see Risks Analysis).

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Necessary materials, not included in the test

- Test tube for sampling of plasma, serum or blue blood.
- Laboratory centrifuge.
- Timer.
- Lancets sterile (only for whole blood from a finger).
- Heparinized capillary pipes and distribution flasks, measuring bolus (only for the whole blood from the finger).
- Control panel of positive biologic samples, containing the cardiac Troponin I (cTnI (Lot No:CTNI180514 -01).
- Control panel of negative biologic samples, not containing the cardiac Troponin I (cTnI (Lot No:170228).

Preparation for work

Before the test set components should be heated to the room temperature.

Before opening the package bring it to the room temperature.

Take the test holder from the sealed package and use it during one hour.

Place the test holder on the clean and even surface.

Analysis method


For serum and plasma samples:

Holding the dropper vertically, place 2 drops of the serum or plasma (about 50 mcl) on the sample support and then add 1 dropper of the buffer solution (about 40 mcl) and switch on the timer. See the picture below.

For whole blood samples from the vein: Holding the dropper vertically, place 3 drops of the whole blood (about 75 mcl) on the sample support and then add 1 drop of the buffer solution (about 40 mcl) and switch on the timer. See the picture below.

For whole blood from the finger: While using the capillary pipe: Fill the capillary pipe and place about 75 mcl of the whole blood from the finger on the sample support in the test holder, add 1 drop of the buffer solution (about 40 mcl) and switch on the timer. See the figure below.

Add when the color line(s) appear(s). Count results in 10 minutes. Do not read results in 20 minutes.

• INTERPRETATION OF RESULTS

(See illustration above)

POSITIVE RESULT:* Two lines appear. One color line should appear in the control line area (C), and another clear color line— in the test line area (T).

* **NOTE:** The intensiveness of coloring in the test line area (T) depends on the Troponin I concentration in the sample. In this juncture any level of coloration in the test line area (T) is considered as positive result.

NEGATIVE RESULT: The appearance of one color line in the control line area (C). No line in the test line area (T).

INVALID RESULT: No control line. Insufficient sample volume or incorrect test procedure are most probable reasons of the control line absence. In such case study the procedure again and repeat the test with the new test holder. If the problem is not solved, immediately stop using the test system and contact the local distributor.

• QUALITY CONTROL

The test procedure comprises the test procedure control with the LFIA method. The color line in the control line area (C) is considered as the internal instrument for the test procedure control. Its appearance confirms the due flowing in the membrane.

Control standards are not included in the test holder set; anyway, according to principles of the good laboratory practice is it recommended to use both positive and negative control to confirm the test performance procedure and to check the due character of the test performance. Samples for external quality control have to be tested according to the standards of quality control established in your laboratory.

• RESTRICTIONS

- Only for *in vitro* diagnostics.
- The set is used for cTnI detection in whole blood, serum or plasma samples (K2EDTA, sodium citrate, sodium heparin, lithium heparin).
- This test set cannot determine the quantitative value of the cTnI or the level of increase of its concentration.
- The set shows only the cTnI presence in the samples and cannot be used as sole criteria for the diagnostics of the myocardium infarction (MI).
- In order to confirm obtained results additional analyses of samples should be performed in conformity with guidelines of local healthcare authorities, for example, the immunoenzyme analysis.
- Results should be interpreted only in combination with other clinical data, being at physician's disposal.
- This test is designed only for screening purposes.
- The negative test result does not exclude the possibility of the myocardium infarction. So for the exact diagnosis the result, obtained with the set, should be used in combination with the clinical data.
- The false positive result can be observed at patients with unstable stenocardia, intra coronary thrombosis or a chronic renal failure and also at some other diseases.
- The false-negative result does not exclude a myocardial infarction since the research can be conducted prior to release of troponin.
-

• Attention!

• This set will provide the correct result only if the application instruction is strictly observed.

• EXPECTED VALUES

Results of set use were compared with results of use of available ones – the leading commercial test cTnI EIA with use of clinical samples (reliable marketed products), and in the course of comparison was pointed out the good test sensibility, making 98,8%, with the specificity of 98,9%.

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The set of reagents "Express test for the troponin I assay in the human plasma, serum and whole blood by the immunochromatography method" (CARDIAC TROPONIN I RAPID TEST) (cat.№. CTNI-C41) was evaluated on the basis of samples, obtained from evidently healthy blood donors, as well as on the basis of samples, obtained from patients. Results show that the test sensibility is 99,8%, and the specificity - 99.5% relatively to other express tests.

Method		Other express test		Final result
The set of reagents “Express test for the troponin I assay in the human plasma, serum and whole blood by the immunochromatography method” (CARDIAC TROPONIN I RAPID TEST) (cat.№. CTNI-C41)	Results	Positive result	Negative result	
	Positive result	158	7	165
	Negative result	2	603	605
Final result		160	610	770

*Confidence intervals

- **Precision**

For results of the same analysis

The between-run precision was determined for 15 series of following 5 samples: samples levels at negative, cTnI 1,0ng / mL, cTnI 5.0ng / mL, cTnI 10ng / mL, cTnI 40ng / mL. Samples were determined correctly in more than 99% cases.

For results of different analyses

The interserial precision was determined for 15 independent analyses with use of same three samples: negative, cTnI 1,0ng / mL, cTnI 5.0ng / mL, cTnI 10ng / mL, cTnI 40ng / mL. Samples were determined correctly in more than 99% cases.

Cross-reactivity

The set of reagents "Express test for the troponin I assay in the human plasma, serum and whole blood by the immunochromatography method" (CARDIAC TROPONIN I RAPID TEST) passed the check on the basis of following positive samples: 10 000 ng / ml Skeletal Troponin I, 2000 ng / ml Troponin T, 20 000 ng / ml cardiac myosin, rheumatoid factor (RF), HAMA, HCV, HBsAg, HBsAb, HBeAg, HBeAb, HBcAb, syphilis, HIV, anti-H.pylori, MONO, anti-CMV, anti-RUBELLA and anti-TOXO. Results did not show any cross-reactivity.

Interfering substances

Following potentially interfering substances were added to both negative and positive samples

Acetaminophen: 20 mg/dl	Caffeine: 20 mg/dl	Bilirubin: 1000 mg/dl
Acetylsalicylic acid: 20 mg/dl	Creatine: 200 mg/dl	Oxalic acid: 600 mg/dl
Gentisic acid: 20 mg/dl	Ascorbic acid: 20 mg/dl	Cholesterol: 800 mg/dl
Hemoglobin: 1,000 mg/dl	Albumen: 10,500 mg/dl	Triglycerides: 1,600 mg/dl

According to analysis results no substance is interfering in given concentrations.

SHELF LIFE

Set shelf-life – 2 years.

The set with expired shelf life cannot be used.

- **STORAGE AND STABILITY**

Storage (stability)

The set should be stored in the manufacturer's package at the room temperature or in the refrigerator (at 2 - 30°C), providing for the regulated temperature mode with daily registration of the temperature) during the whole shelf life. DON'T FREEZE.

The set is stable before the expiration of the storage period, indicated on the sealed package.

Storage (stability) after opening

The opened buffer solution should be stored at the room temperature or in the refrigerator (at 2 - 30°C) within the whole shelf life, indicated on the package. DON'T FREEZE. The content of one individual package is designated for single use. Don't use after the shelf life is expired.

Storage (stability) (imitation of transportation conditions)

Should be transported in the original manufacturer's package, protecting the production from external effects, on all kinds of transport in covered vehicles according to rules of transportation, valid at this kind of transport at 2 - 30 °C during the whole shelf life.

While carrying out the stevedoring and transportation requirements of manipulation signs, put on the package and transportation tare, should be strictly observed.

The medical appliance should not be used if transportation and storage terms of the medical appliance are not observed!

RELEASE TERMS

For healthcare institutions.

MANUFACTURER'S GUARANTEES

The manufacturing site has got the quality management system, certified as corresponding with requirements of international standards of the Directive 98/79/EC, ISO 13485.

On matters, referred to the quality and circulation of the set of reagents **"Express test for the troponin I assay in the human plasma, serum and whole blood by the immunochromatography method" (CARDIAC TROPONIN I RAPID TEST)**, please, contact: ROTANA Limited Liability Company, ROTANA LLC, 117342, Russia, city of Moscow, General Antonova str., 3A, tel./fax +7 (495) 646-28-61, med@rotana-rf.ru.

Manufacturer:

Hangzhou Biotest Biotech Co., Ltd, Китай
17#, Futai Road, Zhongtai Street, Yuhang District, Hangzhou, 311121-P.R. China
Tel.: +86-571-88630899.
Fax: +86-571-88633388.

Authorized representative of the manufacturer in the territory of the RF:

ROTANA Limited Liability Company, ROTANA LLC
117342, Russia, city of Moscow, Generala Antonova str., 3A, tel./fax +7 (495) 646-28-61,
med@rotana-rf.ru

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International Registration Director

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