CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

CIVIL CODE § 1189


A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California)
 County of Orange)
 On June 29, 2020 before me, Steve Oberlander, Notary Public
 Date Here Insert Name and Title of the Officer
 personally appeared Michelle Mirzaee
 Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Steve Oberlander
 Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Appendix For the IFV

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Appendix for the Instruction for use**Onyx Liquid Embolic System (LES) (hereafter referred to as «Onyx LES»)****1. Name of medical device**

Onyx Liquid Embolic System (LES), versions:

I.

1. Onyx Liquid Embolic System (LES), model 105-7000-060, contains:

Onyx-18 LES (6% EVOH) in a sterile vial (1.5 ml) – 1 pc.;

DMSO in a sterile vial (1.5 ml) – 1 pc.;

Onyx LES delivery syringe sterile (1 ml) – 2 pcs.;

DMSO syringe sterile (1 ml) – 1 pc.;

Instruction for use for Onyx Liquid Embolic System (LES).

2. Onyx Syringe-Catheter Interface Adapter, model 103-1207 (if needed), contains:

Adapter – 20 pcs.;

Instruction for use for Onyx Syringe-Catheter Interface Adapter – 1 pc.

II.

1. Onyx Liquid Embolic System (LES), model 105-7000-080, contains:

Onyx-34 LES (8% EVOH) in a sterile vial (1.5 ml) – 1 pc.;

DMSO in a sterile vial (1.5 ml) – 1 pc.;

Onyx LES delivery syringe sterile (1 ml) – 2 pcs.;

DMSO syringe sterile (1 ml) – 1 pc.;

Instruction for use for Onyx Liquid Embolic System (LES).

2. Onyx Syringe-Catheter Interface Adapter, model 103-1207 (if needed), contains:

Adapter – 20 pcs.;

Instruction for use for Onyx Syringe-Catheter Interface Adapter – 1 pc.

III.

1. Onyx Liquid Embolic System (LES), model 105-7000-065, contains:

Onyx-20 LES (6,5% EVOH) in a sterile vial (1.5 ml) – 1 pc.;

DMSO in a sterile vial (1.5 ml) – 1 pc.;

Onyx LES delivery syringe sterile (1 ml) – 2 pcs.;

DMSO syringe sterile (1 ml) – 1 pc.;

Instruction for use for Onyx Liquid Embolic System (LES).

2. Onyx Syringe-Catheter Interface Adapter, model 103-1207 (if needed), contains:

Adapter – 20 pcs.;

Instruction for use for Onyx Syringe-Catheter Interface Adapter – 1 pc.

2. With respect to authorized representative of the manufacturer of the medical device – full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity

Authorized representative of the Manufacturer:

LLC «Medtronic»

Medtronic, LLC, 123112, Moscow, Presnenskaya Nab., b.10, floor 9, office III, room 41

Tel.: (495) 580-73-77,

Fax: (495) 580-73-78,

www.medtronic.ru

3. Technical characteristics of the medical device

The main parameters and characteristics of MD

Name	Characteristic	Units	Value
Onyx LES delivery sterile syringe	nominal volume	ml	1
	maximum "dead" space volume	ml	<0.142
	overall syringe length	mm	161.62±0.254
	size of the syringe at the widest part	mm	29.21±0.254
	nominal syringe outer diameter	mm	9.9±0.13
	nominal syringe inner diameter	mm	4.95-5.05
	value of scale division	ml	0.025
	volume between graduation lines with numbers	ml	0.1
	syringe barrel length	mm	86.487 ±0.254
	plunger length	mm	90.17±0.254
	plunger emergence length	mm	14.656±0.254
	silicon piston length	mm	6.731±0.127
	Luer connector taper rate	%	6
DMSO delivery sterile syringe	nominal volume	ml	1
	maximum "dead" space volume	ml	<0.07
	overall syringe length	mm	116.332±0.254
	size of the syringe at the widest part	mm	29.21±0.254
	nominal syringe outer diameter	mm	9.9±0.25
	nominal syringe inner diameter	mm	4.95-5.05
	value of scale division	ml	0.025
	volume between graduation lines with numbers	ml	0.1
	syringe barrel length	mm	86.614±0.254
	plunger length	mm	90.17±0.254

MEDTRONIC

9775 TOLEDO WAY
IRVINE, CA
92618

949-837-3700 [T]

Page 2 of 8

	plunger emergence length	m m	18.288±0.254
	silicon piston length	m m	6.731±0.127
	Luer connector taper rate	%	6
Onyx liquid embolic in a sterile vial	volume	ml	1.6 +/- 0.1
	nominal liquid viscosity	cSt	16.82±1.68 (Onyx 18), 20.34±2.03 (Onyx 20), 31.46±3.15 (Onyx 34),
	EVOH concentration	%	6.0 (Onyx 18), 6.5 (Onyx 20), 8.0 (Onyx 34)
	Density	g/ ml	1.72 ± 0.1 (Onyx 18), 1.72 ± 0.1 (Onyx 20), 1.71 ± 0.1 (Onyx 34)
	recommended solution injection speed	ml /m in	0.16
	maximum solution injection speed	ml /m in	0.3
	continuous shaking time	mi n	20
	maximum material expansion (37 °C)	%	≤5
	time of remaining of embolized area	mo nt hs	≥12
	The minimum vessel size, in which embolic liquid must be able to penetrate	m m	150
DMSO solution in a sterile vial	volume	ml	1.6 ± 0.1
	DMSO concentration in the solution	%	99.9
	recommended solution injection speed	ml /m in	≤ 0.16
	maximum solution injection speed	ml /m in	0.3
	Viscosity	cSt	1.43 ± 0.01

	Density	g/ ml	1.10
4. Information on the presence in the medical device of a medicinal product for human use, human and/or animal materials			
Not Applicable			
5. Information on the procedure of installation, adjustment, calibration and other actions necessary to put the medical device into operation			
Not Applicable			
6. List of national standards applied by the manufacturer			
Available upon customer request.			
7. Procedure and conditions for the disposal or destruction of a medical device			
This product does not harm the environment or human health if it would be right to exploit, transport, store and utilize. The system should be disposed per institutional guidelines and in accordance with the local regulations. After using the product it refers to the epidemiologically safe waste (Class A according to Russia Sanitary Regulations and Standards 2.1.7.2790). Wastes that do not have contact with the patient's body fluids and infectious patients.			

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© 2016 г. Национальная ассоциация