To whom it may concern

Letter of Declaration

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1. Kamal AH, Tefferi A, Pruthi RK. How to interpret and pursue an abnormal prothrombin time, activated partial thromboplastin time, and bleeding time in adults. Mayo Clin Proc 2007; 82: 864-73.
2. Leung LLK. Perioperative evaluation of bleeding diathesis. Hematology Am Soc Hematol Educ Program 2006; 457-61.
3. Eikelboom JW, Hirsh J. Monitoring unfractionated heparin with the aPTT: Time for a fresh look. Thromb Haemost 2006; 96: 547-52.
4. Bates SM, Johnston M, Hirsh J, Ginsberg JS. Use of a fixed activated partial thromboplastin time ratio to establish a therapeutic range for unfractionated heparin. Arch Intern Med 2001; 161: 385-91.

5. Hirsh J, Raschke R. Heparin and low-molecular-weight heparin. The seventh ACCP Conference on Antithrombotic and Thrombolytic Therapy. Chest 2004; 126: 188S-203S.
6. Olson JD, Arkin CF, Brandt JT, et al. College of American Pathologists Conference XXXI on laboratory monitoring of anticoagulant therapy: laboratory monitoring of unfractionated heparin therapy. Arch Pathol Lab Med. 1998; 122: 782-98.
7. Teruya J, West AG, Suell MN. Lupus Anticoagulant Assays. Questions answered and to be answered. Arch Pathol Lab Med. 2007; 131: 885-9.
8. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2008.
9. Hattersley PG, Hayse D. The effect of increased contact activation time on the activated partial thromboplastin time. Am J Clin Pathol 1976; 66: 479-82.
10. Ambrus JL, Schimert G, Lajos TZ, et al. Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery. J Med 1971; 2: 65-81.
11. Prasad RN, Koh SC, Viegas OA, Ratnam SS. Effects on hemostasis after two-year use of low dose combined oral contraceptives with gestodene or levonorgestrel. Clin Appl Thromb Hemost 1999; 51: 60-70.
12. Solomon GE, Hilgartner MW, Kutt H. Coagulation defects caused by diphenylhydantoin. Neurology 1972; 22: 1165-71.
13. Young DS, ed. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington: AACC Press, 2000.
14. Greinacher A, Warkentin TE. The direct thrombin inhibitor hirudin. Thromb Haemost 2008; 99: 819-29.
15. Soloway HB, Cox SP, Donahoo JV. Sensitivity of the activated partial thromboplastin time to heparin: effect of anticoagulant used for sample collection. Am J Clin Pathol 1973; 59: 760-2.
16. Lippi G, Franchini M, Montagnana M et al. Quality and reliability of routine coagulation testing: can we trust that sample? Blood Coag Fibrinol 2006; 17: 513-9.
17. CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline. CLSI document H47-A [ISBN 1-56238-301-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 1996.
18. CLSI. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline - Second Edition. CLSI document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2000.

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© 2010 Siemens Healthcare Diagnostics Products GmbH.

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Addendum to Instruction for Use
Reagents for in vitro studies of hemostasis:
Dade Actin FS Activated PTT Reagent

1. Indications and contraindications for the use of the medical device

Not applicable for this medical device.

2. Information about potential consumers of the medical device

For in vitro diagnostic use. For Professional Use. Only clinical laboratory professional user, medical laboratory technician can use this medical device.

3. Use, storage and transportation conditions

Shelf life - 24 months

Storage conditions

Temperature 2 - 8 °C. Atmospheric pressure. Humidity 10-80%.

Stored unopened at 2 to 8 °C, Actin FS can be used by the labeled expiry date.

Transportation conditions

Refrigerated reagent ships in storage conditions. Reagent is stable till expired date indicated on package.

Use conditions.

After use, store the closed container at 2 to 8 °C.

Once opened, the reagent is stable for 7 days at 2 to 15 °C. Do not freeze. If the reagent is left to stand, a green deposit may form consisting of ellagic acid and lipids. Before use, mix by inverting. Avoid contamination with plasma.

On-board stability

Information about on-board stability is specified in the Reference Guides (Application Sheets) for the different coagulation analyzers.

Signs of expiry: Deviations from the normal laboratory value in the determination of normal plasma or controls.

4. A procedure for disposal or destruction of the medical device

The generation of waste should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Waste should not be disposed of untreated to the sewer unless fully compliant with the requirements of all authorities with jurisdiction.

The wastes of Calcium Chloride are classified as Class B (epidemiologically dangerous wastes) in accordance with Sanitary and epidemiologic rules and standards SanPiN 2.1.7.2790-10 (Russian regulation).

5. The description of the sterilization method, sterilization validation methods (including bioburden, pyrogen and residual sterilant tests) and packaging validation methods (if the medical device is supplied in a sterile state)

Not applicable for this medical device.

6. The list of materials of animal and (or) human origin with the information on their biological compatibility and safety, source (donor) selection, sampling, processing, storage, and handling of such materials (if any)

Not applicable for this medical device.

7. Procedure precautions

Please read the Instruction for Use before using device. Follow the Instruction for Use (the IFU), disregard of the IFU can involve wrong results. Don't use the medical device with damage of packaging. Don't use samples that are different with those that are indicated in the IFU, this can involve wrong results.

Don't use the medical device with expired date. Check the expiration date before testing. Don't mix reagents from boxes with different lots.

8. The list of national and international regulations / standards, which correspond to the medical device

Manufacturing process, distribution and final quality control conform to requirements of Council Directive 98/79/EC and national standards.

Summary of Reference Standards

Standard Number	Title (International Standard)
EN ISO 13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes
EN 13612:2002 EN 13612:2002/AC:2002	Performance evaluation of in vitro diagnostic medical devices
EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
EN ISO 15223-1:2016	Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
EN ISO 18113-1:2011	In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 1: Terms, definitions and general requirements
EN ISO 18113-2:2011	In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 2: In vitro diagnostic reagents for professional use
EN ISO 23640:2015	In vitro diagnostic medical devices -- Evaluation of stability of in vitro diagnostic reagents
EN 62366-1:2015	Medical devices -- Part 1: Application of usability engineering to medical devices
EN 14136:2004	Use of external quality assessment schemes in the assessment of the performance of in vitro diagnostic examination procedures

9. Warranties

This product was designed, manufactured, tested and packaged under all applicable safety precautions. The manufacturer guarantees stability of the medical device until an expiration date at appropriate conditions of transportation, storage and intended use indicated in the Instruction for Use. The manufacturer guarantees safety of the medical device, lack of unacceptable risk for health and the environment during using of Actin FS assay if all the appropriate conditions provided by the manufacturer for this medical device is kept up.

10. Complaints

For product quality issues, please contact the authorized representative of the manufacturer:

LLC «Siemens Healthcare»

Address: 115093, Moscow, ul. Dubininskaya 96, the Russian Federation

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Dr. Schmölz, Notar

File-No. 1486 / 2019

It is officially certified hereby that **Andrea Güntner and Sabine Eucker**, Emil-von-Behring-Straße 76, 35041 Marburg - both
personally known to me - have signed this document.
Marburg, 07.08.2019

Dr. Schmölz, Notary public

A handwritten signature in blue ink, likely belonging to Dr. Anton Sebastian Schmölz, is written above the text "Dr. Schmölz, Notary public".

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Addendum to Instruction for Use
Reagents for in vitro studies of hemostasis:
Calcium Chloride Solution

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Storage at 2 to 25 °C: The expiry date is indicated on the label.

Atmospheric pressure. Humidity – 10-80%. Use by expired day.

Stability once opened:

8 weeks at 2 to 25 °C

Shipping condition: Refrigerated

Product has been shown to tolerate temperature of at least

30°C (86°F) for 36 hours, at least 35°C (95°F) for 12 hours, plus at least 45°C (113°F) for 4 hours.

Freeze/thaw cycles have NOT been shown to be tolerated.

Product has NOT been shown to tolerate freezing temperature (35°C (0°F))

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The wastes of Calcium Chloride are classified as Class B (epidemiologically safe wastes) in accordance with Sanitary and epidemiologic rules and standards SanPiN 2.1.7.2790-10 (Russian regulation).

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9. Warranties

This product was designed, manufactured, tested and packaged under all applicable safety precautions. The manufacturer guarantees stability of the medical device until an expiration date at appropriate conditions of transportation, storage and intended use indicated in the Instruction for Use. The manufacturer guarantees safety of the medical device, lack of unacceptable risk for health and the environment during using of Calcium Chloride assay if all the appropriate conditions provided by the manufacturer for this medical device is kept up.

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Shelf life - 30 months

Atmospheric pressure. Humidity – 10-80%.

Store Control Plasma N unopened at 2 to 8 °C and use by the expiry date given on the label.

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at 15 to 25 °C 4 hours

at ≤-20 °C 4 weeks

Reconstituted Control N can be frozen and thawed once. The reconstituted control plasma must be frozen as rapidly as possible in a tightly closed container. Thawing should be accomplished at 37 °C within 10 minutes. The reconstituted control plasma should not be exposed to 15 to 25 °C for longer than 2 hours after thawing.

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The generation of waste Control N should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Waste should not be disposed of untreated to the sewer unless fully compliant with the requirements of all authorities with jurisdiction.

The wastes of Control N are classified as Class B (epidemiologically dangerous wastes) in accordance with Sanitary and epidemiologic rules and standards SanPiN 2.1.7.2790-10.

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Control Plasma N contains pooled plasma collected from selected healthy blood donors.

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

Dr. Schmölz, Notar

File-No. 1736 / 2019

It is officially certified hereby that **Andrea Güntner and Dr. Sven Liebler**, Emil-von-Behring-Straße 76, 35041 Marburg - both personally known to me - have signed this document.

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
Handwritten signature of A.A. Milevskaya.

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Handwritten signature of the acting notary.

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To whom it may concern

Letter of Declaration

Herewith we, Siemens Healthcare Diagnostics Products GmbH, Marburg, Germany, in our responsibility as legal manufacturer would like to confirm that the attached documents are true and correct copies of the original documents.

Documents:

- Instruction for Use Control Plasma P
- Addendum to Instruction for Use Control Plasma P

Sincerely,

Siemens Healthcare Diagnostics Products GmbH

By: *A. V. Güntner*
Name: Andrea Güntner
Title: Senior Manager, Regulatory Affairs

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© 2008 Siemens Healthcare Diagnostics Products GmbH.

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Addendum to Instruction for Use

Reagents for in vitro studies of hemostasis:

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Shelf life - 24 months

Atmospheric pressure. Humidity – 10-80%.

Store Control Plasma P unopened at 2 to 8 °C and use by the expiry date given on the label.

The reagent should be transported in accordance with storage condition.

Stability after reconstitution:

at 15 to 25 °C 4 hours

at ≤-20 °C 4 weeks

Reconstituted Control P can be frozen and thawed once. The reconstituted control plasma must be frozen as rapidly as possible in a tightly closed container. Thawing should be accomplished at 37 °C within 10 minutes. The reconstituted control plasma should not be exposed to 15 to 25 °C for longer than 2 hours after thawing.

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The generation of waste Control P should be avoided or minimized wherever possible. Disposal of this product, solutions and any by-products should at all times comply with the requirements of environmental protection and waste disposal legislation and any regional local authority requirements. Dispose of surplus and non-recyclable products via a licensed waste disposal contractor. Waste should not be disposed of untreated to the sewer unless fully compliant with the requirements of all authorities with jurisdiction. The wastes of Control P are classified as Class B (epidemiologically dangerous wastes) in accordance with Sanitary and epidemiologic rules and standards SanPiN 2.1.7.2790-10.

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Control Plasma P contains pooled plasma collected from selected healthy blood donors.

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Summary of Reference Standards

Standard Number	Title (International Standard)
EN 13612: 2002 EN 13612: 2002/AC:2002	Performance evaluation of in vitro diagnostic medical devices
EN 13641:2002	Elimination or reduction of risk of infection related to in vitro diagnostic reagents
EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
EN ISO 15223-1:2016	Medical devices – Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements
EN ISO 17511:2003	In vitro diagnostic medical devices – Measurement of quantities in biological samples – Metrological traceability of values assigned to calibrators and control materials
EN ISO 18113-1:2011	In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 1: Terms, definitions and general requirements
EN ISO 18113-2:2011	In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 2: In vitro diagnostic reagents for professional use
EN ISO 23640:2015	In vitro diagnostic medical devices -- Evaluation of stability of in vitro diagnostic reagents (ISO
EN 62366-1:2015	Medical devices -- Application of usability engineering to medical devices
EN 14136:2004	Use of external quality assessment schemes in the assessment of the performance of in vitro diagnostic examination procedures

9. Warranties

This product was designed, manufactured, tested and packaged under all applicable safety precautions. Siemens Healthcare Diagnostics Products GmbH guarantees stability of the medical device until an expiration date at appropriate conditions of transportation, storage and intended use indicated in the Instruction for Use.

The manufacturer guarantees safety of the medical device, lack of unacceptable risk for health and the environment during using of product if all the appropriate conditions provided by the manufacturer for this medical device is kept up.

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Address: 115093, Moscow, ul. Dubinskaya 96, the Russian Federation

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
Dr. Schmölz, Notar

File-No. 1738 / 2019

It is officially certified hereby that **Andrea Güntner and Dr. Sven Liebler**, Emil-von-Behring-Straße 76, 35041 Marburg - both personally known to me - have signed this document.
Marburg, den 18.09.2019

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To whom it may concern

Letter of Declaration

Herewith we, Siemens Healthcare Diagnostics Products GmbH, Marburg, Germany, in our responsibility as legal manufacturer, would like to confirm that the attached document is a true and correct copy of the original documents.

Document:

Russian Instruction for Use of the product **Test Thrombin Reagent** with additional Addendum for the Russian Authority

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By:

Name: Andrea Güntner
Title: Senior Manager, Regulatory Affairs

Sabine Eucker
Sabine Eucker
Specialist, Regulatory Affairs

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1. Sirridge, M.S. Laboratory evaluation of hemostasis. 2nd Ed. Philadelphia: Lea & Febiger; 1974: 148.
2. Biggs, R. Human blood coagulation, haemostasis and thrombosis. 2nd Ed. Oxford: Blackwell Scientific Publications; 1976: 722.
3. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 1908781898 USA, 2008.

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Addendum to Instruction for Use
Reagents for in vitro studies of hemostasis:
Test Thrombin Reagent

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Storage conditions

Temperature 2 - 8 °C. Atmospheric pressure. Humidity 10-80%.

Stored unopened at 2 to 8 °C, Test Thrombin reagent can be used by the labeled expiry date.

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Once Opened Stability

The once opened stability is applicable for the reconstituted reagent. Confirmed Once Opened Stability:

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The wastes of Calcium Chloride are classified as Class B (epidemiologically dangerous wastes) in accordance with Sanitary and epidemiologic rules and standards SanPIN 2.1.7.2790-10 (Russian regulation).

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Summary of Reference Standards

Standard Number	Title (International Standard)
EN ISO 13485:2016	<i>Medical devices. Quality management systems. Requirements for regulatory purposes</i>
EN 13612:2002 EN 13612:2002/AC:2002	<i>Performance evaluation of in vitro diagnostic medical devices</i>
EN ISO 14971:2012	<i>Medical devices – Application of risk management to medical devices</i>
EN ISO 15223-1:2016	<i>Medical devices -- Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements</i>
EN ISO 18113-1:2011	<i>In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 1: Terms, definitions and general requirements</i>
EN ISO 18113-2:2011	<i>In vitro diagnostic medical devices – Information supplied by the manufacturer (labelling) – Part 2: In vitro diagnostic reagents for professional use</i>
EN ISO 23640:2015	<i>In vitro diagnostic medical devices -- Evaluation of stability of in vitro diagnostic reagents</i>
EN 62366-1:2015	<i>Medical devices -- Part 1: Application of usability engineering to medical devices</i>
EN 14136:2004	<i>Use of external quality assessment schemes in the assessment of the performance of in vitro diagnostic examination procedures</i>

9. Warranties

This product was designed, manufactured, tested and packaged under all applicable safety precautions. The manufacturer guarantees stability of the medical device until an expiration date at appropriate conditions of transportation, storage and intended use indicated in the Instruction for Use. The manufacturer guarantees safety of the medical device, lack of unacceptable risk for health and the environment during using of assay if all the appropriate conditions provided by the manufacturer for this medical device is kept up.

10. Complaints

For product quality issues, please contact the authorized representative of the manufacturer:
LLC «Siemens Healthcare»

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

Dr. Schmölz, Notar

File-No. 1484 / 2019

It is officially certified hereby that **Andrea Güntner and Sabine Eucker**, Emil-von-Behring-Straße 76, 35041 Marburg - both
personally known to me - have signed this document.

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