PRODUCT INFORMATION
Alcon Laboratories, Inc.

ACRY[®] IQ
ReSTOR[®]
MULTIFOCAL IOL

ENGLISH PRODUCT INFORMATION.....	2-7
DANISH PRODUKTINFORMATION.....	8-12
DUTCH PRODUCTINFORMATIE.....	13-17
FRENCH INFORMATIONS SUR LE PRODUIT.....	18-22
GERMAN GEBRAUCHSINFORMATION.....	23-27
GREEK ΠΛΗΡΟΦΟΡΙΕΣ ΓΙΑ ΤΟ ΠΡΟΪΟΝ.....	28-32
ITALIAN INFORMAZIONI SUL PRODOTTO.....	33-37
NORWEGIAN PRODUKTINFORMASJON.....	38-42
PORTUGUESE FOLHETO INFORMATIVO.....	43-47
RUSSIAN ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ.....	48-52
SPANISH INFORMACIÓN DEL PRODUCTO.....	53-57
SWEDISH PRODUKTINFORMATION.....	58-62
UKRAINIAN ІНФОРМАЦІЯ ПРО ВИРІБ.....	63-67

PRODUCT INFORMATION
Alcon Laboratories, Inc.

ACRY*Sof* IQ ReSTOR®

MULTIFOCAL IOL


**STERILE UV and Blue Light Filtering Foldable
Single-piece Apodized Diffractive Aspheric Posterior Chamber Lens**

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

DESCRIPTION

The AcrySof® IQ ReSTOR® +2.5 D Apodized Diffractive Aspheric Posterior Chamber Intraocular Lens (IOL) is an ultraviolet and blue light filtering foldable multifocal intraocular lens. The optical portion consists of a high refractive index material with a proprietary blue light filtering chromophore which filters light in a manner that approximates the human crystalline lens in the 400-475 nm blue light wavelength range (Boettner and Wolter, 1962). The optical portion is biconvex and consists of a soft acrylic material capable of being folded prior to insertion, allowing placement through an incision smaller than the optic diameter of the lens. After surgical insertion into the eye, the lens gently unfolds to restore the optical performance. The biconvex optic contains an aspheric apodized diffractive structure with a central refractive zone on the anterior surface. The apodized diffractive structure divides incoming light to provide a range of functional vision from distance to near. The anterior surface of the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL Model SV25T0 is designed with negative spherical aberration to compensate for the positive spherical aberration of the cornea*. The physical properties of these lenses are described in Figures 1-3 and Table 1.

**Figure 1: Physical Characteristics, AcrySof® IQ ReSTOR®
+2.5 D Multifocal IOL Model SV25T0
(All dimensions in millimeters)**


* The effect(s) of this aspheric design feature have not been clinically assessed.

Table 1: Physical Characteristics of AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL

	SV25T0
Optic Type	Apodized Diffractive Aspheric Optic With a Central Refractive Zone
Optic Material	Ultraviolet and blue light filtering Acrylate/Methacrylate Copolymer
UV cutoff at 10% T	See Figure 2
Index Of Refraction	1.55
Optic Powers	For available base power range see Alcon Product Guide (+ 2.5 diopters of add power for near vision)
Haptic Configuration	STABLEFORCE® Haptic
Haptic Material	See optic material
Haptic Color	Yellow
Optic Diameter (mm)	6.0
Overall Length (mm)	13.0
Haptic Angle	0°


Figure 2: Spectral Transmittance
(measured in air)

NOTE:

- Human lens data from Boettner and Wolter (1962).

Figure 3: Theoretical Percentage of Light Energy at 550 nm Wavelength

MODE OF ACTION

The AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL is intended to be positioned in the posterior chamber of the eye, replacing the natural crystalline lens. This position allows the lens to function as a refractive medium in the correction of aphakia. This IOL has a biconvex optic containing an aspheric apodized diffractive structure with a central refractive zone on the anterior surface. The apodized diffractive structure divides incoming light to provide a range of functional vision from distance to near. Alternative near add powers and energy balances provide clinicians more options to offer their patients.

INDICATIONS

The AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL is intended for primary implantation in the capsular bag of the eye for the visual correction of aphakia secondary to removal of a cataractous lens in adult patients, with and without presbyopia, who desire near, intermediate and distance vision with increased spectacle independence.

WARNINGS

1. Some visual effects may be expected due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or radial lines around point sources of light under nighttime conditions.
2. A reduction in contrast sensitivity as compared to a monofocal IOL may be experienced by some patients and may be more prevalent in low lighting conditions. Therefore, multifocal patients should exercise caution when driving at night or in poor visibility conditions.
3. The physician should consider the following points that are unique to the use of the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL:
 - The surgeon must target emmetropia to achieve optimal visual performance.
 - Patients with significant preoperative (determined by keratometry) or expected postoperative astigmatism > 1.0 D may not achieve optimal visual outcomes.
 - Care should be taken to achieve IOL centration, as lens decentration may result in a patient experiencing visual disturbances under certain lighting conditions.

PRECAUTIONS

1. Prior to surgery, prospective patients must be informed of the possible risks and benefits associated with the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL.
2. Based on experience with other multifocal IOLs, patients may experience glare, haziness, halos around lights, double vision, or blurred vision.
3. As with all multifocal IOLs, spectacle independence rates will vary. Patients may need glasses when reading small print or looking at small objects.
4. Posterior capsule opacification (PCO), when present, develops earlier into clinically significant PCO with the ReSTOR® lenses as compared to the monofocal control.
5. The safety and effectiveness of the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL have not been substantiated in patients with preexisting ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions.

Before Surgery

- Significant irregular corneal aberration
- Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this lens
- Amblyopia
- Clinically severe corneal dystrophy (e.g., Fuchs')
- Rubella, congenital, traumatic or complicated cataracts
- Extremely shallow anterior chamber, not due to swollen cataract
- Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g., iritis or uveitis)
- Aniridia
- Iris neovascularization
- Glaucoma (uncontrolled or controlled with medication)
- Microphthalmos or macrophthalmos
- Optic nerve atrophy
- Previous corneal transplant
- Pre-existing ocular conditions which may negatively impact stability of the implant
- Color vision deficiencies

Studies have shown that the color vision discrimination is not adversely affected in individuals with IOLs made from AcrySof® Natural IOL material and normal color vision. The effect on vision of the AcrySof® Natural IOL in subjects with hereditary color vision defects and acquired color vision defects secondary to ocular disease (e.g., glaucoma, diabetic retinopathy, chronic uveitis, and other retinal or optic nerve diseases) has not been studied.

During Surgery

- Mechanical or surgical manipulation required to enlarge the pupil
- Vitreous loss (significant)

- Anterior chamber bleeding (significant)
- Uncontrollable positive intraocular pressure
- Complications in which the IOL stability could be compromised

Physicians considering lens implantation in such patients should explore the use of alternative methods of aphakic correction and consider lens implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

- 6 Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from lens implantation when such conditions exist.
- 7 A high level of surgical skill is required for intraocular lens implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on intraocular lens implantation before attempting to implant intraocular lenses.
- 8 As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: lens repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
- 9 Care should be taken to remove viscoelastic from the eye at the close of surgery.
10. Do not resterilize these intraocular lenses by any method.
11. Do not store intraocular lenses at temperatures over 45°C (113°F).
12. Use only sterile intraocular irrigating solutions (such as **BSS®** or **BSS PLUS®** solution) to rinse and/or soak lenses.

CALCULATION OF LENS POWER

Accurate biometry is essential for successful visual outcomes. Preoperative calculation of required lens power for the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOL should be determined by the surgeon's experience and preference. A reference SRK/T A-Constant value for optical biometry equipment such as IOLMaster** or LenStar** is listed on the outer label. This reference A-Constant anticipates the use of both corneal power and axial length values from optical biometry equipment with standard settings for a typical patient population and a spectacle far point at 6 meters. IOL power calculation methods are often included with biometry equipment, and they are also described in the references below. In general, lens constants must be "personalized" to compensate for such things as differences in instrumentation, surgical techniques, and IOL power calculation methods that may exist between different clinical sites.

**IOLMaster is a trademark of Carl Zeiss; LenStar is a trademark of HAAG-STREIT

Hoffer, K J. The Hoffer Q formula: A comparison of theoretic and regression formulas. *J. Cataract Refract. Surg.* 19:700-712, 1993.

Holladay, J T., et al. Standardizing constants for ultrasonic biometry, keratometry, and IOL power calculations. *J. Cataract Refract. Surg.* 23:1356-1370, 1997.

Olsen, T. Calculation of intraocular lens power: A review. *Acta Ophthalmol Scand.* 85: 472-485, 2007.

Retzlaff, J A., Sanders, D.R., and Kraff, M. *Lens Implant Power Calculation*, 3rd ed. Slack, Inc., Thorofare, N.J., 1990.

<http://www.augenklinik.uni-wuerzburg.de/ulib/index.htm>

SELECTION OF ADD POWER

Selection of add power for near vision should be determined by the clinician using standard optometric practices.

DIRECTIONS FOR USE

- 1 Examine the label on the unopened package for model, powers (base and add), proper configuration, and expiration date.
- 2 After opening the cardboard storage container, verify lens case information (e.g., model, power, and serial number) is consistent with information on outer package labeling.
- 3 This device is sterile until the inner pouch is opened. Inspect the pouch carefully for tears, cuts, punctures or other signs that the pouch has been opened or damaged. DO NOT implant the IOL if the sterility has been compromised (see RETURNED GOODS POLICY).
- 4 To remove the lens, open the undamaged pouch and transfer the case to a sterile environment. Carefully open the case to expose the lens.
- 5 To minimize the occurrence of marks on the lens due to handling, all instrumentation should be scrupulously clean. Any forceps used for lens handling must have round edges and smooth surfaces.
- 6 When removing the lens from the case, DO NOT grasp the optical area with forceps. The IOL should only be handled by the haptics. Handle lenses carefully to avoid damage to lens surfaces or haptics. DO NOT attempt to reshape haptics in any way.
- 7 Rinse the lens thoroughly using sterile intraocular irrigating solution such as **BSS®** or **BSS PLUS®** solution. Prior to insertion, the lens should be carefully examined to ensure that particles have not adhered during handling.

8. Alcon recommends that the AcrySof® IQ ReSTOR® +2.5 D Multifocal IOLs are to be used with an Alcon approved delivery system.
9. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Current techniques, appropriate instrumentation, and a list of their equivalents for delivery and implantation are available from Alcon. Surgeons should verify that appropriate instrumentation is available prior to surgery.
10. DO NOT reuse this IOL. This device is for single use only.

PATIENT REGISTRATION AND REPORTING

Each patient must be registered with Alcon Laboratories, Inc. immediately following implantation of one of these lenses. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens box and mailing it to Alcon Laboratories, Inc.

Patient registration is essential for Alcon Laboratories, Inc. long-term patient follow-up program and will assist us in responding to adverse event reports.

The Patient Identification Card included in the package is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye care practitioner the patient consults in the future.

Adverse events that may reasonably be regarded as lens-related and that were not previously expected in nature, severity, or degree of incidence should be reported to Alcon Laboratories, Inc.

This information is being requested from all surgeons in order to document potential long-term effects of intraocular lens implantation.

Surgeons should use the following address and telephone number for reporting adverse events involving these intraocular lenses: Alcon Laboratories, Inc., Medical Safety (AB 2-6), 6201 South Freeway, Fort Worth, TX 76134-2099. Call Collect: (817) 551-4445.

Outside the United States, contact local Alcon offices or distributors regarding any reports of adverse events.

HOW SUPPLIED

These apodized diffractive optic posterior chamber intraocular lenses are supplied dry, in a package terminally sterilized with ethylene oxide, and must be opened only under aseptic conditions (see DIRECTIONS FOR USE).

EXPIRATION DATE

Sterility is guaranteed unless the pouch is damaged or opened. The expiration date is clearly indicated on the outside of the lens package. Any lens held after the expiration date should be returned to Alcon Laboratories, Inc. (see RETURNED GOODS POLICY).

RETURNED GOODS POLICY

Outside the United States, contact local Alcon offices or distributors regarding returned goods policy.

REFERENCES

Boettner, E. A. and Wolter, J. R. Transmission of the ocular media. *Invest. Ophthalmol.* 1:776-783, 1962.

SYMBOLS USED ON LABELING

SYMBOL	ENGLISH
IOL	Intraocular lens
PC	Posterior chamber
PCL	Posterior chamber lens
UV	Ultraviolet
D	Diopter
$\varnothing_B$	Body diameter (Optic diameter)
$\varnothing_T$	Overall diameter (Overall length)
	Do not reuse
 or 	Use by (YYYY-MM: year-month)
STERILE EO	Sterilized by ethylene oxide
SN or SN	Serial number
	Attention: See instructions for use
	Upper Limit of Temperature
	Manufacturer
EC REP	Authorized EU Representative

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- Hoffer, K.J. The Hoffer Q formula: A comparison of theoretic and regression formulas. *J. Cataract Refract. Surg.* 19:700-712, 1993.
- Holladay, J.T., et al. Standardizing constants for ultrasonic biometry, keratometry, and IOL power calculations. *J. Cataract Refract. Surg.* 23:1356-1370, 1997.
- Olsen, T. Calculation of intraocular lens power: A review. *Acta Ophthalmol Scand.* 85: 472-485, 2007.
- Retzlaff, J.A., Sanders, D.R., and Kraff, M. *Lens Implant Power Calculation*, 3rd ed. Slack, Inc., Thorofare, N.J., 1990.
- <http://www.augenklinik.uni-wuerzburg.de/ulib/index.htm>

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© 2012 Novartis
Alcon Laboratories (UK) Ltd