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Carl-Schurz-Straße 1
41453 Neuss
Germany

Postal Address:
ESPE Platz
82229 Seefeld
Germany

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To whom it may concern


August 11, 2020

Registration documents for Imprint™ Bite – IFU

We, 3M Deutschland GmbH, Carl-Schurz-Str. 1, 41453 Neuss, herewith confirm that following pages are an excerpt of the Information for Use.

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Instruction for use: Imprint™ Bite

Product Description

Imprint™ Bite is a fast setting addition reaction silicone material with high final hardness used for bite registration. The application form, cartridge, allows application inside the mouth directly from the mixing cannula, even in centric occlusion. The dimensions of the bite registration remain stable. The mixing ratio is 1:1 (by volume).

For details on all products mentioned below, please refer to the respective Instructions for Use.

☞ Instructions for Use should not be discarded for the duration of product use.

Indications

- Bite registration

Application area

- Dentistry. For professional use only.

Contraindications, adverse events and risks

Contraindications

There is no known contraindications.

Adverse events (Potential complications)

No adverse events reported. From September 2012 to September 2019, 3M ESPE has been notified of no reportable adverse event(s) for Imprint™ Bite.

No product recalls have occurred since market introduction by 3M ESPE.

MD Versions

Vinyl Polysiloxane Bite Registration Material Imprint™ Bite

- 4x Imprint™ Bite cartridges; 50 ml,
- 20 Mixing tips Green
- Instruction for use

Technical description

Imprint™ Bite cartridges: length 124 ±2 mm; width diameter 50 ±2 mm; height 28 ±2 mm; weight 108 g ±10%;

Volume – base paste/catalyst: 50ml

Mixing tips: length 88.5 ±1 mm; O.D 16 ±0.1 mm; weight 3.6 g ± 5%.

Procedure and conditions for the disposal or destruction of a medical device.

Disposal methods: Dispose of content/container in accordance with the local/regional/national/international/regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B),.

Precautionary Measures

3M MSDSs can be obtained from www.mmm.com, or contact your local subsidiary.

Preparation of the Dispenser

- Insert the cartridge into the Dispenser.
- Prior to application check that both cartridge openings are not clogged.
- When using a new cartridge, extrude a small amount of paste without mixing tip until base and catalyst paste are extruded evenly.
- Attach a Mixing Tip Green to the cartridge.
- Make sure that base paste and catalyst paste are mixed completely and are extruded in a uniform color.
- Do not extrude material that has set inside the mixing tip by force, as this may damage the cartridge and mixing tip and result in formation of leaks.

Application

Apply directly from the mixing tip onto the row of teeth, or in centric occlusion from a vestibular direction. To simplify application of the material, an intraoral tip Bite may be attached to the mixing tip.

Allow the paste to cure in centric occlusion and remove from the mouth after 1 minute.

Times	min:sec
Mixing and application	0:20
Setting time from start of mixing	1:00

The stated times are valid at 23°C/74°F.

Times depend on the ambient temperature.

Hygiene

Place the bite registration in a standard disinfection solution for the period of time recommended by the manufacturer. Subsequently, rinse under running water for approx. 15 sec.

Sterility

Not sterile.

Notes

Mouth wash solutions may impair the setting process of Imprint™ Bite, and hence, must be removed by thorough rinsing and drying. Disposable latex gloves impair the setting of silicone bite registration materials. Vinyl gloves are more suitable. Avoid contacting the material with any clothing, since the material cannot be removed by chemical cleaning.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Storage and Shelf-life

Store the product at 15-25°C/59-77°F.

Do not use after the expiration date.

Transportation conditions: Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Maintenance and repair: This medical device isn't subject to maintenance and repair.

Drugs substances: Medicinal substances: Medical devices do not contain or deliver medicinal substances.

Human Blood Derivatives: Medical devices do not contain or deliver human blood derivatives.

Animal-Derived Substances: No materials of human origin are intentionally being added to the products.

Labels and packaging symbols

Products are supplied in packaging labeled with the following information:

	Use by date		Do not reuse
	Catalogue number		Lot number
	Caution		Contained items
	conformity sign		Product meets the requirements of the applicable EC directives
	European Authorized Representative		Temperature limit (storage conditions)
	Manufacturer		Green dot
	Recycling code (Non-corrugated fiberboard)		

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

RUS: Authorised representative of the manufacturer in RF:

AO «3M Russia», 108811, Moscow, Poselenie Moskovskiy, Kievskoe shosse, 22 km, 6, bld.1

tel.: +7 495 784 74 74

Manufacturer of the medical device

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
Dental Products Carl-Schurz-Str. 1

41453 Neuss - Germany

I hereby certify that the document
affixed hereto is a true copy
of the original.

Starnberg, August 12, 2020

Carina Kellner, Notarassessorin
(appointed as official substitute of
Michael Volmer,
Notary in Starnberg, Bavaria, Germany)

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