MEDpro medical

Approved

G.J. ten Brinke
Director

MEDpro medical B.V.
Landjuweel 11, 3905 PE Veenendaal
The Netherlands
T: +31318769080 F: +31318769107
www.medpromedical.nl

Instructions for Use

Urethral Stents (single and in set)

MEDpro medical b.v. • Office address: Landjuweel 11, 3905 PE Veenendaal
• Mail address: Handellaan 7, 3906 BX Veenendaal
☎ +31 318769080 • 📠 +31318769107 • @ info@imedpromedical.nl • KvK: 51758644
www.medpromedical.nl

Instruction for Use

Double J Stents

Device Description

- The Double J Stents are provided for use in Urology procedures. The configurations available include:

Family	Size (x0,1Fr.)	Length (cm)	Type	Color/Material	Options
Both Ends Single Loop (DJSBES) One End Multiloop (DJSOEM) Both Ends Multiloop (DJSBEM)	30	12	One End Open (OEO)	Blue (*B) Green (*G) White (*W) Long Duro (*LD)	With Clamp (/C)
	35				With Guide wire (/G)
	38				With Clamp + Guide wire (/CG)
	40				With Suture (/S)
	45				With Suture + Clamp (/SC)
	47				With Suture + Clamp + Guide Wire (/SCG)
	48	24	Both Ends Open (BEO)		Non standard pusher with Length xx cm (\Pxx)
	50				
	55				
	60				
	65				
	70				
	80				

Indications

- The Double J Stents are indicated for use in Urology procedures to ensure the patency of a Ureter.

Contraindications

- No Specific Contraindications.

Precautions, Warnings & Adverse effects

Carefully read all instructions for use. The device shall only be applied for its intended use and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.

- Caution: Each physician is responsible for using the appropriate technique and deciding on the indication for use of this device based on his own experience, training and medical judgment. The doctor must be trained in the proper use of the device.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
- Warning: Possible complications associated with this device are related to insertion and use techniques and may include perforation and infection.
- Warning: All components of the "Double J Stents" are for single use only.
- Caution: Do not use the device if there is any indication that the sterility of the device has been compromised.
- Adverse effects: Possible adverse effects of guide wire placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.

Condition of Use

- This device should be administered to humans only by physicians authorized by the Competent Authority of the country in which the physician is practicing.
- This device is a one- time use only. Re-use may result in clinical complications.

Shelf-life of the Product

- The shelf life of the Double J Stents is 3 years (Indicated on product label with the following use-by symbol and date).

Replacement of damaged or unsuitable components

- If any component is found damaged or unsuitable for use, do not use. Contact "Distributor" for replacement.

- The Double J Stents have been sterilized by Exposure to Ethylene Oxide. Sterility indicators are on each package. The imprinted label (20 x 6 mm or ¾ inch diameter circles) will change color after ethylene oxide exposure. Do not use the product unless the sterility indicators are the correct color.
- If the sterile package is damaged or possibly opened, do not use. Call "Distributor" and we will replace the product.
- The "Double J Stents" are packed single. Product in each pack must be utilized immediately when opened. Product should not be re-sterilized.

Directions for Use

- The involved renal collecting system should be visualized via intravenous or retrograde pyelography.
- Select the appropriate size Stent for patient's anatomy. A well sized Stent should be fully coiled within the renal pelvis with the loops free floating in the bladder.
- Confirm the information on the label and that the product has not reached its expiration date and there is no damage to the packaging or device.
- The retrieval line can be removed prior to placement if desired or the retrieval line can be left intact to extend externally.
- Insert flexible end of guide wire into the renal pelvis for DJS Both Ends Open or the rigid end of guide wire for DJS One End Open, using clinically. For DJS Both Ends Closed use rigid end of Guide wire through the side hole of the stent using appropriate retrograde technique.
- Pass tapered tip of Stent over guide wire and through cystoscope.
- Advance Stent up the ureter.
- Withdraw guide wire slowly to allow distal coil to form in the renal pelvis. Placement should be confirmed under fluoroscopy or x-ray.
- Monitor the Stent as required. Retrieval line indwelling time should not exceed THIRTY (30) days for blue, green and white DJS to avoid possible cord encrustation.
- For Long Duro DJS (yellow) retrieval line indwelling time should not exceed NINETY (90) days it is recommended that when long term use is indicated the Stent should be evaluated every 30 days.

URETRIC STENT (DJ Stent) REMOVAL INSTRUCTIONS

- Retrieve cystoscopically by gently pulling on Stent or retrieval line with grasping forceps or equivalent. If resistance is encountered during removal of the Stent, stop and determine cause of resistance before proceeding.

Explanation of symbols used on label

Symbol	Description	Symbol	Description		
REF	Catalogue reference	CE 2265	CE Mark		Do not Re sterilize
	Name & address of manufacturer		Date of manufacture		Keep Dry
	Lot Number		Use By		Keep Away from Sunlight
	EO sterilized		Read The Documentation		Temperature Limitations
	Warning		Not for second use		

MEDpromedical BV
Landjuweel 11
3905PE Veenendaal
The Netherlands
info@medpromedical.nl

MEDpro medical B.V.
Landjuweel 11, 3905PE Veenendaal
The Netherlands
T: +31318769090 F: +31318769107
www.medpromedical.nl

Instruction for Use

Loop Stents

Device Description

- The Loop Stents are provided for use in Urology procedures. The configurations available include:

Family	Size (Fr.)	Length	Type
LS	4	24	One End Open (OEO) Both Ends Open (BEO)
	5	26	
	6	28	
		30	

Indications

- Loop stents are used to ensure the patency of a ureter.

Contraindications

- No known Contraindications.

Precautions, Warnings & Complications

- Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.
- Caution: All medical staff is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
- Warning: All components of the Loop Stents are for single use only.
- Caution: Do not use the device if there is any indication that the sterility of the device has been compromised.
- Adverse effects: Possible adverse effects of Loop Stents placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.

Condition of Use

- This device should be administered to humans only by physicians authorized by the Competent Authority of the country in which the physician is practicing.
- This device is a one-time use only. Re-use may result in clinical complications.

Shelf-life of the Product

- The shelf life of the Loop Stents is 3 years (Indicated on product label with the use-by symbol and date). Store the Loop Stents at a dark, dry, cool place. Avoid extended exposure to light.

Sterility – This product is Sterile unless the package has been opened or damaged

- The Loop Stents have been sterilized by Exposure to Ethylene Oxide. Sterility indicators are on each package. The imprinted label (20 x 6 mm or ¾ inch diameter circles) will change colour after ethylene oxide exposure. Do not use the product unless the sterility indicators are the correct colour.
- If the sterile package is damaged or possibly opened, do not use. Call “Distributor” and we will replace the product.
- The Loop Stents are packed single. Product in each pack must be utilized immediately when opened. Product should not be re-sterilized.

Replacement of damaged or unsuitable components

- If any component is found damaged or unsuitable for use, do not use. Contact "Distributor" for replacement.

Directions for Use

- Suggested instructions for using Loop Stents:
 - Carefully read all instructions prior to using this device.
 - The involved renal collecting system should be visualized via intravenous or retrograde pyelography.
 - Select the appropriate size stent for patient's anatomy. A well sized stent should be fully coiled within the renal pelvis with the loops free floating in the bladder.
 - Insert flexible end of guide wire into the renal pelvis using clinically appropriate retrograde technique.
 - Pass tapered tip of the stent over guide wire and through cystoscope.
 - Load the loop positioner over the guide wire proximal to the stent and advance stent up the ureter. The black marker on the Loop stent, located just above the bladder loops, should be positioned at the ureteral orifice.
 - Withdraw guide wire slowly to allow distal coil to form in the renal pelvis. Placement should be confirmed under fluoroscopy or X-ray.

Explanation of symbols used on label

Symbol	Description	Symbol	Description
REF	Catalogue reference		CE Mark
	Name & address of manufacturer		Date of manufacture
	Lot Number		Use By
	Eto sterilized		Read The Documentation
	Warning		Not for second use
	Do not Re sterilize		Keep Dry
	Keep Away from Sunlight		Temperature Limitations

MEDpromedical BV
Landjuweel 11
3905PE Veenendaal
The Netherlands
info@medpromedical.nl

MEDpro medical B.V.
Landjuweel 11, 3905PE Veenendaal
The Netherlands
T: +31318769080 F: +31318769107
www.medpromedical.nl

Instruction for Use

Tumor Stents

Device Description

- The Tumor Stents are provided for use in Urology procedures. The configurations available include:

Family	Size (Fr.)	Length (cm)	Type
TS	5/10	24	Both Ends Open (BEO)
	6/12	26	
	7/14	28	
		30	

Indications

- Tumor Stents optimize drainage, even in cases where the ureter is extremely compressed. The soft material of the pigtail section offers high patient comfort, while the special stiffened shaft portion helps to prevent the stent from collapsing. The tumor stent side holes are strategically placed to allow drainage and to prevent the ingrowth of tissue into the stent.

Contraindications

- No Specific Contraindications.

Precautions, Warnings& Complications

- Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.
- Caution: All medical staff is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
- Warning: All components of the Tumor Stents are for single use only.
- Caution: Do not use the device if there is any indication that the sterility of the device has been compromised.
- Caution: Stent must not remain indwelling more than eight (8) weeks. If the patient's status permits, the stent may be replaced with a new stent. These stents are not intended as permanent indwelling devices. Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered. A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements. Individual variations of interaction between stents and the urinary system are unpredictable. Periodic evaluation via cystoscope or radiographic means suggested. The stent must be replace if encrustation hampers drainage
- Adverse effects: Possible adverse effects of Tumor Stents placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.

Condition of Use

- This device should be administered to humans only by physicians authorized by the Competent Authority of the country in which the physician is practicing.
- This device is a one- time use only. Re-use may result in clinical complications.

Shelf-life of the Product

- The shelf life of the Tumor Stents is 3 years (Indicated on product label with the following use-by symbol and date). Store the Tumor Stents at a dark, dry, cool place. Avoid extended exposure to light.

- Sterility – This product is sterile unless the package is damaged.**
- The Tumor Stents have been sterilized by Exposure to Ethylene Oxide. Sterility indicators are on each package. The imprinted label (20 x 6 mm or 3/8 inch diameter circles) will change colour after ethylene oxide exposure. Do not use the product unless the sterility indicators are the correct colour.
 - If the sterile package is damaged or possibly opened, do not use. Call "Distributor" and we will replace the product.
 - The Tumor Stents are packed single. Product in each pack must be utilized immediately when opened. Product should not be re-sterilized.

Replacement of damaged or unsuitable components

- If any component is found damaged or unsuitable for use, do not use. Contact "Distributor" for replacement.

Directions for Use

- Suggested instructions for using Tumor Stents:
Place the stent over a previously placed 0.038 inch diameter guide wire. NOTE:- Guide wire should be coiled in the urinary bladder.
Advance the stent into the collecting system to the bladder. Check position of the larger diameter ureteral portion of the stent. Assure that it is placed within the ureteropelvic junction (UPJ), to prevent post operative structuring.
Remove the guide wire slowly. Check for at least a full coil in the bladder and that the larger stent section is located in the ureteropelvic junction. This should be confirmed fluoroscopically.

Explanation of symbols used on label

Symbol	Description	Symbol	Description
REF	Catalogue reference		CE Mark
	Name & address of manufacturer		Date of manufacture
	Lot Number		Use By
	EIO sterilized		Read The Documentation
	Warning		Not for second use
	Do not Re sterilize		Keep Dry
	Keep Away from Sunlight		Temperature Limitations

MEDpromedical BV
Landjuweel 11
3905PE Veenendaal
The Netherlands
info@medpromedical.nl

MEDpro medical B.V.
Landjuweel 11, 3905PE Veenendaal
The Netherlands
T: +31318769080 F: +31318769107
www.medpromedical.nl

Instruction for Use Endopyelotomy Stents

Device Description

- The Endopyelotomy Stents are provided for use in Urology procedures. The configurations available include:

Family	Size (Fr.)	Length (cm)	Type
ES	5/10	24	Both Ends Open (BEO)
	6/12	26	
	7/14	28	
		30	

Indications

- Endopyelotomy Stents are used for stenting the ureteropelvic junction and the ureter following a percutaneous endopyelotomy while providing both nephrostomy and ureteral drainage.

Contraindications

- No known Contraindications.

Precautions, Warnings & Complications

- Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.
- Caution: All medical staff is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
- Warning: All components of the Endopyelotomy Stents are for single use only.
- Caution: Do not use the device if there is any indication that the sterility of the device has been compromised.
- Caution: Stent must not remain indwelling more than eight (8) weeks. If the patient's status permits, the stent may be replaced with a new stent. These stents are not intended as permanent indwelling devices.
- Caution: Do not force components during removal or replacement. Carefully remove the components if any resistance is encountered.
- Caution: A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Caution: Individual variations of interaction between stents and the urinary system are unpredictable. Periodic evaluation via cystoscope or radiographic means suggested. The stent must be replace if encrustation hampers drainage
- Adverse effects: Possible adverse effects of Endopyelotomy Stents placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.

Condition of Use

- This device should be administered to humans only by physicians authorized by the Competent Authority of the country in which the physician is practicing.
- This device is a one-time use only. Re-use may result in clinical complications.

Shelf-life of the Product

- The shelf life of the Endopyelotomy Stents is 3 years (Indicated on product label with the use-by symbol and date).

Store the Endopyelotomy Stents at a dark, dry, cool place. Avoid extended exposure to light.

Sterility – This product is Sterile unless the package has been opened or damaged

- The Endopyelotomy Stents have been sterilized by Exposure to Ethylene Oxide. Sterility indicators are on each package. The imprinted label (20 x 6 mm or ½ inch diameter circles) will change colour after ethylene oxide exposure. Do not use the product unless the sterility indicators are the correct colour.
- If the sterile package is damaged or possibly opened, do not use. Call "Distributor" and we will replace the product.
- The Endopyelotomy Stents are packed single. Product in each pack must be utilized immediately when opened. Product should not be re-sterilized.

Replacement of damaged or unsuitable components

- If any component is found damaged or unsuitable for use, do not use. Contact "Distributor" for replacement.

Directions for Use

- Suggested instructions for using Endopyelotomy Stents :
 - Place the stent over a previously placed 0.038 inch diameter guide wire. NOTE:- Guide wire should be coiled in the urinary bladder.
 - Advance the stent into the collecting system to the bladder. Check position of the larger diameter ureteral portion of the stent. Assure that it is placed within the ureteropelvic junction (UPJ), to prevent post operative structuring.
 - Remove the guide wire slowly. Check for at least a full coil in the bladder and that the larger stent section is located in the uretero pelvic junction. This should be confirmed fluorooscopically.

Explanation of symbols used on label

Symbol	Description	Symbol	Description
REF	Catalogue reference		CE Mark
	Name & address of manufacturer		Date of manufacture
	Lot Number		Use By
	EIO sterilized		Read The Documentation
	Warning		Not for second use
	Do not Re sterilize		Keep Dry
	Keep Away from Sunlight		Temperature Limitations

MEDpromedical BV
Landjuweel 11
3905PE Veenendaal
The Netherlands
info@medpromedical.nl

MEDpro medical B.V.
Landjuweel 11, 3905PE Veenendaal
The Netherlands
T: +31318769080 F: +31318769107
www.medpromedical.nl

The documents contains 5 (including cover) pages, bound.
G.J. ten Brinke
Director

MEDpro medical B.V.
Landjuweel 11, 3905PE Veenendaal
The Netherlands
+31318769080 F: +31318769107
www.medpromedical.nl

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