Instruction for Use Suprapubic Catheters

Device Description

- The Suprapubic Catheters are provided for use in Urology procedures. Product consists of Malecot Catheter, matching hollow needle obturator and Urine bag connector. The configurations available include:

Family	Size (Fr.)	Length (cm)	Color/Material
SPC	10	25	*PU
	12		
	14		
	16		
	18		

Indications

- Suprapubic Catheters are used to provide bladder drainage by percutaneous placement of a Malecot catheter.

Contraindications

- No known Contraindications.

Precautions, Warnings & Adverse effects

- Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.
- Caution: All medical staff is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.
- Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.
- Warning: All components of the Suprapubic Catheters are for single use only.
- Caution: The indwelling time should not exceed four (4) weeks. Evaluation and replacement beyond four (4) weeks must be considered.
- Caution: Do not use the device if there is any indication that the sterility of the device has been compromised.
- Adverse effects: Possible adverse effects of Suprapubic Catheters placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.

Condition of Use

- This device should be administered to humans only by physicians authorized by the Competent Authority of the country in which the physician is practicing.
- This device is a one-time use only. Re-use may result in clinical complications.

Shelf-life of the Product

- The shelf life of the Suprapubic Catheters is 3 years (Indicated on product label with the following use-by symbol and date). Store the Suprapubic Catheters at a dark, dry, cool place. Avoid extended exposure to light.

Sterility – This product is Sterile unless the package has been opened or damaged

- The Suprapubic Catheters have been sterilized by Exposure to Ethylene Oxide. Sterility indicators are on each package. The imprinted label (20 x 6 mm or ½ inch diameter circles) will change colour after ethylene oxide exposure. Do not use the product unless the sterility indicators are the correct colour.
- If the sterile package is damaged or possibly opened, do not use. Call "Distributor" and we will replace the product.
- The Suprapubic Catheters are packed single. Product in each pack must be utilized immediately when opened. Product should not be re-sterilized.

Replacement of damaged or unsuitable components

- If any component is found damaged or unsuitable for use, do not use. Contact "Distributor" for replacement.

Directions for Use

- Suggested instructions for using Suprapubic Catheters:
 - Straighten the malecot catheter and guide the needle obturator into the catheter tip. Secure its position with the luer lock to close the malecot wings
 - Fill the bladder.
 - Make a small 2mm incision in the skin approximately 1-2 inches above the symphysis pubis (below the palpable bladder dome)
 - Insert the catheter vertically into the bladder on patient's without previous Pelvic surgery or at a 30 degree angle toward the symphysis pubis on patients with previous pelvic surgery. Gradually advance the catheter 4-5cm beyond the first sign of bladder fluid. The catheter midpoint is an approximation.
 - If bladder pressure is low and urine fails to flow spontaneously through the hollow obturator, withdraw urine to ascertain that the catheter is within the bladder.
 - After releasing the luer lock in order to open the malecot wings, carefully remove the hollow needle obturator. The catheter must be stabilized and the needle rotated while pulling or removing. Care must be exercised to avoid jerking the catheter forward.
 - Use the urine bag connector to connect the malecot catheter to a drainage bag.
 - Allow the bladder to empty.
 - Withdraw the malecot catheter until the malecot wings meet the resistance of the bladder dome. Advance the catheter approximately 2cm back into the bladder to allow for movement.

Explanation of symbols used on label

Symbol	Description	Symbol	Description
REF	Catalogue reference	CE 2265	CE Mark
	Name & address of manufacturer		Date of manufacture
	Lot Number		Use By
	Eto sterilized		Read The Documentation
	Warning		Not for second use
	Do not Re sterilize		Keep Dry
	Keep Away from Sunlight		Temperature Limitation

MEDpromedical BV
Landjuweel 11
3905PE Veenendaal
The Netherlands
info@medpromedical.nl

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STATEMENT OF AUTHENTICITY

The identity of the person stated below:

Mr. G.J. ten Brinke

Passport

number BM128FRL9

issued in Veenendaal, on 28th April 2014 ,

has been established by me, *mr.* Johan van Heteren, notary awaiting appointment, as a substitute for *mr.* Roland Ottens , civil-law notary practising in Veenendaal, on the basis of the identity document referred to above.,

The aforesaid person has placed his/her signature in my presence.

The significance of this authentication is strictly limited to ascertaining the identity and confirming the signature of the aforesaid person. Consequently I, the civil-law notary, do not give an opinion on the content and its possible legal effect or other aspects of the documents attached.

This statement of authenticity does not constitute an authorization of the aforesaid person to perform a legal act on behalf of another natural and/or legal person.

Signed in Veenendaal (The Netherlands) on 18th of September 2015

A handwritten signature in blue ink, appearing to be "J. van Heteren", is written over a large, stylized blue 'X' mark.

APOSTILLE

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1 Country: THE NETHERLANDS

This public document

2 has been signed by:

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3 acting in the capacity of:

candidate notary

4 bears the seal/stamp of

Mr R Ottens

Notaris te Veenendaal

Certified

5 in Utrecht

6 on 18-9-2015

7 by the court registrar

8 no 15/2363-2

9 Seal/stamp

10 Signature:

J. van Berkel

MEDpro medical

To the Federal Service on Surveillance
in Healthcare of Russian Federation

Veenendaal 16-09-2015

Explanation letter

Hereby we, the company "MEDpro medical b.v.", Landjuweel 11, 3905PE Veenendaal, the Netherlands, would like to provide the clarification for the registration dossier "Suprapubic catheter, length 25 cm, size 10 Fr; 12 Fr; 14 Fr; 16 Fr; 18 Fr" (file No. 12954 from 26.03.2015) in accordance with the Inquiry under No. 04-19847/15 from 13.07.2015 on preparation of information and data:

- The absence of information intended to familiarize the consumer with the construction of medical devices (drawings, chart)

This medical device (catheter) is a simple device. All medical staff is familiar with this kind of device and is trained to use it. There is no need to indicate drawings in the IFU.

The product is visible in pouch.

- The ensured by manufacturer values of basic parameters, performances of medical device are presented not in full.

The IFU refers to several product codes and it is impossible to insert these data within the IFU.

Moreover the absence of these does not affect features and characteristics having an impact on the quality, efficiency and safety of medical device.

The IFU has been prepared according to the requirements of the Directive 93/42/EEC and to the applicable harmonized standards.

- The diameter of catheter should be measured in mm.

The diameter of catheter is measured in French, because it is international unit of measurement for catheters, sheath, etc. It is common used measurements in Europe, USA and other countries and it is widely known for specialist in medicine. The measurement in mm we indicated in technical file.

- Transport conditions:

The transport and storage conditions are indicated on product label in form of symbols as follows «Keep dry», «Keep away from sunlight», «Temperature limitation 12°C - 35°C». These conditions are considered as normal. Used symbols are common and meet the requirements of standard EN 1041 «Information supplied by the manufacturer with medical devices». In the IFU we explain symbol used on label.

According to the EN 1041 (Annex A, item 13.3 i) the transport and storage conditions should be indicated in the IFU only in particular case if they differ from normal conditions or if they are critical for the safety use of product. In other cases it is not obligatory to indicate the transport and storage conditions in the IFU. Labeling is a part of operational documentation and it is attached to this letter.

MEDpro medical

– **Manufacturer's warranty:**

Suprapubic catheters are supplied sterile and intended for one-time use only, for that reason the manufacturer's warranty is not applied as the manufacturer indicates the shelf-life.

The shelf-life is indicated on product label according MDD 93/42/EEC. Labeling is a part of operational documentation and it is attached to this letter. For that reason we don't indicate this information in IFU.

– **Disposal conditions:** Dispose the used catheters in strict accordance with official and institutional standards for medical waste disposal.

Disposal conditions differ in the various countries, for that reason we don't include this information in IFU.

Substance hazard category according to Russia relevant regulation: Category B- Epidemiological hazardous waste.

– **The materials (marks of materials) used in medical product are not indicated.**

The material of catheter (PU) is indicated on label and in IFU. Labeling is a part of operational documentation and it is attached to this letter. The main material for the catheter is indicated as the catheter is the part that stays in the body for the indicated use. The other items are used as an aid but do not touch or stay in the body.

– **The possible adverse effects are not indicated.**

The IFU include following information: "Possible adverse effects of Catheters placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures." Doctors are familiar already with these effects as they are documented in literature.

It can be added the common adverse effects noticed are: Lower back pain, chills, fever

– **One can see in the item 3 of technical file "The Catheter is intended for temporary drainage of urine at different levels (kidney, bladder)", but the IFU doesn't include the instruction for use at level of kidney).**

The Catheter is intended for temporary drainage of urine at level of bladder. Item 3 of technical file has misprinted intended use. Please find updated technical file (version 2).

– **There are no removal instructions of medical device.**

Gently pull the catheter out by holding the hub.



There are no any specific removal instructions. All medical staff is familiar with this kind of device and is trained to use it.

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