STATE OF NORTH CAROLINA

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Department of The
Secretary of State

APOSTILLE

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1. COUNTRY: **UNITED STATES OF AMERICA**
2. This Public Document
has been signed by MELISSA BOTTOMS
3. acting in the capacity of NOTARY PUBLIC
COMMISSION EXPIRES 1/27/2024
4. bears the seal/stamp of COUNTY OF STOKES, NORTH CAROLINA

CERTIFIED

5. at Raleigh, North Carolina
6. the 28TH DAY OF MAY, 2020
7. by Secretary of State or Deputy Secretary of State, State of North Carolina
8. No. 6813
9. Seal/Stamp
10. Signature

Scan to verify online.

Claine J. Marshall
Claine J. Marshall

Certification# 4564391135

Verify the issuance of this Apostille at https://www.sosnc.gov/online_services/search/by_title_verification

Secretary of State

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document and, where appropriate, the status of the state or states where the public document bears. This Apostille does not certify the content of the document for which it was issued. This Apostille is not valid for use anywhere within the United States of America, its territories or possessions. (Verify the issuance of this Apostille at <https://www.sosnc.gov/apostille>)

INSTRUCTIONS FOR USE

ENDOSCOPIC BIOPSY NEEDLES

INTENDED USE

The product is used to obtain tissue samples for histopathological examination during various endoscopic procedures.

Howell Biliary Aspiration Needle HBAN-1-22

INTENDED USE

This device is used for aspiration biopsy in the biliary system.

NOTES

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

WARNINGS

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

CONTRAINDICATIONS

Those specific to ERCP.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with ERCP include, but are not limited to: pancreatitis, cholangitis, aspiration, perforation, hemorrhage, infection, sepsis, allergic reaction to contrast or medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be fully retracted into sheath prior to introduction into, advancement through or withdrawal from endoscope. Failure to retract needle may result in damage to endoscope.

INSTRUCTIONS FOR USE

1. Uncoil and straighten catheter completely.

2. Prior to inserting device into endoscope, retract needle completely into sheath. Ensure needle fully retracts into sheath and that no sharp tip extends beyond ball tip.
 3. Introduce device into accessory channel of endoscope and advance in short increments to desired biopsy site. Do not straighten pre-bend more than is necessary for insertion down channel. **Note:** A large-bore endoscope is preferred to minimize straightening.
 4. Position endoscope tip away from and slightly above ampulla.
 5. Place ball tip, with needle retracted, on lower lip of biliary orifice. Advance catheter into cannulating position. Proceed with cannulation of biliary system.
 6. Once cannulation is achieved, advance ball tip to appropriate position for aspiration biopsy.
 7. Carefully advance needle out of sheath. With needle in desired location, remove stylet from Luer lock fitting on handle.
 8. Attach a 20 cc syringe to Luer lock fitting, then aspirate. Continue aspiration and gently move needle in small increments back and forth within biopsy site. **Note:** Several biopsy locations may be chosen and aspiration repeated without withdrawing needle from bile duct. Between each placement of needle for biopsy, suction should be released.
 9. When aspiration is complete, withdraw needle from biopsy site and into ball tip by pulling back on Luer slip fitting, then withdraw device from endoscope.
 10. Prepare specimen for examination per institutional guidelines.
- Upon completion of procedure, dispose of device(s) per institutional guidelines for biohazardous medical waste.

EchoTip ProCore High Definition Ultrasound Biopsy Needle ECHO-HD-19-C with 10ml syringe; EchoTip ProCore High Definition Ultrasound Biopsy Needle ECHO-HD-22-C with 10ml syringe; EchoTip ProCore High Definition Ultrasound Biopsy Needle ECHO-HD-25-C with 10ml syringe; EchoTip ProCore HD Ultrasound Biopsy Needle ECHO-HD-3-20-C with 10ml syringe

IMPORTANT INFORMATION

Please review prior to use.

Intended Use

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal lesions, mediastinal masses, lymph nodes and intraperitoneal masses within or adjacent to the gastrointestinal tract.

NOTES

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site. Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage and acute pancreatitis. Those associated with EUS needle biopsy include but are not limited to: pain, death, peritonitis, portal vein gas and thrombosis, pneumoperitoneum and tumor seeding of the needle tract.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the biopsy site.

When targeting multiple sites, replace device for each site.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)

2. Prepare syringe as follows:

a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.

b. Turn stopcock 90° to "closed" position.

c. Pull plunger back until it is locked into place at desired setting, creating suction.

d. Set prepared syringe aside until biopsy is desired.

3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE


1. Identify desired biopsy site by endoscopic ultrasound.

2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)

3. Attach device to accessory channel port by rotating device handle until fittings are connected.

4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.

5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
 6. Advance needle into lesion.
 7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle.
 8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.
 9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to facilitate biopsy. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during biopsy.
 10. Upon completion of biopsy, turn stopcock to "close" position.
 11. Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
 12. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.
 13. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.
 14. Push in plunger to expel specimen, then prepare per institutional guidelines.
 15. Remaining material may be recovered for examination by reinserting stylet or flushing device.
- Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

EchoTip Ultra Endoscopic Ultrasound Needle ECHO-19 with 10ml syringe; EchoTip Ultra Endoscopic Ultrasound Needle ECHO-25 with 10ml syringe; EchoTip Ultra Endoscopic Ultrasound Needle ECHO-1-22 with 10ml syringe

INTENDED USE

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

NOTES

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Do not use this device for any purpose other than stated intended use.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer slip fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)

2. Prepare syringe as follows:

- a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
- b. Turn stopcock 90° to "closed" position.
- c. Pull plunger back until it is locked into place at desired setting, creating a suction.
- d. Set prepared syringe aside until aspiration biopsy is desired.


3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. **Note:** Reference mark


for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)
3. Attach device to accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device has been attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion.
7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle. Preserve stylet for use if additional cell collection is desired.
8. Attach Luer lock fitting of previously prepared syringe securely onto fitting on handle.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection.
10. Upon completion of cell collection, retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
11. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.
12. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe and attach new syringe to Luer lock fitting on handle.
13. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.
14. Aspirate remaining can be recovered for examination by flushing device.
15. For additional cell collection from same lesion, gently reinsert stylet into metal fitting on handle. **Note:** Prior to reinserting stylet, wipe with saline or sterile water. While supporting sheath at Luer lock fitting, advance stylet in small increments until stylet hub is engaged in fitting.
16. Additional samples may be obtained by repeating Step 2 of "System Preparation" through Step 14 of the "Instructions for Use."

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

EchoTip Ultra Endoscopic Ultrasound Needle ECHO-3-22 with 10ml syringe

INTENDED USE

This device is used to sample targeted submucosal gastrointestinal lesions through the accessory channel of an ultrasound endoscope.

This device is indicated for adult use only.

NOTES

Do not use this device for any purpose other than stated intended use.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

Please read all instructions before using this device.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy.

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

WARNINGS

This device is supplied sterile and intended for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.

Not for use in the heart or vascular system.

Not for use with a wire guide. Use with a wire guide may result in damage to the wire guide or stripping of the wire guide.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

When targeting multiple sites, replace device for each site.

Needle must be fully retracted into sheath and thumbscrew on safety ring must be locked at 0 cm mark to hold needle in place prior to introduction, advancement or withdrawal of device.

Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)

2. Prepare syringe as follows:

a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.

b. Turn stopcock 90° to "closed" position.

c. Pull plunger back until it is locked into place at desired setting, creating a suction.

d. Set prepared syringe aside until aspiration biopsy is desired.

3. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained. Note: Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy site by endoscopic ultrasound.

2. With needle retracted into sheath and safety ring locked at 0 cm mark, introduce ultrasound needle into accessory channel of endoscope and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)

3. Attach device to accessory channel port by rotating device handle until fittings are connected.


4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in the window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device has been attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.

5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.

6. Advance needle into lesion. **Caution:** When advancing needle, excessive force may result in unintended needle extension

7. Remove stylet from needle by gently pulling back on plastic hub seated in metal fitting of needle handle. Preserve stylet for use if additional cell collection is desired.
 8. Attach Luer lock fitting of previously prepared syringe securely onto metal fitting of handle.
 9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments back and forth within biopsy site.
- Note:** Do not remove needle from biopsy site during cell collection.
10. Upon completion of cell collection, retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
 11. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope.
 12. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe and attach new syringe to Luer lock fitting on handle.
 13. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.
 14. Aspirate remaining can be recovered for examination by flushing device.
 15. For additional cell collection from same lesion, gently reinsert stylet into metal fitting of handle. Note: Prior to reinserting stylet, wipe with saline or sterile water. While supporting sheath at Luer lock fitting, advance stylet in small increments until stylet hub is engaged in fitting.
 16. Additional samples may be obtained by repeating Step 2 of "System Preparation" through Step 14 of the "Instructions for Use."

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

EchoTip Ultra High Definition Ultrasound Access Needle ECHO-HD-19-A with 10ml syringe

INTENDED USE

This device is used to access or sample submucosal and extramural lesions of the gastrointestinal tract, access of the following: the intra- or extra-hepatic bile ducts, pancreatic ducts, cystic duct, gallbladder or for delivery of injectable materials into tissues through the accessory channel of an ultrasound endoscope.

NOTES

Do not use this device for any purpose other than stated intended use.

This needle is intended to be used with legally marketed drugs and devices.

If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify COOK Medical for return authorization.

This device can be used together with all ECHO compatible products.

Aspiration syringe not intended to be used for injection.

Store in a dry location, away from temperature extremes.

Use of this device restricted to a trained healthcare professional.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and /or transmission of disease.

CONTRAINDICATIONS

Those specific to primary endoscopic procedure to be performed in gaining access to desired site.

Coagulopathy

POTENTIAL COMPLICATIONS

Those associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest, damage to blood vessels, nerve damage and acute pancreatitis.

WARNING

Not for use in the heart or vascular system.

PRECAUTIONS

When performing fine needle injection (FNI), ensure that injectable material does not contact untargeted site.

When targeting multiple sites, replace device for each site.

Refer to package label for minimum channel size required for this device.

Needle must be retracted into sheath and thumbscrew on safety ring must be locked to hold needle in place prior to introduction, advancement or withdrawal of device. Failure to retract needle may result in damage to endoscope.

SYSTEM PREPARATIONS

1. For FNI and/or Access, skip to step 4.

2. For fine needle aspiration (FNA), examine syringe. It has two plunger locks that must be depressed to advance plunger. Tip of syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (See fig. 1)


3. Prepare syringe as follows:

- a. With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b. Turn stopcock 90° to "closed" position.
 - c. Pull plunger back until it is locked into place at desired setting, creating suction.
 - d. Set prepared syringe aside if aspiration biopsy is desired.
4. Advance device into ultrasound endoscope to determine preferred sheath length. To adjust length, loosen thumbscrew lock on sliding sheath adjuster and slide until preferred length is attained.
- Note:** Reference mark for sheath length will appear in sliding sheath adjuster window. (See fig. 2) Tighten thumbscrew on sliding sheath adjuster to maintain preferred sheath length.

INSTRUCTIONS FOR USE

1. Identify desired biopsy/injection site by endoscopic ultrasound.
2. Introduce ultrasound needle into accessory channel and advance device in small increments until Luer lock fitting at base of sliding sheath adjuster meets fitting on accessory channel. (See fig. 3)
3. Attach device to endoscope accessory channel port by rotating device handle until fittings are connected.
4. With ultrasound endoscope and device straight, adjust needle to desired length by loosening thumbscrew on safety ring, and advancing it until desired reference mark for needle advancement appears in window of safety ring. (See fig. 4) Tighten thumbscrew to lock safety ring in place. **Note:** Number in safety lock ring window indicates extension of needle in centimeters. **Caution:** During needle adjustment or extension, ensure device is attached to accessory channel. Failure to attach device prior to needle adjustment or extension may result in damage to endoscope.
5. While maintaining position of ultrasound endoscope, extend needle by advancing needle handle to pre-positioned safety ring.
6. Advance needle into lesion, intra- or extra-hepatic bile ducts, pancreatic duct, cystic duct or gallbladder.
7. Remove stylet from needle by rotating stylet cap anti clockwise and gently pulling it back from metal fitting of handle. Preserve stylet for use if additional punctures are desired later.
8. For maintaining access a 0.035" wire guide (Cook Medical recommended) can be advanced through needle. Skip to step 12. **Note:** a 0.035 wire guide can also be advanced post aspiration. For FNA, attach Luer lock fitting of previously prepared syringe securely onto fitting on handle. For FNI, follow manufacturer's guidelines for use of injectable materials and skip to step 12.
9. Turn stopcock to "open" position aligned with syringe, allowing negative pressure in syringe to aspirate cells. Gently move handle in small increments, back and forth within biopsy site. **Note:** Do not remove needle from biopsy site during cell collection
10. Upon completion of FNA, turn stopcock to "close" position.
11. Retract needle into sheath by pulling back on handle shaft. Lock thumbscrew on safety ring at 0 cm mark.
12. Disconnect Luer lock fitting from accessory channel and withdraw entire device from ultrasound endoscope. **Note:** If a wireguide has been advanced, withdraw needle over wireguide.
13. The remaining steps are for FNA only.
14. Unlock thumbscrew on safety ring and extend needle. Exchange locking syringe for an empty syringe; attach to Luer lock fitting on handle.
15. Push in plunger to expel aspirate specimen, then prepare per institutional guidelines.
16. Remaining aspirate may be recovered for examination by reinserting stylet or flushing device.
17. For additional aspiration from same targeted site, gently reinsert stylet into metal fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer lock fitting, advance the stylet in small increments, until the hub is engaged in the fitting. Rotate stylet cap clockwise until the fittings are connected.
18. Additional samples may be obtained by repeating Step 2 through 16 of the "Instructions for Use".

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

EchoTip Ultra Endobronchial High Definition Ultrasound Needle ECHO-HD-22-EBUS-O with 10ml syringe; EchoTip Ultra Endobronchial High Definition Ultrasound Needle ECHO-HD-25-EBUS-O with 10ml syringe

ECHOTIP ULTRA ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is

attached to the accessory channel of an Olympus endobronchial ultrasound (EBUS) scope using an adapter. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for sample collection, if desired, and an accessory channel adapter. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic ultrasound procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution. This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbcrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with an Olympus EBUS scope.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

1. Examine the adapter and note features. (Fig. 1)
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub. (Fig. 2)
4. Close the adapter ensuring it is locked. (Fig. 3)
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (Fig. 4)
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d) Set prepared syringe aside until aspiration is required.



INSTRUCTIONS FOR USE







1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. (Fig. 5)
Caution: If resistance is encountered on the needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the endoscopes accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 6). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of ultrasound scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for the needle advancement appears in the window of the safety ring. (Fig. 7) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (steps 9-11) or, if desired, other techniques that may or may not incorporate the use of this stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during sample collection.

12. Upon completion of sample collection, if using a syringe, turn the stopcock to the "closed" position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
 13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.
 14. Unlock the thumbscrew on the safety ring and extend the needle.
 15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
 16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
 17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. Note: Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
 18. Additional samples may be obtained by repeating Step 2 through 16 of the "Instructions for Use".
 19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise and withdraw the entire device from the scope.
 20. To remove the adapter from scope, pull back on the locking clip, lift the top and slide off the scope.
- Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes ECHO-HD-22-EBUS-O-C with 10ml syringe; EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Olympus Scopes ECHO-HD-25-EBUS-O-C with 10ml syringe

ECHOTIP PROCORE ENDOBRONCHIAL HD ULTRASOUND BIOPSY NEEDLE FOR OLYMPUS SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree or gastrointestinal tract.

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of an Olympus endobronchial ultrasound (EBUS) scope using an adapter. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for biopsy, if desired, and an accessory channel adapter. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy and gastrointestinal endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution. This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the target site.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with an Olympus EBUS scope.

PLACEMENT OF ACCESSORY CHANNEL ADAPTER

1. Examine the Adapter and note features. (Fig. 1)
2. Place the adapter, in the open position, in line with the scope accessory channel metal hub.
3. Slide the adapter on to the accessory channel ensuring that the base of the adapter slides under the accessory channel metal hub. (Fig. 2)
4. Close the adapter ensuring it is locked. (Fig. 3)
5. The adapter will now function as a Luer lock to which the needle can be connected. **Note:** If irrigation of accessory channel is desired/necessary, irrigate through closed adapter prior to insertion or after removal of needle.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (Fig. 4)
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d) Set prepared syringe aside until aspiration biopsy is required.

INSTRUCTIONS FOR USE

1. Identify the target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. (Fig. 5)
Caution: If resistance is encountered on the needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. Note: The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 6). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining position of ultrasound scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for the needle

advancement appears in the window of the safety ring. (Fig. 7) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.

7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.

8. Standard vacuum syringe techniques may be applied for biopsy (steps 9-11) or, if desired, other techniques that may or may not incorporate the use of this stylet may be used.

9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if an additional biopsy is required later.

10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer the fitting on the needle handle.

11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during for biopsy.

12. Upon completion of for biopsy, if using the syringe, turn the stopcock to the "closed" position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.

13. Disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.

14. Unlock the thumbscrew on the safety ring and extend the needle.

15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.

16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.

17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.

18. Additional samples may be obtained by repeating Step 2 through 16 of the "Instructions for Use".






19. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope.




20. To remove the adapter from the scope, pull back on the locking clip, lift the top and slide off the scope.

Upon completion of the procedure, dispose of the device, including the adapter, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

EchoTip Ultra Endobronchial High Definition Ultrasound Needle ECHO-HD-22-EBUS-P with 10ml syringe; EchoTip Ultra Endobronchial High Definition Ultrasound Needle ECHO-HD-25-EBUS-P with 10ml syringe

ECHOTIP ULTRA ENDOBRONCHIAL HD ULTRASOUND NEEDLE FOR PENTAX SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used to sample targeted submucosal and extramural lesions within or adjacent to the tracheobronchial tree through the accessory channel of an ultrasound endoscope for Fine Needle Aspiration (FNA).

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of a Pentax endobronchial ultrasound (EBUS) scope. The device has a removable sheath adjuster to allow the user to insert the needle into the scope in the same orientation repeatedly. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for sample collection if desired. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution. This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Ensure the stylet is fully inserted when advancing the needle into the target site.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with a Pentax EBUS scope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when the stopcock is in "open" position, aligned with syringe. (Fig. 1)

2. Prepare syringe as follows:

- a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
- b) Turn stopcock to "closed" position.
- c) Pull plunger back until it is locked into place at desired setting, creating a suction.
- d) Set prepared syringe aside until aspiration biopsy is required.

INSTRUCTIONS FOR USE

1. Identify the desired target site by endobronchial ultrasound.

2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.

3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.






4. Attach the device to the endoscopes accessory channel port by rotating the device handle clockwise until the fittings are connected.

5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred sheath length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 2). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining the position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for needle advancement appears in the window of the safety ring. (Fig. 3) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for sample collection (steps 9-11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.
9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting of the needle handle. Preserve the stylet for use if additional sample collection is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to aspirate the sample. Move the needle handle back and forth within the target site. **Note:** Do not remove the needle from the target site during the sample collection.
12. Upon completion of the sample collection, if using the syringe, turn the stopcock to the "closed" position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. To remove the needle from the scope: unlock the thumbscrew on the sliding sheath adjuster while maintaining attachment of the sliding sheath adjuster to the scope and remove needle. (Fig. 4.) Or, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope. **Note:** When reinserting device into scope for additional passes and before advancing it into target site, you must position the handle onto the scope such that the thumbscrews on the handle are aligned in the same direction as they were prior to device removal.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing device.
17. For additional sample collection from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer lock fitting, advance it in small increments, until the hub is engaged in the fitting.
18. If the needle was removed from the scope while maintaining attachment of the sheath adjuster to the scope in Step 13: readvance the sheath through the attached sheath adjuster, ensuring that both thumbscrews are aligned on same side. Tighten the thumbscrew on the sliding sheath adjuster at the preferred sheath length.
19. Additional samples may be obtained by repeating Step 2 through 16 of "Instructions for Use".

20. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope. Upon completion of the procedure, dispose of the device, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Pentax Scopes ECHO-HD-22-EBUS-P-C with 10ml syringe; EchoTip ProCore Endobronchial HD Ultrasound Biopsy Needle for Pentax Scopes ECHO-HD-25-EBUS-P-C with 10ml syringe

ECHOTIP PROCORE ENDOBRONCHIAL HD ULTRASOUND BIOPSY NEEDLE FOR PENTAX SCOPES

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

INTENDED USE

This device is used with an ultrasound endoscope for fine needle biopsy, (FNB), of submucosal and extramural lesions within or adjacent to the tracheobronchial tree.

DEVICE DESCRIPTION

The device consists of a handle with adjustable components to allow the user to adjust the extension of the needle and sheath. The device has a 4.1 Fr outer sheath for protection of the needle and the scope, and it is attached to the accessory channel of a Pentax endobronchial ultrasound (EBUS) scope.

The device has a removable sheath adjuster to allow the user to insert the needle into the scope in the same orientation repeatedly. The device has a maximum needle extension length of 5 cm and sheath adjustment length of 3 cm. The device is supplied with a 10 ml syringe to utilize the standard vacuum technique for biopsy if desired. The device is available in two different needle sizes: 22 and 25 ga.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired position to visualize the targeted site.

Relative contraindications include, but are not limited to: coagulopathy.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with bronchial endoscopy include, but are not limited to: allergic reaction to medication, allergic reaction to nickel, aspiration, cardiac arrhythmia or arrest, damage to blood vessels, death, discomfort, fever, hemorrhage, hypotension, infection, inflammation, nerve damage, pain, perforation, pneumoperitoneum, respiratory depression or arrest, sepsis, septicemia/bacteremia and tumor seeding.

WARNINGS

Not for use in the heart or vascular system.

The tip of the needle and stylet are sharp and could cause injury to the patient or user if not used with caution. This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

If the package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook Medical for return authorization.

PRECAUTIONS

Refer to package label for minimum channel size required for this device.

Ensure the stylet is fully inserted when advancing the needle into the target site.

The needle must be retracted into the sheath and the thumbscrew on the safety ring must be locked at the 0 cm mark to hold the needle in place prior to introduction, advancement or withdrawal of the device. Failure to retract the needle may result in damage to the scope.

Do not use this device for any purpose other than stated intended use.

Use of this device is restricted to a trained healthcare professional (Rx only).

When targeting multiple sites, adequate care should be taken to minimize the risk of tumor seeding.

This device is intended for use with a Pentax EBUS scope.

SYSTEM PREPARATION

1. Examine syringe. It has two plunger locks that must be depressed to advance plunger. The tip of the syringe has a Luer lock fitting with a stopcock on the side port. Air can be exchanged when stopcock is in "open" position, aligned with syringe. (Fig. 1)
2. Prepare syringe as follows:
 - a) With stopcock in "open" position, depress plunger locks and fully advance plunger into syringe.
 - b) Turn stopcock to "closed" position.
 - c) Pull plunger back until it is locked into place at desired setting, creating a vacuum.
 - d) Set prepared syringe aside until biopsy is required.


INSTRUCTIONS FOR USE





1. Identify the desired target site by endobronchial ultrasound.
2. With the needle retracted into the sheath and the thumbscrew on the safety ring locked at the 0 cm mark to hold the needle in place, introduce the ultrasound needle into the accessory channel of the scope. **Caution:** If resistance is encountered on needle introduction, reduce angulation of the scope until smooth passage is allowed.
3. Advance the device, in small increments, until the Luer lock fitting at the base of the sliding sheath adjuster meets the Luer fitting of the accessory channel port.
4. Attach the device to the accessory channel port by rotating the device handle clockwise until the fittings are connected.
5. Adjust the sheath to the desired position, ensuring that it is visible on the endoscopic image, emerging from the working channel of the scope. To adjust length, loosen the thumbscrew lock on the sliding sheath adjuster and slide until the preferred length is attained. **Note:** The reference mark for the sheath length will appear in the sliding sheath adjuster window (Fig. 2). Tighten the thumbscrew on the sliding sheath adjuster to maintain the preferred sheath length.
6. While maintaining the position of the scope, set the needle extension to the desired length by loosening the thumbscrew on the safety ring, and advancing it until the desired reference mark for the needle advancement appears in the window of the safety ring. (Fig. 3) Tighten the thumbscrew to lock the safety ring in place. **Note:** The number in the safety lock ring window indicates the extension of the needle in centimeters. **Caution:** During needle adjustment or extension, ensure the device has been attached to the accessory channel of the scope. Failure to attach the device prior to needle adjustment or extension may result in damage to the scope.
7. Extend the needle, by advancing the needle handle of the device to the pre-positioned safety ring, into the target site. **Caution:** If excessive resistance is encountered on needle advancement, retract the needle into the sheath with the thumbscrew locked at the 0 cm mark, reposition the scope and attempt needle advancement from another angle. Failure to do so may result in needle breakage, device damage or malfunction.
8. Standard vacuum syringe techniques may be applied for biopsy (steps 9-11) or, if desired, other techniques that may or may not incorporate the use of the stylet may be used.

9. Remove the stylet from the ultrasound needle by gently pulling back on the plastic hub seated in the Luer fitting on the needle handle. Preserve the stylet for use if additional biopsy is required later.
10. Attach the Luer lock fitting of the previously prepared syringe securely to the Luer fitting on the needle handle.
11. Turn the stopcock to the "open" position aligned with the syringe, allowing the negative pressure in the syringe to facilitate biopsy. Move the needle handle back and forth within the biopsy site. **Note:** Do not remove the needle from the biopsy site during biopsy.
12. Upon completion of biopsy, if using the syringe, turn the stopcock to the "closed" position. Retract the needle completely into the sheath by pulling back on the needle handle and lock the thumbscrew on the safety ring at the 0 cm mark to hold the needle in place.
13. To remove the needle from the scope: unlock the thumbscrew on the sliding sheath adjuster while maintaining attachment of the sheath adjuster to the scope, and remove needle. (Fig. 4.) Or, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counterclockwise, and withdraw the entire device from the scope. **Note:** When reinserting device into the scope for additional passes and before advancing it into target site, you must position the handle onto the scope such that the thumbscrews on the handle are aligned in the same direction as they were prior to device removal.
14. Unlock the thumbscrew on the safety ring and extend the needle.
15. Use a new syringe or the stylet to expel aspirate the specimen, then prepare per institutional guidelines.
16. Remaining material may be recovered for examination by reinserting the stylet or flushing the device.
17. For additional biopsy from the same target site, gently reinsert the stylet into the Luer fitting on the needle handle. **Note:** Prior to reinserting the stylet, wipe with saline or sterile water. While supporting the stylet at the Luer fitting, advance the stylet in small increments, until the hub is engaged in the fitting.
18. If the needle was removed from the scope while maintaining attachment of the sheath adjuster to the scope in Step 13: readvance the sheath through the attached sheath adjuster, ensuring that both thumbscrews are aligned on same side. Tighten the thumbscrew on the sliding sheath adjuster at the preferred sheath length.
19. Additional samples may be obtained by repeating Step 2 through 16 of "Instructions for Use".
20. Upon completion of the procedure, disconnect the Luer lock fitting of the device from the accessory channel port, by rotating the device handle counter clockwise, and withdraw the entire device from the scope. Upon completion of the procedure, dispose of the device, per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Store in a dry location, away from temperature extremes.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Please notify Cook for return authorization.

Cook Medical Europe Ltd

O'HALLORAN ROAD

NATIONAL TECHNOLOGY PARK

LIMERICK, IRELAND

TEL: (+353 61 334440 FAX: (+353 61 334441

WWW.COOKMEDICAL.COM

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SCHAG".

Russia, 119002, Moscow, Karmanitskiy per., house 9, floor 5, office 501.

phone: (+95) 956-13-09;

Fax: (+95) 956-1310.

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

SHELF LIFE

The "use by" date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

The shelf-life is 3 years

STERILIZATION DETAILS

Endoscopic Biopsy Needle are supplied sterile

STORAGE RULES

The devices are to be stored in a dry location, at room temperature.

Endoscopic biopsy needles are single use only, do not attempt to repair, carry out maintenance or reuse.


DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste, class B.

Signed:

Ayesha Ford
Regulatory Affairs Specialist
Wilson-Cook Medical, Inc.


06 May 2020
Date

State of North Carolina, Forsyth County

I, Melissa Bottoms a Notary Public certify that, Ayesha Ford personally came before me and acknowledged that she is Regulatory Affairs Specialist of Wilson-Cook Medical Inc Cook Endoscopy, and acknowledged on behalf of the corporation, the due execution of the foregoing instrument.

Witness my hand and notarial seal this 16 day of May, 2020.

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