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Warnings

Do not use if package is damaged. Do not use if the product sterilization barrier or its packaging is compromised.

- Contents are sterile unless package is opened or damaged. DO NOT RESTERILIZE. For single use only. Do not use after the expiration date.
- It is the surgeon's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.
- It is extremely important to select the appropriate size and type components. Failure to use the largest possible components or improper positioning may result in loosening, bending, cracking, or fracture of the device or bone or both.
- Use only Ti-6Al-4V screws with Ti-6Al-4V devices. Use only stainless steel screws with stainless steel devices.
- Smith & Nephew Variable-Angle Locking Mini-Fragment Plating System devices and accessory components should not be placed across growth plates in pediatric patients.

Precautions

U.S. Federal law restricts this device to sale by or on the order of a physician.

- Use extreme care handling and storing implant components. Cutting, bending or scratching the surface of metal components may cause internal stresses which may significantly reduce the strength and fatigue resistance.
- Intraoperative fracture or breaking of instruments can occur. Instruments which have experienced extensive use or excessive force are susceptible to fracture. Instruments should be examined for wear and damage prior to surgery. Single use devices should not be reused due to risks of breakage, failure or patient infection.
- Continuous screening with an image intensifier (fluoroscopy) during guide wire insertion and whenever cannulated instruments are advanced over a guide wire is recommended to prevent unintended guide wire advancement and penetration into the surrounding tissues.
- Intraoperative cleaning of cannulated instruments is recommended to prevent accumulation of bone debris in the cannulation.
- If a computer assisted surgery system is used, consult the applicable software and hardware reference manuals provided by the manufacturer to ensure proper operation of this equipment.

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- For computer assisted surgery systems, it is extremely important to correctly select input parameters (e.g., bony landmarks). Operators of this equipment should be familiar with the anatomy relevant to the procedure. Failure to provide proper input could cause problems such as violation of critical anatomical structures and malpositioned implants.
- Postoperative instructions to patients and appropriate nursing care are critical. Early load bearing substantially increases implant loading and increases the risk of loosening, bending or breaking the device. Early load bearing should only be considered where there are stable fractures with good bone-to-bone contact.

- Patients should be directed to seek a medical opinion before entering potentially adverse environments that could affect the performance of the implant, such as electromagnetic or magnetic fields, including a magnetic resonance environment.

- While the surgeon must make the final decision regarding implant removal, wherever possible and practical for the patient, fixation devices should be removed once their service as an aid to healing is complete.

- After use, these components may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.

- Surgical technique information is available on request.

Adverse Effects

- Loosening, bending, cracking or fracture of implant components.

- Loss of anatomic position with malunion may occur.

- Infections, both deep and superficial, have been reported with internal fixation devices.

- Vascular disorders including thrombophlebitis, pulmonary emboli, wound hematomas, and avascular necrosis may result from the surgery and concomitant use of internal fixation devices.

- Metal sensitivity reactions and/or allergic reactions to foreign materials.

- Penetration of a K-wire, screw, or peg into a joint can occur.

- Tissue reactions which include macrophage and foreign body reactions adjacent to implants can occur.

Magnetic Resonance Imaging (MRI) Safety

Smith & Nephew fracture fixation devices have not been reviewed by the FDA for safety and compatibility in the MR environment. Fracture fixation components have not been tested for heating or migration in the MR environment. Known risks of exposing implant devices to the MR environment include displacement, torque, and radio frequency induced heating. Implant devices may also create image artifacts in MR scans.

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Packaging and Labeling

Components should only be accepted if received by the hospital or surgeon with the factory packaging and labeling intact. If the sterile barrier has been broken or is past the expiration date, return the component to Smith & Nephew, Inc. Use of damaged or expired product increases the risk of infection, which may lead to revision surgery.

Sterilization

For components provided sterile, the sterilization method is noted on the label. Sterile implant components are supplied sterile to a Sterility Assurance Level (SAL) of 10^{-6} . Sterile packaged components are supplied in protective sterile barrier packaging. Inspect packages for punctures or other damage prior to surgery. If the sterile barrier has been broken, return the component to Smith & Nephew, Inc. If not specifically labeled sterile, components are supplied non-sterile and must be cleaned and sterilized prior to surgery. For non-sterile trauma implants (i.e. plates and screws), remove all original packaging and labeling inserts prior to sterilization. It is important that adequate cleaning be carried out prior to sterilization. Please see the document, "Recommendations for decontamination and sterilization of Smith & Nephew orthopaedic devices," which is provided with Smith & Nephew instrument sets, for further information on cleaning instructions and validated sterilization procedures. DO NOT REUSE implant components or single use disposable instruments.

Recommended Steam Sterilization Cycle Parameters

- Dynamic Air Removal (Prevacuum) Steam Cycle: 132° C (270° F) for 4 minutes or 135° C (275° F) for 3 minutes. Minimum vacuum drying time:
 - Wrapped devices - 15 minutes
 - Containerized devices - 30 minutes
 - Gravity Displacement Steam Cycle: 132° C (270° F). Exposure time:
 - 15 minutes for instruments not in a containment device
 - 30 minutes* for devices in a containment device
- Minimum vacuum drying time of 30 minutes.
- Immediate Use Steam Sterilization (IUSS): 132° C (270° F). Exposure time:
 - dynamic air removal (prevacuum) - 4 minutes.

For Non-U.S. Customers:

- United Kingdom Steam Cycle: Prevacuum cycle. 134° C (273° F) for 3 minutes. Minimum vacuum drying time of 30 minutes. (Note: Sterilization evacuation and pulsing should be carried out in accordance with HTM 2010).

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- World Health Organization (WHO) Steam Cycle: 134° C (273° F). Exposure time 18 minutes and a minimum vacuum drying time of 30 minutes.

Containment devices should be wrapped with a central supply wrap (CSR) or placed in a reusable rigid container for sterilization.

Note to U.S. Customers: Sterilizers and wraps used in the sterilization process must be FDA-cleared.

*This sterilization cycle is not considered by the United States Food and Drug Administration (US FDA) to be a standard sterilization cycle. Only use sterilizers and accessories (such as sterilization wraps sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been FDA-cleared for the selected sterilization cycle specifications (time and temperature).

Retrieval and Analysis of Removed Implants

The most important part of surgical implant retrieval is preventing damage that would render scientific examination useless. Special care should be given to protect the implant from damage during handling and shipment. Follow internal hospital procedures for the retrieval and analysis of implants removed during surgery. When handling removed implants, use precautions to prevent spread of bloodborne pathogens. If the implant will be returned to Smith & Nephew, Inc. for analysis, contact Customer Service using the phone numbers outlined below.

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Definitions

Biological Indicator (BI): Test system that contains viable microorganisms that provide a defined resistance to a specified sterilization process.

Caddy: A small flip top containment device that contains multiple small implants (i.e. screws, plates, etc.).

Chemical Indicator (CI): Test system that reveals change in one or more pre-defined process variables based on a chemical or physical change resulting from exposure to a process. CIs assist in the detection of potential sterilization failures that could result from incorrect packaging, incorrect loading of the sterilizer, or sterilizer malfunctions. The "pass" response of a CI does not prove that the item monitored by the indicator is sterile.

Cleaning: The removal of contaminants to the extent necessary for further processing or for intended use.

Containment Device: A reusable rigid sterilization container, instrument case/cassette, caddy, or organizing tray and any reusable accessories for use in healthcare facilities for the purpose of containing reusable medical devices for sterilization.

Critical Water: Water that has been extensively treated (usually by a multi-step treatment process that could include a carbon bed, softening, Deionized (DI), and Reverse Osmosis (RO) or distillation) to ensure that the microorganisms and the inorganic and organic materials are removed from the water. A final submicron filtration could also be part of the treatment process. This water is mainly used for the final rinse or steam generation.

Disinfection: A process to reduce the number of viable microorganisms to a level previously specified as being appropriate for a defined purpose.

Immediate Use: The shortest possible time between a sterilized item's removal from the sterilizer and its aseptic transfer to the sterile field. A sterilized item intended for immediate use is not stored for future use or held from one case to another.

Manual Cleaning: The removal of contaminants from an item to the extent necessary for further processing or for intended use without the use of an automated process. Ultrasonic cleaning is considered a manual cleaning step for the purposes of this document.

Manufacturer: The natural or legal person with responsibility for design and/or manufacture of a medical device with the intention of making the medical device available for use, under his name; whether or not such a medical device is designed and/or manufactured by that person himself or on his behalf of another person(s).

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Medical Device: Instrument, apparatus, implant, reagent for in-vitro use, software, material or other similar or related article intended by the manufacturer to be used alone or in combination, for human beings for one or more of the specific medical purpose(s) of:

- diagnosis, prevention, monitoring, treatment or alleviation of disease; diagnosis, monitoring, treatment, alleviation of or compensation for an injury; investigation, replacement, modification or support of the anatomy or of a physiological process; supporting or sustaining of life; control of conception; disinfection of medical devices; providing information by means of in vitro examination of specimens derived from the human body; and does not achieve its primary intended action by pharmacological, immunological or metabolic means, in or on the human body, but which may be assisted in its function by such means.

Processing: Preparation of medical devices or activity to prepare a new or used healthcare products for its intended use. In this document processing includes cleaning, disinfection and sterilization.

Reusable Medical Device: A medical device designated or intended by the manufacturer as suitable for processing and reuse.

Single Use Medical Device (SUD): A medical device that is designated or intended by the manufacturer for one-time use only. For clarity, a single use device that has come into contact with blood, tissue, or bodily fluids is not intended to be further processed and used again.

Sterile: Free from viable microorganisms.

Sterility Assurance Level (SAL): The probability of a single viable microorganism occurring on an item after sterilization, expressed as the negative exponent to the base 10.

Sterilization: The process used to render product free from viable microorganisms.

Utility Water: Water as it comes from the tap that might require further treatment to achieve the specifications. This water is mainly used for flushing, washing, and rinsing.

Validation: Confirmation, through the provision of objective evidence, that the requirements for a specific intended use or application have been fulfilled.

Washer-Disinfector: Equipment that is designed to clean and disinfect product.

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Recommended Cleaning Instructions

Cleaning is the single most important step in preparing a device for reuse. Effective cleaning must be carried out to achieve proper disinfection/sterilization. Cleaning begins at the point of use to prevent soil and contaminants from drying on the medical devices after use. Thorough cleaning and rinsing are vital to processing reusable medical devices. Also, thorough rinsing is important for the removal of any residual cleaning agents from the medical devices. The purpose of cleaning and rinsing is to remove all adherent visible debris and to reduce the number of particulates, microorganisms, and pyrogens. The recommended cleaning instructions in this document include both manual and automatic washing/disinfection procedures. While manual cleaning is the most universal method of cleaning, automatic washing is preferred. The cleaning processes presented in this brochure have been validated. Other methods of cleaning may be suitable but must be validated by the user of the device.

Important Information/Recommendations for Use

Automatic washer/disinfector

Washer-disinfectors are not only used to clean devices, but also to provide intermediate to high level disinfection with a hot water rinse. Cleaning is dependent upon thorough coverage of the devices and the force of the water spray. Therefore, all sections of the device must be accessible for ease of cleaning and penetration of cleaning agents. The automatic washer/disinfector equipment should be operated following the manufacturer's instructions for use.

Detergents

Low foaming detergents with a pH range between 6.0 and 8.0 are recommended. Detergents with a pH outside this range can have an adverse effect or be damaging to some medical devices and containment devices. Enzymatic detergents aid in the removal of organic soil such as blood.

Detergents should be used at the concentration and temperature recommended by the detergent manufacturer.

Some alkaline detergents have been formulated to be safe for processing medical devices. Refer to the manufacturer's limitations and warnings for information concerning specific materials that are adversely affected by the detergents.

Loaner Instrument Set

It is the hospital's responsibility to clean, disinfect, package and sterilize all loaner devices upon receipt. In addition, it

is also the responsibility of the hospital to clean, disinfect, package and sterilize the instrument sets prior to their return to Smith & Nephew.

Manual cleaning tools

Hospital tools necessary for manual cleaning include: surgical scrub brushes, chenille pipe cleaners, soft low linting cloths, cotton tip

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applicators, and several various sizes and lengths of brushes.

Do not use abrasive cleaning tools (i.e. scouring pads or metal brushes).

Cleaning tools must be cleaned and inspected after use. Cloths should be clean and lint free and changed frequently. Brushes should be clean. Discard worn brushes and disposable cleaning tools.

Safety precautions

Personal Protective Equipment (PPE) should be worn when handling or working with contaminated or potentially contaminated devices. PPE includes: gown, mask, goggles or face shield, and shoe covers.

Universal precautions are standards of infection control practices designed to reduce the risk of transmission of bloodborne infections. Universal precautions should be observed by all hospital personnel that work with contaminated or potentially contaminated devices. Exercise caution when handling devices with sharp points or cutting edges.

Water

The quality of water should be carefully considered for use in cleaning reusable devices. The water quality can affect the life of the device. Water hardness is a concern because deposits left on medical devices may result in ineffective cleaning and sterilization. Final rinsing should be carried out using critical water.

Warnings

- It is the responsibility of the user to ensure that the cleaning process is performed following these procedures to achieve the desired result.
- For ultrasonic cleaning, the enzymatic detergent solution should be changed before it becomes heavily soiled so that effective cleaning is not inhibited.
- Most electronic devices cannot be submerged during cleaning and could sustain permanent damage as a result of submersion. Check the "Instructions for Use" for information on submersion of electronics.
- Do not sonicate powered devices.
- All cleaning should be performed in a manner designed to minimize exposure to bloodborne pathogens. Manual cleaning should be done while the instrument is immersed.

Limitations on Processing

Limitations of cleaning instructions

These recommended procedures are intended as a general guide for cleaning of medical devices. Some devices are labeled with more specific instructions. Refer to the Instruction for Use for that device.

Limitations of processing

Repeated processing of reusable medical devices has minimal effect on the devices.

The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful

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Flexible reamers: Includes all devices that have a laser cut flexible reamer design.

Examples: Flexible screw drills and flexible shafts.

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Powered devices: Those that use pneumatic power or have power cords and require the use of electricity to operate.

Examples: ACCURIS handpiece and pulse lavage handpiece.

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Containment devices: Instrument cases, trays, caddies and lids

Non-critical devices: Those that do not contact the patient or may contact intact skin.

Examples: Power supply boxes and power cable surfaces.

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Initial Treatment at the Point of Use

Cleaning begins at the point of use to prevent soil and contaminants from drying on the devices. Keep the devices moistened after use and before cleaning.

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Transportation

Used devices must be transported to the processing location (i.e. Sterile Processing) in closed or covered containers to prevent contamination risks.

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Preparation Before Cleaning

Disassemble each device into its component parts. Containment devices must be cleaned separately from the medical devices. If you have questions concerning the disassembly of any Smith & Nephew Orthopaedics device, contact your Smith & Nephew sales representative.

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Water Quality for Cleaning

Use utility water for flushing, washing and rinsing during cleaning except for the final rinse. For the cleaning final rinse use critical water.

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Cleaning: Manual Devices Without challenging design features

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.
2. Immerse and soak for a minimum of one (1) minute in enzymatic detergent.
3. With the device completely immersed in the cleaning solution, use a surgical scrub brush to

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remove visible debris from the surfaces of the device.

4. Rinse for a minimum of one minute using a sufficient volume of critical water that will completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing.

5. Check for visible debris. Repeat cleaning if debris is visible.

Devices With challenging design features

1. Rinse in cold water <45 C (113° F) to remove visible debris and

to prevent coagulation of blood.

2. Immerse and soak for a minimum of five (5) minutes in enzymatic detergent.

3. With the device completely immersed in the cleaning solution, remove additional soil from the surfaces and from the challenging design features (i.e. holes, lumens, hinged/mating surfaces, interfaces, crevices, serrations) using common hospital cleaning tools.

a. Move and/or retract all moveable parts and remove visible debris using a brush.

b. Scrub lumens or holes with a brush of an appropriate size

to ensure that the full width and depth is accessed. Use a twisting action with the brush. Thoroughly flush lumens with enzymatic detergent using a syringe.

c. Open hinged devices and scrub hinged area with a brush. d. Scrub crevices with a brush.

e. Scrub the surfaces of the device with a toothbrush styled brush.

4. Sonicate the device in enzymatic detergent in its fully opened position for a minimum of 15 minutes in an ultrasonic cleaner.

5. Rinse for a minimum of one minute using a sufficient volume of critical water that will completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing. If the components of the device are movable or can be retracted, it is necessary to do so during the rinsing process. Thoroughly flush lumens with critical water using a syringe.

6. Check devices for visible debris (see "Verifying Cleaning").

Repeat cleaning if debris is visible.

Flexible reamer devices

1. Gently bend the flexible reamer in a U-shape and rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.

2. Immerse and soak for a minimum of ten (10) minutes in enzymatic detergent.

3. With the device completely immersed in the cleaning solution, scrub lumen with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. (The bristles should be stiff enough to remove bone and tissue.) Thoroughly flush lumens with enzymatic detergent using a syringe.

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4. With the device completely immersed in the cleaning solution, scrub the surface with a surgical scrub brush to remove all visible debris from the surface and crevices.

5. With the device completely immersed in the cleaning solution, gently bend the reamer in a U-shape and scrub the surface with a scrub brush.

Note: Gently bend at several locations along the length to access all crevices.

6. Gently bend the reamer in a U-shape and rinse thoroughly with warm water making sure to irrigate the lumens and crevices. Thoroughly flush lumens with water using a syringe.

Note: Use a clean brush during the rinse cycle and move it back and forth through the lumen several times during rinsing.

7. Sonicate in enzymatic detergent for a minimum of 15 minutes in an ultrasonic cleaner.

8. Rinse for a minimum of one minute using a sufficient volume of critical water to completely immerse the instrument. Water should be changed a minimum of two times to ensure thorough rinsing. Thoroughly flush lumens with critical water using a syringe.

9. Check the devices for visible debris (see "Verifying Cleaning"). If debris is visible, repeat cleaning steps 2 – 8.

Powered Devices

Leave the device that directly connects to the motorized device (i.e. the hose that connects to a powered handpiece) in place during cleaning to prevent an excessive amount of water from entering the motor.

Do not sonicate powered devices.

Clean the Motorized Device as Follows:

1. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.

2. With the connecting device still in place, immerse the motorized device, with the connecting instrument still in place, in enzymatic detergent and soak for three minutes.

3. With the device completely immersed in the cleaning solution, scrub the surfaces of the device in the enzymatic detergent using cleaning brushes.

4. Rinse for a minimum of one minute using a sufficient volume of critical water to completely immerse the instrument. Water should be changed a minimum of two times to ensure thorough rinsing.

5. Remove the motorized device from the water before disassembly of the connecting device. Clean the connecting device separately following the cleaning instructions for the correct Cleaning Product Group.

Containment Devices

1. Remove medical devices from the containment device prior to cleaning. Inspect the containment device for visible debris.

2. If dried visible debris is observed, follow the manual instructions for "Devices with challenging design features."

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3. If dried visible debris is not observed, clean following the manual instructions for "Devices **without** challenging design features."

Non-critical Devices

Sterile hospital grade covers can be used whenever possible to prevent contamination of non-critical devices. Porous items like foam that are contaminated must be discarded. Clean surfaces by wiping with a cloth and using detergent and water.

Verifying Cleaning

1. After cleaning, visually inspect devices under normal lighting for any remaining visible debris.

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3. If not visibly clean, repeat cleaning and reinspect the device.

Drying

Use clean filtered compressed air or a clean lint free towel to

remove moisture from the devices.

Cleaning: Automated

Devices Without challenging design features

1. Manual precleaning:

IS NOT REQUIRED if the device does not have dried-on visible debris. Place the device directly into the automatic washer for cleaning.

IS REQUIRED if the device does have dried on visible debris. Follow the manual cleaning steps below prior to placing the device in the automatic washer for cleaning.

a. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.

b. Immerse and soak for a minimum of one (1) minute in enzymatic detergent.

c. With the device completely immersed in the cleaning solution, use a surgical scrub brush to remove visible debris.

d. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Water should be changed a minimum of two times to ensure thorough rinsing.

2. Place the device in the automatic washer and run the recommended automatic washer steps (see section "Automatic Washing Cycle Steps and Parameters.")

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a. Rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.

b. Immerse and soak for a minimum of five (5) minutes in enzymatic detergent.

c. With the device completely immersed in the cleaning solution, remove additional visible

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debris from the device surfaces and from the challenging design features (i.e. holes, lumens, hinged/mating surfaces, interfaces, crevices, serrations) using common hospital cleaning tools.

-- Move and/or retract all moveable parts and use a brush to remove visible debris.

-- Scrub lumens and holes with a brush of an appropriate size to ensure that the full width and depth is accessed.

Use a twisting action with the brush. Thoroughly flush lumens with enzymatic detergent using a syringe.

-- Open hinged devices and scrub hinged area with a brush.

-- Scrub crevices with a brush.

d. Sonicate the device in its fully opened position for a minimum of 15 minutes in an ultrasonic cleaner containing enzymatic detergent.

e. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. If the components of the instrument are movable or can be retracted, do so during the rinsing process. Change the water a minimum of two times to ensure thorough rinsing. Thoroughly flush lumens with water using a syringe.

2. Load the devices into the washer such that all design features of the device are accessible to cleaning and such that design features that might retain liquid can drain (for example, hinges should be open and lumens and holes positioned to drain).

3. Run the recommended automatic washer steps (see section "Automatic Washing Cycle Steps and Parameters.")

4. Check devices for visible debris (see "Verifying cleaning"). Repeat cleaning if debris is visible and reinspect.

Flexible reamer devices

1. Manual precleaning is required for all devices in this product group. Follow the manual cleaning steps below prior to placing the device in the automatic washer for cleaning.

a. Gently bend the flexible reamer in a U-shape and rinse in cold water <45° C (113° F) to remove visible debris and to prevent coagulation of blood.

b. Immerse and soak for a minimum of ten (10) minutes in enzymatic detergent.

c. With the device completely immersed in the cleaning solution, scrub the lumen with a brush of an appropriate size to ensure that the full width and depth is accessed. Use a twisting action with the brush. (The bristles should be stiff enough to remove bone and tissue.)

d. With the device completely immersed in the cleaning solution, scrub the surface with a surgical scrub brush to remove all visible debris from the surface and crevices. Thoroughly flush lumens with enzymatic detergent using a syringe.

e. With the device completely immersed, gently bend the reamer in a U-shape and scrub the surface with a scrub brush. **Note:** Bend at

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several locations along the length to access all crevices.

f. Gently bend the reamer in a U-shape and rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Change the water a minimum of two times to ensure thorough rinsing. Make sure to irrigate the lumen and crevices. Thoroughly flush lumens with enzymatic detergent using a syringe.

g. Sonicate for a minimum of 15 minutes in an ultrasonic cleaner containing enzymatic detergent.

2. Place the device in the automatic washer and run the recommended automatic washer steps (see the section "Automatic Washing Cycle Steps and Parameters").

Powered devices

Leave the device that directly connects to the motorized device (i.e. the hose that connects to a powered handpiece) in place during cleaning to prevent an excessive amount of water from entering the motor. Do not sonicate powered device. Clean the motorized device as follows:

1. Rinse in cold water $<45^{\circ}\text{C}$ (113°F) to remove visible debris and to prevent coagulation of blood.

2. With the device completely immersed in the cleaning solution, use a surgical scrub brush to remove visible soil from the surfaces of the device and connecting device.

3. Rinse for a minimum of one minute using a sufficient volume of water to completely immerse the device. Change the water a minimum of two times to ensure thorough rinsing.

4. Place the motorized device with the connecting device attached to it in the automatic washer and run the recommended automatic washing steps (see the section "Automatic Washing Cycle Steps and Parameters").

Clean the connecting device separately following the cleaning instructions for the correct Cleaning Product Group.

Containment devices

1. Remove all medical devices from the containment device and inspect it for visible debris.

2. If dried visible debris is observed follow the manual cleaning instructions for "Devices with challenging design features".

a. Then place the containment device in the automatic washer and run the recommended automatic washer steps (see the section "Automatic Washing Cycle Steps and Parameters").

3. If dried debris is not observed, place the containment device in the automatic washer and run the automatic washer steps given below (see section "Automatic Washing Cycle Steps and Parameters").

Automatic washing cycle steps and parameters

Load the devices into the automatic washer-disinfector.

-- Connect any device lumen(s) to the rinsing ports of the washer-disinfector. If no direct connection is possible, place the lumen(s) on

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injector jets or in injector sleeves of the injector basket.

-- Place blind holes at a downwards incline for drainage.

-- Place all instruments in the open position whenever possible (i.e. hinged/mating features, articulating surfaces). Minimum cycle parameters: Five (5) minute cold prewash <45°C (113°F) Five (5) minute neutral enzymatic detergent wash

Five (5) minute neutral detergent wash

One (1) minute rinse with critical water.

Disinfection

All Smith & Nephew orthopaedics reusable devices can be thermally disinfected using the parameters below. minimum 91°C for 1 minute 25 minute minimum dry time. Use clean filtered compressed air or a clean lint free towel to remove remaining moisture from the devices.

Verifying cleaning

1. After cleaning, visually inspect devices under normal lighting for the removal of visible debris.

2. For difficult to view design features, apply 3% hydrogen peroxide. Bubbling indicates the presence of blood. Rinse the devices thoroughly with critical water following hydrogen peroxide testing.

3. Repeat cleaning if not visibly clean and reinspect the device.

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Inspection and function testing

All reusable devices: Visually inspect for damage or wear, including components in their disassembled state prior to re-assembly. Ensure components are re-assembled securely.

Hinged devices: Check for smooth movement of the hinge without excessive "play."

Locking mechanisms: Check for action.

Cutting features: Check edges for distortion and/or large nicks. Edges should be continuous.

Trials: Articular surfaces should be smooth and free of cracks and deep nicks.

Mating parts: Check to make sure that mating parts fit together without complications. Ensure components are re-assembled securely.

Reamer/drill bits: Inspect "chuck" end for burrs and distortion that might hinder insertion into a drill.

Hammering surfaces: Inspect for burrs and large nicks.

Driving devices: Inspect plastic ends for cracks and large nicks.

Metal surfaces: Inspect for corrosion and major deformation.

Powered devices: Verify that power is supplied when the device is turned on and ceases when the device is turned off.

Measuring devices (i.e. gauges, calipers): Check for legible measuring markings.

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Implants

Implants are sold as single-use devices. The method of initial sterilization for implants that are supplied sterile is noted on the package label. Do not reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.

Containment Devices

Smith & Nephew Orthopaedics containment devices are categorized into families based on the design, density and material of the inner and outer containment devices. Each family has been validated to a 10^{-6} Sterility Assurance Level (SAL). Smith & Nephew designs the proper placement of devices into each containment device layout. The layout has predetermined brackets and artwork descriptions. If devices are added, the user is responsible for validation of the new layout. When applicable, adhesive labels shipped with the containment device should be applied to the outer containment device. Large labels typically have a designated outlined box with the instructions "PLACE LABEL HERE." Small labels are placed on the end caps under the labels.

CI and BI placement: Some Smith & Nephew containment device designs have multiple layers of instrument caddies or trays. Testing has shown that the following CI and BI placement locations should be used for qualification testing by the hospital or healthcare facility:

-- **For containment devices that do not contain internal instrument case caddies:** At the center and at each corner of the bottom internal containment device or tray.

-- **For containment devices that contain internal instrument case caddies:** At each corner of the caddy on the bottom level of the case.

Containment device weight: Smith & Nephew containment devices are designed to achieve a total case weight (outer case, device, plus inner trays/caddies) of 25 pounds (11.34 kgs) or less. Please be aware that there are some older models that cannot be separated to achieve a weight below 25 pounds (11.34 kgs). These older models have been validated to achieve a 10^{-6} SAL. Only hospital trained personnel should be utilized for inspection and maintenance of containment devices. Modifications to the containment devices should be made only by Smith & Nephew unless the materials and instructions for modifications are

supplied by Smith & Nephew (see "Contact Information"). The sterilization of containment devices is validated with the devices placed and positioned in the predetermined placement locations of the containment device. A single absorbent towel (i.e. a huck towel) can be placed under the containment device to aid in drying. Smith & Nephew is the only authorized service/repair company for Smith & Nephew containment devices. Containment devices in need of repair/replacement must be returned to

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References

- AAMI TIR 12, "Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: a Guide for Device Manufacturers,"
- AAMI TIR 30, "A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices."
- AAMI TIR 34, "Water for reprocessing of medical devices."
- ANSI/AAMI ST79, "Comprehensive guide to steam sterilization and sterility assurance in health care facilities."
- AAMI/ANSI ST81, Sterilization of medical devices- Information to be provided by the manufacturer for the processing of resterilizable medical devices.
- AORN, Perioperative Standards and Recommended Practices.
- Centers for Disease Control, Guideline for Disinfection and Sterilization in Healthcare Facilities, 2008.
- Health Canada June 2011, "Guidance Document on Recognition and Use of Standards under the Medical Devices Regulations."
- ISO 17664, Sterilization of medical devices — Information to be provided by the manufacturer for the processing of resterilizable medical devices.
- German Redbook, Proper Maintenance of Instruments, Instrumenten-Aufbereitung [Instrument Preparation Working Group].
- Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff, March 17, 2015.
- UK Department of Health, Health Technical Memorandum 01-01: Decontamination of reusable medical devices. Parts A and B.
- ANSI/AAMI ST77 "Containment devices for reusable medical device sterilization".
- SO 15883-1, Washers-disinfectors- Part : Requirements for the development, validation and routine control of a sterilization process for medical devices
- ISO 17665-1 Sterilization of healthcare products, most heat- Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.

Contact Information

Please contact authorized representative Smith & Nephew, LLC:
105120, Moscow, 2nd Syromyatnichesky lane, 1, floor 9, prem 1, room 1
Tel: +7 (495) 755-55-03

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Glossary of Symbols

	Medical Device
	Contains hazardous substances
	Sterilized using vaporized hydrogen peroxide
	Single patient – multiple use
	Electronic instructions of use

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Interpretation of symbols used on the package.

	Manufactured by:
	Catalog number
	Authorized representative of the EU
	Batch code
NON-STERILE	Non-sterile
	Caution! Refer to the instructions for use
R ONLY	Caution: U.S. Federal law restricts this device to sale by or on order of a physician
	Conformity to EU standards
QTY: (1)	Quantity

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VLP MINI-MOD fixation system for small fragments instrument components are intended for use together with VLP MINI-MOD system components for small bone fragments fixation.

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The subject Smith & Nephew products do not incorporate non-viable materials of animal origin, nor are they present in the manufacturing process or present in the final device.

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The subject Smith & Nephew products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device, therefore this section does not apply.

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The products are intended for use in the internal environment of the body. Avoid deliberate ingress of moisture during operation. If the rules of operation are observed, this medical device does not affect people with medical products implanted into the human body, pregnant women, and women during breastfeeding, children, adults with chronic diseases.

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Instruments should be kept in a designated place, access to which is restricted, which is well ventilated and protected from dust, moisture, insects, rodents and excessive temperature and humidity.

Reusable devices that will be stored between cleaning and sterilization should be dried with a low-linting, non-abrasive soft cloth to prevent microbial contamination that could result from wet instruments. Containment devices can be stacked for storage.

Recommended transportation and storage ranges are presented below.

Transportation conditions	
Humidity	40 ...90%
Temperature	-20...35°C
Pressure	86.6...106.7kPa
Storage conditions	
Humidity	40 ...90%
Temperature	-20...35°C
Pressure	86.6...106.7kPa

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Blunt or damaged tools should be disposed as epidemiologically hazardous waste (hazard class B). Unused products can be disposed as household waste (hazard class A). If the implant needs to be returned to Smith & Nephew, Inc. for analysis, contact by the phone number listed below.

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Smith & Nephew LLC (S&N) warrants all of its products to the original purchaser (Customer) against defect in workmanship and materials for one hundred eighty (180) days from date of purchase from Smith & Nephew; except for NAVIO products, one (1) year from date of installation (Warranty Period). For any products found to not be in conformance with this warranty during the

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Warranty Period, this warranty provides and is restricted to, as elected by S&N, either (i) repair or replacement of such products without charge and within a reasonable period of time or (ii) a refund or credit in the amount of the purchase price of such products. If S&N repairs or replaces product under this warranty and requests Customer to return such product, Customer must ship such product to S&N freight prepaid by Customer. Customer shall be invoiced for any replacement product if Customer does not return the requested replaced product within thirty (30) days after S&N's shipment of the replacement product. This warranty does not cover and is voided by any of the following: (i) product packaged or labeled by someone other than S&N or its authorized agents; (ii) product not used in compliance with the specifications, instructions or claims for use of the product; (iii) equipment product used in conjunction with disposables or accessories not specified for use with such equipment; (iv) equipment product used in conjunction with reprocessed disposables or accessories; (v) modification of product; (vi) product past its expiration date; (vii) normal wear and tear; (viii) damage due to misuse, reprocessing, alteration, unauthorized repair or negligent handling or damage due to lack of care by the owner, user or handler of the product including but not limited to storage, handling or cleaning; or (ix) any other damage inflicted to products by the owner, user or handler. A 90-day warranty on repairs applies to the defective component repaired. This warranty applies only to the original purchaser from S&N (or its authorized distributor) and is not transferable. EXCEPT TO THE EXTENT PROHIBITED OR OTHERWISE REQUIRED BY APPLICABLE LAW, THIS WARRANTY IS THE SOLE WARRANTY OF S&N, AND ALL OTHER WARRANTIES OF ANY KIND OR DESCRIPTION WHATSOEVER, INCLUDING WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY AND FITNESS FOR A PARTICULAR PURPOSE, EXPRESSED OR IMPLIED, ARE EXCLUDED.

**AUTHORIZED REPRESENTATIVE ON THE
TERRITORY OF THE RUSSIAN FEDERATION:**

Smith & Nephew, LLC, Russia, 105120, Moscow, 2nd Syromyatnichesky lane., 1, floor 9, prem. 1, room 1
Tel: +7 (495) 755-55-03

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State of Tennessee)
County of Shelby)

Subscribed (or affirmed) before me this 24th day
of March (month), 2020, by Maureen Whitson
(name of signer), whose identity was proven to me by personal
knowledge or satisfactory evidence.

Linda Rega (Signature of Notary)

(Seal of Notary)

Total number of pages: 41

Name: Maureen Whitson

Title: Regulatory Affairs Specialist

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