CARESTREAM HEALTH, INC
150 VERONA STREET
ROCHESTER, NEW YORK 14608

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Carestream

INSTRUCTION FOR USE





The RVG 5200, RVG 6200 digital intra-oral X-ray sensors, with accessories

Photo, diagram, drawing medical device

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MANUFACTURER

Carestream Health, Inc., 150 Verona Street, Rochester, New York, 14608, USA

Тел. 585-627-1800, Fax: 585-627-1800, E-mail: health.imaging.tsc@carestreamhealth.com

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Notice

The User & Installation Guide for RVG 5200 and RVG 6200 includes information on the device as well as its installation and usage. We recommend that you familiarize yourself with this guide to make the most effective use of your system.

RVG 5200 and RVG 6200 are intended to produce images of the dento-maxillofacial region of the human anatomy at the direction of health care professionals.

The information contained in this guide may be subject to modification without notice, justification or notification to the persons concerned.

No part of this guide may be reproduced without the express permission of Carestream Health, Inc.

The US Federal law restricts this device to sale by or on the order of a physician.

This document is originally written in English.

Manual Name: *RVG 5200 and RVG 6200 User & Installation Guide.*

Document Code: SM846

Revision Number: 02

Print Date: 2014-03

The brand names and logos reproduced in this guide are copyright.

RVG 5200 and RVG 6200 comply with Directive 93/42/EEC relating to medical equipment.

Important: We recommend that you consult the *RVG 5200 and RVG 6200 Safety, Regulatory and the Technical Specification User Guide* before using RVG 5200 and RVG 6200.

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INDICATION

The RVG 5200 and RVG 6200 digital intra-oral X-ray sensors are intended to produce an image of the dental area at the direction of dentists, and orthodontists for X-ray imaging of the dento-maxillofacial area.

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KNOWN CONTRAINDICATIONS

Allergies, destruction of the jaw, bleeding.

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MEDICAL RISKS AND UNDESIRE SIDE-EFFECTS

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SAFETY INSTRUCTIONS

The following special messages emphasize information or indicate potential risk to persons or equipment:

WARNING: Warns you to avoid injury to yourself or others by following the safety instructions precisely.

CAUTION: Alerts you to a condition that might cause serious damage

Important: Alerts you to a condition that might cause problems.

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Note: Emphasizes important information.

Tip: Provides extra information and hints.

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Read and understand this Safety Information before using RVG 5200 and RVG 6200.

You are responsible for the operation and maintenance of this device. Only legally qualified persons can operate this device. They MUST have training to use the device. When necessary, have a trained authorized service technician carry out inspection and maintenance operations.

Install this device in an X-ray room that complies with current installation standards. From this location, you must be able to maintain visual or audio communication with the patient and be able to access the Acquisition Interface module during exposure.

Do NOT operate the device if there is the threat of an earthquake. Following an earthquake, ensure that the device is operating satisfactorily before using it again. Failure to observe this precaution may expose patients to hazards.

X-ray equipment is hazardous to patients and the operator if you do not observe the exposure safety factors and operating instructions.

To dispose of the device or its components, contact a service technician.

No modification of this device is allowed

Do NOT use this device in conjunction with oxygen-rich environments. This device is not intended for use with flammable anesthetics or flammable agents.

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It is not recommended to use accessories other than those specified in this document and sold by Carestream Health.

Single-use disposable protective hygienic sleeves must cover the RVG sensor before placing it in the patient's mouth.
To prevent cross-contamination, use a new hygienic barrier for each new patient and disinfect the RVG sensor.

To ensure the best quality images, you must connect the RVG sensor to a USB 2.0 port on the BACK of the workstation.

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Functional Components Overview

Types of RVG Sensor

The RVG sensor is radio-sensitive. The active surface of the RVG sensor is the flat surface marked with #1 or #2 indicating the size.

RVG Sensor Size	Use
Size 1(universal RVG sensor)	Regular periapical and retro-coronary procedures.
Size 2	Bitewings and peri-apical procedures.

The non-reactive surface of the RVG sensor is rounded and contains the cable attachment.

RVG Sensor Overview

Figure 1: RVG Sensor

- 1.Active surface of the RVG sensor.
- 1.Non-reactive surface of the RVG sensor.
- 3.USB 2.0 connector.

Sharing the RVG Sensor Between Workstations

Single RVG Sensor/Multiple Workstations

You can share the RVG sensor between several workstations to provide access for several practitioners based on an agreed-upon arrangement. Each workstation must have the Carestream Dental Imaging Software and corresponding drivers installed. To share the RVG sensor between several workstations, move it from workstation to workstation. The RVG sensor is automatically recognized and operational when you connect it to a USB 2.0 port on the back of the workstation.

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Sharing Images Between Workstations

To share images between workstations, you can connect the workstations to a network without having to change the configuration described above. The Carestream Dental Imaging Software only needs to access a shared database on the same workstation or on a remote workstation.

You can print images either on a printer attached to each workstation or to a printer shared on the network.

Using the Different Positioning Systems

Apply the same rules for positioning the RVG sensor in the mouth that you use in classic radiology. You can use different systems for positioning the RVG sensor in the mouth.

A starter kit is provided with each RVG sensor.

X-ray Generator Compatibility

The RVG sensor is compatible with all X-ray generators that meet the current standards of intraoral radiology. We recommend a high-frequency X-ray generator. The X-ray generator must operate with a voltage of 60 to 70 kV.

X-ray generators from Carestream meet the requirements.

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Important: The RVG sensor is NOT compatible with X-ray generators that have a voltage LESS THAN 60 kV.

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Imaging Software Overview
Computer System Requirements

For the minimum computer system requirements for RVG 5200 and RVG 6200, see the RVG 5200 and RVG 6200 Safety, Regulatory and Technical Specifications User Guide.
If necessary you must update your computer system configuration.
RVG 5200 and RVG 6200 operate with the Carestream Dental Imaging Software.

You can acquire:

- Single images.
- Full mouth series (FMS).

Single Image Acquisition Overview

You can acquire individual images using the Carestream Dental Imaging Software.
Figure 2 Carestream Dental Imaging Software with One Active RVG Sensor Connected.

The Carestream Dental Imaging Software can display up to three RVG sensors connected to the workstation.

Icon	Explanation
	RVG sensor is connected to the workstation and is ready for image acquisition.
	RVG sensor is connected to the workstation but is not ready for image acquisition.

Full Mouth Series (FMS) Image Acquisition Overview

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
Click  in the Imaging Window to access the RVG FMS Acquisition interface.

The Full Mouth Series (FMS) is a static representation of the mouth of the patient using a series of intraoral images.

The images are placed in fixed numbered frames.

Figure 3: RVG FMS Acquisition Interface Home Page

1. Dental arch: Highlights the acquisition zone.

2. FMS template: Displays frame templates for acquisition.

- Green highlight: Frame ready for new acquisition

- Blue highlight: Frame in the revue and retake mode. This mode interrupts the automatic acquisition sequence. The retake images are displayed in the Retake Image gallery.

3. Retake Image gallery: Displays all the retake images acquired for a specific frame.

4. Preview screen: Displays the current acquired image.

5. Preferences: Displays the preferences for FMS template selection.

6. Available RVG sensors: Displays a maximum of three RVG sensors with their name.

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

- Green:RVG sensor is connected to the workstation and is ready for image acquisition.
 - Blue:RVG sensor is connected to the workstation but is not ready for image acquisition.
7. Refresh: Restarts the timer.
 8. Timer: Displays the timer for acquisition.

Figure 4: FMS Preferences


Preferencesenable you to select:
 You can select your preferences before starting to acquire images.
 If you try to change the FMS template after you have finished your image acquisitions you are prompted with a warning that indicates that you risk loosing some of the images.

Enhancement applied at acquisition

Enhancement type applied to acquired images:

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- Perio: Optimizes the display of periodontal tissues.
 - Endo: Optimizes the contrast values over the entire range (by default).
 - Dentin-to-enamel Junction: Optimizes the values at the crown, the amelo-dentinal junction and the roots.
- Sharpness filter Image filtering to increase image contrast applied to acquired images.
- FMS templates FMS template options to select for acquisition.
- You can modify existing templates and create new templates (see the Carestream Dental Imaging Software Online Help).
- Sensor activation duration (minutes): N/A for RVG 5200 and RVG 6200.

You can select your preferences before starting to acquire images. If you try to change the FMS template after you have finished your image acquisitions you are prompted with a warning that indicates that you risk losing some of the images.


Figure 5: FMS Retake Image Gallery

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
A. Is the FMS retake image gallery and only displays the images acquired for the frame highlighted in blue in the FMS template .

B. A blue circle in the corner of the FMS frame indicates that there are retake images for this specific frame. The images are automatically saved unless you want to select and delete them.

C. is the FMS image enhancement toolbar and applies either to a single selected frame (highlighted in blue) or to the entire FMS template.

-  Perio:Optimizes the display of periodontal tissues.
-  Endo:Optimizes the contrast values over the entire range.
-  Dentin-enamel Junction:Optimizes the values at the crown, the amelo-

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dentinal
junction and roots.

Sharpness enhancement: Optimizes the clarity of the image.

Brightness enhancement: Optimizes the brightness of the acquired image.

Contrast enhancement: Optimizes the contrast of the acquired image.

Refresh: Resets to the initial state of the current image.

Important: All enhancements applied to images as well as all retake images are transferred to the Carestream Dental Imaging Software when you close the FMS Acquisition interface.

Setting Up

1. Insert the Carestream Dental Imaging Software DVD-ROM (1/2) in the DVD-ROM drive and install the software (see the Carestream Dental Imaging Software documentation).
2. Insert the Drivers DVD-ROM (2/2) in the DVD-ROM drive. The Installer Language dialog box is displayed.

3. Select the installer language and click OK. The Welcome to the Drivers Setup Wizard dialog box is displayed.

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

5. Select RVG Family and click Install. The RVG Driver Setup Welcome dialog box is displayed.

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

6. Click Next. The Choose Components dialog box is displayed.

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

7. Select RVG 5200/6200 and click Install. The Installing dialog box is displayed.

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
When the installation is complete, the Completing the RVG Setup Wizard dialog box is displayed.

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
8. Click Reboot now, then click Finish. The computer shuts down and restarts.
9. Connect the RVG sensor to a USB 2.0 port on the back of the workstation.


The Found New Hardware Wizard is displayed.

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
10. Click Install the software automatically (Recommended), then click Next.


11. Click  on your desktop to start the RVG 5200/6200 Service Tools. The Service Tools home page is displayed.

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

2. Click Post Installation Procedure then follow the instructions on screen.
You only need to run the Post Installation Procedure on one workstation for each RVG sensor.

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Important: Do NOT connect the RVG sensor before installing the driver software.

Important: Do NOT uncheck any boxes that are already checked.

Important: To ensure the best quality images, you must connect the RVG sensor to a USB 2.0 port on the BACK of the workstation.

WARNING: You must run the Post Installation Procedure to check that the installation of the sensor is correct.

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If your workstation is running Microsoft Windows XP, a warning is displayed and you need to click Continue Anyway.

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
Acquiring Single Images

To prepare the acquisition of a single image, follow these steps:

1. Select the appropriate size of RVG sensor.
2. Connect the RVG sensor to a USB 2.0 port on the back of the workstation.


3. Access the Imaging Window from the patient file. The icon  is displayed in the Imaging Window toolbar indicating that an RVG sensor is connected to the workstation and is ready for Acquisition.

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

4. Select an appropriate positioner for the region of interest and the size of the sensor.
5. Cover the RVG sensor with a disposable hygienic sleeve specifically designed for each size of RVG sensor.


6. Place the protected RVG sensor in the biteblock of the RVG sensor positioner.


7. Position the RVG sensor in the mouth of the patient depending on the region of interest.

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

8. Move the X-ray generator tube head close to the patient and align it with the tooth of the patient and the RVG sensor.

9. Select the X-ray exposure time according to the region of interest and the patient type. Follow the user instructions of your X-ray generator. The following tables provide guidelines for exposure times for an X-ray generator at 70 kV and 7mA. Add your values for the exposure time in seconds in the column on the right.

Table 1 ADULT exposure times

Acquisition Mode	Suggested Exposure Time in Seconds	Your Exposure Time in Seconds
Upper incisor/canine	0,18	
Upper premolar	0,24	

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To acquire a single image, follow these steps:

1. Tell the patient to remain still.
2. Position yourself either two meters behind the X-ray generator or outside the door.

3. Make sure that the  is displayed in the Dental Imaging Software toolbar indicating that the RVG sensor is connected to the workstation and is ready for acquisition.

4. Trigger the X-ray with the remote control of the X-ray generator. The image is immediately displayed in the Imaging Window.

5. Check the image and if the quality is:

- Not satisfactory, for example, if the exposure quality indicator is red, retake the X-ray.

- Satisfactory, remove the X-ray generator tube head.

The ideal image quality is when the indicator of the quality of exposure is a full green bar in the

Control Panel. Avoid under-exposed or over-exposed images indicated by a partial or full red bar.

Under-exposed Image	Ideal Image Quality	Over-exposed Image
---------------------	---------------------	--------------------

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

6. Remove the RVG sensor from the mouth of the patient.
7. Remove the hygienic sensor protection and throw it away.
8. Clean and disinfect the RVG sensor after each patient .

-  Important: To ensure the best quality images, you must connect the RVG sensor to a USB 2.0 port on the BACK of the workstation.
-  Important: Use a NEW hygienic sleeve for each new patient to prevent cross-contamination.
-  Important: Always insert the RVG sensor holding it horizontally for the comfort of the patient.
-  Important: Make sure that the tube head is not shaking.
-  Important: These are suggested exposure times and need to be adjusted for your specific X-ray generator. For dark images, reduce the exposure time and for grainy images, increase the exposure time.

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Important: Make sure you can keep visual contact with the patient during the X-ray.


Important: Do NOT pull the RVG sensor by its cable when you remove the hygienic protection.

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
Acquiring FMS Images

To prepare the acquisition of FMS images, follow these steps:

1. Select the appropriate size of RVG sensor .
2. Connect the RVG sensor to a USB 2.0 port on the back of the workstation.


3. Access the Imaging Window from the patient file.

4. Click  in the Imaging Window to access the RVG FMS Acquisition interface.  is displayed in the RVG FMS Acquisition interface indicating that the RVG sensor is connected to the workstation and is ready for acquisition.

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


5. Select an appropriate positioner for the region of interest and the size of the RVG sensor.
6. Cover the RVG sensor with a disposable hygienic sleeve specifically designed for each size of RVG sensor.

7. Place the protected RVG sensor in the biteblock of the RVG sensor positioner.
8. Position the RVG sensor in the mouth of the patient depending on the region

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

9. Move the X-ray generator tube head to the patient and align it with the tooth of the patient and the RVG sensor.

10. Select the X-ray exposure time according to the region of interest and the patient type. Follow the user instructions of your X-ray generator. The following tables provide guidelines for exposure times for an X-ray generator at 70 kV and 7mA. Add your values for the exposure time in seconds in the column on the right.

Table 3 ADULT exposure times

Acquisition Mode	Suggested Exposure Time in Seconds	Your Exposure in Seconds
------------------	------------------------------------	--------------------------

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 Important: To ensure the best quality images, you must connect the RVG sensor to a USB 2.0 port on the BACK of the workstation.

 Important: Use a NEW hygienic sleeve for each new patient to prevent cross-contamination.

 Important: Always insert the RVG sensor holding it horizontally for the comfort of the patient.

 Important: Make sure that the tube head is not shaking.

 Important: These are suggested exposure times and need to be adjusted for your specific X-ray generator. For dark images, reduce the exposure time and for grainy images, increase the exposure time.

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
To acquire FMS images, follow these steps:

1. Tell the patient to remain still.
2. Position yourself either two meters behind the X-ray generator or outside the door.

3. Make sure that  is displayed in the RVG FMS Acquisition interface indicating that the RVG sensor is connected to the workstation and is ready for acquisition. The selected frame for the image is highlighted in green.

4. Trigger the X-ray with the remote control of the X-ray generator.
 - The image is displayed in the preview screen of the RVG FMS Acquisition interface.
 - The next frame is automatically highlighted in green, ready for the next acquisition.
5. Continue acquiring images until the whole FMS template is completed.
6. Check the image and if the quality is:
 - Unsatisfactory, for example, if the exposure quality indicator is red, retake the X-ray (see

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

“Retaking Images” on page 24).

- Satisfactory, remove the X-ray generator tube head.

The ideal image quality is when the indicator of the quality of exposure is a full green bar in the Control Panel. Avoid under-exposed or over-exposed images indicated by a partial or full red bar.

Under-exposed Image	Ideal Image Quality	Over-exposed
		

7. Remove the RVG sensor from the mouth of the patient.
8. Remove the hygienic sensor protection and throw it away.
9. Clean and disinfect the RVG sensor after each patient.

 Important: Make sure you can keep visual contact with the patient during the X-ray




 Important: Do NOT pull the RVG sensor by its cable when you remove the hygienic protection.

Retaking Images

If you need to retake images either while you are going through the FMS template acquisition sequence or after the FMS template acquisition is finished, follow these steps:

1. Click the frame you want to retake images.
2. Check the image quality in the preview screen.
3. Retake another image or images if you are not satisfied with the image quality.

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

The FMS retake image gallery «A» only displays the images acquired for the frame highlighted in blue in the FMS template «B». A blue circle in the corner of the FMS frame indicates that there are retake images for this specific frame. The images are automatically saved unless you want to select and delete them.

4. Select an image and apply image enhancement. The image enhancements are automatically saved.

5. Exit the RVG FMS Acquisition interface when you have finished all the acquisition sequences.

The FMS template with the acquired images and the applied image enhancements are saved and

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displayed in the Imaging Window.
The retake images are also saved in the Imaging Window but not as a part of the FMS template.

- 6. Remove the X-ray generator tube head.
- 7. Remove the RVG sensor from the mouth of the patient.
- 8. Remove the hygienic sensor protection and throw it away.
- 9. Clean and disinfect the RVG sensor after each patient.

If you retake images before finishing all of the FMS template, you have interrupted the automatic acquisition sequence. To relaunch the automatic acquisition, click the next frame in the acquisition sequence.

Important: Do NOT pull the RVG sensor by its cable when you remove the hygienic protection.

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Troubleshooting

Malfunction	Possible Cause and Action
After triggering the X-rays, no image is displayed.	- Make sure that  is displayed in the Dental Imaging Software toolbar indicating that an RVG sensor is connected to the workstation and is ready for acquisition. - Make sure that  is displayed in the RVG FMS Acquisition interface indicating that an RVG sensor is connected to the workstation and is ready for acquisition. - Make sure the RVG sensor is correctly aligned with the X-ray generator. - Make sure the settings of the X-ray generator are correct. - Make sure that the RVG sensor is connected to a back of the workstation.

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The image is pale and grainy.	• The exposure time is too short; increase it . • The X-ray generator is too far from the patient with respect to the selected dose. • Check the contrast and brightness settings of the monitor and ensure there are no reflections on the screen. • The X-ray generator voltage is too low; have the X-ray generator checked.
The image is too dark.	• The exposure time is too high; lower it. • Check the monitor settings (contrast and brightness) and ensure there are no reflections on the screen.
The image is blurred.	• Patient moved during exposure. • X-ray generator head was not stable. • Use an image filter to enhance the contrast.
The image is white.	• The image is white. • X-ray dose is insufficient. • Ensure the X-ray generator is producing X-rays; have it checked by a certified technician.

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APPLICATION

Install the unit in the X-ray study, which corresponds to the current standards for installation.

With the installation location you should be able to maintain visual and audio contact with the patient

and to have access to the interface module receiving images during the exposure.

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LEGEND OF LABEL

Symbols	Means
	Type BF device symbol complying with the IEC 60601-1 standard.
	In the European Union, this symbol indicates: Do NOT discard this product in a trash receptacle; use an appropriate recovery and recycling facility. Contact your local sales representative for additional information on the collection and recovery programs available for this product.
	Warning
	Refer to instruction manual/booklet
	Manufactured date.
	Manufacturer's address.
	Note

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**THE MATERIALS FROM WHICH IT IS MADE OR ITS MEDICAL
DEVICE BASIC PARTS**

Component	Material
RVG 5200 size 1 sensor; RVG 5200 size 2 sensor; RVG 6200 size 1 sensor; RVG 6200 size 2 sensor,	Polysulfone UDEL P-1700 BK 937
Toothbrush type holders for size 1, vertical applications for left side positioning; Toothbrush type holders for size 1, horizontal applications for right side positioning; Toothbrush type holders for size 2, horizontal bitewing; Toothbrush type holders for size 2, PA right side;	Polysulfone UDEL P-1700 BK 937 Iron rod A316L

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- DO NOT cold sterilize

* Note: Plastic parts have a limited life and should be replaced periodically. Any method of sterilization will shorten the life of plastic parts.

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TECHNICAL SUPPORT AND REPAIR

Toothbrush holders	1. Remove any residue with hot water and soap
Bite blocks	
RINN Arm & ring	1. Disassemble the metal arm and the plastic ring. 2. Remove any residue from components with hot water and soap.

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
Robert C. Meagher
Senior Director
International Regulatory Affairs
Carestream Health, Inc.
150 Verona Street
Rochester, New York 14608
Telephone 585-627-6528

STATE OF NEW YORK
COUNTY OF MONROE

Signed and attested before me on March 10, 2015 by

Notary Public

CARESTREAM HEALTH, INC.
150 VERONA STREET
ROCHESTER, NEW YORK 14608

RACHEL MARSHALL
NOTARY PUBLIC-STATE OF NEW YORK
No. 01MA6240276
Qualified in Monroe County
My Commission Expires May 02, 2015

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