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
PIP and Swap

SCROLL   

Selecting a secondary input will display a small image of the secondary source in the upper right corner of the display.

Secondary Image Size Control

The size of the secondary image is controlled by pressing  or  buttons. Pressing the  button will cycle through:

- Small PIP = Width of secondary image is 25% of total screen width.
- Large PIP = Width of secondary image is 40% of total screen width.
- Split-Screen = Primary and secondary have equal width side by side (half of screen width).
- Split-Screen Overcan = Primary and secondary have equal width with overcan applied.
- Full Screen Primary = No secondary image displayed.

The images at the bottom of the page illustrate the above sequence.

Pressing the  button will cycle through the above sizes in reverse order.

Pressing the **SCROLL** / **SWAP** button will exchange the primary and secondary inputs, and exchange their respective locations on the display. Pressing the  button a second time will restore the inputs to their original primary/secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

Small PIP

Large PIP

Split Screen

Split-Screen Overcan

Full Screen Primary


Primary and Secondary input selection is described on the front of this card.

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


PIP and Swap

Selecting a secondary input will display a small image of the secondary source in the upper right corner of the display.

Secondary Image Size Control

The size of the secondary image is controlled by pressing ◀ or ▶ buttons. Pressing the ▶ button will cycle through:

Small PIP = Width of secondary image is 25% of total screen width.

Large PIP = Width of secondary image is 40% of total screen width.

Split-Screen = Primary and Secondary have equal width side by side (half of screen width).

Split-Screen Overscan = Primary and Secondary have equal width with overscan applied.

Full Screen Primary = No secondary image displayed.

The images at the bottom of the page illustrate the above sequence.

Pressing the ◀ button will cycle through the above sizes in reverse order.

Pressing the SCROLL / SWAP button will exchange the primary and secondary inputs, and exchange their respective locations on the display. Pressing the button a second time will restore the inputs to their original primary/secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

Primary and Secondary input selection is described on the front of this card.

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Split Screen = Primary and Secondary have equal width side by side (half of screen width).

Full Screen Primary = Secondary image is equal width with primary image.

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Pressing the SCROLL / SWAP button will exchange the primary and secondary inputs, and exchange their respective locations on the display. Pressing the button a second time will return the inputs to their original primary/secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

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
Secondary Image Size Control
The size of the secondary image is controlled by pressing the **ROLL** button. Pressing the **ROLL** button will cycle through:

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Full Screen = Primary and Secondary have equal width with overscan applied.

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Pressing the **ROLL** button will cycle through the above sizes in reverse order.

Pressing the **SWAP** button will exchange the primary and secondary inputs, and exchange their respective locations on the display. Pressing the **ROLL** button a second time will restore the inputs to their original primary/secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

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Pressing the ◀ button will cycle through the above states in reverse order.

Pressing the SCROLL / SWAP button will exchange the primary and secondary inputs, and exchange their respective status. Pressing the button will return the primary and secondary inputs to their original primary/secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

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The size of the secondary image is controlled by pressing the $\leftarrow$ or $\rightarrow$ buttons. Pressing the $\rightarrow$ button will cycle through:

Small PIP = Width of secondary image is 25% of total screen width.
 Large PIP = Width of secondary image is 40% of total screen width.
 Split Screen = Primary and Secondary have equal width side by side (half of screen width).
 Full Screen = Primary and Secondary have equal width with secondary image displayed.

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Pressing the $\leftarrow$ button will cycle through the above sizes in reverse order.

Pressing the SCROLL / SWAP button will exchange the primary and secondary inputs, and exchange the respective locations on the display. Pressing the button a second time will restore the inputs to the original Primary/Secondary status. It is not necessary for both images to be displayed in order to swap primary and secondary images.

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
Secondary Image Size Control
The size of the secondary image is controlled by pressing $\leftarrow$ or $\rightarrow$ buttons. Pressing the $\rightarrow$ button will cycle through:

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Split-Screen Overscan = Primary and Secondary have equal width with overscan around

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Pressing the **SWAP** button will exchange the primary and secondary inputs, and exchange their respective locations on the display. Pressing the button a second time will restore the inputs to their original positions. Pressing the button a third time will swap the primary and secondary images.

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