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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SIEMENS SHENZHEN
MAGNETIC RESONANCE LTD. on the annexed DOCUMENT is
genuine.

China Council for the Promotion
of International Trade

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证明 查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Siemens Shenzhen Magnetic Resonance Ltd.,
Siemens MRI Center,
Gaoxin C. Ave, 2nd
Hi-Tech Industrial Park
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

User manual Supplement for Russia

The magnetic-resonance system MAGNETOM C! with accessories

Siemens Shenzhen Magnetic Resonance Ltd.

Name: Mai Zhi Hua

Position: Regulatory Affairs Specialist

1. The name of the Medical Device.

The magnetic-resonance system MAGNETOM C! with accessories.

2. Manufacturing information of the Medical Device.

Legal manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421
www.siemens.com.cn/ssmr

Manufacturing sites:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421
www.siemens.com.cn/ssmr

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja street 96,
115093 Moscow, Russia
Tel.: +7 495 737 1252
www.siemens-healthineers.com/ru/

3. Designation of the medical device.

3.1. Indication:

MAGNETOM C! system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal and oblique cross sectional images, spectroscopic images and/or spectra, and that displays the internal structure and/or function of the head, body, or extremities.

Other physical parameters derived from the images and/or spectra may also be produced. Depending on the region of interest, contrast agents may be used. These images and/or spectra and the physical parameters derived from the images and/or spectra when interpreted by a trained physician yield information that may assist in diagnosis.

MAGNETOM C! system may also be used for imaging during interventional procedures when performed with MR compatible devices such as in-room displays and MR-safe biopsy needles.

3.2. Authorized operating personnel:

MAGNETOM C! system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to the necessary user training. This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual.

3.3. Contraindications and Possible side effects:

An MR examination is contraindicated for patients with electronic or electronically conducting implants or metals, especially those containing ferromagnetic foreign bodies. These may also be patients engaged in occupations or activities which may cause accidental lodging of ferromagnetic materials or metal fragments from military activities.

Contraindications for MR examinations are:

- *Electronic implants: Pacemakers, insulin pumps*
- *Artificial heart valves*
- *Aneurysm clips*
- *Metal splinters in the eye (danger of retinal detachment)*
- *Artificial anus (anus praeter) with magnetic closure*
- *Transdermal drug patches with metallic backings*
- *Electrically conductive implants and prostheses*
- *Metallic spirals for contraception (IUDs = Intra Uterin Devices)*
- *Transdermal or other similar implants for body piercings as well as magnetic piercings*

3.4. Physiological effects:

Due to the presence of strong alternating electromagnetic fields, patients may experience various physiological effects during MR examinations:

- *Peripheral nerve stimulation through low-frequency fields of the gradient coils*
- *Warming of body tissue through RF fields of the RF-transmitter coil*


It is generally accepted that no published evidence supporting the occurrence of cumulative and/or long-term effects after exposure to EMF emitted by the MR equipment exists.

3.5. Conditions for use:

The MAGNETOM C! MR system must be operated according to the indication and only qualified persons with the necessary knowledge in accordance with country-specific regulations, e.g. physicians, trained radiological technicians or technologists, subsequent to necessary user training.

This user training must include basics in MR technology as well as safe handling of MR system. The user must be familiar with potential hazard and safety guidelines the same way the user is familiar with emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the Operator Manual.

3.6. System's overview

4. Class

4.1.

Classific

Class

Classifi

135130

4.2.

The
Class


4.3.

Acc

4.4.

Th

4.

4. Classification of the Medical Device:

4.1. Type of the medical device in accordance with the nomenclature for classification:

Classification according to Annex IX, MDD

Class

Ila

Classification in accordance with Russia Nomenclature classification:

135130

4.2. Class of potential Risk of the medical device application:

The Medical Device is intended for transmission of energy in the form of ionizing radiation, it belong to the Class 2a.

4.3. Class of the Medical Device in accordance with IEC 60601-1:

According to IEC 60601-1:2005, the Medical Device is the Class 1.

4.4. Type on contact with the human body:

The Medical Device is non-sterile, nonexpendable.

4.5. Software classification:

The Medical Device's Software safety Class according to IEC 62304:2006 is **Class B**.

4.6. **Applied parts:**

According to IEC 60601-1 Cl. 7.2.10 and 8.3 the following parts are applied parts of following types:

- Patient table >>> **Type B** applied part;
- Transmit/receive Body Coil >>> **Type B** applied part;
- Head Array Coil >>> **Type BF** applied part;
- Neck Coil medium >>> **Type BF** applied part;
- Neck Coil large >>> **Type BF** applied part;
- Extremity Array Coil small >>> **Type BF** applied part;
- Extremity Array Coil large >>> **Type BF** applied part;
- Extremity Array Coil XL >>> **Type BF** applied part;
- Shoulder Array Coil >>> **Type BF** applied part;
- 4-channel Wrist Array Coil >>> **Type BF** applied part;
- Breast Array Coil >>> **Type BF** applied part;

The system is not classified as defibrillator-proof applied part (IEC 60601-1:2012 Clause 8.5.5).

4.7. **SanPiN 2.1.7.2790-10 class A**

4.8. **IP protection class:**

IP protection class according to IEC 60529 >>> IP X₀

4.9. **Laser positioner:**

Laser light Class: **2M** (according to IEC 60825 - 1:2007).

4.10. **Operating mode:**

Operating conditions according to IEC 60601-1.

System: Continuous mode of operation with short-term load

5. Technical Data: Specifications and Characteristics.

5.1. **Technical Data:**

MAGNET SYSTEM

Magnet Parameters	
C-shaped permanent magnet with vertical magnetic field	
Operating field strength	0.35 Tesla (-0.0094 T ~ -0.0047 T)

Patient Aperture	
The open design provides free patient access from three sides	
Upper magnet width	137 cm (± 1 cm)
Vertical gap distance (cover)	41 cm (± 5 mm)
Table top upper pole cover	38.5 cm (± 5 mm)

Hybrid Shim System	
Provides maximum homogeneity of the magnetic field	
Preshim	Factory preshim for short installation time
Passive shim	Computer based passive shim to adapt to environmental static interferences
Active shim	The 3 gradient coils are used for shimming the linear terms

Volume Selective 3D Shim	
Patient-specific automated online shim, for optimization of the magnetic field homogeneity within the imaging volume	
Time to shim	approx. 20 sec

Isocenter Imaging	
All extremity joints are scanned in the magnet isocenter for optimum results	

Homogeneity ¹		
DSV	Typical	Guaranteed
36 cm	± 2.5 ppm	± 4 ppm (FWHM)
20 cm	± 1.2 ppm	± 2.5 ppm (FWHM)

Fringe Field ²		
X-axis: horizontal	(in front of the magnet)	2.1 m (-0.2m/+0)
	(behind the magnet)	3.0 m (-0.2m/+0)
Z-axis: horizontal	(left / right)	2.3 m (-0.2m/+0)
Y-axis: vertical	(top of the magnet)	3.2 m (-0.2m/+0)

	(bottom of the magnet)	2.7 m (-0.2m/+0)
--	------------------------	------------------

1. DSV=Diameter spherical volume (x, y, and z direction) FWHM=Full Width Half Maximum of the proton spectrum measured in 13 planes with 12 points per plane.
2. Distances of the 0.5 mT line (pacemaker safety limit) from the magnet isocenter.

The respective

P
A
T
T
V

PATIENT HANDLING

The respective configuration is set by Siemens Service during MR system installation.

Patient Table	
Allows fast and flexible patient positioning outside the magnet gap	
Table dimensions (LxW)	210x79 cm (±1mm)
Table height (fixed)	81.5 cm (±5mm)
Longitudinal travel (manual)	± 83 cm (±5mm)
Lateral travel (manual)	- 81.8/+ 21.2 cm (±5mm)
Position accuracy	± 1 mm
Max. patient weight	200 kg

Patient Positioning Aids
Standard set of cushions, sandbags and restraining straps for comfortable and stable patient positioning
Arm rest for comfortable positioning of the arms
Optional set of vacuum cushions, 'Comfort Kit' (three sizes)

Patient Positioning Tools
Laser light localizer for positioning of center slice
Manual transfer to center of imaging volume
Positioning tools for comfortable and stable positioning of the patient

Patient Communication
Assistance call squeeze-bulb for patient
Ergonomically designed patient communication unit - may be placed at any convenient location on the desktop
- Intercom communication with the patient - Response to the patient's activation of the squeeze-bulb - Volume of speaker in control room - Connection to hi-fi stereo system - Independent volume control of voice

• Loudspeaker
• Microphone
Automatic voice command for breathhold examinations

NOISE-RELATED DATA

Noise-related data, patient

- Max. noise measured: 114.8 dB(A)
- Measurement reproducibility: ± 3 dB
- Legal limit: 99 dB(A)

Hearing protection specified for noise attenuation of at least 19 dB is required for patients (114.8 dB(A) + 3 dB – 99 dB(A)).

Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4

Noise for personal in the examination room

Noise measured : 84.4 dB(A) (± 3 dB)

Noise in the operator room (at the syngo Acquisition Workplace)

The operating personnel in the operator room (syngo Acquisition Workplace) does not require hearing protection

GRADIENT SYSTEM

General features	
Actively shielded (AS), Z-axis gradient coil system	
Integrated with tx coil	
Shortest TE and TR for high image contrast and best signal-to-noise ratio	
Water-cooled coil and amplifier for maximum performance	
Gradient coil diameter	974 mm (± 1 mm)
Patient-specific shimming	

Active shim via 3 linear channels

Actively shielded flat gradient coil system

Integrated into the magnet pole shoes

Gradient Amplifier	
Ultrafast solid-state technology with very low switching losses	
Max. gradient field strength	24 mT/m
Max. slew rate	55 T/m/s
Duty cycle	100%
Effective Performance	
Max. eff. gradient field strength	39.4 mT/m
Max. eff. slew rate	92.5 T/m/s
Gradient Cooling System	
Gradient coils	closed loop water cooling
Power electronics	closed loop water cooling
C! Gradient	
Max. output voltage*	400 V
Max. output current*	150 A

* Values for each of the 3 gradient axes

RF SYSTEM

Transmit amplifier	
Bandwidth	800 kHz
Peak power	2.5 kW
Gain instability	0.15 dB

Receive path	
Receiver signal resolution	32 bits
Receiver signal dynamic	up to 24 bits
Number of independent receiver channels	4

Max. output sampling rate (complex)	1 MHz
Preamplifier noise figure	<0.8 dB
Preamplifier ADC sampling rate	10 MHz
Preamplifier gain	29 dB/53 dB with automatic control
Dynamic range	128 dB with automatic control

RF COILS

Array coils (standard)	Applications
Head Array Coil	• Head imaging • MR Angiography • Combined head/neck imaging • TMJ (temporomandibular joints)
Neck coil, Large	• Cervical spine • Neck • Larynx/ Esophagus • MR Angiography • Combined head/neck imaging
Body/Spine Array, Large	• Thoracic (incl. heart) • Abdomen • Pelvis • Thighs • MR Angiography • Spine imaging

Array coils (optional)	Applications
Body/Spine Array, Small Body/Spine Array, Medium Body/Spine Array, XL	• Thoracic (incl. heart) • Abdomen • Pelvis • Thighs • MR Angiography • Spine imaging
Breast Array coil	• Simultaneous high-resolution imaging of both breasts in coronal and transverse directions

	• Uni- or bi-lateral imaging of the breasts in sagittal direction
Extremity Array Coil, Large Extremity Array Coil, Small Extremity Array Coil, XL	• High resolution imaging of the knee and ankle
Shoulder Array Coil	• High resolution imaging of the right or left shoulder with excellent visualization of small anatomical structures
Wrist Array Coil	• High resolution imaging of the wrist and adjacent region with excellent visualization of small anatomical structures
Neck coil, Medium	• Cervical spine • Neck • Larynx/ Esophagus • MR Angiography • Combined head/ neck imaging

Flexible coils (optional)	Applications
MultiPurpose coil 16	• Extremities • Neck • Interventional procedures
MultiPurpose coil 21	• Extremities • Neck • Interventional procedures
BodySpine_Array_XXL	• Thoracic (incl. heart) • Abdomen • Pelvis • Thighs • MR Angiography • Spine imaging

Dedicated coils (optional)	Applications
Cordless Coil	• Interventional procedures

SEQUENCES

Spin Echo family of sequences

- Spin Echo (SE)—Single, Double and Multi Echo (up to 32 echoes); Inversion Recovery (IR)

- 2D/3D Turbo Spin Echo (TSE)—Restore technique for shorter TR times while maintaining excellent T2 contrast; TurboIR: Inversion Recovery for STIR, DarkFluid T1 and T2, TrueIR; Echo Sharing for dual contrast TSE
- 2D/3D HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo)— Inversion Recovery for STIR and DarkFluid contrast

Gradient Echo family of sequences

- 2D/3D FLASH (spoiled GRE)—dual echo for in-/opposed phase imaging
- 3D VIBE (Volume Interpolated Breath-hold Examination)—quick fat saturation; double echo for in-phase/opposed phase 3D imaging; Inline Breast Evaluation
- 2D/3D MEDIC (Multi Echo Data Image Combination) for high resolution T2 weighted orthopedic imaging and excellent contrast
- 2D/3D TurboFLASH—3D MPRAGE; single shot T1 weighted imaging e.g. for abdominal imaging during free breathing
- 2D/3D FISP (Fast Imaging with Steady State Precession)
- 2D/3D PSIF—PSIF Diffusion
- Echo Planar Imaging (EPI)—diffusion weighted; single shot FID e.g. for Perfusion weighted imaging; 2D/3D Segmented EPI (SE)
- 2D/3D TrueFISP—Shared Phases Real-time TrueFISP (w/o ECG gating)
- ce-MRA sequence with Inline subtraction and Inline MIP
- 2D/3D Time-of-Flight (ToF) Angiography—single slab and multi slab; triggered and segmented
- 2D/3D Phase Contrast Angiography—multiple velocity encoding; triggered

COMPUTER SYSTEM

Host computer (<i>syngo</i> Acquisition Workplace)	
New state-of-the-art CPU generation	
Intuitive Windows interface supports multi-tasking for fast operation	
Double processor	Intel Xeon™
Clock rate	≥ 2.6 GHz
Main memory (RAM)	4 GB
1 st hard disk (system SW)	≥ 146 GB
2 nd hard disk (data base)	≥ 146 GB
3 rd hard disk (images)	≥ 146 GB (approx. 220000 images, matrix 256 ² and 512 ² uncompressed)
CD-R writer	approx. 4000 images 256 ² ; DICOM Standard, ISO 9660

DVD-R writer	approx. 25000 images 256 ² ; DICOM Standard, ISO 9660
Media drives	CD/DVD-ROM drive
Color LCD monitor	Screen matrix 1280 × 1024 pixels

Measurement and Reconstruction System (MARS)	
New state-of-the-art CPU generation	
CPU	1 × Dual Core Xeon™
Clock rate	≥ 2.33 GHz
Main memory (RAM)	≥ 8 GB
1 hard disk	≥ 500 GB
Reconstruction speed	≥ 983 recons per second (256 ² FFT, full FoV) ≥ 2773 recons per second (256 ² FFT, 25% recFoV)
Parallel Scan and Recon	Simultaneous scan and reconstruction of up to 4 data sets

NETWORK/DATA COUPLING

The MR system provides the level of safety according IEC 60950-1 outside patient environment.

All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

IMAGE VIEWING AND FILMING

Image recall

Images are stored in DICOM format, allowing fast image access and recall

- Combines images in a new series, e.g. post-processed images
- Loads images automatically
- Inline Display

Image display

- Various display layouts selectable
- Up to 3 patients can be simultaneously active in the viewer

- Image annotation and labeling
- Non-interpolated display
- Fast paging through up to 500 images with 15 images/s for full screen display

Windowing

- Freely selectable window width and center
- Windowing on succeeding images
- Auto-windowing for optimized contrast
- Saves and sends window values

Automatic movie for cine display

Interactive movie paging by dragging the mouse or automatic movie mode by clicking the icon

Evaluation

Parallel evaluation of up to 40 regions of interest

- Circle
- Rectangle
- Freehand ROI
- Pixel lens with position marker
- Statistical evaluation
- Area
- Standard deviation
- Mean value
- Min/max values
- Image scrolling
- Magnification
- Distance
- Angle

2D Post-processing

Image manipulations

- Reversal of gray-scale values
- Image rotation by 90° or by user-defined angle
- Flip horizontally/vertically
- Image zoom and pan

- Shutter
- Annotation

Position display

Displays measurement

Argus viewer

Viewing software

- Efficiency
- Multiple
- Single
- Rapid
- Creation

Dynamic analysis

Arithmetic operations

- Addition
- Arithmetic
- Calculation
- Difference
- Calculation
- Time
- Analysis

Several evaluation

Image Filter

The Image

Mean Curve

Time-int

Filming

- Shutter
- Annotation

Position display

Displays measured slice positions on localizer image and selected series

Argus viewer

Viewing software for cardiac MR studies and large data sets

- Efficient cine review of cardiac and other dynamic data sets
- Multiple sorting options
- Single movie as well as 2, 4, or 8 simultaneous slices together in movie mode
- Rapid avi creation of 1 to 8 slices simultaneously
- Creates and edits DICOM structured reports

Dynamic analysis

Arithmetic operations on images and series (e.g. for evaluation of contrast media studies)

- Addition, subtraction, multiplication, division of single images and whole series
- Arithmetic mean and standard deviation across a range of selected images
- Calculation of T₁ and T₂, and logarithmic images
- Differentiation/integration of selected images
- Calculation of a mean slope image from a range of selected images
- Time-to-peak evaluation (TTP)
- ADC maps

Several evaluation functions may be started consecutively in the background

Image Filter

The Image Filter software uses adaptive filtering and adjusts to the local image content automatically

Mean Curve

Time-intensity analysis for contrast-enhanced examinations

- Creates and edits DICOM structured reports

Filming

- Connection via DICOM Basic Print
- Interactive filming

- Filming parallel to other activities
- Independent scanning and documentation—no wait time due to camera delays
- Freely selectable positioning of images onto virtual film sheet
- Selectable various film layouts
- Mother-in-Child display
- Windowing, image zoom and pan on film sheet
- Configurable image text
- Simultaneous handling of multiple film jobs
- Up to 100 virtual film sheets

Printing on paper

Interface and software for printing images on paper (laser printer not included)

Grey levels and color printing supported

Data format Postscript Level 2

3D POST-PROCESSING

MPR—Multiplanar Reconstruction

Real-time multiplanar reformatting of secondary views

- Viewing perspectives: sagittal; coronal; axial; oblique; double oblique; curved (freehand)
- Reconstruction along polygon and/or curved (freehand) cut lines
- Reconstruction based on reconstructed planes possible
- Reconstruction of user-defined ranges of parallel, radial or freehand cuts
- Selectable slice thickness and slice increment of reconstructed images
- Storing of post-processing protocols
- Annotations and 2D evaluations such as distance and ROI

MIP—Maximum Intensity Projection

3D reconstructions of vessels from a 3D data set, or a

2D sequential slice data set (acquired with dedicated MR Angiography sequences)

- Volume of Interest (Vol) defined to increase reconstruction speed and to improve image quality
- Freehand MIP: virtually eliminating interfering overlays by manually drawn, irregular contours in all dimensions
- Arbitrary views along any direction defined interactively with mouse-driven virtual trackball
- Multiple view angles around any orthogonal axis

- Projection
- MIP thin

MinIP—Minimum
Similar to MIP b

SSD—Shaded
Three-dimensi

- Select
- Multi
- Rect

AMBIENT CO

SPECIFIC

The spe
to the r

INTEG

- Projections displayed as single images, as interactive movie or by fast paging
- MIP thin/MIP thick

MinIP—Minimum Intensity Projection

Similar to MIP but reconstructs the minimum intensity (e.g. for Dark Blood techniques)

SSD—Shaded Surface Display

Three-dimensional display of surfaces, such as contrast-enhanced vessels

- Selectable variable threshold values
- Multiple view angles around any orthogonal axis
- Rectangular and irregular Volumes of Interest (Vol) can be defined to improve image quality

AMBIENT CONDITIONS

Control room	
Temperature	15–30 °C
Rel. humidity	40 to 80%

Electronics room	
Temperature	18–24 °C (with integrated cabinet temperature control)
Rel. humidity	40 to 80%

Examination room	
Temperature	21–24 °C
Rel. humidity	40 to 60%

SPECIFIC ABSORPTION RATES

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

INTEGRATED COOLING SYSTEM

Heat dissipation	≤ 8.9 kW (400 V - 50 Hz)
	≤ 9.4 kW (460 V - 60 Hz)
Cooling capacity	≥ 8 kW
Water setting temperature	24.5 °C

• **LASER**

Class	2M
Radiant exposure	≤ 25 W/m ²
Wavelength	650 nm

LINE POWER SUPPLY

Variable	Values	Tolerance
Line voltage	380 V, 400 V, 420 V, 440 V, 460 V, 480 V	± 10%
Line frequency	50 Hz, 60 Hz	± 1 Hz
Connection value	23 kVA	
On-site protection	at least 63 A	

For all products, line power has to be supplied via an on-site system contactor or via another multipole shut-down mechanism. Room installation has to be in compliance with VDE 0100-710.¹

¹ In all countries outside the Federal Republic of Germany, compliance with local and national legal regulations is required. However, we strongly recommend compliance with the regulations described herein—to the extent permitted by relevant local and national laws—in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate

380V/400V/420V/440V/460V/480V~	
3PHASE	50/60 Hz
SHORT-TIME	kVA: 21
LONG-TIME	kVA: 6

POWER AND ENERGY CONSUMPTION

System off ¹	0.9 kW
Stand-by	3.5 kW
Typical measurement	4.5 kW
Highest average power	13 kW

¹ Magnet heating always on

GENERAL CLASSIFICATIONS

Protection class	I
Components used • Patient table • Local coils	Type B applied part Type B or BF applied part
IP protection class according to IEC 60529	IP Xo
Explosion protection	The MR system is not intended for operation in areas prone to explosion (e.g., highly flammable mixtures of anaesthesia gases with air or oxygen or nitrous oxide)
Operating mode	Continuous operation with short-term load
Disinfection Receptors/ system components	Disinfectants without alcohol, ether
Degree of safety in the presence of a flammable anesthetics mixture with air or with oxygen or with nitrous oxide	No AP or APG category equipment

DIMENSIONS

Component	Width [cm]	Depth [cm]	Height [cm]	Weight [kg]	Heat dissipation [kW]
Magnet (incl. cover)	137 ±1cm	209 ±1cm	189 ±1cm	16,000±100kg	1.0
Patient table	79 ±1mm	210±1mm	82*** ±5mm		
Required min. room height			225		
syngo Acquisition Workplace (table+monitor)	120±5mm	80±5mm	72+45=117±5mm		2
Host computer	22±5mm	80±5mm	49±5mm		
Electronics cabinet, incl. system control, RF system, gradient power system, image processor	95.4±1cm	64.8±1cm	197.5±1cm	602.5±10kg	10
Cooling system	80±1cm	64.8±1cm	193±1cm	310±10kg	

*)without electronics, table and table bracket

**)without attachments

***)Above floor level

5.2. Software:

MAGNETOM C! basic software information	
Name and Version:	syngo MR C15 (VC15D)
Date of release	Nov 26, 2013
Size in GB	1.40 GB

6. The list of components and accessories.

The magnetic-resonance system MAGNETOM C! with accessories			
#	Название принадлежности	Name of component	Description
BASIC BLOCK:			
1.	Магнит.	Assembly #C C-shaped permanent magnet	Open, C-shaped permanent magnet with a field strength of 0.35 Tesla with whole body RF transmit coil and gradient coil.
2.	Градиенты.	Assembly #C Gradient system engine	The whole-body gradient system engine.
3.	Интегрированная в магнит катушка для тела.	Assembly #C Transmit/receive Body Coil	Integrated no tune transmit/receive Body Coil
4.	Приемо-передающий радиочастотный тракт.	MAGNETOM C! RF system	Digital RF Transmit and Receive System for computer controlled signal processing
5.	Кожух магнита.	MAGNETOM C! Magnet cover	Magnet cover
6.	Стол пациента.	Patient table #C	Patient table freely movable in longitudinal and transverse direction for effective patient positioning and optimal

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s ents to	Катушка для тела/позвоночника - 1 шт.	Special coils kit #C	Coil set for head/neck and body imaging. Automatic coil tuning reduces preparation time to a minimum.
System: for	Катушка для конечностей (при необходимости, не более 3 шт.)	Extremity Array Coil, small #C Extremity Array Coil, XL #C Extremity Array Coil, large #C Extremity Array Coil, large #C	Phased Array Coil for imaging of the lower extremities
splays om of	Катушка для плечевого сустава (при необходимости, не более 2 шт.)	Shoulder Array Coil #C Shoulder Array Coil #C	Phased Array Coil for imaging of both Shoulders.
tings: ctors,	Катушка для отделов позвоночника/тела (п ри необходимости, не более 5 шт.)	Body/Spine Array Coil, small #C Body/Spine Array Coil, medium #C Body/Spine Array Coil, medium #C Body/Spine Array Coil, XL #C Body/Spine Array Coil, XXL #C	Receive array coil for body imaging
host display eed	Катушка для рук и запястья (при необходимости).	4-channel Wrist Array Coil #C	Auto tune phased Array receive coil, anatomically adapted to the shape of the wrist for imaging of the hand and proximal regions.
ge of	Набор беспроводных поверхностных катушек (при необходимости, не более 3 шт.)	MP2 Coil Set, small #C	Coil kit consisting of two semiflexible ring coils and an adapter for imaging of the peripheral joints and the neck region.
and ies	Катушка для молочных желез (при необходимости).	Breast Array Coil #C	iPAT compatible coil for high SNR imaging of the breasts. Array coil with 4

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			of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.
4.	Программный пакет трехмерной ортовизуализации.	CISS & DESS #C CISS & DESS	Advanced 3D imaging sequences and protocols which are unique to Siemens. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.
5.	Программный пакет расширенной васкулярной визуализации.	Advanced Angio Package #C	Advanced MR Angiography sequences and protocols optimized for the head/neck, body and peripheral vasculature. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.
6.	Программный пакет динамического просмотра.	Argus viewer	Argus viewer can be used to display cine studies. Argus Viewer allows users to load a large list of dynamic data sets and view them comfortably. This is a feature that greatly reduces the reading and review time for MR cardiac studies. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number

version number , are basic syngo MR 11/02/2013			and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
ng protocols to nal. n installed ware and ne request on number are ic syngo MR 11/02/2013 raphy ocols	Программное обеспечение для оптимизации использования контрастного вещества.	Care Bolus syngo #C	Care bolus technique for real-time visualization of contrast medium enhancement on the integrated online display with the possibility of switching from a 2D sequence to 3D MRA sequence interactively with a single mouse click. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
ure. l. talled re and request n number re go MR 02/2013 used to . Argus to load nic nem a nd rdiac lled and request number	Программный пакет кардиовизуализации.	Advanced Cardiac Package #C	The package contains special sequences and scan protocols for MR studies of the heart. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
	Программный пакет интерактивной визуализации в реальном масштабе времени со специальной электронной «МЫШЬЮ».	MAGNETOM C! Inline software	Inline and offline image filter software for noise suppression in the MR images. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR

10.	Программный пакет адаптивной фильтрации.	MAGNETOM C! Image Filter	C15 (VC15D) from 11/02/2000 The MAGNETOM C! integrates all aspects of the clinical MR exam into an efficient, cost-effective environment. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo C15 (VC15D) from 11/02/2000
11.	Стол оператора, не более 4 шт.	Table syngo 1,2m #C	The table is fixed to the magnet system and has ergonomically designed holding grips for moving the table top in longitudinal or transverse direction. The table comes with a comprehensive set of patient positioning cushions.
12.	Контейнер к столу оператора, не более 3 шт.	Office Container syngo #C	614 mm wide extra case for the syngo host computer with sliding front door to allow change of storage media.
13.	Подушки с песком для комфортной укладки пациента - 4 шт.	Comfort Kit #C	Vacuum cushions for positioning and stabilization of the patient during the examination. The comfort kit may significantly reduce patient set-up times and optimize image quality by minimizing the occurrence of motion artifacts caused by involuntary patient motion.
14.	Жидкокристаллический монитор.	LCD Monitor	High-resolution LCD color flat screen monitor, with integrated gamma correction for optimum display of MR grayscale images and automatic backlight control.

11/02/2011			for long-term brightness stability.
integrates the clinical, efficient, environment onal. installed software and the request sion number are basic syngo MR 11/02/2011 d to the has igned moving ransverse ble comes ve set of g case for puter door to orage r ne e s and ity by rence of used by motion. color with rection of MR d control	Блок измерения физиологических параметров.	PMU Wireless Physio Control #C PMU Wireless Physio Control #C	Physiological Measurement Unit (PMU) - Wireless Physio Control for wireless triggering, synchronizes the measurement with the physiological cycles of cardiac and/or respiratory motion. Wireless technology for all sensors allows fast and easy patient set-up and comfort, and robust cardiac or respiratory signal transmission as it eliminates the need to attach cables to the patient. The Wireless Physio Control contains wireless VCG, respiration and pulse sensors and a charging station as all sensors are powered by rechargeable batteries.
	Специальный шкаф с полками для хранения катушек.	Coil storage cart	Specially designed non-ferromagnetic cart for easy storage of the most commonly used coils and accessories.
	Модуль для теленаблюдения за пациентом.	Patient supervision	Module for monitoring the patient during an MR examination, conveniently mounted on the wall of the examination room.
	Примагнитная консоль оператора (компьютерный блок, программное обеспечение, монитор, клавиатура, электронная «мышь» и тележка).	syngo in-Room Acquisition Workplace	syngo in-Room Acquisition Workplace: including host computer, LCD color display and separate high speed image processor
	Программный пакет функциональной	Argus 4D Ventr. Function syngo	Argus 4D Ventricular Function:

	оценки сердца.	#C	- Calculation of volume cardiac data of a given patient very quickly and easily. - Parametric results and volume-time curves are calculated upon automatic creation and adaptation of a 4D model of the left ventricle. The resulting 4D model of the patient's heart can be visualized superimposed on anatomical images as a reference. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.	Програ трехме цветов отобра анатом объем
20.	Программный пакет анализа кровотока.	Argus 4D Ventr. Function syngo #C	Argus Function: - Automatic, semi-automatic, or manual segmentation of the left and right ventricle. - Volumetric analysis and wall thickness analysis. - Output of parametric results, volume-time curves and bull's-eye plots. - DICOM Structured Reporting. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.	Прог трех мани изоб реал вре
21.	Программный пакет оценки динамических обследований.	Dynamic Analysis	Standard software permits advanced evaluations: addition/subtraction.	24. Пр

			ction, division/ multiplication, T1 and T2 calculations, and standard deviation. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
2.	Программный пакет трехмерного цветового отображения анатомических объемов.	3D VRT syngo	This application provides a dedicated evaluation software for colored, volume rendered (VRT) visualization of MR and CT data as a complement to traditional MIP display. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
3.	Программный пакет трехмерных манипуляций с изображениями в реальном масштабе времени.	3D Post- Processing Card	- Manual MPR (Multi-Planar Reconstruction) - MIP (Maximum Intensity Projection) - MinIP (Minimum Intensity Projection) - SSD images (Shaded Surface Display) Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
4.	Программный пакет	Fly Through	Integrated syngo image

	виртуальной эндоскопии.	syngo	analysis is a vital component of the MR examination.\nThis application provides dedicated evaluation software for simulation virtual endoscopy or bronchoscopy and for Fly-Through in vascular structures, where real endoscopic procedures are impossible. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.	оценки
25.	Программный пакет наложения серий изображений.	Image Fusion syngo	This application provides dedicated evaluation software for spatial alignment (matching) and visualization of image data either from different modalities (CT, MR, NM, PET) or from the same modality but from multiple examinations of the same patient. It supports optimal diagnostic outcome (fusion of morphological and functional information and therapy planning. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2011.	7. Прог комп спин изоб
26.	Программный пакет количественной	Vessel View syngo	Integrated syngo image analysis is a vital	28. Пр ав к

ital e MR his vides tion ulation of y or for ascular e real dures onal. ninstalled ware and he request ion number are sic syngo MR 11/02/2013 provides a on al g) and image data ent (NM, PET) odality e same ts outcome ogical ormation) g. l. talled re and request n number re C go MR 02/2013 ge	оценки сосудов.	Vessel View AVS	component of the vascular MR examination. syngo Vessel View provides a dedicated evaluation software for automated image analysis of MR and CT angiography data for small and large vessels. The combination of automated segmentation and easy-to-use editing tools, provides users with a rapid way to extract vessels and quantify and grade vessel lesions. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
	Программный пакет композиции спинальных изображений.	Spine Composing syngo #C Spine Composing syngo #C	This application provides dedicated evaluation software for creating fullformat images from overlapping MR volume data sets of the spine acquired at multiple stages. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
	Программный пакет авторизации доступа к системе и данным.	syngo Security Package	Software option for general regulatory security rules, providing enhanced security features including user management and audit trail functionality. Software is optional.

			Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2015.
29.	Программное обеспечение просмотрной станции СИНЕТ.	Software viewing workstation	Software viewing workstation Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2015.
30.	Программное обеспечение для ДИКОМ-документирования СИНЕТ.	Software DIVOM documentation workstation	Software DIVOM documentation workstation Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2015.
31.	Программное обеспечение для записи на компакт-диски СИНЕТ.	Software CD recording workstation	Software CD recording workstation Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2015.
32.	Программный пакет трехмерной обработки и отображения поверхностей.	Multi Channel Application Suite #C	The Application Suite offers a complete range of clinically optimized sequences, protocols and workflow functionalities for all body regions.

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с power and moving components in the room.	минераловолоконные вставки, электровыключатели)		
Probe for internal field to move the for the M	Панель шумопоглощающая для отделки радиочастотной кабины.	Filter plate for RF-Cabin	Filterplate in low leakage circuitry for shielded rooms (MRI) 6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A
cher	Роллеры для сейсмозащиты, не более 8 шт.	Seismic Kit	Seismic Kit comprises: - tufnol spacers which have hole offsets increments and washers - radial gap checking tool - vertical gap setting tool - tufnol spacer kit - loctite - washers
n for			
GXT4 regarding	Контур для охлаждения магнита.	Contour of cooling system.	Contour of cooling system.
	Климатотехническая установка.	Cooling system	Cooling system
t	Дополнительные подкладки для пациента, не более 20 шт.	Comfort Kit #C	Vacuum cushions for positioning and stabilization of the patient during the examination. The comfort kit may significantly reduce patient set-up times and optimize image quality by minimizing the occurrence of motion artifacts caused by involuntary patient motion.
cabin, building. for the RF			
er			
late in.	Опора для позиционирования рук пациента.	Patient hands positioning support	Patient hands positioning support
cor g for	Специальные подушки для улучшения качества визуализации, не более 5 шт.	Comfort Kit #C	Vacuum cushions for positioning and stabilization of the patient during the examination. The comfort kit may significantly reduce patient set-up times and optimize image quality by minimizing the occurrence of

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caused by ent motion	расширенной количественной оценки перфузии.	Evaluation syngo	syngo provides a task card for detailed post- processing of brain perfusion data sets. Color display of the relative Mean Transit Time (relMTT), relative Cerebral Blood Volume (relCBV), and relative Cerebral Blood Flow (relCBF) is supported. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013
er d for and nsensitive pin Echo body onal al. stalled are and e request on number are ic syngo MR 1/02/2013 PU and aphy cols	Программный пакет удаленного доступа к управлению сканером.	syngo Expert-i #C	This software application enables support by "remote staff" for planning and processing. Software is optional. Software is pre-installed on the basic software and is activated at the request of the user. Version number and release date, are similar to the basic software, namely syngo MR C15 (VC15D) from 11/02/2013

7. Capabilities of the Medical Device and ways of integration with other medical devices.

Name of the device	Manufacturer
Color display HL1916	Shenyang Torch-Bigtide
Host R640	Fujitsu Siemens
Intercom	SASSE
Cordless Coil	RAPID
Cushions	PolyForm
Cushion Vakuum Satz	Pi-Products
Vakuum Cushion	Pi-Concept
Respirationskissen 005	Pi-Products

Sprechanlage	Sasse
C2! ICS	Airsys
Host xw4600	HP
Remote Color Display 19"	Siemens A&D SE DT/ EIZO GmbH
z400 Computer Workstation	HP
Color display DSC1910-D/19" Monitor(10656056)	Shenyang Torch-Bigtide
DVI_Splitter_DS01	Shenyang Torch-Bigtide
Host Z420(10928295/10864003)	HP

8. Labeling information.

8.1. At electronic cabinets.


Name plate label:

Gradient system label:


Gradient Coil label:

Magnet label:

PERU ECG Respiratory module label:

Pulse module label:

Pulse sensor label:

Body Spine A

Body Spine


Body Spine Array L label:


Body Spine Array M label:

Body Spine Array S label:

Body Spine Array XL label:

Body Spine Array XXL label:

Body Spine Array XXL C035

(01) 04056869038766	GTIN
(240) 08635505	MODEL
(422) 156	ORIGIN CN
(21) 1153	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, HI - Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Body spine lower mr coil label:

BODY SPINE LOWER MR COIL

(01) 04056869049854	GTIN
(240) 10352808	MODEL
(422) 156	ORIGIN CN
(21) 5633	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, HI - Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Breast coil label:

Breast coil c035

(01) 04056869049564	GTIN
(240) 08635653	MODEL
(422) 156	ORIGIN CN
(21) 1141	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, HI - Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Extremity Array Coil L label:

Extremity Array Coil L

	(01) 04056869038827	GTIN
	(240) 08635588	MODEL
	(422) 156	ORIGIN CN
	(21) 1978	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, HI - Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Extremity Array Coil S label:

Extremity Array Coil S

	(01) 04056869038834	GTIN
	(240) 08635596	MODEL
	(422) 156	ORIGIN CN
	(21) 1598	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, HI - Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Extremity Array XL label:


Ext


Head Array l

Multipur


Multit

Head Array label:

Multipurpose Coil 16S label:

Multipurpose Coil 21 label:

Multipurpose Coil 21 C035

	(01) 04056869038858	GTIN
	(240) 08635612	MODEL
	(422) 156	ORIGIN CN
	(21) 1613	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Neck Coil L label:

Neck Coil L C035

	(01) 04056869038803	GTIN
	(240) 08635562	MODEL
	(422) 156	ORIGIN CN
	(21) 1373	SERIAL


IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Neck Coil M label:

Neck Coil M C035

	(01) 04056869038872	GTIN
	(240) 08635679	MODEL
	(422) 156	ORIGIN CN
	(21) 2106	SERIAL

IVK

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Shoulder Array Coil label:

Wrist Array

GTIN
Model
ORIGIN
SERIAL
Made in
Made in

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
AIRSYS Refrigeration Engineering Technology (Beijing) Co., Ltd.
East to LuGu Street ShijingShan District, Beijing 100040, P. R. China
Tel: +86-10-68656161 Fax: +86-10-68652453

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AIRSYS Refrigeration Engineering Technology (Beijing) Co., Ltd

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AIRSYS Refrigeration Engineering
Technology(Beijing)Co., Ltd.
East to LuGu Street ShiJingShan
District, Beijing 100040, P.R. China
Tel:+86-10-68656161 Fax:+86-10-68652453

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

Labels according to EN 50419:2004:

Within the EU, products identified with this symbol are subject to guidelines 2002/96/EG for old electrical or electronics system, modified by guidelines 2003/108/EG. Please contact Siemens

Service in case of questions about returning and disposing the MR system and/or its components and accessories. **Safety certificate**

Revision label:

Revision:														Model - No.:																		
00	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98

Patent label:

This product is covered by one or more of the following US-patents:		
4587504	5541514	6073042
4714887	5545990	6084409
4743853	5546081	6084410
4769603	5546299	6100693
4777807	5550471	6101644
4825159	5557204	6107799
4879538	5574372	6111412
4922204	5581184	6114854
4945321	5581187	6118337
4972852	5592091	6118681
4991586	5614827	6128523
4992769	5623207	6133737
	5627470	6145026
5034692	5633586	6148225
5111378	5642264	6150880
5128615	5662112	6154068
5138259	5663646	6157191
5144241	5668474	6157193
5153517	5678549	6160400
5189370	5680045	6160445
5200701	5681327	6160453
5206591	5691678	6168403
5208537	5692508	6170019
5210512	5695872	6172558
5230090	5708361	6188923
5235283	5712567	6195031
5245282	5777475	6195578
5258718	5828215	6198287
5291741	5850143	6205349
5294886	5857970	6215911
5304929	5872500	6218839
5307014	5879298	6232548
5309107	5884489	6236204
5329266	5886524	6236209
5332990	5901036	6240310
5345178	5913863	6246239
5363078	5932936	6253101
5396174	5944663	6265872
5442292	5948521	6295465
5451877	5969568	6297633
5453866	5973527	6297637
5459401	5990625	6300761
5459540	5991179	6307374
5467017		6335620
5471142	6009341	6339332
5474067	6023799	6342785
5492124	6025716	6342786
5512828	6025720	6351123
5515002	6026315	6366090
5519321	6043851	6369569
5523689	6064204	6385480

For laser light localizer:

The EAC label:

The Eurasian Conformity mark is a product marking to indicate that the product conforms to all appropriate Technical Regulations of the Customs Union assessment procedures.

Artwork of the MD labeling in Russian*:

(*please note that Labeling mentioned below is applied only for the MAGNETOM C! systems which are delivered to Russia not earlier than May 2019)

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Depth

Height (t

8.2. Additional system characteristics

Magnet parameters	
C-shaped permanent magnet with vertical magnetic field	
Operating field strength	0.35 Tesla (-0.0094 T ~ -0.0047 T)

RF system	
RF transmit coils	Body coil
Peak power of transmitter amplifier	2.5 kW
Transmitter bandwidth	250 kHz
Max. applied RF field B ₁₊	35.25 μT
Max. specified B ₁₊ rms	12.1 μT
Distance to isocenter, where RF transmit field is reduced by	
3 dB	0.18 m (±0.02m)
10 dB	0.60 m (±0.08m)

Gradient system	
Actively shielded flat gradient coil system (Whole body gradient system)	
Max. gradient field strength	24 mT/m
Min. rise time	400 μs
Max. slew rate	55 T/m/s
Effective performance	
Max. eff. gradient field strength	39.4 mT/m
Max. eff. slew rate	92.5 T/m/s

Patient table	
Width	79 cm (31.1 in.) (±1mm)
Depth	210 cm (82.7 in.) (±1mm)
Height (fixed)	81.5 cm* (32.1 in.) (±5mm)

Max. patient weight	200 kg (440 lbs)
Longitudinal travel (manual)	± 83 cm (± 32.6 in.) (± 5 mm)
Lateral travel (manual)	$- 81.8/+ 21.2$ cm ($- 32.2/+ 8.3$ in.) (± 5 mm)
Positioning accuracy	± 1 mm

* above floor level

Patient comfort	
The open design provides free patient access from three sides	
Distance between table top and upper pole cover	38.5 cm (15.2 in.) (± 5 mm)
Intercom system	with loudspeaker and microphone

9. The IEC 60601-1:2007 Deviation

IEC 60601-1:2007
CORR.2(2007)

IEC 60601-1-1

IEC 60601-1-2

IEC 60601-2-1

IEC 62464-1

IEC 62304-1

IEC 62366-1

ISO 14971

ISO 13485

ISO 10993

ISO 15223

10.

The Me

11

On dis
observ

9. The list of International Standards and Directives applied to the Medical Device.

IEC 60601-1:2005+CORR.1(2006)+ CORR.2(2007)+AM1(2012)	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2:2007	Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
IEC 60601-1-6:2010	Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 60601-2-33: 2010 +A1: 2013	Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 62464-1:2007	Magnetic resonance equipment for medical imaging – Part 1: Determination of essential image quality parameters
IEC 62304:2006	Medical device software - Software life-cycle processes
IEC 62366:2007+A1:2014	Medical devices - Application of usability engineering to medical devices
ISO 14971:2007	Medical devices — Application of risk management to medical devices
ISO 13485:2016	Medical devices Quality management systems Requirements for regulatory purposes
ISO 10993-1:2009	Biological evaluation of medical devices: Part 1 Guidance on selection of tests
ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

10. Sterility.

The Medical device is non-sterile.

11. Disposal.

On disposing of the system or parts thereof, currently valid environmental legislation must be observed.

Examples of environmentally relevant components are:

- Accumulators and batteries
- Transformers
- Capacitors
- Monitor picture tubes
- Phantoms

For details contact your local customer service representative or your Siemens regional office.

NOTICE: System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

After usage the product can be classified as a medical waste of class "A".

All medical waste shall be collected, used, decontaminated, placed, stored, transported, listed and disposed of in accordance with the laws on public health and epidemiological welfare.

12. Transportation requirements and conditions.

12.1. Mechanical environmental conditions

Transport and handling loads according to DIN EN 60721-3-2,
Mechanical load class 2M2

Transport position any X restricted

Define transport position: standing vertically on the pallet.

12.2. Climatic environmental conditions

Temperature, Humidity – and pressure according to DIN EN 60721-3-2, Climate class 2K4

Temperature -20 °C to 60 °C

Humidity 10% to 90%

Air pressure 70 to 106kPa

To be fulfilled specifications/ regulations where checked in accordance with quality check list for designed packings.

The specifications in accordance with packed good passport are fulfilled.

13. Storage requirements and conditions.

Storage conditions based on system part localization:

Climatic conditions:

Room	Room temperature	Relative air humidity	Temperature gradient	Air pressure
Examination room	24 ± 3 °C	40–80%	</=3 ° C/h	70kPa-106kPa
Evaluation room/ Control room	15°C – 30 °C	40–80%	n.a.	

Equipment room	18°C – 24 °C	40–80%	</=1° C/5min	
----------------	--------------	--------	--------------	--

If electrical cabinet and cooling system are installed in the control room, equipment room's climate conditions are to be applied.

Heat dissipation:

Room	Component	Heat dissipation to air
Examination room	Magnet	max. 2kW
Evaluation room/ Control room	Host computer	max. 1.8kW
	Monitor	
Equipment room	ECA	max. 11.4kW
	ICA	

14. Service contacts

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract.

It is also possible to take out a service contract based on system usage.

Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubinskaja street 96,
115093 Moscow, Russia
Tel.: +7 495 737 12 52
www.siemens-healthineers.com/ru/

Customers support center (24 hours):

8800 550 5022

Service.ru@siemens-healthineers.com

15. Service life.

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of **10** (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and operating safety of these products.

Please contact Siemens about these measures early enough!

16. Warranty.

Siemens guarantees that the Equipment has no defects caused by installation and, as of the date of transfer, is compliant with the approved specifications and has no defects in materials and products within 14 (fourteen) months as of the Equipment delivery date. Warranty terms and conditions for components made by the third parties and supplied along with the Equipment shall be detailed in the respective Order.

Warranty period for operable Equipment components supplied to replace the defected ones shall be (six) months as of their delivery date. Total warranty period for defective components repaired and installed by Siemens shall not exceed 14 months as of the Equipment delivery date.

Defective Equipment (its part), Spare parts, high-vacuum components, other types of Products replaced within the warranty period shall be handed over to Siemens to be further returned to the manufacturer no later than 14 (fourteen) days as of the date of signing the documents (certificate/order form) provided by Siemens, confirming that the respective defect was eliminated. Title to such spare part shall be transferred to Siemens.

17. Contacts for reclamation.

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC
Dubinskaya street 96
115093 Moscow,
Russia
Tel.: +7 495 737 10 00
Fax: +7 495 737 13 20

Contacts for questions and reclamations in the country of Legal manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center,
Gaoxin C. Ave. 2nd, Hi-Tech Industrial Park,
518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA
Tel.: +86 0755 26525421

Glossary

Name of the component
RF system
syngo Acquisition
PC Keyboard
PMU Wireles
syngo in-Ro
Head Array
CP Head/Ne
Neck coil, L
Neck coil, M
Body/Spine
Body/Spine
Body/Spine
Body/Spine
Body/Spine
Breast Arr
Extremity
Extremity
Extremity
Shoulder
Wrist Arr
4-channe
MultiPurp
MultiPur
Cordless

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