ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

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(sign and seal)

27 February 2020

COMPRESSION DEVICE PreludeSYNC
INSTRUCTIONS FOR USE

PreludeSYNC

English

Radial Compression Device

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The Prelude Sync™ Radial Compression Device is a sterile, single use disposable device used to assist in gaining and maintaining hemostasis of the radial and ulnar artery following catheterization procedures. It has a soft wristband with a secure hook and loop fastener. The band delivers adjustable compression of the puncture site with an inflatable bulb, and a check valve for easy inflation and deflation with a syringe. A clear curved backer plate provides optimal visualization of the puncture site and ease of placement.

INDICATIONS FOR USE

The Prelude Sync is a compression device used to assist in gaining hemostasis of arterial percutaneous access sites.

CONTRAINDICATIONS

- Patients hypersensitive to the materials of the compression device.
- Patients with infection or other serious skin diseases at the site of puncture.
- Pertaining to radial artery catheterization only: patients with an abnormal Allens test, radial pulse, or insufficient dual artery supply.
- Not indicated for femoral artery compression.

WARNINGS

- Prior to inflation of bulb, confirm that air is being injected into the Prelude Sync and NOT the side port of the sheath or other device.
- Ensure the band is fastened securely around the wrist without slack.
- Over-inflation of bulb (above 20mL of air) may cause pain, numbness, artery occlusion, or damage to the device.
- Under-inflation of bulb, or failure to secure band without slack around the wrist, may compromise the ability of the device to assist hemostasis of the artery.
- Arterial pulse distal to the compression device should be monitored to ensure the artery is not completely occluded as arterial damage or thrombosis may occur.
- Patients should not be left unattended while the Prelude Sync is in use.
- Do not leave the Prelude Sync on for an inappropriately long period of time as tissue damage may occur.
- Do not expose the Prelude Sync to organic solvents, as they may cause damage to the device.

CAUTIONS

- Maintain sterile field during application.
- This device should be used by clinicians with adequate training in the use of the device.
- Sterile if package is unopened and undamaged.

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POTENTIAL COMPLICATIONS

Possible complications that may result from use of this device include, but are not limited to: hematoma, recurrent bleeding, local venous thrombosis, nerve damage, pain or numbness, complex regional pain syndrome, allergic reaction, and artery occlusion.

INSTRUCTIONS FOR USE

Device Placement

1. Ensure site is clean and dry.

NOTE: This device needs to be positioned differently for the following uses:

- Radial artery use: Ensure the CURVED section of the clear plate is on the THUMB side of the wrist.
- Ulnar artery use: Ensure the CURVED section of the clear plate is on the LITTLE FINGER (fifth digit) side of the wrist.

Figure 1.

2. Aspirate the sheath, then withdraw the sheath approximately one inch (2-3 cm).


Figure 2.

3. Place the center of the "crosshairs" over the arteriotomy (location where the sheath entered the artery, approximately 1-2mm proximal to the skin puncture site).


Figure 3.

4. Fasten the band securely around the wrist without any slack, but do not overtighten.

5. Fill the Prelude Sync syringe (included) with 20mL of air.
NOTE: Nominal air inflation: 15mL
 Maximum fill volume: 20mL
6. Attach and completely engage the Prelude Sync syringe to the tubing line labeled "AIR".
7. Slowly inflate the bulb with air while simultaneously removing the sheath. Once the sheath is completely removed, continue to inject air into the bulb until bleeding has stopped.

Figure 4.

- NOTE:** While inflating, maintain pressure on the syringe plunger to avoid any inadvertent release of air.
8. Slowly withdraw air from the bulb until there is oozing from the access site. Once oozing is observed, re-inject up to 4mLs of air into the bulb until hemostasis is achieved. If bleeding or hematoma are present, additional air may be injected as needed to achieve patent hemostasis.


Figure 5.

9. Remove syringe.
NOTE: If bleeding is observed at any time, inject additional air (not exceeding the max fill volume of 20mL) until bleeding stops.
10. Optional tubing clip: You may snap the loose tubing into the clip to secure it to the device.
11. Per hospital protocol, record the patient's vital signs and ensure adequate distal perfusion is maintained (patent hemostasis). If necessary, adjust air volume in bulb.
NOTE: Air volume and compression time may differ according to patient's condition, anticoagulant dosage, and size of puncture site.
12. Optional device/chart stickers are included to facilitate recording the time of inflation and mLs, as well as the time of deflation and mLs removed.


Figure 6.

Device Removal

1. At the recommended device removal time, withdraw 2mL of air from the Prelude Sync and observe the access site for bleeding.

Figure 7.

2. If bleeding is present, inject air until bleeding stops. Wait 30 min. and repeat step 1. of the device removal instructions.
3. If no bleeding is present, continue to remove 2mL of air every 15 minutes x3, or until pressure is fully released.
4. Once the air is removed and hemostasis is confirmed, carefully remove the Prelude Sync. Place a sterile dressing per hospital protocol.

Rx Only	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult accompanying documents. Read instructions prior to use.
	Single use.
	Do not use if package is damaged.
STERILE EO	Sterilized using ethylene oxide.
REF	Catalog Number
LOT	Lot Number
	Used By

© 2017 Merit Medical Systems, Inc. All rights reserved.

CE⁰⁰⁸⁶

www.merit.com

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
EC Customer Service +31 43 3588222

401766001_004 2017-06-13

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The PreludeSYNC DISTAL Radial Compression Device is a sterile, single use disposable device used to assist in gaining and maintaining hemostasis of the radial artery following catheterization procedures. It has a soft wristband and thumb saddle strap with secure hook and loop fasteners. The band delivers adjustable compression of the puncture site with an inflatable balloon, and a check valve for easy inflation and deflation with a syringe. A clear curved backer plate provides optimal visualization of the puncture site and ease of placement.

INDICATIONS FOR USE

The PreludeSYNC is a compression device used to assist in gaining hemostasis of arterial percutaneous access sites.

CONTRAINDICATIONS

- Patients hypersensitive to the materials of the compression device.
- Patients with infection or other serious skin diseases at the site of puncture.
- Not indicated for femoral artery compression.
- Patients with an abnormal Allens test, radial pulse, or insufficient dual artery supply.

WARNINGS

- Prior to inflation of balloon, confirm that air is being injected into the PreludeSYNC DISTAL and NOT the side port of the sheath or other device.
- Ensure the band is fastened tightly around the wrist and hand without slack.
- Prolonged over-inflation of balloon may cause pain, numbness, artery occlusion, nerve damage or damage to the device.
- Under-inflation of balloon, or failure to secure band without slack around the wrist and hand, may compromise the ability of the device to assist hemostasis of the artery, resulting in bleeding, and/or hematoma.
- Arterial pulse distal to the compression device should be monitored to ensure the artery is not completely occluded as arterial damage or thrombosis may occur.
- Patients should not be left unattended while the PreludeSYNC DISTAL is in use.
- Do not leave the PreludeSYNC DISTAL on for an inappropriately long period of time as tissue damage, or arterial occlusion may occur.
- Do not expose the PreludeSYNC DISTAL to organic solvents, as they may cause damage to the device.

CAUTIONS

- Maintain sterile field during application.
- This device should be used by clinicians with adequate training in the use of the device.
- Sterile if package is unopened and undamaged.

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
POTENTIAL COMPLICATIONS

Possible complications that may result from use of this device include, but are not limited to: hematoma, recurrent bleeding, local venous thrombosis, nerve damage, pain or numbness, complex regional pain syndrome, allergic reaction, and artery occlusion.


INSTRUCTIONS FOR USE

Device Placement

1. Ensure site is clean and dry.
2. For left hand placements use catalog code SDRB-REG-LT. For right hand placements use catalog code SDRB-REG-RT.
3. For a large patient wrist/hand, you may add the extension band (included) prior to placing the PreludeSYNC DISTAL (FIGURE 1).

4. Aspirate the sheath, then withdraw the sheath approximately one inch (2-3 cm).
5. Place the band on the hand/wrist with band #1 on the little finger side, band #2 on the thumb side, and band #3 facing up toward the thumb (FIGURE 2).

- Place the center of the balloon "target" over the arteriotomy (location where the sheath entered the artery, approximately 1-2mm proximal to the skin puncture site) (FIGURE 3).


FIGURE 3

- Fasten band #1 and #2 tightly around the wrist, then wrap band #3 around the saddle of the thumb without any slack (FIGURE 4).
NOTE: If the sheath is in the way of securing band section #3, then band section #3 may be secured after removal of sheath.

FIGURE 4

Fasten band securely around the wrist

- Fill the syringe (included) with 10mL of air.
NOTE: Maximum balloon fill volume: 10mL.
- Engage the syringe into the valve on the top of the device.
- Slowly inflate the balloon with air while simultaneously removing the sheath. Once the sheath is completely removed, continue to inject air into the balloon until bleeding has stopped.
NOTE: While inflating, maintain pressure on syringe plunger to avoid any inadvertent release of air.
- Adjust air volume in balloon to achieve hemostasis. If bleeding or hematoma are present, additional air may be injected until bleeding stops.
- Remove the syringe.
NOTE: The patient should be able to move their hand/wrist freely without any bleeding. If desired; ask the patient to move their hand and wrist and check the site for any bleeding. If bleeding is observed, inject additional air (not exceeding the max fill volume of 10mL) until bleeding stops.
NOTE: If bleeding is observed at any time, inject additional air (not exceeding the max fill volume of 10mL) until bleeding stops.
NOTE: Keep included syringe accessible for device deflation and removal.
- Per hospital protocol, record the patient's vital signs and ensure adequate distal perfusion is maintained. If necessary, adjust air volume in balloon.
NOTE: Air volume and compression time may differ according to patient's condition, anticoagulant dosage, and size of puncture site.
- Device/chart stickers are included to facilitate recording the time of inflation and mLs, as well as the time of deflations and mLs removed.

DEVICE DEFLATION AND REMOVAL

- If desired, all air may be removed slowly at the appropriate amount of time to achieve complete hemostasis; or, a small amount of air (such as 1-2 mL) may be removed from the balloon periodically based on the hospital's preferred clinical technique/hospital protocol to aid in maintaining hemostasis.
- Upon reaching complete hemostasis, slowly withdraw remaining air from the PreludeSYNC DISTAL and observe the access site for bleeding.
NOTE: Complete hemostasis time may differ according to patient's condition, anticoagulant dosage, and size of puncture site.
- If bleeding is present, inject air until bleeding stops. Wait sufficient time and repeat the device deflation and removal instructions.
- Once the air is removed and hemostasis is confirmed, carefully remove the PreludeSYNC DISTAL.
- Place a sterile dressing per hospital protocol.

Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.		Caution: Consult accompanying documents. Read instructions prior to use.
	Single use.		Do not re-sterilize
	Do not use if package is damaged.		Sterilized using ethylene oxide.
	Left hand device		Right hand device
	For electronic copy scan QR code or go to www.merit.com/ifu and enter IFU ID Number. For printed copy, call U.S.A or E.U Customer Service.		

CE 0086

www.merit.com

Manufacturer:
Merit Medical Systems, Inc.
1600 West Merit Parkway,
South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A Customer Service 1-800-356-3748

Authorized Representative:
Merit Medical Ireland Ltd,
Parkmore Business Park West,
Galway, Ireland
EU Customer Service +31 43 358 82 22

403907001EN_002 ID 2018-03-12

INDICATION FOR USE

Compression devices are designed to achieve and maintain hemostasis during vascular catheterization procedures. The PreludeSync and PreludeSYNC DISTAL are a compression devices used to assist in gaining hemostasis of arterial percutaneous access sites.

DEVICE SPECIFICATIONConfigurations:

1. Compression device PreludeSYNC, composed of:

- compression device, length: 24 cm, 29 cm.
- syringe - 20 ml.
- the informational sticker.
- instruction for use.

2. Compression device PreludeSYNC DISTAL, composed of:

- compression device configuration: left, right;
- syringe - 10 ml;
- the informational sticker;
- instruction for use.

Table 1: PreludeSYNC specification

Device	Dimensional specifications	Mass, g
Compression device PreludeSYNC, length: 24 cm	240mm x 44.2mm	≤ 12.2
Compression device PreludeSYNC, length: 29 cm	290mm x 44.2mm	≤ 12.5
Compression device PreludeSYNC DISTAL, configuration: left, right	Loop 178mm x 25.4mm Short Hook 39.4mm x 25.4mm Long Hook 147.5mm x 25.4mm	≤ 11.5

Table 2: PreludeSYNC Syringe Inflator specification

Specification	10ml Syringe	20ml Syringe
Fill volume	10ml	20ml
Total graduated capacity	10ml	20ml
Maximum usable capacity	11ml	20ml
Mass (weight)	≤16g	≤21g
Graduation intervals	at least 1.0 ml	at least 5.0 ml




The main materials used in the manufacture of the PreludeSYNC and PreludeSYNC Distal are identified below:

Table 3: PreludeSYNC material specification

PRODUCTS COMPONENT	MATERIAL
Fastener	Nylon/PVC
Bladder	Polyurethane
Backer Plate	Thermoplastic Polyurethane
Fill tube with check valve	Connector tube: PVC Check valve: Polycarbonate Stem: Silicone
Adapter Cap	Polycarbonate
Adhesive	UV Adhesive
"AIR" Label	Polypropylene film, blue

PACKAGING

One (1) PreludeSYNC DISTAL/ PreludeSYNC Radial and one (1) syringe inflator are placed in a thermoformed PETG blister tray. The formed blister tray is sealed with a Tyvek lid using validated equipment in a controlled environment. Five (5) blister trays are placed inside a shelf carton. Four (4) shelf cartons are placed into a shipper carton.

Applicable Symbols	
	This symbol means that the device fully complies with European Directive 93/42/EEC.
	Manufacturer
	Authorized Representative
	Single use
	Do not resterilize
	Reference Number
	Lot Number
	Sterilization: Ethylene-Oxide Gas
	Caution: U. S. federal law restricts this device to sale by or on the order of a physician.
	Use by
	Caution: Consult Accompanying Documents
	For electronic copy scan QR code, or go to www.merit.com.lfu
	Do not use if package is damaged.
	Non-pyrogenic
	Manufacturing date

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Ethylene Oxide

MAINTENANCE:

NOT APPLICABLE. The product is intended for single use only.

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WARRANTY**Limited Warranty and Disclaimer**

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

Instructions for Use
Compression Device, SafeGuard
Dated February 27, 2020
Consisting of 11 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 27th day of February 2020, before me, Jennifer Fordham, a notary public, personally appeared Robert Wolfarth, Senior Manager, Regulatory Affairs proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public
My Commission Expires: June 27, 2021

Robert Wolfarth
Senior Manager, Regulatory Affairs
Merit Medical Systems, Inc. USA

(sign and seal)

27 February 2020

COMPRESSION DEVICE SafeGuard
INSTRUCTIONS FOR USE

SAFEGUARD

24cm PRESSURE ASSISTED DEVICE

INSTRUCTIONS FOR USE

Read instructions prior to use.
Product not made with natural rubber latex.

DEVICE DESCRIPTION

The SAFEGUARD 24 cm is a single use disposable device. SAFEGUARD has a clear medical grade polyurethane window and bladder, a clear medical grade PVC flexible fill tube, and a pressure sensitive, self-adhesive peel backing. A valve on the end of the fill tube enables a syringe to be connected to inflate the central bladder with air to provide pressure to the puncture site. The SAFEGUARD pressure assisted device has a sterile dressing with a clear window that facilitates visibility of the access site without removal or manipulation of the device.

INDICATIONS

The indications for use for the SAFEGUARD 24cm pressure assisted device are to assist in obtaining and maintaining hemostasis. The device is also indicated in the reduction of active compression time in femoral artery cannulation following diagnostic and interventional procedures with an ACT of 140 seconds or less, using a 6 Fr. and smaller sheath size.

CONTRAINDICATIONS

The adhesive portion of the SAFEGUARD device should not be used over excoriated skin.

CAUTIONS

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: With over-inflation i.e., above 40 mL of air, the bulb may begin to expand radially and could compromise the adhesive properties of the device.

CAUTION: Under-inflation of device could compromise the ability of the device to assist in obtaining and maintaining hemostasis.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

- Use proper aseptic techniques while handling product.
- Do not use if package is damaged
- Inspect device prior to use to verify that no damage has occurred during shipping.

POTENTIAL ADVERSE EFFECTS

Possible adverse effects that may result from the use of this device:

- Hematoma
- Local bleeding
- Arterio-venous fistula or pseudoaneurysm

	SAFEGUARD MAT N=100	Manual Compression Historical Control N=85
Avg. AC/TTH Manual Compression (minutes)	7.7 ± 3.3	28

	SAFEGUARD MAT N=101		Manual Compression Historical Control N=85	
	Diagnostic	Interventional	Diagnostic	Interventional
Major Complications				

Major Complications

Access site bleeding requiring transfusion	0	0	0	0
Vascular repair or the need for vascular repair (via surgery, ultrasound-guided intervention, catheter embolization, or stent-graft)	0	1	0	2
Any new ipsilateral lower extremity ischemia	0	0	0	0
Surgery for access site related nerve injury	0	0	0	0
Access site-related infection requiring IV antibiotics or extending hospitalization	0	0	0	0

Minor Complications

Access site hematoma >6 cm (after sheath pull)	1	1	-	-
Bleeding requiring > 30 minutes to re-establish hemostasis	0	0	-	-
Non-treated pseudoaneurysm	0	0	-	-
Non-treated arteriovenous (AV) fistula documented by ultrasound	0	0	-	-
Skin excoriation at site of Safeguard after removal of dressing	0	0	-	-
Skin erythema	0	0	-	-
Allergic reaction to adhesive	0	0	-	-
Ipsilateral lower extremity arterial emboli, transient loss of lower extremity pulse, or deep vein thrombosis	0	0	-	-
Access site-related vessel laceration	0	0	-	-
Transient access-site related nerve injury	0	0	-	-
Access site wound dehiscence	0	0	-	-
Localized access site infection treated with intramuscular or oral antibiotics	0	0	-	-

*None Reported

There is no statistical difference in the rate of major complications between the control group and the treated group.

SAFETY AND EFFECTIVENESS RESULTS

A clinical trial was conducted to evaluate the safety and effectiveness of the SAFEGUARD Manual Assist Technique (MAT) in reducing active compression time compared to historical manual compression data. (Reduced Vascular Complications After Percutaneous Coronary Interventions with a Non-mechanical Suture Device: Results from the Randomized RACE Study, Sanborn, TA)

PRE-HEMOSTASIS, or MANUAL ASSIST TECHNIQUE (MAT)

PLACEMENT OF SAFEGUARD

1. Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port and to allow for easy sheath removal.

Note: Placement may require adjustment based on the patient's anatomy, angle of the puncture site, and the presence or absence of a procedural sheath.

2. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.
3. Pull the procedural sheath back approximately 1" (2.5 cm) so that when SAFEGUARD is adhered to the skin the sheath hub is outside the area of the SAFEGUARD adhesive.





Note: It is recommended that you aspirate the sheath prior to removal to prevent distal embolization from residual clot in sheath.

4. Remove the adhesive backing and place the bulb where you would position your fingers to hold manual compression (for example, in femoral artery procedures, typically the point of maximum femoral pulse). Make sure SAFEGUARD is completely adhered to the skin.

5. Attach and completely engage the appropriately configured syringe to inflate the desired volume a maximum volume of 40 mL's of air into the bulb to apply pressure on the arteriotomy site. Syringe must be completely engaged in the valve to inflate/deflate the bulb. Remove syringe.

Note: Maintain pressure on the plunger while detaching syringe from the SAFEGUARD valve. Observe that the desired pressure is achieved and maintained.

6. Remove sheath, then immediately apply manual compression directly over inflated bulb.
7. Hold manual compression until hemostasis has been achieved.*
- Slowly release manual compression.
- Check distal/proximal pulses to assure flow is maintained.
- Confirm hemostasis by viewing the site through the inflated bulb window.

* Recommendations (MAT only):

Diagnostic patients - minimum 5 minutes

Interventional patients - minimum 10 minutes

8. Per hospital protocol, periodically check the site through the bulb window to confirm hemostasis and to manage the bulb volume and resultant pressure as needed. Continue to check distal/proximal blood flow to assure patency.
9. Deflate bulb every two hours to allow for capillary refill and to assess the site. Re-inflate the bulb if necessary.
10. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

11. Prior to discharge of the patient, remove SAFEGUARD and apply sterile dressing per hospital protocol.

POST-HEMOSTASIS TECHNIQUE

1. When hemostasis at the access site has been achieved, apply the SAFEGUARD device with the access site visible under the bulb window of the SAFEGUARD device. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.

Note: Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to valve inflate/deflate port.

2. Attach an appropriately sized and configured syringe to the valve of the SAFEGUARD device.

Note: Syringe must be completely engaged valve to inflate/deflate the bulb.



3. Inflate the bulb of the SAFEGUARD device with air to the desired volume of air (24cm maximum of 40 mL's) to apply pressure on the arteriotomy site and remove the syringe. Check distal/proximal pulses to assure flow is maintained.

Note: Maintain pressure on the plunger while detaching syringe from the valve of the SAFEGUARD device. Observe that the desired volume is achieved and maintained.

4. Per hospital protocol, periodically check the site through the bulb window to assure hemostasis is maintained and the bulb maintains pressure.
5. Deflate bulb every two hours and assess the site. Re-inflate the bulb if necessary.
6. Deflate the bulb by attaching an appropriately sized and configured syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove the plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

7. Prior to discharge of the patient, remove the SAFEGUARD device and apply sterile dressing per hospital protocol.

-  Maintain sterile field during application.
-  Keep away from sunlight
-  Do not Resterilize
-  Do not Reuse
-  Do not use if package is damaged
-  Catalog Number
-  Lot Number
-  Use By
-  Sterilized Using Gamma
-  Caution: Consult accompanying document
-  Manufacturer
-  Authorized Representative

MERITMEDICAL®

CE 0086

Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

EC REP

Authorized Representative: Merit Medical Ireland Ltd.
Parkmore Business Park West, Galway, Ireland

www.merit.com

ID 050115 400373004/A

-  Maintain sterile field during application.
-  Keep away from sunlight
-  Do not Resterilize
-  Do not Reuse
-  Do not use if package is damaged
-  Catalog Number
-  Lot Number
-  Use By
-  Sterilized Using Gamma
-  Caution: Consult accompanying document
-  Manufacturer
-  Authorized Representative

MERITMEDICAL®

CE 0086

Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A.
1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

EC REP

Authorized Representative: Merit Medical Ireland Ltd.
Parkmore Business Park West, Galway, Ireland

www.merit.com

ID 050115 400373004/A

INDICATION FOR USE

Compression devices are designed to achieve and maintain hemostasis during vascular catheterization procedures.

SAFEGUARD 12 cm: The indications for use for the SAFEGUARD 12 cm pressure assisted device are to assist in obtaining and maintaining hemostasis.

SAFEGUARD 24cm: The indications for use for the SAFEGUARD 24cm pressure assisted device are to assist in obtaining and maintaining hemostasis. The device is also indicated in the reduction of active compression time in femoral artery cannulation following diagnostic and interventional procedures with an ACT of 140 seconds or less, using a 6 Fr. and smaller sheath size.

DEVICE SPECIFICATION

Configurations:

Compression device SafeGuard, versions:

1. Compression SafeGuard device, length 12 cm composed of:

- compression device;
- syringe - 7 ml (if necessary);
- instruction for use.

2. Compression SafeGuard device, length 24 cm composed of:

- compression device;
- syringe - 40 ml (if necessary);
- instruction for use.

Table 1: Safeguard specification

Device	Dimensional specifications	Mass, g
Compression device SafeGuard device, length 12 cm	120mm x 71.1mm	<4.5
Compression SafeGuard device, length 24 cm	240mm x 142.2mm	<18

Table 2: Safeguard Syringe Inflator specification

Specification	7ml Syringe	40ml Syringe
Fill volume	7ml	40ml
Total graduated capacity	7ml	40ml
Maximum usable capacity	11ml	60ml
Mass (weight)	≤16g	≤36
Graduation intervals	at least 1.0 ml	at least 5.0 ml

The main materials used in the manufacture of the Safeguard 12&24cm are identified below:

Table 3: Safeguard material specification

PRODUCTS COMPONENT	MATERIAL
Non-woven 12cm Sub-assembly	Polyester non-woven medical tape; blue ink
Bladder	Polyurethane
Peel backing	HDPE, white
Fill tube with Luer check valve	Connector Tube: PVC Check valve: PVC
Adapter Cap	Polycarbonate Valve Adapter Cap
O-ring	Silicone
Molded Backer w/ Doghouse	Polyurethane resin
Dymax Adhesive	Medical Device UV Adhesive

PACKAGING

The single device is loaded into a pouch and so that the device graphics are visible when viewed through the clear side of the pouch. A label is placed on the paper side of the pouch manufactured from uncoated Tyvek sealed to polyethylene film. Each shelf carton contains ten (10) labeled pouches. The carton shall also contain an appropriate IFU. Each corrugated shipping case shall contain ten (10) labeled shelf cartons.

The Safeguard device composed with the Safeguard Syringe shall be placed on a die-cut card and loaded into a pouch and so that the device graphics are visible when viewed through the clear side of the pouch.

SAFEGUARD 12CM: Each shelf carton for shall contain ten (10) labeled pouches oriented in alternating directions with the pouch labels facing up, toward the tear-away window.

SAFEGUARD 24CM: Each shelf carton for shall contain five (5) labeled pouches oriented in alternating directions with the pouch labels facing up, toward the tear-away window.

The carton shall also contain an appropriate IFU. Each corrugated shipping case shall contain four (4) labeled shelf cartons.

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS

Transportation temperature	-35 °C to +55 °C
Storage temperature	+5 °C to +27 °C
Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

Use Conditions - in medical and preventive treatment facility.

SHELF LIFE: Shelf life – 3 year.

STERILIZATION: Gamma sterilization

MAINTENANCE: NOT APPLICABLE. The product is intended for single use only.

DISPOSAL

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability


UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

Customer Service

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
+ 1 801 253 1600

Merit Technologies LLC
119017 Moscow,
1st Kazachiy pereulok,
bld 7, 1 floor, room 6, Russian Federation
+ 7 495 221 89 02

INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION
Merit Medical's SAFEGUARD Radial compression device is a sterile, single use disposable device. It has a clear medical grade polyurethane window and bulb (balloon) that facilitates visualization of the puncture site, a clear medical grade PVC flexible tube, and a pressure sensitive, self-adhesive peel backing. A valve on the end of the tube enables a standard Luer or proprietary syringe to be connected to inflate and deflate the bulb with air to provide compression of the transradial puncture site. Product not made with natural rubber latex.

SAFEGUARD Radial is a compression device to assist hemostasis of the radial artery after a transradial procedure.

CONTRAINDICATIONS
The adhesive portion of SAFEGUARD Radial should not be used over excoriated skin.

- Prior to inflation of bulb, confirm that you are **ONLY** injecting air into **SAFEGUARD** radial compression device and **NOT** to side port of sheath or other device.
- **Over-inflation** of bulb (above 7 mL of air) could cause bulb to excessively expand and cause pain, numbness, or radial artery occlusion.
- **Under-inflation** of bulb could compromise ability of device to assist hemostasis of radial artery.
- **Maintain sterile field** during application.

- Read instructions prior to use.
- **RE Only:** Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
- This device is intended for single use only. Do not reuse or resterilize.
- Do not use if package is damaged.
- This device should be used by clinicians with adequate training in the use of the device.

REUSE PRECAUTION STATEMENT

For single-patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

POTENTIAL COMPLICATIONS
Possible complications that may result from use of this device include, but are not limited to: hematoma, bleeding, pain or numbness, radial artery occlusion.

INSTRUCTIONS FOR USE

1. Ensure site is clean and dry.
2. After procedure, aspirate sheath, then withdraw sheath one inch (2.5 cm). [Figure 1]
3. Prior to use, the bulb of SAFEGUARD Radial compression device approximately halfway (from distal), then apply center of bulb over antecubital site. [Figure 2]

Note: The arrow on SAFEGUARD Radial should be pointing up the arm of the patient.

4. Remove remainder of adhesive backing (on side of a time) and completely secure around wrist. [Figure 3]
5. Attach and completely engage a standard syringe or proprietary syringe to the tubing line labeled "AIR"

WARNING: Prior to inflation of bulb, confirm that you are **ONLY** injecting air into SAFEGUARD Radial compression device and **NOT** to side port of sheath or other device.

6. Inflate the bulb to (difficult) with a maximum volume of 7 ml.

WARNING: Over-inflation of bulb above 7 mL of air could cause bulb to excessively expand and cause pain, numbness, or radial artery occlusion.

WARNING: Under-inflation of bulb could compromise ability of device to assist hemostasis of radial artery.

puncture site by viewing site through bulb window. If bleeding is observed at any time, inject more air (not exceeding the

8. Per hospital protocol, ensure adequate distal perfusion is

Note: Air volume and compression time may differ according to patient's condition, anthropometric data, and size of puncture.

9. Before removing device, confirm that bleeding has stopped.

Note. Be sure not to create a vacuum by pulling back on the

10. Carefully remove SAFEGUARD Radial compression device and

apply sterile dressing per hospital protocol

INDICATION FOR USE

Compression devices are designed to achieve and maintain hemostasis during vascular catheterization procedures. **SAFEGUARD Radial** is a compression device to assist hemostasis of the radial artery after a transradial procedure.

DEVICE SPECIFICATION

Configurations:

Compression SafeGuard Radial device, length 26 cm composed of:

- compression device;
- syringe - 7 ml (if necessary);
- instruction for use.

Table 1: Safeguard Radial specification

Device	Dimensional specifications	Mass, g
Compression device SafeGuard Radial	260mm x 41.3mm	<10

Table 2: Safeguard Radial Syringe Inflator specification

Specification	7ml Syringe
Fill volume	7ml
Total graduated capacity	7ml
Maximum usable capacity	11ml
Mass (weight)	≤16g
Graduation intervals	at least 1.0 ml

The main materials used in the manufacture of the Safeguard Radial are identified below:

Table 3: Safeguard Radial material specification

PRODUCTS COMPONENT	MATERIAL
Non-woven Sub-assembly	Polyester non-woven medical tape; blue ink
Bladder	Polyurethane
Peel Backing	HDPE, white
Fill tube with Luer check valve	Connector tube: Tygon PVC Check valve: PVC
Molded Backer w/ Doghouse	Polyurethane resin
Dymax Adhesive	Medical Device UV Adhesive

PACKAGING

The device is placed on a die-cut card and loaded into a pouch and so that the device graphics are visible when viewed through the clear side of the pouch. A label is placed on the paper side of the pouch manufactured from uncoated Tyvek sealed to polyethylene film.



Each shelf carton contains ten (10) labeled pouches. The carton shall also contain an appropriate IFU.

Each corrugated shipping case shall contain ten (10) labeled shelf cartons for single device.

Each corrugated shipping case shall contain four (4) labeled shelf cartons for Safeguard with Syringe.

Applicable Symbols

	This symbol means that the device fully complies with European Directive 93/42/EEC.
	Manufacturer
	Authorized Representative
	Single use

	Do not re-sterilize
	Reference Number
	Lot Number
	Sterilization: Gamma Irradiation
	Caution: U. S. federal law restricts this device to sale by or on the order of a physician.
	Use by
	Caution: Consult Accompanying Documents
	For electronic copy scan QR code, or go to www.merit.com.lfu
	Do not use if package is damaged.
	Non-pyrogenic
	Manufacturing date

TRANSPORTATION, STORAGE AND OPERATION CONDITIONS

Transportation temperature	-35 °C to +55 °C
Storage temperature	+5 °C to +27 °C
Usage	15 °C to 35 °C, relative humidity ≤80%, atm. Pressure 84.0 – 106.7 kPa (630-800 mm Hg).

Use Conditions - in medical and preventive treatment facility.

SHELF LIFE: Shelf life – 3 year.

STERILIZATION: Gamma sterilization

MAINTENANCE: NOT APPLICABLE. The product is intended for single use only.

DISPOSAL

Dispose the devices according to local environmental requirements.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

Customer Service

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A.
+ 1 801 253 1600

Merit Technologies LLC
119017 Moscow,
1st Kazachiy pereulok,
bld 7, 1 floor, room 6, Russian Federation
+ 7 495 221 89 02

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