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19-729017.09.15

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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of CHENGDU PUTH MEDICAL PLASTICS PACKAGING CO., LTD on the annexed INSTRUCTIONS FOR USE is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Sun Wenhui

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(Date: Sep. 06, 2019)

APPROVED
General Manager at
Chengdu Puth Medical Plastics Packaging Co., Ltd.
Hua Liao
2018

INSTRUCTIONS FOR USE

Medical device

«PUTH Holder for Vacuum Blood Collection Systems»

Produced by

«Chengdu Puth Medical Plastics Packaging Co., Ltd.», P.R.China

This Instructions for use apply to PUTH Holder for Vacuum Blood Collection Systems manufactured by Chengdu Puth Medical Plastics Packaging Co., Ltd., P.R. China, (hereinafter - holder).

The holder is available in the following versions:

1. PUTH Holder for Vacuum Blood Collection Systems, colorless, 50 pcs.;
2. PUTH Holder for Vacuum Blood Collection Systems, colorless, 100 pcs.;
3. PUTH Holder for Vacuum Blood Collection Systems, colorless, 250 pcs.;
4. PUTH Holder for Vacuum Blood Collection Systems, green, 50 pcs.
5. PUTH Holder for Vacuum Blood Collection Systems, green, 100 pcs.;
6. PUTH Holder for Vacuum Blood Collection Systems, green, 250 pcs.;
7. PUTH Holder for Vacuum Blood Collection Systems, pink, 50 pcs.;
8. PUTH Holder for Vacuum Blood Collection Systems, pink, 100 pcs.;
9. PUTH Holder for Vacuum Blood Collection Systems, pink, 250 pcs.;
10. PUTH Holder for Vacuum Blood Collection Systems, yellow, 50 pcs.;
11. PUTH Holder for Vacuum Blood Collection Systems, yellow, 100 pcs.;
12. PUTH Holder for Vacuum Blood Collection Systems, yellow, 250 pcs.

MANUFACTURER OF THE MEDICAL DEVICE

Chengdu Puth Medical Plastics Packaging Co., Ltd.

Manufacturer address:

No.6, Kexin Road, Western Hi-Tech Zone, Chengdu, P.R. China, 611731

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER

«REAMED» Co.Ltd, Russia

115419, Russia, Moscow, Shabolovka 34, I, room 6.

Tel.: (495) 933-34-36, 935-72-89;

E-mail: office@reamed.ru

INTENDED USE AND OPERATING PRINCIPLE

Intended Use

PUTH Holder for Vacuum Blood Collection Systems is part of the system for vacuum collection of venous blood, designed to safely position the needle for blood collection and its connection to a vacuum tube.

The medical device is intended for professional use by medical professionals engaged in the procedure for taking venous blood.

Indications: use for the venous blood collection procedure.

Contradictions: no.

Side effect: no (no contact with a patient).

This is a single use medical device.

Operating Principle

The principle of action of the medical device is based on the ability to hold the needle at the time of drawing blood into a vacuum tube.

When puncturing the vein with a special double-sided needle (special needle or butterfly needle) on the one side of the needle and the rubber tube stopper on the other side of the needle, the tube is filled with the required blood volume under vacuum, the holder is a holding lock, allowing manipulation of the held needle and as a guide when the vacuum tube is connected to the needle at the time of blood collection.

DESCRIPTION AND SPECIFICATIONS

The holder is made of plastic (propathene) in the form of a truncated cone with transparent walls finger flange on one side and a central branch with internal thread for screwing in the needle for blood collection on the other side.

The holder can have the following color options: colorless, green, pink, yellow.

Transparent walls of the holder allow to control the procedure for blood collection.

Nominal sizes and related tolerances are in the table:

Hight	46±3 mm
Diameter of the end to the needle	20±2 mm
Diameter of the part with finger flange	30±2 mm/26±2 mm (See Fig.. 1)
Diameter of the outlet with internal thread	8±2 mm
Mass	Not more than 3 g

The thread of the holder is compatible with the needles for blood collection and should provide a secure attachment of the needle at the time of blood collection. The internal diameter of the holder allows to use it with vacuum tubes of nominal diameter of up to 16 mm.

Packaging

Holders are packaged in plastic bags of 50, 100 or 250 pcs (middle packaging).

Holders in middle packaging are placed in shipping cardboard boxes.

RISKS of USE and PRECAUTIONS

This medical device is of risk class 1 (low individual risk and low public health risk).

The type of contact with patient - no.

When use as intended for the medical device, the risks are associated with the procedure for blood collection and possible contamination of the product with blood.

For safety it is necessary:

1. Use the product strictly in accordance with the intended use and according to the instructions for use.
2. Use disposable rubber or plastic gloves during the blood collection procedure. Human blood samples are potentially infectious material capable of long-term preservation and transmission of the viruses of HIV, hepatitis or any other infectious agents.
3. Dispose of the medical device in accordance with the instructions for use in accordance with the requirements adopted by health care facilities.
4. Do not use the product after the expiration date.
5. Do not reuse the product.

DIRECTIONS for the USE of MEDICAL DEVICE

Before the blood collection procedure prepare all the necessary items, check the expiration date, the absence of mechanical damage, a sufficient amount. Arrange the necessary items on the desktop so that they are convenient to use.

Perform blood collection procedure only with protective gloves.

1. The needle for venous blood collection is screwed into the needle holder until it stops, the needle cap is previously removed to pierce the tube cover.

2. Carry out the procedure of blood collection with required number of tubes.
3. After removing the needle from the vein, the needle can only be detached if a special container is used. In the absence of such a container, the disposal of the holder is carried out with the needle; it is prohibited to put protective cap on the used needle.
4. Used tools are placed in a special container for subsequent safe disposal. The holder is a single-use medical device.

MAINTENANCE AND REPAIR

Disposable medical device. The concepts of maintenance and repair for a medical device are not applicable.

DISPOSAL AND DESTRUCTION

Used holders must be disposed of in accordance with the rules and regulations adopted by hospitals. Disposal must take into account the risk of possible blood contamination of the product.

Unused holders by their characteristics correspond to solid domestic waste. Special measures for their disposal are not required.

STORAGE CONDITIONS

Shelf life

The shelf life when storage and transportation conditions are observed is 5 years.

Storage conditions

In the manufacturer's packaging at T from -50 to + 45°C.

Without direct sunlight.

Avoid deformation of the product and / or damage to the packaging.

Transportation conditions

Transportation by any kind of covered transport at T from -50 to + 50 °C.

LIST of STANDARDS for the MEDICAL DEVICE

EN ISO 13845:2016 «Medical devices. Quality management systems. Requirements for regulatory purposes».

ISO 15223-1:2016 (GOST R ISO 15223-1-2014) «Medical devices. Symbols to be used with medical device labels, labeling, and information to be supplied. Part 1. General requirements».

MANUFACTURER WARRANTIES

Chengdu Puth Medical Plastics Packaging Co.,Ltd., P.R. China, guarantees the compliance of the medical device «PUTH Holder for Vacuum Blood Collection Systems» with the technical documentation requirements provided that the transportation and storage conditions established by this Instructions for use.

Quality claims are sent to the authorized representative of the manufacturer.

The EXPLANATION of SYMBOLS WHICH CAN BE USED for LABELLING

Symbol	Symbol name	Symbol	Symbol name
	Do not reuse		Limit for storage horizons
	LOT number		Manufacturer
	catalogue number		fragile
	manufacturing date3a		Avoid exposure to sunlight
	Upper side		Use until
	Keep dry		temperature range

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