IFU Butterfly

DIC1

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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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Signature:

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(Date: Oct. 21, 2020)

证明书查询 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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DESCRIPTION

Butterfly needle:

1 vein needle tube cover; 2 vein needle tube; 3 wing; 4 plastic tube; 5 connector (Luer Adapter); 6 puncture needle tube; 7 rubber shield; 8 hub.

Needle specification is marked or recognized by according to wing and hub color.

Butterfly needles with a flexible catheter and a luer adapter for drawing blood into vacuum tubes are used as a system for taking venous blood mainly in children and those with "difficult" veins.

Devices for blood sampling with a butterfly needle are particularly relevant in pediatrics, gerontology, veterinary medicine (especially in small animals) and in the intensive care unit.

Each butterfly needle is sterile and individually wrapped.

The dimensions of the butterfly needle with a catheter correspond to the size of the injection needle. The color coding of the wing is made in accordance with the ISO standard 6009 (depending on the diameter).

Page 1 of 9

Butterfly needles are disposable. All the needles are siliconized (by silicon oil) and have a double oblique cut with a V-shaped laser sharpening. Butterfly needles have projections in the form of wings for puncture of hard-to-reach and thin veins. Butterfly needles for drawing blood also differ in the length of the connecting tube-catheter.

Butterfly needles are equipped with a flexible catheter made of pyrogen-free material that does not contain latex; at the end of the catheter there is a luer adapter for connection with the holder. Available versions of devices for blood collection with holder, in the latter, you can immediately insert the tube without attaching an additionally purchased holder. After removing the luer adapter from the catheter, you can connect the butterfly needle to the syringe or to the solution infusion system and perform both blood collection and intravenous injection or infusion "for one injection".

A transparent catheter provides visual control of the needle hitting the vein, since with the success of venipuncture, venous blood is immediately visualized in the catheter tube.

At the end of the catheter there is Luer adapter for screwing into a standard holder. The needle of the luer adapter is facing the test tube and has a protective rubber membrane (valve) with nipple properties.

The membrane prevents blood from leaking out, which allows blood to be taken in several tubes.

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After package is opened, please check the expiry date and integrity of products. Do not use the product if there is any damage, deformation or caps falling.

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INSTRUCTION FOR USE

1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.

2) Offer the patient or help him to take a comfortable position: sitting or lying down.

3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample).

4) Treat your hands hygienically, dry them.

5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.

6) Wear non-sterile gloves.

7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.

8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.

9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.

2) Insert the needle into the needle holder and screw in until it stops.

3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.

4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.

5) Remove the cap on the other side of the needle.

6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site, to tighten the skin.

7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.

8) Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.

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SYMBOLS ON PACKAGE

Table 1. Symbols description

	Batch code
	Date of manufacturer
	Used-by date
	Manufacturer
	Do not reuse
	Sterilized by radiation
	Do not use if pack is damaged
	Keep away from sunlight
	Keep dry
	Temperature range
	Fragile, handle with care
	Up
	Stacking layer limit

	Rostest symbol
	EC conformity

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Table 2. Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device. Values of the parameters, characteristics (features) of the medical device are presented in the tables below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	size	Needle size* (mm)	Catheter tube size	Color
Butterfly needle	20G	0,9 mm x 19 mm**	24 cm, 30 cm	Yellow
Butterfly needle	21G	0,8 mm x 19 mm	24 cm, 30 cm	Dark-green
Butterfly needle	22G	0,7 mm x 19 mm	24 cm, 30 cm	black
Butterfly needle	23G	0,6 mm x 19 mm	24 cm, 30 cm	Dark-blue
Butterfly needle	25G	0,5 mm x 19 mm	24 cm, 30 cm	Orange
Butterfly needle	27G	0,4 mm x 19 mm	30 cm	Light-Grey

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Needle cathogory – thin walled

The needle sharpening angle has cut at an angle of $(11 \pm 2)^\circ$ or a "short" cut at an angle of $(17 \pm 2)^\circ$.

Connection dimensions

Connection length of Luer and catheter tube (for butterfly needle) -8 mm $\pm 10\%$

Connection diameter of Luer and catheter tube (for butterfly needle) -5,5 mm $\pm 10\%$

Connection length of screw connection (for all needles) -5 mm $\pm 10\%$

Connection diameter of screw connection (for all needles) -4,5 mm $\pm 10\%$

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Connection strength for butterfly needle details

Connection strength of the butterfly needle tube and Luer adapter is not less than 15N

Connection strength of the butterfly needle tube and hub is not less than 15N

Connection strength of the butterfly needle hub and cap is 01-05 N

Diameter inner and outer of butterfly needle catheter	1.5 $\pm$ 0.1mm 2.2 $\pm$ 0.1mm
--	------------------------------------

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Plastic Tube (catheter) Polypropylene (PP)

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Transportation of products should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the temperature from -25 to +25 ° C. During transportation to provide protection from the sun, rain and snow, do not mix with harmful or toxic substances.

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PACKAGE. SCOPE OF DELIVERY

Butterfly needle and individually packed by PE pouch, packed in bag (50 or 100 pouches per bag) with IFU and transport carton (20 or 50 bags per carton).

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THIS IS TO CERTIFY THAT: the seal of HEBEI XINLE SCI & TECH CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

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(Date: Oct. 21, 2020)

书查询网 site for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Multi-sample needle

Flashback needle

Registration Certificate #XXX from XX.XX.XXXX

1. Multi-sample needle

1,2 mm x 25 mm (18G*1");
1,2 mm x 32 mm (18G*1¼");
1,2 mm x 38 mm (18G*1½");

1,1 mm x 25 mm (19G*1");
1,1 mm x 32 mm (19G*1¼");
1,1 mm x 38 mm (19G*1½");

0,9 mm x 25 mm (20G*1");
0,9 mm x 32 mm (20G*1¼");
0,9 mm x 38 mm (20G*1½");

0,8 mm x 25 mm (21G*1");
0,8 mm x 32 mm (21G*1¼");
0,8 mm x 38 mm (21G*1½");

0,7 mm x 25 mm (22G*1");
0,7 mm x 32 mm (22G*1¼");
0,7 mm x 38 mm (22G*1½");

0,6 mm x 25 mm (23G*1");
0,6 mm x 32 mm (23G*1¼");
0,6 mm x 38 mm (23G*1½");

2. Flashback needle

0,8 mm x 25 mm (21G*1");
0,8 mm x 38 mm (21G*1½");


0,7 mm x 25 mm (22G*1");
0,7 mm x 38 mm (22G*1½");

0,6 mm x 25 mm (23G*1");
0,6 mm x 38 mm (23G*1½");

DESCRIPTION

Multi-sample needle:


1 Back cap; 2 Puncture needle tube; 3 Hub; 4 Rubber Shield; 5 Front cap 6. Vein needle tube

Needle specification is marked or recognized by according to hub and (or) front cap color.

Flashback needle:

1 Back cap; 2 Vein needle tube; 3 Hub with transparent camera; 4 Puncture needle tube; 5 Rubber shield; 6 Front cap

Needle specification is marked or recognized by according to back cap color.

Multi-sample needle

It is a sterile needle with sharpening at both ends. One side of the needle is designed to be inserted into the patient's vein, the other side, closed with a rubber shield facing the tube in order to puncture the rubber part of the tube cap. The rubber shield keeps the system tight during a tube change when blood samples are taken into several tubes.

The surface of the needle facing the patient is treated with silicone, Double needle ensures the closeness of the system, in which blood enters the tube without contact with the external environment.

needle case is provided with a perforated label to visually monitor the integrity of the package.
The case has an elongated shape.
The needle tube is color-coded depending on the needle diameter in accordance with ISO 6009.
A flexible rubber shield at the end of the needle, designed to prevent backflow of blood;
Double oblique cut from the sagittal end of the needle; triple spear sharpening laser;
There is a thread on the needle tube for screwing the needle into the needle holder.
Case for color coding of needle diameter.

Flashback needle

Flashback needle has transparent camera for visualizing blood flow.
Flashback needle - one part of the needle is intended for insertion into the patient's vein, the other - closed with a rubber valve, in order to puncture the rubber part of the tube. Between them there is a transparent elongated modified visual camera (1.5 cm), due to the increased visualization of blood flow (control of needle penetration into the patient's vein). The rubber valve keeps the system tight during a tube change when blood samples are taken into several tubes. A double-sided needle ensures the closeness of the system, in which blood enters the tube without contact with the external environment.
The needle case is provided with a perforated label to visually monitor the integrity of the package.
The case has an elongated shape to prevent blunting of the needle.
Case for color coding of needle diameter.

DESIGNATION and INDICATION. FIELD OF ACTION

The product is intended for use in venous blood collection for analysis and / or other research.

Field of application.

The product is intended for use in the outpatient and inpatient hospitals.

CONTRAINDICATIONS

There are no contraindications at this time.

SIDE EFFECTS

There are no side effect at this time.

PRECAUTIONS. WARNINGS

Precautions:

Before and after the procedure, hygienic treatment of hands.

Gloves are required during the procedure.

The use of a non-punctured container for used needles is mandatory.

With the threat of blood splashing, the use of protective equipment is mandatory (mask, goggles, etc.)

Matters need attention (Warnings):

The needle cap shall be pulled out uprightly and make sure it will not touch rubber shield.

Check if the needle and holder match tightly enough. Or else, it may lead to collection problem.

When the package is opened, please check the expiry date and integrity of products. Do not use the product if there is any damage, deformation or caps falling.

The products are for single-use only. Once the package is open, please use the product immediately. Recover the cap after use to prevent pollution. The products should be used and monitored by professional staff.

The tourniquet shall not be released before the needle is pulled out. (Or else, the blood pressure change due to release tourniquet may lead to blood flowing back to vein).

INSTRUCTION FOR USE

- 1) Identify the patient, introduce himself, and explain the course and purpose of the procedure. Ensure that the patient has informed consent for the upcoming blood collection procedure. In the absence of such, clarify further actions with a doctor.
- 2) Offer the patient or help him to take a comfortable position: sitting or lying down.
- 3) Mark the tubes, indicating the patient's full name, department (in order to eliminate errors in the identification of the biomaterial sample).
- 4) Treat your hands hygienically, dry them.
- 5) Treat your hands with an antiseptic. Do not dry, wait until the antiseptic is completely dry.
- 6) Wear non-sterile gloves.
- 7) Select tubes corresponding to the declared tests or laboratory tests required by the patient, prepare a needle, holder, alcohol wipes, and a plaster.
- 8) Place a tourniquet on a shirt or diaper 7-10 cm above the venipuncture site. Harness need to impose no more than one minute. A longer vein squeeze time may affect the test results due to changes in the concentration of certain parameters in the blood.
- 9) Ask the patient to make a fist. It is impossible to set for a hand an exercise stress (vigorous "squeezing and unclenching a fist"), since this may lead to changes in blood concentrations of some indicators. Choose venipuncture site. The most common are the ulnar and saphenous veins, but smaller and fuller veins of the dorsal surface of the wrist.

Performing the procedure:

- 1) Take the needle and remove the protective cap from it. If a double-sided needle is used, remove the protective cap in gray or white.
- 2) Insert the needle into the needle holder and screw in until it stops.
- 3) Disinfect the venipuncture site with a gauze napkin or swab moistened with an antiseptic solution, in a circular motion from the center to the periphery.
- 4) Wait until the antiseptic solution is completely dry (30-60 s). You cannot wipe and blow the puncture site, so as not to bring microorganisms. It is also impossible to palpate the vein after disinfection. If difficulties arise during venipuncture and the vein is palpated again, this area must be disinfected again.
- 5) Remove the cap on the other side of the needle.
- 6) To clasp the patient's forearm with the left hand so that the thumb is 3-5 cm below the venipuncture site, to tighten the skin.
- 7) Position the needle along the line with the vein bevel up and puncture the vein at an angle 15-30 ° to the skin.

Insert the previously prepared tube into the needle holder until it stops and hold it until the blood stops flowing into the tube. The harness must be removed immediately after the start of blood flow into the tube. Make sure the patient opens his fist. The blood passes into the tube until it fully compensates for the vacuum created in it. If the blood does not flow, it means that the needle went through the vein - in this case, you need to slightly pull out the needle, but do not remove it until the blood goes to the tube. The accuracy of filling the tube is $\pm 10\%$ of the nominal volume.

9) Remove the tube from the holder.

10) Immediately after filling, the tube must be gently inverted to mix the sample with the filler: a tube without anticoagulants - 5-6 times; test tube with citrate - 3-4 times; a tube with heparin, EDTA and other additives - 8-10 times. The tubes cannot be shaken - it can cause foaming and hemolysis, as well as lead to mechanical lysis of red blood cells.

If necessary, a number of other tubes are inserted into the needle holder to obtain the desired volume of blood for various studies. Re-insert the needle into the vein for this is not necessary.

11) After all the required tubes are filled, attach a dry sterile cloth to the venipuncture site and remove the needle.

III Termination of the procedure:

1) Place the used needle together with the disposable holder in a sharps container. In order to prevent contact with blood, it is forbidden to disassemble the needle and the holder in your hands!

2) Ensure that the patient does not have external bleeding in the venipuncture area.

3) Apply a pressure bandage on the arm or a bactericidal patch.

4) Handle gloved hands with a disinfectant.

5) Disinfect all used material. Remove gloves, place in a disinfection tank or a waterproof bag / container for disposal of class B waste.

6) Treat your hands hygienically, dry them.

7) Refine patient well-being.

8) Make an appropriate record of the results of the service in the medical records or arrange a referral.

9) Arrange the delivery of test tubes with the obtained laboratory material to the laboratory

Be sure to put on the label on the tube, the name of the medical organization, the patient's name and initials, his gender, age, date and time of blood collection and send the tube to the laboratory. Labels must always be filled in the presence of the patient.

LABELING

The following information is on the sealing label of multi sample and flashback needles:

-needle name

-manufacturer name and logo

-needle size (needle tube outer diameter in mm x needle tube length in mm)

- or code
- thin walled needle
- symbol "H"
- batch code
- manufacturing date
- sterilization date
- shelf life
- words «non pyrogenic», «non toxic»
- Sterility symbol
- Do not reuse symbol
- Do not use if pack is damaged

For group/ transport package of multi sample and flashback needle the following information is:

- device name
- quantity
- needle size (needle tube outer diameter in mm x needle tube length in mm), catheter length (for butterfly needle)
- sterility symbol
- thin walled needle
- batch code
- manufacturing date
- sterilization date
- shelf life
- precautions and storage requirements including:
 - «do not use if package are damaged»
 - «Do not use after expiration date»
 - "Used in conjunction with a vacuum tube"
- words «Non pyrogenic», «Non toxic»
- Symbol "H"
- Do not reuse symbol
- Storage conditions symbols
- Manufacturer trademark symbol;
- Authorized representative name and address in Russia;
- manufacturer name and address;
- Registration Certificate number and date
- Rostest symbol
- transportation conditions symbol
- CE Mark

SYMBOLS ON PACKAGE

Table 1. Symbols description

LOT	Batch code
------------	------------

	Date of manufacturer
	Used-by date
	Manufacturer
	Do not reuse
	Sterilized by radiation
	Do not use if pack is damaged
	Keep away from sunlight
	Keep dry
	Temperature range
	Fragile, handle with care
	Up
	Stacking layer limit
	Rostest symbol
	EC conformity

GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.

Table 2. Guaranteed by manufacturer values of main parameters, characteristics (features) of medical device.

Values of the parameters, characteristics (features) of the medical device are presented in the tables below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	Size	Needle size* (mm)	Color
Multi-sample needle	18G	1,2 mm x 25 mm 1,2 mm x 32 mm 1,2 mm x 38 mm	pink
Multi-sample needle	19G	1,1 mm x 25 mm 1,1 mm x 32 mm 1,1 mm x 38 mm	ivory
Multi-sample needle	20G	0,9 mm x 25 mm 0,9 mm x 32 mm 0,9 mm x 38 mm	Yellow
Multi-sample needle	21G	0,8 mm x 25 mm 0,8 mm x 32 mm 0,8 mm x 38 mm	Dark-green
Multi-sample needle	22G	0,7 mm x 25 mm 0,7 mm x 32 mm 0,7 mm x 38 mm	black
Multi-sample needle	23G	0,6 mm x 25 mm 0,6 mm x 32 mm 0,6 mm x 38 mm	Dark blue
Flashback needle	21G	0,8 mm x 25 mm 0,8 mm x 38 mm	Dark-green
Flashback needle	22G	0,7 mm x 25 mm 0,7 mm x 38 mm	black
Flashback needle	23G	0,6 mm x 25 mm 0,6 mm x 38 mm	Dark blue

*Nominal outer needle tube diameter in mm X nominal needle tube length in mm.

Needle category – thin walled

The needle sharpening angle has cut at an angle of $(11 \pm 2)^\circ$ or a “short” cut at an angle of $(17 \pm 2)^\circ$.

Connection dimensions

Connection length of screw connection (for all needles) - 5 mm $\pm 10\%$

Connection diameter of screw connection (for all needles) - 4,5 mm $\pm 10\%$

The number of turns of the screw connection - 2.5 turns

MAIN RAW MATERIALS

Vein needle tube - Austenitic SS

Hub and wings - ABS

Vein needle tube cap (back) polyethylene PE

needle tube cap (front) Polypropylene
shield - Halogenated butyl rubber

CONDITIONS FOR TRANSPORTATION, STORAGE AND OPERATION. SHELF LIFE

Store this product in a ventilated and dry place at room temperature.

The product should be protected from sunlight, rain, snow during storage. Store separately from toxic or harmful items.

Storage

Temperature	+4 - +25°C
Atmospheric pressure	700 - 1060 hPa
Relative humidity	Not more 75%

Operation

Temperature	15-30° C
Atmospheric pressure	700 - 1060 hPa
Relative humidity	35- 75%

Transportation of products should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the temperature from -25 to +25 °C. During transportation to provide protection from the sun, rain and snow, do not mix with harmful or toxic substances.

Transportation

Temperature	From -25 to +25 °C.
Air pressure	700 - 1060 hPa
Relative humidity	Not more than 80%

The shelf life of the product is 3 years from the date of sterilization. After the expiration date use of the product is prohibited.

WARRANTY OBLIGATION

All questions related to the operation and quality of the product can be addressed to the Authorized representative of the manufacturer in the territory of the Russian Federation.

The manufacturer guarantees that the product meets the requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation.

The manufacturer shall not be liable for any damage or damage caused by the use of the product until proven guilty.

Warranty shelf life is -3 years

The manufacturer shall not be liable for any damage or damage caused by the use of the product after the expiry date indicated on the packaging.

For questions about the quality of a medical product, manufactured by Xiamen Sciendox Biological Technology Co., Ltd, China, please contact:

REAL CARE LLC

Askaya st. 1 building 10 floor 2 office VII, room 3, Moscow Russia 111524
phones: +7(495)234-57-31

POTENTIAL USERS

A patient without age limit.

Specialist with a diploma of a standard pattern on graduating from a secondary vocational medical educational institution in the following specialties: "General Medicine", "Nursing", "Obstetrics". There are skills to perform this simple medical service.

Conditions:

Outpatient

Stationary

STERILITY

The devices are sterile. Sterilization method – radiation.

PACKAGE. SCOPE OF DELIVERY

Multi-sample and flashback needles are packed sealed with label sealing which will keep sterile effectively, box is ivory board with IFU. Boxes are packed in transport carton.

Multi-sample needle package: 100pcs per box with IFU, 50 boxes per transport carton.

Flashback needle package: 50pcs or 100 pcs per box with IFU, 50 boxes per transport carton.

REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not applicable.

REQUIREMENTS FOR THE THE DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

Used needles are potentially hazardous medical waste to the spread of infectious diseases transmitted through blood and are class B medical waste — epidemiologically hazardous waste. Disinfection and disposal of used tubes should be carried out in accordance with the requirements of sanitary regulations of SanPin 2.1.7.2790-10 dated 12.12.2010 «Sanitation and epidemiological requirements for medical waste circulation» and MU 287-113 "Guidelines for disinfection, presterilizing cleaning and sterilization of medical products".

Unused devices, including expired ones, should be disposed of as Class A waste and can be disposed of as household waste.

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

THE NAME AND LEGAL ADDRESS OF THE MANUFACTURER OF MEDICAL DEVICES AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN RF:

MANUFACTURER:

MANUFACTURER:

Hebei Xinle Sci&Tech Co., Ltd, China

No. 189, Nanhuan Road 050700 Xinle City, Hebei Province People's Republic of China

Manufacturing site

Xinle Sci&Tech Co., Ltd, China

189, Nanhuan Road 050700 Xinle City, Hebei Province People's Republic of China

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

REAL CARE LLC

Perovskaya st. 1 building 10 floor 2 office VII, room 3, Moscow Russia 111524

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