SRI TRANG
GLOVES THAILAND

INSTRUCTION FOR USE FOR MEDICAL DEVICE:

**Medical diagnostic (examination)
disposable non-sterile gloves "BENOVY"**

Director

(signature)

/Mr. Chia Chee Ping/

Director

(signature)

/Ms. Thanawan Sa-ngiamsak/

Initial issue of operational documentation (instruction for use for medical device "Medical diagnostic (examination) disposable non-sterile gloves "BENOVY""): 01/07/2019.

2019

Signed before me on 7 May 2020 by
Mr. Chia Chee Ping and Ms. Thanawan Sa-ngiamsak

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TABLE OF CONTENTS

1 Data on the medical device	3
1.1 Name of the medical device.....	3
1.2 Manufacturer (producer) of the medical device.....	4
1.3 Purpose of the medical device	4
1.4 Scope of the medical device.....	4
1.5 Intended users of the medical device	4
1.6 Application conditions of the medical device.....	4
1.7 Indications for use of the medical device.....	4
1.8 Contradictions of the medical device.....	5
1.9 Instruction for use of the medical device	5
1.10 Warning and precautions	6
1.11 Shelf life of the medical device.....	6
2 Transportation, storage and disposal conditions of the medical device.....	7
2.1 Transportation conditions of the medical device	7
2.2 Storage conditions of the medical device	7
2.3 Disposal conditions of the medical device.....	7
3 Disinfection and cleaning methods & equipment of the medical device.....	7
4 Specification of the medical device	8
4.1 Description of the medical device.....	8
4.2 Physical specification of the medical device	10
4.2.1 Physical and mechanical performance of the medical device.....	10
4.2.2 Inspection levels and acceptable quality level of the medical device.....	11
5 Adverse reactions.....	12
6 Circumstances, in which the consumer must consult healthcare professional	12
7 Manufacturer's warranties	12
8 Data on Authorized representative of the manufacturer	12

1 Data on the medical device

1.1 Name of the medical device

Medical diagnostic (examination) disposable non-sterile gloves “BENOVY”, in hardware versions:

1. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex smooth powdered (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 5-P. Type 122560.
2. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex finger textured powdered (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 5-P-FT. Type 122560.
3. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex finger textured powdered (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 58-P-FT. Type 122560.
4. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex finger textured powder free chlorinated (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 5-1C. Type 122540.
5. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex finger textured High Density powder free with inner polymer coating with double chlorinated outer surface (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 62-1P/2C. Type 122540.
6. Medical diagnostic (examination) disposable non-sterile gloves “BENOVY” natural latex textured powder free with inner polymer coating (beaded cuff), color: natural, size: XS, S, M, L, XL. Series: 58-1P. Type 122540.

1.2 Manufacturer (producer) of the medical device

Manufacturer (producer) of the medical device is:

Sri Trang Gloves (Thailand) Public Company Limited
110 Kanjanavanit Road, Pahtong, Hat Yai, Songkhla 90230, Thailand

Addresses of the medical device manufacture:

110 Kanjanavanit Road, Pahtong, Hat Yai, Songkhla 90230, Thailand
109/2 Kanjanavanit Road, Pahtong, Hat Yai, Songkhla 90230, Thailand
352 Kanjanavanit Road, Pahtong, Hat Yai, Songkhla 90230, Thailand
189, Moo 7, Phlai Wat, Kanchanadit, Surat Thani 84160, Thailand
85 Moo 6, Khuan Thani, Kantang, Trang 92110, Thailand

1.3 Purpose of the medical device

The gloves are intended for single use as protective barrier in the process of wearing on healthcare professional's hands during examination/treatment of patient or for other medical goals.

1.4 Scope of the medical device

Medical gloves are used during medical examinations, diagnostic and therapeutic procedures, as well as for handling contaminated medical materials.

1.5 Intended users of the medical device

Medical personnel, who perform medical examinations, diagnostic and therapeutic procedures, as well as personnel, who work with contaminated medical materials.

1.6 Application conditions of the medical device

Medical examination disposable gloves are used in medical or healthcare facilities indoors in a dry place, away from direct sunlight, at temperature from 10°C to 30°C and relative humidity of not more than 80% in accordance with the instruction for use on medical device.

Gloves for external use only (the device is non-invasive).

1.7 Indications for use of the medical device

The indications for use of the medical device are the necessity of conduction of medical examinations, diagnostic and therapeutic procedures.

1.8 Contradictions of the medical device

Contradictions of the medical device are as follows:

1. The patient or healthcare professional has an allergic reaction to natural latex proteins;
2. The patient or healthcare professional has an allergic reaction to the powder (applicable only to powdered gloves);
3. The patient or healthcare professional has an allergic reaction to residual process chemicals in the medical device.

1.9 Instruction for use of the medical device

How to put the gloves on:

- 1) Make sure to wash your hands before putting the gloves on or changing them.
- 2) Take a glove from the box. It is recommended to touch the bounded surface of the glove corresponding to the wrist (upper edge of the cuff).
- 3) Pull on the first glove.
- 4) Take the second glove with your bare hand, touching the bounded surface of the glove corresponding to the wrist.
- 5) To avoid touching the skin of the forearm with a gloved hand, rotate the outer surface of the glove to be put on with folded fingers in the glove so that you could pull the glove over the other hand.
- 6) Having put the gloves on, touch only those objects provided for by the indications for use and conditions of the use of the gloves.
- 7) To avoid contamination, change gloves before proceeding with the new procedure.

How to take the gloves off:

- 1) To take a glove off, pick it up at the level of the wrist without touching the skin of the forearm and pull it off your hand, turning it inside out.
- 2) Holding the removed glove with the gloved hand, it is necessary to place the fingers of your bare hand between the glove and wrist. Remove the second glove by rolling it down the arm and put into the first glove.
- 3) Discard the gloves and wash your hands thoroughly.

1.10 Warning and precautions

- For single use only;
- For individual use only;
- Every time it is necessary to perform hygienic treatment of your hands before putting on or changing gloves, as well as after removing the gloves;
- It is necessary to choose the gloves of the right size – they should be easily put on and should not impede movement or squeeze the hands. Do not use in case of size mismatch;
- It is necessary to choose the gloves depending on the type of the performed manipulations;
- Touching is permissible only to the bounded surface of the glove corresponding to the wrist (upper edge of the cuff);
- Depending on the type of the performed manipulation, to keep an optimal level of protection, it is recommended to change gloves every 60-120 minutes;
- It is necessary to change the gloves in case of any defect, puncture, damage, contamination, fluid contact with skin under the glove, as well as stickiness. If one glove is damaged, change both ones;
- It is necessary to put the gloves on dry hands only;
- Do not wet and keep away from moisture;
- Avoid exposure to direct sunlight;
- The composition of the medical device includes natural rubber latex that may cause an allergic reaction;
- Powder may cause an allergic reaction (applicable only to powdered gloves);
- Gloves contain residual process chemicals that may cause an allergic reaction;
- Gloves do not guarantee 100% hand protection against punctures, cuts and damage caused by exposure to chemicals.

1.11 Shelf life of the medical device

The established shelf life of the examination disposable non-sterile powdered natural latex gloves is 5 years.

The established shelf life of the examination disposable non-sterile powder free natural latex gloves is 3 years.

2 Transportation, storage and disposal conditions of the medical device

2.1 Transportation conditions of the medical device

There are no special requirements for transportation other than requirements for carriage by land, sea and air. During transportation, the storage conditions should not be violated.

2.2 Storage conditions of the medical device

Examination gloves are stored indoors in a dry place, away from direct sunlight, heat, ignition, ultraviolet and ozone sources, at temperature of 10-30°C and relative humidity of not more than 80%.

2.3 Disposal conditions of the medical device

Disposal of the medical device “Medical diagnostic (examination) disposable non-sterile gloves “BENOVY”” in the Russian Federation is carried out by the medical organization as class B waste in accordance with SanPiN 2.1.7.2790-10 “Sanitary-epidemiological requirements for medical waste management” and approved guideline of the organization.

3 Disinfection and cleaning methods & equipment of the medical device

Disposable examination non-sterile gloves are disposable medical device. Reuse of the medical device is a gross violation of application methods.

4 Specification of the medical device

4.1 Description of the medical device

Form. Straight fingers, thumb and other fingers in the same plane, can be used for either hand (ambidextrous). Rolled cuff.

Surface. Powdered gloves – modified corn starch is used; powder-free gloves are treated with chlorine – chlorinated and/or polymer coating treated.

Texture of the outer surface. Gloves can be smooth, textured or finger textured.

Sizes. Depending on the size, table 1 gives the values of width, length and thickness of the medical device.

Table 1 – Sizes of the medical device:

Value name	Size					
	XS	S	M	L	XL	
Width W (figure 1), mm	≤ 80	80±10	95±10	110±10	≥110	
Length l (figure 1), mm, min	220		230			
Thickness (in points, shown on figure 2), mm, min						
- smooth surface						0.08
- textured surface						0.11
Thickness of rolled cuff, mm, max	2.50					


Figure 1 – Points for width and length measurement of the glove


Figure 2 – Points for wall thickness measurement of the glove

4.2 Physical specification of the medical device

4.2.1 Physical and mechanical performance of the medical device

Physical and mechanical parameters of the medical device are represented in table 2.

Table 2 – Values of physical and mechanical parameters of gloves:

DESCRIPTION	SPECIFICATION
1. Before accelerated aging:	
- force at break, N, min	7.0
- elongation at break, %, min	650
- tensile strength, MPa, min	18
2. After accelerated aging:	
- force at break, N, min	6.0
- elongation at break, %, min	500
- tensile strength, MPa, min	14

4.2.2 Inspection levels and acceptable quality level of the medical device

For the purposes of quality control medical gloves are selected and examined in accordance with standard ISO 2859-1.

Hermeticity. Hermeticity is determined in accordance with EN 455-1. Medical gloves are hermetic. Acceptable quality level (AQL) is ≤ 1.5 .

Powder. The content of residual powder is determined in accordance with the requirements of ASTM D 6124.

- For powdered gloves, the value obtained does not exceed 10 mg/dm^2 per glove.
- For powder-free gloves, the value obtained does not exceed 2 mg per glove.

The protein content in latex gloves.

- The content of water extractable proteins is not more than $200 \text{ } \mu\text{g/dm}^2$ when tested in accordance with ASTM D 5712.
- The content of antigenic proteins is not more than $10 \text{ } \mu\text{g/dm}^2$ when tested in accordance with ASTM D 6499.

5 Adverse reactions

Manifestation of dermal hypersensitivity reactions, when using the device is possible.

6 Circumstances, in which the consumer must consult healthcare professional

The patient is recommended to consult healthcare professional in the event of dermal hypersensitivity reactions, when using the device.

7 Manufacturer's warranties

The manufacturer (producer) guarantees preservation of all the parameters and characteristics of the medical device taking into account acceptable quality level (AQL) ≤ 1.5 , as well as the quality of the medical gloves within all shelf life from the date of manufacture upon the condition of package integrity, observance of transportation, storage and operation conditions in accordance with the instruction.

8 Data on Authorized representative of the manufacturer

Please contact the Authorized representative of the manufacturer for any and all purposes, related to circulation of the medical device on the territory of the Russian Federation:

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