160472

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Shanghai MicroPort Medical (Group) Co., Ltd. on the annexed INSTRUCTIONS FOR USE is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Zhou Pingfan

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(Date: Aug. 17, 2020)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

I affirm
Dr. Zhaohua Chang
Chairman & CEO
Shanghai MicroPort Medical (Group) Co., Ltd

10.08.2020

INSTRUCTIONS FOR USE

WALTZ™ CoCr Coronary Stent System

Shanghai MicroPort Medical (Group) Co., Ltd.

NOTICE: The information contained herein is a **CONFIDENTIAL AND PROTECTED PROPERTY RIGHT**. The contents of this document may not be reproduced or made available to third parties without the prior consent of Shanghai MicroPort Medical (Group) Co., Ltd and should not be used in any unauthorized way.

This operational documentation has been prepared to provide information for the registration of a medical device in the territory of the Russian Federation. Information is presented in the amount required in accordance with the Rules for the registration of medical devices (Decree of the Government of the Russian Federation dated December 27, 2012 No. 1416), Order of the Ministry of Health of the Russian Federation dated January 19, 2017 No. 11n "On approval of the requirements for the content of technical and operational documentation of the manufacturer (manufacturer) a medical device."

INSTRUCTIONS FOR USE

WALTZ™ CoCr Coronary Stent System

Before use, read these Instructions and pay special attention to all points from the sections Warnings and Precautions!

1 Introduction

The WALTZ™ CoCr Coronary Stent System (hereinafter WALTZ™ Stent System) includes:

- WALTZ™ stent – a pre-mounted L605 Cobalt Chromium alloy stent
- Delivery system – Rapid exchange (RX) Coronary Catheter

The delivery system is a quick change catheter (RX). The catheter has a semi-compliant balloon on which the stent is mounted. The catheter has two radiopaque labels for marking the ends of the stent to facilitate proper placement of the stent during PCI. Two radiopaque marks located under the balloon and marking the working length of the balloon and the length of the stent when opened.

Two proximal axial marks on the delivery system (at a distance of 95 cm and 105 cm from the distal end), showing the relative position of the delivery system to the end of the shoulder or femoral guide catheter. The working length of the catheter is 140 cm.

The sketch of WALTZ™ Stent System can refer to Fig.1.


Fig.1 WALTZ™ Stent System.

User

The WALTZ™ stent can be used by well-trained personnel under the guidance of a qualified physician who knows the general benefits and risks of working with stent systems.

Type of contact with the body

Stent - lifelong contact with the internal environment of the body (with blood or lymph).
Quick change balloon catheter for stent delivery - short-term contact with blood and lymph.

The characteristics of the device components are given in Table 1.1, Table 1.2.

Table 1.1: Product Description WALTZ™ Stent

Component	Characteristic
Delivery system	Rapid Replacement Balloon Catheter (RX)
Stent delivery system balloon	Semi-compliant balloon with two radiopaque marks for stent marking
Stent Length Options (mm)	8; 13; 16; 18; 21; 23; 26; 29 31; 33

Component		Characteristic
Stent Diameter Options (mm)		2,25; 2,5; 2,75; 3,0; 3,5; 4,0 2,5; 2,75; 3,0; 3,5; 4,0
Stent material		Medical Cobalt Chrome Alloy L-605
Balloon material		Polymer brand Pebax
X-ray contrast strip		Pt-Ir Alloy
Hydrophilic coating		Polyvinylpyrrolidone (DSM and JMEDTECH)
Delivery system working length		1400±30 mm
Distances to the distal tip of the catheter	Guide conductor outlet	280±15 mm
	The first mark on the rod	950±20 mm
	The second mark on the rod	1050±20 mm
Stent Platform Thickness		86.4±12.7 MKM
Soft tip length, mm		3,0 ±0,2 mm
Outer diameter of soft tip, mm		0.5334 ± 0.0127 mm
Delivery system design		Single port for access to the inflation hole. The outlet of the conductor is located at a distance of 28 cm from the end. For conductors of size ≤ 0.014" (0.356 mm)
Stent delivery system balloon		The length and position of the installed stent is determined by the boundaries of 2 radiopaque platinum-iridium marks.
Cylinder inflation pressure		Nominal inflation pressure: 9 atm (912 kPa) for all sizes; Estimated Burst Pressure (RDR): 16 atm (1621 kPa) for diameters from 2.25 to 3.5 mm; 14 atm (1419 kPa) for a diameter of 4.0 mm.
Guide Catheter Inner Diameter		≥ 6F (0,070"/1,78mm)

Table 1.2: Specifications

Parameter	Value
Stent	
Surface area without stent and external stent area	For stents with the specified diameter, minimum and maximum diameter, the surface area without a stent should be more than 80% of the total surface area of the cylinders with the same diameter and length.

Parameter	Value														
Compression Resistance at Radially Applied Load	The table below shows the values of compression resistance for stents of different diameters at 12% deformation.														
	<table><tr><th>Nominal diameter (mm)</th><th>Compression resistance at 12% stent deformation (kPa)</th></tr><tr><td>Ø2,25</td><td>≥ 110</td></tr><tr><td>Ø2,5</td><td>≥ 105</td></tr><tr><td>Ø2,75</td><td>≥ 100</td></tr><tr><td>Ø3,0</td><td>≥ 95</td></tr><tr><td>Ø3,5</td><td>≥ 90</td></tr><tr><td>Ø4,0</td><td>≥ 85</td></tr></table>	Nominal diameter (mm)	Compression resistance at 12% stent deformation (kPa)	Ø2,25	≥ 110	Ø2,5	≥ 105	Ø2,75	≥ 100	Ø3,0	≥ 95	Ø3,5	≥ 90	Ø4,0	≥ 85
	Nominal diameter (mm)	Compression resistance at 12% stent deformation (kPa)													
	Ø2,25	≥ 110													
	Ø2,5	≥ 105													
	Ø2,75	≥ 100													
	Ø3,0	≥ 95													
	Ø3,5	≥ 90													
Ø4,0	≥ 85														
Corrosion	Corrosion susceptibility: breakdown potential should be ≥ 600 mV, protection potential should be ≥ 300 mV Corrosion abrasion (carried out after fatigue test): the stent's galvanic potential must be ≥ 300 mV														
Stent length to diameter ratio	When expanding the stent from minimum to maximum diameter, the change in stent length should be ≤ 20%.														
Effect / Extension Profile	1) After the control process, the delivery system has no visible damage (for example, kinks), stents do not detach and do not move; 2) The distance between the outer diameter of the balloon and the outer diameter of the corrugated stent profile should be ≤ 5% both before and after tracking, and more than 8 mm of curvature.														
Adherence to the vessel wall	The outer surface of the stent should fit snugly against the vessel wall; deviation and separation of the stent from the vessel wall is not allowed.														
Recoil	≤ 10%														
Stress / Strain Analysis Determination of safety factor fatigue strength Fatigue Test	The stent should maintain its integrity without any obvious defects or fractures after a fatigue test under physiological modeling for 10 years (380,000,000 cycles).														
Stent diameter at balloon pressure	Nominal pressure 9 atm														
	<table><tr><th>Nominal diameter (mm)</th><th>Diameter of the expanded stent (mm)</th></tr><tr><td>Ø2,25</td><td>2,25±0,22</td></tr><tr><td>Ø2,5</td><td>2,50±0,25</td></tr><tr><td>Ø2,75</td><td>2,75±0,27</td></tr><tr><td>Ø3,0</td><td>3,00±0,30</td></tr><tr><td>Ø3,5</td><td>3,50±0,35</td></tr><tr><td>Ø4,0</td><td>4,00±0,40</td></tr></table>	Nominal diameter (mm)	Diameter of the expanded stent (mm)	Ø2,25	2,25±0,22	Ø2,5	2,50±0,25	Ø2,75	2,75±0,27	Ø3,0	3,00±0,30	Ø3,5	3,50±0,35	Ø4,0	4,00±0,40
	Nominal diameter (mm)	Diameter of the expanded stent (mm)													
	Ø2,25	2,25±0,22													
	Ø2,5	2,50±0,25													
	Ø2,75	2,75±0,27													
	Ø3,0	3,00±0,30													
	Ø3,5	3,50±0,35													
Ø4,0	4,00±0,40														

Parameter	Value	
Bending / twisting	The minimum radius of curvature of the expanded stent should be ≤ 8.0 mm.	
Visibility	Stents should be visible with digital subtraction angiography (CSA) or fluoroscopy.	
Bond strength Tensile Strength	Validation items	Minimum fracture force (N)
	Connection of the distal part of the cylinder	2
	Connection of the proximal part of the balloon	5
	The distal part of the shaft	5
	Adapter tube	5
	Inside	3
	Conductor outlet	5
	Sheath compound	5
	Connection between Luer tip and strain relief tube	10
	The connection between the proximal part of the rod and the tube to relieve tension	10
Appearance: stent system	The surface of the delivery system should be smooth without any visible foreign particles; There should be no twisting (especially around the guide hole). There should be no damage at the end. The stent frame elements do not stagger or overlap, the frame does not tip over, there are no cracks, and there are no visible foreign substances.	
Component Size Compatibility	The maximum outer diameter of the stent system.	See profile / diameter dimensions
	Lumen diameter and conductor smoothness	The conductor should move smoothly without any obstacles through the clearance of the guide; maximum pushing force must not exceed 0.3 N.
	Working length	1400 $\pm$ 30 mm
	Other sizes	See sizing
	Compatibility with Recommended Accessories	A simulated use test must pass successfully. There should be no leakage of contrast medium at the connection points. The contrast medium should easily enter the stent system balloon and effectively expand the stent. Verification of hemostasis should be successfully completed. After the product has passed the three-way hemostatic valve, there should be no kink in the delivery system or deformation of the stent, and

Parameter	Value		
	it should be possible to successfully remove the product through the three-way hemostatic valve.		
Stent Profile / Diameter	Corrugated Stent Profile: Stent Size (mm)	Profile (inches)	Profile (mm)
	Ø2,25	≤ 0,042"	1,0668
	Ø2,5	≤ 0,043"	1,0922
	Ø2,75	≤ 0,045"	1,143
	Ø3,0	≤ 0,045"	1,143
	Ø3,5	≤ 0,049"	1,2446
	Ø4,0 mm	≤ 0,050"	1,27
Displacement force	① The initial displacement force must be at least 2.0 N with an initial movement of 0-5 mm (the initial displacement force is the force required to create the displacement). ② Record the maximum displacement force (maximum displacement force is the force required to completely separate the stent from the delivery system).		
A leak	When filling the system with a contrast agent, there should be no leaks or damage, such as tears or protrusions on the delivery system.		
	There should be no leaks in the Luer tip, in any area of the joint or any other elements.		
	Air must not enter through the assembly points of the delivery system during the suction process.		
Bending / twisting	The stent / delivery system should not have kinks during delivery and retrieval in a simulated use test.		
Pushing ability	During the simulation of the delivery process, the stent system is combined with a guiding catheter, a conductor barrel and can smoothly deliver the stent to the desired position.		
Imitation use	1) Delivery process: The delivery system matches the wire guide very well, smoothly delivers the stent to its destination and crosses the winding sections of the model. The delivery system can move along the conductor smoothly, without any significant pauses. The system has no kinks in the delivery process. 2) The process of expansion: The stent is easily stretched and straightened at nominal pressure and has a good visible shape and does not fall off after expansion. There are no leaks during expansion. 3) The extraction process: The catheter can be removed smoothly while the stent does not move, and the delivery system does not have kinks in the extraction process. 4) Compatibility with the walls of the vessel: The external surface of the stent after deployment corresponds to the wall of the vessel.		

Parameter	Value		
	the stent does not protrude and does not separate from the walls of the vessel.		
Traceability	During simulated use, the stent and delivery system can smoothly move along the guide conductor, cross the winding sections of the model, ultimately reaching the target area.		
Stent shortening	Stent shortening	Nominal Length (mm)	
	≤12%	8≤L<20	
	≤10%	20≤L<30	
	≤8%	L≥30	
Delivery system and stent system			
Visibility	X-ray contrast strips of a balloon catheter should be visible during fluoroscopy.		
Balloon filling time	Tank filling time ≤ 20 seconds		
Balloon Duration	Nominal diameter (mm)	Stent length (mm)	
		8, 13, 16, 18, 21, 23, 26, 29	31, 33
	2,25	≤ 9 сек	—
	2,5	≤ 9 сек	≤ 12 сек
	2,75	≤ 9 сек	≤ 12 сек
	3,0	≤ 9 сек	≤ 12 сек
	3,5	≤ 9 сек	≤ 15 сек
	4,0	≤ 11 сек	≤ 15 сек
Estimated burst pressure (RDR) of the cylinder	Estimated burst pressure (RDR) is based on test results that show that statistically 99.9% of cylinders will not burst at this or lower pressure with a confidence of at least 95%. For stents with diameters from Ø2.25 to Ø3.5 mm, the balloon RDR should be ≥16 atm. For stents with a diameter of Ø4.0 mm, the balloon RDR should be ≥14 atm.		
Estimated fatigue characteristics of the cylinder	The delivery system must not leak or break after 10 cycles of filling and emptying.		
Excessive balloon expansion at both ends of the stent	The diameter of the balloon extending beyond the ends of the stent should be no more than 5% of the outer diameter of the stent at the design burst pressure.		
Hemostasis	At a pressure of 25 kPa and a pushing force of 2 ± 0.2 N, the volume of water flowing out of the hemostatic valve in 20 minutes should not exceed 10 ml.		
Outer diameter of the catheter shaft (nominal)	Distal part of the shaft: Ø 2.25~2.5 mm: 2.7F (0.035″[0,00138mm])		

Parameter	Value
	$\varnothing$ 2.75~3.0 мм: 2.8F (0.037''[0,00146мм]) $\varnothing$ 3.50~4.0 мм: 3.0F (0.039''[0,00153мм]) Proximal part of the shaft: Other: 2.3F (0.030''[0,00118мм]) $\varnothing$ 3.50~4.0 мм с длиной 29, 31 и 33: 2.6F (0.034''[0,00134мм])

Stent Functioning Factors

- The degree of probability of occurrence of contact, galvanic, pitting and fretting corrosion, the possibility of oxidation: the probability is low, provided by the use of bioinert materials
- Danger of tissue perforation: Perforation or rupture of tissues can occur if a stent is used that is not suitable for the diameter of the vessel, or if precautions regarding the use of the stent are not followed.
- The level of force transmitted by the stent to the surrounding tissue: 0.10-0.15 N / mm
- The period of expected stent ingrowth, penetration into tissues, tissue perforation, stent displacement and migration - about 4 weeks
- Penetration in tissue 50%
- Information on biodegradation: the probability is low, provided by the use of bioinert materials.
- Minimum force to remove the system: at least 2N.

DELIVERY FORM

- **STERILE: ONLY FOR ONE TIME USE.** This product is sterilized with ethylene oxide (EtO). Pyrogen-free.
- Do not use if the packaging has been opened or damaged.
- Do not use if label is incomplete or illegible.
- Do not re-sterilize.

CONTENTS OF DELIVERY:

One (1) WALTZ™ Stent.

One (1) Directions for Use

PERFORMANCE OPTIONS:

Diameter (mm)	Length (mm)									
	8	13	16	18	21	23	26	29	31	33
$\varnothing$ 2.25	TG2208	TG2213	TG2216	TG2218	TG2221	TG2223	TG2226	TG2229	×	×
$\varnothing$ 2.5	TG2508	TG2513	TG2516	TG2518	TG2521	TG2523	TG2526	TG2529	TG2531	TG2533
$\varnothing$ 2.75	TG2708	TG2713	TG2716	TG2718	TG2721	TG2723	TG2726	TG2729	TG2731	TG2733
$\varnothing$ 3.0	TG3008	TG3013	TG3016	TG3018	TG3021	TG3023	TG3026	TG3029	TG3031	TG3033
$\varnothing$ 3.5	TG3508	TG3513	TG3516	TG3518	TG3521	TG3523	TG3526	TG3529	TG3531	TG3533
$\varnothing$ 4.0	TG4008	TG4013	TG4016	TG4018	TG4021	TG4023	TG4026	TG4029	TG4031	TG4033

«×»: Not available

STORAGE: Store in a dry, dark, well-ventilated place.

2 PURPOSE AND SCOPE

Purpose: To increase the diameter of the lumen of the coronary vessels in patients with symptomatic coronary artery disease when the size of the patient's coronary artery lesion is not more than 30 mm long with a vessel diameter from 2.25 to 4.0 mm.

Scope: medical and medical-preventive medical institutions.

3 INDICATIONS

The WALTZ™ Stent System is indicated for improving coronary luminal diameter in the following:

- Patients with symptomatic ischemic heart disease due to discrete *de novo* or restenotic native coronary artery lesions (length < 30 mm) with reference vessel diameters ranging from 2.25 mm to 4.0 mm.
- Patients with symptomatic ischemic heart disease due to lesions in saphenous vein bypass grafts (length < 30 mm) with reference vessel diameters ranging from 2.25 mm to 4.0 mm.
- Restoring coronary flow in patients experiencing acute myocardial infarction, as confirmed by ST segment elevation or angiographic findings, who present within 12 hours of symptom onset with native coronary artery lesions of length < 30 mm with a reference vessel diameter of 2.25 mm to 4.0 mm.

4 Contraindications

The WALTZ™ Stent System is generally contraindicated for use in:

- Patients with a hypersensitivity to cobalt chromium;
- Patients with showing active full or local serious infection focus at beginning;
- Patients with stomach haemorrhage or active gastric ulcer recently;
- Patients with brain stroke in six months;
- Patients who have aggressive focus, cavitory or pulmonary consolidation were found in chest X-ray during patient choice of OPS;
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated;
- Patients with coronary artery reference vessel diameter < 2.25 mm or > 4.0mm;
- Patients with unresolved vessel thrombus at the lesion site;
- Patients with the distal vessel lesions or a few of vascular bed in the end of lesions;
- Patients with coronary artery myocardial bridge oppressed;
- Patients with simple coronary artery cramp;
- Patients with diffuse disease to the identified lesions;
- Patients with serious calcified lesions;
- Patients with tortuous vessels in the lesion;
- Other patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or stent delivery.

5 WARNINGS

- Before use, check the integrity and tightness of the sterile inner packaging.
- Careful selection of patients is necessary, as the use of the product may be associated with stent thrombosis, vascular complications and / or bleeding.
- Patients with hypersensitivity to the cobalt-chromium alloy L-605 (including cobalt, chromium, tungsten, nickel) may experience an allergic reaction caused by this implant.

- The use of this product is unacceptable in patients who are not able to receive the recommended therapy with antiplatelet drugs.

6 PRECAUTIONS

6.1 General Precautions

- Stent implantation should only be performed by physicians who have received appropriate training.
- Stent placement should be performed at hospitals where emergency coronary artery bypass graft surgery is accessible.
- Subsequent stent blockage may require repeat dilatation of the arterial segment containing the stent. Long-term outcomes following repeat dilatation of the stent are presently unknown.
- Risks and benefits should be considered in patients with severe contrast agent allergies.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and delivery system withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.

6.2 Use of Multiple Stents

When multiple stents are required, stent materials should be of similar composition. Placing multiple stents of different materials in contact with each other may increase potential for corrosion. To avoid the possibility of dissimilar metal corrosion, do not implant stents of different materials in tandem where overlap or contact is possible.

Potential interactions of the WALTZ™ stent with other bare metal stents have not been evaluated. When multiple bare metal stents are required, use only WALTZ™ stents in order to avoid potential interactions with other stents.

6.3 Brachytherapy

The safety and effectiveness of the WALTZ™ stent in patients with prior brachytherapy of the target lesion have not been established. The safety and effectiveness of the use of brachytherapy to treat in-stent restenosis in a WALTZ™ stent has not been established. Both vascular brachytherapy and the WALTZ™ stent alter arterial remodeling. The synergy between these two treatments has not been determined.

6.4 Use in Conjunction with Other Procedures

The safety and effectiveness of using mechanical atherectomy devices (such as directional atherectomy catheters and rotational atherectomy catheters), or laser angioplasty catheters in conjunction with WALTZ™ stent implantation have not been established.

6.5 Use in Special Populations

6.5.1 Pregnancy

The WALTZ™ stent has not been tested in pregnant women or in men intending to father children. Effects on the developing fetus have not been studied. Effective contraception should

be initiated before implanting a WALTZ™ stent and continued for one year after implantation. While there is no contraindication, the risks and reproductive effects are unknown at this time.

6.5.2 Lactation

A decision should be made whether to discontinue nursing prior to stent implantation considering the importance of the stent to the mother.

6.5.3 Pediatric Use

The safety and effectiveness of the WALTZ™ stent in pediatric patients below the age of 18 years have not been established.

6.5.4 Geriatric Use

Clinical studies of the WALTZ™ stent did not suggest that patients age 65 years and over differed with regard to safety and efficacy compared to younger patients.

6.6 Magnetic Resonance Imaging (MRI)

Non-clinical testing has demonstrated that the WALTZ™ stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 1.5 or 3 Tesla
- Spatial gradient field of 810 Gauss/cm or less
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode) for 15 minutes of scanning or less.

The WALTZ™ stent should not migrate in this MRI environment. Non-clinical testing at field strengths greater than 3 Tesla has not been performed to evaluate stent migration or heating. MRI at 1.5 or 3 Tesla may be performed immediately following the implantation of the WALTZ™ stent.

Stent heating was derived by using the measured non-clinical, in vitro temperature rises in a GE 3 Tesla, Signa HDxt scanner in combination with the local specific absorption rates (SARs) in a digitized human heart model. The temperature rise was derived by validated calculation. In this testing, single WALTZ™ stent up to 33 mm in length produced a temperature rise of less than 2.8 °C at a SAR of 2.0 W/kg for 15 minutes of MR scanning. These calculations do not take into consideration the cooling effects of blood flow.

The effects of MRI on stents with fractured struts are unknown.

As demonstrated in non-clinical testing, an image artifact can be present when scanning the WALTZ™ stent. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

6.7 Stent Handling Precautions

- **Each stent is for single use only.** Do not resterilize or reuse this device. Note the "Use by"

(expiration) date on the product label.

- **Do not remove the stent from the delivery system. Removal may damage the stent and/or lead to stent embolization. These components are intended to perform together as a system.**
- **Do not induce a vacuum on the delivery system prior to reaching the target lesion. This may cause dislodgement of the stent from the balloon.**
- **The delivery system should not be used in conjunction with other stents.**
- **Special care must be taken not to handle or disrupt the stent on the balloon especially during delivery system removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.**
- **Do not manipulate, touch, or handle the stent with your fingers, which may cause contamination or stent dislodgement from the delivery balloon.**
- **Use only the appropriate balloon inflation media (see Section 9.3.3 – Operator’s Instructions, Delivery System Preparation). Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in stent deployment.**

6.8 Stent Placement Precautions

6.8.1 Stent Preparation

- **Do not prepare or pre-inflate the delivery system prior to stent deployment other than as directed. Use the balloon purging technique described in Section 9.3.3 – Operator’s Instructions, Delivery System Preparation.**
- **Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system (see Section 1.0 – Product Description).**

6.8.2 Stent Implantation

- **The vessel should be pre-dilated with an appropriate sized balloon. Failure to do so may increase the difficulty of stent placement and cause procedural complications.**
- **Do not expand the stent if it is not properly positioned in the vessel (see Section 6.9 – Precautions, Stent System Removal).**
- **Implanting a stent may lead to vessel dissection and acute closure requiring additional intervention (CABG, further dilatation, placement of additional stents, or other).**
- **Although the safety and effectiveness of treating same coronary artery with more than one WALTZ™ stents has not been established, if this is performed, place the stent in the distal lesion before the proximal lesion in order to minimize dislodgement risk incurred by traversing through deployed stents.**
- **Stent placement may compromise side branch patency.**
- **Do not exceed Rated Burst Pressure (RBP) as indicated on product label. See Table 9.1, Typical WALTZ™ Compliance. Balloon pressures should be monitored during inflation. Applying pressures higher than specified on the product label may result in a ruptured balloon with possible arterial damage and dissection.**
- **An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded**

stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed as the stent may be damaged when retracting the undeployed stent back into the guiding catheter.

- Should any resistance be felt at any time during coronary stent system withdrawal, the stent delivery system and guiding catheter should be removed as a single unit (see Section 6.9 – Precautions, Stent System Removal).
- Stent retrieval methods (i.e., using additional wires, snares, and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.
- Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion, could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and, if necessary, a new dilatation catheter exchanged over the guide wire to complete the expansion of the stent.

6.9 Stent/System Removal Precautions

Should any resistance be felt at any time during either lesion access or removing the delivery system post-stent implantation, the stent delivery system and the guiding catheter should be removed as a single unit.

When removing the delivery system and guiding catheter as a single unit, the following steps should be executed under direct visualization using fluoroscopy:

- Confirm complete balloon deflation. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to the guiding catheter position. In some cases it may be necessary to slightly retract the guiding catheter in order to prevent unplanned guiding catheter movement and subsequent vessel damage. In cases where unplanned guiding catheter movement has occurred, a coronary tree angiographic assessment should be undertaken to ensure that there is no damage to the coronary vasculature.
- DO NOT retract the delivery system into the guiding catheter.
- Position the proximal balloon marker just distal to the guiding catheter tip.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, and remove the guiding catheter and delivery system as a single unit.

Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in loss or damage to the stent and/or delivery system components and/or vasculature.

If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

6.10 Post-Procedure Precautions

- Care should be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, or balloon catheter to avoid disrupting the stent placement, apposition, geometry.
- Post-dilatation: All efforts should be made to assure that the stent is not under-dilated. If the deployed stent is not fully apposed to the vessel wall, the stent may be expanded further with a larger diameter balloon. The post-dilatation can be done using a low-profile, high pressure, non-compliant balloon catheter. The balloon should not extend outside of the stented region (see Operator's Manual – 9.7 Post-Deployment Dilatation of Stent Segments).

Non-clinical testing on single and overlapped stents has demonstrated that the WALTZ™ stent is MR Conditional (see Section 6.6 Magnetic Resonance Imaging (MRI)).

7 ADVERSE EVENTS

Adverse events may be associated with the use of a coronary stent in native coronary arteries:

- Acute myocardial infarction
- Arrhythmias, including VF and VT
- Death
- Dissection
- Emboli, distal (air, tissue or thrombotic emboli)
- Emergent coronary artery bypass surgery
- Hemorrhage, requiring transfusion
- Hypotension / hypertension
- Infection and / or pain at insertion site
- Ischemia, myocardial
- Perforation
- Pseudoaneurysm, femoral
- Restenosis of stented segment
- Spasm
- Stent embolization
- Stent thrombosis / occlusion
- Stroke / cerebrovascular accident
- Total occlusion of coronary artery

7. INDIVIDUALIZATION OF TREATMENT

The risks and benefits described above should be considered carefully for each patient before use of the WALTZ™ stent. Patient selection factors to be assessed should include a judgment regarding the risk of antiplatelet therapy. Stenting is generally avoided in those patients at

heightened risk of bleeding (e.g., those patients with recently active gastritis or peptic ulcer disease, see **Section 4 – Contraindications**).

Premorbid conditions that increase the risk of a poor initial result and the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed. The relation of baseline and procedural variables to Major Adverse Cardiac Events (MACE) was examined. Multivariable modeling suggested that treatment assignment remained an independent predictor of clinical and angiographic outcomes even after adjusting for other baseline and procedural confounding variables.

8. OPERATOR'S INSTRUCTIONS

9.1 Inspection Prior to Use

- Before opening the product, carefully inspect the sterile inner pouch. Do not use if the integrity of the sterile inner pouch has been compromised.
- Do not use after the "Use by" date.
- Peel open the inner pouch and pass or drop the product into the sterile field using an aseptic technique.
- Carefully remove the system from the coiled hoop and inspect it for bends, kinks, and other damage. Do not use the product if any damage is noted.
- Carefully remove the protective stent sheath and mandrel together, visually inspect the stent to ensure that it has not been damaged or displaced from its original position (between proximal and distal marker bands) on the balloon.

Note: do not manipulate, touch, or handle the stent with your fingers, which may cause contamination, or stent dislodgement from the delivery balloon.

9.2 Materials Required

Quantity	Material
1	Guiding catheter [≥ 6 F (1.8 mm, 0.070 inch) inner diameter]
2-3	20 cc syringe
1,000 u/500 cc	Sterile Heparinized Normal Saline (HepNS)
1	Guidewire [≤ 0.014 inch (0.36 mm) diameter]
1	Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
N/A	Contrast medium diluted 1:1 with heparinized normal saline
1	Inflation device
1	Stopcock (3-way minimum)

1

Torque device

9.3 Preparation

9.3.1 Packaging Removal

- ♦ Carefully remove the delivery system from its protective tubing for preparation of the delivery system. Do not bend or kink the hypo tube during removal.
- ♦ Remove the product mandrel and protective stent sheath together by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally.

9.3.2 Guide Wire Lumen Flush

Flush the guide wire lumen with HepNS using the flushing tool. Insert the flushing tool into the tip of the catheter and flush until fluid exits the guide wire exit notch.

Note: Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent on the balloon.

9.3.3 Delivery System Preparation

- 1) Fill a 20 cc syringe with 5 cc of contrast/heparinized normal saline mixture (1:1).
- 2) Attach the syringe to the stopcock; attach the stopcock to the inflation Luer port of delivery system.
- 3) With the tip down, orient the delivery system vertically.
- 4) Open the stopcock to delivery system; pull negative for 30 seconds; release to neutral for contrast fill.
- 5) Close the stopcock to the delivery system; purge the syringe of all air.
- 6) Repeat steps 2 through 4 until all air is expelled. If bubbles persist, do not use the product.
- 7) Attach a prepared inflation device (only partially filled with contrast media) to stopcock.
- 8) Open the stopcock to the delivery system.
- 9) Leave on neutral

Note: While introducing the delivery system into the vessel, **do not induce negative pressure on the delivery system.** This may cause dislodgement of the stent from the balloon.

9.4 Delivery Procedure

- 1) Prepare the vascular access site according to standard practice.
- 2) Pre-dilate the lesion with a PTCA catheter of appropriate length and diameter for the vessel/lesion to be treated.
- 3) Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as widely as possible.
- 4) Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion.
- 5) Advance the stent delivery system over the guide wire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent

across the lesion; perform angiography to confirm the position of the stent. If the position of the stent is not optimal, it should be carefully repositioned or removed (see Precautions – 6.9 Stent/System Removal Precautions). The inside edges of two balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion segment of the vessel.

- 6) Sufficiently tighten the rotating hemostatic valve. The stent is now ready to be deployed.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, do not force passage.

Remove the entire system as a single unit. See Section 6.9 – Precautions, Stent System Removal for specific delivery system removal instructions.

9.5 Deployment Procedure

- 1) Prior to stent expansion, reconfirm the correct position of the stent relative to the target lesion via the radiopaque balloon markers.
- 2) Open the stopcock to delivery system and apply negative pressure to purge the balloon of air with inflation device.
- 3) Turn the stopcock on the delivery system to the off position and purge the inflation device of air. Close the side port of the stopcock.
- 4) Under fluoroscopic visualization, inflate the balloon to deploy the stent. Deploy the stent slowly by pressurizing the delivery system until stent is completely expanded. Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall. **Do not exceed the labeled rated burst pressure (RBP) of 16 atm (1621 KPa) for the 2.25 mm – 3.5 mm stents and 14 atm (1419 kPa) for the 4.0 mm stents sizes.**
- 5) Fully cover the entire lesion and balloon treated area (including dissections) with the WALTZ™ stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
- 6) Confirm that the stent is adequately expanded by angiographic injection through the guiding catheter. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).

9.6 Removal Procedure

- 1) Deflate the balloon by pulling negative pressure on the inflation device. Confirm complete balloon deflation before attempting to remove the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to the guiding catheter position.

- 2) Fully open the rotating hemostatic valve.
- 3) While maintaining the guide wire position and negative pressure on the inflation device, withdraw the delivery system.

Note: Should any resistance be felt at any time during either lesion access or removal of the delivery system post stent implantation, the entire system should be removed as a single unit. See Section 6.9– Precautions, Stent System Removal for specific delivery system removal instructions.

- 4) Tighten the rotating hemostatic valve.
- 5) Repeat angiography to assess the stented area. If post-dilatation is necessary, ensure that the final stent diameter matches the reference vessel diameter. **Assure that the stent is not under-dilated.**

9.7 Post-Deployment Dilatation of Stent Segments

All efforts should be taken to assure that the stent is not under-dilated. If the deployed stent size is still inadequate with respect to vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guidewire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the following limits.

Nominal Stent Diameter	Dilatation Limit
<u>2.25 mm</u> to 2.5 mm	2.76 mm
2.75 mm to 3.0 mm	3.50 mm
3.5 mm to 4.0 mm	4.25 mm

Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

9.8 In Vitro Compliance Information

Table 9.1: Typical WALTZ™ Stent Compliance

Pressure			Stent					
			Inner Diameter (mm)					
atm	kPa		<u>2.25</u>	2.5	2.75	3.0	3.5	4.0
6	608		<u>2.14</u>	2.10	2.42	2.65	3.08	3.55
8	811		<u>2.21</u>	2.35	2.60	2.88	3.30	3.75

9	912	Nominal	2.25	2.50	2.75	3.00	3.50	4.00
10	1013		2.29	2.58	2.82	3.03	3.58	4.08
12	1216		2.36	2.65	2.93	3.10	3.67	4.18
14	1419		2.44	2.72	3.00	3.18	3.73	4.25**
16	1621	RBP*	2.52	2.76	3.05	3.27	3.80	4.35
18	1824		2.59	2.83	3.18	3.35	3.90	4.45

***Do not exceed the rated burst pressure (RBP)**

****Rated Burst Pressure (RBP): 14 atm (1419 kPa) for 4.0 mm diameters**

Note: The Stent Diameter values listed are actual average stent diameters at the specific balloon inflation pressures obtained during in vitro testing at 37°C and **do not take into account lesion resistance. The stent sizing should be confirmed angiographically.**

10 PACKAGING

The cobalt-chromium stent on the WALTZTM delivery system (hereinafter referred to as the WALTZTM stent system) consists of four layers (from inside to outside): rolled tube, inner bag (or inner bag, sterile barrier system), outer box (protective packaging) and transport box (transport packaging).

Layer 1: Rolled Tube

Each WALTZTM stent system has a protective sheath that covers the stent / balloon and is inserted into a welded coiled tubing.

Layer 2: Inner Pack

The rolled tube is placed in an internal heat sealable bag. The inner bag is a sterile barrier system that prevents microbiological contamination of the product. The inner bags are placed in a corrugated cardboard box for sterilization with ethylene oxide.

Layer 3: External carton

The inner bag is placed in a labeled cardboard box with instructions for use (PPI). The box is designed to protect packaged products from damage from the moment of their packaging to the time of use.

Layer 4: Shipping Box

Individual items are placed in suitable boxes designed to protect packaged products from damage during transport.

11 TRANSPORTATION CONDITIONS

The product should be transported at a temperature of 10 ° C to 30 ° C. Keep away from moisture. Keep away from sunlight.

12 WORKING TEMPERATURE

Use at a temperature of + 32 ° C to + 42 ° C (human body, blood flow)

13 DISPOSAL INSTRUCTIONS

On the territory of the Russian Federation, the product should be disposed of in accordance with SanPiN 2.1.7.2790-10 for class B waste (epidemiologically hazardous waste). After hardware-based disinfection methods using physical methods and changing the appearance of the waste, which

excludes the possibility of their reuse, class B waste can be accumulated, temporarily stored, transported, disposed of and disposed of together with class waste (epidemiologically safe waste similar in composition to solid waste). Packaging of Class B medical disinfected waste should be labeled to indicate that the waste has been decontaminated.

14 SHELF LIFE

Under store requirements conditions, the "Use by" of WALTZ™ CoCr coronary stent system is 2 years.

15 WARRANTY AND LIMITATIONS

Shanghai MicroPort Medical (Group) Co., Ltd. warrants that all the product shall be free of defects in materials and workmanship upon receipt. Shanghai MicroPort Medical (Group) Co., Ltd. will not be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from use of its product.


Shanghai MicroPort Medical (Group) Co., Ltd. indicates definitely that this device is intended for one time use only and makes no warranty, representation or condition of any kind, whether expressed or implied (including any warranty of merchantability, suitability or fitness for a particular purpose) respecting the reuse of this catheter. Moreover, Shanghai MicroPort Medical (Group) Co., Ltd. assumes no responsibility or liability for incidental or consequential damages which may result from such reuse. Ask us directly any other questions related to the products.

16 COMPLIANCE

The medical device complies with the provisions of Directive 93/42 / EEC (Appendix II) and the following standards: EN ISO 13485, EN ISO 14971, EN ISO 14630, EN ISO 25539-2, EN ISO 10555-4, series standards EN ISO 10993, EN ISO 11135, EN ISO 11737-1, EN ISO 11607-1, EN ISO 15223-1.

17 GRAPHIC SYMBOLS FOR LABELS FOR MEDICAL PRODUCTS

The primary packaging label and the outer box label are provided with three (3) stickers — adhesive metering labels for attaching to the medical history of the patient receiving the implant. Stickers contain the following information:

- manufacturer identification; MicroPort manufacturer logo;
- product information: trademark "WALTZ", stent size (diameter x length) mm;
- the symbol "series number";
- the symbol "catalog number"
- the symbol "expiration date"

Stickers are part of the labels (perforation is applied on the edge of each sticker - punching holes along the tear line, which makes it easy to separate the sticker from the label and paste it into the patient's medical history).

17.1 Explanation of symbols on the marking

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18 MANUFACTURER

Shanghai MicroPort Medical (Group) Co., Ltd.

Add: 501 Newton Rd., ZJ Hi-Tech Park, Shanghai 201203, P.R.China

Tel: (86) (21) 38954600

Fax: (86) (21) 50801305

Email: customerservice@microport.com

Website: www.microport.com

19 OFFICIAL REPRESENTATIVE IN THE RUSSIAN FEDERATION

ETS-PRO Limited Liability Company (ETS-PRO LLC)


Address: Russia, 105523, Moscow, Schelkovo highway, 100, bldg. 1, of. 5051

Phone: (495) 781-63-95

E-mail: info@ets-pro.com

Registration certificate No. _____ from XX.XX.20XX

 I affirm
Dr. Zhao Hua Chang
Chairman & CEO
Shanghai MicroPort Medical (Group) Co., Ltd

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