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案號: E 0929	日期: AUG. 05 2020
Case No.:	Date:
本文件之簽名或蓋章, 在臺灣臺北地方法院所屬民間公證人陳建源事務所認證。 公證人 陳 建 源	
Attested at the Cheng, Chien-Yuan Notary Public Office, Taiwan Taipei District Court, R.O.C., that the signature(s)/seal(s) in this document is/are authentic.	
Notary Public	✓ ✓ 4
Cheng, Chien-Yuan	

9F-2, No. 52, Sec. 2, Chang An E. Rd, Taipei, Taiwan
Tel: (02)25074587

APPROVE

Sarah/ Lu

(name)

Vice President

(position)

(signature)

« 04 » 08 2020
«day» month (numerals)

Stamp

Besmed

Instructions for use of the medical device: BREATHING FILTER in different variants

Manufacturer: Besmed Health Business Corp.

Address: No.

5, Lane 116, Wu-Kong 2nd Rd, Wu-Ku District, New Taipei City 248, Taiwan

Telephone: 886-2-2290-3959

Fax: 886-2-2299-9076

E-mail address: info@besmed.com

1. Medical device name

BREATHING FILTER (hereinafter, filter or device)

Design variants:

Breathing filter with Luer-lock port FL-70311

Breathing filter with Luer-lock port FL-70312

Breathing filter with Luer-lock port FL-70313

Breathing filter with Luer-lock port FL-70314

Breathing filter with Luer-lock port FL-70114

Breathing filter with Luer-lock port FL-70211

Breathing filter FL-70320

Breathing filter with oxygen port FL-70125

Breathing filter with oxygen port FL-70126

Breathing filter with oxygen port and sanitation port FL-70128

Breathing filter with oxygen port and sanitation port FL-70129

Breathing filter with oxygen port and sanitation port FL-70130

Breathing filter FL-70124

2. Purpose of the medical device

Filters are medical devices designed to be placed in the breathing circuit proximally relative to the patient to remove and retain microbes while capturing the heat and moisture generated by the patient's exhalation for use in heating and humidifying therapeutic gases inhaled by the patient.

3. Scope of the medical device

Surgery. Designed for use in intensive care units.

4. Medical device usage conditions

at temperatures from 0 to 40 °C

air humidity from 20 to 80%

5. Indications for use

for removal and retention of microbes when capturing the heat and moisture generated by the patient's exhalation

6. Contraindications

- This device should not be used in conjunction with heated humidifiers or sprayers.

7. Warnings

All filter components are sterile. Do not use if the package is damaged or opened. Sterilized with ethylene oxide. For single use. Re-sterilization or re-use are prohibited. Does not contain medicinal substances. Biologically safe for humans. Store under the conditions described in the instructions, failure to comply with these conditions may cause damage to the filters. It is prohibited to use in conjunction with heated humidifiers or sprayers. The filter is used only by physician's prescription and only under the supervision of a medical professional.

Store under conditions described in the operating instructions, otherwise the filter will not be suitable for use

Excessive discharge/bleeding may cause the filter to become clogged and make it difficult to breathe; in this case, the filter should be replaced immediately. Make sure that there are no obstacles to the passage of air inside the filter.

Make sure that the aspiration/inspection holes are properly closed during use.

To reduce the risk of damage and contamination by small objects, remove the device from its package only immediately before use.

8. Safety instructions

- Regularly check for increased flow resistance and replace the filter if necessary.
- Avoid contact with acetone or similar solvents that may disrupt its structure and cause the catheter to leak or damage it.
- Maximum filter operation time: 24 hours.
- Do not restore, re-use, soak, rinse, wash or sterilize the medical device, and do not use any disinfectants for its treatment (especially avoid contact with phenol and alcohols).
- The device may be used only with devices whose connectors comply with ISO standards.
- Replace the filter for patients at least once a day to reduce the risk of obstruction and the appearance of a colony of bacteria or fungus in the device.

POSSIBLE ADVERSE REACTIONS

Filter use may cause the following adverse reactions. Complications are listed in an order that does not reflect their frequency or severity: the consequence of blocking the filter with sputum and complications due to a violation of the filter's conductivity – breathing difficulty, dyspnoea, hypercapnia and hypoxia.

10. Medical device manufacturer

Manufacturer: Besmed Health Business Corp.

Address: No.

5, Lane 116, Wu-Kong 2nd Rd, Wu-Ku District, New Taipei City 248, Taiwan

Telephone: 886-2-2290-3959

Fax: 886-2-2299-9076


E-mail address: info@besmed.com

11. Composition of the medical device, the list of design variants

Breathing filter with oxygen port (FL-70126): consists of a filter body with an oxygen port, a PP filter and a sanitation port


Breathing filter with oxygen port (FL-70125): consists of a filter body with oxygen port and a PP filter

Breathing filter with oxygen port and sanitation port (FL-70128, FL-70129, FL-70130): consists of a filter body with oxygen port, PP filter and sanitation port (see table below)


Breathing filter with Luer-lock port FL-70114
Consists of a filter body, a PP filter and a Luer-lock port

Breathing filter with Luer-lock port FL-70314
Consists of a filter body, a PP filter and a Luer-lock port

Breathing filter with Luer-lock port FL-70313 and FL-70312
Consists of a filter body, a PP filter and a Luer-lock port

Breathing filter with Luer-lock port FL-70311
Consists of a filter body, a PP filter and a Luer-lock port

Breathing filter FL-70211
Consists of a filter body and a PP filter

22M/15F

22F

Breathing filter FL-70130
Consists of a filter body, a PP filter and a sanitation port

Breathing filter FL-70320
Consists of a filter body, a PP filter

Breathing filter FL-70124
Consists of a filter body, a PP filter and a tube

The delivery package includes:

1. Breathing filter in the design variant - 1 piece
2. Instructions for use - 1 piece

List of accessories

Not provided

12. Technical characteristics

Breathing filter design variant	FL-70125	FL-70126	FL-70211	FL-70114	FL-70128	FL-70129	FL-70130
Range of breathing volume (BV) (ml)	250-1000	15-50	32	250-1000	allows enrichment of inhaled air with oxygen with a concentration up to 50%	Allows enrichment of inhaled air with oxygen with a concentration up to 50%	Allows enrichment of inhaled air with oxygen with a concentration up to 50%
Dead space (ml)	12.5	23	-	58	10	10	10
Moisture output	9 mg of water/L at BV 500 ml	-	29 mg of water/L at BV 500 ml	33 mg of water/L at BV 500 ml	9,5 mg of water/L BV 500 ml	9,5 mg of water/L BV 500 ml	9,5 mg of water/L BV 500 ml
Flow resistance (cm H ₂ O)	30 L/min 0,62 cm H ₂ O	-	30 L/min 0,75 cm H ₂ O 60 L/min 1,48 cm H ₂ O	30 L/min 0,8 cm H ₂ O 60 L/min 1,5 cm H ₂ O	30 L/min 0,25 cm H ₂ O	30 L/min 0,25 cm H ₂ O	30 L/min 0,25 cm H ₂ O
Weight (g)	6.8	27.60	23.60	29.51	4.0	4.0	6.3
Connectors patient side (mm) apparatus side	15F/22 M	22F/15M 22M/15F	22F/15M 22M/15F	22F/15M 22M/15F	-	-	-


(mm)							
Bacterial filtration efficiency (BFE)	>99.9992%	>99.9992%	>99.9992%	>99.9992%	>99.9992%	>99.9992%	>99.9992%
Viral filtration efficiency (VFE)	>99.9970%	>99.9970%	>99.9970%	>99.9970%	>99.9970%	>99.9970%	>99.9970%

Model	FL-70311	FL-70312	FL-70313	FL-70314	FL-70320	FL-70124
Range of breathing volume (ml)	50-250	250-1000	250-1500	250-1000	250-1000	more than 1500
Dead space (ml)	19	45	60	71	35	90
Moisture output	29 mg/L at BV 500 ml	28 mg/L at BV 500 ml	29 mg/L at BV 500 ml	28 mg/L at BV 500 ml	29 mg/L at BV 500 ml	29 mg/L at BV 500 ml
Flow resistance (cm H ₂ O)	30 L/min 2,6 cm H ₂ O 60 L/min 5,3 cm H ₂ O	30 L/min 1,6 cm H ₂ O 60 L/min 2,6 cm H ₂ O	30 L/min 1,1 cm H ₂ O 60 L/min 2,5 cm H ₂ O	30 L/min 1,6 cm H ₂ O 60 L/min 2,6 cm H ₂ O	7 L/2.5 cm H ₂ O	0.5 L/s 0.5 cm H ₂ O 1.05 L/s 1.4 cm H ₂ O
Weight (g)	11.30	27.90	34.50	38.50	31.0	35.0
Bacterial filtration efficiency (BFE)	>99,9992 %	>99,9992 %	>99,9992 %	>99,9992 %	>99,9992%	>99,9992%
Viral filtration efficiency (VFE)	>99,9970 %	>99,9970 %	>99,9970 %	>99,9970 %	>99,9970%	>99,9970%
Size (height x width)	7.16 cm x 6.00 cm	7.27 cm x 5.91 cm	7.5 cm x 6.23 cm	7.19 cm x 5.26 cm	-	-
HEPA class	-	13	13	13	-	-
Connectors patient side (mm) apparatus side (mm)	15M 22M/15F	22F/15M 22M/15F	22F/15M 22M/15F	22F/15M 22M/15F	22F/15M 22M/15F	22F/15F

Model	FL-70128	FL-70129	FL-70130
Sanitation port	5.0 mm	3.0 mm	2.0 mm
Color	White	Blue	Blue


Conical connectors should meet the following requirements:

Connector size	A	B	C	D	E	Taper F	Tolerance per unit length
15	15.525±0.005	15.165±0.005	14.5±0.005	4.3±0.005	2.2±0.005	1:40	0.025±0.00
22	22.425±0.005	21.94±0.005	19.5±0.005	5.2±0.005	2.2±0.005	1:40	0.025±0.00

Dimensions of conical connectors

Male parts of the 22 mm conical connections, except for those intended for attachment to the face mask, should have a groove in accordance with the figure below).

The breathing volume was set at the rate of 8-9 ml/kg of ideal body weight. The maximum period of filter use is limited to 48 hours.

The fitting on the patient side is a covered conical connector 22 mm in size and, coaxially, a female conical connector 15 mm in size;

both connectors should comply with ISO 5356-1.

The sanitation port should not get disconnected from the filter body. Gas leakage from the sanitation port, brought to a temperature of 20 °C and a pressure of 101.3 kPa, should not exceed 5 ml/min.

13. Description of materials contained in key functional elements / blocks / parts of the medical device, including materials in contact with the human body

Name of the breathing filter component	Material, grade	Supplier
Filter body	Polyvinylchloride EV08EFP02 green pigment Pigment No. Green 7 :CAS NO.1328-53-6 (Chloronated copper phthalocyanine)	Green Field Chemtech co. Ltd.
Filter	Polypropylene PP 5090T	Formosa Plastics Corporation
Sanitation port	TEP Dryflex® blue pigment Pigment No. Blue 15:0:CAS No,147-14-8 (Copper phthalocyanine)	Hotai Industrial co. ltd
Luer-lock port	Polypropylene PP 5090T blue pigment Pigment No. Blue 15:0:CAS No,147-14-8 (Copper phthalocyanine) green pigment Pigment No. Green 7 :CAS NO.1328-53-6 (Chloronated copper phthalocyanine)	Formosa Plastics Corporation
Tube	EVA EV 302	ASIA polymer Corporation
Package	HDPE polyethylene 8050	Formosa Plastics Corporation

14. Preparation for work and operational procedure

14.1 Preparation for work

ATTENTION: Strictly follow aseptic technique and general precautions and procedures established by the protocols of the medical institution.

The procedure is performed according to the guidelines for medical professionals approved in a medical institution.

14.2 Operational procedure

1. Choose a model of the breathing filter or a filter with a heat and humidity exchange that is appropriate for the patient's age and clinical situation.
2. It is recommended to install breathing filters with a heat and humidity exchanger between the end of the endotracheal (tracheostomy tube, laryngeal mask, etc.) and the Y-shaped element of the breathing circuit. Breathing filters with heat and humidity exchange are not recommended for use with active humidifier or aerosol systems.
3. Breathing filters without a heat-and-humidity exchanger may be installed between the end of the endotracheal (tracheostomy tube, laryngeal mask, etc.) and the Y-shaped element of the breathing circuit, or along the inspiration and expiration line. Respiratory filters without heat-and-humidity exchanger are not recommended for use with active humidifier or aerosol systems.
4. Open the package and remove the device.
5. Connect the device to the breathing circuit and the patient's endotracheal tube (tracheostomy tube, laryngeal mask etc.), make sure that the breathing circuit is tight (there is a plug on the Luer-lock port).
6. Set the necessary ventilation parameters taking into account the "dead space".
7. When using a breathing filter or a breathing filter with a heat transfer unit equipped with a Luer-lock connector for gas analysis/capnography connect the analyzer line to the filter port, and make sure that the connection is tight.
8. Replace the device every 24 hours of ventilation or more often as contamination occurs.
9. It is not allowed to use one device in several patients.
10. Disinfection, sterilization and reuse of the device are not provided, dispose of the product after use.

Note

1. Device is used under the physician's supervision.
2. If damage of the package is detected, the use of the filter is not recommended.
3. Obstruction of the device by the discharge from the patient's airways, which disrupts ventilation and requires immediate replacement of the filter, is possible.
4. It is necessary to make sure that the breathing circuit is tight, check the connection of the device with the breathing circuit and the patient's endotracheal tube (tracheostomy tube, laryngeal mask, etc.), check the connection of the gas analysis/capnography line with the filter port, analyze the ventilation monitoring data for leaks.

15. Maintenance

- 15.1 The controlled conditions and environment are constantly monitored for bioburden using the sedimentation plate method. Bioburden control method, the frequency of checks, and the limit value of colony-forming units (CFU) are provided in GTP/QC/14.
- 15.2 Periodical monitoring of the bioburden on equipment, structures, and personnel clothing, is also conducted. Bioburden control method, the frequency of checks, and the limit value of colony-forming units (CFU) are provided in GTP/QC/14.

16. Cleaning and disinfection

Not provided.

17. Safe use conditions

17.1 The device is manufactured under controlled conditions: a class 7 clean room in accordance with ISO 14644-1. The room is equipped with a high-efficiency particulate air filter (HEPA), humidity and temperature are controlled. The room environment is periodically validated to check the effectiveness of the HEPA filter. The room conditions correspond to WI/CVC/02.

17.2 Possible complications of excessive humidification:

- Increased risk of nosocomial infection
- Increased mucosal discharge
- Increased need for suction procedures
- The narrowing and occlusion of the tracheal tube
- Water condensation may block the air duct, causing obturation

Typical characteristics of an unheated medical gas mixture:

Temp. 20°C

Abs. hum. 1 mg H₂O/L

Possible complications of insufficient humidification:

- Restriction and occlusion of the tracheal tube lumen
- Mucosal and ciliary dysfunctions
- Atelectasis
- Increased possibility of postoperative pulmonary complications
- The change in lung mechanics that leads to hypoxemia

Typical characteristics of an isothermal saturation zone:

Temp. 37°C

Abs. hum. 44 mg H₂O/L

17.3 The device is used according to its instructions, failure to follow the instructions may result in damage to the filters and termination of the manufacturer's warranty.

The device does not contain medicinal products.

18. Environmental protection requirements, including requirements for safe destruction and disposal. Environmental protection.

18.1 Products during storage, transportation and operation do not emit toxic substances into the environment and do not (in direct contact) have a harmful effect on the human body and the environment. When working with the devices, no special measures to protect the environment from harmful effects are required.

18.2 All filter manufacturing operations are carried out in rooms equipped with local extract and input ventilation.

Open flame in these rooms is prohibited.

18.3 Environmental protection in accordance with the applicable sanitary and epidemiological rules and regulations of the manufacturer.

18.4 The controlled conditions and environment are constantly monitored for bioburden using the sedimentation plate method. Bioburden control method, the frequency of checks, and the limit value of colony-forming units (CFU) are provided in section

18.5 Periodical monitoring of the bioburden on equipment, structures, and personnel clothing, is also conducted. Bioburden control method, the frequency of checks, and the limit value of colony-forming units (CFU) are provided in section

18.6 Disposed of as class b waste – epidemiologically hazardous waste.

Unused devices are disposed of in accordance with the standards adopted in Medical Practice, as well as the relevant local, state and federal legislative acts and instructions, as class A waste SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for medical waste handling".

19. Service life data

The finished product remains in compliance with technical and functional requirements for at least 5 years from the date of manufacture. Service life – 24 h.

20. Transportation and storage requirements

Transportation of products is carried out by all types of covered vehicles in accordance with the rules of transportation of goods in force for this type of transport, subject to temperature conditions.

Transportation at temperatures from 0 to 40 °C.

Store in a dry place, protected from direct sunlight, at temperatures from 0 to 40 °C.

21. Manufacturer's warranty

Twelve month warranty. The manufacturer guarantees absence of workmanship or material defects if proper storage and use conditions are met. This warranty is void in the event of improper storage or use, as well as in the event of factors related to the patient, his/her diagnosis, treatment, surgery and other aspects that are beyond manufacturer's control and have a direct impact on the manufacturer's products and the results of its use.


22. Medical device labeling

22.1 Description of the MD labeling

The labeling includes the following particulars:

- manufacturer's trademark
- manufacturer's name and address
- name and designation of the product type (series)
- product number according to the manufacturer's numbering system
- other data depending on product requirements (see symbols used for labeling).

Example of labeling for RF:

O-FL-70000-2 Rev 0 201908


Labeling of the outer cardboard box and transport cardboard box is performed with additional symbols regarding correct transportation and storage of devices, such as:

- Temperature limitations


- Protect from moisture

- Protect from direct sunlight

Layout of secondary and transport packaging of the device:

22.2 Description of the symbols, printed on the MD label

	Do not re-use
	Caution! Please refer to the accompanying documentation!
	Use before
	Batch code
	Reference number
	Sterile, sterilized with ethylene oxide.
	Do not re-sterilize
	Date of manufacture
	Manufacturer, address
	European conformity

22.3. Device packaging process

Devices are packed in air-bubble wrap.

Ethylene oxide sterilization procedure according to EN ISO 11607-1.

Information about the sterilization procedure:

1. Exposure time: 3 hours / half-cycle, 6 hours / full cycle.
2. Temperature: $65 \pm 5^\circ\text{C}$
3. Humidity: RH 55~90%

4. Pressure: 500 mm Hg $\pm$ 100 mm Hg

23. List of applicable national and international standards to ensure the safety, efficiency and quality of the medical device

- ☐ EN ISO 14971: 2012 Medical devices. Application of risk management to medical devices
- ☐ EN ISO 13485: 2012 Medical devices. Quality Management Systems. Requirements for regulatory purposes
- ☐ EN ISO 9360-1: 2009 Anaesthetic and respiratory equipment. Heat and moisture exchangers (HMEs) for humidifying respired gases in humans. Part 1. HMEs for use with minimum tidal volumes of 250 ml
- ☐ EN ISO 9360-2: 2009 Anaesthetic and respiratory equipment. Heat and moisture exchangers (HMEs) for humidifying respired gases in humans. Part 2: HMEs for use with tracheostomized patients having minimum tidal volumes of 250 ml
- ☐ EN ISO 1041: 2008 Information supplied by the manufacturer of medical devices
- ☐ EN ISO 5356-1: 2015 Anaesthetic and respiratory equipment. Conical connectors Part 1. Cones and sockets
- ☐ EN ISO 11135-1: Sterilization of health care products. Ethylene oxide. Part 1. Requirements for development, validation and routine control of a sterilization process for medical devices
- ☐ EN ISO 11607-1: 2014 Packaging for terminally sterilized medical devices. Part 1. Requirements for materials, sterile barrier systems and packaging systems

Authorized representative of the manufacturer in the RF (name, address, telephone, mail)

Authorized representative of the manufacturer:

Full or abbreviated company name (in accordance with the Constituent documents)	Odyssey Medical Limited liability company (Odyssey Medical LLC)
Legal address (in accordance with the Constituent documents)	Apt. 4, 36, Zoyi i Aleksandra Kosmodemyanskikh Str., Moscow, 125130
Telephone (at the actual address)	8-925-145-61-34
E-mail	info@odyssey-medical.ru

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