Instructions for use

Intra-articular gel OArt+

ALBOMED GmbH,
Hildebrandstrasse 11,
90592 Schwarzenbruck, Germany
Tel.: +49 (0)9183-956982-0
Fax: +49 (0)9183-956982-9
E-Mail: info@albomed.eu

Responsible person Peter Halfwassen
Position CEO, ALBOMED GmbH

ALBOMED® GmbH
Hildebrandstrasse 11
D-90592 Schwarzenbruck
Tel. +49 (0) 9183 - 95 69 82-0
Fax +49 (0) 9183 - 95 69 82-9

Mma , 1. April 2020
(place) (date)

Alb Half
(signature)

INSTRUCTIONS FOR USE For Russian Federation

Trade name of medical device: Intra-articular gel OArt+

Purpose: For viscosupplementation of intra-articular synovial fluid in damaged joints and for improving of viscoelastic properties.

Field of application:

The gel is used in specialized medical and medical-prophylactic institutions.
Only for use by physicians who have been trained to inject implantable gels.

Information about active ingredient (sodium hyaluronate solution):

Intra-articular gel OArt+ is a clear sterile solution of sodium hyaluronate with different molecular weights and concentration in phosphate buffer solution depending on dosage.

The active ingredient of the intra-articular gel OArt+ (OArt+ 10, OArt+ 44 and OArt+ 120) is sodium hyaluronate. Sodium hyaluronate is obtained by bacterial fermentation of the natural bacteria *Streptococcus equi*, no animal tissue/cells are used. Intra-articular gel OArt+ is completely absorbed, serves as an analogue (implant) of synovial fluid produced in the human body.

The gel gradually dissolves over time. The time of complete biodegradation and excretion from the body is 6 - 12 months and depends on the individual characteristics of the patient.

Indication for application of medical device:

- The gel is indicated for supplementation of synovial fluid in joints and to improve viscoelastic properties of synovial fluid in case of OA.

Contra-Indications:

The product should not be used for patients

- with a known hypersensitivity reaction to any of the components of the product;
- with septic arthritis;
- with infections of the skin or dermatological diseases at the injection side;
- given coagulation inhibitors.

Precautions for use:

- patient should be informed about particularities of gel application, contra-indications and possible side effects;
- before use it's necessary to check that the packaging is not damaged and the expiry date of the product;
- application is not recommended during pregnancy and lactation;
- application is not recommended in patients younger 18 years.

The product must not be used if:

- it's expired;
- the packaging is damaged.

Notice! During primary cleaning of site for injection do not use antiseptics containing quaternary ammonium compounds:

It's recommended to avoid contact with products containing quaternary ammonium compounds.

Interaction with other medical devices:

Simultaneous application of other injectable gels is not recommended since no clinical data related to compatibility of the product with analogue medical devices are available.

Biocompatibility:

Toxicological data show that the product is well tolerated and has low allergenic potential.





Guidelines:

- the application of the medical device should be performed by physicians who have been trained to inject implantable gels;
- for a reliable assessment of the feasibility and maximum efficiency of the application of an intradermal implant, a preliminary assessment of the patient's condition is necessary;
- the application should be performed in accordance with rules of aseptic and antiseptic conditions;
- the injection should be performed without efforts;
- for injection of the gel following size of syringe needle for single use with Luer Lock adaptor is recommended:
 - Needle with size 21 (Ø 0.8 mm) is suitable for injection of Intra-articular gel OArt+ 10 and OArt+ 44;
 - Needle with size 19 (Ø 1.1 mm) is suitable for injection of Intra-articular gel OArt+ 120.
- The product should be injected into the articulation. Intra-vascular or intra-tissular injection should be avoided.
- The local application of ice for 5 - 10 minutes may reduce side effects after the gel injection.
- Simultaneous intake of analgesics or anti-inflammatory drugs by oral route may be useful during the first days after treatment.

Method of administration:

- Remove air bubbles from the syringe before injection.
- The skin is disinfected with local antiseptic.
- The product is injected in the joint, once a week, 1 - 5 times. The beneficial effects of the treatments are maintained during at least 6 months. In the case of intra articular swelling, fluid should be removed by aspiration before the hyaluronic acid injection.
- Depending on the size of the joint: 1, 3 or 5 injections of 1 ml, 2 ml or 4.8 ml may be administered intra-articular.
 - The prefilled syringe of 1 ml is the volume adapted for small joints.
 - The prefilled syringe of 2 ml is the volume adapted for the knee, hip, shoulder and ankle.
 - The prefilled syringe of 4.8 ml is the volume adapted for the knee only.
- The recommended treatment regimen with the mono dose of 4.8 ml for the knee is one injection, which may be repeated 6 months later based on the clinical experience or the pain relief obtained by the patient.
- For additional injections the consultation of the physician is needed.

Guidelines for needle setting:

Side effects:

Before treatment the patients should be clearly informed about possible side effects.

In rare cases immediately or during some period after injection following local side effects may occur:

- inflammatory reactions (redness, swelling) which may be associated with pruritus, pain if pressing. These reactions occur after injection and the duration is till 1 week;
- mild bleeding, especially in the presence of bleeding disorders or taking drugs that reduce blood viscosity;
- hardening or compaction at the injection site;
- skin discoloration at the injection site;
- skin pigmentation at the site of implant placement;
- hypersensitivity to the components of the gel, especially sodium hyaluronate.

The patient should be clearly informed to provide their doctor the information of any side effects lasting more than one week. In this case, the doctor may prescribe the appropriate treatment to the patient.

The information about all adverse and side effects related to gel injection should be provided to the doctor.

Strongly prohibited!

- Intra-vascular injection and also off-label use;
- Dilution of the gel and injection with other substance in one syringe;
- Usage if packaging is damaged;
- Re-sterilization;
- Use after expiry date.

Storage and transporting conditions:

The device in packaging from the manufacturer should be stored at temperature + 2 °C to + 25 °C and relative humidity of 40 to 70 %, in ventilated warehouse, protected from light and atmosphere effects, at distance not less of 1 m from heater in dry place. Avoid from organic solvents, ionizing radiation or UV since it can lead to destruction of plastic. Avoid impact effects.

Keep out of the reach of children!

Attention! Do not freeze!

The device is transported by all types of covered vehicles, except for unheated compartments of aircraft, in accordance with the rules approved in the established manner for each type of

transport. Transportation by sea must be carried out in accordance with the "Rules for the safety of sea transportation of general cargo."

Type of shipment - in containers or by other means as agreed with the consumer, at a temperature of + 2 ° C to + 25 ° C.

Shelf-life:

42 (forty-two) months from date of manufacturing (finish sterilization).

Manufacturer guarantees the stability of gel components until expiry date in case if transporting, storage and application conditions were kept in accordance mentioned in the Instructions for use.

Packaging (set of medical device):

1.0 ml or 2.0 ml of gel contained in glass syringe with nominal volume 2.25 ml and 4.8 ml contained in glass syringe with nominal volume 5.0 ml in blister and consumer packaging – carton box with instructions for use. Graphic symbols applied on consumer packaging are common for international circulation and allow explanation of necessary requirements for medical device by visual demonstration – storage conditions, presentation, etc.

Chemical composition of Intra-articular gel OArt+:

Ingredient	Specification acc. to current edition of European Pharmacopoeia	Unit formula for 1 ml	Product and filling volume syringe
Sodium hyaluronate	Eur.Ph. 1472 current edition	10.0 mg (1.0 %)	OArt+ 10 (1.0 ml)
		22.0 mg (2.2 %)	OArt+ 44 (2.0 ml)
		25.0 mg (2.5 %)	OArt+ 120 (4.8ml)
Sodium chloride	Eur.Ph. 0193 current edition	8.50 mg	OArt+ 10 (1.0 ml) OArt+ 44 (2.0 ml)
		6.50 mg	OArt+ 120 (4.8ml)
Disodium hydrogen phosphate dihydrate	Eur.Ph. 0602 current edition	0.563 mg	
Sodium dihydrogen phosphate dihydrate	Eur.Ph. 0194 current edition	0.045 mg	
Water for injection	Eur.Ph. 0169 current edition	q. s. ad 1 ml	

Technical characteristics of intra-articular gel OArt+:

Tests	Specifications
Appearance	Clear, colorless, particle free viscoelastic gel
Osmolality	270 – 400 mOsmol/kg
pH	6.8 – 7.4
Extractable volume	> to nominal volume / syringe
Sodium hyaluronate content	9.0 – 11.0 mg/ml (OArt+ 10, 1.0 %) 19.8 – 24.2 mg/ml (OArt+ 44, 2.2 %)

	22.5 – 27.5 mg/ml (OArt+ 120, 2.5 %)
Particles	> 10 µm: 25000 per syringe > 25 µm: 5000 per syringe
Sodium Hyaluronate Molecular Weight	0.8 – 2.5 million Dalton (OArt+ 10, OArt+ 44) 1.2 – 3.2 million Dalton (OArt+ 120)
Viscosity at zero shear rate	> 10 000 mPas (OArt+ 10) > 50 000 mPas (OArt+ 44) > 200 000 mPas (OArt+ 120)
U.V. Transparency	Conform
Bacterial endotoxins	< 0.5 U.I. / ml

Gel is sterile. Sterilization is performed by moist heat (autoclaving).

Presentation:

Intra-articular gel OArt+ is presented:

1. Intra-articular gel OArt+ 10 as set:

- 1.0 % sodium hyaluronate solution, 1.0 ml in syringe.
- instruction for use -1 pcs.

2. Intra-articular gel OArt+ 44 as set:

- 2.2 % sodium hyaluronate solution, 2.0 ml in syringe.
- instruction for use -1 pcs.

2. Intra-articular gel OArt+ 120 as set:

- 2.5 % sodium hyaluronate solution, 4.8 ml in syringe.
- instruction for use -1 pcs.

Procedure for disposal:

Used syringes should be disposed as potentially infected devices related to Class B of medical wastes. Product with damaged packaging or expired products should be disposed as medical waste Class A in accordance with the Sanitary Rules and Norms 2.1.7.2790-10

Manufacturer: Albomed GmbH

Hildebrandstrasse 11, 90592 Schwarzenbruck, Germany

Tel. +49 (0) 9183-95 69 82-0, e-mail: info@albomed.eu

Manufacturing site: Tronjestrasse 14, 44139 Dortmund, Germany

Tel. +49 (0) 9183-95 69 82-0, e-mail: info@albomed.eu

Distributor: Mercure Pharma Consulting, C/o A.A.C.I. 6 Avenue Du Pre Felin PAE Les Glaisins - Annecy-le-Vieux 74000 Annecy France

Tel. + 33 673 29 31 89, +33 613 54 73 25, e-mail: vreje.iliozer@orange.fr, stephan.chianea@yahoo.fr

Authorised representative on the Territory of Russian Federation: LLC Tirupharm

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