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Mr. Bu NingBin

General Director

We appreciate you chose vein finder VeinSight, version VS400 with accessories (hereinafter equipment). Please read this Manual carefully before operation to use it correctly. Repair and troubleshooting may be performed by trained personnel only certified by the Manufacturer.

BLZ Technology (Wuhan) Co., Ltd (hereinafter "Company") is the holder of all intellectual property rights. Neither part of this Manual may be amended or translated in either way without written consent of BLZ Technology (Wuhan) Co., Ltd.

Warranty

The Company gives one year workmanship warranty and six months battery warranty. The warranty begins on the date of purchase. This warranty doesn't cover consumables.

Immediately contact shipping company if damage occurred during transportation. All supporting documents related to equipment shall include information about models and serial numbers.

Claims related to appearance, mechanical damages and content of delivery shall become null and void after acceptance of equipment!

Package shall be retained during the warranty period!

To start the warranty a customer shall provide supporting purchase document containing serial number and purchase date and shall send it by post, e-mail or fax to the Manufacturer or the Manufacturer's authorized representative:

Manufacturer:

BLZ TECHNOLOGY (WUHAN) CO., LTD.

Room 2166, 2176 Building E, Cyberport Dongxing Road, East Lake New Technology Development Zone, 430074, Wuhan, China

tel.: +86-27-65520255

web: www.veinsight.com

Manufacturer's representative in the Russian Federation:

LIMITED LIABILITY COMPANY "GENERAL MEDITECH RUSCOMPANY"

22 Delovaya street, suite 2, office 7, Nizhny Novgorod city, Nizhny Novgorod region, 603193

Phone/Fax: +7(831) 438-44-98

e-mail: gm-rus@inbox.ru

Warranty shall not cover the following:

- Mechanical damages;
- Damages resulted from natural disaster, fire or social factors;
- Damages related to power mains instability;
- Damages related to contact of inner parts with foreign objects, substances, liquids or insects;
- Damages resulted from incompliance of power, telecommunication or cable mains with state standard or similar damages;
- Damages related to use of incompliant consumables and spare parts;
- Damages of equipment or accessories related to incompliance with operating conditions;
- Incompliance of transportation and storage terms;
- Damages resulted from changes or update of equipment by third parties without the Manufacturer's or the Manufacturer's representative's consent;
- Repair of equipment during warranty period by third parties.

BLZ Technology (Wuhan) Co., Ltd and/or its official representative in the Russian Federation shall bear safety, reliability and operation liability if:

- Maintenance, repair, changes or service was performed by BLZ Technology (Wuhan) Co., Ltd or its official representative;
- Equipment is used in compliance with this Manual;
- Power complies with respective standards if a charging unit was connected to external power mains.

No	Edition	Name of the document
1	0	User Manual


Table of content

1. Medical product intended use	5
1.1 Warnings and precautions	6
1.2 Regulatory documents	7
1.3 Symbols	8
1.4 Proposed use and contraindications	9
1.5 Environment and indoor requirements	9
1.6 Battery	10
2. Functional characteristics	11
2.1 Control board	13
2.2. Focus control	14
2.3 Projected image description	14
3. Risks related to use of medical device	15
4. Specifications	16
5. Description of accessories	18
6. Maintenance and repair	19
6.1 Troubleshooting	19
7. Disinfection and cleaning	19
8. Electromagnetic compatibility	20
9. Transportation and storage	25
10. Marking and package	25
11. Disposal	27
Annex 1	28

1. Medical product intended use


Equipment is used to get projection of superficial vein pattern. It is a portable device that helps medical staff to detect respective superficial veins of any age and any physical peculiarities patients.

Vein finder VeinSight, model VS400

1. Vein finder
2. Charger
3. Network adapter
4. Calibration sheet

Accessories:

1. Mobile support
2. Fixed support

- Equipment shall not be treated as the only way to determine superficial veins and is designated for medical staff only.

- When using the equipment users shall always comply with respective valid hospital medical practice rules and procedures to prepare accurate medical conclusion.
- It shall be used when it is necessary to determine vein pattern. If used correctly, the equipment enables determination of superficial veins together with medical procedures, for example, venipuncture.
- It detects superficial veins by infrared light and projects the pattern on the skin surface right above veins.
- Projected vascular system enables medical staff choosing a vein of necessary size and position for venipuncture and other procedures.
- It illuminates only superficial vascular system. Maximal vein illumination depth depends on patients' category.
- Reasons for different illumination depth are skin conditions (for example, eczema, tattoo), hair, scratches or very rough skin and lubricous (greasy) skin.
- If positioned directly above studied area the equipment determines the middle part of a vein.
- Width of projected vein may differ from actual width because of physical diversities of patients, but vein middle line remains precise.
- Equipment power supply is implemented from a battery that ensures prolonged operation.
- No training is required to work with the equipment.
- The equipment requires no regular or preventive maintenance but regular cleaning.

1.1 Warnings and precautions

Please learn the below operational safety rules to avoid damages:

- Equipment is mainly used to find superficial veins but it shall not replace accurate medical conclusion based on visual and tactile vein finding and professional skills.
- Equipment is an auxiliary unit for medical staff to find veins.
- Equipment only projects a pattern of superficial veins. Vein finding depth is limited due to a patient's peculiarities.
- Equipment only reflects superficial veins patter and provides no diagnostic information.
- Do not place the equipment in fire because it will lead to explosion of battery.
- Equipment contains no parts to be maintained by an operator. All parts shall be purchased at the Manufacturer or its official representative.
- Equipment shall be charged only by a charger supplied. This will help to avoid eventual risk.
- Reliability of equipment may be infringed if foreign accessories are used.
- Use only supplied power mains charger to charge the battery.

- Repair of equipment and replacement of battery shall be performed by an attested professional certified by the Manufacturer.
- Do not expose eyes to the projector light for a prolonged time.
- Keep equipment clean and avoid vibrations.
- Equipment shall be dry. Avoid moisture in equipment and network adapter.
- High temperature sterilization is not required. Do not sterilize by electron beam or by γ radiation.
- Avoid electromagnetic or high frequency interferences near equipment because strong electromagnetic interferences emitted by such sources as wireless sensors or mobile phones and etc. aggravate normal operation.
- Inspect equipment each two months. If evident damage is revealed, please contact official representative of the Manufacturer.
- Check safety labels. Check operative condition of equipment and its compliance with description provided in User Manual.

Notice: Equipment makes no clinical diagnosis. Equipment contains no medications.

Equipment may be used by emergency departments, oncology, pediatrics, laboratory, vascular and plastic surgery and with patients of any age, weight and health status as specified in medical practice rules.




Note: Check equipment and its parts before use. Do not use equipment in case of any damages.

1.2 Regulatory documents

Equipment corresponds to below regulatory and technical requirements:

EN 60601-1:2006+A11+A1+A12	Medical electrical equipment and systems. General safety requirements providing for main functional characteristics.
EN 60601-1-2:2014	General safety requirements providing for main functional characteristics. Parallel standard. Electromagnetic compatibility. Tests and requirements.
EN 60601-1-2:2015	General safety requirements providing for main functional characteristics. Parallel standard. Electromagnetic compatibility. Tests and requirements.
ISO13485:2012	Medical devices – Quality management systems – Regulatory requirements.
ISO13485:2012/AC:2012	AC edition – Medical devices – Quality management systems – Regulatory requirements.
EN ISO 14971:2012	Medical devices – Risk management application to medical devices.

1.3 Symbols

	EU compliance. Equipment totally complies with European Directive 93/42/EEC
	Recycling material
	Direct current voltage
	Warning
	Manufacturer
	Date of manufacture
	Disposal as household waste is prohibited (EU Directive on identification of disposed electric and electronic equipment)
	Manufacturer's contact phone
	See User Manual
	See Operating Manual
	Serial number
	European representative
	Equipment interface
	Charger interface
	Fragile, handle with care
	Avoid rain, keep dry

	Place vertically
	Do not turn over
	Vertical stack limitation: 5
	Package may be recycled

1.4 Proposed use and contraindications

Equipment may be used both at hospital and outside because it is indispensable when dealing with oncology patients, drug-dependents, tuberculosis and ambustial patients, full-figured people, children, dark skin patients, old people and patients exposed to prolonged intravenous infusions. These categories usually have acute and chronic venous insufficiency thus making vein finding for diagnostics and treatment difficult.

Equipment will help to:

- Accelerate vein finding because sometimes time is critical;
- Accelerate vein finding with patients having poor venous bed condition;
- Decrease injuries during venipuncture thus decreasing patient’s discomfort and promoting faster tissue healing;
- Conduct medical procedures in a very short time; it is life essential in some situations.
- Minimize tissue damage and promote healing.

Do not use the equipment for eyes veins or near eye veins!!!

1.5 Environment and indoor requirements

Operating environment shall be free from noise, vibration, dust, corrosive, explosive or flammable materials.

Rooms where equipment will be used shall meet the following requirements:

- No thermal shocks;
- Ambient temperature, moisture and pressure shall comply with the following values:

Temperature	(10-30) °C
-------------	------------

Relative humidity	(10÷90) % (no condensing)
Pressure	(70-106) kPa

- No corrosive or flammable gases;
- Severe impacts or vibrations are prohibited;
- No sprinkling of oil or chemical reagents.

Note: Duration of self-sufficient operation may be sufficiently decreased if equipment is used in conditions other than recommended by the Manufacturer (i.e. extremely hot or cold ambient temperature).

1.6 Battery

- Internal battery in equipment itself.
- The battery shall be fully charged before use of equipment.
- Avoid battery depletion; charge the battery before discharge.
- Charge the battery if it is in protection mode after short circuit or overcharging.
- Check charging and discharge every two months and contact official representative of the Manufacture if any problem appears.
- Run-down batteries shall not be disposed with household waste. Contact waste disposal company for details.
- Charge the battery with network adapter supplied.

Note: If battery power is less 3%, the equipment emits continuous sound. Push any button to stop sounding (about 0.5 seconds). Immediately charge the battery.

Charge indicator is in upper right corner of projected image. Charge the equipment is battery power is low. Observe the following charging procedure:


1. Connect charger with network adapter to alternating current line.
2. Shift power switch to OFF position.
3. Connect equipment to the charger.
4. The equipment is charging if the light is blue.

Note: Temperature of network adapter may be higher 48 °C and a person who contacts it may feel discomfort or pain. Do not touch the adapter for more than 1 minute.

State of battery charge may be detected by light indicator (1):

- Blue – battery is not fully charged, continue charging.
- Green – battery is fully charged. We recommend to disconnect charger.
- Blinking light – battery temperature may be outside the required one or battery is not charging. Check equipment, charger and network adapter, and connect network adapter again to remedy the problem.

Note: Operational life of the battery is 2 years. Replace the battery after this time expired.


Note: Wrong connection or opening of network adapter is prohibited. Life threatening! Network adapter has voltage-carrying electronic parts.

2. Functional characteristics


Equipment projects near infrared beams on a patient's skin. Blood hemoglobin is responsive to this and this vein pattern is reflected on skin. The process is immediate and requires no preparation. We use digital full-featured visualization technology that produces good resolution image. This image may be displayed and color and size may be changed.

No installation or adjustment is required before use.

Key functional elements:


Imaging Unit

Charging Cradle

2.1 Control board

1. ON/OFF: switching power on/off: the equipment switches off when the switch is in OFF position. Shift the switch to ON position to switch the equipment on.

2. Power: power button.

Power on: hold button “Power” for about 0.5 seconds until the equipment beeps twice. The equipment will be powered on.

Power off: hold button “Power” for about 2 seconds until the equipment beeps twice. LED indicators will turn off and the equipment will be powered off.

Standby mode: press “Power” button (less 0.5 seconds) with the equipment powered on and LED indicators lighting. The equipment will go to standby mode and switch off.

Active mode: press “Power” button (less 0.5 seconds) in standby mode. The equipment will beep and LED indicators will light.

3. Color: Button for changing the color of an image projected.

Five colors of image projection are available: light blue, green, natural, red and blue. Light blue color is set as default.

4. Size: Button for choosing the size of an image projected

Three sizes of image projected are available: full size, 2/3 and 1/2. Full size is set as default.

5. Inverse: Button for changing brightness and contrast.

Projected image of a blood vessel is brighter and surrounding tissues are dimmed. Blood vessel projection is set dimmed as default.

6. Mode: Button for changing the mode of an image projected.

Hold “Mode” button (about 0.5 seconds) to change brightness. Three brightness degrees are available. Shift UNIVERSAL mode to FINE detailed mode by holding the button (about 2 seconds). Second brightness degree and universal mode are set as default.

Note: Inversion and size of projected image may be chosen both in universal and fine modes. Previous mode immediately cancels after activation of another mode.

Note: Universal: standard mode of image projection for all patients. Fine mode: to find small and deep veins.


2.2. Focus control

- Place the equipment exactly above investigated area at the distance about 200 mm.
- Switch the equipment on and ensure the projection port is perpendicular to investigated area.
- Slightly adjust the height to make text on the side of an image clear. Focusing is completed.

Note: Image of a vein projected may be thicker or thinner real vein due to the distance between equipment and skin, physical peculiarities of a patient and different investigation modes. We recommend medical conclusion made under standard procedures and experienced staff in administration of medications shall prevail.

2.3 Projected image description

Projected image is given below. Text frame contains information about name of equipment, current mode and battery condition

3. Risks related to use of medical device

Possible risk	Risk reduction measures
Power factors	
High temperature	Network adapter temperature may exceed 48 °C and thus it may cause discomfort or pain. Contact time limitation is less 1 minute.
Stored energy	Battery capacity is 15.48 watt-hours. Thus the battery may explode if exposed to high temperature. Do not place the equipment in fire because it may lead to battery explosion. Battery may be replaced only by attested specialist certified by an authorized manufacturer.
Operation risks	
Wrong measurement	Assess adjustment by checking division between projected line and printed line on calibration sheet. If obvious division between projected and printed lines persists, stop using the equipment and contact authorized representative of the Manufacturer.
Adverse feedback due to output resulted from electromagnetic disturbance	Strong electromagnetic disturbance may cause equipment reset. The equipment will restore after reboot.
Information risks	
Wrong description of operating characteristics	Operating characteristics comply with clause 1 “Medical device intended use” and clause 2 “Functional characteristics” of this User Manual
Wrong specification of intended use	Equipment may be used by emergency departments, oncology, pediatrics, laboratory, vascular and plastic surgery and with patients of any age, weight and health status as specified in medical practice rules.
Wrong specification of accessories designated for use with medical device	Accessories include mobile support and fixed support, see clause 5 “Description of Accessories” of User Manual

4. Specifications



Item	Class
Electric shock protection	Equipment with embedded power source without working parts
Insulation class	Class II
Water and solid particles protection	IPX0
Software safety	Class A
Time rating	Continuous
Operational life of equipment	5 years
Operational life of battery	2 years

Parameter	Characteristics
Weight	480 g
Dimensions	60x60x240 mm
Battery	Embedded, li-ion battery, 3.6V, 4.3 A/h
Storage energy time requirement	5 hours
Charging time	3 hours
Network adapter characteristics	Rated input: 100-240V AC, 50-60 Hz, 0.4 A Rated output: 5V DC, 2 A
Light source for vein finding	Near infrared light
Peak length of infrared wave	~850 nm
Visible light wave length	~530 nm
Image processing speed	≤25 frame / s
Vein image accuracy	0,25 mm
Vein recognition depth	10 mm
Image resolution	854x480 pixels
Optical emission	None
Focus distance	About 20 cm
Color of image projection	5 colors
Image size	3 sizes
Brightness degree	3 degrees
Inverse mode	Yes
Quality improvement/fine image adjustment mode	Yes
Standby mode	Yes
Low battery signal	Yes

Ambient pressure (for operation, storage and transportation)	700 hPa ~ 1060 hPa
Moisture during operation	10-90 % (without condensate)
Temperature during operation	10 °C ~ 30 °C
Moisture during storage and transportation	0-90 % (without condensate)
Temperature during storage and transportation	0 °C ~ 60 °C
Mobile support	
Frame, dimensions	ABS plastic 5 rays, 600*600*190 mm
The length of the support bar	645 mm
Diameter of the support bar	38 mm
The length of the flexible piston	700 mm
Diameter of the flexible piston	18 mm
Wheels	5 pcs. diameter 75 mm
Weight	4600 gm
Fixed support	
Bracket	Aluminium alloy, 50*70*95 mm
Capture bracket	Rectangular 3 ÷ 55 mm
The length of the flexible piston	50 mm
Diameter of the flexible piston	18 mm
Weight	850 gm

5. Description of accessories

Accessories include mobile support and fixed support.

Mobile support	
	Light, mobile and easy to move wheeled support that fully complies with staff convenient working space demands. Advantages of wheeled support – mobility, light weight and prompt movement. Structural size and stability are perfect for equipment placement and transportation. The support has convenient moving handle. Metal frame ensures necessary stability with small area occupied. It is easy for assembling and requires no additional maintenance.
Fixed support	
	It is needed for fixing on horizontal or vertical surface, precise positioning and holding of equipment. No additional maintenance is required.

6. Maintenance and repair

Note: The equipment requires no regular and preventive maintenance.

The equipment was designed in such a way that no electric or mechanical parts shall be repaired by a buyer and/or operator. Do not open the device and do not replace electric or mechanical parts. Always contact official representative of the Manufacturer.

Even slightest interference to the equipment’s integrity will void manufacturer’s warrantee.

Unauthorized interference may affect compliance of the equipment with valid rules and may reduce operators’ or patient’s safety. Only trained specialists certified by the Manufacturer may conduct repair and troubleshooting.

- Do not put the equipment in liquid. Ensure no liquid contacts with the equipment.
- Avoid pressure or heat sterilization of the equipment.
- Do not clean the equipment by organic solvents and rough cloth to avoid damages.

Note: Do not use the equipment in case of any abrasions or dirt on the lens window. Contact authorized representative of the Manufacturer if window is damaged.

6.1 Troubleshooting

Trouble	Reason	Resolution
No image projected after equipment is loaded	Battery low	Charge battery
	Power source failure	Reset after repeated switching on
	Operational failure	Reset after repeated switching on
Loading failure	System failure	Reset after switching off
		Reset after repeated switching on
No vein pattern	Operational failure	Reset after repeated switching on

Note: In case of any operational abnormalities outside listed above please contact authorized representative of the Manufacturer.

7. Disinfection and cleaning

We supply non sterile equipment and it shall be cleaned. If a working space is severely soiled or there is a lot of dust or sand, clean the equipment more frequently. Consult equipment cleaning, disinfection and sterilization rules of the hospital prior cleaning.

- Do not clean the equipment on mobile of fixed support.

- Disconnect network adapter and ensure the equipment is switched off prior cleaning.
- Clean lens window every month. Use clean and soft cloth impregnated with 95% spirit. Air dry everything.
- Outer surface of the equipment may be cleaned by clean and soft cloth, sponge or cotton wool impregnated with 70% spirit twice a month. We recommend to twist excess detergent before cleaning.

Note: Clean carefully and avoid scratches on lens window.

Note: Regular check is needed. Comply with cleaning and sterilization rules of a hospital before each use.

8. Electromagnetic compatibility

Medical devices shall be operated with care due to electromagnetic effect. Thus, equipment shall be placed and/or used under information contained in supporting documents (in this case see tables below). Hand-held or mobile radio equipment may affect the equipment operation.

Risk of cross interference with another equipment:

Equipment causes no interference for other nearby medical devices used during medical examinations and treatment. The equipment shall not be connected to another device because it contains embedded power source.

Risk of electromagnetic interferences and remedial measures:

Equipment shall not be used in immediate neighborhood, in contact or in crossed electromagnetic fields with other electric or radio equipment. If such placement is necessary or may not be avoided, take special precautions for ensure correct operation at basic configuration (for example, constant and visual control of no abnormal operation signs or faults).

Below tables contain information about main characteristic of electromagnetic compatibility of this equipment. The equipment is in CISPR 11 category of medical devices (Industrial, research and medical equipment. Radio interference characteristics. Limit values and methods of group 1 measurements) and complies with class A requirement.

If electromagnetic interference of the equipment or other nearby electric devices is possible we recommend to change location of interfering equipment or devices (phone, radio, mobile antennas). Place the equipment away the working area if possible or switch off not needed equipment.

Electromagnetic radiation:

The equipment is designated to be used in below electromagnetic conditions. Operator shall guarantee compliance of equipment with requirements listed.

Radiation control	Compliance	Electromagnetic conditions - Guideline
Radiation radio frequency CISPR11	Group 1	Equipment uses radio-frequency energy for internal operation only. Thus, radio frequency radiation of equipment is very low and create almost no interference for nearby electronic devices.
Radiation radio frequency CISPR 11	Class A	Equipment may be used in all non-household institutions, in household institutions and in institutions directly connected to public low-voltage mains of household buildings provided the following warning will be communicated: Warning: for medical staff only. May cause radio interference or affect operation of nearby equipment. The following measures may be taken to reduce effect: realignment, moving or shielding the area.
Harmonic radiation IEC 61000-3-2	Class A	
Power fluctuations and oscillations IEC 61000-3-3	Complies	

Electromagnetic interference resistance test:

The equipment is designated to be used in below electromagnetic conditions. Operator shall guarantee compliance of equipment with requirements listed.

Resistance test	Registered level IEC 60601	Minimal level of compliance	Electromagnetic conditions – Guideline
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV at contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	±8 kV at contact ±2 kV, ±4 kV, ±8 kV, ±15 kV in air	Floor shall be wooden, concrete or ceramic. If floor is coated with synthetic materials, relative air humidity shall not be less 30%.

Conductivity/pulses IEC 61000-4-4	± 2 kV for power line ± 1 kV for input/output line	± 2 kV for power line	Network characteristics shall meet commercial or steady conditions.
Overvoltage (surges) IEC 61000-4-5	± 1 kV between phases ± 2 kV between phase and earth	± 1 kV between phases	Network characteristics shall meet commercial or steady conditions.
Voltage drop, Short interruptions and change of power supply on power mains IEC 61000-4-11	$< 5\% U_T$ (break $> 95\%$ in U_T) for half-cycle $40\% U_T$ (break 60% in U_T) for 5 cycles $70\% U_T$ (break 30% in U_T) for 25 cycles $< 5\% U_T$ (break $> 95\%$ in U_T) for 5 seconds	$< 5\% U_T$ (break $> 95\%$ in U_T) for half-cycle $40\% U_T$ (break 60% in U_T) for 5 cycles $70\% U_T$ (break 30% in U_T) for 25 cycles $< 5\% U_T$ (break $> 95\%$ in U_T) for 5 seconds	Network characteristics shall meet commercial or steady conditions. If the equipment shall be used during power failures, we recommend to connect it to a UPS or a battery.
Magnetic field with power mains frequency (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	In case of abnormal operation place the equipment away sources of magnetic field with power mains frequency or create a magnetic shield. Measure magnetic field with power mains frequency prior the equipment installation. Use equipment in areas with low magnetic field only.
NOTE: U_T means accrual current voltage in power line before test			

Test for resistance to induced and radiated interference:

The equipment is designated to be used in below electromagnetic conditions. Operator shall guarantee compliance of equipment with requirements listed.

Interference resistance test	Test registered level IEC 60601	Minimal level of compliance	Electromagnetic conditions – Guideline
Conductive radio frequency interference under IEC61000-4-6	3 Vrms from 150 kHz to 80 MHz	3 Vrms from 150 kHz to 80 MHz	Mobile and hand-held communication tools that use radio frequency range including cables shall not be used in close neighborhood with the equipment. Observe minimal recommended distance calculated under equation for such equipment transmitter frequency.
Radiated radio frequency interference under IEC61000-4-3	3 V/m from 80 MHz to 2.7 GHz Resistance test of enclosure port for radio frequency wireless equipment (ref. Table 9 of EC 60601-1-2:2014)	3 V/m from 80 MHz to 2.7 GHz Resistance test of enclosure port for radio frequency wireless equipment (ref. Table 9 of EC 60601-1-2:2014)	Recommended distance: $d=[3,5/V_1] \times P^{1/2}$ $d=1.2 \times P^{1/2}$ from 80 MHz to 800 MHz $d=2.3 \times P^{1/2}$ from 800 MHz to 2,5 GHz where <P> is maximal rated output power in watt (W) specified by the transmitter manufacturer and <d> is recommended distance in meters (m). Field density value from fixed radio frequency transmitters measured on site (a) shall be less compatibility level in each frequency range (b). Interference may appear close to equipment marked with the following symbol:

NOTE 1: The highest frequency ranges are applied for 80 MHz and 800 MHz

NOTE 2: It is possible that specified guidelines may not be applied to all situations. Buildings, objects and people absorption and reflection influence electromagnetic field propagation.

a) Electromagnetic field strength value radiated by fixed transmitters such as cell towers for mobile phones and mobile radio systems, amateur transmitters, AM and FM radio transmitters and TV radio transmitters may not be precisely specified theoretically. Electromagnetic measurements on site are necessary to define influence of such transmitter. Please ensure that the equipment works correctly if field strength on site where equipment is used is higher the above level of radio frequency compliance. Additional measures such as change of equipment orientation or positioning may be needed if equipment works incorrectly.

b) Field strength shall not exceed 3 V/m away the range from 150 kHz to 80 MHz.

Recommended distance from hand-held and mobile radio frequency equipment to the equipment:

The equipment is designated for use in electromagnetic environment with controlled radio frequency oscillations. An operator may limit electromagnetic interference by minimal distance from hand-held or mobile radio frequency communication equipment (transmitters) in compliance with below recommendations which are based on maximal output power of communication equipment.			
Maximal rated output power of transmitter, W	Distance to equipment based on transmitter frequency, m		
	from 150 kHz to 80 MHz $d=1.2\times P^{1/2}$	from 80 MHz to 800 MHz $d=1.2\times P^{1/2}$	from 800 MHz to 2,5 GHz $d=2.3\times P^{1/2}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
As for transmitters with maximal output power not specified in the table, we recommend to define distance d (in meters) under equation applied for transmitter frequency where P is maximal rated output power in watt (W) as specified by transmitter manufacturer.			

Note 1 Use distance applied to higher frequency range if frequency is 80 and 800 MHz.

Note 2 It is possible that specified recommendations may not be applied to all situations. Buildings, objects and people absorption and reflection influence electromagnetic field propagation.

9. Transportation and storage

Avoid extremely high or cold temperatures and moisture or dust when storing and transporting the equipment.

Store and transport the equipment in compliance with the following requirements: 0 ~ 90% (without condensate), temperature: 0 °C ~ 60 °C., ambient pressure 700 hPa ~ 1060 hPa.

Store away sources of flammable and explosive gases. Avoid strokes and vibrations during transportation.

10. Marking and package

Label of equipment shall contain the following:

- Name of the Manufacturer;
- Serial number and date of production;
- Name of the equipment and model;
- Power supply voltage, frequency and intake power (on power input);
- Address of manufacture and contact phone;
- International mark under IEC 60601-1 and ISO 10079-1 related to the type of equipment, category and CE marking (Directive 93/42/EEC);
- Symbol – ref. User Manual if needed;
- Symbol – Do not dispose as household waste.

Label of charger shall contain the following:


- Name and model of product;
- Serial number and date of manufacture (year, month);
- Power supply voltage, frequency and intake power (on power input and output);
- Symbol – ref. User Manual if needed;
- Symbol – Charger's interface;

- Symbol – Do not dispose as household waste (EU Directive on identification of disposed electrical and electronic equipment).


Label of Mobile support shall contain the following:

- Name of the equipment and model;
- Serial number and date of production;
- Name of the Manufacturer.

Label of Fixed support shall contain the following:

- Name of the equipment and model;
- Serial number and date of production;
- Name of the Manufacturer.

Specially created cardboard box is designated to protect the equipment during transportation or storage. On the one side there is information about storage and transportation conditions and on another side there is a graphical image about package content.

Meaning of six symbols on the package:

- 

Place vertically
- 

Fragile, handle with care
- 

Avoid rain, keep in dry place
- 

Vertical stack limitation: 5
- 

Limit on the number of stacked tiers: 5
- 

Do not turn over
- 

Package may be recycled

11. Disposal

Equipment is made of plastics and electronic components. No specific disposal recommendations are required. Dispose in accordance with applicable medical waste disposal rules and requirements for disposal of class A waste.



Do not dispose together with household waste. Не утилизируйте аппарат вместе с бытовым мусором. Contact certified electronic waste acceptance and processing factory to clarify disposal procedures. Addresses of such factories may be found at an authorized environment authority or at municipality.





Package (cardboard and inserts) may be disposed as paper waste.




Batteries may not be disposed as household waste. Contact waste disposal company for further information.





Annex 1

Accessories

Mobile support	Fixed support
	
	

Mobile support assembly stages	
1  Wheel Frame	Mount wheels by pushing to respective opening in the frame
2  Support bar Spacer A1 Spacer A2 Screw A1	Connect support bar to the frame by screw A1 and spacer A1 and A2 using hex-nut wrench.
3 	Place metal ring of the handle on support bar.
4  Screw A2	Connect the handle to support bar by screw A2 using hex-nut wrench A2.
5  Screw A3	Screw flexible piston to support bar and fix it with screw A3 using hex-nut wrench A3.

6 	Put ring lock on upper part of flexible piston, insert metal ball to rear part of the support and screw the ring lock on the rear part of the support, Ensure that support may keep the equipment at any angle.
7 	Click the locking block, insert the equipment in the support and turn the locking block for fixation.
8 	Assembly is done.

Fixed support assembly stages	
1 	Put ring lock on upper part of flexible piston, insert metal ball to rear part of the support and screw the ring lock on the rear part of the support, Ensure that support may keep the equipment at any angle.
2 	Click the locking block, insert the equipment in the support and turn the locking block for fixation.
3 	Turn locking screw and fix the bracket on the side of a table.
4 	Assembly is done.

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