**AUTHENTICATION OF
SIGNATURE OF A PERSON
SIGNING ON BEHALF OF A
BODY CORPORATE OR IN THE
NAME OF ANOTHER PERSON
ON A DECLARATION**

I, the undersigned, Sigalit Tzipora Rubin Flor Notary holding license no. 2063442 hereby certify that on 26.12.2019 appeared before me at my office located at 86 Hadkalim st. Kanionit Hakikar Pardes hanna Karkur, Israel.

Mr. **Levin Mark Meir** whose identity was proved to me by Identity Booklet No. 014140941 issued on 02.09.2007

And Mr. **Topaz Tal** whose identity was proved to me by Identity Booklet No. 040901399 issued on 02.08.2015

and after I was satisfied that they understand the language of the declaration-English, that is known to me as well, they read in my presence the declaration attached and marked with the letter A1-13

And after having inquired and being satisfied that Mr. **Levin Mark Meir** and Mr. **Topaz Tal** understood the contents of the affidavit and after I warned them that they must tell the truth and that they will be subject to the penalties prescribed by law if they fail to do so, they confirmed the correctness of the affidavit and signed it voluntarily in my presence on behalf of corporation

Paltop Advanced Dental Solutions Ltd company number 514458744 located at Hashita 5, Industrial Park Caesarea, Israel . I confirm that I have been submitted the document letter of signature rights issued by Company's lawyer on 3.7.2019 as written evidence, to my satisfaction, for the purpose of proving their competence to sign as aforesaid.

In witness whereof, I hereby authenticate the signature of Mr. **Levin Mark Meir** and Mr. **Topaz Tal** by my own signature and seal this day 26.12.2019

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To whom it may concern,

The Company Paltop Dental Advanced Solution Ltd, company no. 514458744, Hashita 5 Industrial Park IL- P.O Box 3568- 3088900 Caesarea, Israel. represented by CEO, Tal Topaz, and by CFO, Meir Levin acting on the basis of the Charter, declare that the information provided here below is accurate and complete in all respects.

Instruction for Use (11-50026R Rev.02)

Signature: _____

Tal Topaz
CEO

Date: 26/12/2019

Signature: _____

Meir Levin
CFO

Date: 26.12.19

Instruction for Use - Paltop Dental Implant System

11-50026R Rev 2 dated: 28/10/2019

Warning

The following descriptions are insufficient to allow immediate use of the Paltop Dental Implant System. Knowledge of dental implantology and instruction in the handling of the Paltop Dental Implant System provided by an operator with the relevant experience are always necessary.

The Paltop Dental Implant System has not been evaluated for safety and compatibility in the MR environment. The Paltop Dental Implant System has not been tested for heating or migration in the MR environment.

Do not use if package is damaged.

Disclaimer of liability

This product is part of the Paltop Dental Implant System and may only be used in conjunction with the corresponding original components and instruments according to Paltop LTD. instructions and recommendations. Use of products and tools made by third parties detracts from the efficient functioning of the Paltop Dental Implant System and will void any warranty or other obligation, expressed or implied by Paltop LTD. Advice on the use of our products is provided verbally, in writing, by electronic media or in demonstrations. It does not relieve the user of Paltop LTD. products of the responsibility of determining whether or not any product is suitable for the intended purpose, indications and procedures. Use of this product is not within the control of Paltop LTD. it is the responsibility of the user. All liability for loss or damage attributable to the use of this product is excluded. Within the framework of the Conditions at Sale and Delivery of Paltop LTD, we warrant the quality of the Paltop Dental Implant System.

Warranty

All Paltop LTD. products are guaranteed to be free from defects in material and workmanship at the time of shipping. All of our products are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals, for any failure due to misuse or application which is not in accordance with the designer's intentions. Products will lose their right to carry the CE sign in case they will be changed or repaired at another place than Paltop LTD or a company authorized by Paltop LTD to do the change or repair.

Availability

Some items in the Paltop Dental Implant System are not available in all countries.

Copyright and trademarks

Paltop and/or other products and logos from Paltop that are mentioned here are trademarks or registered trademarks of Paltop LTD. and/or its affiliates

Advanced Dental Solutions Ltd.

PALTOP®

With you to the TOP

A3K

PALTOP DENTAL IMPLANTS

Intended Use

The Paltop Dental Implant System is indicated for use in surgical and restorative applications for placement in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, in order to restore the patient's chewing function. The Paltop Dental Implant System is indicated also for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The Paltop Dental Implant

Paltop dental implant system is essentially a substitute for a natural root. Once fitted, these components provide the foundation for long-term support of crowns, bridges or dentures.

The Paltop dental implant is a screw type implant with internal hexagonal or conical connection, which is intended for single use.

The Narrow Implant is indicated for use in surgical and restorative applications for placement in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws where the horizontal space is limited by the adjacent teeth and roots, to provide support for prosthetic devices, such as artificial teeth, in order to restore the patient's chewing function. The Paltop Narrow Implant is indicated also for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The Paltop Dental Prosthetic component is a pre-manufactured component directly connected to the endosseous dental implant and is intended for use as an aid in prosthetic rehabilitation.

The Paltop Surgical Instrument is essentially, as other dental surgical instruments are available in the market, a stainless steel or titanium device intended for use in reconstructive oral surgery to drill or cut into the upper or lower jaw and may be used to prepare the bone to insert a dental implant.

Conditions for use, transport and storage

Products can be used, transported and stored in temperature between -18°C and 55°C and RH between 0 to 100%.

Disposal

After the use the device has to be disposed as a contaminated healthcare waste

1. Description of the Paltop surfaces

Paltop dental implants are made of biocompatible TI 6AL 4V ELI [Titanium Grade 5] and have a large grit sandblasted, acid-etched surface.

2. Indications and Contraindications

The Paltop Dental Implant System is indicated for use in surgical and restorative applications for placement in the bone of the upper or lower jaw to provide support for prosthetic devices, such as artificial teeth, in order to restore the patient's chewing function. The Paltop Dental Implant System is indicated also for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The Narrow Implant is indicated for use in surgical and restorative applications for placement in the mandibular central, lateral incisor and maxillary lateral incisor regions of partially edentulous jaws where the horizontal space is limited by the adjacent teeth and roots, to provide support for prosthetic devices, such as artificial teeth, in order to restore the patient's chewing function. The Paltop Narrow Implant is indicated also for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

Contraindications

Serious internal medical problems, bone metabolism disturbances, uncontrolled bleeding disorders, inadequate wound healing capacity, poor oral hygiene, maxillary and mandibular growth not completed, poor general state of health, uncooperative, unmotivated patient, drug or alcohol abuse, psychoses, prolonged therapy-resistant functional disorders, xerostomia, weakened immune system, illnesses requiring periodic use of steroids, titanium allergy, uncontrollable endocrine disorders.

Relative contraindications

Previously irradiated bone, diabetes mellitus, anticoagulation drugs / hemorrhagic diatheses, bruxism, parafunctional habits, unfavorable anatomic bone conditions, tobacco abuse, uncontrolled periodontitis, temporomandibular joint disorders, treatable pathologic diseases of the jaw and changes in the oral mucosa, pregnancy, inadequate oral hygiene.

Local contraindications

Inadequate bone volume and/or quality, local root remnants, in whom adequate sizes, numbers or desirable position of implants are not reachable to achieve safe support of functional or eventually parafunctional loads.

3. Principles of treatment planning

The surgical part of implant treatment must be preceded by a comprehensive patient evaluation, preoperative diagnostics and therapy planning. Inadequate treatment planning may cause implant loss.

General precautions

One hundred percent implant success can never be guaranteed.

With respect to pediatric patients, routine treatment is not recommended until the end of the jaw bone growth has been properly documented. Pre-operative hard tissue or soft tissue deficits may yield a compromised esthetic result or unfavorable

implant angulation. The device is intended for single use only. reuse of the device may cause microbial contamination and loss of performance.

In case Short implants are used, clinician should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal.

Precautions Pre-surgery

A careful clinical and radiological examination of the patient has to be performed prior to surgery to determine the psychological and physical status of the patient. Special attention has to be given to patients who have localized or systemic factors that could interfere with the healing process of either bone or soft tissue or the osseointegration process (e.g., cigarette smoking, uncontrolled diabetes, oro-facial radiotherapy, steroid therapy, infections in the neighboring bone).

In general, implant placement and prosthetic design must accommodate individual patient conditions. In case of bruxism or unfavorable jaw relationships reappraisal of the treatment option may be considered.

Precautions at Surgery

Particular caution should be used when placing narrow platform implants in the posterior region due to risk of prosthetic overload. Lack of adequate quantity and/or quality of remaining bone, infection and generalized diseases may be potential causes for failure of osseointegration both immediately after surgery, or after osseointegration is initially achieved.

Besides the mandatory precautions for any surgery such as of asepsia, during drilling in the jaw bone one must avoid damage to nerves and vessels by referring to anatomical knowledge and preoperative radiographs. Failure to recognize actual lengths of step/twist drills relative to radiographic measurements or drilling beyond the depth intended can result in permanent injury to nerves or other vital structures, potentially resulting in permanent numbness to the lower lip and chin or leading to hemorrhage in the floor of the mouth.

4. Side-effects and interactions, complications with Paltop implants

Immediately after insertion of dental implants, activities that demand considerable physical exertion should be avoided. Possible complications following insertion of dental implants are:

Temporary symptoms

Pain, swelling, phonetic difficulties, gingivitis

More persistent symptoms

Chronic pain in connection with the dental implant, permanent paresthesia, dysesthesia, loss of maxillary/mandibular ridge bone, localized or systemic infection, oroantral or oronasal fistulae, unfavorably affected adjacent teeth, irreversible damage to adjacent teeth, fractures of implant, jaw, bone or prosthesis, esthetic problems, nerve damage, exfoliation, hyperplasia.

It is recommended that the dental personnel brief the patients on the precautions to be taken, such as: abstinence from smoking and maintenance of dental hygiene.

Adverse events

Treatment by means of implants may lead to loss of bone, biologic or mechanical failures including fatigue fracture of implants.

5. Placement of an implant after preparing the implant bed

- 5.1 Before the procedure, general contra-indication of the patients medical status should be evaluated clinically and radio graphically.
- 5.2 Bone quality of the patient should be examined.
- 5.3 Implant length and diameter should be selected.
- 5.4 The surgical field should be exposed.
- 5.5 Using a marking drill an initial drill should be performed.
- 5.6 The drilling should gradually create a suitable hole into which the implant will be inserted.
- 5.7 The implant will be inserted into the hole.
- 5.8 Using hex drives the screw should be removed.
- 5.9 Cover screw, healing screw or abutments should be applied according to the clinician's judgment.

Do not use a high torque or high speed when inserting the implant (do not exceed 15 rpm).

6. Healing phase for Paltop dental implants

6.1 Immediate restoration of Paltop implants

Paltop implants allow, within the scope of indications, immediate restoration in single tooth gaps and in an edentulous or partially dentate jaw. Good primary stability and an appropriate occlusal load are essential. Two or more adjacent implants should be prosthetically connected together.

6.2 Delayed restoration of Paltop implants

Delayed restoration is possible after 6 weeks:

if bone quality is good and bone volume adequate

using $\varnothing 4.2$ mm or more implants with a length of 8mm or more

Delayed restoration is possible after 8 week:

with cancellous bone quality

using $\varnothing 3.75$ mm or more implants

Advanced Dental Solutions Ltd.

PALTOP®

With you to the TOP

A ≠ X

using implants with a length of 6 mm or more

In situations where the implant surface is not completely in contact with the bone, or where bone augmentation measures are necessary, a healing phase appropriate to the situation must be planned. A radiographic examination is recommended after a healing phase of 6-12 weeks, before prosthetic restoration.

7. Packaging, Sterility, and Sterilization Method

Paltop dental implants are delivered sterile. The intact sterile packaging protects the gamma-sterilized implant from outside influences and, if stored correctly, ensures sterility up to the expiration date. When removing the implant from the sterile packaging, the rules of asepsis must be observed. The sterile packaging must not be opened until immediately prior to insertion of the implant. Implants with damaged sterile packaging must not be used due to risk of contamination. It is recommended to have a replacement implant on hand. Paltop Prosthetic Components and Surgical Tools are provided either sterile or non-sterile according to the indication on the labeling. Paltop Surgical Kits are provided Non-Sterile. Products provided as non-sterile must be sterilized prior to use. Standard autoclave steam sterilization is recommended using a pre-vacuum method with a minimum exposure time of 4 minutes (four pulses) at 135 °C / 275 °F in an FDA cleared sterilization pouch. While drying times may vary according to load size, a minimum of 16 minutes is recommended. This cycle has been validated by the overkill method to a sterility assurance level (SAL) of 10⁻⁶ according to requirements or guidances outlined in ISO 17665-1 Sterilization of health care products -- Moist heat -- Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical device, ANSI/AAMI ST79 and ISO 17665-2 Sterilization of health care products -- Moist heat -- Part 2: Guidance on the application of ANSI/AAMI/ISO 17665-1. Prosthetic Components are to be used immediately after sterilization.

8. Shelf Life

For sterile products only, the shelf life of the product has been validated for 5 years. There is no shelf life for products which are not provided sterile.

PALTOP DENTAL PROSTHETICS

1. Description of Raw Material

Paltop dental prosthetics are made of biocompatible Ti 6AL 4V ELI, Gold parts (Au60%Pt19%Pd20%) and Zirconia (ZrO₂).

2. Indications

Paltop Dental Prosthetic component is a pre-manufactured component directly connected to the endosseous dental implant and is intended for use as an aid in prosthetic rehabilitation. Paltop Dental Prosthetic component are intended for single use. Reuse of the device may cause microbial contamination and loss of performance.

3. Contraindications

Do not use in cases of hypersensitivity to one or more of the metals contained in the alloy.

4. General precautions

Collaboration between the surgeon, restorative dentist and dental laboratory technician is essential for a successful implant treatment.

Temporary PEEK abutments are recommended for use for up to 30 days.

5. Cleaning

We recommend cleaning and/or disinfecting abutment before sterilization. We recommend using chlorhexidine gluconate 0.1~0.5w/v% solution, or alcohol 70% for 10~30 min, at room temperature. Do not use such disinfectants as Chlorine-based disinfectant (sodium hypochlorite, etc.), Peroxide-based disinfectant (oxydol, etc.), non-rust-proofing disinfecting surfactant. Place abutment in ultrasonic cleaner for 5 minutes at room temperature. If any cleaner is used, be sure to inject sterile water and completely remove the cleaner, running water can be used instead of sterile water. Completely dry and inspect abutment for integrity and flaws.

6. Sterility

The device may be delivered sterile or non-sterile. Please see indication on sticker. The device must be sterile before use. it may be steam-sterilized using Standard Sterilization Methods: steam sterilization for 10 minutes at 275° F (135° C). Steam sterilization is necessary only in case of NON-STERILE products or to re-sterilize "Not Used Devices".

7. Clinical Procedure

1. Modifications of abutments could be performed using appropriate tools. Use water irrigation if adjustments are performed intra-orally.
2. Prior to placement inspect the implant-abutment interface. Damage to the interface can lead to improper fitting. Do not process the interface area.
3. Fasten the abutment with the provided screw using a recommended torque value as per table below.
4. After cementing the prosthesis on the abutment excess cement should be removed immediately.

8. Storage:

No special condition required.

SURGICAL INSTRUMENTS

1. Caution

Please read all information contained in this insert. Incorrect handling and care as well as misuse can lead to premature wear of surgical/dental instruments. All persons using this device should have knowledge in the use and handling surgical instruments, accessories and related equipment. Paltop Surgical Instruments are reusable and should be replaced after being used up to 20 times, or prior to that if cutting efficiency declines.

2. Initial Use of New Instrument

Unless removed from the sterile packaging, every instrument must be cleaned and sterilized before it is used for the first time. All instruments were developed for sterilization by autoclave and have shown good results using this method.

3. Inspection and Functional Check

It is very important to carefully examine each surgical/dental instrument for breaks, cracks or malfunctions before use. It is especially essential to check areas such as blades, points, ends, stops and snaps as well as all movable parts. Do not use damaged instruments. Never attempt to make repairs yourself. Service and repairs should be referred to trained qualified persons only. Refer questions about repair to the manufacturer or your biomedical engineering department.

4. Cleaning and Maintenance

Every surgical/dental instrument should be disinfected and thoroughly cleaned after each use. Proper cleaning, inspection and maintenance will help ensure correct function of the surgical instrument. Clean, inspect and test each instrument carefully. Sterilize all instruments before surgery. A good cleaning and maintenance procedure will extend the useful life of the instrument.

Special attention should be paid to slots, stops, ends hollow tubes and other highly inaccessible areas. Do not use damaged instruments. If an instrument can be disassembled, it should be cleaned in the disassembled state.

(NOTE: Keep all parts together and protected from being misplaced, do not intermix with other or similar devices of other manufacture).

Cleaning and rinsing must take place immediately after each use for best effect. Failure to clean promptly may result in adherent particles or dried secretions that may resist cleaning and complicate or resist future sterilization. Instruments must be completely cleaned and rinsed of all foreign matter. Use warm water and a commercially available instrument presoak or cleaning agent. Enzymatic cleaners should be used to remove protein deposits.

•

Do not use corrosive cleaning agents (i.e. bleach). Cleaning solutions and rinses at or near a neutral pH (7.0) are best.

•

Do not use abrasive cleaners.

•

Only a soft bristle brush should be used.

•

Rinse and maintain all parts and inaccessible areas like inside channels etc.

•

Ultrasonic cleaners can be used. Check and retighten any fittings that may have vibrated loose.

•

Can be disinfected in the washing machine up to 203° F (95° C).

•

Rinse thoroughly with distilled water.

•

Prepare for storage and / or sterilization.

After cleaning and rinsing, dry instruments completely and carefully with compressed air (Highly inaccessible areas like inside channels etc. have to be blown out).

NOTE: after cleaning and before sterilization, treat all instruments with an oil which is considered as being physiologically safe (paraffin oil according to DAB 8 of Ph. Eur. USP XX), especially their blades, ends, stops, snaps and all movable parts.

5. Storage and Re-use

Instruments should be stored in a clean, dry moisture free area. The instruments should be stored individually in their shipping carton or in a protective tray with partitions. Protect tips with cloth, gauze or tubing if stored in drawers. Steam autoclave sterilization is recommended. Thoroughly clean instruments of all debris,

Advanced Dental Solutions Ltd.

PALTOP®

With you to the TOP

A 104

tissue and foreign matter prior to sterilization. Follow the sterilizer manufacturer's instructions for operation and loading of steam autoclaves. There must be direct steam exposure to all surfaces of the instruments being sterilized including the internal surface and tubes channels. Allow instrument to air cool to room temperature before use.

6. Standard Sterilization Methods

The device may be delivered sterile or non-sterile. Please see indication on sticker. The device must be sterile before use. it may be steam-sterilized using Standard Sterilization Methods: steam sterilization for 10 minutes at 275° F (135° C).

7. Handling

All surgical/dental instruments should be handled with the greatest care when being transported, cleaned, treated, sterilized and stored. This is especially true for blades, fine points and other sensitive areas. Surgical/dental instruments corrode and their functions are impaired if they come into contact with aggressive materials. The instruments should not be exposed to acids or other aggressive cleaning agents.

8. Warranty

All Paltop LTD. products are guaranteed to be free from defects in material and workmanship at the time of shipping. All of our products are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals, for any failure due to misuse or application which is not in accordance with the designer's intentions. Products will lose their right to carry the CE sign in case they will be changed or repaired at another place than Paltop LTD or a company authorized by Paltop LTD to do the change or repair.

Manufactured by

Paltop Advanced Dental Solutions LTD

Hashita 5 PO Box 3568

Casearea 30889

Israel

Tel +972-4-6271711, Fax +972-4-6275363

MedNet GmbH

Borkstrasse 10

48163 Muenster

Germany

Tel +49-251-32266-60, Fax: +49-251-32266-22

Authorized representative in Russia:

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CE 1936

A 2 K

REF	Catalog Number
LOT	Batch Code
	Caution, consult accompanying documents
	Do not reuse
	Non-Sterile
	Do not use if package is damaged
STERILE R	Sterilized using gamma radiation
Rx only	By Prescription Only
	Use by date
	Do not resterilize
	Manufacturer
EC REP	EC representative

	Spade	Step Drills				Final Drills			
		2.0 1.5	2.4 2.0	3.2 2.4	3.8 3.2	3.25	3.75	4.2	5.0
Speed RPM	1200- 1500	900- 1100	700- 900	500- 700	400- 600	200 - 400			

Prosthetic / Prothèse / Prothese / Protesi / רכיב	Torque Value / Valeur de couple / Drehmomentwert / Coppia di serraggio / ערך מומנט הידוק (N/Cm)
Healing Cap / Capuchon de cicatrisation / Heilabdeckung / Cappetta di guarigione / כיפת ריפוי	15
Screw / Vis / Schraube / Vite / בורג 80-70002	32-35
Screw / Vis / Schraube / Vite / בורג 80-70005	25-28
Screw / Vis / Schraube / Vite / בורג 80-70010	28-30
PEEK Abutment / PEEK / Piliers PEEK / Monconi in PEEK / מבנה PEEK	15
Multi-Unit Abutments / Prothèse Vissée / Mehreinheit / Monconi Multi Unit / מולטי יוניט	32-35
Ball Abutments / Prothèse Amovible / Herausnehmbarer Zahnersatz / Monconi A Sfera / מבנה כדורי	32-35

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