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
« June 24 2019


INSTRUCTION FOR USE FOR MEDICAL DEVICE

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Sworn to and subscribed
before me this
24 day of June, 2019

**LIESA S. DeNARDO
NOTARY PUBLIC OF NEW JERSEY
COMMISSION EXPIRES 06/12/2022**

Straumann® Membrane Flex™

Resorbable Collagen Membrane For Dental Surgery

INSTRUCTIONS FOR USE

Description:

Straumann® Membrane Flex™ Resorbable Collagen Membrane is a white, nonfriable, resorbable, single layered, conformable collagen membrane matrix manufactured from purified porcine peritoneum. Straumann® Membrane Flex™ is sterilized by gamma irradiation and is supplied sterile, non-pyrogenic, and for single use only.

Intended use:

Straumann® Membrane Flex™ is a bioresorbable, implantable collagen material that is intended for use in oral surgical procedures as a resorbable membrane material for use in: simultaneous use of guided bone regeneration (GBR)-membrane and implants;

- augmentation around implants placed in immediate extraction sockets;
- augmentation around implants placed in delayed extraction sockets;
- localized ridge augmentation for later implantation;
- alveolar ridge reconstruction for prosthetic treatment;
- filling of bone defects after root resection, cystectomy or removal of retained teeth;
- guided bone regeneration in dehiscence defects; and
- guided tissue regeneration procedures in periodontal defects.

Contraindications:

Straumann® Membrane Flex™ is contraindicated in patients who have:

- acute infections or contaminated wound in the oral cavity;
- a known allergy to collagen; and/or
- clinically significant renal, hepatic, cardiac, endocrine, hematologic, autoimmune or systemic disease, which in the physician's judgment, will prevent safe implantation or likely healing.

Administration:

Straumann® Membrane Flex™ is packaged in a double sterile pouch. The outer pouch should be opened carefully, allowing the inner pouch to be placed onto a sterile field. The membrane should be removed from the inner pouch with sterile gloves or instruments.

The bone defect is exposed by a mucoperiosteal flap and basic surgical procedures are performed (e.g. curettage). Space-making material such as autologous bone, demineralized bone matrix and ceramic materials may be used to fill the defect.

Straumann® Membrane Flex™ can be placed either dry or hydrated. If the clinician prefers the handling characteristics of the hydrated collagen, the membrane can be hydrated in sterile water or saline solution for approximately five minutes prior to final placement.

Straumann® Membrane Flex™ can be trimmed to the size and shape of the defect in the dry or wet state using sharp, sterile scissors.

Straumann® Membrane Flex™ should overlap the walls of the defect by at least 2 mm to allow complete bone contact and to prevent gingival connective tissue invasion below the material.

Fixation of the membrane may be indicated to avoid displacement due to loading or mobilization. The membrane can be sutured in place using absorbable sutures and a noncutting needle.

Resorbable tacks can also be used to affix the membrane. The mucoperiosteal flap is sutured over the collagen membrane and the wound should be closed completely to avoid accelerated resorption due to membrane exposure. The membrane is expected to be resorbed in approximately 12 to 16 weeks.

Post-operative procedures:

Straumann® Membrane Flex™ is completely resorbable and should not be removed. Patients should rinse with an antimicrobial agent such as chlorhexidine gluconate (Peridex) twice daily for four weeks following surgery. Beginning 24 hours after surgery, the wound site may be additionally swabbed with a cotton-tipped applicator dipped in the antimicrobial agent.

The patient should refrain from brushing the treated area for two weeks following the surgery. After this period, the patient may be instructed to gently brush the area with a soft toothbrush. Dental floss should not be used prior to four weeks following surgery. Coronal scaling and prophylaxis can be performed at follow-up visits, if indicated.

The patient should be seen seven to ten days following surgery for wound evaluation and removal of any closing sutures or periodontal packing. These follow-up visits should be repeated every two weeks thereafter, up to eight weeks following surgery. The patient may return to normal oral hygiene routine.

Probing and subgingival scaling should not be performed prior to six months following surgery to prevent damage to immature tissues. Other assessments of clinical health may be repeated, including plaque, bleeding and tooth mobility indices.

Warning:

Clinicians should use care in screening their patients for any known allergies to collagen. Hypersensitivity reactions have been noted with the use of other products containing collagen; therefore, the possibility exists of developing a local sensitivity response to *Straumann® Membrane Flex™*.

Precautions:

As with all surgical procedures, caution should be exercised when treating medically compromised patients such as patients receiving long-term steroidal therapy or currently taking anticoagulants. Patients with clinically significant systemic diseases, indicating a history of anaphylactic reactions, autoimmune diseases, uncontrolled diabetes or severe hypertension have not been implanted with the membrane; therefore, the safety and effectiveness for those patients have not been determined. Nor has it been evaluated in pregnant women, children and/or in patients with conditions involving extremely severe defects with little periodontium or bone.

Straumann® Membrane Flex™ cannot be resterilized or reused. Open, unused *Straumann® Membrane Flex™* must be discarded. In vivo stability may be adversely affected if resterilized. Cross-contamination and infection may occur if re-used.

Do not use if the product sterilization barrier or its packaging is compromised.

Adverse reactions:

Possible complications that can occur with any dental surgery include infection, swelling of the intraoral tissue, thermal sensitivity, gingival recession, excessive gingival bleeding, flap sloughing, resorption or ankylosis, with loss of crestal bone height, pain, or complications associated with the use of anesthesia. Minor discomfort may occur for a few days.

Storage:

The product should be stored at room temperature (15°C/59°F to 30°C/86°F). Avoid excessive heat and humidity.

How Supplied:

Straumann® Membrane Flex™ is supplied sterile, non-pyrogenic and for single use only. One (1) membrane per package, size 15 x 20 mm, 20 x 30 mm, or 30 x 40 mm.

Caution:

Federal (U.S.A.) law restricts this device to sale by or on the order of a physician or dentist.

Labeling Symbols:

Symbols may be used on some international package labeling for easy identification.

	Caution, Consult accompanying documents
	Use by date
	Do not re-use
	Do not use if package is damaged

LOT	Batch code
STERILE R	Sterilized using irradiation
	Temperature limitation
REF	Catalog number
R_x only	Caution: Federal (U.S.A.) law restricts this device to the sale by or on the order of a physician or dentist
	Manufacturer
EC REP	Authorized European Representative

Distributed by (USA):

Straumann USA, LLC
60 Minuteman Road
Andover, Massachusetts 01810 USA

Distributed by:

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland

Manufacturer:

Collagen Matrix, Inc.
15 Thornton Road
Oakland, New Jersey 07436 USA

Authorized European Representative:

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Russia specific Storage and Product Information:

Technical characteristics:

Straumann® Membrane Flex™	Variant 15 mm x 20 mm	Variant 20 mm x 30 mm	Variant 30 mm x 40 mm
Acidity (extractable pH)	6 - 9		

The product is a white, nonfriable, semi-permeable, resorbable, single layered, conformable collagen membrane matrix and meets the requirements of the international standard BS EN 22803 "Dentistry. Membrane materials for guided tissue regeneration in oral and maxillofacial surgery. Contents of technical File".

The product is biocompatible (nontoxic) and meets the requirements of international standards ISO 10993 "Biological evaluation of medical devices" and ISO 7405 "Dentistry. Evaluation of biocompatibility of medical devices used in dentistry".

Application condition:

Medical institutions of dental profile.

Disposal:

Shelf (sterility) life of the product is 3 years from the date of manufacture under observance of storage and transportation regulations.

The product shall be disposed of in accordance with regulatory documentation and requirements applicable to this type of medical product in the country of distribution.

Warranty:

Straumann® Membrane Flex™ do not contain any hazardous, non-toxic substances that would require special handling.

Shelf (sterility) life of the product is 3 years from the date of manufacture under observance of storage and transportation regulations.

The manufacturer guarantees the operability and maintenance of the sterility of the products during the entire shelf life of 3 years, subject to storage and transportation conditions. The product can be transported by any means of transport, such as covered railway wagons, containers, freight vehicles, ship holds, aviation modules in accordance with the rules for these to modes of transport, in compliance with the conditions of transportation without risk of damage to the package. The product does not affect the environment during operation. Minor excursions from the recommended temperature range during storage, shipping, and handling are acceptable. The manufacturer cannot guarantee the product to meet specifications, if stored at temperatures outside of the recommended temperature range for more than 3 days.

Complaints:

Issues on application of the medical product, undesirable events identified in course of product application, and quality complaints from the customers in the territory of the Russian Federation shall be directed to authorized representative at below contact details:

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