Sworn to and subscribed
before me this
13th day of February, 2020.

MATEUSZ F. LESZCZAK
NOTARY PUBLIC OF NEW JERSEY
Commission # 50065913
My Commission Expires 08/10/2022

«УТВЕРЖДАЮ» / «APPROVE»
ЭТИКОН, Эл-Эл-Си (ETHICON, LLC)
475 C Street Los Frailes Industrial Park
Suite 401, Guaynabo, Puerto Rico
00969, USA

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1

INSTRUCTION FOR USE

Absorbable surgical suture unidirectional with antibacterial
coating with needle for knotless tissue fixation STRATAFIX™
Spiral MONOCRYL™ Plus in configurations

Version for use on the territory of Russian Federation

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NAME OF PRODUCT

Absorbable surgical suture unidirectional with antibacterial coating with needle for knotless tissue fixation STRATAFIX™ Spiral MONOCRYL™ Plus in configurations.

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DESIGNER

ETHICON, INC.,
Route 22 West, P.O. Box 151, Somerville, NJ 08876, USA

PHYSICAL MANUFACTURER

Ethicon, LLC,
Highway 183 Km 8.3, San Lorenzo, Puerto Rico, 00754, USA

INTENDED USE

The STRATAFIX™ Spiral MONOCRYL™ Plus Device is indicated for use in soft tissue approximation where the use of absorbable sutures is appropriate.

USAGE CONDITIONS

This medical device is intended for professional use in the medical and health care institutions.

AREA OF USE

Surgery

POTENTIAL CONSUMER

Surgical Suture STRATAFIX™ Spiral MONOCRYL™ Plus is intended to be used by professional doctor, who is aware about surgical procedures, techniques and materials, as well as how to work with absorbable surgical materials.

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INDICATIONS

The STRATAFIX™ Spiral MONOCRYL™ Plus Device is indicated for use in soft tissue approximation where the use of absorbable sutures is appropriate.

CONTRAINDICATIONS

The STRATAFIX™ Spiral MONOCRYL™ Plus is not to be used where extended approximation of tissues under stress is required and is not to be used in conjunction with or for fixation of prosthetic devices (e.g., heart valves or synthetic grafts) that are nonabsorbable in nature. This device, being absorbable, should not be used where extended approximation of tissue under stress is required, such as in fascia. The STRATAFIX™ Spiral MONOCRYL™ Plus Device should not be used in patients with known allergic reactions to IRGACARE* MP (triclosan).

WARNINGS

- The fixation loop of the STRATAFIX™ Spiral MONOCRYL™ Plus Device should be deployed by users familiar with surgical procedures, techniques and tissues and the use of absorbable devices before employing the STRATAFIX™ Spiral MONOCRYL™ Plus Device for wound closure, as risk of wound separation may vary with the site of application and the suture material used.
- Physicians should consider the quantity and quality of tissue in which the loop will be anchored.
- The use of this suture may be inappropriate in highly vascularized tissue or fragile tissue that cannot withstand potential further tightening/constriction within the loop.
- Care should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.
- Avoid contacting STRATAFIX™ Device and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.
- Physicians should consider the in vivo performance (under ACTIONS section) when selecting a device for use in patients. The use of this device may be inappropriate in elderly, malnourished or debilitated

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patients, or in patients suffering from conditions which may delay wound healing.

- As with any foreign body, prolonged contact of any device with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.
- As an absorbable device, STRATAFIX™ Spiral MONOCRYL™ Plus may act transiently as a foreign body.
- Acceptable surgical practice should be followed for the management of contaminated or infected wounds.
- Small bowel obstruction (SBO): including volvulus, bowel infarction, and significant morbidity, have been reported due to barbs or barbed suture ends hooking onto adjacent small bowel and/or mesentery, such as in peritoneal closure.
- As this is an absorbable device, the use of supplemental nonabsorbable devices sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention or which may require additional support.
- The safety and effectiveness of the STRATAFIX™ Spiral MONOCRYL™ Plus Device have not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, osseous tissue, tendinous tissue, ophthalmic surgery or for use in microsurgery.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.
- Discard open packages and unused devices. Discard packaging that is damaged or opened.
- Do not use after expiration date.

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PRECAUTIONS

- The STRATAFIX™ Spiral MONOCRYL™ Plus Device contains a loop end and barbs to anchor tissues and does not require knots to secure the device. Tying of knots on the barbed section of the device will damage the barbs and potentially reduce their effectiveness. For the device to function properly, the STRATAFIX™ Spiral MONOCRYL™ Plus Device must first be anchored by using the fixation loop in robust tissue. Then subsequently engage the barbs into the tissue in a standard closure pattern. When completing placement of the barbed segment in subcuticular tissue, take at least one additional pass in the reverse direction to lock the STRATAFIX™ Spiral MONOCRYL™ Plus Device in place.
- Skin devices that remain in place longer than 7 days may cause localized irritation and should be snipped off or removed if indicated.
- Subcuticular devices should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption.
- Under some circumstances, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Consideration should be taken in the use of absorbable devices in tissue with poor blood supply as device extrusion and delayed absorption may occur. In handling this or any other device material, care should be taken to avoid damage from handling.
- Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.
- Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound separation are typical or foreseeable risks associated with any suture and hence are also potential complications associated with the STRATAFIX™ Spiral MONOCRYL™ Plus Device.
- Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may

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cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks.

- Discard used needles in “sharps” containers.

ADVERSE REACTIONS

Adverse effects associated with the use of synthetic absorbable devices include tissue separation, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distention occurs, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, device extrusion and delayed absorption in tissue with poor blood supply and transitory local irritation at the wound site, as well as allergic reaction to IRGACARE® MP (triclosan). Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

DESCRIPTION

The STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Device is an antibacterial monofilament, synthetic absorbable device prepared from a copolymer of glycolide and ε-caprolactone. The device contains IRGACARE® MP (triclosan), a broad spectrum antibacterial agent, at no more than 2360 µg/m. Poliglecaprone 25 copolymer has been found to be nonpyrogenic and elicits only a slight tissue reaction during absorption. The STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control

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Device, Variable Loop Design consists of barbed suture material, armed with a surgical needle on one end and a fixation loop at the opposite end. The STRATAFIX™ Spiral MONOCRYL™ Plus Device is designed to anchor with a closed loop at one end and a unidirectional barbed section on the other end. The STRATAFIX™ Spiral MONOCRYL™ Plus Device barbs are oriented in one direction to allow tissue approximation without the need to tie surgical knots. While the formation of barbs in the STRATAFIX™ Spiral MONOCRYL™ Plus Device reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strengths. For this reason, the strength of the STRATAFIX™ Spiral MONOCRYL™ Plus Device can be compared to USP knot strength of non-barbed sutures. USP designations for diameter are used to describe the STRATAFIX™ Spiral MONOCRYL™ Plus Device suture material after barbing, except for minor variation in suture diameter with a maximum overage of 0.1 mm.

The USP and EU Pharmacopoeia sizes of the STRATAFIX™ Spiral MONOCRYL™ Plus Device are further defined in Table 1.

Table 1.

USP DEVICE SIZE DESIGNATION	EU PHARMACOPOEIA DEVICE SIZE (Metric/Ph.Eur.) DESIGNATION	STRATAFIX™ Spiral MONOCRYL™ Plus	
		USP	Metric/Ph.Eur.
2-0	3	2-0	3
3-0	2	3-0	2
4-0	1,5	4-0	1,5

TENSILE STRENGTH

STRATAFIX™ Spiral MONOCRYL™ Plus Device straight tensile strength is stronger than the knot tensile strength for a USP and EU Pharmacopoeia poliglecaprone 25 device of the equivalent size as shown in Table 2.

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ACTIONS

The results of implantation studies in animals using the STRATAFIX™ Spiral MONOCRYL™ Plus Device indicate that approximately 62% of its original strength remains after 7 days implantation and approximately 27% of its original tensile strength at 14 days post implantation. All of the original tensile strength is lost by 21 days post implantation. The absolute strength remaining 14 days post implantation meets or exceeds that historically observed with plain or chromic surgical gut devices. Data obtained from implantation studies show that absorption of STRATAFIX™ Spiral MONOCRYL™ Plus Devices is essentially complete by 91 days (Table 3).

Table 3.

DAYS IMPLANTATION	APPROXIMATE % ORIGINAL STRENGTH REMAINING
7 days	62% of USP knot tensile strength
14 days	27 % of USP knot tensile strength

Using *in vitro* studies, STRATAFIX™ Spiral MONOCRYL™ Plus Device with IRGACARE® MP (triclosan) has been shown to inhibit colonization of the device by Staphylococcus aureus, Staphylococcus epidermidis, methicillin-resistant S. aureus, methicillin-resistant S. epidermidis, Escherichia coli, and Klebsiella pneumoniae, which are microorganisms known to be associated with surgical site infections. Animal studies have demonstrated that STRATAFIX™ Spiral MONOCRYL™ Plus Device

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inhibits bacterial colonization of the device after direct in vivo challenge with bacteria.

RISK ANALYSIS

The evaluation of risk management was conducted in accordance with the International Standards Organization EN ISO 14971:2012, entitled Medical Devices – Application of risk management to medical devices.

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INSTRUCTIONS FOR USE

STRATAFIX™ Spiral MONOCRYL™ Plus Device is designed to be used in continuous suture patterns. STRATAFIX™ Spiral MONOCRYL™ Plus Device is intended to be used without anchoring knots to begin or terminate the device line.

Use as required per surgical procedure.

1. Secure the fixation loop portion to robust tissue by taking a bite in the designated tissue, passing the needle through the loop and pulling tautly around the anchoring tissue.

2. After STRATAFIX™ Spiral MONOCRYL™ Plus Device is anchored into the tissue, grasp the strand and approximate tissue to the desired tension. Bites, or passes through the tissue, can be taken in a continuous suturing technique to approximate the wound. Care should be taken to utilize the STRATAFIX™ Spiral MONOCRYL™ Plus Device on the barbed segments only. Do not attempt to approximate wounds using the non-barbed segment near the needle, as the barbs are required for successful wound approximation with the STRATAFIX™ Spiral MONOCRYL™ Plus Device.

3. a. To complete the closure in subcuticular closure, take at least one pass in the reverse direction. Then pass the needle across the incision and take a split thickness bite perpendicular to the incision and exit the skin.

b. To complete the closure in all other tissue layers, take at least two passes in the reverse direction to complete the closure.

4. Then gently pull on the free end of the device and cut flush with the surface of the tissue.

To remove the entire STRATAFIX™ Spiral MONOCRYL™ Plus Device: Cut the STRATAFIX™ Spiral MONOCRYL™ Plus Device at the end anchored with the loop, and then pull the distal end to remove the device.

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MATERIALS

The STRATAFIX™ Spiral MONOCRYL™ Plus Suture:

Suture- Poliglecaprone 25 copolymer, a copolymer of glycolide (75%) and ε-caprolactone (25%). Manufactured Ethicon Inc, Cornelia, state Georgia, USA.

Cover- antibacterial agent IRGACARE® MP (triclosan), manufactured BASF Grenzach GmbH, Germany, GMP certificate №DE_BW_01_GMP_2018_0014

Needle

Stainless steel AISI grade 420*, Supplied by Carpenter Technology Corp. Reading, PA or

Stainless steel AISI grade 4310*, Supplied by Zapp Medical Alloys GmbH or

Ethalloy Stainless steel*, manufactured by ATI Allvac, USA.

Needle coverage – Silicone (CAS Number 01330-20-7

63148-57-2), Supplied by Nusil Corporation or

Multipass (polydimethylsiloxane silicone with the addition of wax powders), Supplied by MicroPowders Inc and Nusil Corporation

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ACCESSORIES

The device in this file does not have accessories included.

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STERILITY

STRATAFIX™ Spiral MONOCRYL™ Plus Suture is a single-use device provided sterile and sterilized by Ethylene Oxide (EO). Sterilization method is validated according to ISO 11135-1:2014.

The ethylene oxide sterilization process for the Device is validated to provide a sterility assurance level (SAL) of 10⁻⁶.

The remaining concentration of EO and ethylene chlorohydrin is in accordance with ISO 10993-7:2008.

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HOW SUPPLIED

The STRATAFIX™ Spiral MONOCRYL™ Plus Device is available sterile, undyed. Each STRATAFIX™ Spiral MONOCRYL™ Plus Device is single armed with various needle sizes. Device is placed into a plastic tray packaging which is placed into peel open foil package. Twelve (12) pouches and an IFU are placed into a carton box. The boxes are wrapped with a film and tear tape is placed during the process to ease the opening of the boxes.

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CRITERIA OF UNSUITABLE CASES OF MEDICAL DEVICES USE

Do not use of the packaging is open or damaged. Discard open packages and unused devices.

EXPIRATION DATE

Shelf life of the product is 24 months.

LIFE CYCLE/ACTIONS

The results of implantation studies in animals using the STRATAFIX™ Spiral MONOCRYL™ Plus Device indicate that approximately 62% of its original strength remains after 7 days implantation and approximately 27% of its original tensile strength at 14 days post implantation.

Data obtained from implantation studies show that absorption of STRATAFIX™ Spiral MONOCRYL™ Plus Devices is essentially complete by 91 days.

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STORAGE AND TRANSPORTATION CONDITIONS

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Storage: -60°C до +30°C

transportation conditions: -60°C до +50°C

air humidity is N/A (packaged in moisture proof foil).

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WARRANTY

Manufacturer guarantees to the end user absence of defects and conformity to product specifications for the entire shelf life of the product under the conditions of adherence to storage, transportation and package integrity requirements. There are no warranties that extend beyond the description on the face hereof. Suitability for use of this medical device for any particular surgical procedure should be determined by the user in conformance with the manufacturer's instructions for use.

RECYCLING

Discard open packaging and unused devices.

Discard used needles in "sharps" containers.

Dispose the product and its residual elements or waste items in accordance with local laws and regulations.

After being used the products are potentially infectious waste.

Dispose of open packages and unused products.

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ENVIRONMENTAL REQUIREMENTS IN THE APPLICATION OF THE MEDICAL DEVICE

This medical device is in use, transportation and storage, has no negative effects on humans and the environment.

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