Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation

Extract from the Technical Documents

Support Device with Water Level Adjustment for Swimming Pool Rehabilitation, versions: 4x2 m, 4x3 m, 4x4 m, 4x5 m, 4x6 m, 4x7 m, 4x8 m, 4x9 m, 4x10 m, 5x2 m, 5x3 m, 5x5 m, 5x6 m, 5x7 m, 5x8 m, 5x9 m, 6x2 m, 6x3 m, 6x6 m, 6x7 m, 6x8 m, 7x2 m, 7x3 m, 8x2 m, 8x3 m, 9x2 m, 9x3 m, with accessories.


Manufacturer

EWAC Medical

The Netherlands

www.ewac.com

EWAC Medical

Contents Page

1. Preface	3
2. Intended use	3
3. Manufacturer	3
4. Design information	4
5. Operation, methods of use	7
6. Indications and contra-indications	8
6.1 Indications	8
6.2 Contra-indications	8
7. Technical data	9
8. Materials	11
9. Storage / Transportation	12
10. Operating Conditions	12
10.1 General	12
10.2 Requirements for pool water	13
10.3 Requirements for the water in the compression unit	13
10.4 Requirements for the environment	13
10.5. Maintenance (by the user)	13
10.6 Daily maintenance	13
10.7 Weekly maintenance	14
10.8 Monthly maintenance	14
10.9 Annual maintenance	15
11. Trouble shooting	16
12. Service and maintenance	16
13. Transport packing, package contents.....	17
14. Disposal	17
15. List of the international regulations/standards the medical device complies with	17
16. Service life.....	18
17. Use of risk management for medical devices.....	18
18. Requirements to environment protection	18

1. Preface

We very much appreciate your choice of the Support Device with Water Level Adjustment for Swimming Pool Rehabilitation¹. This unit incorporates the very latest technologies, and particular attention has been paid to reliability, accuracy, safety and a long working life.

In order to operate the Support Device with Water Level Adjustment for Swimming Pool Rehabilitation correctly and safely, it is very important that you should read these operating instructions carefully before use and keep it with the unit.

Keep this user manual in range for all involved personnel at all time.

The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation has been designed according to up-to-date safety and quality requirements and complies with the following standards:

- The Certificate of Conformity № RQA662680 valid until 20.12.2016 confirming that EWAC B.V. quality management system conforms to EN ISO 13485:2012;
- The EWAC B.V. Declaration confirming that EWAC B.V. quality management system meets the requirements of Annex II, Article 3 of the Directive EC 93/42/EEC.

2. Intended Use

The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation (hereinafter referred to as the "support device") is used for early post-stroke, post-traumatic or post-surgical rehabilitation of patients of any age group.

The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation is a modern, safe and easy-to-handle apparatus, which has been developed, produced, installed and calibrated with greatest care and attention. The support device is built in conformity to the European guidelines for medical electrical equipment: IEC 601-1 and for swimming pool equipment: EN 13451-11 and is CE marked. EWAC B.V. can not be held responsible for any discomfort to the operator or to the patient due to faulty diagnosis, misuse or mishandling of the equipment and or accessories, or due to misunderstanding of the operating procedures, or maintenance errors.

3. Manufacturer

EWAC B.V., Dulleweg 16-18 1721 CT Broekop Langedijk, the Netherlands

This support device shall be operated according to the valid extract from the technical documents and user manual.

4. Design information

Support Device with Water Level Adjustment for Swimming Pool Rehabilitation, versions: 4x2 m, 4x3 m, 4x4 m, 4x5 m, 4x6 m, 4x7 m, 4x8 m, 4x9 m, 4x10 m, 5x2 m, 5x3 m, 5x5 m, 5x6 m, 5x7 m, 5x8 m, 5x9 m, 6x2 m, 6x3 m, 6x6 m, 6x7 m, 6x8 m, 7x2 m, 7x3 m, 8x2 m, 8x3 m, 9x2 m, 9x3 m, with accessories (hereafter "MD" or "Support Device with Water Level Adjustment for Swimming Pool Rehabilitation, with accessories") includes:

I. The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation, variants:

1. 4 x 2 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 54 pcs;
- Tile clips, 108 pcs;
- Floor frame, dimensions 4X2m, 1 pc;
- Stainless steel cables, 20 m long, 4 pcs;
- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 120 pcs;
- Nuts M12, 120 pcs.

2). Hydraulic system:

- Compressor unit;
- High-pressure hose, 10 m long;
- Piston, 3,500 mm long, diameter is 23 cm.

3). Control box.

4). Step-down transformer.

5). Connecting cable, 15 m long.

6). User manual.

2. 4 x 3 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 80 pcs;
- Tile clips, 160 pcs;
- Floor frame, dimensions 4X3m, 1 pc;
- Stainless steel cables, 20 m long, 4 pcs;
- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 180 pcs;
- Nuts M12, 180 pcs.

2). Hydraulic system:

- Compressor unit;
- High-pressure hose, 10 m long;
- Piston, 3,500 mm long, diameter is 23 cm.

3). Control box.

4). Step-down transformer.

5). Connecting cable, 15 m long.

6). User manual.

3. 4 x 4 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 106 pcs;
- Tile clips, 212 pcs;

- Floor frame, dimensions 4X4m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 240 pcs;
 - Nuts M12, 240 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

4. 4 x 5 m

- 1). The height-adjustable artificial floor:
- Tiles, dimensions 30x50 cm, 134 pcs;
 - Tile clips, 268 pcs;
 - Floor frame, dimensions 4X5m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 300 pcs;
 - Nuts M12, 300 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

5. 4 x 6 m

- 1). The height-adjustable artificial floor:
- Tiles, dimensions 30x50 cm, 160 pcs;
 - Tile clips, 320 pcs;
 - Floor frame, dimensions 4X6m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 360 pcs;
 - Nuts M12, 360 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

6. 4 x 7 m

- 1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 187 pcs;
 - Tile clips, 374 pcs;
 - Floor frame, dimensions 4X7m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 360 pcs;
 - Nuts M12, 360 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

7. 4 x 8 m

- 1). The height-adjustable artificial floor:
- Tiles, dimensions 30x50 cm, 214 pcs;
 - Tile clips, 428 pcs;
 - Floor frame, dimensions 4X8m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 480 pcs;
 - Nuts M12, 480 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

1. 4 x 9 m

- 1). The height-adjustable artificial floor:
- Tiles, dimensions 30x50 cm, 240 pcs;
 - Tile clips, 480 pcs;
 - Floor frame, dimensions 4X9m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 540 pcs;
 - Nuts M12, 540 pcs.
- 2). Hydraulic system:
- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

9. 4 x 10 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 267 pcs;
 - Tile clips, 534 pcs;
 - Floor frame, dimensions 4X10m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 600 pcs;
 - Nuts M12, 600 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

10. 5 x 2 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 67 pcs;
 - Tile clips, 134 pcs;
 - Floor frame, dimensions 5X2m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 150 pcs;
 - Nuts M12, 150 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

11. 5 x 3 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 100 pcs;
 - Tile clips, 200 pcs;
 - Floor frame, dimensions 5X3m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 225 pcs;
 - Nuts M12, 225 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.

6). User manual.

12. 5 x 5 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 167 pcs;
- Tile clips, 334 pcs;
- Floor frame, dimensions 5X5m, 1 pc;
- Stainless steel cables, 20 m long, 4 pcs;
- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 375 pcs;
- Nuts M12, 375 pcs.

2). Hydraulic system:

- Compressor unit;
- High-pressure hose, 10 m long;
- Piston, 3,500 mm long, diameter is 23 cm.

3). Control box.

4). Step-down transformer.

5). Connecting cable, 15 m long.

6). User manual.

14. 5 x 7 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 234 pcs;
- Tile clips, 468 pcs;
- Floor frame, dimensions 5X7m, 1 pc;
- Stainless steel cables, 20 m long, 4 pcs;
- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 525 pcs;
- Nuts M12, 525 pcs.

2). Hydraulic system:

- Compressor unit;
- High-pressure hose, 10 m long;
- Piston, 3,500 mm long, diameter is 23 cm.

3). Control box.

4). Step-down transformer.

5). Connecting cable, 15 m long.

6). User manual.

15. 5 x 8 m

1). The height-adjustable artificial floor:

- Tiles, dimensions 30x50 cm, 267 pcs;
- Tile clips, 534 pcs;
- Floor frame, dimensions 5X8m, 1 pc;
- Stainless steel cables, 20 m long, 4 pcs;
- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 600 pcs;
- Nuts M12, 600 pcs.

2). Hydraulic system:

- Compressor unit;
- High-pressure hose, 10 m long;
- Piston, 3,500 mm long, diameter is 23 cm.

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18. 6 x 3 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 120 pcs;
 - Tile clips, 240 pcs;
 - Floor frame, dimensions 6X3m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 270 pcs;
 - Nuts M12, 270 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;

- Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

19. 6 x 6 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 240 pcs;
 - Tile clips, 480 pcs;
 - Floor frame, dimensions 6X6m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 540 pcs;
 - Nuts M12, 540 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

20. 6 x 7 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 280 pcs;
 - Tile clips, 560 pcs;
 - Floor frame, dimensions 6X7m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 630 pcs;
 - Nuts M12, 630 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

21. 6 x 8 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 320 pcs;
 - Tile clips, 640 pcs;
 - Floor frame, dimensions 6X8m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 720 pcs;
 - Nuts M12, 720 pcs.
- 2). Hydraulic system:

- Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
 - 4). Step-down transformer.
 - 5). Connecting cable, 15 m long.
 - 6). User manual.

22. 7 x 2 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 94 pcs;
 - Tile clips, 188 pcs;
 - Floor frame, dimensions 7X2m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 210 pcs;
 - Nuts M12, 210 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

23. 7 x 3 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 140 pcs;
 - Tile clips, 280 pcs;
 - Floor frame, dimensions 7X3m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 315 pcs;
 - Nuts M12, 315 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

24. 8 x 2 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 107 pcs;
 - Tile clips, 214 pcs;
 - Floor frame, dimensions 8X2m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 240 pcs;

- Nuts M12, 240 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

25. 8 x 3 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 160 pcs;
 - Tile clips, 320 pcs;
 - Floor frame, dimensions 8X3m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 360 pcs;
 - Nuts M12, 360 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

26. 9 x 2 m

- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 120 pcs;
 - Tile clips, 240 pcs;
 - Floor frame, dimensions 9X2m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;
 - Anchors, dimensions 200x8 cm, 4 pcs;
 - Bolts M12x60, 270 pcs;
 - Nuts M12, 270 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

27. 9 x 3 m


- 1). The height-adjustable artificial floor:
 - Tiles, dimensions 30x50 cm, 180 pcs;
 - Tile clips, 360 pcs;
 - Floor frame, dimensions 9X3m, 1 pc;
 - Stainless steel cables, 20 m long, 4 pcs;

Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation


- Anchors, dimensions 200x8 cm, 4 pcs;
- Bolts M12x60, 405 pcs;
- Nuts M12, 405 pcs.
- 2). Hydraulic system:
 - Compressor unit;
 - High-pressure hose, 10 m long;
 - Piston, 3,500 mm long, diameter is 23 cm.
- 3). Control box.
- 4). Step-down transformer.
- 5). Connecting cable, 15 m long.
- 6). User manual.

II. Accessories:

1. Easily-removable bars, 4 m long, 2 pcs.
2. Easily-removable wall bars, dimensions 200x120 cm, 2-5 pcs.
3. Balance board, dimensions 50x50 cm, 2-5 pcs.
4. Rings bar, dimensions 40x60x140-160 cm, 2-5 pcs.
5. Obstacle course, 60x30x30 cm, 2-5 pcs.
6. Stairs, 55x75x140cm, 1 pc.
7. Depth indicator, dimensions 25x10x10 cm, 1 pc.

Picture 1 (Draft). Height-adjustable artificial floor: 1 – tiles, 2 – frame, 3 – cable, 4 – high-pressure hose, 5 – piston.

Picture 2. Compressor unit

The support device (picture 1) is in height adjustable, without fixed steps, between a predetermined minimum and maximum depth. The height is 0 cm to 170 cm and all given dimensions are project depending.

The walking surface of the floor is made of polyester tiles with a non-skid surface fitted onto a stainless steel frame. These tiles can be removed individually for inspecting the floor or for cleaning the swimming pool section under the floor.

Stainless steel cables from four anchors fitted in the pool wall suspend the floor. The raising and lowering movement of the floor is produced by a hydraulic system, which operates on water pressure and is installed underneath the support device. This hydraulic system is made in full from stainless steel.

The hydraulic system consists of a cylinder unit with piston, rod guiding, cable wheels, travel limitation, etc., as well as the hose fracture valve mounted directly on the cylinder. If the high-pressure hosing is accidentally broken, this valve prevents the piston rod from unintentionally moving.

The hydraulic unit is connected to the floor conduit on the inside of the pool by a chemical resistant high-pressure rubber pipe (suitable for an operational pressure of 16 Bar and a burst pressure of 40 bars, with threaded couplings).

A compression unit is installed outside the pool (Picture 2), in the technical room, and

EWAC Medical

Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation

supplies the required water pressure for the hydraulic unit. The pressure piping for this unit is connected externally underneath the pool to the floor conduit, fitted for this purpose.

A control cabinet with 2 splash-proof switches for controlling the support device is fitted either in the wall or in the pool area.

In order not to affect the quality of the swimming pool water, the water circulation in the pool must not be hampered by the installation of an adjustable false floor.

With the EWAC construction and design, this hindrance is kept at a minimum. The corner points of the polyester tiles have notches so that, when assembled, a flow opening of 6 mm is automatically achieved which, together with the open supporting frame construction, prevents "dead" water occurring anywhere under the floor. The tile dimensions are 30 x 50 cm, the tiles are made from hot-pressed polyester. They are light blue and have a rough surface. These tiles are attached to the supporting frame by stainless steel retaining clips.

5. Operation, Methods of Use

The support device is operated by means of a wall mounted control box in the vicinity of the support device. The controls of the support device operate on a low voltage (12 V) and are therefore safe to use in a humid and wet environment. The control box itself is splash waterproof protected (IP 54).

The operator of the floor needs to constantly monitor the pool area while operating the support device.

WARNING: It is strictly forbidden to dive into a pool with a support device!

Diving is extremely dangerous, as the water depth cannot be estimated properly. The manufacturer strongly recommends indicating this in the hydrotherapy room where the pool with the support device is located.

EWAC Medical also recommends raising the support device to just below water level after use. This prevents persons from trying to access the pool in the undesirable circumstance in which the pool area is unattended.

Methods of use:

The first option of the support device is to set zero pool depth to wheel non-mobile patients sitting in their cart (wheelchair) up into the pool without any additional arrangements.

- Depth of 40 to 60 cm is to be set for preschool children. Such water level makes children feel confident, safe and free.
- Depth of 60 to 90 cm is to be set for children of junior to middle school age (depending on their age) supervised by an instructor.
- Depth of 100 to 130 cm is to be set for joint rehabilitation of parents and their children. Such depth allows children to reach the pool floor while parents can lessen most of their own weight. Type of exercises: swimming, games, walking with obstacle courses, ring game, ball games.
- Depth of 140 to 150 cm is suitable for different water exercises. Most often it is used for rehabilitation of post-stroke and orthopedic patients. Type of exercises: walking, exercising lower limbs while sitting in the underwater chair, underwater obstacles.
- Maximum submersion of the pool floor is set for swimming, aqua aerobics. No special procedures are used.

Patients perform disease-specific exercises – walking, running, leg raising or abduction, sitting-up.

6. Indications and contra-indications

6.1. Indications

- Traumas (fractures, ligamentous injury, muscle damage) and consequences
- Musculoskeletal diseases
- Neurological diseases, consequences of strokes and cerebro-cranial traumas associated with defect of dynamic stereotype
- Vascular abnormalities (including endarteritis), conditions like vegetovascular dystonia, migraines
- Emotional impairments: neurosis, chronic fatigue syndrome, depressed mood and qualm
- In pediatrics rehabilitation activities in swimming pool are useful in the presence of neurological disturbance, arrested development, scoliosis, cerebral palsy, hyperthymia and even autism.

6.2. Contra-indications

- Acute disease period;
- Catarrhal diseases;
- Skin diseases.

7. Technical Data

Name	Value	Tolerance
Support Device with Water Level Adjustment for Swimming Pool Rehabilitation		
Maximum allowed patient's weight	200 kg	-
Dimensions	min. 9 m ² and max. 48 m ²	±5%
Tile dimensions	30x50 cm	±5%
Length of the connecting cable	15 m	±5%
Floor frame dimensions	4 x 2 m; 4x3 m; 4x4 m; 4x5 m; 4x6 m; 4x7 m; 4x8 m; 4x9 m; 4x10 m; 5x2 m; 5x3m; 5x5 m; 5x6 m; 5x7 m; 5x8 m; 5x9 m; 6x2 m; 6x3 m; 6x6 m; 6x7 m; 6x8 m; 7x2 m; 7x3 m; 8x2 m; 8x3 m; 9x2 m; 9x3 m	±5%
Length of the stainless steel cables	20 m	±5%
Anchor dimensions	200 x 8 cm	±5%
Bolts	M12x60	±5%
Nuts	M12	±5%
Tile weight	2 kg	±5%
Total weight of 1 m ² artificial floor (frame + tiles)	60 kg	±5%
Weight of the transformer assembled	17 kg	±5%
Weight of the compression unit assembled	34 kg	±5%
Piston, dimensions	Length is 350 cm, diameter is 23 cm	±5%
Hydraulic medium	Potable water	-
Mains water supply	Cold water 15 mm	-
High-pressure hose	Diameter is 19 mm, length is 10 m	±5%
Drain connection	PVC 20 mm (open connection)	±5%
Mains voltage	3 x 230/400 V 50 Hz	-
Power consumption	max. 2,2 KVA (project depending)	±10%
Travel speed	approx. 700 mm/min. (project depending)	±10%
Maximum height adjustment	170 cm (project depending)	±5%
Burst pressure of the rubber pipe	40 bar	±5%

Extract from the Technical Documents. Support Device with Water Level Adjustment for
Swimming Pool Rehabilitation

Maximum pressure	16 bar	±5%
Control voltage	12 V	-
Easily-removable bars	Length is 4 m	±5%
Easily-removable wall bars	200x120 cm	±5%
Balance board	50x50 cm	±5%
Rings bar	40x60x140-160 cm	±5%
Obstacle course	60x30x30 cm	±5%
Stairs	55x75x140 cm	±5%
Depth indicator	25x10x10 cm The indicator is powered by a transformer down 5V, 250 mA power consumption	±5%
Number of rungs on the wall bars (pcs.)	2	-
The diameter of the bars and other components	50 mm	±5%
The strength of the cables breaking	2450 kg	±5%
Maximum withstand load rack	150 kg	±5%
Maximum withstand load rings	150 kg	±5%
The maximum withstand load of the Swedish wall	150 kg	±5%
Maximum withstand load beams	150 kg	±5%
Power consumption of the hydraulic system	2.5 kW	-
The range of depth measurements (Depth indicator)	from 0 to 170 cm	±5%
Discrete (Depth indicator)	0.01 m	-
Color display (Depth indicator)	red	-
The size of the control panel	3x16x26 cm	±5%
The control panel is powered from the transformer down	5B	-

Safety class: type B, applied part class I according to IEC EN 60601-1.

8. Materials:

Component	Material
Support device	Tiles – high-impact antistatic fiber-reinforced polyester Yapiser Dye Penn Color International B.V.
	Clips – AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%.
Easily-removable bars	AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Control box	Polyvinyl chloride of SG3 type
Easily-removable wall bars	AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Balance board	High-impact antistatic fiber-reinforced polyester Yapiser AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Rings bar	AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Obstacle course	AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Stairs	AISI 316L stainless steel (material number is 1.4401) Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Floor frame	AISI 316L stainless steel (material number is 1.4401)
Stainless steel cables	Alloy: C: 0.03%, Cr: 16-18%, Ni: 10-14%, Si: up to 1%, Mn: up to 2%, P: up to 0.045%, S: up to 0.03%, Mo: 2-3%
Anchors	
Bolts M12x60	
Nuts M12	
Piston	
(rod guiding, cable wheels,	

Extract from the Technical Documents. Support Device with Water Level Adjustment for
Swimming Pool Rehabilitation

travel limitation, hose fracture valve)	
High-pressure hose length 10 m	Styrene-butadiene rubber (SBS, SBS) Globalprene 3411

9. Storage / Transportation

The support device may be shipped by all modes of transport in covered vehicles in accordance with the current shipment rules and regulations for every type of transport.

When transported the support device is considered weather-proof at the temperature of -20° up to $+40^{\circ}$, relative air humidity of up to 90%, no mechanical stress, vibration load of 10 to 55 Hz, displacement amplitude of 0.35 mm and shock loads of 10g peak acceleration and 10 ms shock acceleration duration. Permitted angle of inclination is up to 10° .

When stored the support device is considered weather-proof at the temperature of -20° up to $+40^{\circ}$ and relative air humidity of up to 90%.

10. Operating Conditions

10.1 General

The frame of the support device is made from stainless steel quality AISI 316 L (material number 1.4401). Unfortunately, even stainless steel can corrode under certain conditions. Hereunder a list with possible causes for corrosion:

- a. Iron particles in the water can cause galvanic corrosion of the stainless steel.
- b. Too high concentration of acid and/or salt content can attack the stainless steel.
- c. Water, which is too hard or too soft, can cause deposits on the stainless steel.

10.2 Requirements for pool water

To avoid corrosion of the stainless steel, EWAC has drawn up the following guidelines for the pool water:

1. The hardness of the water should be lower than 8 German hardness degrees, that is equivalent to 6,4 English or Clark degrees.
2. The acid content should be between $\text{pH} = 7,2$ and $\text{pH} = 7,6$.
3. The free chlorine content should be between 0,4 and 1 mg/l.
4. The chloride (salt) content should not exceed 100 mg/l.
5. A 5-micron filter has to be installed in the water supply pipe of the bath. This prevents iron particles from entering the bath where they can cause galvanic corrosion. This is also the reason why the water supply pipe should be made of a non-ferrous material.
6. The pump capacity of the filtration system should be more than a quarter of the pool volume per hour.

As already stated: be sure that no iron particles, other than stainless steel AISI 316, are present in the bath because they can cause galvanic corrosion.

10.3 Requirements for the water in the compression unit

- The hardness of the water should be lower than 8 German hardness degrees, that is equivalent to 6,4 English or Clark degrees.
- The acid content should be between $\text{pH} = 7,2$ and $\text{pH} = 7,6$.

Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation

- A 5-micron filter has to be installed in the water supply pipe of the unit. This prevents small iron particles from entering the hydraulic cylinder, where they can cause galvanic corrosion. This is also the reason why the water supply pipe should be made of a non-ferrous material.
- The chloride (salt) content should not exceed 100 mg/l.

10.4 Requirements for environment

- The hydrotherapy area must be ventilated and conditioned. The room temperature should not exceed 30° C. The humidity must not exceed the 70 %.

The technical area must be ventilated in 4 x room contents per hour, because of the storage of chemicals, like chlorine and acid.

10.5 Maintenance (by the user)

Note: Before cleaning or maintenance starts the support device must be raised to the highest level. Please switch off the control unit by turning the key switch to the OFF position and remove the key in order to prevent any unexpected interference by others.

Caution: Opening of the equipment and/or attempts to repair by unauthorized persons can adversely affect the safety of the installation and is therefore not permitted.

10.6 Daily maintenance

The daily maintenance of the support device is limited: only cleaning of the walking surface of the support device with a cloth or soft broom is required. Avoid any use of detergents. Never use acid cleaning agents in the hydrotherapy section, because even the vapour of such cleaning agent might corrode the stainless steel parts of the support device and/or other parts.

Cleaning agents should also never be allowed to come into contact with the pool water.

If the stainless steel of the support device shows signs of corrosion, the following actions should be taken immediately:

- a. Check the chloride (salt) content of the pool water. This must not exceed 100 mg/l.
- b. Check the chlorine and pH level of the water.
- c. Check if iron particles are in the pool water.
- d. Try to remove the corroded spots. Please see chapter 4.6.
- e. Inform your local authorized distributor.


10.7 Weekly maintenance

The support device works on water hydraulic pressure. This compression unit is installed in the technical area. The level of the water in the reservoir of the compression unit should be checked every week. Normally the technical department does this, as they are also in charge of the maintenance of the filtration system, which is installed in the same area.

10.8 Monthly maintenance

EWAC recommends a monthly check-up for the presence of the stainless steel clips, which are located between the polyester tiles of the support device. See drawing 3.

In the mains water supply to the compression a special 5-micron filter has been installed to prevent the appearance of iron particles inside the water hydraulic system. This filter should be kept in a good condition and must be checked and cleaned. However, if no pollution is shown every month the time for checking can be set to a longer period. In most cases the technical department will check this filter as well.

Picture 3: Attachment of the clips

10.9 Annual maintenance

Once per year the support device should be cleaned entirely. Furthermore all moving parts must be checked for wear. Also the water reservoir of the compression unit must be cleaned.

Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation
For this annual maintenance it is necessary to empty the pool and all tiles should be removed from the stainless steel.

Warning: Before starting the annual maintenance the special service supports must be installed under the support device in order to avoid accidents.

Persons and/or companies who are authorized by EWAC B.V. should only perform the annual maintenance.

11. Trouble shooting

All abnormal performances or problems that occur can indicate that the support device is not working properly anymore. The support device has to be stopped immediately. Before putting the system back into operation, the failure has to be corrected. The operator himself can solve some of the problems easily. These are listed below.

Problem	Possible cause	Solution
Support device does not sink or raise.	Key switch is not turned on.	Switch on the key switch.
Green indicator light does not light.		Have the power supply check by authorized personnel.
Support device is moving irregularly.	Hydraulic system contains air.	Check the water level in the reservoir of the compression unit.

If it does not work out to solve the above-mentioned problems, please contact your local distributor or an authorized service engineer.

If small corrossions on the stainless steel are visible, the best way of removing them is to do the following: Brush the corroded spot using a simple nailbrush or dish-mop with a normal household cleaning polish. The brush must be moved into the direction of the grinding of the stainless steel surface. Never make turns or use aggressive detergents.

If the corrosion cannot be removed according to the above mentioned, it is possible to use a scotch-brite sponge to be used into the direction of the grinding of the stainless steel.

If there are any doubts how to clean the stainless steel, please contact your local dealer.

12. Service and maintenance

It is necessary to perform an annual inspection to control the good working and safety of the unit. We advise to keep your equipment under service-maintenance by your local qualified service office. Further we recommend filing all activities with respect to carried out maintenance. It is not allowed for unqualified personnel to open the equipment for any reason whatsoever. Qualified personnel can request for a service manual at the local

Extract from the Technical Documents. Support Device with Water Level Adjustment for Swimming Pool Rehabilitation

dealer, this to perform a good and adequate maintenance to the unit. EWAC Medical can not be held responsible for any discomfort to the operator or to the patient due to careless maintenance.

Any claims regarding unit operation please send to the authorized representative of EWAC Medical: OOO Oktopus, 1, Yeniseyskaya Street, bldg 1, office 325, Moscow, 129344, Russia.

According to approved safety standards, always use original components if, for any reason, the unit has to be repaired.

13. Transport packing, package contents, marking

Package contents of the Support Device with Water Level Adjustment for Swimming Pool Rehabilitation with accessories include: height-adjustable floor; hydraulic system; control box; hydraulic system; step-down transformer; connecting cable; User Manual; and the following accessories: easily-removable bars; easily-removable wall bars; balance board; rings bar; obstacle course; stairs, depth indicator.

The support device with accessories is packed into polyethylene film and wooden boxes.

Marking of this medical device includes the following information:

- address of the manufacturer;
- manufacture date;
- device reference number according to the manufacturer's numeration system;
- name of the medical device;
- CE certificate marking.

Transportation marking includes such handling instructions as "Fragile", "Be Cautious", "Avoid water", "Avoid sunrays".

Dulleweg 18, 1721 CT Langedijk, The Netherlands
www.ewacmedical.com

Adjustable Floor

S / N : 4431

Year of manufacture : 2014

14. Disposal

The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation must be disposed of by local public law institutions according to the local requirements and the local disposal regulations.

EWAC Medical

15. List of internal regulations / standards the medical device complies with

EN 980:2008. Medical devices. Symbols for use in the labeling of medical devices.

ISO 10993-1:2009. Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing.

ISO 10993-5:2009. Medical Devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity.

ISO 10993-10:2010. Medical Devices. Biological evaluation of medical devices. Part 10. Tests for irritation and skin sensitization.

ISO 13485:2003. Medical devices. Quality management systems. Requirements for regulatory purposes.

ISO 15223-1:2012. Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied. Part 1. General requirements.

IEC 60601-1. Medical electrical equipment. Part 1. General requirements for basic safety and essential performance.

IEC 60601-1-2. Medical electrical equipment. Part 1-2. General requirements for basic safety and main operational characteristics. Collateral standard. Electromagnetic compatibility. Requirements and tests.

16. Service life

Service life of the device is 10 years since the manufacture date provided it is operated according to the Manual.

17. Use of risk management for medical devices

The Support Device with Water Level Adjustment for Swimming Pool Rehabilitation has been studied in order to analyze any risks associated with such devices. This analysis was scheduled for the initial stage of product development and the report was ready by the beginning of production of the first lot. In case the risks and technologies associated with the devices change, the analysis will be held simultaneously. The main methods and procedures were chosen according to ISO 14971:2007. This standard specifies the process of defining hazards associated with medical devices.

18. Requirements to environment protection

The support device when operated, transported or shipped has no negative effect on humans or environment.

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