State of Minnesota

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SECRETARY OF STATE

APOSTILLE

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1. Country: United States of America

This Public Document: Instruction for Use

2. has been signed by Irma Vargas

3. acting in the capacity of Notary Public, State of Minnesota

4. bears the seal/stamp of Irma Vargas, Notary Public, State of Minnesota

CERTIFIED

5. at: St. Paul, Minnesota

6. Date: February 17, 2016

7. by: Secretary of State, State of Minnesota

8. File No: 8734034-11

9. Seal/Stamp:

10. Signature:

Steve Simon

Steve Simon
Secretary of State

To Whom It May Concern:

The attached is a copy of the Instruction for Use (IFU) for 3M™ Littmann® Stethoscopes.


Scott S. Sardeson
International Regulatory Affairs and
Quality Compliance Leader
3M Health Care

Feb 16, 2016
Date

State of Minnesota)
County of Ramsey)SS

Subscribed and sworn to before me this 16th day of February 2016.

Irma Vargas, Notary Public
My Commission Expires: January 31, 2020

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2016

Littmann® Stethoscopes (warranty period):

Littmann® Classic II Infant (3 years), Littmann® Master Classic II (3 years), Littmann® Select (2 years), Littmann® Lightweight II S.E (2 years), Littmann® Master Cardiology (7 years), Littmann® Cardiology III (5 years), Littmann® Cardiology S.T.C (6 years), Littmann® Classic II Pediatric (3 years), Littmann® Classic III (5 years)

INFORMATION FOR 3M LITTMANN STETHOSCOPES USERS:

3M™ Littmann™ Stethoscopes have become a recognized standard in medical industry because of their high sensitivity, durability and usability. Having regard to their unique features, including variable diaphragm, "warm" rim, comfortable headset and durable dual eartubes, it is not difficult to understand why professional doctors all over the world prefer working with stethoscopes of this brand.

Headset Positioning

Littmann Stethoscopes are equipped with comfortable and efficient headset. Eartubes are oriented in accordance with the anatomy of the ear - the eartubes are permanently set at an angle, in the ear canal direction

- To reduce spring tension in the headset, hold each eartube at the bend near the eartip and gradually pull apart to the desired width
- To increase spring tension, grasp the headset with one hand where the eartubes enter the plastic tubing and squeeze repeatedly

The stethoscope kit also includes soft-sealing eartips for eartubes to improve the acoustic properties.

Eartips

The Littmann stethoscopes use new mechanism of eartips installation to eartubes, where the eartips are tipped and latched at the eartube end. For eartip cleaning or replacing, pull it off the eartube while rotating.

For your safety and stethoscope efficiency, make sure that eartips are tightly fixed to eartubes.

Spare eartips for Littmann stethoscopes are available in 3M representative offices.

Chestpieces

Littmann stethoscopes are equipped with a combined two-way chestpieces with an open bell at the one side and floating diaphragm at the other side. To switch the mode hold the chestpiece by the stem with one hand, and with the other hand turn it until it clicks. Stethoscopes of MC and CIISE models are also available in versions with unilateral chestpieces with adjustable diaphragm.

Low-frequency tones auscultation:

Low-frequency tones auscultation is available with a conventional open bell or a floating diaphragm - in both cases it is necessary to slightly press the chestpiece to the patient's chest. The diaphragm will transmit all tones in perceived diapason, but low frequency tones will dominate and suppress high-frequency tones.

High-frequency tones auscultation:

For high-frequency tones auscultation, you can use a normal floating diaphragm and an adjustable diaphragm, strongly pressing the chestpiece. The adjustable diaphragm is moved inside and is in

contact with an absorbent ring positioned underneath. Thus, at pressing the chestpiece the low frequency amplitude will drop significantly, enhancing the high-frequency tones legibility.

Rim

To clean the stethoscope diaphragm it is usually enough to wipe it with alcohol or wash with soap and water, without removing the rim.

If, however, you need to remove the rim and the diaphragm rim from the chestpiece, it is necessary to operate as follows:

Chestpiece disassembling: Turn the chestpiece with diaphragm upward, push with thumbs the rim at one point and remove it from the chestpiece.

Chestpiece assembling: Place the diaphragm on the chestpiece surface. Insert the chestpiece with the diaphragm in the rim groove, press them with thumbs to the ring at one point and, starting from this point, push ring over the entire circumference. After assembling, check that the diaphragm is completely seated into the rim groove.

STETHOSCOPE CARE

- Keep your stethoscope away from high and low temperatures; avoid stethoscope contact with solvents and oils.
- Eartubes and plastic tube should be wiped with alcohol or soap and water.
- Eartips can be removed from the eartubes for thorough cleaning.
- Plastic tube or other non-metallic parts lifetime can be improved by using special protective solutions (ex. Armor All®).

Note: Do not immerse your stethoscope in any liquid or subject it to steam sterilization process. To sterilize stethoscope use the gas sterilizer 3M™ Steri-Vac™ in cold mode, followed by neutralization for 36 hours.

IMPORTANT

Littmann® stethoscopes tubing is made of polyvinyl chloride (PVC). PVC becomes hard and may crack when exposed to skin for long periods of time. If you like to wear your stethoscope around your neck, wear it over a collar to keep the tubing flexible.

Littmann® stethoscopes tubing damage by cracks formation when exposed to skin is not covered by warranty.

Stored and transport at room temperature.

UTILIZATION

Dispose of contents/container in accordance with applicable local/regional/national/international regulations.

STORAGE AND TRANSPORT

At room temperature

Manufacturer: 3M Health Care, USA

Authorized representative of the manufacturer in the Russian Federation: ZAO «3M Russia», 121614, Moscow, Krylatskaya str. bldg. 17/3
Tel.: +7 495 784 7474

SYMBOLS

 CE Mark, conformity to European Directives

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