
Certification for Information in Dossier

The submitted dossier for the Russian application contains the following:


- **Instruction for Use for Polaris 5.5, 6.35 and Array Instruments and accessories for spinal surgery**

I certify that the documents which are submitted are true and accurate to the best of my knowledge. In addition, Zimmer Biomet Spine, Inc. takes responsibility for any legal obligations related to this submission.

Brendon Heyerdahl
Regulatory Affairs Specialist
Zimmer Biomet Spine, Inc.

Subscribed and sworn before me on this 20 day of December, 2019, in the County of Jefferson, State of Colorado.

Notary Public, State of Colorado

Zimmer Biomet instruments are manufactured from high quality stainless steel and other quality materials that conform to ASTM specifications during F899 Stainless Steel Billet, Bar, and Wire for Surgical Instruments. Instruments are provided non-sterile and should be sterilized prior to use.

CLEANING INSTRUCTIONS

- All instrument components should be disassembled and rinsed with tap water.
- Prepare solution of enzymatic surgical detergent and tap water by adding 2 oz. Of Enzol (Enzymatic Detergent, Johnson & Johnson) to 1 gallon of warm tap water (72°/72° C to 109°/43° C). Instrument components should be immersed for 5 minutes.
- Scrub components with soft brushes and rinse with a brisk stream of tap water (72°/72° C to 109°/43° C) until all visual soil is removed.
- Using an ultrasonic cleaner, sonicate all components in the same enzymatic solution (see step 2) for 10 minutes.
- Manually clean all components with soft brushes and repeat step 3 until all visual soil is removed. If more soil should be repeated until there is no gross soil or residue remains.
- Dry components with a soft gauze cloth.
- Inspect cleaned instrument for loose screws and pins, clamp alignment, cracks and other mechanical irregularities.
- Lubricate instruments to protect instruments during sterilization and storage. This should be done with a water soluble, preservative lubricant after each cleaning. The lubricant should contain a chemical preservative to prevent bacterial growth and be made with distilled water. Excess lubricant should be wiped off prior to storage and sterilization.

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- All torque wrenches should be re-calibrated every 6 months.
- Please refer to ASTM Standard F1744, "Standard Guide for Care and Handling of Stainless Steel Surgical Instruments" for additional information.
- Prior to use, instruments should be visually inspected and function should be tested to assure instruments are functioning properly. If instruments are discolored, have loose screws/pins, are out of alignment, are cracked or have other irregularities, Do Not Use.

STERILIZATION

All packaging materials must be removed prior to sterilization. All instrument components should be sterilized in a hospital setting and that components may move freely. The following steam sterilization parameters are recommended:

Cycle:	Pre-vacuum Steam
Temperature:	270°/132° C
Time:	4 minutes
Drying Time:	20 minutes
Note:	Allow for cooling

Individuals or hospitals not using the recommended method, temperature, and time are advised to utilize any alternative method or cycles using an approved method or standard.

INFORMATION

Please see cover for Customer Service contact information.

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
Zimmer GmbH
Sulzerallee 8
CH-8404 Winterthur
Switzerland
+41 058.854.80.00

Printed in U.S.A.
© 2017 Zimmer Biomet Spine, Inc. All rights reserved.
RIN 060195.C04 Rev 03 2017-08

EN - Manufacturer
CS - Výrobce
DA - Producent
DE - Hersteller
ET - Tootja
EL - Κατασκευαστής
ES - Fabricante
FR - Fabricant
IT - Produttore
LV - Ražotājs
LT - Gamintojas
HU - Gyártó
NL - Fabrikant
NO - Produsent
PL - Producent
PT - Fabricante
RO - Producător
SK - Výrobca
SL - Proizvajalec
FI - Valmistaja
SV - Tillverkare
TR - Üretici
ZH - 製造商
KO - 제조업체
RU - Производители

EN - Authorized Representative in the European Community
CS - Autorizovaný zástupce pro Evropskou unii
DA - Autoriseret repræsentant i EU
DE - Autorisierter Vertreter in der Europäischen Gemeinschaft
ET - Valitud asutaja Euroopa Ühenduses
EL - Εξουσιοδοτημένος εκπρόσωπος της Ευρωπαϊκής Κοινότητας
ES - Representante autorizado en la Comunidad Europea
FR - Représentant autorisé dans la Communauté européenne
IT - Rappresentante autorizzato per la Comunità Europea
LV - Pārvietotais pārstāvis Eiropas Kopienā
LT - Įgalintas atstovas Europos Bendrijai
HU - Hivatalos képviselő az Európai Közösségekben
NL - Gedelegeerdde vertegenwoordiger in de Europese Gemeenschap
NO - Autorisert europeisk representant
PL - Autoryzowany przedstawiciel w Unii Europejskiej
PT - Representante autorizado na Comunidade Europeia
RO - Reprezentant autorizat în Comunitatea europeană
SK - Autorizovaný zástupca pre európske spoločenstvo
SL - Pooblaščen predstavnik v Evropski skupnosti
FI - Valtuutettu edustaja Euroopan yhteisössä
SV - Autoriserad representant inom Europeiska gemenskapen
TR - Avrupa Topluluğu Yetkili Temsilcisi
ZH - 授權經銷代表
KO - 인가된 공식 영업점
RU - Уполномоченный представитель в ЕС

EN - Date of manufacture
CS - Datum výroby
DA - Fremstillingsdato
DE - Herstellungsdatum
ET - Valmistamise kuupäev
EL - Ημερομηνία κατασκευής
ES - Fecha de fabricación
FR - Date de fabrication
IT - Data di fabbricazione
LV - Ražošanas datums
LT - Pagaminimo data
HU - Gyártási dátum
NL - Productiedatum
NO - Produksjonsdato
PL - Data produkcji
PT - Data de fabrico
RO - Data fabricației
SK - Datum výroby
SL - Datum izdelave
FI - Valmistuspäivämäärä
SV - Tillveringsdatum
TR - Üretim tarihi
ZH - 製造日期
KO - 제조일
RU - Дата производства

EN - Use by date
CS - Datum expirace
DA - Udløbsdato
DE - Verwendbar bis
ET - Kasutada enne
EL - Ημερομηνία λήξης
ES - Fecha de caducidad
FR - Utiliser avant le
IT - Usare entro
LV - Izmantoj līdz
LT - Sunaudoti iki
HU - Felhasználható
NL - Uiterste gebruiksdatum
NO - Bruk innen
PL - Stosować przed upływem
PT - Prazo de validade
RO - Data expirării
SK - Spotrebovať do dátumu
SL - Uporabiti do datuma na embalaži
FI - Viimeinen käyttöpäivä
SV - Sista förbrukningsdag
TR - Son kullanma tarihi
ZH - 使用期限
KO - 유효기간
RU - Дата окончания срока

EN - Catalog number
CS - Katalogové číslo
DA - Katalognummer
DE - Katalognummer
ET - Kataloogi number
EL - Αριθμός καταλόγου
ES - Número de catálogo
FR - Numéro de catalogue
IT - Numero di catalogo
LV - Kataloga numurs
LT - Katalogo numeris
HU - Katalógus szám
NL - Catalogusnummer
NO - Katalognummer
PL - Numer katalogowy
PT - Número de catálogo
RO - Numar catalog
SK - Katalógové číslo
SL - Katalogski številka
FI - Tuotenumero
SV - Katalognummer
TR - Katalog numarası
ZH - 目錄號
KO - 카탈로그 번호
RU - Номер каталога

EN - Rx only
CS - Pouze na lékařský předpis
DA - Kun på recept
DE - Verschreibungspflichtig
ET - Saadaval ainult retseptiga
EL - Μόνο με ιατρική συνταγή
ES - Exclusivamente bajo receta médica
FR - Sur ordonnance uniquement
IT - Solo su prescrizione
LV - Tik ar recepti
LT - Tik pagal receptą
HU - Kizárólag rendelvényre kapható
NL - Alleen op recept
NO - Kun på resept
PL - Wyłącznie na receptę
PT - Apenas mediante receita médica
RO - Numai pe bază de rețetă
SK - Len na lékařský předpis
SL - Samo na recept
FI - Vain lääkärin määräyksellä
SV - Receptbelagt
TR - Sadece reçete ile satılır
ZH - 仅限处方药
KO - 의사처방전만 사용 가능
RU - Только по рецепту

EN - Batch code
CS - Označení šarže
DA - Batchkode
DE - Chargenbezeichnung
ET - Partii kood
EL - Καθίσκος παρτίδας
ES - Número de lote
FR - Numéro de lot
IT - Codice del lotto
LV - Partijas kods
LT - Partijos kodas
HU - Tételkód
NL - Particode
NO - Seriekode
PL - Kod serii
PT - Número de lote
RO - Cod lot
SK - Číslo šarže
SL - Serijska številka
FI - Erän tunnus
SV - Batchkod
TR - Parti kodu
ZH - 批次代码
KO - 배치 코드
RU - Код партии

EN - Do not use if package is damaged
CS - V případě poškození obalu, výrobek nepoužívejte
DA - Må ikke benyttes, hvis emballagen er beskadiget
DE - Nicht bei beschädigter Verpackung
ET - Ära kasutada, kui pakend on kahjustatud
EL - Μην χρησιμοποιείτε αν η συσκευασία έχει υποστεί ζημιά
ES - No utilizar si el envase está dañado
FR - Ne pas utiliser si l'emballage est endommagé
IT - Non usare se la confezione è danneggiata
LV - Nelietojiet, ja iepakojums ir bojāts
LT - Nenaudokite, jei pažeistas pakuotė
HU - Ha a csomagolás sérült, ne használja fel
NL - Niet gebruiken indien verpakking is beschadigd
NO - Må ikke brukes hvis forpakningen er skadet
PL - Nie używać, jeżeli opakowanie jest uszkodzone
PT - Não usar se a embalagem estiver danificada
RO - Nu utilizați dacă ambalajul este deteriorat
SK - Nepoužívať, ak je poškodený obal
SL - Ne uporabljajte, če je embalaža poškodovana
FI - Älä käyttää mikäli pakkaus on vaurioitunut
SV - Använd inte om förpackningen är skadad
TR - Ambalaj hasarlıysa kullanmayın
ZH - 包裝受損切勿使用
KO - 포장에 손상이 경우에는 사용하지 마십시오
RU - Не использовать, если упаковка повреждена

STERILE

EN - Sterilized using ethylene oxide
CS - Sterilizováno ethylenoxidem
DA - Steriliseret med ethylenoxid
DE - Mit Ethylenoxid sterilisiert
ET - Steriliseend etüleenoksidi
EL - Αποστειρωμένο με οξείδιο του αιθυλενίου
ES - Esterilizado utilizando óxido de etileno
FR - Stérilisé à l'oxyde d'éthylène
IT - Metodo di sterilizzazione con ossido di etilene
LV - Sterilizēts, izmantojot etilēna oksīdu
LT - Sterilizuota etileno oksidu
HU - Etilén-oxidál sterilizálva
NL - Gesteriliseerd door het gebruik van ethylenoxide
NO - Sterilisert ved hjelp av etylenoksyd
PL - Produkt wykalony jonkiem etylenu
PT - Esterilizado com óxido de etileno
RO - Sterilizat cu oxid de etilenă
SK - Sterilizované s použitím etylénoxidu
SL - Sterilizirano s etilen oksidom
FI - Steriloitu käyttämällä etylenioksidia
SV - Steriliserad med etylenoxid
TR - Etilen oksit kullanılarak sterilize edilmiştir
ZH - 环氧乙烷灭菌
KO - 산화이틸렌으로 멸균처리
RU - Стерилизуют с использованием этиленоксида

EN - Consult instructions for use
CS - Přečtěte si návod k použití
DA - Læs før du bruger
DE - Gebrauchsanleitung beachten
ET - Vi kasutusjuhendit
EL - Συμβουλευτείτε τις οδηγίες χρήσης
ES - Consulte las instrucciones de uso
FR - Consultez le mode d'emploi
IT - Consultare le istruzioni per l'uso
LV - Skatīt lietošanas instrukcijas
LT - Skaitlykite naudojimo instrukcijas
HU - Olvassa el a használati utasítást
NL - Lees voor het gebruik de gebruiksaanwijzing
NO - Kontroller instruksjonene for bruk
PL - Zapoznaj się z instrukcją użycia
PT - Consultar as instruções de utilização
RO - Consultați instrucțiunile de utilizare
SK - Pozrite si návod na použitie
SL - Uporabite navodila za uporabo
FI - Katso käyttöohjeesta
SV - Se bruksanvisningen
TR - Kullanma talimatına bakınız
ZH - 请参考使用说明
KO - 사용 설명서 참조
RU - Обратитесь к инструкциям по использованию

EN - Non-sterile
CS - Nesterilní
DA - Usteril
DE - Nicht steril
ET - Mittestereinne
EL - Μη στείρο
ES - Sin esterilizar
FR - Non stérile
IT - Non sterile
LV - Nesteriāls
LT - Nesterilus
HU - Nem steril
NL - Nietsteriel
NO - Usteril
PL - Wyrób nisteryczny
PT - Não esterilizado
RO - Nesteril
SK - Nesterilné
SL - Nesterilna
FI - Epästerilinen
SV - Icke-steril
TR - Steril Değildir
ZH - 非无菌
KO - 비살균
RU - Нестерильно

EN - Caution, consult accompanying documents
CS - Pozor! Přečtěte si přírodní dokumenty
DA - Forsigtig, læs før du medfølgende dokumenter
DE - Vorsicht, Begleithunterlagen beachten
ET - Ettevaatust, lugeda kaasnevat dokumente
EL - Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
ES - Precaución, consulte los documentos adjuntos
FR - Attention, consultez la documentation jointe
IT - Attenzione, consultare la documentazione allegata
LV - Uzmanību, uzmanī skatīt pievienotos dokumentus
LT - Dėmesio, skaitlykite pridėtus dokumentus
HU - Figyelem, olvassa el a mellékelt dokumentumokat
NL - Waarschuwing, lees de bijbehorende documenten
NO - Forsiktig, se i medfølgende dokumentasjon
PL - Przestroga, należy zapoznać się z dołączoną dokumentacją
PT - Cuidado, consultar os documentos em anexo
RO - Atenție, consultați documentația însoțitoare
SK - Výstraha, pozrite si sprievodnú dokumentáciu
SL - Pozor, uporabite priložene dokumente
FI - Varoitus, katso mukana olevista asiakirjoista
SV - Vaktigt, se medföljande dokumentation
TR - Dikkat, birlikte gelen dokümanlara bakınız
ZH - 注意, 参阅附件文件
KO - 주의, 첨부 문서 참조
RU - Осторожно, обратитесь к сопроводительным документам

EN - Do not reuse
CS - Nepoužívejte opakovaně
DA - Må ikke genbruges
DE - Nicht zur Wiederverwendung
ET - Mitte korduvkasutada
EL - Μην επαναχρησιμοποιείτε
ES - No reutilizar
FR - Ne pas réutiliser
IT - Non riutilizzare
LV - Neļietojiet atkārtoti
LT - Nenaudoti pakartotinai
HU - Ne használja fel újra
NL - Niet hergebruiken
NO - Kun til engangsbruk
PL - Nie używać ponownie
PT - Não reutilizar
RO - Nu reutilizați
SK - Nepoužívať opakovane
SL - Ponovna uporaba ni dovoljena
FI - Älä käyttää uudelleen
SV - Får inte återanvändas
TR - Bir kereden fazla kullanılmayın
ZH - 切勿重复使用
KO - 재사용하지 마십시오
RU - Не использовать повторно

EN - Do not re-sterilize
CS - Výrobek znovu nesterilizujte
DA - Må ikke resteriliseres
DE - Nicht erneut sterilisieren
ET - Mitte steriiseerida korduvalt
EL - Μην επαναποστείρωστε
ES - No reesterilizar
FR - Ne pas restériliser
IT - Non riesterilizzare
LV - Nesterilizēt atkārtoti
LT - Nesterilizuoti pakartotinai
HU - Ne sterilizálja újra
NL - Niet opnieuw steriliseren
NO - Må ikke resteriliseres
PL - Nie wyjątkować ponownie
PT - Não voltar a esterilizar
RO - Nu reesterilizați
SK - Nesterilizujte opakovane
SL - Ponovna sterilizacija ni dovoljena
FI - Älä steriloitaa uudelleen
SV - Får inte omsteriliseras
TR - Tekrar sterilize edilmeyin
ZH - 切勿重复灭菌
KO - 재멸균 처리하지 마십시오
RU - Не переиспользуйте

(English)

NON-STERILE INSTRUMENTS

INTRODUCTION

Zimmer Biomet Spine provides instruments and instrument trays that are generally manufactured from stainless steel, aluminum and polymeric materials. Zimmer Biomet Spine instruments and instrument trays are not intended for use with other manufacturer's instruments, implants or instrument trays. Instruments should be placed in the intended location as indicated on the instrument tray.

Reusable surgical instruments and instrument trays must be cleaned and sterilized prior to use. The following steps apply to the cleaning of instruments and instrument trays.

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CLEANING INSTRUCTIONS

1. All instrument components should be disassembled and rinsed with tap water. The inner shaft of the inserter and trial holders should be disassembled prior to cleaning.
2. Prepare solution of enzymatic surgical detergent and tap water by adding 2 oz./59 ml of Enzol (Enzymatic Detergent, Johnson & Johnson) to 1 gal./3.8 L of warm tap water (72°F/22°C to 109°F/43°C). Instrument components should be immersed in solution for 5 minutes.
3. Scrub components with soft brushes and rinse with a brisk stream of tap water (72°F/22°C to 109°F/43°C) until all visual soil is removed.
4. Using an ultrasonic cleaner, sonicate all components in the same enzymatic solution (see step # 2) for 10 minutes.
5. Manually clean all components with soft brushes and repeat step # 3 until all visual soil is removed. This step should be repeated until there are no signs of soil or residue remaining. Instruments and trays should be visually clean prior to further processing.
6. Dry components with a soft gauze cloth.
7. Inspect cleaned instrument for wear, loose screws and pins, damp alignment, cracks and other irregularities. If instruments are discolored, show evidence of corrosion, have loose screws/pins, are out of alignment, are cracked or have other irregularities, Do Not Use.

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U.S. Sterilization Parameters:	
Cycle:	Pre-vacuum Steam
Temperature:	270°F/132°C
Time:	4 minutes
Drying Time:	20 minutes
Note:	Allow for cooling

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Sterilization Parameters for Use Outside of the U.S.:

Cycle:	Pre-vacuum Steam
Temperature:	275°F/135°C
Time:	3 minutes
Drying Time:	20 minutes
Note:	Allow for cooling

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
Zimmer Biomet Non-Sterile Instruments

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Zimmer Biomet Spine, Inc.
10225 Westmoor Dr.
Westminster, CO 80021 USA
+1 800.447.3625
zimmerbiomet.com


Zimmer GmbH
Sutzelallee 8
CH-8404 Winterthur
Switzerland
+41 056 854.80.00

Printed in U.S.A.
© 2017 Zimmer Biomet Spine, Inc. All rights reserved.
PIN 060195.C05 Rev 05 2017-08

EN - Manufacturer
CS - Výrobce
DA - Producent
DE - Hersteller
ET - Tootja
EL - Κατασκευαστής
ES - Fabricante
FR - Fabricant
IT - Produttore
LV - Ražotājs
LT - Gamintojas
HU - Gyártó
NL - Fabrikant
NO - Produsent
PL - Producent
PT - Fabricante
RO - Producător
SK - Výrobca
SL - Proizvajalec
FI - Valmistaja
SV - Tillverkare
TR - Üretici
ZH - 制造商
KO - 제조업체
RU - Производитель

EN - Authorized Representative in the European Community
CS - Autorizovaný zástupce pro Evropskou unii
DA - Autoriseret repræsentant i EU
DE - Autorisierter Vertreter in der Europäischen Gemeinschaft
ET - Volitatud esindaja Euroopa Ühenduses
EL - Εξουσιοδοτημένος εκπρόσωπος στην Ευρωπαϊκή Κοινότητα
ES - Representante autorizado en la Comunidad Europea
FR - Représentant autorisé dans la Communauté européenne
IT - Rappresentante autorizzato per la Comunità Europea
LV - Pārvietotais pārstāvis Eiropas Kopienā
LT - Įgaliojasis atstovas Europos Bendrijoje
HU - Hivatalos képviselő az Európai Közösségekben
NL - Gedelegeerdde vertegenwoordiger in de Europese Gemeenschap
NO - Autorisert europeisk representant
PL - Autoryzowany przedstawiciel w Unii Europejskiej
PT - Representante autorizado na Comunidade Europeia
RO - Reprezentant autorizat în Comunitatea europeană
SK - Autorizovaný zástupca pre Európske spoločenstvo
SL - Pooblaščen predstavnik v Evropski skupnosti
FI - Valtuutettu edustaja Euroopan yhteisössä
SV - Auktoriserad representant inom Europeiska gemenskapen
TR - Avrupa Topluluğu Yetkili Temsilcisi
ZH - 歐羅巴區代表
KO - 유럽 공동체 공식 영업점
RU - Уполномоченный представитель в ЕС

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ATTENTION OPERATING SURGEON

DESCRIPTION

Zimmer Biomet Spine provides instruments and instrument trays that are generally manufactured from nitinol, titanium alloy, stainless steel, aluminum and polymeric materials, e.g., PEEK, Radel and silicone. Zimmer Biomet Spine instruments and instrument trays are not intended for use with other manufacturer's instruments, implants or instrument trays. Instruments should be placed in the intended location as indicated on the instrument tray.

Reusable surgical instruments and instrument trays must be cleaned and sterilized prior to use. The following steps apply to the cleaning of instruments and instrument trays:

PRECLEANING AT THE POINT OF USE

1. Soil should be wiped from the device with a sponge or towel moistened with water.
2. Cannulated instruments should be flushed to prevent drying of debris inside.
3. To prevent blood and/or debris from drying, devices should be placed in a container and covered with a towel that has been moistened with water.

CLEANING INSTRUCTIONS

1. All instrument components should be disassembled and rinsed with tap water. The inner shaft of the instruments should be disassembled prior to cleaning.
2. Prepare solution of enzymatic surgical detergent and tap water by adding 2 oz./59ml of Enzoi (Enzymatic Detergent, Johnson & Johnson) to 1 gal./3.8L of warm tap water (72°F/22°C to 109°F/43°C). Instrument components should be immersed in solution for 5 minutes.
3. Scrub components with soft brushes and rinse with a brisk stream of tap water (72°F/22°C to 109°F/43°C) until all visual soil is removed.
4. Using an ultrasonic cleaner, sonicate all components in the same enzymatic solution (see step #2) for 10 minutes.
5. Manually clean all components with soft brushes and repeat step #3 until all visual soil is removed. This step should be repeated until there are no signs of soil or residue remaining. Instruments and trays should be visually clean prior to further processing.
6. Conduct a final rinse with a brisk stream of tap water (72°F/22°C to 109°F/43°C) for the entire instrument including any inner shafts.
7. Dry components with a soft gauze cloth.
8. Inspect cleaned instrument for wear, loose screws and pins, clamp alignment, cracks and other irregularities. If instruments are discolored, show evidence of corrosion, have loose screws/pins, are out of alignment, are cracked or have other irregularities, Do Not Use.

CARE AND HANDLING

- All torque wrenches should be re-calibrated every 6 months.
- Please refer to ASTM standard F1744, "Standard Guide for Care and Handling of Stainless Steel Surgical Instruments" for additional information.
- Prior to use, instruments should be visually inspected for wear and function should be tested to assure instruments are functioning properly. If instruments are discolored, show evidence of corrosion, have loose screws/pins, are out of alignment, are cracked or have other irregularities, Do Not Use. Instrumentation that appears damaged should be returned to the manufacturer.

STERILIZATION

All packaging materials must be removed prior to sterilization. All instrument components should be sterilized in a loosened state such that components may move freely. The following steam sterilization parameters are recommended:

U.S. Sterilization Parameters:

Cycle:	Pre-vacuum Steam
Temperature:	270°F/132°C
Time:	4 minutes
Drying Time:	30 minutes
Note:	Allow for cooling

FDA cleared sterilization wraps should be used to maintain sterility after processing.

Sterilization Parameters for Use Outside of the U.S.:

Cycle:	Pre-vacuum Steam
Temperature:	275°F/135°C
Time:	3 minutes
Drying Time:	30 minutes
Note:	Allow for cooling

Zimmer Biomet does not recommend stacking of trays during the sterilization process.

Individuals or hospitals not using the recommended method, temperature, and time are advised to validate any alternative methods or cycles using an approved method or standard.

Additional information about instrumentation (including assembly/disassembly) may be found in the system surgical technique, and/or reference guide.

INFORMATION

For further information, please contact the Customer Service Department at:

Zimmer Biomet Spine, Inc.
10225 Westmoor Dr.
Westminster, CO 80021 USA
+1 800.447.3625
zimmerbiomet.com

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© 2010 Zimmer Biomet Spine, Inc. All rights reserved. Zimmer Biomet Spine, Inc. is a registered trademark of Zimmer Biomet Spine, Inc. in the United States and other countries. Zimmer Biomet Spine, Inc. is not responsible for the use of the product in any country where it is not registered.