State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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1. Country (Pays): United States of America

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2. has been signed by (a été signé par) Ryan J. Bartz**3. acting in the capacity of (agissant en qualité de) Notary Public, State of Minnesota****4. bears the seal / stamp of (est revêtu du sceau / timbre de)
Ryan J. Bartz, Notary Public, State of Minnesota****Certified
Attesté****5. at (à) St. Paul, Minnesota****6. Date (Date) 11/26/2019****7. by (par) Secretary of State, State of Minnesota****8. File No (Sous n°) 1120679200076****9. Seal / Stamp (Sceau / Timbre):****10. Signature (Signature):**

Steve Simon
Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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November 25, 2019

To whom it may concern

In my capacity as Senior Regulatory Affairs Associate, 3M ESPE Dental Products,
2510 Conway Avenue, St. Paul, MN, 55144-1000, USA, I certify that attached IFU for:

3M ESPE RelyX™ Try-In Paste

in English language is truthful and accurate

Angie Draper

Angie Draper
Sr. Regulatory Affairs Associate
3M ESPE Dental Products

25 Nov 2019

Date


State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 25th day of November 2019.

R. J. Bartz

Ryan Bartz, Notary Public
My Commission Expires: January 31, 2024

Instructions for Use

3M ESPE RelyX™ Try-In Paste

3M ESPE RelyX™ Try-In Pastes are used to guide the dentist in the selection of the shade required for the final cementation with 3M ESPE RelyX Veneer cement and RelyX Ultimate adhesive resin cement. The shades of RelyX Try-In paste specifically match the final cured shade of RelyX Veneer cement and RelyX Ultimate cement. RelyX Try-In paste contains a polyethylene glycol resin, a zirconia/silica filler and pigments. The RelyX Try-In pastes are water soluble for easy clean up.

Relyx Try-In Paste can be supplied as refills:

1) 1 Relyx Try-In Paste (A1T shade); 2) 1 Relyx Try-In Paste (A3T shade), 3) 1 Relyx Try-In Paste (B0,5T shade); 4) 1 Relyx Try-In Paste (TRT shade); 5) 1 Relyx Try-In Paste (WOT shade).

Indications

- Color matching guide for RelyX Veneer cement and RelyX Ultimate adhesive resin cement.

Instructions for use

Preparation

1. To prepare the veneer, the internal bonding surface of a ceramic veneer surface should be micro-etched with a hydrofluoric acid etchant. The dental laboratory typically does this. The internal surface of a composite veneer should be roughened by sandblasting or air abrasion.
2. Remove the provisional restoration, if present, from the tooth. Clean tooth with pumice; rinse thoroughly and lightly air dry.

Directions

1. Check the fit and esthetics of the restoration by trying in the veneer with water or the Translucent shade of RelyX Try-In paste. If the fit and esthetics are acceptable, proceed directly to the bonding procedure using the Translucent shade of RelyX Veneer cement.
2. If color adjustment is needed, select the appropriate shade of RelyX Try-In paste. Apply RelyX Try-In paste directly onto the internal surface of the veneer. Gently seat the veneer onto the tooth preparation. With multiple veneers, start at the midline and seat veneers sequentially working back. The shade of RelyX Try-In paste is formulated to precisely match the corresponding cured shade of the cement.
3. Once the fit and esthetics are verified, remove the veneers and thoroughly rinse the RelyX

Try-In paste from the tooth and restoration with a water spray and then dry.
The veneer is now ready for final bonding and cementation. Refer to the appropriate product's instruction for use.

Storage and Use

This product is designed to be used at room temperature. The humidity and pressure are not regulated. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.
Do not expose the paste to elevated temperatures.

Technical characteristics:

Relyx Veneer Try-in paste syringe: Length (mm): $117,5 \pm 11$ (filled full assembly), Width (mm): $31,75 \pm 3$; Height (mm): $12,7 \pm 1,6$; Diameter (mm): $8,9 \pm 1$; Mass (g): $6,8 \pm 5\%$ (filled assembly).
Relyx Try-in Paste cap: Length (mm): $12,7 \pm 1,6$, Diameter (mm): $4,8 \pm 1$; Mass (g): $0,24 \pm 0,1$ (cap only).

List of standards: ISO 10993-1: 2003(E), ISO 10993-3: 2003, ISO 10993- 5: 1999, ISO 10993-6: 2007, ISO 10993-10: 2002, ISO 10993-11: 2006, ISO 7405: 1997, ISO 13485: 2003 +Cor, 1:2009, ISO 14971: 2000, ISO 4049: 2009, ISO 15223-1: 2012, EN 980: 2008, EN 1041: 2008/ Amd 1:2013

Labelling information

Products are supplied in packaging labeled with the following information: Product name and REF number; name and address of the manufacturer; Do not reuse; European Authorized Representative sign; Contents; LOT number; Storage conditions; Shelf life; Information about trademarks and symbols - Caution, See accompanying document, Symbol CE; Flammable; Corrosive substance; Ecological packaging, Recycling sign.

Cleaning and Disinfection

Disinfect the products using an intermediate level disinfection process (liquid contact) as recommended

Medicinal substances: Medical devices do not contain or deliver medicinal substances.

Human Blood Derivatives: Medical devices do not contain or deliver human blood derivatives.

Animal-Derived Substances: No materials of human origin are intentionally being added to the products.

Disposal methods: Dispose of content/container in accordance with the local/regional/national/international regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological requirements for the handling of medical waste" to the waste class B.

Maintenance and repair

This medical device isn't subject to maintenance and repair.

Manufacturer of the medical device

– Developer

3M ESPE Dental Products

Address: 2510 Conway Avenue, St, Paul, Minnesota, 55144-1000, USA

Phone: 651-7334365

– Manufacturer

3M ESPE Dental Products

Address: 2510 Conway Avenue, St, Paul, Minnesota, 55144-1000, USA

Phone: 651-7334365

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Information valid as of November 2019

Customer Service

Authorized representative of the manufacturer in the Russian Federation: 
OAO "3M Russia", 108811, Moscow, Poselenie Moskovskiy, Kievskoe shosse, 22 km, 6, bld.1,
Russia /phone +7 (495) 578-74-74

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