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To whom it may concern

Letter of Declaration

Herewith we, Siemens Healthcare Diagnostics Products GmbH, Marburg, Germany, in our responsibility as legal manufacturer would like to confirm that the attached documents are true and correct copies of the original documents.

Documents:

- Instruction for Use Dade Owren's Veronal Buffer
- Addendum to Instruction for Use Dade Owren's Veronal Buffer

Sincerely,

Siemens Healthcare Diagnostics Products GmbH

By:

Name: Andreas Wiegand

Title: Senior Director, Regulatory Affairs

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© 2008 Siemens Healthcare Diagnostics Products GmbH.

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**SUPPLEMENT TO THE INSTRUCTIONS FOR USE
OF MEDICAL DEVICE**

**“Kits and separate packages of in vitro reagents for studying hemostasis, including:
- sample diluent Dade Owren's Veronal Buffer
by Siemens Healthcare Diagnostics Products GmbH,
Germany**

2019

1. Purpose of the medical device

Diluting buffer for coagulation tests. Sample diluent Dade Owren's Veronal Buffer is designed specially for studying hemostasis, but can also be used in any laboratory analysis where isotonic barbitol buffer with pH 7.35 is used.

2. Indications and contraindications for use of medical device, possible side effects associated with the use of medical device

Not applicable for this medical device.

3. Information about possible consumers of medical device, specific requirements for proficiency level of possible users

For Professional Use.

This medical device shall be used only by qualified, trained (in the field of clinical laboratory diagnostics) personnel.

4. Information in accordance with data contained in the State Register of Medicinal Products for Human Use on the medicinal product and (or) pharmaceutical substance contained in the medical device

Not applicable for this medical device.

5. Safety requirements and precautions

Precautionary statements

For in Vitro Diagnostic Use.

To avoid contamination, perform pipetting in view of precautionary statements. Visible turbidity and flocculation are the signs of microbiological contamination.

It is not intended for internal and outward application by people or animals.

For professional use only.

Procedure-related precautions

Before you start, read the instruction in its entirety. Carefully follow the instruction. Failure to follow the instruction may lead to irregular test results.

Do not use if the package is not intact.

Using a reagent to analyze the samples different from those that can be used with this assay may lead to irregular test results.

Do not use reagents beyond the expiration date. Always check the expiration date before testing.

Do not mix the reagents from batches with different numbers.

6. Transport conditions

During the transportation in refrigerator, the product remains stable until the expiration of shelf life indicated on the labeling under the following conditions:

Temperature: from +2°C to +8°C.

Relative humidity: from 10 to 80 % without condensation.

Transportation in the user's premises must be performed only in closed securely reinforced containers in an upright position. Persons responsible for the transportation must be aware of the actions to be taken in case of accidents or leakage.

7. A procedure and conditions for disposal or destruction of the medical device

During the use of medical device, waste is generated, which must be categorized and destroyed (disposed) as an epidemiologically hazardous waste (SanPiN 2.1.7.2790-10: B class).

If possible, the generation of waste shall be avoided or minimized.

When disposing this product, solutions and derivatives, the legislative requirements for environment protection and waste disposal, as well as any requirements of local governing bodies must be complied with. The surplus of product and products not subject for recycling shall be disposed by a contracting organization having a license for waste disposal. Do not release the crude wastewater into sewerage, unless it conflicts with all applicable requirements of competent authorities.

Used package shall be recycled. If recycling is impossible, dispose the product only in the waste incineration plants or at the landfill site.

The remaining product and containers must be disposed in a safe way. Empty containers can have some amount of the remaining product. Avoid the spread of spilled material and its penetration into soil, drains and sewerage.

8. Warranty commitments

The manufacturer guarantees stability of the medical device until an expiration date subject to the conditions of transportation, storage and use indicated in the Instruction.

The manufacturer guarantees the safety of medical device, lack of unacceptable risk of causing harm to the human health and environment during the use of reagent kit as intended under the conditions indicated by the manufacturer.

9. The list of national and international regulatory documents/standards to which a MD corresponds

Process of production and distribution of products and final quality control fully correspond to the requirements set out in Directive 98/79/EC and U.S. national standards of quality and technical safety.

List of applicable standards:

Standard	Description
EN 980	Medical devices. Graphical symbols for use in the labelling of medical devices
EN 13612	Performance evaluation of in vitro diagnostic medical devices
EN 13640	Stability Testing of in vitro diagnostic reagents
EN 13641	Elimination or reduction of risk of infection related to in vitro diagnostic reagents
EN ISO 14971	Medical devices. Application of risk management to medical devices
DIN ISO 15223	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied.
DIN EN ISO 18113	In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling).

10. Reclamations

For product quality issues, please contact:

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File-No. 2272 / 2019

It is officially certified hereby that **Andreas Wiegand** and **Dr. Sven Liebler**, Emil-von-Behring-Straße 76, 35041 Marburg - both personally known to me - have signed this document.

Marburg, November 25, 2019

Dr. Adler, Notary Public
as officially appointed representative of
Dr. Anton Sebastian Schmölz, notary public in Marburg

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