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Statement

ProteomeTech Inc. declares that we issued Maintenance Documentation as attached.

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APPOINTMENT

The reagent kit «Immunochromatographic assay for semi-quantitative detection of the level of total immunoglobulin class G in human serum, plasma or whole blood samples» «ImmuneCheck» for semi-quantitative determination of human immunoglobulin G in human serum, plasma or whole blood in order to determine the immunological status of a person by Immunochromatography method.

Application area

Clinical Laboratory Services.

Medical device for diagnostics in-vitro only.

Potential consumers of products

Health facilities of any type.

The professional level of the potential users

Doctor of clinical laboratory diagnostics, medical laboratory technician.

Possible side effects

Side effects are not revealed.

Contraindications

If the patient has a severe form of haemophilia.

Description of the target analyte, the information about its scientific validity.

Human IgG is most of all immunoglobulins in the blood. These antibodies captured antigens in the body and thus protect the organism. In particular, Fc-fragments of IgG molecules bind to specific receptors on neutrophils and macrophages which effectively destroy hostile bacteria. These actions are important in the antibody- dependent cell-mediated cytotoxicity (ADCC). Deficiency of IgG may be linked to sinus infections, ear infections and other respiratory infections, pneumonia and bronchitis and so forth. Therefore, quantification of IgG in the blood can be used as an important method of assessing the state of human health.

OPTIONS AND ACCESSORIES

The set is available in two basic versions:

- I. ImmuneCheck IgG.
- II. ImmuneCheck IgG-F.

Title	Description	ImmuneCheck IgG	ImmuneCheck IgG-F
The test - cassette	Container of plastic size 70x19 mm, with holes for sample introduction and an analytical zone comprising a test strip size 60x4 mm. Packaged in individual package size 120x80 mm with aluminium foil along with the desiccant. The test strip includes a membrane on which there are: - a control line C: goat anti-mouse IgG; - test line T1: mouse antibody to human	20 pcs.	1 pc.

	IgG; - test line T2: human IgG; - on gold conjugated goat anti-human IgG; - on gold conjugated mouse IgG		
The vial with 1 mL of buffer solution	Plastic bottle (high density polyethylene) with a lid dropper (polypropylene, polyethylene, silicone) 2 mL volume with an adhesive label (polypropylene film, size: 38.1×19.2 mm) containing 1 mL of a buffer solution without foreign particles or contaminants, comprising: - sodium carbonate; - casein; - sodium chloride.	20 pcs.	1 pc.
Capillary tube	Material: Polycarbonate Volume: 5 uL	20 pcs.	1 pc.
Lancet-scarifier	Material: low-density polyethylene Size: diameter of the needle 28G	-	1 pc.
Alcohol wipes	Material: the nonwoven fabric is saturated 70% isopropyl alcohol Size: 30x35 mm (napkin), 35x55 mm (packing)	-	1 pc.
Instructions for use	Material: 75 g/m ² paper with a high degree of whiteness Size: 210x297 mm	1 pc.	1 pc.
Table for determining of the concentration of IgG	Material: 75 g/m ² paper with a high degree of whiteness. Size: 80 x 120 mm	2 pcs.	1 pc.

Notes.

1. The kit includes all the reagents necessary for the production ICA except timer and disposable gloves.
 2. In the composition of the product there are no medicines and pharmaceutical substances.
 3. As part of the product contains materials (lancet) that come into direct contact with the short body of the patient, and there are no materials that come into indirect contact with the patient's body, and personnel using the product (instructions for use) in the performance of operational documentation requirements.
 4. The use of the product has no population and demographic aspects.
- The kit components are packed in a box, in which the enclosed instructions for use, Certificate of Analysis.

BASIC CONSUMER CHARACTERISTICS

The reagent kit ImmuneCheck IgG is designed for the study of 20 samples, and is designed for use in health care facilities.

The reagent kit ImmuneCheck IgG-F is designed for the testing of only 1 sample and is designed for use in health care facilities.

Methods of sterilization of products

The product does not require sterilization.

The software work products

Is absent.

Maintenance and repair products

The product requires no maintenance and cannot be repaired.

OPERATING PRINCIPLE

The test work is based on a qualitative immunochromatographic analysis.

In conventional tests, using a one-step sandwich assay when hIgG is present in an excessive amount, there is an erroneous diagnosis (false negative), no color on the test line. All test strains comprising hIgG-antibodies saturated free immunoglobulin which blocks the formation of a "sandwich complex" type antibody-immunoglobulin antibody ("prozone effect" or high-dose effect (hook effect)).

ImmuneCheck is a quick semi-quantitative diagnostic test which determines the amount of human IgG using two different mechanisms: a sandwich immunoassay (Test Line 1) and competitive immunoassay (Test Line 2). On the test line of human IgG 1 from the sample and antibody to human IgG, conjugated to gold nanoparticles bind with antibodies to IgG, fixed on the test line by a sandwich method, and coloring occurs line. At the test line 2 human IgG in the sample and the human IgG, fixed line, competes with the anti-human IgG, conjugated gold nanoparticles staining line. To determine concentration of human IgG, the color intensity of the test lines 1 and 2 is compared with the intensity of the Control line color (C), which is not affected by the change of human IgG concentration, and it has a constant color intensity. Solving the problem of the hook effect, these products provide quantitative diagnostics in POCT format and applied to a sample of serum, plasma or whole blood.

ANALYTICAL AND DIAGNOSTIC FEATURES

No	Index	meaning
1	Analysis Range / Sensitivity	0.12~4,000 mg / dL or 0.0012~40 mg / mL
2	The imprecision (CV)	Deviations in assay: less than 20% Deviations between analyses: at least 30%
3	Range of accuracy	70~140%
4	Compliance with another method of quantitative analysis	74.2 ~ 132.1 % agreement and correlation coefficient of 0.9243 (R^2) compared with the turbidimetric immunoassay
5	Positive control	Positive standards are prepared according to SOP.
6	Negative control	Negative standards are prepared according to SOP.
7	The time to reach stable results	20 minutes

THE TESTING SAMPLES

Serum, plasma or whole blood of human.

Prior to analysis, storage of serum or plasma samples at a temperature of 2-8°C for no more than 2 weeks; if necessary, longer storage (up to 2-3 months) - at a temperature of -20°C or lower. Blood cells and other solids should be removed from taken after storage at 2-8°C or frozen at a temperature of -20°C and thawed serum or plasma sample. Repeated freezing and thawing of serum or plasma is not allowed.

Store whole venous blood at 2-8°C for no more than 72 hours.

Freezing whole blood is not allowed.

Whole finger capillary blood is tested immediately.

Serum, plasma or whole human blood

1. Preparation of a sample of whole venous blood.

1) Whole blood is drawn into syringes or into drainage tubes that are filled with an anticoagulant (heparin, sodium citrate, or EDTA).

2) Whole blood samples should be tested immediately after collection. In the case of storage at 2~8°C, they should be tested within 72 hours.

2. Preparation of serum or plasma samples

1) Samples of serum or containing anticoagulant (heparin, sodium citrate or EDTA) plasma should be tested immediately after sampling.

2) Do not store the samples at 2~8°C over 14 days. Prolonged storage of samples should be at a temperature below -20°C.

Haemolysed samples of serum, plasma, and whole blood cannot be used.

INTERFERENCE AND CROSS-REACTIONS

Hemoglobin up to a concentration of 5 g / l, triglycerides up to a concentration of 4.5 mmol / l, bilirubin up to a concentration of 360 μ mol / l and ascorbic acid up to a concentration of 60 mg / dl do not affect the correct determination of IgG concentration by the Reagent kit "Immunochromatographic test system for semi-quantitative determination of the concentration of total class G immunoglobulin in human serum, plasma or whole blood «ImmuneCheck».

PRECAUTIONS

- 1) For in vitro diagnostic use only.
- 2) The test results are used to determine the clinical diagnosis, and after evaluating all the clinical and laboratory results, the doctor must make a specific clinical diagnosis or treatment decision.
- 3) Do not smoke, eat or drink in places where samples or test reagents are used.
- 4) Avoid direct hand contact when handling reagents.
- 5) Repeated freezing and thawing of serum / plasma samples can lead to incorrect results.
- 6) All samples may contain unknown infectious materials. When handling samples, wear disposable gloves and wash your hands after testing.
- 7) Do not put needles, knives or any other objects that could result in personal injury when handling human samples and reagents to prevent accidents.
- 8) All patient samples and kit components used should be considered biohazard. They should be disposed of in accordance with the relevant instructions.
- 9) Do not use kits if the cassette is damaged or expired.
- 10) Cartridges, capillary tubes and lancets included in the kits are disposable. Do not reuse.
- 11) Use the cassette immediately after opening, as it is vulnerable to moisture.

TESTING PROCEDURE

Equipment and materials.


Gloves.
Timer.
Refrigerator with freezer.

The samples used.

Whole blood, serum and plasma can be used as samples. Use the following methods of sample preparation:

(1) Capillary whole blood.

- 1) Disinfect the skin for blood collection (e.g., fingertip) with an alcohol swab, allow the skin to dry.
- 2) Pierce the skin area with a lancet-scarifier. For this:
 - rotate the cap along arrow (1) Fig.1
 - remove the protective cap from the lancet-scarifier (2) Fig.1.

(Fig. 1)

- attach the lancet-scarifier to the skin at the future site of the puncture, and then click on the scarifier to activate the device (1) Fig.2.

3) Collect 5 uLo of a sample of capillary blood by filling its capillary tube¹ by contacting the tip of a disposable capillary tube with blood at the puncture site (2) Fig. 2.

NOTE If the blood sample was not completed before the end of the capillary tube, the result may not be correct due to insufficient sample.


4) Unscrew and remove the dropper cap from the vial with the buffer solution and put the capillary tube filled with blood into the vial. Reinstall the removed dropper cap and tighten the thread securely.

5) Gently tilt the bottle from side to side for 10-15 seconds so that the blood sample in the capillary tube mixes well with the buffer solution (3) Fig. 2.

WARNING If the blood is poorly mixed with a buffer solution, the result may be incorrect.

6) Open aluminum bag, remove the test cassette and place on a flat surface.

7) Break off spout in the lid-dropper bottle with diluted sample (4) Fig. 2.

(Fig. 2)

8) Holding the bottle vertically drip 4 drops (about 100 uL) of the diluted sample in the sample well (S) of test cassette, slowly squeezing the vial by fingers (5) Fig. 2.

NOTE If the amount of the sample is too high or too low, the result may not be correct. Drip 4 drops smoothly, vertically holding the vial.

CAUTION Stability test cassette can be reduced if a strip within the cassette is exposed to moisture.

9) Analyze the result after 20 minutes after addition of the sample (6) Fig. 2.

¹ The capillary tube is not used for invasive purposes and is used only for the introduction of the sample (letter from the manufacturer), so sterilization is not required.

WARNING Do not analyze results in 40 minutes or more after the introduction of the sample.

(2) Whole venous blood

- 1) Collect blood from the venous venipuncture into the tube containing anticoagulant (heparin, sodium citrate or EDTA). The collected blood must be used immediately or within three days when stored at 2~8°C.
- 2) Collect 5 uL of sample of venous blood filling capillary tube by contacting the tip with a disposable capillary venous blood sample (1) Fig. 3.

(Fig. 3)

NOTE If the blood sample was not completed before the end of the capillary tube, the result may not be correct due to insufficient sample.

- 3) Unscrew and remove the dropper cap from the vial with the buffer solution and put the capillary tube filled with venous blood into the vial. Reinstall the removed dropper cap and tighten the thread securely.
- 4) Gently incline the vial for 10-15 seconds so that the venous blood sample in the capillary tube mixes well with the buffer solution (2) Fig. 3.

WARNING If the blood is poorly mixed with a buffer solution, the result may be incorrect.

- 5) Open the aluminum package, remove the test cassette and place on a flat surface.
- 6) Break off the spout at the dropper cap of the diluted sample bottle (3) Fig. 3.
- 7) While holding the bottle vertically, drip 4 drops (about 100 uL) of the diluted sample into the sample well (S) of the test cassette, slowly squeezing the bottle with your fingers (4) Fig. 3.

NOTE If the amount of the sample is too high or too low, the result may not be correct. Drip 4 drops smoothly, vertically holding the vial.

CAUTION. Stability test cassette can be reduced if a strip within the cassette is exposed to moisture.

- 8) Analyze the result after 20 minutes after addition of the sample (5) Fig. 3.

WARNING. Do not analyze results in 40 minutes or more after the introduction of the sample.

(3) Plasma

- 1) Collect blood from the venous venipuncture into a tube containing anticoagulant (heparin, sodium citrate or EDTA).
- 2) Separate the plasma from a sample of venous blood taken standard method. The separated plasma can be used for two weeks when stored at 2~8°C or must be stored at -20°C for long term use.

NOTE. Before taking the sample with a capillary tube, blood cells and other solids should be removed from the sample taken after storage at 2~8°C or frozen at -20°C and thawed serum / plasma sample.

Hemolyzed samples of serum, plasma, and whole blood cannot be used. Centrifuge serum / plasma samples at 2,000~3,000 rpm.

Repeated freezing and thawing of serum / plasma samples should be avoided.

3) Collect 5uL of plasma sample filling the capillary tube by contacting the tip with a disposable capillary plasma sample (1) Fig. 3.

NOTE If the plasma sample has not been completed before the end of the capillary tube, the result may not be correct due to insufficient sample.

4) Unscrew and remove the dropper cap from the bottle with buffer solution and put the capillary tube filled with plasma in the bottle. Reinstall the removed dropper cap and tighten the thread securely.

5) Carefully turn the bottle for 10-15 seconds so that the plasma sample in the capillary tube mixes well with the buffer solution (2) Fig. 3.

CAUTION If the plasma is poorly mixed with a buffer solution, the result may be incorrect.

6) Open the aluminum bag, remove and place the test cassette on a flat surface.

7) Break off the spout at the dropper cap of the diluted sample bottle (3) Fig. 3.

8) While holding the bottle vertically, drip 4 drops (about 100 uL) of the diluted sample into the sample well (S) of the test cassette, slowly squeezing the bottle with your fingers (4) Fig. 3.

NOTE If the amount of the sample is too high or too low, the result may not be correct. Drip 4 drops smoothly, vertically holding the vial.

CAUTION Stability test cassette can be reduced if a strip within the cassette is exposed to moisture.

9) Analyze the result after 20 minutes after addition of the sample (5) Fig. 3.

WARNING Do not analyze results in 40 minutes or more after the introduction of the sample.

(4) Serum

1) Collect blood from the venous venipuncture into a tube containing no anticoagulant.

2) Leave the sample at room temperature. After 30 minutes, separate the serum by centrifugation, and collect the supernatant. The separated serum may be used for two weeks when stored at 2~8°C or stored at -20°C for longer storage.

NOTE. Before sampling capillary tube, blood cells and other solids must be removed from taken after storage at 2~8°C or frozen at -20°C and thawed serum / plasma sample.

You cannot use hemolyzed samples of serum, plasma and whole blood. Centrifuge serum / plasma samples for 2,000~3,000 rpm.

Refreezing serum / plasma samples should be avoided.

3) Collect 5 uL of serum sample, filling the capillary tube by contacting the tip with a disposable capillary serum sample (1) Fig. 3.

NOTE If serum sample was not completed before the end of the capillary tube, the result may not be correct due to insufficient sample.

4) Unscrew and remove the dropper cap from the bottle with buffer solution and put the capillary tube filled with plasma in the bottle. Reinstall the removed dropper cap and tighten the thread securely.

5) Carefully turn the bottle for 10-15 seconds so that the plasma sample in the capillary tube mixes well with the buffer solution (2) Fig. 3.

CAUTION If the serum is poorly mixed with a buffer solution, the result may be incorrect.

- 6) Open the aluminum bag, remove and place the test cassette on a flat surface.
- 7) Break off the spout at the dropper cap of the diluted sample bottle (3) Fig. 3.
- 8) While holding the bottle vertically, drip 4 drops (about 100 μ L) of the diluted sample into the sample well (S) of the test cassette, slowly squeezing the bottle with your fingers (4) Fig. 3.

NOTE If the amount of the sample is too high or too low, the result may not be correct. Drip 4 drops smoothly, vertically holding the vial.

CAUTION Stability test cassette can be reduced if a strip within the cassette is exposed to moisture.

- 9) Analyze the result after 20 minutes after addition of the sample (5) Fig. 3.

WARNING Do not analyze results in 40 minutes or more after the introduction of the sample.

5.3. Determination of IgG concentrations

- 1) Determine the concentration of human IgG, visually comparing relative colour intensity of the test line T1 and T2 with the control line (C), which is independent of the concentration of human IgG via <Tables to determine the concentration>.

The colour intensity of the control line designated as 1. Find the closest result in <Table to determine the concentration> and read the corresponding concentration (mg / mL or mg / dL).

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www.biofocus.ru

- (1) Perform the following steps to make a decision on the concentration.
 - a) Determine the concentration range based on the test line T2 colour intensity.
 - b) Check the position in <Table determining the concentration> by comparing the relative ratio of colour intensity of the test line T1 with a control line (C).
 - c) Determine the result of <Table determining the concentration>.
 - d) The corresponding value of the concentration (mg / mL or mg / dL) is the IgG concentration in the sample.

2) Background data

- When it is considered that the actual result is somewhere between the two adjacent strips at <Table solutions of concentrations> read the concentration (mg / mL or mg / dL) among these two bands to that described in the upper part.
- When used for the analysis of whole blood sample, the concentration of IgG in the serum can be calculated by dividing the "concentration (mg / mL or mg / dL)" whole blood to the value (100-hematocrit (HCT)), and multiplying the result by 100

- graphic symbol "Person";
- graphic symbol of "the number of panels per pack"
- graphic symbol «Invitro diagnostic (for in vitro diagnostic)» IVD;
- graphic symbol "Keep at a temperature 1~30°C";
- graphic symbol "Caution! Refer to the instructions for use";
- graphic symbol "Read the instructions for use";
- graphic symbol "The ban on the re-use";
- graphic symbol "biological risk";
- RU number and date;
- barcode;
- sign NeT (NEW EXCELLENT TECHNOLOGY)

Transportation (Group) packing:

- the manufacturer's name and address;
- the name of the authorized person in the Russian Federation, the address, phone number, website, e-mail;
- full name of the product;
- the abbreviated name of the product;
- the number of products;
- graphic symbol "best before";
- graphic symbol "date of issue";
- graphic symbol "series number" (LOT);
- graphic symbol "manufacturer";
- graphic symbol "Keep at a temperature 1~30°C";
- graphic symbol "Caution! Refer to the instructions for use";
- graphic symbol "Read the instructions for use";
- Net weight;
- circulation mark in the market.

Symbols used on packaging and labels

IVD	Only for invitro diagnostics
	Caution! Refer to the instructions for use
	Store at a temperature of 1~30°C
	The number of tests in the package
	Do not re-use
	Refer to the instructions
	Biological risk
	Manufacturer

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appeared before me and admitted said
principal's subscription to the
STATEMENT

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This is hereby attested on this 30 th
day of Jul. , 2020 at this office.

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HANLIM LAW AND NOTARY OFFICE

BELONG TO SEOUL SOUTHERN
DISTRICT PROSECUTOR'S OFFICE

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KYU-DAI KWON

Attorney at law acting as Notary Public

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public since
February 7, 2020 under Law No. 15150.

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