laserocom
medical and cosmetic laser systems

tflTM
ThermaFractionaLase
Professional TFLs

OPERATING INSTRUCTION

Version 01st March 2013

Title – Operating Instruction THERMAFRACTIONALASE Professional TFLs

The MANUFACTURER

LASER.COM S.r.l.
Galleria Spagna, 35
35127 Padova – Italy.
Tel. +39 049 87 00 242
Fax +39 049 76 28 339
e-mail info@2000laser.com

The PRODUCER FACTORY

ELTECH S.R.L.,
Via Castagnole 20/H-interno 2
31100 Treviso, Italy

The device is manufacture in compliance with the provisions of
Council Directive 93/42/EEC and 2007/47/CE concerning medical devices.
Its compliance is based on the following standards : CEI EN 60601-1.

1. Warning and safety information	5
1.1. Highlighting of warning and safety information	5
1.2. Intended use	5
1.3. Wireless phone interference	6
1.4. Disposal of the device	7
2. Glossary, symbols e abbreviations	8
2.1. Symbols on the THERMAFRACTIONALASE Professional TFLs	8
2.2. Glossary	9
2.3. Abbreviations	9
3. Introduction	10
3.1. Classification	10
3.2. Safety precautions	10
3.3. Transport and storage	14
4. Installation	15
4.1. Packaging contents	15
4.2. Interlock	16
4.3. Connection of the foot switch	17
4.4. Connection of net and switching the laser unit on	18
4.5. Connector of the Scanner handpiece	19
4.6. Connection of the Scanner handpiece	20
4.7. How to use correctly the Scanner handpiece	23
4.8. Scanner handpiece	26
5. Operating Instructions	31
5.1. Switching the device on and off	31
5.2. Electronic access key	32
5.3. Use of the device	32
5.4. List of treatments	36
5.4.1. AREA	37
5.5. Warning and error messages	38
6. Indications and medical precautions	40
6.1. Indications	40
6.2. Contraindications	41
7. Post-treatment cleaning	42
7.1. Cleaning of the device	42
7.2. Cleaning of the scanner handpiece	43

8. Maintenance and assistance	44
8.1. Maintenance	44
8.2. Safety check	45
8.3. Basic assistance	45
8.4. Technical support	46
 9. Manufacturer's declaration regarding electromagnetic compatibility.....	47
9.1. Definitions	47
9.1.1. Emission (electromagnetic)	47
9.1.2. Interference immunity	47
9.1.3. Immunity level	47
 10. Labeling	50

1. Warning and safety information

1.1. Highlighting of warning and safety information

To prevent any personal injury or material damage, please observe the warning and safety information provided in the present operating instructions.

All such information is highlighted as follows :

NOTICE

For additional information,

CAUTION

If there is any risk of damage to the laser unit,

WARNING

If there is any hazard to the life or health of persons.

This symbol indicates that you have to take action.

This symbol indicates that a certain result will occur

1.2. Intended use

THERMAFRACTIONALASE Professional TFLs (Therma**Fractiona**Lase) is a table laser diode device, usable in, dermatology, aesthetic

In the equipment some treatment programs are already stored, the parameters are taken from the current medical literature (see paragraph 5.4). The operator can set up the parameters manually in order to be able to use the equipment with his specific requirements.

The equipment THERMAFRACTIONALASE Professional TFLs has a double emission of 1540 nm +/- 10nm with the maximum power output of 15 W CW and 980 nm +/- 15 nm on 2 W.

This laser unit must only be used by trained personnel in compliance with the applicable occupational safety regulations and accident prevention measures as well as the present operating instructions.

NOTICE

The user is obliged to use only faultless materials, to observe the correct application as well as to protect himself or herself, the patient and other persons against hazards.

WARNING

This laser unit is not intended for operation in areas subject to explosion hazards or in the vicinity of flammable materials or substances.

WARNING

Public legal provisions may include special safety regulations for the protection of persons against laser radiation. These regulations must be complied with.

WARNING

Using controls or settings or performing procedures other than those specified in this manual may result in hazardous radiation exposure.

1.3. **Wireless phone interference**

CAUTION

To ensure the operational safety of medical electrical equipment, the use of mobile wireless phones in practice or hospital environments must be prohibited.

1.4. Disposal of the device

NOTICE

TO ALL OPERATORS "Inkeeping with the 2002/95/CE, 2002/96/CE e 2003/108/CE Directive, with respect to the reduction and disposal of dangerous substances used in electronic and electric devices."

The sign of the barred trash bin displayed on equipment or packaging thereof indicates that upon completion of its life cycle, it should be disposed of separately from other waste.

The separate disposal of dead equipment is the producer's responsibility. Therefore, the user must contact the producer and follow their established procedure for its disposal.

The proper disposal of this equipment will automatically provide for the recycling and proper processing and disposal of the same which helps to prevent possible negative effects on our environment and health, and encourages the recycling of its parts. The improper or illegal disposal by the user will entail the application of administrative sanctions according to the laws and regulations in force.

2. Glossary, symbols e abbreviations

2.1. Symbols on the THERMAFRACTIONALASE Professional TFLs

CE mark in accordance with Council Directive 93/42/CEE and 2007/47/CE, stating the manufacturer's Notified Body

0476

Type B applied part according to IEC 60601-1.

Please refer to manual first

Do not dispose with domestic waste


Laser radiation warning.


Specification of laser output power and wavelength of IR and aiming beam

States the compliance of the THERMAFRACTIONALASE Professional TFLs with IEC 60825-1, edition 1.2.

Warns of potential laser radiation hazards when opening the laser unit

Warns of class 4 laser radiation hazards when using the laser unit

"LASER STOP" button: Press this button in case of an emergency

Use only PROTEK PMP 130-12-Q12

2.2. Glossary

CONTINUOS EMISSION	Continuous laser emission
PULSED EMISSION	Pulsed laser emission (Chopped Mode)
FREQUENCY	Number of laser pulses per second
HERTZ	Measuring unit for frequency
INTERLOCK	Safety device that stops laser radiation when the door of the treatment room is opened.
JOULE	Unit of measure for emitted energy
WATT	Unit of measure for laser power
STOP	End of treatment or treatment break
TIME	Treatment time setting mode

2.3. Abbreviations

Cm ²	Square centimeter
Hz	Hertz
S	Seconds
W	Watt
mW	Milliwatt (one thousandth of a Watt)
J	Joule
nm	Nanometer
V	Volt
IR	Infrared diode
NOHD	Nominal Ocular Hazard Distance according to EN60825-1 : 2003

3. Introduction

3.1. Classification

According to the applicable standards, THERMAFRACTIONALASE Professional TFLs is classified as follows

- Class I, type B according to CEI EN 60601-1:2007
- Class IIb according to Council Directive 93/42/EEC
- Class A according to CEI EN 60601-1-2:2007 of the electromagnetic compatibility
- Class IV laser according to CEI EN 60825-1:2009
- Degree of protection IP 20 (enclosure not water-proof)

CAUTION

The laser unit itself and the accessories cannot be sterilized.

WARNING

The laser unit is not suitable for use in presence of anaesthetics that are flammable when in contact with air, oxygen or nitrogen monoxide.

3.2. Safety precautions

THERMAFRACTIONALASE Professional TFLs is manufactured in compliance with the provisions of Council Directive 93/42/EEC and 2007/47/CE concerning medical devices. Always observe the following precautions :

WARNING

It is recommended to verify that the part to treatment not present neo, tattoo or dark zones that could absorb excessive amounts of energy.

Do not treat absolutely these parts.

WARNING

Before starting the treatment with THERMAFRACTIONALASE Professional TFLs, verify that the patient does not have any cutaneous anomalies in the interesting zone (for example: herpes) or pregnant. In these cases the use of the device is not recommended.

WARNING

Before starting the treatment with THERMAFRACTIONALASE Professional TFLs, verify that the patient is not wearing metallic jewel, for example rings, earrings or piercing.

WARNING

It is recommended to verify that the part to treatment not present skin cancer. Do not treat absolutely this part.

WARNING

After the treatment avoid the direct sun exposure on the treated zone and feed the skin with moisturizing and protective creams, it is recommended also the use of a hat.

WARNING

The execution of operational passages, setup or tests of functioning not described in the present manual can cause a harmful exposure to the laser radiations.

Laser.com Srl is not responsible for damages caused to the non correct use of the device or to the disrespect of the instructions and caution indicated in the present manual

WARNING

Never insert your fingers or any objects into the exit ports. This is important in order to avoid damage to the optical system.

WARNING

In case of an emergency, switch the laser unit OFF immediately. To do this, press the "LASER STOP" button on the control panel.

WARNING

Pay attention to the labels attached on the unit.

WARNING

Do not cover the bottom of the device. The cooling system is located.

WARNING

Operation of this laser unit by unauthorized persons must be prohibited in order to prevent incorrect or improper use. Laser equipment not in use must be protected against unauthorized access. This can be achieved for example by switching off the laser unit after use, so that the electronic access key must be entered prior to further operation.

WARNING

Set up THERMAFRACTIONALASE Professional TFLs properly and completely before putting it into operation (see chapter 4).

WARNING

Verify that the line voltage corresponds to the voltage indicated on the rating plate of the adapter or in the technical specifications.

WARNING

Do not use the laser unit if a visual inspection reveals damage.

WARNING

If you accidentally spill liquid on the unit, immediately stop treatment, disconnect the power supply cable and contact the authorized service center for assistance.

WARNING

Do not attempt to disassemble the laser unit in any case. This is limited exclusively to trained personnel authorized by Laser.com S.r.l.

WARNING

Never use optical instruments such as microscopes, e loupes or magnifiers together with the original protective goggles. Adequate eye protection is no longer ensured in this case.

WARNING

Always cover the connector of the optical fiber with the protection caps when the fiber is not inserted, in order to avoid the damaging of the optical part.

Always verify that powder, dirt or other bodies strangers do not penetrate in the opening for connection of the optical fiber.

WARNING

The unit is not suitable for use in the presence of anaesthetics that are flammable when in contact with air, oxygen or protoxide.

WARNING

Oxygen-saturated materials such as cotton wool can catch fire due to the high temperature that the unit reaches during operation. Label removers and flammable solutions used for cleaning and disinfecting the laser unit should be allowed to evaporate before using the device. Be aware of the fire risk caused by flammable gases.

WARNING

Never direct the laser beam toward paper, plastics or dark surfaces. They may catch fire or be damaged due to the high temperatures produced by the laser beam.

WARNING

In order to avoid damages in the optical part of the TFL scanner handpiece do not ever direct the laser beam toward paper, plastics or dark surfaces. To test the laser, only thermal paper (for example POS device) should be used..

CAUTION

Avoid interference of the laser emission with optical sensors of devices operated in the vicinity of THERMAFRACTIONALASE Professional TFLs.

3.3. **Transport and storage**

3.3.1. **Transport and storage**

THERMAFRACTIONALASE Professional TFLs is supplied in a case that ensures proper and easy transport.

However please observe the following:

To ensure appropriate storage, the laser unit must always be kept in the case supplied. When stored in the box provided by Laser.com S.r.l., the device withstands the following environmental conditions:

- Temperatures from -40°C to +70°C
- Relative humidity from 10% to 95%
- Atmospheric pressure from 800 hPa e 1060 hPa

CAUTION

Do not leave the laser unit in a car parked in the sun, because the temperature inside the car can reach levels that may damage its components.

3.3.2. **Operating conditions**

THERMAFRACTIONALASE Professional TFLs may be operated in the following environmental conditions:

- Temperatures from +15°C to +35°C
- Relative humidity from 20% to 75%
- Atmospheric pressure from 800 hPa to 1060 hPa

CAUTION

Following transport and storage, allow the laser unit to acclimate for about 1 hour in order to reduce the risk of malfunction due to condensation.

4. Installation

NOTICE

If national or local legal regulations require that the installation of the THERMAFRACTIONALASE Professional TFLs be performed by specialized personnel, these regulations must be complied with.

4.1. Packaging contents

THERMAFRACTIONALASE Professional TFLs is delivered with the following:

Quantity	Description
1	Device THERMAFRACTIONALASE Professional TFLs (main unit)
1	Operating Instruction
1	Power Cord
1	Goggles for operator
1	Goggles for patient
1	Scanner handpiece TFL
1	Foot switch
1	Power adapter
1	Interlock connector

4.2. Interlock

As foreseen by the normative the device

THERMAFRACTIONALASE Professional TFLs is supplied with connector interlock for the automatic interruption of the emission laser at distance. Such connector is normally connected to a switch which interrupts the functioning of the laser in case a person enters in the room in which it is used.

Connect the interlock in the socket correspondent.

NOTICE

The installation must be performed by a qualified electrician who is also responsible for the installation and maintenance of the electrical system to which the device is connected.

NOTICE

Please request the technical data sheet with wiring diagram for installation of the interlock device.


NOTICE

If national or local legal regulations require additional or different measures to be taken for the safety of doctor or patient, these regulations must be complied with.

4.3. Connection of the foot switch

➤ To connect the cable of foot switch with the socket correspondent.

Foot switch

4.4. Connection of net and switching the laser unit on

➤ Connect the plug of the PROTEK PMP130-12-Q12 power supply provided in endowment, to the plug placed on the back of THERMAFRACTIONALASE Professional TFLs.

➤ Switch on the device by pressing the adequate switch, placed over the display.

CAUTION

*Verify that the rating corresponds to the power supply rating of the adapter.
An incorrect power supply can compromised the operation of the device and its safety.*


CAUTION


If the display doesn't show any image to make sure that the key of "Laser STOP" has not been pressed.

4.5. Connector of the scanner handpiece


The following instructions explain how to insert the fiber in the connector:


➤ This is the connector.

➤ Unscrew the two parts of the connector

➤ Insert the optic fiber into the connector through the groove.

➤ Rescrew the parts of the connector.

4.6. Connection of the scanner handpiece

In case the optical fibers used are different from the ones supplied, with the device or just as spare parts, it is not possible to guarantee the correct calibration of the equipment, and therefore its security and effectiveness.

WARNING

Before connecting the fiber optic to the device verify the full integrity.

CAUTION

Before connecting the fiber into the device, always verify that the two components of the connector are firmly screwed.

➤ Remove the protection caps from the fiber optic.

➤ Remove the protection caps from the device.

➤ Connect the fiber optic to the device.

CAUTION

Unscrew the connector from the device pressing only on the nearest part to the laser equipment. On the contrary, there is the risk that two components of the connector unscrew them, not allowing to release the fiber by device. If this happened, screw strongly again the connector.

WARNING

Never direct the laser beam or the aiming beam toward the eyes of a person. Patient and operators must always wear the laser protective goggles delivered along with the laser unit.

WARNING

Always protect the optical fiber socket as well as the optical fiber itself with the special protection caps. Do not allow dust, dirt or other foreign bodies to enter the optical fiber socket. Always make sure that the optical system is clean before connecting the optical fiber.

WARNING

Before and during the laser emission verify that the aiming beam is present and that it projects a regular shape. If it does not come projected some shape, or if it is much irregular one, the optic fiber or the device could be damaged. In this case to proceed as follows:

- *Switch off the device and verify that fiber and the relative connections don't show any kind of damage. Verify the correct connection of the fiber.*
- *Verify the condition of the optical fiber, and if it's possible, try to use a new fiber.*
- *If you cannot detect any damage on the optical fiber and the signal of the aiming laser is not visible with the new optical fiber either, switch off the laser and contact Laser.com Srl or your authorized service center.*

WARNING

Using the laser unit without a properly working aiming beam may cause injuries to operators, assistants and doctors.

**CAUTION**

Do not touch the end surfaces of connectors and protect against damage.

**CAUTION**

Do not kink or pinch the optical fiber, as this would cause a breaking hazard.

**CAUTION**

Do not pull on the optical fiber.

**CAUTION**


Note the maximum permissible bending radius of the optical fiber!

- *Short term (during use): 100 x radius of optical fiber*
- *Long term (during storage): 600 x radius of optical fiber*

4.7. How to use correctly the Scanner handpiece

Before starting the treatment on patients, it is important to execute tests on thermic paper to obtain handiness and trust with the scanner handpiece


To do this, follow the instructions :


Get a piece of thermal paper and place it on a plain surface. Be aware that the thermal paper may stain with ink and damage (for overheating) the plain surface.

i NOTICE

The prominence part of the thermal paper must be the one without thermal ink, to be more precise the part normally marked by the advertisement.

Set from the menu the energy at about 21 mJ.


Lay on the paper the scanner handpiece TFL, begin to move it and at the same time activate the laser using the footswitch.

Perform a back straight movement at a uniform speed.


After the first transit, turn the paper and check the result.


The correct result has to be a uniform distribution of the points on the paper.

If the performance is not correct, the distribution of the points is not homogeneous. Here under are described the common errors a not practiced hand can perform.

➤ **ERROR 1:** There is a high compactness of points in a little zone. This means that the footswitch has been kept pressed but the handpiece is not in movement.


SOLUTION : to avoid this mistake, it is necessary, first of all, to release the footswitch and then to stop the movement of the hand.

➤ **ERROR 2:** The distribution of the points on the paper is not constant.

This means the hand speed has not been the same during the application.

SOLUTION: practice the hand to perform a movement the most regular as possible. To do this perform straight movements and at constant speed.

➤ **ERROR 3:** The points are too drawn near, sometimes overlapped, the hand speed is too low. This means that the set parameters are not compatible with the hand speed

SOLUTION: in this case the solution is to increase the hand speed.

This solution allows to increase the distance among the lines created by the scanner handpiece TFL

NOTICE

It is important to perform several tests. In this way the operator acquires experience and learn to estimate, modifying the parameter of speed of the TFL area, the correct speed of his own hand.

WARNING

Only after having acquired certainty and a good hand with the scanner handpiece TFL on the thermal paper, the treatments on the patient skin can begin.

4.8. Scanner handpiece

WARNING

Before and during the laser emission verify that the aiming beam of the scanner handpiece is present. If it does not come projected some shape, or if it is much irregular one, the optic fiber or the device could be damaged.

WARNING

Before using the scanner handpiece verify its correct operating, in particular the lens vibration.

WARNING

It is recommended to verify that the part to treatment not present neo, tattoo or dark zones that could absorb excessive amounts of energy.

Do not treat absolutely these parts.

WARNING

Before starting the treatment with THERMAFRACTIONALASE Professional TFLs, verify that the patient does not have any cutaneous anomalies in the interesting zone (for example: herpes) or pregnant. In these cases the use of the device is not recommended.

WARNING

Before starting the treatment with THERMAFRACTIONALASE Professional TFLs, verify that the patient is not wearing any metallic jewel, for example rings, earrings or piercing.

WARNING

It is recommended to verify that the part to treatment not present skin cancer. Do not treat absolutely this part.

WARNING

Before starting the treatment it is recommended to cleanse the skin. Do not use inflammable products.

WARNING

Before starting the treatment it is recommended to test the ache tolerance of the patient. Do not exceed 70 mJ of energy dosage.

**WARNING**

Before starting the treatment verify the energy dosage to be irradiated. The device stores in its memory the energy dosage irradiated in the previous treatments.

**WARNING**

We recommend clients not to sunbath before and after THERMAFRACTIONALASE Professional TFLs treatments.

**WARNING**

After the treatment avoid the direct sun exposure on the treated zone and feed the skin with moisturizing and protective creams, it is recommended also the use of a hat.

**WARNING**

It is recommended not to use the scanner handpiece TFL within 1 cm of the bony orbit of the eye.

**WARNING**

It is recommended not to perform the treatment with THERMAFRACTIONALASE Professional TFLs in case of active infections, blood disorders and immunosuppression of the patient. These ones can compromise the healing ability of the skin. Also clients with medical conditions or receiving treatment that significantly compromise healing response should not be treated.

**WARNING**

Perform the treatment with THERMAFRACTIONALASE Professional TFLs only on intact and healthy skin, with the exception of acne-affected skin.

**WARNING**

Diabetes. Clients must provide clearance from their primary care physician before being treated with THERMAFRACTIONALASE Professional TFLs.

**WARNING**

THERMAFRACTIONALASE Professional TFLs is not yet evaluated on darker skin types including types V and VI.

WARNING


Any hair that lies on the skin surface during treatment will absorb laser emission, limit full contact and cause unnecessary discomfort or present other complication to the skin. Therefore, all treatment areas should be reasonably clear of hair prior the treatment.


WARNING

Although little scabs may rarely appear after treatments, it is recommended not to scratch them or take them off the skin. Scabs must fall by themselves and naturally. Patients with the tendency of removing their scabs and keloids should not be treated or warned not to do that.


The scanner handpiece is already supplied complete of optical fiber strenghtened and the connectors. Therefore it is not necessary no particular operation before its use.

Scanner handpiece


The connectors of the scanner handpiece


Screw the connectors.


Scanner handpiece shall be carefully drawn near the skin until light touch with surface, the scanner handpiece shall be held with an angle 90° to the treated area surface. Begin the movement from top downward (see the arrows on the figure) and simultaneously activate the laser radiation by pressing the foot switch. It is important not to press the skin surface with the scanner handpiece as the scanner handpiece lens may scrape the skin. Scanner handpiece shall touch the treated area. Perform slow, smooth and vertical (from top downward) movement over the surface.


Having finished the movement, release the foot switch, following which remove the scanner handpiece. Prior to another scan with the scanner handpiece, wait for a few minutes in order to check if the patient has any possible pain of the treated area.

Having performed the vertical scan, turn the scanner handpiece by the angle of 90 degrees in order to change the direction of next scan. Perform next scans in the like manner by moving horizontally (from left to right).

Moving the scanner handpiece, the operator must give particular emphasis to regular movements. It is impermissible to make curved or angled lines, to move the scanner handpiece in the wrong direction, to linger in some position. It is necessary to strictly observe information displayed.

! WARNING

The operator must not move the scanner handpiece in the same direction of the microbeam generation.

! WARNING


When moving the handpiece on the surface, the operator must pay particular attention on the correct vertical movement. Do not move the handpiece in the wrong way and strictly follow what the image on the tool shows. When using it wrongly, the high energy may damage the patient's skin.

! WARNING

The TFL scanner handpiece must be gently put on the skin and lightly touch the surface while moving up and down as shown by the arrows on the label. It is also important not to press on the skin since the lens on the handpiece may scratch the surface.

! WARNING

Perform a constant straight movement on the skin, as shown below. No curves or bows must be delineated with the movement because the high energy may be harmful to the patient.

! WARNING

While treating the surface do not ever halt in one single position because the energy would be too much focused in that single area and therefore it may hurt the patient.


5. Operating Instructions

5.1. Switching the device on and off

In order to switch on the device use the switch over the display. The control of the device utilizes a comfortable interface with touch-screen.

➤ Touch the screen to continue, or wait the next screen.

In order to switch off the device use follow the instruction:


➤ Enter the Power off area

➤ Touch OK on the display.

5.2. Electronic access key

We remind you that THERMAFRACTIONALASE Professional TFLs has to be used by authorized personnel only. For this purpose when the unit is started, once the auto-test has run, an image requires the introduction of an access code which works as an electronic key. This access code is represented by the sequence of four number.

➤ Insert the code 0000 and press Ok.

WARNING

Do not write down the electronic access key somewhere else. Risk of misuse of the laser by unqualified persons

5.3. Use of the device

You can select the desired application from the areas of application available by pressing the arrow key on the touch screen.

**REFRESHING
AREA**

**SCAR/ACNE
SCAR AREA**

**DARK SPOTS
AREA**

**TFL MANUAL
SETTING AREA**


It allows the customer to set up ten various programs manually. The programs set up to you will remain stored for forward us.

Confirm the desired application area by pressing the "OK" key and select the treatment from the list that you see on the screen.

See paragraph 5.4 for the complete list of the treatments.

The key "SETUP" allows to enter to a menu that allows the customer to personalize the device:

Once selected the area it is visualized a picture like the following, which allows to select the desired application.

In the inner part of the screen the application to be selected are shown, eight simultaneously.


The single arrows allow to scroll application one at a time.

Double arrows allow to scroll the complete page. In this way operator can all the list faster.

Once selected a treatment program it is visualized a warning message on the screen in order to remember the obligation to wear the protective goggle during the treatment.

➤ *Respect the warning and wait for some second.*


In the case that the screen visualizes an error message consult the paragraph 5.5. "Wrong message".

Description of the parameters in the start screen:

- **Speed:** speed of the scanner (from 1 to 10). The **speed** should be set depending on the duration of the treatment that the operator wish to perform. At maximum **speed** (10) the operator should move the handpiece at 1 cm/second.
- **Ton:** duration of a single microbeam
- **ENERGY:** explain the energy of the treatment
- **SHOT:** number of spot forming the width of one single line (3-5-7-9). The shot must be set depending on the duration of the treatment the operator wish to perform.

For example, if you set the shot at 3, the operator should move the handpiece at 1cm /second.

This the example of SHOT 3: 

This the example of SHOT 5: 

After having chosen the selected program, it is sufficient to press the foot switch in order to activate the emission of the laser beam. A written "LASER ON" is visualized from the device in order to inform that the emission is activated.

An acoustic signal informs of the laser emission, if it is not deactivated from the customer. In order to interrupt the emission momentarily, and to return to the phase of ready, enough to release the foot switch.

Pressing key "STOP" during the phase of ready it is possible to exit the treatment and to return to the begin menu.

WARNING

Anytime you can interrupt the laser emission by pressing the key of "Laser STOP" on the device.

5.4. List of treatments

The parameters of the treatments are loaded following the indications of medical literature. Operator can modify the parameters of 10 personal programs, and can do it in the specific manual settings. So he can create personal programs for his own particular necessities.

Follow these warnings carefully:

WARNING

It is recommended to verify that the part to treatment not present neo, tattoo or dark zones that could absorb excessive amounts of energy.

Do not absolutely treat these parts.

WARNING

Please run the treated zone with the scanner handpiece, keeping in contact with the skin of the patient,, in order to check the sensibility of the patient and the energy density, especially if the patient's skin is particularly brown.

WARNING

Please start to perform the treatment in LOW modality. This fact allows to test the pain sensibility of the patient. After it, it is possible to perform the treatment in MID modality, and then in HIGH modality.

WARNING

Before a new scansion with the scanner handpiece TFL, wait some minutes, in order to test eventual pain on the treated area of the patient.

5.4.1. AREA

Area		Ton (ms)	mJ x microbeam	Skin Types	Range of recommended age
REFRESHING	LOW	1.4	21	(I,II,III,IV)	< 35
	MID	2.0	30	(I,II,III)	36 – 45
	HIGH	2.7	45	(I,II)	> 45
ACNE SCARS/SCARS	LOW	2.7	40	(I,II,III,IV)	< 40
	MID	3.4	51	(I,II,III)	41 – 50
	HIGH	4.0	60	(I,II)	> 50
DARK SPOT	LOW	2.7	40	(I,II,III,IV)	< 35
	MID	3.4	51	(I,II,III)	36 – 50
	HIGH	4.0	60	(I,II)	> 50

WARNING

Only acne and dark spot, like melasma, cloasma, lentigos (age spots and sun spots) can be considered pathologies and marked by CE 0476.

WARNING

The parameters above must be selected according to the age of the patient too. Young skins have a greater absorption than aged, therefore it is recommended to use LOW and MID for the first ones and MID and HIGH for the second ones.

WARNING

Consider the age-range shown in the column "Range of recommended age" just indicative. It is necessary the doctor values the proper dosage of energy to plan for each single patient.

If the beam density is 600 microbeam/cm², the speed of the hand of the operator must be, at least, 1 cm/sec.


NOTICE

The diameter of each microbeam delivered by THERMAFRACTIONALASE Professional TFL is less than 100 µm.

5.5. Warning and error messages

If this signal appears, it means that the connector of the scanner handpiece TFL has not been inserted in correct way or that the scanner handpiece TFL does not operate correctly. Control the connector or check the right operating of the device. Push the touch screen.


The scanner handpiece TFL needs periodical calibration. This signal remembers that the handpiece must be calibrated. Push Yes to continue, but contact immediately the technical support.


WARNING

A missed calibration can cause a bad functioning of the scanner handpiece TFL and it can provoke cutaneous damages to the patient.

If this signal appears, it means that the optic fibre has not been inserted in correct way. Insert the optic fiber to continue.

ATTENTION!

Shut the
access door
and insert
the interlock

interlock sign

The interlock alarm is on.

Verify the connection of the connector
interlock

CAUTION!

Both operator
and patient must wear
protection glasses

Use the protection glasses.

Respect the warning

ATTENTION!

Please wait
excessive temperature!

Equipment has been overheated
and its operation has been blocked in
order to prevent the damaging. Wait
the cooling of the device for some
minutes for being able to resume the
treatment.

6. Indications and medical precautions

6.1. Indications


Laser therapy has now become a multi-use treatment method. Particularly in dermatology and aesthetic medicine.

Thermafractionalase employs a 15W maximum power and a double wavelength of 1540nm and 980nm. The device works in a fractional mode, this means that the aiming beam hits the localized sectors producing a feeling of slight itching but absolutely painless. This is due to the exclusive TFL scanner system.

The scanner makes some small thermal lesions on about the 20% of the treated surface, hundreds of invisible microbeams through a slight punctiform redness. Such a cutaneous event lessens itself in 24 hours.

The fractionation of the laser affected area is of basic importance during the cure and new derma activation procedure, this because the fulfilled microlesions are separated among them from intact tissue which baits an intense process of tissue regeneration and of fibroblasti production which create themselves collagen and elastina.

The following figure shows the penetration of the electro-magnetic radiation in connection with the wavelength.

The combined use of the 980nm wavelength allows to create a further tissue biostimulation effect, accelerating in this way the regeneration period of the treated area.

The emission at 980 nm has a good absorption from the water, and therefore it raises the thermal effects associated to the use of the laser in all that treatments in which the thermal effect is not a contraindication but a beneficent effect.

The use of the combination 980 nm, is therefore very effective when the thermal effect can be associated to the effect of biostimulating in therapy, as for instance in the pain therapy.

The main biological effects induced by this radiation are the following:

- Increase in the flow of blood due to capillary and artery vasodilatation;
- Modification of the hydrostatic pressure in the capillaries;
- Stimulation of the electrolyte exchange;
- Rise in the perceptive threshold of sensitive nerve endings whose result is a lower sensitivity to pain

6.2. Contraindications

- Pregnancy, although there are no studies that prove any damage to the fetus deriving from this type of treatment, its results may be affected by the hormonal fluctuation common to this state.
- Herpes, must have completely healed from infection prior to treatment.
- Open injuries, the treatment can be performed only on healthy skin.
- Medications, do not perform the treatment on patients who have taken anticoagulant or photosensitive drugs.
- Allergies, verify that the patient does not have any allergic reaction to any substance used during or after the treatment.
- Epilepsy, all devices which emit light could trigger epileptic crisis.
- Cheloids, although crusting phenomena is rare, the voluntary removal of these crusts may give rise to cheloids. Subjects prone to cheloid formation should be excluded.
- Active infections, immunosuppression, immunodeficiency and active infections can depress the skin's capability to regenerate. Therefore treatment must restart only after complete recovery.
- Metabolic disorders (for example diabetes). Metabolic disorders can alter the natural healing process of the skin. These patients must obtain prior written consent from their personal physician.
- The high temperature of the patient (above 36,6)

7. Post-treatment cleaning

WARNING

After finishing the treatment, switch off the laser and disconnect the power cord from the power supply.

WARNING

The THERMAFRACTIONALASE Professional TFLs and all the accessories are not sterilizable.

7.1. Cleaning of the device

Dust can be removed from THERMAFRACTIONALASE Professional TFLs using a soft and dry cloth. Hard spot can be removed using a wet cloth.

To disinfect the THERMAFRACTIONALASE Professional TFLs use products that have been tested and approved by the competent national bodies, or have been proven to have bactericidal, fungicidal and virucidal properties.

Follow carefully the indications of disinfectant.

7.2. Cleaning of the scanner handpiece

WARNING

Do not use solvents or any aggressive solutions on the crystal and the scanner handpiece in general.

Dust can be removed from scanner handpiece using a soft and dry cloth. Hard spot can be removed using a wet cloth.

To disinfect the scanner handpiece use products that have been tested and approved by the competent national bodies, or have been proven to have bactericidal, fungicidal and virucidal properties.

Follow carefully the indications of disinfectant.

NOTICE

Optical fiber can be damaged easily if it suffers an excessive flexions or it is used in a not correct way. This fact can compromise the usefulness of the treatment.

8. Maintenance and assistance

8.1. Maintenance

THERMAFRACTIONALASE Professional TFLs does not require particular maintenance. In case of any operating failure, please see paragraph 8.4.

The manufacturer recommends to do all these operations regularly for the correct working of THERMAFRACTIONALASE Professional TFLs:

ACTION	TIME	PERSON IN CHARGE
Check the handpiece	Before each treatment session	Operator
Check the calibration of scanner handpiece TFL	When the warning message of calibration appears	LASER.COM S.r.l., authorized service center
Check the power of emission of laser	Six months	Operator
Technical safety check	Every 2 years	LASER.COM S.r.l., authorized service center

NOTICE

If national or local regulations require additional safety checks for your laser unit, these regulations must be complied with and the corresponding checks must be performed.

The manufacturer accepts responsibility for the safety of the laser unit only if the following prerequisites are fulfilled:

- Modifications of the laser unit or repair work must be performed by authorized personnel.
- The electrical system at the premises where the THERMAFRACTIONALASE Professional TFLs is installed must comply with the applicable legal requirements for electrical installations.
- The laser unit must be used in compliance with the instructions provided in the present manual.
- The check of the scanner handpiece TFL must be performed when the calibration message appears.

8.2. Safety check

The following safety checks must be performed every 2 years by a qualified service engineer:

- Visual inspection of the unit and its accessories for mechanical damage that might impair operation
- General function check
- Check of the audible and visual indicators
- Check of electric safe
- Check of power of laser

8.3. Basic assistance

In case of defective operation, consider the following points before calling the service authorized by LASER.COM S.r.l.:

- Check the connection of power cord.
- Check the connection of Interlock.
- Check the connection of the optic fibre.
- Check the connection of the two parts which compose the fiber connector
- Check the correct operating of the scanner handpiece
- Be sure that all the operations have been correctly carried out.
- Check the connection of the foot switch.
- Check the emergency stop button.

8.4. Technical support

If technical assistance is needed, please contact our company at the following address:

LASER.COM srl
Galleria Spagna, 35
35127 Padova – Italy.
Tel. +39 049 87 00 242
Fax 39 049 76 28 339
e-mail info@2000laser.com

In case it is necessary to send back THERMAFRACTIONALASE Professional TFLs for assistance, please insert the device into the supplied case and pack it safely into a carton.

Laser.com s.r.l. is the only provider of technical service for its products, unless otherwise indicated by the company. This is to guarantee the product's safety and for the preservation of the product's characteristics.

9. Manufacturer's declaration regarding electromagnetic compatibility

NOTICE

THERMAFRACTIONALASE Professional TFLs complies with all requirements for electromagnetic compatibility (EMC) according to CEI EN 60601-1-2: 2003

9.1. Definitions

9.1.1. Emission (electromagnetic)

When electromagnetic energy is emitted by a source

9.1.2. Interference Immunity

The ability of a device or system to work without errors even if there is electromagnetic interference

9.1.3. Immunity level

The maximum level of a certain electromagnetic interference that affects a particular device or system, where the device or system remains operative with a certain level of performance

ELECTROMAGNETIC EMISSION		
EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIROMENT- GUIDANCE
RF Emission CISPR 11	Group 1	The model THERMAFRACTIONALASE Professional TFLs uses RF energy only for its internal function. Therefore, its RF emission are very low and are not likely to cause any interference in nearly electronic equipments The model THERMAFRACTIONALASE Professional TFLs is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emission CISPR 11	Class A	
Harmonic emission IEC 61000-3-2	Class A	
Voltage fluctuation / flicker emission IEC 61000-3-3	Complies	

Interference immunity			
THERMAFRACTIONALASE Professional TFLs is intended for use in the electromagnetic environment specified below. The customer or the user of the THERMAFRACTIONALASE Professional TFLs should assure that is used in such an environment.			
Immunity test	EN 60601-1-2 Test level	Compliance level	Electromagnetic environmental- guidance
Electrostatic discharge (ESD) EN 61000 - 4 - 2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electric fast transient/ burst EN 61000-4-4	+/- 2 kV for power supply	+/- 2 kV for power supply	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	+/- 1 kV differential mode +/- 2 kV common mode	+/- 1 kV differential mode +/- 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptios and voltage variation on power supply input lines EN 61000-4-11	< 5% Ut (> 95% dip in Ut) for 0,5 cycle 40% of Ut (60% dip in Ut) for 5 cycle 70% of Ut (30% in Ut) for 25 cycle < 5% Ut (> 95% dip in Ut) for 5 cycle	< 5% Ut (> 95% dip in Ut) for 0,5 cycle 40% of Ut (60% dip in Ut) for 5 cycle 70% of Ut (30% dip in Ut) for 25 cycle < 5% Ut (> 95% dip in Ut) for 5 cycle	Main power quality should be that of a typical commercial or hospital environment. If the user of the THERMAFRACTIONALASE Professional TFLs requires continued operation during power main interruption, it is recommended that the THERMAFRACTIONALASE Professional TFLs be powered from an interruptible power supply or a battery.
Power frequency(50/60 Hz) magnetic field EN 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the monitor of the THERMAFRACTIONALASE Professional TFLs further from source of power frequency magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is

			sufficiently low.
--	--	--	-------------------

RF immunity

THERMAFRACTIONALASE Professional TFLs is intended for use in the magnetic environment specified below. The customer or the user of the THERMAFRACTIONALASE Professional TFLs should assure that is used in such an environment.

Immunity test	EN 60601-1-2 Test level	Compliance level	Electromagnetic environmetal-guidance
Conducted RF EN 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and communications equipment should be used no closet to any parts of THERMAFRACTIONALASE Professional TFLs, included cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d=1,2 \cdot \sqrt{P}$ 150 kHz to 80 MHz $d=1,2 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d=2,3 \cdot \sqrt{P}$ 800 MHz to 2,5 GHz Where "P" is the maximum output power rating of the transmitter manufacture and "d" is the recommended separation distance in meters (m).
Interference RF irradiate EN 61000-4-3	3 Vrms Da 80 MHz a 2,5 GHz	3 Vrms	

Field strengths from fixed transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range

Interference may occur in the vicinity of equipment marked with the symbol.

Recommended separation distances between portable and mobile RF communications equipment and the THERMAFRACTIONALASE Professional TFLs

THERMAFRACTIONALASE Professional TFLs model is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and THERMAFRACTIONALASE Professional TFLs as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	From 150 kHz to 80 MHz $D = 1,2 \cdot \sqrt{P}$	From 80 MHz to 800 MHz $D = 1,2 \cdot \sqrt{P}$	From 800 MHz to 2 GHz $D = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance "d" in meters can be determined using the equation applicable to the frequency of the transmitter, where "P" is the maximum output power rating of the transmitter in watts according to the transmitter manufacturer


NOTICE


At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

10. Labeling

Label 2
Label 3


Label 5

Label 1: Warns of class 4 laser radiation hazards when using the laser unit


Label 2: Warns of potential laser radiation hazards when opening the laser unit

Label 3: States the compliance of the THERMAFRACTIONALASE Professional TFLs with IEC 60825-1, edition 1.2


Label 4: Laser radiation warning

Label 5: Specification of laser output power and wavelength of diode and aiming beam

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

THERMAFRACTIONALASE Professional TFLs

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