Notariat II Tuttlingen

Notar Haller

Bahnhofstraße 103 ♦ 78532 Tuttlingen
Tel.: 07461/98372 ♦ Fax: 07461/98358

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Important information for users of KARL STORZ devices

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1 Important information for users of KARL STORZ equipment

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your confidence in the brand name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

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E.A.S.I. =
Endoscopic Automatic System for Irrigation

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This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HYSTEROMAT E.A.S.I.®. All required details and all actions on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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General warning

Users of the HYSTEROMAT E.A.S.I.® must be thoroughly trained in its use for the applicable medical procedures.

The device and accessories should not be used until you have familiarized yourself fully with the instruction manual, assembly, and the care instructions. Instruction manuals should be carefully studied and be available to the surgical team during the procedure. It is essential to follow the instructions contained in the instruction manual, paying particular attention to the Cautions and Warnings.

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
- 2 Images of the equipment
2.1 'Stand alone' operation

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
- 2.2 Combined operation with
UNIDRIVE® S III

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

3 **Controls, displays,
connectors and their uses**

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Controls, displays, connectors and their uses

- 1 Power switch, illuminated
0 = off
- 2 Display and touch screen
- 3 Irrigation pump with locking lever
- 4 Suction pump with locking lever
- 5 Pump roller
- 6 Pressure sensor for the pressure measuring chamber of the irrigation tube
- 7 Tube holder for the suction tube
- 8 Footswitch connection socket
- 9 Connection for external systems
(e.g., UNIDRIVE® S III)
- 10 Connections for SCB*
- 11 Potential equalization connector
- 12 Power cord socket
- 13 Fuse holder

*The SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions and remote display of equipment parameters.

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Potential equalization connector

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Applied part type BF in accordance with IEC 60601-1

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Electronic information product pollution control (China RoHS)

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This device has been labelled in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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RxONLY

Federal law restricts this device to sale or use by or on the order of a physician

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1	Important information for users of KARL STORZ equipment	III
2	Images of the equipment	IV
2.1	'Stand alone' operation	IV
2.2	Combined operation with UNIDRIVE® S III	IV
3	Controls, displays, connectors and their uses	V
4	Symbols employed	VII
5	General information	4
5.1	The HYSTEROMAT E.A.S.I.® device	4
5.2	Intended use	4
5.3	User qualification	5
5.4	Operating modes	5
5.5	Operation	5
6	Warnings and cautions	6
6.1	Definitions	6
6.2	Training in the operation and function of the device	12
6.3	Patient profile	13
6.4	User profile (specialist physician)	13
6.5	User profile (assistant)	14
6.6	Use	14
6.7	User position	15
6.8	Safety precautions at the site of installation	16
6.9	Safety precautions when operating the device	16
7	Explanation of terms (in alphabetical order)	17
8	Installation and operating instructions	22
8.1	Unpacking the equipment	22
8.2	Basic equipment	22
8.3	Installing and connecting up the device (preparations)	22
8.4	Operating the unit	25
8.5	Preparing connection of the irrigation tube	26
8.6	Preparing connection of the suction tube	29
8.7	Combined operation with UNIDRIVE® S III	32
9	System description	33
10	Menu description	37
10.1	Components of the header/footer	37
10.1.1	Components of the display	37
10.2	Start screen	38
10.3	Operation Selection menu	39
10.3.1	Function and use	39
10.4	Procedure List	40
10.4.1	Function/use	40

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5 General information

5.1 The HYSTEROMAT E.A.S.I.® device

The HYSTEROMAT E.A.S.I.® is a combined pressurized suction and irrigation pump with a continuous fluid flow for hysteroscopic interventions. Alternatively, the two pumps can be operated independently of each other with flow regulation, e.g. for laparoscopic interventions.

In pressurized mode, the control compensates for the flow-related fall in pressure between the supply and the opening of the hysteroscope both vis-a-vis the display and setpoint. This allows the user to work directly with intrauterine pressure.

The device accumulates the irrigation liquid supplied and can inform the user when preset limit values have been exceeded (monitoring of irrigation fluid supply).

The HYSTEROMAT E.A.S.I.® is operated directly via a color screen with a touch-sensitive surface (touch screen) or by an external prime control unit via the SCB interface. The operation displays the information required in the specific context to enable a comprehensive assessment of the situation in a clear manner, and guides and aids the user with all device-related work.

5.2 Intended use

Suction/irrigation pumps and their accessories are used for the introduction of irrigation fluids into organs and operating fields as well as the suctioning off of irrigation and bodily fluids, secretions and tissue during endoscopic interventions.

For safety reasons, the device must only be used as intended. Use of the HYSTEROMAT E.A.S.I.® is contraindicated if, in the opinion of the responsible physician, the patient's health is endangered through its use.

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The HYSTEROMAT E.A.S.I.® may only be operated with tubing sets, sheaths and accessories which have been designated as suitable for the device by KARL STORZ.

This applies in particular for tubes and sheaths, since the device uses the predetermined characteristics of these elements for determining the intrauterine pressure. No responsibility for the safe operation of the device can be accepted if any accessories other than those specified are used.

The device can be used in combination with a HF surgical device.

5.3 User qualification

The HYSTEROMAT E.A.S.I.® may only be used by physicians and trained medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the device.

5.4 Operating modes

The HYSTEROMAT E.A.S.I.® has 4 operating modes for surgical and diagnostic interventions/operations. The operating modes are 'HYS' (hysteroscopy) for gynecological diagnostics, 'RES' ('resection') for surgical hysteroscopic interventions, 'LAP' ('laparoscopy') for laparoscopic interventions and 'SHA' ('shaver') for surgical hysteroscopic interventions with the shaver.

5.5 Operation

The device features a color monitor with a touch-sensitive screen for making entries (touch screen). The information displayed and the possible entries are adapted to the specific work flow status.

Both pump beds feature a rotating tensioning lever with which the user can open or close the pump bed for insertion and removal of the tubes.

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6 Warnings and cautions

Please read this manual carefully and follow its instructions precisely. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this instruction manual, they should be carefully reviewed to ensure the safe and effective operation of this device. To make the signal words stand out more clearly, they are accompanied by a pictogram.

6.1 Definitions

WARNING: A Warning indicates that the personal safety of the patient or physician may be endangered. Failure to observe a Warning could result in injury to the patient or physician.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: Notes contain additional information on the safe operation of the device.

WARNING: Read this instruction manual thoroughly and be familiar with its contents prior to using this device. In particular, read the chapter on safety instructions to avoid putting your patients, your personnel or yourself at risk.

WARNING: The electrical installations in the operating room in which the device is connected and operated must comply with the applicable IEC standards.

WARNING: Keep the unit out of reach of patients.

WARNING: For safety reasons, when the system is in use the output sockets of the device and the patient must not be touched at the same time.

WARNING: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely.

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Safety instructions

Warnings and cautions

WARNING: Combinations of medical devices are only assured to be safe if:

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

WARNING: To avoid the risk of an electric shock, this unit may only be connected to a power supply network with a protective conductor

WARNING: When used in combination with HF surgical devices: Only HF operating modes with a maximum output voltage of 5.5 kVpp may be selected.

WARNING: Grounding reliability can only be achieved when the equipment is connected to 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

WARNING: Test this device prior to each surgical procedure to ensure that it is functioning correctly. It should not be used if any damage is evident.

WARNING: Always unplug the device before performing any maintenance work on it.

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WARNING: Do not open the device! Danger of an electrical shock. Arrange for service work to be carried out by the manufacturer or by personnel authorized by the manufacturer only. Any opening of the device by unqualified persons will result in the forfeiting of all warranty claims.

WARNING: Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's wellbeing is maintained in the event of such interference.

WARNING: Risk of explosion if used in the presence of flammable anesthetics.

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Safety instructions

Warnings and cautions

WARNING: Reel pump (pressure-controlled pump motor). Do not touch the pump reels 5 when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

WARNING: The irrigation tubing set must be filled with irrigation fluid before use on a patient.

WARNING: The fluid level in the irrigation bag must be monitored carefully to ensure that no air is pumped into the cavity. Another bag of irrigant should be started prior to emptying the bag being used.

WARNING: The intracavitary pressure is always equal to or lower than the pressure at the outlet of the HYSTEROMAT E.A.S.I.® which is used to regulate the device.

WARNING: The user should always use the lowest possible pressure for sufficient distention and/or irrigation.

WARNING: Do not reuse disposable tubing sets. Discard after use.

WARNING: If a footswitch no longer reacts when pressed, the device must not be operated.

WARNING: Before you activate the pump, make sure that the tubing sets have been inserted properly and that the pressure measuring chamber is correctly positioned.

WARNING: The time display showing the seconds in the header indicates that the device control is working properly. If the changing seconds are not visible, the device must not be operated.

WARNING: Before you start the intervention, make sure that the tubing sets used are not leaking or kinked and that the ends of the tubes are properly secured.

WARNING: The device regulates based on the pressure measured at the outlet of the irrigation pump. The pressure in the cavity depends, among other things, on the height difference between the device and cavity. Try to keep the difference as small as possible.

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Safety instructions

Warnings and cautions

WARNING: Any treatment may only be performed if visual observation of the action of the device is ensured.

WARNING: Attaching a pressurized measuring chamber (cambered diaphragm) can result in an erroneous pressure measurement. The irrigation tube should always be unpressurized when installed.

WARNING: Observe the applicable local regulations when disposing of disposable tubing.

WARNING: Observe the availability signal during power-up. If it can't be heard the device must not be used.

WARNING: Screw the patient tube's connector firmly into the pump tube's coupling to avoid leakages.

WARNING: Servicing may only be performed by specialists authorized by KARL STORZ.

WARNING: The device displays the intrauterine pressure. The correct functionality depends on correct selection of the used sheath. Make sure to use matching sheath and procedure.

WARNING: The HYSTEROMAT E.A.S.I.® may be subject to sudden failures. Therefore you should ensure that a replacement unit is available or be prepared to apply alternative procedures if the intervention requires permanent fluid transportation/distension.

WARNING: The HYSTEROMAT E.A.S.I.® relies on the characteristics of the used sheath for correct pressure measurement. Never use the device with sheaths other than those specified in the instruction manual.

WARNING: The one-day tubing set for irrigation can remain in the reel pump for a maximum of one surgery day (24 hours) during which it can be used for several operations. Afterwards it must be disposed of.

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WARNING: Valves and tube clamps should be either completely opened or completely closed.

WARNING: Make sure that the set volume of the loudspeakers allows for perception of acoustic signals over normal background noise.

WARNING: If the touch screen no longer reacts when pressed, the device must not be operated.

WARNING: The device must not be operated if the touch screen or display are defective.

WARNING: It is not permitted to press more than one touch point simultaneously.

WARNING: The device may only be operated under the environmental conditions specified.

WARNING: Please ensure that there is sufficient air circulation for your device.

WARNING: The signal input and signal output parts of this device are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

WARNING: The device may only be operated at the voltage stated on the rating plate.

WARNING: When cleaning the device, always observe the instructions in the instruction manual (see section 'Cleaning and care').

WARNING: When replacing fuses, only use fuses of the correct rating.

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WARNING: Avoid allowing liquids to enter the housing. Do not store liquids on or above the unit. If, despite precautions, liquid enters the device, sufficient time should be planned for evaporation (this also applies to the formation of condensation).

WARNING: This unit may only be operated with the Irrigation Tubing Set 031717-01 intended for this purpose and which is indicated by the bar on the packaging.

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Safety instructions

Warnings and cautions

WARNING: Do not use in a magnetic resonance (MR) environment. This device is MR Unsafe.

WARNING: The HYSTEROMAT E.A.S.I.® is not to be used where MRI, CT, Diathermy, Electrocautery or RFID equipment may cause Electromagnetic interference (EMI). If the user is unsure or suspects that such equipment may be causing technical problems, the clinical/biomedical engineer of the facility should be consulted.

WARNING: To prevent electrical medical device equipment from being adversely effected by unknown sources of electromagnetic interference, it is highly recommended that Ad hoc radiated RF testing be conducted by a Clinical/Biomedical engineer or an EMC professional.

WARNING: – Continuous Flow Hysteroscopy – If a liquid distension medium is used, strict fluid intake and output surveillance should be maintained to ensure that fluid deficit is known at all times. Depending on whether non-electrolyte or electrolyte solution is being used, when excessive fluid deficit occurs, consideration should be given to stopping further infusion and concluding the procedure. Potential Complications of Continuous Flow Hysteroscopy:

- Hyponatremia
- Hypothermia
- Uterine perforation resulting in potential injury to bowel, bladder, major blood vessels, and ureter.
- Pulmonary edema
- Cerebral edema
- Death

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PRECAUTION:

Gravity fed intrauterine fluid distension can usually be accomplished with pressures in the range of 35-75 mmHg. Hanging the fluid distension medium 42" above the patient can generate intrauterine pressure of approximately 80 mmHg. Unless the systemic blood pressure is excessive, it is seldom necessary to use pressures greater than 75-80 mmHg.

CAUTION: When cleaning the device surface, the area around the pressure sensor ⑤ should be wiped carefully due to its sensitive surface.

NOTE: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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Before using the unit on the patient, it is imperative that you be acquainted with how the unit operates and is controlled.

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6.2 Training in the operation and function of the device

Training should only be provided by persons equipped to do so on the basis of their skills. Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises and to give you information about further training options.

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6.3 Patient profile

Species: Human
Ethnic group, race, skin color: irrelevant
Nationality, language: irrelevant
Gender:
LAP: irrelevant; HYS: female; RES: female;
SHA: female
Age: irrelevant
Weight: age appropriate
Health: suitable for endoscopic interventions
Mental health:
irrelevant, patient is never the user
Affected body part or tissue type
HYS:
Body part: uterus
Access port: vagina
State: suitable for hysteroscopic interventions
RES:
Body part: uterus
Access port: vagina
State: suitable for hysteroscopic interventions
SHA:
Body part: uterus
Access port: vagina
State: suitable for hysteroscopic interventions
LAP:
Body part: abdominal cavity
Access port: abdominal wall
State: suitable for laparoscopic interventions

6.4 User profile (specialist physician)

- Recognized medical skills in the area of gynecology/laparoscopy (specialist physician).
- Adequate powers of comprehension to assess the situation rationally at all times during the application.
- Adequate language skills in at least one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation and use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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6.5 User profile (assistant)

- Completed training as medical assistant.
- Adequate language skills in at least one of the languages used on the device and in the instruction manual.
- Be thoroughly trained in the operation and use of the device.
- Have been trained thoroughly at least once in the use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

6.6 Use

The device may only be used in hospitals and medical practices in the ambient conditions given in the technical data.

- Not in intensive care units.
- Not in mobile facilities (cars or trucks, ambulances, aircraft, ships or trains).
- Position it on a vibration-free, fixed base.
- Visual distance from the physician to the device when operated: max. 2 m.
- Visual distance for operating the device: 30 – 70 cm.
- Viewing angle to the front of the device: max. $\pm 45^\circ$.
- The use of the device requires the division of work.
- Frequency of use: up to several times a day.
- Length of use: from a few minutes to several hours a day.
- Can only be moved by installing it on an equipment trolley.
- Can only be swiveled by installing it on ceiling supply units.

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- Patient must not necessarily be under general anesthesia during the intervention (local anesthesia).
- May be used with other devices (light source, HF surgery etc.).
- To be operated in ambient light between 100 Lux and 1500 Lux.
- Ambient temperature in immediate proximity of the device during operation: max. 40°C.
- Sound pressure from background noise < 60 dBA.
- The device is placed outside the sterile area.
- Must only be operated with parts for single use.

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6.7 User position

- For operation, the user should always be standing in front of the device within a 'cone of vision' with an opening angle of $\pm 45^\circ$ at a distance of approx. 30 – 70 cm to the front panel.
- For observation of the actual values and for operation of the footswitch during the application, a visual distance to the device of 2 m is assumed (the tube length is 2.5 m).

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
6.8 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

The device is equipped with a connector for potential equalization. It should be connected up in accordance with the applicable national regulations.

WARNING: The device must not be used in oxygenated environments. It is not intended for use in hazardous zones.

This means, among other things:

WARNING: When using easily combustible and explosive inhalation anesthetics or mixtures thereof, the device may not be operated in the demarcated hazard zone.

This also applies to easily combustible and explosive chemicals, e.g. skin disinfectants and surface disinfectants.

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6.9 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operates properly before use.

During treatment with the HYSTEROMAT E.A.S.I.®, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring vital signs and the anesthetic.

Any treatment may only be performed if visual observation of the action of the device is ensured.

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7 Explanation of terms (in alphabetical order)

 **NOTE:** The terms are in alphabetical order in all languages. As such, it is possible that the German text on one page does not correspond to the English or Spanish text in terms of content.

Activated: When the line voltage is switched on, the user can move the rollers of the irrigation and suction pump via a button on the screen. This is referred to as 'activation of the pumps'.

 **NOTE:** Further preconditions must be satisfied prior to activation. These are described on page 22 (Installation and operating instructions).

Bipolar Instrument Mode: The Bipolar Instrument Mode is selected by the user if an increased irrigation pressure is necessary, e.g. if an HF electrode is inserted into the working channel of the instrument (endoscope), resulting in a reduction of the instrument's cross section for irrigation liquid.

 **WARNING:** The Bipolar Instrument Mode is designed for the compensation of the HF instruments being used. Do not use while no HF instruments are in use.

Blood/Mucus Mode: In Blood/Mucus Mode, the flow is increased to enable improved suction of blood and mucous particles, thus providing the user with a better view. The Blood/Mucus Mode can be set in the operating menu 'Pressure/Flow Settings' (only with Continuous Flow procedures).

In the Blood/Mucus Mode, the pressure can be increased in percentages by the user. An increase in pressure can, for instance, prevent the cavity from collapsing or push back bleeding. With a 50% increase, the pressure climbs by way of example from 60 to 90 mmHg.

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Continuous Flow: Continuous Flow CF: In the Continuous Flow Mode the pump firstly irrigates in order to build up pressure. The suction pump is inactive until the pressure has been achieved. Once the pressure in the uterus has been achieved, the Continuous Flow Mode is activated. Both pumps then irrigate and suck with the same flow rate (e.g. 50 ml/min). As such, there is no fall in pressure.

The Continuous Flow Mode offers the advantage that clear fluid is continuously conveyed. In contrast to this, with the Single Flow Mode the fluid in the uterus tends to be cloudy or bloody.

Day set/day tube set: The day set concerns only the irrigation tube and not the suction tube.

Purpose of the day set: There is an interface in the irrigation tube set; this is kept sterile. Before each intervention a new patient tube (length approx. 2 m) is connected to this interface. After each operation, rather than replacing the entire irrigation tube only the patient tube of the irrigation tube set is replaced.

This means that the costly sterile parts of the irrigation tube set (e.g. pressure measuring chamber) only have to be replaced once a day; only the non-sterile (e.g. contaminated with blood) patient tube is replaced after every intervention.

Deactivated: Opposite of → activated.

Double Flow Mode: This mode is available for laparoscopy. Here, the suction and irrigation flow rates can be changed and the respective pumps activated and deactivated independently of each other.

Flow modes: There are three different flow modes:

- → Continuous Flow Mode
- → Single Flow Mode
- → Double Flow Mode (for laparoscopy)

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LAP: The LAP procedures are exclusively suction and irrigation procedures.

A pneumoperitoneum for laparoscopic procedures is created via insufflation, e.g., with carbon dioxide gas.

Menus: A distinction is made between

- Operating menus (the user works with these during the intervention – e.g. Normal Mode operating menu)
- 'Settings' menu: this is used to set all basic settings (e.g. prior to the intervention), e.g. pressure/flow settings, monitor language settings etc.

Normal Operation: During Normal Operation, regulation is based on the intrauterine pressure. This is done using characteristics which are saved in the device for the specific sheaths. These characteristics differ depending on the type of intervention (hysteroscopy or resection) and the associated sheath.

The intrauterine pressure is calculated automatically in Normal Operation. This means that the user does not have to worry about falls in pressure from the irrigation liquid bag via the pressure measuring chamber to the instrument. The pressure shown on the monitor in the operating menu 'Normal Operation' is the actual intrauterine pressure.

In the operating menu 'Normal Operation', the user can only change the pressure and not the flow. However, using the menu 'Pressure/Flow Settings', he can set a specific flow for a specific pressure. The flow is then automatically supplied by the system during irrigation depending on the pressure selected.

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Operating state: There are three different pump operating states for hysteroscopy and resection interventions.

- Normal Operation (when neither the Bipolar Instrument Mode nor the Blood/Mucus Mode has been activated)

- → 'Bipolar Instrument' Mode

- → Blood/Mucus Mode

Only the Normal Operation Mode is available for laparoscopic interventions.

Pressure increase: In the → Blood/Mucus Mode, the pressure can be increased in percentages by the user. An increase in pressure can, for instance, prevent the cavum from collapsing or push back bleeding. With a 50% increase, the pressure climbs by way of example from 60 to 90 mmHg.

Procedures: With the HYSTEROMAT E.A.S.I.®, the user works exclusively with procedures. On delivery, the device features predefined procedures for GYN (HYS, RES and SHA) and LAP.

The user cannot change, rename or delete these predefined procedures.

They can, however, copy a predefined procedure and save this under another name and then change the parameters individually (e.g. the desired pressure or flow).

The following applies: a GYN procedure (HYS or RES) with Continuous Flow sheath can only be changed into a GYN procedure with the same Continuous Flow sheath.

A GYN procedure (HYS) with Single Flow sheath can only be changed into a GYN procedure with the same Single Flow sheath.

A LAP procedure can only be changed into a LAP procedure.

A SHA procedure can only be changed into a SHA procedure.

Shaver: Surgical mode for hysteroscopic interventions with the aid of a shaver system (UNIDRIVE® S III).

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**8 Installation and
operating instructions**

8.1 Unpacking the equipment

Carefully unpack the HYSTEROMAT E.A.S.I.® and its accessories from their packing.

Check the delivery for missing items and possible damage.

Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packaging for later use; this can come in handy for transportation purposes.

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8.2 Basic equipment

- 1 HYSTEROMAT E.A.S.I.® 26340020-1
- 1 Power cord 400 A (400 B)
- 1 SCB connecting cord 200901 70
- 1 Instruction manual 96126014DER
- 1 Basic Set Tubing Sets 031917-10 consists of 5 Irrigation Tubing Sets 031717-01 and 5 Suction Tubing Sets 031217-01

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**8.3 Installing and connecting up
the device (preparations)**


WARNING: Observe the section "Safety precautions at the site of installation" on page 16.

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
Place the device on a flat surface.
The difference in height between the device and the operating field should be as small as possible (see warning on page 8).

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Potential equalization

Connect the potential equalization cable to the potential equalization connector ⑪.

! WARNING: The device's ground line should be installed by a qualified electrician.

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Connecting the power cord

Connect the power cord when the device is switched off. Insert the power cord into the socket ⑫ as far as it will go.

The device may only be operated at the voltage stated on the rating plate.

! WARNING: The instrument may only be operated with the power cord delivered by KARL STORZ or a similar power cord which has a national inspection seal.

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! WARNING: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

i NOTE: Ensure that access to the power plug is ensured at the site of installation.

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Connecting the footswitch

Insert the connecting cable of the footswitch into the corresponding socket ⑬.

i NOTE: The footswitch is an optional accessory.

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Connecting the shaver


Insert the control cable 20701070 of the shaver into the corresponding socket ⑭.

i NOTE: The shaver and the control cable are optional accessories.

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Connecting up the SCB cable (optional)


In order to use this function, insert one end of the SCB cable (both ends are identical) into the connection socket ⑩ of the HYSTEROMAT E.A.S.I.®. Connect the other end to an available SCB connection socket on a device in the existing SCB network or a device with which you wish to establish a network.

① **NOTE:** To prevent the SCB connecting cable from being pulled out accidentally, the SCB connector possesses a protection device. Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑩. Connect the other end of the cable to the KARL STORZ SCB control device (KARL STORZ Communication Bus) or other SCB devices (see 'KARL STORZ SCB control NEO System' instruction manual).

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
8.4 Operating the unit

CAUTION: Before switching the HYSTEROMAT E.A.S.I.® on, make sure that no tubing set has been inserted into the device and that the cables are properly connected.

Press the power switch ●.

CAUTION: Do not touch the touch screen while you switch the device on.

Once the device has been switched on, the first thing that is displayed on the screen is the KARL STORZ logo. After this, a white screen appears.

WARNING: Check the white screen for color errors. In the event of any errors the device must not be used.

During this period, the device runs a self test. Once the test has been successfully completed, an acoustic availability signal sounds. An animation appears on the screen which assists you in inserting the tubing set.

If the self test fails, the message 'Safe State' is displayed and an alarm sounds. Verify that the device is properly connected and that the preparations have been made according to the instructions. For further information about troubleshooting please refer to the section 'Technical description: Information and alarm signals' from page 83.

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Animation


The animation shows the order in which the tubing set should be inserted. The animation runs in an endless loop. This menu is exited when the device detects that:

- the irrigation and suction pump levers have been opened and closed and the pressure measuring chamber on the sensor (S) has been adapted (Continuous Flow Mode).
- the irrigation pump lever has been opened and closed and the pressure measuring chamber on the sensor (F) has been adapted (Single Flow Mode).

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**8.5 Preparing connection of the
irrigation tube**

! WARNING: The packaging for the tubing set must not be damaged (sterility!). Check the tubing set for damage etc. before use.

Remove the irrigation tube (blue, with pressure measuring chamber) from the sterile packaging and connect it.

Pass the irrigation tube on for installation in the pump while observing aseptic procedures (safety precautions).

! WARNING: The connection of the irrigation tubing set (sheath) to the patient must remain in the sterile area.

! WARNING: Always connect the patient's end of the irrigation tube (blue) to the irrigation connector of the shaft. If a suction tube is being used, the patient's end (red) must be connected to the suction connector of the shaft (see section 8.6).

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Unlocking the irrigation pump

WARNING: Do not touch the pump reels ⑤ when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.


Unlock the irrigation pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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
Positioning the pressure measuring chamber

Position the pressure measuring chamber and turn it clockwise until fixed.

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Inserting the irrigation tube

Insert the irrigation tube into the pump bed.


WARNING: Make sure that the tube is completely inserted into the pump bed.

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Locking the irrigation pump

Lock the irrigation pump by turning the lever counterclockwise until it clicks into place.

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
Installation and operating instructions

Connect the irrigation tube to the irrigation bag(s). Make sure that no liquids are stored on or above the device (see the section 'Safety instructions', p. 10).

NOTE: The tubes for the irrigation bag are fitted with a blue clamp which can be used to pinch off the tubes if necessary.

WARNING: The irrigation pump is not self-priming. As such, the irrigation bag must always be located to the side and above the device.

Day tube sets

If the device is approved for use with day tube sets, the steps for inserting the irrigation tube set 031767-01 must be performed as described on page 26-27. The separately packed patient tube 031162-01 is first connected with the sheath and then with the irrigation tube via the coupler.

WARNING: The patient tube must be disposed of after each intervention, while the tube inserted in the device may be reused.

WARNING: A new patient tube must be used for each patient.

WARNING: To guarantee the sterility of the connected tube set during patient change-over, the patient tubes are equipped with a separately packed, sterile plug. This must be screwed into the tube coupler from removal of the patient tube from the tube set until the next patient tube is connected.

Handling of the day tube set:

- Following the end of the intervention, the patient tube must be removed from the irrigation tube and disposed of correctly.

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- The plug is then removed from the small sterile packaging and used to seal off the open end of the pump tube.
- The new patient tube is removed from the tube set for the next patient.
- The plug in the pump's irrigation tube is now removed and disposed of, and the new patient tube is connected to the open pump tube under aseptic conditions (safety measures).
- The new plug in its sterile packaging is then connected from the sterile area.

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8.6 Preparing connection of the suction tube

WARNING: The packaging for the tubing set must not be damaged (sterility!). Check the tubing set for damage etc. before use.

Remove the suction tube (red) from the sterile packaging and connect this to the instrument's sheath.

Pass the suction tube on for installation in the pump while observing aseptic procedures (safety precautions).

WARNING: The connection of the suction tubing set to the patient must remain in the sterile area.

NOTE: When using the pump in Single Flow Mode, the steps for connection of the suction tube are not required. If relevant, observe the instruction manual of the suction system used.

NOTE: The tube to the sheath is fitted with a red clamp which can be used to pinch it off if necessary.

WARNING: Always connect the patient's end of the suction tube (red) to the suction connector of the shaft.

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Unlocking the irrigation pump


! WARNING: Do not touch the pump reels ⑤ when the pump is active (even if it is not moving). An activated pump may start moving at any time. Danger of crushing and/or of clothes (or hair) getting pulled into the device.

Unlock the suction pump by turning the lever clockwise as far as it will go. As soon as the lever is actuated, the pump can no longer be activated.

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Inserting the suction tube


Insert the suction tube into the pump bed. Make sure that the tube is completely inserted into the pump bed.

The tube connector at the patient end should be inserted into the suction tube bed at the top in such a way that the locking lever can be turned counterclockwise to the end position without meeting obstacles. The outflow end of the suction tube is inserted in the tube bed at the bottom.

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Locking the suction pump

Lock the suction pump by turning the lever counterclockwise until it clicks into place.

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



Connecting the waste container

Connect the end of the suction tube which is not
close to the patient (with the text 'Waste') to the
waste container.

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8.7 Combined operation with UNIDRIVE® S III

- NOTE:** Observe the instruction manual for the UNIDRIVE® S III.
- NOTE:** The UNIDRIVE® S III is suitable for combined operation with the HYSTEROMAT E.A.S.I.®. The instruction manual describes the two devices as a functional unit. Regardless of this, the UNIDRIVE® S III can also be operated with another pump. If third-party pumps are used, no liability can be accepted for potential problems or risks resulting from a combination with such devices.

Connection of the connecting cable

Connect the devices with the accompanying control cable 20701070 (see Fig. A). The HYSTEROMAT E.A.S.I.® can now be operated using the UNIDRIVE® S III footswitch in 'SHA' mode (UNIDRIVE® S III) (see Chapter 10.13).

Connection of the irrigation tube

Connect the blue colored patient tube to the irrigation sheath (see Fig. B).

- NOTE:** Additional irrigation can be connected if necessary.

Connection of the suction tube

Connect the red colored suction tube to the shaver handpiece (see Fig. C).

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9 System description

Reference

This chapter describes the IEC 60601-1 requirements for medical electrical systems: Medical electrical equipment, section 16.

Definition

A medical electrical system is a combination of individual devices, as specified by the manufacturer, of which at least one is a medical electrical device, and which are inter-connected by a functional connection or the use of a portable multiple socket outlet.

Consisting of the following components:

- UNIDRIVE® S III 20701020-1
- HYSTEROMAT E.A.S.I.® 26340020-1
- Control cable 20701070

Additionally includes the accessory components, handpieces and shaver blades specified in the respective current catalog.

Further system components are not covered by this instruction manual. All modifications or additions to this system require a complete new assessment and documentation for the purposes of IEC 60601-1:2005. Risk management must be observed as stipulated by this standard.

WARNING: All instruction manuals included with the delivery of the individual devices remain valid and must be observed by the user of the system.

WARNING: Components and equipment which are not part of this system must not be connected to the system without additional follow-up documentation, assessment and additional notes for the user.

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Field of application

Aside from the applied parts, no other parts of the system are suitable for use within the patient environment.

Explanation

The system must not be connected to other non-medical devices. All deviations from this instruction require a new assessment of risks and additional documentation with additional warnings and notes for the user if applicable. All medical electrical components are considered devices for the purposes of IEC 60601-1.

Multiple socket outlets

According to IEC 60601:2005, Section 3, 3.67, a multiple socket outlet is two or more sockets, which are intended to be connected to flexible cables, leads or ME devices, or to be integrated into them, for connection to a power supply network carrying a comparable voltage.

WARNING: Multiple socket outlets must not be placed on the floor.

There are two permissible options for connecting to the line voltage supply:

- Each device is supplied from a separate wall socket outlet.
- Both devices are operated in combination via a multiple socket outlet in connection with an isolation transformer (e.g., in a video trolley).

If using any other connection method or combination of the entire system, new measurements must be taken to ensure that the maximum leakage current as per IEC 60601-1:2005 is not exceeded. A multiple socket outlet must be employed, which can only be accessed using a tool, in order to ensure that such a system cannot be modified subsequently.

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WARNING: The use of freely accessible multiple socket outlets is not permitted.

WARNING: The use of additional multiple socket outlets or an extension cable with the system is not permitted.

WARNING: Do not simultaneously touch the patient and non-ME equipment in the proximity of the patient.

Maximum load permitted by the system

For the maximum load permitted by the system, please refer to the instruction manual supplied. The maximum load of a multiple socket outlet in connection with a KARL STORZ isolation transformer depends on the model used and must therefore be checked.

NOTE: If an isolation transformer is connected to the multiple socket outlet, as is described on page 34 of this manual, the isolation transformer is then considered a part of the medical electrical system referred to here.

Reprocessing

Instructions for reprocessing of the system components can be found in the instruction manuals for the system components. Further instructions are not intended for the system.

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Preventative maintenance measures

Any maintenance work which may be necessary to carry out is detailed in the instruction manual for the individual system components. This includes, for example, the timely replacement of wearing parts. In addition, it is advisable to inspect the supply line of any video trolley employed (where applicable) for mechanical damage, and to arrange for a replacement to be carried out by a specialist if damage is detected, before every use.

Once the system has been assembled for the first time with all data lines, signal lines, electrical and non-electrical components, an electrical safety inspection must be carried out. This must be carried out on site or included within the scope of the delivery, if the system has already been inspected for the first time at the factory.

This inspection must be carried out at least once a year. Each alteration to the system (connections, signal lines, cables and components) requires documentation and a new repeat inspection according to the following criteria:

These measurements must be carried out according to EN 62353:2008 for equipment of type BF.

For the permissible limit values for protective ground resistance, earth leakage current, housing leakage current and patient leakage current, please refer to the current standard.

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- Header
- Footer
- Content area

Buttons

- ① Shortcut to the 'Settings' submenu. Here, amongst other things, the following can be set
 - the desired screen language
 - the desired volume of the acoustic signals
 - the 'Fluid Supply Monitor'
- ② Switch acoustics signals ON/OFF
Desired volume: The desired signal volume can be set in the 'Settings' menu (see ①). The state 'switched off' is retained until the user deactivates this function or switches the device off and on again.
- ③ Fluid Supply Monitor. When this button is pressed, the menu 'Fluid Supply Monitor' (monitoring of the amount of fluid in the irrigation bag) is called up, but only if this function has been activated in the menu 'Extended Settings' (see ①).
- ④ Call up the Procedure List.

Displays

- ③ Texts describing the current menu level
- ⑥ Time

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10.2 Start screen


After switching on the supply voltage, a KARL STORZ logo with the device name is displayed for approx. 5 s (operating system loading time) followed by a white screen for at least 2 s.

This should allow the user to detect any pixel errors on the display.

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10.3 Operation Selection menu

10.3.1 Function and use

After starting, the operation type to be performed can be selected.


The operation types 'Hysteroscopy', 'Resection', 'Laparoscopy' and 'Shaver' are available, and can be activated using the buttons ① or ② or ③ or ④.

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10.4 Procedure List

10.4.1 Function/use

After selecting the operation type, the Procedure List with the procedures saved for this kind of intervention is displayed. On the first page, the last four procedures to have been used are shown, with the most recently used one first. The top procedure is shown as 'selected' (orange).

As of page 2, the other procedures are listed according to the procedure number.

If the user wishes to create a new procedure, this can be done using the 'Show Sheaths' button or by copying an existing procedure.

The HYSTEROMATE E.A.S.I.® contains numerous predefined procedures. For hysteroscopy, these begin at 150, for resection at 250, for laparoscopy at 399 and for SHA at 550. Specific key values are saved in each of these procedures.

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

10.4.2 Buttons

- ① Jump to the previous page of the Procedure List
- ② Jump to the next page of the Procedure List
- ③ Shows the current page number/total number of pages
- ④ Jump back to Operation Selection (Hysteroscopy/Resection/Laparoscopy)
- ⑤ Procedure selection is saved. Following this, the system goes to the menu level 'Tube Set Insertion' or into the operating menu 'Normal Operation GYN' depending on whether the tube has already been inserted or not
- ⑥ Name of procedures. When you select a procedure, this is highlighted in orange
- ⑦ Show Sheaths (only for hysteroscopy, resection and Shaver operations). If you press this button, you can create a new procedure (starting with the selection of a sheath, thus the name 'Show Sheaths')

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10.5 Sheath List

10.5.1 Function/use

WARNING: The HYSTEROMAT E.A.S.I.® relies on the characteristics of the sheath used for correct pressure measurement. Only use the device with gynecology sheaths from KARL STORZ. KARL STORZ cannot accept any liability for the use of sheaths from other manufacturers.

Creation of one's own procedure via sheath selection.

Select one of the sheaths saved for this type of intervention (hysteroscopy/resection) and confirm this selection.

On the first page, the last four sheaths to have been used are shown, with the most recently used one first. The top sheath is shown as 'selected'.

As of page 2, the other sheaths are listed according to the sheath number.

By pressing the button 'Show Procedures' you return to the Procedure List.

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
10.5.2 Buttons

- ① List of sheaths. When you select a sheath, this is highlighted in orange
- ② Shortcut to the 'Procedure List'
- ③ Jump to the previous page of the Sheath List
- ④ Jump to the next page of the Sheath List
- ⑤ Shows the current page number/total number of pages
- ⑥ Jump back to the 'Procedure List'
- ⑦ Confirmation of sheath selection, shortcut to the menu level 'Flow Mode Selection'

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
- ① Specification of setpoint for intrauterine pressure.
Example: In this menu, a value of 40 mmHg is specified. This value of 40 is then saved and displayed in the operating menu 'Normal Operation GYN'.
- ② Specification of setpoint for normal flow.
Example: In this menu, a value of 40 ml/min is specified. Irrigation with a flow of 40 ml/min is then achieved in the operating menu 'Extended Display GYN'.
Please note that in the Normal Operation GYN operating menu only the pressure values can be changed and not the flow values. The flow values can only be changed in the menu 'Extended Display GYN'. If the menu 'Extended Display GYN' is to be shown, the button 'Extended Display' in the menu 'Extended Settings' must be pressed as a precondition.
- ③ Specification of setpoint for pressure increase in the Blood/Mucus Mode.
Example: If a value of 80% is specified in this menu, a pressure increase of 40% in addition to the set intrauterine pressure (e.g. 100 mmHg) will be realized in the operating menu 'Normal Operation GYN', but only if the button 'Blood/Mucus Mode' is pressed. Instead of 100 mmHg the user would then work with a value of 140 mmHg.
- ④ Specification of setpoint for flow in the Blood/Mucus Mode
Example: If a value of 150 ml/min is specified in this menu, the irrigation flow rate will be 150 ml/min during pump operation, but only if the 'Blood/Mucus Mode' button was pressed in the operating menu 'Extended Display GYN'.
- ③ Specification of setpoint for SHA suction flow: 140 ml/min

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Origin of the values: The values which are displayed in this menu originate from the procedure which was selected after starting.

Example: If an intrauterine pressure of 40 mmHg was set in the procedure selected after starting, the value 40 mmHg is shown in the menu 'Pressure/Flow Settings'.

Saving changes permanently: Changes which are made here are only retained for as long as no new sheath is selected. If the settings should be retained for longer, the values must be saved in a new procedure. To do this, the button 'Save as New Procedure' (5) must be pressed.

Overwriting the values of the operating menu: Any previously changed settings in the operating menu 'Normal Operation GYN' are overwritten upon confirmation of the settings with the button (7).

Example: You have manually set an intrauterine pressure of 35 mmHg in the menu 'Normal Operation GYN'. If you now set a value of 40 mmHg in the menu 'Pressure/Flow Settings' and save this by pressing the button (7), from now on a value of 40 mmHg is shown in the menu 'Normal Operation GYN'. This value applies until the device is switched off.

Pressures of 65 mmHg upwards: If the intrauterine pressure is set to more than/equal to 65 mmHg, a menu is called up once in which it must be expressly confirmed that pressures greater than 65 mmHg are permitted.

Standard specification values for each new procedure (for HYS or RES):
Intrauterine pressure: 45 mmHg
Normal flow: 50 ml/min
Blood/Mucus Mode Pressure increase: 30 %
Blood/Mucus Mode flow: 150 ml/min

Standard specification values for each new procedure (for SHA):
Intrauterine pressure: 45 mmHg
Suction flow: 140 ml/min

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10.6.2 Buttons


- ⑥ Return to the previous menu; the changed values are discarded.
- ⑦ Return to the previous menu; the changed values are saved.
- ⑧ Increase the specified setpoint for the intrauterine flow.
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) which is to be displayed in the 'Intrauterine pressure' field of the operating menu 'Normal operation GYN' is specified.
- ⑨ Increase the specified setpoint for the intrauterine pressure.
- ⑩ Decrease the specified setpoint for the normal flow.
'Normal flow' refers to the flow associated with Normal Operation and not the flow seen in the Blood/Mucus Mode.
- ⑪ Increase the specified setpoint for the normal flow.
- ⑫ Decrease the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.
- ⑬ Increase the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.
- ⑭ Decrease the specified setpoint for the flow in the Blood/Mucus Mode.
- ⑮ Increase the specified setpoint for the flow in Blood/Mucus Mode.

- ⑫ Decrease the specified setpoint for the suction flow (SHA).
- ⑬ Increase the specified setpoint for the suction flow (SHA).

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10.7 Save as New Procedure

10.7.1 Function/use

The current settings are saved as a new procedure together with a designation, name and number. The designation and name are restricted to 20 alpha-numeric characters or a total width of 200 pixels.

10.7.2 Buttons

- ① Tab – switches between the input windows (display elements 1, 2 and 3)
- ② Return to the menu level 'Pressure/Flow Settings'; the entries are discarded
- ③ Proceed to the menu level 'Normal Operation GYN'; the entries are saved

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10.8 Tube insertion/animation

10.8.1 Function/use

Animation to aid the assistant with insertion of the tube.

This menu level is exited as soon as it is detected that all of the components necessary for the relevant operation (irrigation function/suction function) have been connected.

Example 1: The Continuous Flow Mode is decisive for the selected procedure 'Cont.-Flow HYS B.I.O.H.'. For this mode, both an irrigation and suction tube must be inserted. The menu level is only exited once this has been done.

Example 2: For a Single Flow procedure, only the irrigation tube has to be inserted in order to leave this menu level.

Within a fixed period of 24 hours, operation with a day set is permitted if this has been selected accordingly in the 'Extended Settings'.

This means that the irrigation tube (excluding the patient tube element) can remain in the pump.

After switching on the device, a state change in terms of the position detection of the pressure measuring chamber and irrigation pump is not necessary. The relevant animations cease.

The animation is comprised of the following images (see next page).

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10.9.1 Function/use

Explanation of terms: see pages 17 – 21

Use: Working level for Normal Operation for therapeutic and diagnostic hysteroscopies for the flow modes → 'Continuous Flow' and → 'Single Flow'.

Setting the desired intrauterine pressure: Using the '+/-' buttons, the setpoint of the displayed cavity pressure can be changed. This can be done independently of the pump state (→ activated/→ deactivated).

Blood/Mucus Mode ④: When the pumps are both activated and deactivated, it is possible to switch between Normal Operation and → Blood/Mucus Mode within Continuous Flow procedure using the button.

Once the Blood/Mucus Mode is set, the values which were set in the Pressure/Flow Settings are displayed.

Using the '+/-' buttons, the percentage pressure increase of the intrauterine pressure can be changed.

The pump cannot be switched on whilst the Blood/Mucus Mode button is pressed.

Bipolar Instrument Mode ③: The → Bipolar Instrument Mode button can only be actuated when the pumps are activated. If a HF electrode or similar instrument is used, this function should maintain the pressure in the cavity.

If no electrode or an instrument which is too small are operated with this function, the control system causes this Mode to be exited automatically to prevent the patient from being injured.

① NOTE: HYS sheaths (without working channel) and RES sheaths have no Bipolar Instrument Mode. This means that this button is not enabled for these sheaths.

⚠ WARNING: The Bipolar Instrument Mode is designed for the compensation of the HF instruments being used. Do not use while no HF instruments are in use.

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Footswitch The special modes 'Bipolar Instrument' and 'Blood/Mucus Mode' can be executed either via the touch screen and/or a connected footswitch (5). The following applies for the footswitch:

Grey pedal for Bipolar Instrument, blue pedal for Blood/Mucus Mode.

NOTE: The footswitch symbol is only displayed if a footswitch is actually connected.

Idle Mode: Idle mode is another fully automated additional function for reducing liquid consumption when changing instruments or when removing tissue samples from the body. The Idle Mode is activated when the flow exceeds a value of 320 ml/min for more than a second.

To ensure patient safety, on activation of the Idle Mode the Instrument Mode is automatically deactivated.

In Idle Mode, the irrigation pump works continuously at a rate of 100 ml/min; the suction pump is motionless to prevent air being sucked in.

The device pressure is shown in this phase. After being in Idle Mode for at least five seconds, once a certain pressure has been reached (depending on the sheath), the display switches from Idle Mode to normal automatic operation.

The changed parameters are retained even after deactivation of the pumps, the Blood/Mucus Mode and/or the 'Bipolar Instrument' Mode switch to Normal Operation.

After switching the line voltage off and on again or by jumping back to the Procedure List, these are then replaced with the values of the selected procedure.

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10.9.2 Buttons

Explanation of terms: see pages 17 – 21

- ① By pressing the + button, the desired intrauterine pressure in mmHg is increased. The number '50', which appears here in the content area, was specified by the manufacturer for this specific procedure. The predefined number can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer. Predefinitions: from number 150 for HYS, from number 250 for RES, from number 399 for LAP and from number 550 for SHA.
- ② By pressing the – button the desired intrauterine pressure in mmHg is reduced.
- ③ Bipolar Instr.: By pressing this button, the → Instrument Mode (Cut Mode) is switched on. The button is then shown in orange.
- ④ Blood/Mucus: By pressing this button, the → Blood/Mucus Mode is switched on. The button is then shown in orange.
- ⑤ Footswitch symbol, if a footswitch is connected.
- ⑥ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

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
10.9.3 Displays

- ⑦ Procedure names are comprised of the following 4 elements:
- Procedure icon
 - Procedure number
 - Designation 1
 - Designation 2
- ⑧ Designation of the pressure/increase in pressure. During Normal Operation: Intrauterine pressure. When switching on the Blood/Mucus Mode: Increase in pressure (change in percent)
- ⑨ Symbol for pressure
- ⑩ Pressure (during Normal Operation three digits, in Blood/Mucus Mode two digits with percentage sign)

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10.10 Extended Display GYN

10.10.1 Function and use

The display mode 'Extended Display' which can be enabled via the menu 'Extended Settings' has the same properties and functions as the Normal Mode.

The only difference is that the flow is additionally displayed. This can be adjusted at any time, except with Single Flow sheaths.

10.10.2 Buttons

Explanation of terms: see page 17 – 21

- ① By pressing the + button, the desired intrauterine pressure in mmHg is increased. The number '50', which appears here in the content area, was specified by the manufacturer for this specific procedure. The predefined number can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer (from number 100).
- ② By pressing the – button, the desired intrauterine pressure in mmHg is decreased.
- ③ Bipolar Instr.: By pressing this button, the → 'Bipolar Instrument Mode' is switched on. The button is then shown in orange.

WARNING: The Bipolar Instrument Mode is designed for the compensation of the HF instruments being used. Do not use while no HF instruments are in use.

- ④ Blood/Mucus: By pressing this button, the → Blood/Mucus Mode is switched on. The button is then shown in orange.
- ⑤ Footswitch symbol, if a footswitch is connected.
- ⑥ By pressing this button, the pump is activated (the pump reels move).
- ⑦ By pressing the + button, the desired flow in ml/min is increased.
- ⑧ By pressing the – button, the desired flow in ml/min is decreased.

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10.10.3 Displays

- ⑨ Procedure names are comprised of the following 4 elements:
- Procedure icon
 - Procedure number
 - Designation 1
 - Designation 2
- ⑩ Designation of the pressure/increase in pressure. During Normal Operation: Intrauterine pressure. When switching on the Blood/Mucus Mode: Increase in pressure (change in percent)
- ⑪ Symbol for pressure
- ⑫ Pressure (during Normal Operation three digits, in Blood/Mucus Mode two digits with percentage sign)
- ⑬ Designation of flow
- ⑭ Symbol for flow
- ⑮ Flow in ml/min

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
2E: Confirm sheath!

To start the pump confirm the following settings:

"Bottocchi 4 20152B/BO B4.2"

Operation Mode	Hysteroscopy
Flow Mode	Continuous Flow
Intrauterine Pressure	45 mmHg
Normal Flow	50 ml/min
Blood/Mucus Press., Intr.	80 %
Blood/Mucus Flow	150 ml/min

Do not withdraw instrument from the sheath during operation! This can lead to an increased intrauterine pressure!

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10.11 Sheath and settings confirmation

On pressing the start button, you are asked whether the selected sheath is connected. The set or preset pressure and flow values continue to be displayed.

The user is requested to check the settings and sheath and confirm with 'Yes' or 'No'.

On pressing the 'Yes' button ②, the screen jumps to the operating menu ('Normal Operation' GYN' or 'Extended Display GYN') and the pump is automatically started.

On pressing the 'No' button ①, the screen jumps to the procedures and sheaths ('Procedure List' menu).

WARNING: Do not withdraw instrument from the shaft during operation! This can lead to an increased intrauterine pressure!

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10.12 LAP operating menu

10.12.1 Function/use

Use: Normal working level for the laparoscopic mode. Using the '+/-' buttons, the setpoints for the suction and irrigation flow parameters can be changed. This can be done independently of the pump state (activated/deactivated). Each individual pump can be activated and deactivated independently with its own Start/Stop button. If the footswitch is connected, the symbol is shown on the display. The suction pump can also be started by pressing the gray (left) footswitch pedal without having to press the relevant Start/Stop button on the display (only if the footswitch is connected). The suction pump remains activated until the footswitch or Start/Stop button is pressed again.

10.12.2 Buttons

- ① Decrease the setpoint for irrigation flow, 100–1300 ml/min
- ② Increase the setpoint for irrigation flow, 100–1300 ml/min
- ③ Activate/deactivate the irrigation pump status display (locking function) with On/Off
- ④ Decrease the setpoint for suction flow, 100–1300 ml/min
- ⑤ Increase the setpoint for suction flow, 100–1300 ml/min
- ⑥ Activate/deactivate the suction pump status display (locking function) with On/Off

10.12.3 Displays

- ⑦ Laparoscopy mode
- ⑧ Procedure icon, number, designation and name
- ⑨ Text 'Irrigation flow [ml/min]'
- ⑩ Symbol (irrigation) flow, is always displayed
- ⑪ Flow (four digits, in ml/min)
- ⑫ Text 'Suction flow [ml/min]'
- ⑬ Symbol (suction) flow, is always displayed
- ⑭ Flow (four digits, in ml/min)
- ⑮ Footswitch symbol, if a footswitch is connected

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10.13 Shaver operating menu

① **NOTE:** This mode of operating a KARL STORZ pump together with a KARL STORZ morcellator is not authorized for users inside the U.S.

10.13.1 Function/use

Explanation of terms: see pages 17 – 21

Use: This working level is selected for operation with the UNIDRIVE® S III for work with the shaver.

Setting the desired pressure ①: Using the '+/-' buttons, the setpoint of the displayed cavity pressure can be changed. This can be done independently of the pump state (→activated/→ deactivated).

Setting the desired (suction) flow ②: Using the '+/-' buttons, the setpoint of the suction flow rate can be changed. This can be done independently of the pump state (activated/deactivated).

Footswitch ③: is displayed if the connecting cable between the UNIDRIVE® S III and the HYSTEROMAT E.A.S.I.® has been inserted and the UNIDRIVE® S III is switched on. --> If there is no UNIDRIVE® S III connected, a message appears on the display and the pump cannot be started.

① **NOTE:** The footswitch symbol ③ is only displayed if a footswitch is actually connected.

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10.13.2 Buttons

Explanation of terms: see pages 17 – 21

- ③ By pressing the – button the desired pressure in mmHg is reduced.
- ④ By pressing the + button, the desired pressure in mmHg is increased.
- ⑤ By pressing the – button, the desired (suction) flow in ml/min is reduced.
- ⑥ By pressing the + button, the desired (suction) flow in ml/min is increased.

These predefined parameters can be changed in the 'Settings' menu for self-defined procedures, but not for procedures which have been defined by the manufacturer.

- ⑧ By pressing this button, the pump is activated (the pump reels remove). When the pump is activated, the button is shown in orange. Pressing this button again serves to deactivate the pump.

10.13.3 Displays

- ⑦ Footswitch symbol, if a footswitch is connected
- ⑨ Procedure names are comprised of the following 4 elements:
 - Procedure icon
 - Procedure number
 - Designation 1
 - Designation 2
- ⑩ Designation of the pressure.
During Normal Operation: Pressure.
- ⑪ Symbol for pressure
- ⑫ Pressure (three digits during Normal Operation)
- ⑬ Designation for (suction) flow during Normal Operation: Suction Pump Flow.
- ⑭ Symbol for (suction) flow
- ⑮ (Suction) flow (three digits during Normal Operation)

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10.14 Fluid supply monitoring of the irrigation bag

10.14.1 Function/use

Accessing the menu: This menu is accessed by pressing the button 'Fluid Supply Monitor' (see the orange highlighted irrigation bag at the bottom of the screen).

This menu is only displayed if the 'Fluid Supply Monitor' is activated (enabled) in the menu 'Extended Settings'.

Use: This menu monitors the level of fluid in the irrigation bag. This menu is available regardless of whether the pumps are activated or deactivated. For safety reasons, when the pump is activated the display returns to the previous menu if no button is pressed within 3 seconds. Values which were modified will then be discarded.

The selected irrigation bag size and remaining available fluid volume are displayed. The value for the irrigation bag size can be changed with the +/- buttons (buttons ① and ②). Moreover, the remaining volume value can be modified, if required (buttons ③ and ④), for example if the procedure is started with a partly-used irrigation bag. Using the ⑤ button, the available fluid supply can also be reset to the size of the bag, for example when a new full irrigation bag is connected. All functions can be used regardless of the status of the irrigation and suction pumps (activated/deactivated).

- ⑥ Discard the values and exit the menu
- ⑦ Confirm and save the values and exit the menu

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The modified parameters are retained even if the pumps are deactivated. If the power is switched off and on, however:

- the irrigation bag size is reset to the original value which was set in the menu 'Extended Settings'
- the remaining fluid volume is set to the same value as the irrigation bag size (it is assumed that a new irrigation bag is used for each new intervention)

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10.14.2 Buttons

NOTE: The values for the irrigation bag size are overwritten by any modifications made in the section 'Irrigation bag size' in the 'Extended Settings' menu.

NOTE: Please do not forget to save your settings ⑦.

- ① ② Decrease/increase irrigation bag size setpoint. Possible setpoints: 1 – 3 – 5 – 10.
- ③ ④ Decrease/increase the actual value for available fluid supply.
Range: 'Irrigation bag size' up to 0.0
- ⑤ Reset the available fluid supply to the value of the selected bag size (see input field above) – for example when changing the irrigation bag.
- ⑥ Return to the previous menu level, changes are discarded.
- ⑦ Return to the previous menu level, changes are saved.

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11.1 Basic 'Settings' menu

Accessing the menu

This menu is reached by pressing the 'Settings' button (2).

Explanation of terms: see pages 17 – 21

11.1.1 Function/use

This menu is the basic menu (main menu) in the menu structure 'Settings'. It is the starting point for the submenus in which settings for the HYSTEROMAT E.A.S.I.® can be made and displayed.

11.1.2 Buttons

- ① Shortcut to the menu 'Pressure/Flow Settings'. This is where you enter the initial values for the 'operating menu'.
- ② Shortcut to the menu 'Extended Settings'. In this menu, you can:
 - activate the Fluid Supply Monitor
 - set the standard value for the irrigation bag size
 - activate the use of → Day sets
 - activate the 'Extended Display GYN' to allow not only the pressure but also the flow to be set
- ③ Shortcut to the 'Device Settings' menu. In this menu, you can:
 - set the language
 - set the display brightness
 - set the volume for the information and alarm signals
 - set the volume for button actuation
 - proceed to the menu 'Date and Time Settings'

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④ Shortcut to the menu 'Device Configuration'. In this menu, the following information is displayed:

- serial number
- operating hours
- user interface software version
- software process control

⑤ Shortcut to the menu 'System Log'. This is where the information and alarm signals are recorded.

⑥ Shortcut to the menu 'Service Programs'. Access only with password for authorized service personnel.

⑦ Shortcut to the menu 'Procedure Administration'. In this menu, you can:

- save predefined procedures (e.g. 100 Cont. Flow HYS B.I.O.H.) under another name and set different values for the newly created procedure (e.g. different pressure, different sheath)
- rename procedures (excluding the procedures predefined on delivery)
- delete procedures (excluding the procedures predefined on delivery)
- overwrite procedures with other values (excluding the procedures predefined on delivery)

⑧ Return to the previous menu level.

⑨ Return to the previous menu level. Saving the settings.

① **NOTE:** When you switch to the menu level 'Operation Selection' the buttons ① and ⑦ are not available.

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11.2 Pressure/flow settings in GYN

11.2.1 Function/use

Input menu for the pressure and flow parameters:


- Specification of setpoint for intrauterine pressure.
Example: In this menu, a value of 40 mmHg is specified. This value of 40 is then saved and displayed in the 'GYN operating menu'. Here, the value for the current operation can be changed at any time.
- Specification of setpoint for normal flow.
Example: In this menu, a value of 40 ml/min is specified. Irrigation with a flow of 40 ml/min is then achieved in the 'GYN operating menu'. Here, the value for the current operation can be changed at any time.
- Specification of setpoint for pressure increase in Blood/Mucus mode.
Example: If a value of 20% is specified in this menu, a pressure increase of 20% in addition to the set intrauterine pressure (e.g., 120 mmHg) will be realized in the 'GYN operating menu', but only if the 'Blood/Mucus mode' button is pressed (this is only possible in the flow mode 'Continuous Flow', not in 'Single Flow'). Instead of 120 mmHg the user would then work with a value of 144 mmHg ($120 + 20\% \text{ of } 120 = 120 + 24$). The percentage increase in pressure which is specified in this menu can still be changed at any time for the current operation in the GYN operating menu.
- Specification of setpoint for Blood/Mucus mode flow.
Example: If a value of 150 ml/min is specified in this menu, the irrigation flow rate will be 150 ml/min during pump operation, but only if the 'Blood/Mucus mode' button was pressed in the 'GYN operating menu' (this is only possible in the 'Continuous Flow' flow mode).

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

11.2.2 Buttons

- ① Decrease the specified setpoint for the intrauterine flow.
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) which is to be displayed in the 'Intrauterine pressure' field of the operating menu 'GYN' is specified. In the menu 'GYN', this value can still be changed specifically for the operation at any time.
- ② Increase the specified setpoint for the intrauterine flow.
- ③ Decrease the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.
- ④ Increase the specified setpoint for the percentage increase in pressure in the Blood/Mucus Mode.
- ⑤ Decrease the specified setpoint for the normal flow.
Here 'normal flow' refers to the flow associated with Normal Operation and not the flow seen in the Blood/Mucus Mode.
- ⑥ Increase the specified setpoint for the normal flow.
- ⑦ Decrease the specified setpoint for the flow in the Blood/Mucus Mode.
- ⑧ Increase the specified setpoint for the flow in the Blood/Mucus Mode.
- ⑨ Return to the previous menu level; the changed values are discarded.
- ⑩ Return to the previous menu level; the changed values are saved.

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11.3.1 Function/use

Input menu for the flow parameters:

- Specification of setpoint for irrigation flow.
Example: In this menu, a value of 500 ml/min is specified. This value of 500 is then saved and displayed in the operating menu 'LAP'. The specified value shown here can be changed specifically for the operation at any time.
- Specification of setpoint for suction flow.
Example: In this menu, a value of 500 ml/min is specified. A suction flow of 500 ml/min is then achieved in the operating menu 'LAP'. The specified value shown here can be changed specifically for the operation at any time.

11.3.2 Buttons

- ① Decrease the specified setpoint for the irrigation flow.
Explanation of 'specified setpoint': In this 'Pressure/Flow Settings' menu, the desired value (setpoint) which is to be displayed in the 'Irrigation flow' field of the operating menu 'LAP' is specified.
- ② Increase the specified setpoint for the irrigation flow.
- ③ Decrease the specified setpoint for the suction flow.
- ④ Increase the specified setpoint for the suction flow.
- ⑤ Return to the previous menu level; the changed values are discarded.
- ⑥ Return to the previous menu level; the changed values are saved.

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
11.4 Pressure/flow settings in SHA

11.4.1 Function/use

Input menu for the pressure and flow parameters:
Depending on the operating mode selected, different parameters can be selected individually.

- Specification of setpoint for irrigation pressure ①.
Example: In this menu, a pressure value of 75 mmHg is specified. This setpoint is saved and displayed in the 'SHA operating menu'.
- Specification of setpoint for suction flow ②.
In this menu, a flow value of 140 ml/min is specified. This setpoint is saved and displayed in the 'SHA operating menu'.

Saving changes permanently: Changes which are made here are only retained for as long as no new procedure is selected. If the settings should be retained for longer, the values must be saved in a new procedure (see section 10.7).

Overwriting the values of the operating menu: Any previously changed settings in this operating menu are overwritten upon confirmation of the settings with the button ④.

11.4.2 Buttons

- ⑤ Decrease the specified setpoint for pressure.
Explanation of 'specified setpoint': In this 'Pressure/flow settings' menu, the desired value (setpoint) for the pressure/flow which is to be displayed in the 'Pressure'/'Flow' field of the 'SHA' operating menu is specified.
- ⑥ Increase the specified setpoint for pressure
- ⑦ Decrease the specified setpoint for flow
- ⑧ Increase the specified setpoint for flow
- ③ Return to the previous menu; the changed values are discarded
- ④ Return to the previous menu; the changed values are saved

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11.5 'Extended Settings' menu

11.5.1 Function/use

This is an entry menu for further settings.

① **NOTE:** Do not forget to confirm your entries by pressing button ①.

11.5.2 Buttons (orange)

- ① Here you can call up the menu which serves to monitor the irrigation bag fluid supply.
- ② ③ Decrease/increase the initial value of the irrigation bag size.
The irrigation bag size is shown in the 'Fluid Supply Monitor' menu if this function has been selected with the button ①.

① **NOTE:** If the initial values are confirmed by pressing the button ⑥, the values set in the 'Fluid Supply Monitor' menu are overwritten!

- ④ If you press this button, the menu 'Extended Display GYN' is enabled. For a GYN procedure, this allows you not only to specifically set the pressure but also the flow for the current operation.
- ⑤ Enable the 'Day Set' function.
- ⑥ Return to the previous menu level and discard the selected values.
- ⑦ Return to the previous menu level and save the selected values.

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11.6.1 Function/use

This is the menu for setting

- ① ② the language displayed on-screen,
- ③ ⑥ monitor brightness,
- ③ ④ the volume of information/alarm signals

WARNING: Make sure that the set volume of the loudspeakers allows for perception of the acoustic signals despite normal background noise.

- ⑦ ⑧ volume of button actuation

In addition, pressing the ⑨ button allows you to switch to the settings menu for date/time. Any modified values that have not been saved are then discarded.

- ⑩ Return to the basic 'Settings' menu. The modified values are discarded.
- ⑪ Return to the basic 'Settings' menu. The modified values are saved.

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11.7 'Date and Time Settings' menu

11.7.1 Function/use

The setting of date and time is self-explanatory.

- ① Return to the basic 'Settings' menu.
The modified values are discarded.
- ② Return to the basic 'Settings' menu.
The modified values are saved.

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11.8 'Procedure Administration' menu

11.8.1 Function/use

- ① By pressing this button, you can overwrite an existing procedure that you defined yourself. Procedures which were defined on delivery cannot be overwritten. This means that the button does not react when you click it for a predefined procedure.
- ② By pressing this button, you can delete an existing procedure that you defined yourself. Procedures which were defined on delivery cannot be deleted. This means that the button does not react when you click it for a predefined procedure.
- ③ Save as New Procedure. If you have changed the values of a predefined procedure and want to have these changed values available permanently, we recommend creating a new procedure. To do so, press this button.
- ④ Rename Procedure. Self-defined procedures can be renamed.

Example: Defining a new procedure and setting other values yourself

Our example has the following objectives:

- A procedure with the sample name '100 Cont. Flow HYS B.I.O.H.' is to be copied
- The following value should be changed:
– Pressure: from 50 mmHg to 40 mmHg
- The new procedure should be saved under the name '160 Cont. Flow HYS B.I.O.H. 40mmHg' under the number '160'.

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Method for saving a new procedure

- Switch the device on and wait a few seconds until the 'Procedure List' appears
- Select an operation type
- On the touch screen, click on the procedure '150 Cont. Flow HYS B.I.O.H.'. This procedure is now highlighted in orange. This means the procedure is now activated
- Click on the tool symbol in the bottom left of the screen
- (The 'Settings' menu appears)
- Click on the 'Procedure Administration' submenu
- Click on 'Save as New Procedure'
- In the orange highlighted field, enter 'Cont. Flow HYS B.I.O.H. 40 mmHg'
- Press the TAB key (the large button on the far bottom left). The second input field is then highlighted orange
- Leave this field empty
- Press the TAB key again (the large button on the far bottom left). The third input field is then highlighted orange
- Enter the identification number for the new procedure here, i.e., '160'
- Click on the check mark in the far bottom right to save your settings
- (The updated procedure list appears. The new procedure 'Cont. Flow HYS B.I.O.H. 40 mmHg' appears at the top and is highlighted in orange, i.e. it is activated)

Method for changing values

- Click on the tool symbol in the bottom left of the screen
- (The 'Settings' menu appears)
- Click on 'Pressure/Flow Settings'
- Under the header 'Pressure', click twice on the minus sign to reduce the value from 50 to 40

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- Click on the check mark in the bottom right to save your settings
- (You are once again in the 'Settings' menu)
- (The procedure '160 Cont. Flow HYS B.I.O.H. 40 mmHg' has now been permanently changed)
- (The updated procedure list appears. The new procedure '160 Cont. Flow HYS B.I.O.H. 40 mmHg' appears at the top and is highlighted in orange, i.e. it is activated)

Method for renaming the procedure

- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Rename Procedure'
- The flashing character (cursor) is now on the left, next to the '1' of '160'
- Click on the 'arrow-to-the-right button' in the bottom, right-hand area until the cursor is to the right of '40'
- Enter a space. The procedure is now called '160 Cont. Flow HYS B.I.O.H. 40 mmHg'
- Save your entries by clicking on the check mark

Method for deleting the procedure

- (You are now in the Procedure List)
- Click on the tool symbol on the bottom left
- (The 'Settings' menu appears)
- Click on 'Procedure Administration'
- Click on 'Delete Procedure'
- Click on the check mark on the bottom right to delete the procedure permanently

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
11.9 'Device Configuration' menu

Function/use

The following elements are displayed:


- Device serial number
- Operating hours
- User interface software version
- Process control software version

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11.10 'System Log' menu

Function/use

Presentation of the saved information and alarms.

Each entry occupies one row and consists of a sequential number, date, time and alarm/information ID.

The most recent entry is allocated the number 0.

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
11.11 'Service Programs' menu

Function/use

The service programs are only available to service personnel. Access is password-protected.

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Confirm increase of intrauterine pressure!

You are going to increase intrauterine pressure over 60 mmHg.

Confirm to increase intrauterine pressure over 60 mmHg.

Click X to continue with intrauterine pressure of max 60 mmHg.

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
11.12 'Confirm increase of pressure' menu

- If the pressure in the activated or deactivated pump status, i.e., also in the 'Pressure/Flow Settings' is increased for the first time to 65 mmHg, this results in a jump to this menu level, where a reference to the increase in pressure appears.
- This increased pressure is then confirmed once for the relevant procedure using the confirm button (2) before the screen jumps back to the previous menu level.
- With the discard button (1), the pressure remains at 60 mmHg and the screen jumps back to the previous menu level. If the pressure is increased again in excess of 60 mmHg, the confirmation process is repeated.

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11.13 'Safe State' menu

- The Safe State is set for certain causes of risks. In this state, the device can no longer be operated normally.
- However, in the 'Settings' area you can access another examination.
- This menu level is a special case. Switching off the line voltage ends the Safe State.

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



12 Service and repair

12.1 Fuse replacement

Turn the device off and remove the power cord from the electrical outlet.

WARNING: The line voltage to the device is only definitively disconnected once the power plug has been unplugged.

Loosen the line fuse holder with a screwdriver or other suitable tool. Remove the line fuse holder.

WARNING: Only use fuses of the correct rating. Replace defective fuses and screw the line fuse holder back on.

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12.2 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

WARNING: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

WARNING: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

WARNING: Always pull out power plug before cleaning!

WARNING: Avoid allowing liquids to enter the housing at all costs.

The design allows only for manual wipe-down disinfection. The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used.

CAUTION: When cleaning the device surface, the area around the pressure sensor ⑥ should be wiped carefully due to its sensitive surface.

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12.3 Maintenance and safety check

12.3.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life. Information regarding maintenance services can be obtained from your local representative or from the manufacturer.

12.3.2 Safety checks/ repeat inspections as per IEC 62353

WARNING: Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Patient leakage current measured according to IEC 62353: The limit values can be found in the current standard.

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Test for proper functioning

A functional check must be carried out in accordance with the instruction manual (see Section 8.4 'Operating the unit').

Documentation

The safety check/repeat inspection and results must be documented.

WARNING: If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

NOTE: For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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12.4 Service and repair

WARNING: Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ genuine KARL STORZ parts only.

12.5 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

WARNING: At the end of its service life, dispose of the device as electronic waste.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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12.6 Repair program

During the repair period, you will generally be loaned a unit, which must then be returned to KARL STORZ as soon as the repaired unit is returned to you.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

In other countries please contact your local KARL STORZ branch or authorized dealer.

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12.7 Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, modifications, or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no guarantee of the correct functioning of devices or instruments which have been repaired by unauthorized third parties.

When shipping the device, use either the original packaging or equivalent packaging marked with the symbols illustrated here.

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12.8 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet the applicable laws and standards, and
- the device has been used in accordance with its operating instructions at all times.

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12.9 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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13.1 Information signals

13.1.1 Configuration

All the information signals created, with the exception of those which display the cause for an inoperable device, (entry 'Safe State' in the 'Reaction' column of the table from page 91) are not self-retaining. This means that the signal is only issued whilst the signal conditions are present (unless something different is defined under 'Reaction').

If the signal conditions are only very brief, there are regulations as to the minimum visual/acoustic signal duration to which the described information signals also have to comply. This is explained in more detail in the following chapters.

Signals are issued depending on priority. Any low priority signal will be overwritten by an alarm signal with a higher priority or a low priority information signal is suppressed so long as there is a signal with a higher priority. If there are several signal conditions of the same priority, only the condition which was recognized first will be signaled. Should signals deviate from this rule, this is noted in the corresponding entry from page 91.

13.1.2 Visual signals

Information signals are displayed inversely (blue font on a white background) compared to normal displays and do not flash.

The visual signal is issued for as long as the signal conditions apply. In the case of signal conditions of a very short duration, it may be difficult for the user to recognize which condition triggered the signal. For this reason, a minimum display duration of 5s is specified.

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Technical description

Specification for acoustic information signals

13.1.3 Acoustic signals


The acoustic signals described below satisfy the requirements in terms of distinctiveness.

The acoustic signal is issued for as long as the signal conditions apply. An entire sequence is played through as a minimum.

The acoustic signal can be switched on or off using a button in the footer of the device display.

The acoustic signal, however, can only be switched off temporarily by the user. After switching on the line voltage, the setting is always 'Audio Alarm On'.

Information signals (general)

The pitch is modulated by ± 2 Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics. In the case of messages which are related to the Safe State, the sequence is repeated every 15 s, but otherwise just played through once.

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
Availability signal


After completion of the self-test after switching on, an availability signal sounds. This signal, which is issued following a successful self-test, has the following characteristics:

The pitch is modulated by +/- 2 Hz with a modulation frequency of 1.5 Hz. There are 5 different harmonics.

Blood/Mucus period

The following periodic signal is issued on starting the Blood/Mucus Mode. This signal is repeated every 15 seconds for the duration of the Blood/Mucus Mode.

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

On ending the Blood/Mucus Mode, the following signal sounds:

Low irrigation fluid supply

If the Fluid Supply Monitor is activated and the irrigation fluid supply is low, a double gong is produced. It lasts for 1.5 seconds.

Key click

Pressing a key on the touch screen is confirmed by a short beep. The volume of this beep can be set separately from the volume of all other information signals in the 'Device Settings' menu and can also be permanently deactivated.

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13.2 Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
20: Dual-pedal footswitch recognized Dual-pedal footswitch was plugged in	3	None
21: Dual-pedal footswitch disconnected Dual-pedal footswitch was disconnected	3	End footswitch special function
22: Shaver system recognized The shaver system is connected and switched on	3	None
23: Shaver system disconnected The shaver system was switched off or disconnected	3	End shaver special function
24: Footswitch/Shaver not supported The current operating context has not assigned a function to the actuated foot- or shaver switch	3	None
25: Blood/Mucus mode period expired The Blood/Mucus Mode ended after the set period of time, the footswitch pedal is still actuated	3	None
26: Low irrigation fluid supply The irrigation fluid supply is coming to an end (< 200 ml)	1	None. Message remains until the irrigation fluid supply is above the limit value or the fluid supply monitor is deactivated.
27: Pressure dome detached Pressure dome was detached	2	Cancel activation (if activated)
28: Irrigation pump unlocked Irrigation pump locking lever was opened	2	Cancel activation (if activated)
29: Suction pump unlocked Suction pump locking lever was opened	2	Cancel activation (if activated) Single Flow: no reaction
2A: Illegal password Incorrect password entered	3	None
2B: Password entry mismatch On entering/confirming a new password, different entries were made	3	None
2C: Pressure dome attached On turning on the line voltage, the pressure dome was attached but the day set was not enabled	2	Menu change not executed, message remains until the pressure dome is removed
2D: Invalid procedure number The procedure number which has been selected for saving has already been used	3	None

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
40: Contact KARL STORZ Service Voltage of the SRAM back-up battery on powering up < 2.7 V	2	None
41: Contact KARL STORZ Service Voltage of the SRAM back-up battery on powering up < 2.3 V	1	None
42: Reset of settings/data The settings/data were reset manually or due to an error (one or more saving areas)	2	None
43: Blood/Mucus mode duration exceeded Current Blood/Mucus mode duration exceeds the setpoint by more than 5 s	2	Cancel activation
44: Perform service On switching on the line voltage incompatible software versions were found by the User Interface and Process Control	1	Safe State
45: Activated ext. switch during power-up On switching on the line voltage a pedal on the footswitch or a button on the shaver handpiece was actuated	1	Switch function is not executed, message remains until the switch element is released
46: Inserted external switch is activated On insertion a pedal on the footswitch or a button on the shaver handpiece was actuated	1	Pedal function is not executed, message remains until the switch element is released
47: Touch screen short circuit On switching on the line voltage the touch screen was pressed	1	Safe State
48: Touch screen short circuit Maximum touch screen actuation time of 40 s exceeded	1	Safe State
4A: Contact KARL STORZ Service On switching on the line voltage a faulty component in the drive system was detected	1	Safe State
4B: Contact KARL STORZ Service On switching on the line voltage a faulty component in the drive system was detected	1	Safe State

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
50: Attach pressure dome again The value of the exposed sensor is outside the tolerance of 24.4 mV – 427 mV at the sensor output (approx. 4–70 mmHg)	1	Cancel activation. Pressure dome animation Renewed activation is only possible once the pressure dome has been removed and replaced
53: Attach pressure dome again A pressure is measured by the sensor which could result in the latter's destruction	1	Cancel activation. Pressure dome animation Renewed activation is only possible once the pressure dome has been removed and replaced
60: Internal temperature too high Temperature inside the device is greater than 55°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 61
61: Critical internal temperature Temperature inside the device is outside the 10°C – 70°C range	2	Cancel activation, message remains until the temperature falls below the limit value
62: Irrigation motor temperature too high Temperature > 70°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 63
63: Irrigation motor temperature critical Temperature > 100°C or disconnection of a sensor line	2	Cancel activation, message remains until the temperature falls below the limit value
64: Suction motor temperature too high Temperature > 70°C	3	None, message remains until the temperature falls below the limit value or is superseded by message 65
65: Suction motor temperature critical Temperature > 100°C or disconnection of a sensor line	2	Cancel activation, message remains until the temperature falls below the limit value

Information messages

Brief description (1 line with approx. 40 characters possible) Possible cause	Priority	Reaction
70: Reposition pressure dome in < 4 h More than 20 hours have elapsed since the last sensor test	2	None, message remains until the pressure dome is detached
71: Attach pressure dome again The offset value of the pressure sensor is no longer up-to-date (pressure dome must be detached and reattached)	1	Activated: Activation remains Message remains Not activated (or immediately after deactivation): Pressure dome animation Activation and deletion of the message is only possible once the pressure dome has been removed and replaced
72: Restart device in < 4 h Since the last self-test over 20 hours ago, the device has been in continuous operation	2	None, message remains until the device is switched off and on again
73: Restart device Since the last self-test over 24 hours ago, the device has been in continuous operation and a new self-test is now necessary	1	Activated: Activation remains Message remains Not activated (or immediately after deactivation): Safe State
80: Contact KARL STORZ Service Supply voltage(s) outside the specification	1	Safe State
81: Contact KARL STORZ Service Supply voltage below 15 V	1	Safe State
90: Unable to download a procedure An attempt was made to download a procedure from the RUI whilst the device was not ready-to-receive.	2	None
91: Invalid procedure download data The data sent from the RUI are faulty	2	None
FD: Internal system error I An unexpected system state is detected	3	None
FE: Internal system error II An unexpected system state is detected	2	Cancel activation
FF: Internal system error III An unexpected system state is detected	1	Safe State

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- Type of protection against electric shock: Protection Class I
- Degree of protection against electric shock: Applied part type BF
- Type of moisture protection: Drip-proof as per IPX1

**13.5 Directive compliance
(for 26340020-1)**

This medical product belongs to Class IIb.

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

NOTE: The code number after the CE mark indicates the responsible notified body.

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13.6 Technical documentation

On request, the manufacturer will provide available circuit diagrams, detailed parts lists, descriptions, instructions for adjustments, and other documentation deemed useful for suitably qualified user personnel authorized by the manufacturer when repairing device components that have been designated as repairable by the manufacturer.

Supply of such technical documentation relating to the device shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the device.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

NOTE: We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
13.7 Blockschaltbild

13.7 Block diagram

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**Spare parts,
recommended accessories**

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**14 Spare parts,
recommended accessories**

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Item	Order no.
Tube sets* — for suction* — for irrigation*	031217-10 031717-10
Day tube set for irrigation* Pump tube set with two puncture needles. Sterile, pack of 10.	031767-10
Patient tube for single use* For use with pump tube sets Sterile, pack of 10.	031162-10
Footswitch	26 3403 30
UNIDRIVE® SIII	20 7010 20-1
Control cable	20 7010 70

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Electromagnetic compatibility (EMC) information

15 Electromagnetic compatibility (EMC) information

WARNING: Medical electrical equipment needs special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this section during installation and operation.

The HYSTEROMAT E.A.S.I.® model 26340020-1 corresponds to IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HYSTEROMAT E.A.S.I.® model 26340020-1 is a Group 1 device (as per CISPR 11).

Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this section provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the device or system to perform its intended use without causing interferences in other devices or systems or non-medical electrical devices. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the devices to different electrical circuits.

If you have any further questions, please contact your local representative or our service department.

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Electromagnetic compatibility (EMC) information

WARNING: The HYSTEROMAT E.A.S.I.® model 26.3400.20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this or other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the electronic HYSTEROMAT E.A.S.I.® model 26.3400.20-1. The accessories and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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Accessories and cables which have been shown to comply with IEC 60601-1-2:

Cable type	Shielded	Length [m]	Ferrite	Used for
MKF 4 E-4396	Yes	4	No	Footswitch
Power Cord	No	2	No	Power Supply
SCB Cord	Yes	1.0 m	No	Connecting cable
PE	No	> 3.0 m	No	Potential equalization

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Table 1 Guidance and manufacturer's declaration – electromagnetic emissions		
HYSTEROMAT E.A.S.I.® model 26340020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the HYSTEROMAT E.A.S.I.® model 26340020-1 should make sure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidelines
RF emissions CISPR 11	Group 1	The HYSTEROMAT E.A.S.I.® model 26340020-1 uses only RF energy for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The HYSTEROMAT E.A.S.I.® model 26340020-1 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2 Guidance and manufacturer's declaration – electromagnetic immunity			
HYSTEROMAT E.A.S.I.® model 26340020-1 is intended for use in the electromagnetic environment specified below. The customer or user of the HYSTEROMAT E.A.S.I.® model 26340020-1 should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidelines
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst acc. to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5 % U_T * (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Complies <5 % U_T * (>95 % dip in U_T) for 0.5 cycle Complies 40 % U_T (60 % dip in U_T) for 5 cycles Complies 70 % U_T (30 % dip in U_T) for 25 cycles Complies <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 4

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

HYSTEROMAT E.A.S.I.® model 26340020-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the HYSTEROMAT E.A.S.I.® model 26340020-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the HYSTEROMAT E.A.S.I.® model 26340020-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 6 Recommended separation distances between portable and mobile HF communications equipment and the HYSTEROMAT E.A.S.I.® model 26340020-1			
The HYSTEROMAT E.A.S.I.® model 26340020-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile HF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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KARL STORZ Endoscopy Canada Ltd
7171 Millcreek Drive, Mississauga
Ontario, L5N 3K3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseia.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseia.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuha@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 Mexico D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

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KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Lingang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang1, 102dong, 8F
239 Jeongjae-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp
www.karlstorz.com

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
E-Mail: info@karlstorz.ca

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100,
Toll free: 800 421-0837 (USA only)
Fax: +1 424 218-8526
Toll free: 800 321-1304 (USA only)
E-Mail: info@ksa.com

KARL STORZ Veterinary Endoscopy America, Inc.
175 Cremona Drive
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia
Miramar Trade Center
Edificio Jerusalem, Oficina 108,
La Habana, Cuba
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia Mexico S.A. de C.V.
Lago Constanza No 326, Col. Chapultec Morales
D.F.C.P. 11520, Mexico, Mexico
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing America do Sul Ltda.
R. Joaquim Floriano, 413, 20º andar -
Itaim Bibi, Sao Paulo, -SP, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufiategui 627 6º Piso
B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
P.O. Box 153, Rolf Olsenvei 28
N-2007, (illegible), Norway
Phone: (illegible)

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12724 Skårholmen, Sweden
Postal address: Po Box 8013,
14108 Kungens Kurva, Sweden
Phone: +46 8 50 56 4800, Fax: +46 8 50 56 4848
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Valimotie 1B, 4. Kerros, 00380 Helsinki, Finland
Phone: +358 96824774, Fax: +358 96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
392 Edinburgh Avenue, Slough
Berkshire, SL1 4UF, Great Britain
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: customerservice@karlstorz-uk.com

KARL STORZ Endoscopia Nederland B.V.
Phone: +31 651 938 738, Fax: +31 135 302 231

KARL STORZ Endoscopy Belgium N.V.
Phone: +32 473 810 451, E-Mail: info@stopler.be

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S.A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: marketing@karlstorz.es

KARL STORZ Endoscopia Italia S.r.l.
Via dell' Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info@karlstorz.it

KARL STORZ Adria Eos d.o.o.
Zadarska 80
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

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KARL STORZ ENDOSKOPE East Mediterranean and Gulf
Building 25 (Al-Hayyan) - Units 205-206
Dubai Health Care City, P.O. Box 118069
Dubai United Arab Emirates
Further - illegible

KARL STORZ Endoscopy India Private Ltd.
D-181, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F.10 - Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
PO Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch Company,
Room 610, China Life Tower No. 6, Chaowai Street,
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 2000300, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch Company
F-5, 24/F, Chuanxing Mansion,
No. 18 Renmin Road South
Chengdu, Sichuan, 610016, People's Republic of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch Company,
Rm 2225, Tower B, City Plaza,
No. 83, Zhongshan Road, Heping District,
Shenyang, Liaoning, 110001,
People's Republic of China
Phone: +86 24 6258 9911, Fax: +86 24 6258 9922
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch Company
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte. Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
Jeljon Tower 1, Suite 905
16 Neuti-Ro, Bundang-Gu
Seongnam, 463-847, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: Stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
6F-1, No. 10, Sec. 1, Beism Rd.
Sindian, Taipei County 231, Taiwan (R.O.C.)
Phone: +886 933 014 160
Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 (0)3 6380-8633
E-Mail: info@karlstorz.co.jp

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