User Manual

PicoStar

Dr. Thomas Weikart
Civil Law Notary

18 Juli 2020

Jene.

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Asclepion Laser Technologies GmbH

Bruesseler Str. 10
D - 07747 JENA
Germany

Phone: +49 (0) 3641 / 7700 - 401

Fax: +49 (0) 3641 / 7700 - 402

E-mail: service@asclepion.com

CE 0123

1	INTRODUCTION	1
1.1	GENERAL INFORMATION	1
1.2	SYMBOLS IN THE USER MANUAL	1
1.3	COPYRIGHT	2
2	INTENDED USE	3
2.1	PROPER USE	3
2.2	CONTRAINDICATIONS	3
2.3	POSSIBLE SIDE EFFECTS	4
2.4	REQUIREMENTS FOR THE USER	4
2.5	INTENDED USE ENVIRONMENT	4
2.6	MEDICAL DEVICE DIRECTIVE / MEDICAL PRODUCT ACT	4
3	SAFETY	5
3.1	GENERAL SAFETY INFORMATION	5
3.2	ELECTRICAL HAZARDS	5
3.3	BIOLOGICAL HAZARDS	6
3.4	RADIO INTERFERENCE	6
3.5	LASER	7
3.5.1	<i>Optical hazards caused by laser emission</i>	7
3.5.2	<i>Laser-induced fire hazard</i>	9
4	DEVICE DESCRIPTION	10
4.1	DEVICE DESIGN	10
4.2	SAFETY EQUIPMENT	11
4.3	LABELS AT THE DEVICE	13
5	INSTALLATION	17
5.1	SCOPE OF DELIVERY	17
5.2	GENERAL INFORMATION ABOUT INSTALLATION	18
5.3	DOOR INTERLOCK CONNECTION	19
5.4	POWER GRID CONNECTION	20
5.5	POTENTIAL EQUALIZATION LINE	20
5.6	FOOT SWITCH CONNECTION	20
5.7	SMARTCARD	20
5.8	HANDPIECE CONNECTION	21
6	TECHNICAL DATA	23
6.1	DEVICE MODELS AND APPLIED PARTS	23
6.2	SPECIFICATIONS	23
7	OPERATION OF THE DEVICE	27
7.1	SWITCHING ON	27
7.2	GENERAL INFORMATION ON TOUCH SCREEN OPERATION	28
7.3	MAIN MENU	29
7.4	SETUP	29
7.4.1	<i>Import/Export data</i>	30
7.4.2	<i>Software/ System information</i>	31
7.5	STARTING THE TREATMENT	31
7.6	STARTING VIA LASER	32
7.6.1	<i>Saving of parameters</i>	32
7.7	STARTING VIA INDICATION	33
7.7.1	<i>Indication list: Tattoo</i>	34

7.7.2	Indication list: Make-up	34
7.7.3	Indication list: Pigment	34
7.7.4	Emission of laser beam (READY mode)	35
7.7.5	Check of the beam delivery system	36
7.8	SWITCHING OFF THE DEVICE	37
8	CLEANING, DISINFECTION AND STERILIZATION	38
8.1	DEVICE	38
8.2	HANDPIECE	38
9	FAULT MESSAGES AND TROUBLESHOOTING	39
9.1	GENERAL INFORMATION	39
10	CUSTOMER SERVICE	41
11	REGULAR MAINTENANCE, SAFETY CHECKS AND CALIBRATION	42
11.1	ROUTINE MAINTENANCE	42
11.2	REGULAR SAFETY CHECK	43
12	DISPOSAL	44
13	EC DECLARATION OF CONFORMITY	45
	APPENDIX A: EMC MANUFACTURER'S DECLARATION	46

1 Introduction

1.1 General information

- This user manual contains the information required for the intended use of the device.
- It is an integral part of the complete product according to the applicable national, European and international directives and laws on product liability.
- The most important objective is to protect persons against risk situations and the device against damages caused by improper use.
- Regular maintenance according to the instructions given in this manual increases the precision and operability of the device during its entire service life.
- It contains the contact data of the Technical Services of Asclepion Laser Technologies GmbH for faults and errors, which go beyond the capabilities of the local representative.

Knowledge of this user manual is required for the operation of the devices. Therefore, please make yourself familiar with the contents and pay special attention to information concerning the safe operation of the device.

1.2 Symbols in the user manual

The following symbols are used to inform of the residual risks due to any shortcomings of the protection measures adopted or to other possible dangers.

WARNING

This symbol indicates a possible danger for an individual's life and health. Non-compliance with these instructions can cause serious danger to health, including life-threatening injuries.

CAUTION

This symbol indicates a possible dangerous situation. Failure to obey these instructions can cause injuries.

The following symbols are used to give instructions on operating the device:

Important information

This symbol indicates important information about the proper use of the device. Failure to obey these instructions can cause damage to property or malfunctions of the device or its surroundings.

Tip

This symbol indicates practical tips and particularly useful information that enables the user to use all the features of the device to best effect.

1.3 Copyright

The specifications are subject to change due to further technical development. Contact the Technical Service of Asclepion Laser Technologies GmbH or your local distributor of Asclepion Laser Technologies GmbH to get latest information.

© Unless expressly authorized, the forwarding and duplication of this document, and the utilization and disclosure of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

2 Intended use

2.1 Proper use

WARNING

The incorrect using of this device can cause serious danger to health. Therefore, all warnings, instructions and information of this user manual are to be observed under all circumstances! The non-compliance is considered as unintended use of the device!

No modification of this equipment is allowed. The non-compliance is considered as unintended use of the device! Thereby the warranty, ensured by Asclepion Laser Technologies GmbH, expires.

The PicoStar is applied for vaporization and ablation of soft tissue, particularly for a selective damage of pigments in the epidermis and dermis and thus the gentle treatment of benign pigmented lesions and artificial lesions. Typical applications also can be:

- Age spots
- Aging skin
- Becker's nevus
- Benign pigmentation
- Café au lait
- Dyschromia
- Hori's macules
- Lentigines
- Melasma
- Nevus Ota
- Nevus spilus
- Permanent Make-up
- Pigmented Lesions
- Scars/Acne
- Striae
- Tattoo
- Wrinkles

2.2 Contraindications

The contraindications listed below stem from many years of clinical experience with laser devices for this intended use. They do not claim to be complete or to have unlimited validity.

- topical therapy with retinoids, imiquimod, 5-fluorouracil, ingenol mebutate, diclofenac, alpha-hydroxy acids, or salicylic acids within the previous month (for rejuvenation)
- pregnant and breastfeeding women (due to a lack of studies, available for this study group)
- active infections in treatment area
- Hypersensitivity to light / Photosensitivity
- History of skin cancer, keloid scarring
- Immunodeficiency
- History of resurfacing procedures such as laser, medium, or deep chemical peel within the previous 6 months

- intense pulsed light or superficial chemical peel within the previous 1 month
- active medications for which sunlight is a contraindication
- history of squamous cell carcinoma or melanoma
- inflammation of the tattooed area

2.3 Possible side effects

- Little discomfort and pain
- Pinpoint bleeding
- Swelling
- Erythema/ Redness/ Edema/ Blisters / Bulla formation
- Hypopigmentation/ Dyspigmentation
- Purpura/ Bruising
- Burning and stinging
- Urticarial reaction
- Crusting/ textural changes in treatment area/ slight fibrosis in treated area
- Scabbing
- Scarring
- Post inflammatory pigmentary changes
- Itching

2.4 Requirements for the user

WARNING

Treatment must only be carried out by a doctor or under the supervision of a doctor (simply referred to as the user in the following).

Active medical products may only be applied by users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience.


It is expressly pointed out that the use of the device is only permitted for users that have been instructed to operate it and confirmed the training course by their signature in the Medical Device Logbook if it is prescribed by a state law.

It is NOT PERMITTED to use laser treatment for pregnant women or children because clinical studies do not exist for these patients.

2.5 Intended use environment

It can be used both in the treatment rooms of medical practitioners and in the treatment rooms of hospitals. Use of the device in an operating room is not permitted (and not meaningful).

2.6 Medical Device Directive / Medical Product Act

This laser device meets the EU Medical Devices Directive 93/42/EEC.

Device class in accordance with Directive 93/42/EEC: II b

UMDNS: 17-775

GMDN: 47877

3 Safety

3.1 General safety Information

As with any medical device, operation of the PicoStar lead to potential dangers, of which the user should be aware before using the device. These risks include optical, electrical and biological dangers and fire hazards.

Observe the national regulations for using a medical device. Applicable law may require that device operators and users follow a specific series of safety precautions when using class IIb medical devices. Such devices may be operated only in accordance with the generally recognized rules of the technology and occupational safety and accident prevention regulations. Furthermore, you may be required to maintain a medical device book. A corresponding copy is provided.

WARNING

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure. Furthermore, there is also a risk of electrical and biological hazards as well as hazards caused by radio interference.

The instructions specified below must be observed:

- No modification of this equipment is allowed.
- The equipment and procedures for operation, configuration, maintenance, inspection and calibration specified in this user manual must be observed.
- Note that regular (annual) safety checks must be carried out for this device in order to ensure the safety of personnel and patients. Failure to regularly carry out the checks can lead to severe health damage.
The results of the checks are to be documented in the medical device book.
- Note that, in order to ensure the safety of personnel and patients, service and repairs on or changes to this device may only be carried out by Asclepion Laser Technologies GmbH or by those authorized to do so.
Service and repairs on or changes to the device carried out by non-authorized persons can lead to severe health damage.

3.2 Electrical hazards

WARNING

The device is operated using line voltage.

To avoid the risk of electric shock, the following must be observed:

- This equipment must only be connected to a mains supply with protective earth.
- Never open the device if you are not trained or authorized to do so.
- Never set up the device such that disconnection from the power supply may be difficult.
- The device adheres to the specified limit values for leakage currents. Nevertheless, ensure that contact with the patient and device does not occur simultaneously when the device is switched on.
- The device is designed to prevent the ingress of liquids over the course of normal use. If, in exceptional cases, liquid enters into the device housing or hand-piece, you must turn the power off and unplug the main power plug. Please contact our Technical Service.

3.3 Biological hazards

The device is developed and manufactured in a way that a safe use is ensured also with regard to compatibility of materials. The contact with materials of the device or with outgoing media (liquids, gases etc.) in normal use is not dangerous.

CAUTION

The spacer (distance holder) of the handpiece comes into contact with the skin of the patient. There is the risk of transmitting microorganisms.

The handpiece used must be thoroughly cleaned and disinfected immediately after each application; if necessary the spacer can be sterilized.

CAUTION

Laser fumes may contain viable tissue particles if tissue is vaporized. There is the risk of microorganisms being transferred and damage of the respiratory tracts due to the smoke itself. If tissue is vaporized, you are advised to use a smoke evacuation device with suitable filters as well as a suitable respiratory protection for the present people.

3.4 Radio interference

The device fulfills the requirements of the IEC 60601-1-2 standard. The system is not affected by electromagnetic radiation from other devices that conform to the same standard. In addition, the system does not generate any electromagnetic radiation that exceeds the limit values of IEC 60601-1-2. The manufacturer's declaration about electromagnetic compatibility in accordance with IEC 60601-1-2 is included as an appendix to this user manual.

WARNING

Medical electrical equipment needs specific precautions regarding the electromagnetic compatibility (EMC). The instructions specified below must be observed:

- The equipment needs to be installed and put into service according to the EMC information provided in the Appendix.
- Portable and mobile RF communications equipment can affect medical electrical equipment.
The same applies to other equipment that does not meet the standards listed above.

This may result in unintended settings.

Switch off cellular phones and similar devices before using the laser device!

- Use of accessories, handpieces and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper function.

Only use accessories, handpieces and cables that are specified for use with this device by Asclepion Laser Technologies.

- Use of this device adjacent or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment should be observed to verify that they are operating normally.
- The device emits small amounts of electromagnetic radiation that can potentially influence devices that do not fulfill the standards listed above. This may result in unintended settings on those devices.

3.5 Laser

3.5.1 Optical hazards caused by laser emission

The lasers are classified according to their potential for danger. The most dangerous laser class is class 4 (the least dangerous is class 1). This device is a class 4 laser.

WARNING

The use of equipment and procedures for operation, maintenance, inspection or calibration other than those specified in this user manual can lead to hazardous irradiation by laser radiation.

Looking into the aiming beam is dangerous and can cause eye damage.

The protective measures specified below must be strictly observed. In particular, everyone in the room where the laser is being used ("laser area") must wear suitable protective eyewear while the laser is active.

The following protective measures are a few essential requirements that must be supplemented further for the specific circumstances with the help of a laser safety responsible.

1. A **LASER SAFETY RESPONSIBLE** with expertise regarding laser applications in medicine must be commissioned in writing and carry out at least the following tasks:
 - Supporting the operator with respect to safe operation and necessary protective measures
 - Monitoring safe operation of the laser
 - Collaborating with specialists for workplace safety for fulfilling their duties, including training for important matters concerning protection from laser beams.

2. LASER AREA

During operation, the area where the maximum permitted radiation level can be exceeded—called the "laser area"—must be cordoned off and marked by a laser warning sign. The laser area is the room in which the laser is being used. Use of the laser must be indicated by **warning lights** and the triangular yellow **laser warning sign** at all entrance points of the LASER AREA.

Figure 1: Laser warning sign

3. PERSONAL EYE PROTECTION

Everyone present in the laser room during a treatment session must wear laser safety glasses or shields for patients!

The laser safety glasses and eye shields for patients must comply with the specifications defined in the technical data section!

The use of safety glasses with defects of any kind is not permitted.

Even when using laser safety goggles the direct view into the laser beam is not permitted.

4. ADDITIONAL SAFETY PRECAUTIONS

- Cover up windows and other openings to the treatment room to prevent emission of laser radiation.

- Restrict access to the treatment room so that only people assisting with treatment and those trained to handle the laser device have access (besides the patient to be treated).
- Make sure that trained personnel assisting with treatment know how to switch off the laser in an emergency.
- Remove all metallic objects such as watches, rings, necklaces and similar objects from the working area and do not use any reflective tools or other reflective materials if at all possible.
(Reflective objects could interrupt the laser beam and deflect it towards areas other than the intended treatment area. Many surfaces, even those that appear to be dull, can easily reflect the emission wavelength of the laser.)
- Only point the activated laser at the intended area to be treated.
- Never look directly into the outlet opening of the handpiece, even if you are wearing laser safety glasses.
- Switch the laser into the STANDBY state when not in use (the laser cannot be activated accidentally while in the STANDBY state).
- Always remove the key from the key switch when the device is switched off and store the key in a secure place.

5. RESPONSIBILITY

The user who triggers the laser emission is responsible for the safety of the laser. This user must ensure that all safety measures are being followed **before** activating the laser. Specifically, everyone present in the room—including the patient—must wear proper protective eyewear or eye protection for patient.

3.5.2 Laser-induced fire hazard

WARNING

A risk of fire and/or explosion exists when the laser output is used in the presence of flammable materials, solutions or gases, or in an oxygen enriched environment

The device must, therefore, never be used in potentially explosive areas (in accordance with the AP and APG classification defined in IEC 60 601-1).

The following protective measures must be followed for this reason.

- Some materials, for example cotton wool when saturated with oxygen may be ignited by the high temperatures produced in normal use of the laser equipment.
Keep clothing as far away from the treatment area as possible.
- The solvents of adhesives and flammable solutions used for cleaning and disinfecting should be allowed to evaporate before the laser equipment is used.
- Attention should also be drawn to the danger of ignition of endogenous gases.
- Do not use any inflammable substances to cool the skin (such as inflammable cooling spray, but e.g. CoolPacks) directly before or during treatment.
- If possible, do not use any inflammable substances for anesthesia if at all necessary.
- Avoid the use of oxidizing gases such as nitrogen oxide (N₂O) and oxygen. Be especially careful when using oxygen. Oxygen increases both the intensity and extent of fire.
- Keep combustible materials in the treatment room to a minimum. If treatment requires the use of combustible materials, moisten the materials.
- Only expose the surfaces of the body required for treatment to the laser emission.
- Always have a small amount of water READY in the treatment room and have a fire extinguisher for electric devices nearby.

4 Device Description

4.1 Device design

The PicoStar laser is a pulsed Nd:Yag solid state laser. The emission wavelength is 1064 nm. The laser radiation is transmitted through an articulated arm with handpiece to the tissue. By a flip mirror it is possible to redirect the laser beam through KTP-crystal to generate the 532 nm wavelength.


Figure 2: Front of the device


Figure 3: Display of the device


Figure 4: Rear of the device

4.2 Safety Equipment

The device is equipped with different safety units that have been provided to prevent incorrect operation and unintentional activation of the system. All persons who operate the laser or assist during the treatment should be familiar with these units.

POWER SWITCH	The power switch separates the device from all poles of the supply voltage. When the device is not being used, the switch should be in the OFF (O) position.
KEY SWITCH	The system is activated with a key switch. Only authorized persons who have access to the key can start the system. The key switch can only switch on the system if the laser emergency STOP switch is not activated. Always take out the key after you have switched off the device and only allow authorized personnel to store it.
LASER EMERGENCY STOP	The purpose of the red laser emergency stop switch on the display is to immediately switch off the laser in an emergency. Do not use the laser emergency stop switch to switch the device on and off in its normal state.
FOOT SWITCH	The foot switch is an electrical button, which activates the laser emission, if the device is in the READY operating state. You can connect the foot switch to its designated jack at the back of the device. It is equipped internally with two redundant switch elements to ensure safety. Always set up the foot switch near the treatment area.

REMOTE INTERLOCK CONTACT	The system is equipped with a remote interlock contact that can be connected at the entryways to the laser room. If the remote interlock contact is open (e.g. if the door is open), laser emission is automatically interrupted and the device goes into STANDBY.
ADDITIONAL POTENTIAL EQUALIZATION	The additional potential equalization forms a basic protective measure in accordance with IEC 60601-1 in the areas of group 1 used for medical purposes (e.g. doctor's offices) and 2 (e.g. operating rooms; additional connections for the additional potential equalization are attached in the rooms for group 2). This fulfills the purpose of equalizing potential differences between different devices in the patient's environment, to prevent the dropout of the ground wire of devices of protection class I (single-fault protection) and to ensure compliance of the required contact voltage with the housing of protection class I devices. The potential equalization line can be connected to the potential equalization connection at the back of the device, to the corresponding connector of other devices and to the fixed installed system (of the room). The potential equalization does not replace the ground wire of the device.
LED RING	The device is equipped with a LED ring on the display that visually indicates laser emission by a pulse-dependent orange flashing. In STANDBY mode the LED ring lights up white. After switching to the READY mode the LED ring lights up orange.
ACOUSTIC SIGNAL	Laser emission is indicated by an acoustic signal during emission.
WARNING SIGNS	The device has various warning labels. These labels must be clearly visible at all times and must be immediately replaced if damaged.
OPERATING STATE STANDBY	The STANDBY operating state prevents the laser from being unintentionally or accidentally activated. Laser emission is not possible if the system is in STANDBY. The operator can only activate emission after the READY button has been pressed. The system goes into STANDBY in the following situations: - After starting the system for the first time - If the system has been in the READY state for a long period of time without the laser being initiated - If the operator has pressed the STANDBY button while the system was in the READY mode.
AUTOMATIC SHUTOFF	The device has an automatic shutoff. If a specific problem occurs, the system switches automatically to the following safe state: The electronic laser lock is closed, the laser emission is immediately interrupted and the foot switch is deactivated. An error message appears on the display that identifies the specific error.

The manufacturer is only liable for safety, reliability and power of the device if:

1. The device is used in agreement with all instructions in this handbook (on safety measures and system application).
2. Installation, assembly, expansions, modifications, repair and maintenance work carried out by qualified and skilled persons.
3. The electrical system at the installation location that fulfills the requirements of IEC 60601-1 and of local regulations.

4.3 Labels at the device

WARNING

The laser device emits strong laser radiation. This radiation can lead to severe eye damage if the device is used improperly.

Therefore, there are warning and information labels on the appropriate places on the device. The entire personnel must know these labels and their meaning to ensure proper use of the device.

All labels must remain in their respective positions and in good condition. Replace damaged labels immediately to ensure proper use of the device.


Figure 5: Labels on the display

Label 1: Laser Emergency STOP

This sign marks the laser emergency STOP button for switching off the laser in an emergency.


Figure 6: Labels on the handpiece

	Label 2: Beam exit aperture This label indicates the exit point of the laser beam out of the transmission system.
	Label 3: Type B Applied Part Degree of protection B of the Application part.
REF 8392300000 SN 8392300100 Asclepion Laser Technologies GmbH	Label 4: Reference No. and serial No. of the handpiece


Figure 7: Label on top of the device

	Label 5: Laser Warning Sign
--	---




Figure 8: Labels on the rear side of the device

	Label 6: Pushing prohibited The device must not be pushed; it must be pulled by the handle provided for this purpose.
	Label 7: Caution: Observe the accompanying documents!
	Label 8: Equipotentiality This contact is used to connect the device additionally to an equipotential system.
	Label 9: Connection for remote interlock contact
	Label 10: Connection for foot switch


Label 11: Laser warning and explanatory label IEC 60 825-1

This label informs about the dangers, the maximum values of laser radiation and the classification of the laser source.

Label 12: Water in for coolant

Label 13: Type plate

Label 14: Water out for coolant

5 Installation

WARNING

The device is operated using line voltage and emits laser radiation. There is a risk of electric shock and injury to the eyes!

The installation and the mandatory functional test may not be performed by anyone other than our service technicians or other technicians trained and authorized by us.

5.1 Scope of delivery

Important information

Check together with the carrier that the transport packaging is intact (no visual damage, no shock or tilt sensor were triggered) and whether all components that are listed below and required for the safe operation of the device are present and not damaged. Incompleteness and/or damages in transit are to be claimed immediately at the carrier and inform the supplier.

Asclepion Laser Technologies GmbH does not accept any liability unless another agreement has been made expressly.

Please, keep the packing for the future safe consignment of the device.

CAUTION

Refrain from start-up action of any kind if the device or its accessories are found to have suffered mechanical damage. There is risk of physical injury if this advice is disregarded. Contact the local representatives of the Technical Service of Asclepion Laser Technologies GmbH in such cases.

The scope of delivery can be different. It depends on the agreement with the dealer.

The basic equipment comprises the following items:

- Base unit
- Power cable
- Foot switch
- Handpiece
- Door interlock connector
- 2 Keys (for key switch)
- Laser safety glasses, eye protection for patient
- User Manual
- Laser warning sign
- SmartCard

5.2 General information about installation

Important information

After installation and any relocation of the device from a cold to a warm environment with a temperature difference of more than 5 °C, allow the device to adjust to the room temperature in an unpacked condition before use. Allow it to adjust for at least:

2 hours at a temperature difference of up to 10 °C

4 hours at a temperature difference of up to 15 °C

8 hours at a temperature difference of more than 20 °C

Afterward, switch on the device without connecting a handpiece and allow the device to warm up for at least 30 minutes.

Make sure that the environment in which the device is to be installed and operated complies with the permitted ambient conditions stated in chapter *Technical data*.

Failure to observe this information can lead to destruction of the device.

CAUTION

Fix the wheels of the device to ensure that it cannot be moved unintentionally. Non-observance of this instruction may cause risks of injuries.

WARNING

The laser device emits high levels of radiation in the infrared range that may cause serious damage to the eye if the device is used improperly.

Every entrance to the laser area (usually the laser treatment room) must be equipped with warning lamps that are lighting (or flashing) as long as the laser emits radiation to warn persons against entering the room.

AIR QUALITY

- The device should be operated in a non-corrosive atmosphere.
- Corrosive materials such as acids can damage electrical wiring and electronic components.
- Dust is destructive to electrical equipment and will reduce the cooling capacity.

OPERATION CONDITIONS

Refer to chapter *Technical data*.

WARNING

Risk of explosion. The device is not protected against explosions. Do not operate the device in explosive atmospheres.

5.3 Door interlock connection

At the rear panel of the device, a connector is provided for the connection of a door interlock contact to the device.

As port a round plug has to be used which can be screwed, type (x)V 40 according to IEC 60130-9, e.g. Lumberg WSV 40 (as delivered). (x): dependent on type, e.g. WS for an angled plug as delivered.


Figure 9: Connections for remote door interlock contact and laser warning lamp

By connecting the two connections remote door interlock contact and external door contact you can prevent laser radiation from being emitted when the door of the laser room is opened.

This safety system serves to shut off the laser beam automatically when the door is being opened thus precluding any risk to the person entering.

The contact connector delivered is factory-fitted with a jumper. To connect a door interlock contact, remove this jumper and connect your door interlock contact to the same terminals. The external contact must be potential-free and designed for at least 24 V DC and 1 A. If a door-interlock contact is not used, make sure to keep the plug with the factory-fitted jumper connected.

Important information

The contact plug may only be connected with the device being turned off and disconnected from the power supply. Non-observance of this instruction can cause damages of the device.

The contacts for the remote interlock and the foot switch must never be connected with line power because otherwise the device will be seriously damaged. Connect these two connectors only in the way described in this section.

Tip

If you have connected the door interlock contact with the entrance door, please check every day before starting a treatment that the laser cannot be switched to the READY mode with this door being open. In case of malfunction please inform the Technical Service.

An external voltage source for a low-voltage warning lamp can be connected using two additional contacts of this plug connector. The switching relay inside the device is designed for the maximum voltage or maximum current as shown in the figure above. The safety requirements below must be followed when connecting an external warning lamp:

WARNING

To avoid the risk of an electric shock with consequential serious health problems as well as fire hazard of the transformer, if you use an external transformer to power the lamp, please note:

A safety transformer complying with IEC 61558-2-6 must be selected.

The internal relay of the device is closed, if the device is in the READY mode.

The prescribed laser warning lamp may, of course, also be turned on and off via a normal light switch if the laser is operated.

5.4 Power grid connection

Important information

Observe all of the technical data for the local power supply and the device when selecting a power socket. First check that your power supply satisfies the previously mentioned connection requirements (supply voltage corresponding to type plate). Connecting the device to an unsuitable power supply can lead to damage to the device.

Please observe the requirements from regulations/standards in your country in their respective applicable versions.

Before connecting the device to the power grid, the shift paddle of the main switch must be on 0 (device de-energized). First plug the power cable into the power input on the rear side of the device before you connect the other side to the power socket. Operation of the laser requires that the device be connected to a single-phase power outlet protected by a **separate** fuse (delay-action fuse preferable).

5.5 Potential equalization line

The potential equalization line can be connected to the potential equalization connection at the back of the device, to the corresponding connector of other devices and to the fixed installed system. The potential equalization does not replace the ground wire of the device.

Important information

Please observe the requirements of the IEC 60601-1 regarding medical electric systems.

5.6 Foot switch connection

You can now connect the foot-operated switch to its designated jack at the back of the device. Insert the connector until mechanical stop position, then turn the union nut clockwise for locking. The footswitch should be placed directly beside the device.


Figure 10: Foot switch

5.7 SmartCard

The SmartCard is an unlimited time card that has to be inserted into the device before the treatment begins. The PicoStar SmartCard can be used for each handpiece.

Figure 11: Front and back of the SMARTCARD

Using the SMARTCARD

1. Insert the SMARTCARD into the card slot on the side of the display.

Figure 12: SmartCard slot

2. Switch on the device.
3. Wait until the card has been read completely (this takes approx. 3 sec).
4. The device is now ready for operation.

Important information

- The SMARTCARD only works on the respective PicoStar and cannot be transferred to another device.
- Simple exchange of the handpiece on the device is possible with the same card.
- In the PicoStar setup menu, you can find all important information about the SmartCard.

5.8 Handpiece connection

The handpiece is screwed to the distal end of the articulated arm. The device is provided with an automatic detection system for the handpiece and thus for the spot size of the treatment beam. This information is sent through an electrical cable to the control module of the device. The plug of the spot detection cable must be put into the corresponding socket, which is fixed at the articulated arm.

Figure 13: Handpiece connection (Example Zoom handpiece ST)

In order to change the handpiece, the plug of the spot detection cable has to be removed from the articulated arm. Then the handpiece can be screwed off.

The spot size of the Zoom handpiece can be changed manually by turning the setting wheel of the handpiece to the desired marking. The device software automatically detects and displays the spot size.

Zoom handpiece ST

Figure 14: Handpiece

Important information

Make always sure that all connections are properly engaged and the handpiece is screwed to the articulated mirror arm.

Make sure that the spot size of the handpiece is identical to the set spot size on the display of the device.

Always bring the handpiece in a safe position, as long as it is not used.

The selected handpiece must be disinfected before each application and cleaned and disinfected after each application as described in chapter *Cleaning, disinfection and sterilization*.

6 Technical data

6.1 Device models and applied parts

Device	Description
PicoStar Item No. M16-9400-000	Device with Nd:YAG laser (including frequency doubling) up to 800mJ 1064nm; mains supply with 200/230 VAC
PicoStar Item No. M16-9500-000	Device with Nd:YAG laser (including frequency doubling) up to 800mJ 1064nm; mains supply with 100/110/120 VAC

Applied parts	Description
Zoom handpiece ST Item No. 839-2400-000	Adjustable spot sizes: 2, 3, 4, 5 mm Laser beam is <u>not</u> collimated It must be worked with a distance holder (spacer)

WARNING

Accessories and disposables other than those specified in this user manual may not be used with this device.

All removable accessories and disposables are approved for this device only and may be used only together with this device.

Only use accessories and disposables provided by Asclepion Laser Technologies GmbH for this device. The use of accessories not tested and released by Asclepion Laser Technologies GmbH for the use with this laser device can cause severe damage to the eye or electrical shocks.

Important information

We strongly advise against the use of accessories made by other manufacturers. Even if an official testing authority should have certified that a specific accessory unit could be used safely, Asclepion Laser Technologies GmbH cannot assume any liability for these products and the damages eventually resulting in their usage with this device.

6.2 Specifications

Basic specification	
Device model	Moveable unit
Display	LCD display
Operator guidance	Touchscreen
Cooling	Internal cooling system (closed internal cooling circuit filled by the manufacturer)
Door interlock connector	5 V / 10 mA

Basic specification	
Laser warning lamp	Potential free relay contact max. 24 V / 1 A (make contact)
Permissible ambient conditions	Temperature: 15°C to 30°C rel. air humidity: max. 85 % (not condensing) Height: max. 2000 m above sea level
Permissible transport and storage conditions	Upright position (not to be laid down at any time!) Temperature: 5°C to +50°C rel. air humidity: 10 % to 90 % (no condensation) Air pressure: 700 – 1070 hPa
Dimensions device with handle	1170 x 450 x 1060 (D x W x H)
Weight	around 160 kg
Power requirements	230 V AC, 50/60 Hz, max. 2500 VA 200 V AC, 50/60 Hz, max. 2500 VA 120 V AC, 50/60 Hz, max. 2500 VA 110 V AC, 50/60 Hz, max. 2500 VA 100 V AC, 50/60 Hz, max. 2500 VA
Main fuse	Over current release 16 A, slow-blow for ≥ 200 V AC Over current release 25 A, slow-blow for ≤ 200 V AC
Overvoltage category	II (IEC 60664-1)
Max. peak of mains voltage	V MT peak = 2500 V peak
Max. nominal mains voltage	V MN r.m.s. = 300 V r.m.s.
Classification according to directive 93/42/EEC	II b
Operating mode	The device has been designed for continuous operation.
Protection class of the application part	B
IP protection class	IP 20

Laser specifications	Nd:YAG laser
Laser model	pulsed Nd:Yag solid state laser
Wavelength	1064 nm 532 nm
Laser class	Class 4

Laser specifications	Nd:YAG laser
Pulse Energy	1064 nm → 20 mJ – 800 mJ 532 nm → 20 mJ – 300 mJ
Pulse duration	250 ± 150 ps
Pulse frequency	0.5, 1, 2, 5, 10 Hz @ 532 & 1064 nm
Spot size	ST Handpiece: ø 2, 3, 4, 5 mm
Tolerance of output power	± 20 %
Laser beam mode	Gaussian mode
Beam divergence at the handpiece output (round angle, 1/e ²)	Zoom Handpiece ST (532 nm/1064 nm): • Typ. 20 mrad @ 2 mm Spot • Typ. 15 mrad @ 3 mm Spot • Typ. 11 mrad @ 4 mm Spot • Typ. 14 mrad @ 5 mm Spot
NOHD - Nominal ocular hazard distance (with handpiece)	Zoom Handpiece ST: • 963 m @ 532 nm • 446 m @ 1064 nm
MPE Limit Aperture	7 mm

Laser specifications	Pilot (aiming) beam
Laser model	Laser Diode
Wavelength	635 nm
Laser class	class 2
Operation mode	Pulsed or CW selectable
Pulse duration	≥ 25 ns

Beam divergence at the handpiece output (round angle, 1/e ²)	Zoom Handpiece ST (532 nm/1064 nm): • Typ. 20 mrad @ 2 mm Spot • Typ. 15 mrad @ 3 mm Spot • Typ. 11 mrad @ 4 mm Spot • Typ. 14 mrad @ 5 mm Spot
Power	≤ 0,9 mW, selectable brightness
MPE Limit Aperture	7 mm

Required laser safety goggles acc. to DIN EN 207:2009	Safety glasses for doctor
	1. Protective properties: Light beams over 480nm are completely blocked 2. Spectral transmittance: at lengths 380-480 nm - 0.45-0.75 3. Refractive index: 1.62 4. Weight (0.07 ± 0.01) kg
	Patient goggles 1. Protective properties: Light beams are completely blocked 2. Spectral transmittance: Light beams are completely blocked 3. Refractive index: Light rays are completely delayed 4. Weight - (0.05 ± 0.01) kg
General requirements	
5. Requirements for the quality of materials and cleanliness of the surface of the light filters: the surfaces of the glasses must be smooth, without scratches, burrs and other damages.	
6. Requirements for the resistance of goggles to ultraviolet radiation and elevated temperatures: - When exposed to ultraviolet radiation, safety glasses must meet the requirements of paragraphs. 1-4. - After testing goggles for 5 hours in a climatic chamber at a temperature of (55 ± 2) ° C and a relative humidity of at least 95%, as well as after holding for at least 2 hours at room temperature, goggles must meet the requirements of paragraphs. 1-4. The permissible relative change in the light transmittance is not more than ± 10%.	
7. Requirements for the resistance of protective goggles to ignition: no part burns or continues to smolder after removing the steel rod.	
8. Requirements for the field of view of safety glasses: safety glasses must have a free field of view in the vertical and horizontal directions of at least 40 ° for each eye	
9. Requirements for the design of safety glasses: Parts of safety glasses that come into contact with human skin are made of materials that do not cause skin irritation.	
10. Requirements for goggle frames: The frame is designed to prevent accidental laser radiation from the sides.	
11. Requirements for the mechanical strength of goggles: The goggles withstand the impact of a steel ball with a nominal diameter of 22 mm and a force of (100 ± 2) N.	

7 Operation of the device

WARNING

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

The equipment and procedures for operation, maintenance, inspection and calibration specified in this user manual must be observed.

The device is operated using line voltage. There is risk of electric shock!

The device adheres to the specified limit values for leakage currents. Nevertheless, ensure that contact with the patient and device does not occur simultaneously when the device is switched on.

7.1 Switching on

Before switching on the device, make sure that the following requirements are fulfilled:

- The device is in proper condition (no damage).
- The voltage specified on the type plate is identical to the local supply voltage.
- The power plug of the device is connected to the correct power supply.
- The footswitch plug is properly connected and screwed to the footswitch socket on the rear side of the device.
- The plug for remote interlock contact is properly connected and screwed to the rear side of the device.
- The handpiece is connected (cleaned and disinfected).
- Warning lamps are switched on at all entrance points.
- Secure the device against unintended moving by fixing the wheels. Otherwise there is the danger of injuries.
- Damaged parts must be exchanged immediately to avoid injuries.
- The patient should be positioned in a way that the areas to be treated are easily accessible and enable patient comfort and safety.

Ensure that all safety measures have been met. Then proceed as follows:

1. Switch the power switch on the rear side of the device to I (On).
2. Then insert the key into the key switch and turn it clockwise as far as it will go.

The starting screen appears on the monitor as the system runs an automatic test routine of assemblies relevant to safety. Afterwards the main menu will be displayed.


Figure 15: Start screen and main menu screen

7.2 General information on touch screen operation

The treatment parameters are to be set via the touch screen on the display at the front. Configuration is very easy and self-explanatory.

- To select a function press the corresponding icon in any screen. A signal is audible.
- Parameters can be selected by pressing the corresponding parameter field. The selected parameter field will be highlighted.
- The laser can be activated by pressing the READY button 
- To switch to the STANDBY mode please press the STANDBY button 

Overview of important functions of the device (icons)

Function/ Parameter	Display	Remarks
Setting keys		Scrolling through lists
Home		Return to main menu
Back		Return to the last screen
Up/Down button		Select values of different parameters
READY mode (start laser emission)		• Activate laser emission (READY mode) • Pressing the footswitch releases laser radiation from the handpiece.
STANDBY mode (stop laser emission)		• Stop laser emission (STANDBY mode) • Unintentional activation is not possible
Cancel		Exit function
Error Button OK		Confirm error message
Delete		Delete treatment setting
Save		Save settings for your treatment

7.3 Main menu

In the main menu you have the option to choose between three different settings: Laser, Indication and Setup.


Figure 16: Main menu screen

Function/ Parameter	Display	Remarks
Laser		Laser parameters can be retrieved and stored
Indication		Shows you a list of available indications
Setup		Choose general display features and show software information

7.4 Setup

You have the option to adjust general settings for your display, such as language and volume. Call up the Setup menu by pressing the *Setup* button on the main menu. The screen changes to the Setup screen.


Figure 17: Setup screen

The following functions can be carried out by pressing the corresponding buttons:

Function/ Parameter	Display	Remarks
Info		Displays software version and information about counter and TimeLock
Sound		Volume and sound settings. Change system level and warning level and type

Import/Export		Export/Import data to or from a USB stick
Language	Aあ	Select between different languages

7.4.1 Import/Export data

To import or export data:

- Select the *Setup menu* 
- Press the *Import/Export* button 

You can either import/export favorites (e.g. saved parameter settings; sound or language settings) or Log files from the system. Log files will be automatically copied during the operation of the device and are used to detect technical errors. Please export these Log files if our technical service asks for it.

Function/ Parameter	Display	Remarks
Export Log files		Export technical data that are relevant for our technical service
Export		Export favorites settings to a USB stick
Import		Import favorites settings from a USB stick

Tip

The use of the USB interface is only allowed in STANDBY mode and only for ex- and import of data which are stored at a suitable USB stick.

Data of parameter settings should be saved routinely at a suitable USB stick, in order to avoid loss of data due to a defective device. The data is only legible on the display of the PicoStar.

7.4.2 Software/ System information

You will find general information about the device and software in the *Setup menu*:

- Select the *Setup menu*
- Press the *Information* button

You can display information about Software, Counter and TimeLock:

Function/ Parameter	Display	Remarks
Counter		Find information about applied laser impulses and foot switch stop stopwatch
TimeLock		TimeLock manager
SmartCard		Find information about the SmartCard type
Software		Find information about the software

7.5 Starting the treatment

The device is suitable for a variety of applications. In order to begin a therapeutic procedure, you have to switch to the treatment screen.

Figure 18: Treatment screen

There are two options to access the treatment screen from the main menu:

- Starting via *Indication*
- Starting via *Laser*

Access to the treatment screen will only be allowed if the correct handpiece is connected. If no handpiece or the wrong handpiece is connected, an information will be displayed on the screen.

The treatment screen allows you to select between spot size (*Spot*), energy density (*Fluence*), treatment speed (*Speed*), wavelength (*Color*), operating mode (*Mode*) and pilot beam (*Pilot*).

The selected parameter will be highlighted. With the up/down  buttons you can manually choose between different value settings. In case of a red highlighted value, the value is not in our recommended applicative range.

7.6 Starting via Laser

By pressing the *Laser* button  in the main menu you will reach to the treatment screen, where you can select and set different parameters.


Figure 19: Treatment screen

Function/ Parameter	Display	Remarks
Spot size	Spot	Spot size currently set on the handpiece
Energy density	Fluence	Fluence settings
Treatment speed	Speed	Speed rate of therapy procedure
Wavelength	Color	Select required wavelength
Operating mode	Mode	PICO pulse mode
Pilot (aiming) beam	Pilot	You can select between pilot intensity (level) and pilot mode (continuous wave (CW) or blink)

Tip

Only the modes feasible for a specific handpiece can be activated. This means that not every handpiece is configured for every mode.

7.6.1 Saving of parameters

- You can save modified settings by pressing the *Save* button 

After pressing this button, a keyboard appears on the display where you can type in the name of your settings.

Figure 20: Keyboard to save parameter settings

To write initial letters in uppercase, press the *Shift* button once; do not hold the button down. After typing a letter, the keyboard automatically shifts back to lowercase letters. You can toggle the keyboard to writing uppercase letters by pressing the *Caps* button.

To call up parameter savings:

- Select the *Main menu*
- Press the *Indication* button
- Select *Favorites* to open stored parameters

7.7 Starting via Indication

By pressing the *Indication* button in the main menu, a list of available indications will appear. Select your indication by pressing the corresponding field. The treatment screen of that indication will appear.



Figure 21: Indication list screen

Function/ Parameter	Display	Remarks
Tattoo		Removal of tattoos
Make-up		Removal of permanent Make-up
Pigment		Pigment removal

7.7.1 Indication list: Tattoo

You have the option to choose between different colors settings.


Figure 22: Indication Tattoo

7.7.2 Indication list: Make-up

In order to remove permanent Make-up you can choose between black and red color.


Figure 23: Indication Make-up

7.7.3 Indication list: Pigment

With the PicoStar you can perform laser treatment of pigmented skin, such as Becker Nevus, Café au lait, Freckles, Hyperpigmented Scars, Lentigo, Nevus Ota and Pigmented Nevus.


Figure 24: Indication Pigment

7.7.4 Emission of laser beam (READY mode)

In order to start laser emission make sure that

- All laser protection measures have been met
- Check the function of the door interlock if you have it installed: Open the door of the treatment room; it must not be possible to switch the laser to the READY mode: (an information is displayed and must be acknowledged)
- The handpiece is cleaned, disinfected and properly connected
- Make sure that the connected handpiece matches the information on your display (e.g. the spot size of the handpiece is identical to the set spot size on the display). If not, switch off the device and contact your customer service.

WARNING

Direct and scattered laser radiation can cause severe eye injuries.

Everyone in the room where the laser is being used must wear protective eyewear or patient blackout glasses suitable for this laser while the laser is active. Other forms of protective eyewear may not ensure sufficient protection!

The user must confirm this before pressing the READY button.

In the READY state, you can also activate a laser pulse by unintentionally pressing the foot switch! Therefore, the laser device must be switched to the STANDBY state immediately during breaks in treatment and at the end of treatment.

The laser beam may be triggered only if the handpiece is precisely directed onto the skin area to be treated.

In the case of any dangerous faults, the laser emergency stop button has to be activated immediately.

WARNING

The intensive laser beam of this device can cause severe damage of the skin.

It is the user's responsibility to correctly determine the zone of treatment and related treatment parameters. He/she is also responsible for correct positioning of the handpiece.

The following functions can be carried out by pressing the corresponding buttons:

Function/ Parameter	Display	Remarks
READY mode (start laser emission)		• Activate laser emission (READY mode) • After starting the READY mode, the yellow laser warning sign appears on the display. • Before activating the treatment beam, the beam delivery system has to be checked • Pressing the foot switch releases the laser radiation from the laser tip of the handpiece with the set parameters. Releasing the foot switch immediately interrupts the laser emission. • The emission of the laser radiation is indicated audibly by a warning tone and visually by flashing of the LED ring on the display in orange.

Function/ Parameter	Display	Remarks
STANDBY mode (stop laser emission)		- Stop laser emission (STANDBY mode) - Unintentional activation is not possible


Figure 25: READY mode (with laser warning sign)

7.7.5 Check of the beam delivery system

The beam delivery system consist of the articulated mirror arm and the handpiece.

Important information

As the pilot (aiming) beam takes the same path through the laser beam delivery system as the working beam, it offers a good method for checking the integrity of the beam delivery system. If the pilot (aiming) beam does not appear at the distal end of the beam delivery system, its intensity is low or if it looks diffuse this may indicate that the beam delivery system is either damaged or misaligned.

If the spot of the aiming beam shows any imperfections at the distal end of the beam delivery system, you must not activate the treatment beam to avoid damage to the device.

Assess the spot of the aiming beam only by directing the handpiece to a white sheet of paper (spacer in contact). Never look directly into the beam (into the handpiece).

The therapy beam may only be activated after the successful check of the beam delivery system. Make sure that you do not press the footswitch unintentionally.

To change the settings of the pilot (aiming) beam, select the parameter *Pilot* in your treatment screen. Choose the Pilot level and the mode (continuous wave (CW) or blink).

Figure 26: Pilot mode and level (aiming beam)

If the aiming (pilot) beam is not seen perfectly, please proceed as follows:

- If the intensity is low, please check the settings in the treatment screen (Pilot mode)
- Switch off the device
- Check once more, that the handpiece is correctly screwed to the articulated mirror arm. Please screw on again and make sure, that the handpiece is straight and not tilted.
- Check once more, that the optics of the handpiece is clean.
- If all is in best order and the pilot beam is not or hardly visible, the device must not be used any longer. In these cases, you should stop using the device and instantly contact the Technical Service.

7.8 Switching off the device

STOP LASER EMISSION

- Exit the *READY* state and switch to the *STANDBY* state by pressing the *STANDBY* button
- Switch off the device by turning the key switch counterclockwise as far as it will go
- Only then switch off the power switch at the rear panel of the device. Switching the device power switch to "0" (OFF) completely disconnects the device from the line voltage

STOP LASER EMISSION IN AN EMERGENCY

- To immediately switch off the laser emission system in an emergency, press the red laser emergency STOP switch on front of the device

WARNING

Direct and scattered laser radiation can cause severe eye injuries.

Never leave the laser device unattended while it is switched on. Completely disconnect the device from the line voltage with the power switch. After switching off the device, ensure that the key is removed and kept in a secure place to prevent unauthorized use of the device.

8 Cleaning, disinfection and sterilization

8.1 Device

WARNING

The device is operated using line voltage. There is risk of electric shock!

Before cleaning and disinfecting, disconnect the device from the line voltage by unplugging the power cable from the power socket.

Cleaning or disinfecting using wet clothes or spraying products such as disinfectants is prohibited.

The device must be cleaned and disinfected at least daily after the last use of the device. The device should be covered with a clean cloth when it is not being used every day.

PROCEDURE:

- Switch off the device by the key switch and power switch. Disconnect the device from the line voltage by unplugging the power cable from the power socket.
- Put on single use gloves.
- All accessible outside surfaces of the device are cleaned using a slightly moistened, soft wipe. Doing this, the connections at the rear side of the device are left out.
- Use Kodan® or Bacillol® wipes to remove stubborn dirt. Abrasive agents must not be used!
- All accessible outside surfaces are carefully disinfected using Kodan® or Bacillol® wipes. Observe the instructions for disinfection given by the manufacturers.
- The connections at the rear side of the device are left out.
- Aggressive disinfectants must not be used!

8.2 Handpiece

CAUTION

The contact surface of the handpiece can come into contact with the patient's skin during treatment. There is a risk of transfer of microorganisms. The handpiece used should be thoroughly cleaned and disinfected according to the reprocessing procedure of the PicoStar handpieces.

If it is necessary to clean the handpiece cover during treatment, the device should be set to STANDBY modus, before you begin cleaning! In this case do not use flammable substances such as e.g. alcohol or isopropyl alcohol near the point of treatment, as there is a, however slight, fire hazard, when subsequently lasering!

Dirt deposits on the optical surfaces can burn in due to the laser radiation and cause side effects such as redness, potential formation of blisters or scabs, and scars in rare cases!

Make sure that accessible optical surfaces are not damaged and remain clean after cleaning and during treatment.

A procedure for cleaning, disinfection of the handpiece and cleaning, disinfection and sterilization of the spacer is described in a separate reprocessing manual *Reprocessing of PicoStar handpieces*.

9 Fault messages and troubleshooting

WARNING

The laser is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

You are strongly advised against trying to do repair works at the device yourself that go beyond the instructions given on the display!

Important information

Improper use or maintenance of the PicoStar device may void the warranty granted by Asclepion Laser Technologies GmbH. Before starting any repair or if you have questions that have not been answered by this manual, contact your service partner.

9.1 General information


If the device does not respond after turning the key switch, check that:

1. The device is connected to the power supply
2. The main power switch on the rear side of the device is switched on.

There are three types of fault messages. General errors (green and yellow) can be remedied by the user. To do so, follow the instructions on the display.


Message that has to be confirmed.

Message where you have to make a selection.

For safety-relevant errors (red), the device shuts off within milliseconds. If this is the case or an error occurs that is not described in this user manual, disconnect the device by unplugging the power plug from the power supply at the power outlet.

Error message without any possibility of action. The device turns off.

In case of any fault the device will switch to the safe mode and stop emission of laser radiation. Please try to remove the fault yourself by following the instructions given at the display. If the error persists, contact Asclepion Laser Technologies GmbH or a service partner expressly authorized by Asclepion Laser Technologies GmbH.

Before you contact the Technical Service of Asclepion Laser Technologies GmbH, **please make a note of the fault message and the error code displayed.**

9.2 Error Messages and Trouble Shooting

You are requested to consult the table below for an initial attempt to remove a fault situation by yourself. If the error persists, you should contact Asclepion Laser Technologies GmbH or a person expressly authorized to handle fault removal responsibilities on behalf of Asclepion Laser Technologies GmbH.

Please take down the name of a given error message and the related error code before you contact the Service Department of Asclepion Laser Technologies GmbH.

Error message	Explanation / Removal
Temperature error (every 15X error)	Operating temperature too low or too high. Check for adequate air supply to casing. Perform system restart. You may resume operation. Contact Service if error occurs again.
Foot switch (40X error)	Unlock the footswitch.
Door Interlock (60X error)	Door contacts are broken or door contact connector is unplugged. Door is open if door was connected to door contact. Connect door contacts conformingly or activate door contact switch. Close door if door was connected to door contact.
Laser error (94X error)	Restart the device. Laser source error.
Smart Card (100X error)	Check that the correct SmartCard is insert correctly into the SmartCard reader on the front side of the device.
Ready mode error (205X error)	Error in switching ready.

10 Customer service

WARNING

The device is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

Never open the device and attempt to repair the device yourself if the laser does not work properly.

Only specialists appropriately qualified and authorized are allowed to open the device.

The device includes no parts of a kind that is intended for maintenance by personnel of the owner/operator. Only service engineers employed with or trained by Asclepion Laser Technologies GmbH are authorized to service the device. This also applies to the prescribed annual safety tests.

If you need assistance by the Technical Service, please contact the local partner of Asclepion Laser Technologies GmbH in your country. We can help you to find your local partner via the address or phone number given below. When contacting your service partner of Asclepion Laser Technologies GmbH, always keep the model and the serial number of the device at hand. Model and serial number are specified on the type plate located at the rear panel of the laser device.

Asclepion Laser Technologies GmbH

Bruesseler Str. 10
D - 07747 JENA
Germany

Phone: +49 (0) 3641 / 7700 - 401
Fax: +49 (0) 3641 / 7700 - 402
E-mail: service@asclepion.com

Please note that the manufacturer, installer and importer will not consider themselves responsible for the effects of servicing on safety, reliability and performance of the device unless the following requirements are met:

- installation, extensions, readjustments, modifications or repairs have been performed by authorized persons,
- the electric installation conditions of the laser room meet the requirements of IEC 64, as amended and
- the device has been used in accordance with the instructions given in this manual.

11 Regular maintenance, safety checks and calibration

The PicoStar does not contain any user-serviceable components. After the regular maintenance, safety checks, calibration and use of the device/handpieces under normal conditions, the PicoStar has a lifetime of 10 years and Zoom handpieces of 12 months.

WARNING

The device is operated with line voltage. There is the risk of an electrical shock and damage to the eye!

The operation, maintenance, test and calibration instructions given in this manually are to be followed.

Do not perform any maintenance work while the device is in operation.

Only service engineers employed with or trained by Asclepion Laser Technologies GmbH are authorized to service the device. This also applies to the prescribed annual safety tests, calibration included.

CAUTION

The use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure.

CAUTION

The PicoStar must be checked by a service technician of Asclepion Laser Technologies GmbH or a qualified service technician of the distributor within the technical safety check at least **every 12 months**, the device is to be calibrated, too if required.

Important information

Asclepion Laser Technologies GmbH granted for the PicoStar a warranty for one year.

Improper procedures for operation, maintenance, safety check and calibration can annul the warranty. If you have any questions that were not included in this manual, please contact your service representative.

11.1 Routine maintenance

Inspection of the outer surfaces of the device

Check if all cables are intact and connected firmly.

Frequency: weekly

To be checked by: personnel of clinic or hospital

Check of the cooling system

Check if all ventilation grids of the cooling system are free.

Frequency: weekly

To be checked by: personnel of clinic or hospital

Check of the accessories

Check if the accessories are in perfect state.

Frequency: before each use

To be checked by: operator of the device

11.2 Regular safety check

SCOPE OF THE ANNUAL SAFETY CHECKS

At least the checks listed below have to be performed according to the legal requirements in their currently valid version as well as to the international standard for Medical Electrical Equipment – re-current test and test after repair of medical electrical equipment (IEC 62353).

The tests should be run in the specified order:

- Visual inspection of laser device and accessories
- Measurement of protective earth resistance
- Measurement of leakage currents
- Measurement of insulation resistance
- Functional test

Additionally to be performed:

- Measurement of the actual output power of the laser modules at the handpiece; adjustment of laser power, if necessary

Tip

Contact your dealer or Asclepion Laser Technologies GmbH directly for detailed specifications on the scope and implementation of the visual inspection.

12 Disposal

The device is to be disposed of according to regulation 2012/19/EU of the European Council for Waste Electric and Electronic Equipment [WEEE] or the specific national regulations on waste electric and electronic devices, as amended.

Please consult our Technical Service if you have any questions regarding disposal.

13 EC declaration of conformity

Appendix A: EMC manufacturer's declaration

Device	Model number
PicoStar Item No. M16-9400-000	Device with Nd:YAG laser (including frequency doubling) up to 800mJ 1064nm; mains supply with 200/230 VAC
PicoStar Item No. M16-9500-000	Device with Nd:YAG laser (including frequency doubling) up to 800mJ 1064nm; mains supply with 100/110/120 VAC

Accessory / designation	Model number	Length / dimensions
Foot switch with cable	866-0806-008	< 2.5 m
Power-supply cord	855-0704-013	< 3.0 m
Zoom handpiece ST	839-2400-000	

Important information

The device is subject to special terms of network connections. The maximally allowed impedance of the power grid is: $I_{Z_{maxl}} = 0.0634 \text{ Ohm}$.

If necessary, the local energy supplier is to be contacted to ensure that the device will be plugged in a supply which has an impedance lower or equal $I_{Z_{maxl}}$.

Important information

For the device, the following essential performances are established:

- No independent switch of STANDBY to READY and conversely
- No independent (unintentional) release of output energy
- No independent changing of the preselected output energy
- No output of pulses in other pulse length than preselected

A.1 Electromagnetic emissions

The PicoStar is suitable for use in the specified electromagnetic environment. The customer and/or the user of the PicoStar should ensure that it is used in an electromagnetic environment as described below.

Emissions measurement	Compliance	Electromagnetic environment guidelines
High-frequent interference emission to CISPR 11	Group 1	The PicoStar uses HF energy only for its internal function. Therefore, its emissions of HF interferences are very low and are not likely to cause any interference in nearby electronic equipment.
High-frequent interference emission according to CISPR 11	Class B	The PicoStar is suitable for use in all establishments, including domestic establishments and those directly connected to the public power supply network that supplies buildings used for domestic purposes. For harmonics :
Emitted interference of harmonics according to IEC 61000-3-2	NA 1)	The PicoStar is intended for use in all establishments, including domestic establishments and those directly connected to the public power supply network that also supplies buildings used for domestic purposes For flickers:
Emitted interference of voltage fluctuations / flickers according to IEC 61000-3-3	Fulfilled	The PicoStar is intended for use in all establishments, including domestic establishments and those directly connected to the public power supply network that also supplies buildings used for domestic purposes. $ Z_{max} = 0.0634 \Omega$

1) for professional use only

A.2 Electromagnetic immunity

Immunity test	Test level acc. IEC 60601-1-2	Compliance level	Electromagnetic environment guidelines
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	+/- 8 kV contact charge +/- 15 kV air discharge	+/- 8kV +/- 15kV	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Fast transient electrical interferences/bursts acc. to IEC 61000-4-4	+/- 2 kV for power lines +/- 1 kV for input/output lines	+/- 2kV +/- 1kV	Mains power quality should be that of a typical commercial and/or hospital environment.
Surges acc. to IEC 61000-4-5	+/- 1 kV series-mode voltage +/- 2 kV common mode voltage	+/- 1kV +/- 2kV	Mains power quality should be that of a typical commercial and/or hospital environment.
Voltage dips, short interruptions and voltage fluctuations of power supply input lines acc. IEC 61000-4-11	0 % UT for ½ cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°)	0%	Mains power quality should be that of a typical commercial and/or hospital environment. If the user of the device requires continued operation during power mains interruptions, too, it is recommended to power the device by an uninterruptible power supply or a battery.
	0 % UT for 1 cycles	0%	
	70 % UT for 25/30 cycles single phase: at 0°	70%	
Magnetic field for the supply voltages (50/60Hz) acc. to IEC 61000-4-8	30 A/m	30 A/m	The magnetic field strength should be at levels typical for commercial and/or hospital environment.

Immunity test	Test level acc. IEC 60601-1-2	Compliance level	Electromagnetic environment guidelines
Conducted			Portable and mobile radio equipment, including cables, should be used not closer to any part of the device than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{6}{3} \right] \sqrt{P}$
HF interferences acc. to IEC 61000-4-3	3 V _{eff} 150 kHz to 80 MHz	3V	
Radiated			where P is the nominal power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength of fixed radio transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
HF interferences acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

^a Field strengths of fixed transmitters, such as base stations for radio telephones (cellular/cordless) and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance levels given above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range from 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

A.3 Recommended protection distances

Transmitter frequency	150 KHz to 6 GHz
Equation	$d = \left[\frac{6}{3} \right] \sqrt{P}$
Rated output of the transmitter (Watt)	Protection distance (meter)
0.01	0.2
0.1	0.63
1	2
10	6.32
100	20

For transmitters rated at a maximum output power not listed above, the protection distance can be estimated by using the equation in the corresponding column, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

Note:
These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Immunity test	Test level acc. IEC 60601-1-2	Compliance level	Electromagnetic environment guidelines
Conducted HF interferences acc. to IEC 61000-4-3	3 V _{eff} 150 kHz to 80 MHz	3V	Portable and mobile radio equipment, including cables, should be used not closer to any part of the device than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distance: $d = \left[\frac{6}{3} \right] \sqrt{P}$ where P is the nominal power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength of fixed radio transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated HF interferences acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

^a Field strengths of fixed transmitters, such as base stations for radio telephones (cellular/cordless) and land mobile radios, amateur radios, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance levels given above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range from 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

A.3 Recommended protection distances

Transmitter frequency	150 KHz to 6 GHz
Equation	$d = \left[\frac{6}{3} \right] \sqrt{P}$
Rated output of the transmitter (Watt)	Protection distance (meter)
0.01	0.2
0.1	0.63
1	2
10	6.32
100	20

For transmitters rated at a maximum output power not listed above, the protection distance can be estimated by using the equation in the corresponding column, where P is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer.

Note:
These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

4.4 Enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	300 – 390	TETRA 400	Pulse Modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745	704 - 787					
780	704 - 787					
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870	800 - 960					
930	800 - 960					
1720	1700 - 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1845	1700 - 1990					
1970	1700 - 1990					
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5500	5100 - 5800					
5785	5100 - 5800					

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

A.4 Enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	300 – 390	TETRA 400	Pulse Modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745	704 - 787					
780	704 - 787					
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870	800 - 960					
930	800 - 960					
1720	1700 - 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1845	1700 - 1990					
1970	1700 - 1990					
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5500	5100 - 5800					
5785	5100 - 5800					

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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File No. W - 988 / 2020

I hereby certify the signature, having been carried out before me, of


Mrs. Dania Di Pietro Paolo,
born on **21.08.1978,**
residing at **Mechelroda,**
business address: **Brüsseler Str. 10, 07747 Jena,**
who identified herself by means of her identity papers (Italian Identity Card).

Mrs. Dania Di Pietro Paolo acts as authorized representative of the limited company:

Asclepion Laser Technologies GmbH in Jena,
business address: **Brüsseler Str. 10, 07747 Jena,**
registered under **HR B 209648** of the commercial register of the local court of Jena,

based on the written power of disposal, which was given to Mrs. Dania Di Pietro Paolo on 22.08.2019/01.10.2019. I, the undersigned civil law notary, confirm, that the original of the named power of disposal lay before me, when Mrs. Dania Di Pietro Paolo signed. I hereby certify, that the copy of the written power of disposal is a true copy of the original, which has been produced to me.

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.....
Dr. Dania Di Pietro Paolo

Power of Attorney

Asclepion Laser Technologies GmbH, Brüsseler Straße 10, 07747 Jena
represented by Managing Director Dr. Danilo Leggieri
hereby grant

Dr. Dania Di Pietro Paolo

the authority for
conduct the necessary actions and prepare necessary documents related to international approvals, registration and licenses of medical devices
for our company with immediately effect.

The power of attorney is limited to the regularly occurring business.

However, the power of attorney granted does not include the power to negotiate, commitments or make dispositions in financial matters of our Company (except for typically incurred registration, approval and / or licence costs).

This power of attorney does not authorize the conduct of the case.

Legally binding declarations made for our company against third parties or internally require the approval and countersignature by the managing director.

Signed letters and other documents are marked with the addition "with full authority" to draw.

The power of attorney can be revoked at any time without giving reasons.

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