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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HANGZHOU TECHDERM
BIOLOGICAL PRODUCTS CO., LTD. on the annexed DOCUMENT
is genuine.

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验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTIONS FOR USE MEDICAL DEVICE

Drafted by: Jane Chen
Approved by: Chai Xiaohua
Approver: Chai Xiaohua
Revision: A/3
Status: Secrecy
Дата: 08.06.2020

2019

HANGZHOU TECHDERM BIOLOGICAL PRODUCTS CO., LTD
China, Hangzhou, Binjiang District, Jiangling Road, 88, East Building 6

Phone/Fax + 86-571 871 655 26

Web: info@techderm.com

1. Name of the medical device.

Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles «Sofiderm» (hereinafter «Sofiderm»)

Manufacturer:

"Biological products Tehderm Hangzhou" / «Hangzhou Techderm Biological Products Co., Ltd.


China, Hangzhou, Binjiang District, Jiangling Road 88, East Building 6, Floor 2

+ 86-571 871 655 26, info@techderm.com


2. The appearance and structure of the medical device.

2.1 Appearance «Sofiderm» Finelines 1.0 with needles.


PHOTO 1.


2.2 Appearance «Sofiderm» Finelines 2.0 with needles. PHOTO 2.


2.3 Appearance «Sofiderm» Derm 1.0 with needles. **PHOTO 3.**


2.4 Appearance «Sofiderm» Derm 2.0 with needles. **PHOTO 4.**

2.5 Appearance «Sofiderm» Deep 1.0 with needles. **PHOTO 5.**


2.6 Appearance «Sofiderm» Deep 2.0 with needles. PHOTO 6.


2.7 Appearance «Sofiderm» Deep 10.0 without needles. PHOTO 7.


2.8 Appearance «Sofiderm» Derm Plus 1.0 with needles. PHOTO 8.


2.9 Appearance «Sofiderm» Derm Plus 2.0 with needles. **PHOTO 9.**


2.10 Appearance «Sofiderm» Derm Plus 10.0 without needles. **PHOTO 10.**


2.11 Appearance «Sofiderm» Derm Plus 20.0 without needles. **PHOTO 11.**

2.12 Appearance «Sofiderm» Derm Subskin 10.0 without needles. **PHOTO 12.**

2.13 Appearance «Sofiderm» Derm Subskin 20.0 without needles. **PHOTO 13.**

Embodiments:

I. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Finelines 1.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in 1ml syringe - 1 ps .;
2. The single-use sterile needle «BD» 30G 1/2"(0,3mm x 13mm)- 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

II. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Finelines 2.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 2 ml - 1 ps .;
2. The single-use sterile needle «BD» 30G 1/2"(0,3mm x 13mm) - 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

III. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm 1.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 1ml- 1 pc .;
2. The single-use sterile needle «BD» 27G 1/2"(0,4mm x 13mm) - 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

IV. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm 2.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in 2 ml-syringe 1 pc .;
2. The single-use sterile needle «BD» 27G 1/2"(0,4mm x 13mm) - 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

V. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Deep 1.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 1ml- 1 pc .;
2. The single-use sterile needle «BD» 27G 1/2"(0,4mm x 13mm)- 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

VI. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Deep 2.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in 2 ml-syringe 1 pc .;
2. The single-use sterile needle «BD» 27G 1/2"(0,4mm x 13mm)- 2 units .;
3. Operating instructions - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece

VII. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Deep 10.0 without needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in 10 ml syringe - 1 ps .;
2. Instructions for use - 1 copy .;

3. The label information for the patient card - 2 pcs.
4. Consumer Packaging - 1 piece.

VIII. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Plus 1.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 1ml- 1 pc .;
2. The single-use sterile needle «BD» size 25G 5 / 8 "(0,5mm x 16mm) - 2 units .;
3. Instructions for use - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

IX. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Plus 2.0 with needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in 2 ml-syringe 1 pc .;
2. The single-use sterile needle «BD» size 25G 5 / 8 "(0,5mm x 16mm) - 2 units .;
3. Instructions for use - 1 copy .;
4. The label information for the patient card - 2 pcs.
5. Consumer Packaging - 1 piece.

X. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Plus 10.0 without needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 10ml- 1 pc .;
2. Instructions for use - 1 copy .;
3. The label information for the patient card - 2 pcs.
4. Consumer Packaging - 1 piece.

XI. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Plus 20.0 without needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 20ml- 1 pc .;
2. Instructions for use - 1 copy .;
3. The label information for the patient card - 2 pcs.
4. Consumer Packaging - 1 piece.

XII. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Subskin 10.0 without needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 10ml- 1 pc .;
2. Instructions for use - 1 copy .;
3. The label information for the patient card - 2 pcs.
4. Consumer Packaging - 1 piece.

XIII. Implant based on sodium hyaluronate "Sofiderm" sterile in syringes with needles or without needles Derm Subskin 20.0 without needles, composition:

1. Gel Sodium hyaluronate 20mg / ml in the syringe 20ml- 1 pc .;
2. Instructions for use - 1 copy .;
3. The label information for the patient card - 2 pcs.
4. Consumer Packaging - 1 piece.

2.3 Description of the medical device.

«Sofiderm»- a gel of hyaluronic acid of non-animal origin in syringes with needles and without needles. It is a sterile, pyrogen-free, nontoxic gel having a composition of chemically cross-linked polymer of sodium hyaluronate of non-animal origin (biosynthesized) in phosphate-buffered saline, placed in a sterile glass syringe.

Types and composition.

Sterile intradermal implant based on hyaluronic acid "Sofiderm» has five different variants: Finelines, Derm, Deep, Derm Plus, Derm Sub- Skin. Table 1 lists the composition of these products.

Type of	Finelines	Derm	Deep	DermPlus	DermSubSkin
Concentration (mg / ml)	20	20	20	20	20
Filling (ml)	1; 2	1; 2	1, 2, 10	1; 2; 10; 20	10; 20
Grain diametr (μm)	100-150	150-280	280-500	500-1250	1250-2000
The gel content in 1 ml					
Sodium hyaluronate (mg)	20	20	20	20	20
Sodium Chloride (mg)	8.5	8.5	8.5	8.5	8.5
Disodium hydrogen phosphate (mg)	0.03	0.03	0.03	0.03	0.03
Sodium dihydrogen phosphate (mg)	0,126	0,126	0,126	0,126	0,126
Sterile water for injection	AD 1 mL	AD 1 mL	AD 1 mL	AD 1 mL	AD 1 mL
1,4-butanediol-diglycidyl ether (ppm)	≤ 2	≤ 2	≤ 2	≤ 2	≤ 2

Table 1. The composition of the five types of products Sofiderm

2.4. Scope of supplying.

«Sofiderm» supplied in a disposable glass syringe with sterile disposable needles or without needles in the individual, consumer packaged with instructions for use and the label information to the patient card. Table 2 shows the specifications of syringes, needles.

Type of	Filling (ml)	Number of syringes (pcs.)	needle size	Number of needles	Information label for the patient card	Instructions for use	consumer packaging
Finelines	1	1	30 G1 / 2" (0,3mm x 13mm)	2	2	one	one
	2	1	30 G1 / 2" (0,3mm x 13mm)	2	2	one	one
Derm	1	1	27 G1 / 2" (0,4mm x 13mm)	2	2	one	one
	2	1	27 G1 / 2" (0,4mm x 13mm)	2	2	one	one
Deep	1	1	27 G1 / 2" (0,4mm x 13mm)	2	2	one	one
	2	1	27 G1 / 2" (0,4mm x 13mm)	2	2	one	one
	10	1	Not	not	2	one	one
Derm Plus	1	1	25 G5 / 8" (0,5mm x 16mm)	2	2	one	one
	2	1	25 G5 / 8" (0,5mm x 16mm)	2	2	one	one
	10	1	Not	not	2	one	one
	20	1	Not	not	2	one	one
Derm Sub Skin	10	1	Not	not	2	one	one
	20	1	not	not	2	one	one

Table 2. Technical characteristics of syringes and needles

2.5. Technical requirements

Analyzed indicator	Specs
Description	Implant is sterile, viscous, transparent, colorless, homogeneous. No visible particles.
Nominal volume of the syringe	95-105% of the specified filling volume
Size distribution	Qualified particles of at least 80%
Degree of reticulation	Finelines – 7%-9%
	Derm – 7%-9%
	Deep – 8%-10%
	DermPlus – 9%-11%
	DermSub-Skin – 10%-11%
Osmotic concentration	270 - 350 mOsm/kg
Sodium hyaluronate concentration	90-120% of the specified value
Dynamic viscosity	Finelines - 35-50 Pa*s
	Derm – 40-60 Pa*s
	Deep - 55-75 Pa*s
	DermPlus - 70-85 Pa*s
	DermSub-Skin - 85-100 Pa*s
BDDE content (no more than)	2 mkg/g (2ppm)
pH (at 25°C ± 2°C)	6,8-7,5
Heavy metals content (no more than)	5 mkg/g
Particles	No

Sterility	Needle – Radiation (Gamma Radiation); finished product - autoclaving
The content of endotoxins in finished products (not more than)	0,5 U/mg
Protein content (not more than)	0,1%
Hemolytic Streptococcal Hemolysin	Without Hemolysis Ring
The presence of materials derived from animal or human tissue, as well as blood and human blood components	Doesn't content
Peptides, amino acids, vitamins	Doesn't content
Medicaments content	Sodium Hyaluronate No. CAS 9067-32-7, manufactured by Bloomage Freda Biopharm Co., Ltd. (CEP certificate R0-CEP 2015-201-red. 00; address: Jinan, Shandon, China), registry number ΦC - 001305 dated January 11, 2016, containing 1,4-butanediol diglycidyl ether (BDDE) No. CAS 2425- 79-8, manufactured by Sigma-Aldrich Ltd. - Sofiderm gel
Plants extract content	Doesn't content
Molecular weight of sodium hyaluronate (for all embodiments)	2 500 000
The mode of biodegradation	FineLine is up to 5 months, Derm- up to 6 months. Deep- up to 9 months. DermPlus- up to 11 months. DermSubskin- up to 13 months
The retention time of the therapeutic effect	FineLine is up to 7 months, Derm- up to 10 months. Deep- up to 12 months. DermPlus- up to 14 months. DermSubskin- up to 16 months

3. Appointment of a medical device.

3.1. Appointment.

“Sofiderm” is intended for the correction of age-related changes of the skin of the face and body (wrinkles and folds), atrophic defects (scars and cicatrix) after surgical interventions and injuries, as well as for use to change the shape of the face, lips and cheekbones by creating additional volume.

After use, products cannot be reused. Products are sterile, pyrogen-free, non-toxic.

“Sofiderm” is used by trained and qualified doctors in the conditions of medical and medical institutions.

“Sofiderm” is 100% biodegradable in human tissues. Products of biodegradation are not toxic, are involved in biochemical reactions taking place in the body. The duration of the correction depends on the following factors: the type of defect (surface, medium or deep); skin tension and stress at the site of implantation; depth of implant insertion inside the tissue; injection technique, etc.

3.2. Indications. Introduction technique

To achieve the effect of a single procedure. The number of procedures is determined by a specialist to achieve the desired effect. It must be remembered that with repeated corrections, the duration of the aesthetic result, as a rule, increases as the newly introduced implant is layered on the remnants of the previous implant. It is recommended to repeat the procedure no earlier than 2-4 weeks after the previous one. A linear-retrograde or bolus implant insertion technique is provided.

3.2.1 "Sofiderm" Finelines.

Sofiderm Finelines can be injected into the skin to effectively remove superficial wrinkles and moisturize the skin, especially the crow's feet, superficial forehead and neck wrinkles ("rings of Venus").

Use mainly linear-retrograde technique.

3.2.2. "Sofiderm" Derm.

Sofiderm Derm is indicated for the correction of atrophy associated with aging; Correction of ptosis of facial tissues, decolletage, omission of the corners of the lips.

Use linear-retrograde and bolus technique.

3.2.3: "Sofiderm" Deep.

"Sofiderm" Deep is used in injection cosmetology and plastic surgery with the purpose of augmentation of the zygomatic region and chin, the formation of facial contours, filling of the nasolacrimal sulcus, replenishment of the lost volume in the temporal region. It is also used to correct scars caused by acne and chickenpox, atrophic scars caused by surgery, trauma and infections, as well as congenital or involutive atrophy of soft tissues.

Use mainly bolus technique.

3.2.4. "Sofiderm" Derm Plus and Derm Sub Skin.

"Sofiderm" Derm Plus and Derm Sub Skin is intended for the correction of age-related changes in the skin of the face and body (wrinkles and folds), atrophic defects (scars and scars) after surgical interventions and injuries, as well as for applying to change the shape of the face, lips, cheekbones by creating additional volume. Use the bolus injection technique.

3.2.5. The maximum possible dose of administration in one procedure for each performance;

On practice:

Fineline, Derm, - up to 10 ml; Deep, DermPlus - up to 30-50ml; Derm Subskin - up to 100 ml

3.2.6. The ability to remove, replace the implant.

Not provided

3.2.7. Interaction with other implants or injected drugs

There is no data on the incompatibility of the hyaluronic acid gel with permanent implants with a correctly performed procedure (i.e., if the integrity of the permanent implant is not violated).

The incompatibility between hyaluronic acid and quaternary ammonium compound, such as, for example, benzylalkonium chloride, is well known. With this in mind, the implant should never be allowed to come into contact with such compounds, as well as medical surgical instruments treated by it.

Incompatibility with other absorbable implants (for example, thread for lifting) is absent.

3.3. Contraindications.

Sterile intradermal implant based on hyaluronic acid "Sofiderm» is contraindicated in:

- For patients with a history of hypertrophic or keloid scars
- For patients with inflammatory diseases of the skin, the active streptococcal infection.
- For patients with a history of anaphylactic reactions to.
- For patients taking anticoagulants and platelet aggregation inhibitors.
- For use in areas containing foreign matter or nondegradable implants (e.g., silicone).
- for use in inflammatory and / or infectious process (acne, herpes and etc.).
- For use in conjunction with laser treatment, deep chemical peeling and dermabrasion.
- For patients with autoimmune diseases receiving immune therapy.
- For patients with known sensitivity to hyaluronic acid.
- For patients with porphyria.
- For patients with epilepsy.
- For pregnant women and during lactation.
- For use in patients younger than 18 years.

3.4. Warning.

- Check the integrity of the syringe and the expiry date before use
- Do not use after the expiration date.
- Do not use if packaging is damaged
- Do not reuse
- Do not inject into blood vessels and muscles.
- Expose resterilized.
- Do not modify the product; changes may affect the sterility and safety.
- Do not use if the product may be damaged.
- Do not use if the gel particles visible, changes in color or turbidity.
- Do not mix with other products.
- Do not use in the presence of symptoms of local inflammation of the skin

3.5. Possible side effects.

Patients should be informed about the side effects and possible after implants injection with hyaluronic acid, which may occur immediately or after a delay.

These include, but are not limited to:

- Redness, erythema, swelling at the injection site that can be associated with response to the injection and will stop spontaneously without intervention within hours / days.
- Hematoma
- Induration in site of injection.
- Bad effect or little effect filling and loss of sensation at the injection site.
- There is a potential risk that the product may be introduced inadvertently entered the blood vessels of the dermis that may cause occlusion of arteries, with the wave characteristic of consequences. No cases have been reported on this product.
- Treatment reactions generally occur within 24-48 hours. Other side effects are also possible within 2-4 weeks.
- Patients should be advised to immediately report any evidence of any adverse effects not mentioned in this leaflet doctor.
- Any other undesirable side effects associated with the injection and fixed specialist should be communicated to the authorized representative and / or directly to the manufacturer.

3.6 Procedure for the application.

➤ recommendations

It is necessary to strictly observe the rules of asepsis. The product is ready for use. The implant should be injected intradermally into the dermis under wrinkles, increasing the volume, thus filling it. When the correction of wrinkles line can remain visible even after its completion.

The result of the procedure depends on the use of appropriate injection method, which can only choose a qualified medical professional.

Expected results of the procedure should be evaluated for each patient individually. Patients with unrealistic expectations should not undergo this procedure.

Before the procedure should always check the patient's medical history to make sure there are no contraindications.

The patient should be informed about all of the indications, contraindications, warnings and precautions.


It should also identify the need for analgesics.

To obtain and maintain the maximum desired result the patient should be informed of the need for additional procedures or procedures.

Needle and syringe assembly

To avoid separation of the needle from the syringe during insertion of the implant is necessary to make the following manipulations:

- remove the cap from the tip, as shown in Figure A, and attach the needle from the kit as shown in Figure B and C.

	Gently in the direction of counter-clockwise motion unscrew the cap covering the Luer-Lock syringe adapter.
	Align the connector needle and the syringe Luer-Lock adapter push and rotate the needle threading operation (at clockwise direction) by a quarter turn
	Make sure that the needle is securely locked in the adapter. Remove the protective cap from the needle.

Preparing for using

The medical device can be used only by a medical personnel with appropriate qualifications and possessing technique of intradermal injection.

In severe pain response in a patient may use a local anesthetic cream.

Skin injection site twice disinfect alcohol or other antiseptic solution to prevent the risk of infection in the course of administration of the drug.

Do not enter the implant too deeply, as this will shorten the resorption of the implant, and thus - term effect of the procedure.

Do not enter the implant too close to the surface, as this may cause noticeable bumps and (or) a change in color. Slowly inject the drug into the dermis, pulling the needle back.

Complete injection for a couple of millimeters to a needle insertion point to the skin to prevent leakage of the implant from the injection point.

For each injection site recommended maximum dose in a 2 ml injection per session.

After insertion of the implant can be easily massage the skin, to evenly distribute the product.

If injectable implant was stored below room temperature, it is necessary to use the product stand at room temperature for at least 30 minutes.

ATTENTION! Any adverse events should be reported to the authorized representative of the manufacturer or directly to the manufacturer

3.7 Conditions of use.

Sterile intradermal implant based on hyaluronic acid "Sofiderm» used in a medical health care institutions.

3.8 Precautions for use.

- Sterile intradermal implant based on hyaluronic acid "Sofiderm» forbidden to use for patients with contraindications.
- Carefully check the integrity of the packaging of the syringe and needle, the shelf life and sterility marker before using the product. Don't use Sterile intradermal implant based on hyaluronic acid "Sofiderm» if the packaging is damaged.
- Prohibited use medical device if the integrity of the blister pack is broken, if any visible signs of damage, if the tip cap to the syringe Luer adapter loosely twisted, or damaged. The gel in the syringe is colorless and transparent, so when you change the color or transparency of the contents of the syringe to apply the product is prohibited.
- Use immediately after opening the package aseptically and precautionary measures applicable to traditional methods of intradermal injections.
- Sterile intradermal implant based on hyaluronic acid "Sofiderm» should be used only as a means of intradermal and should only be administered by qualified trained physicians or other medical personnel, in accordance with applicable laws.
- Never mix Sterile intradermal implant based on hyaluronic acid "Sofiderm» » with other drugs.
- Prohibited administered Sterile intradermal implant based on hyaluronic acid "Sofiderm» into the blood stream.
- Prohibited administered Sterile intradermal implant based on hyaluronic acid «Sofiderm» directly next to the permanent implant, as well as concurrently with other ongoing or planned procedures (laser, micro-current treatment etc.).
- Do Sterile intradermal implant based on hyaluronic acid «Sofiderm» patients who at the injection site there is an infection or inflammation.
- Do not enter the sterile intradermal implant based on hyaluronic acid «Sofiderm» pregnant or lactating women and children.
- Note: The sterile intradermal implant based on hyaluronic acid «Sofiderm» can reactivate latent herpes virus or cause angioedema shock.
- Warn patient that using them drugs acting on platelet function, such as antiinflammatory drugs or aspirin, the patient may experience significant bleeding or bruising at the injection site.
- Studies on the application of sterile intradermal implant based on hyaluronic acid «Sofiderm» patients already implanted medical products based on hyaluronic acid was conducted.
- The impact of the product on the vision and the reaction rate has not been studied, so no special instructions for drivers or people working with complex technological devices, does not exist.
- The product is used in invasive procedures involving the risk of infection. You must follow standard safety precautions associated with using sterile materials / products.

Information for patients.

- Avoid taking aspirin, anti-inflammatory agent or vitamin C during the week preceding the injection.
- There may be a little bleeding after the injection.
- Tell your doctor if you are undergoing anticoagulant treatment or if you are suffering from hemostatic disorders (increased risk of hematoma).
- Within ten (10) days, avoid overheating and damp places (swimming pool, gym, Turkish bath, Jacuzzi and so on.).
- Within a few days, avoid any activity that puts pressure on the treated area.
- If there is pain at the injection site that lasts longer than a week, see your doctor.

Care after procedure.

- Patients are advised to avoid overwork, sunlight, burns or extreme temperature conditions for 24 hours.
- Patients were prohibited to apply makeup for 12 hours before and after treatments.
- Patients should impose ice or cold compress in the area of treatment in 24 hours after the procedure.
- Patients are advised to limit conversations, movements of the face and a laugh in a week after the procedure.

Procedure:

The method of application of the medical device

Preparing for use

The medical device can be used only by a medical personnel with appropriate qualifications and possessing technique of intradermal injection.

In severe pain response in a patient may use a local anesthetic cream.

Skin injection site twice disinfect alcohol or other antiseptic solution to prevent the risk of infection in the course of administration of the drug.

Do not enter the drug too deeply, as this will shorten Sterile intradermal implants based on hyaluronic acid «Sofiderm», and therefore - term effect of the procedure.

Do not enter the implant too close to the surface, as this may cause noticeable bumps and (or) a change in color.

Slowly inject the implant into the dermis, pulling the needle back.

Complete injection for a couple of millimeters to a needle insertion point to the skin to prevent leakage of the drug from the injection point.

For each injection site recommended maximum dose in a 2 ml injection per session.

After the introduction of Sterile intradermal implant based on hyaluronic acid «Sofiderm» can easily massage the skin, to evenly distribute the product.

If Sterile intradermal implants based on hyaluronic acid «Sofiderm» stored at below room temperature, it is necessary to use the product stand at room temperature for at least 30 minutes.

ATTENTION! Any adverse events should be reported to the authorized representative of the manufacturer or directly to the manufacturer.

Stop the procedure immediately if there is a possibility of vascular puncture.

If during the injection process there is considerable resistance to move the injection needle deeper implant to facilitate precipitation. If the resistance persists, try the injection of another entry point. In the case of continuous resistance may require replacement of the injection needle and / or syringe.

Do not apply excessive pressure on the syringe at any time.

If there is an immediate blanching stop injection and massage the area as long as the color is corrected. If normal skin color is not returned or renewed injection process.

Surface injection or introduction of a large amount of gel can result in discoloration and the formation of nodes.

Stop injection before pulling the needle out of the skin to avoid leakage of the gel from the injection site.

Repeat the procedure if you need more adjustment, but only after a careful assessment of the treatment and the patient's condition.

After completion of the injection region is important massaging treated after injection, to ensure uniform distribution of the gel.

If there was overcorrection strongly massage the area to obtain optimal results.

Facials after injection :

Patients should be advised to avoid extreme activities, exposure to sunlight and sunburn or extreme weather conditions within 24 hours.

Patients should be instructed to use a mask or ice cold compress to the treated area for 24 hours after the procedure.

Patients should be advised to limit talking, facial expressions and laughter within one week after the procedure.

Patients who have the procedure in the mouth, it is recommended to rinse your mouth with saline every 3-4 hours was carried out within one week after treatment.

If necessary, you can use oral analgesics.

4. Labeling

4.1 Each consumer (individual) package is the following:

- product name and number and date of the registration certificate;
- Embodiment;
- Sterilization method;
- Special work instructions and warnings;
- the name and address of the manufacturer «Sofiderm»
- the name and address of the manufacturer of needles;
- the name and address of the authorized representative of the manufacturer;

- date of manufacture;
- the date of the expiration date of the product;
- Lot code;
- storage conditions;
- the number of syringes and volume of «Sofiderm» in each syringe;
- the number of needles and their characteristics;
- sodium hyaluronate concentration in mg/ml;
- the inscription "sterile", "Non-toxic";
- symbols in 4.8

4.2 Each Stickers for patient cards applied the following information:

- the name of embodiment
- the date of the expiration date of the product;
- Lot code;
- catalog number;
- the number of syringes and volume of «Sofiderm» in each syringe;
- symbols in 4.6

4.3. The needle marking includes the following information / symbols:

- model name;
- needle size (diameter x length);
- CE marking;
- marking "Sterile. Radiation sterilization";
- number according to the catalog, series number, expiration date;
- "Use only once" symbol;
- Symbols in 4.7

4.4. The following information is applied to each syringe:

- name and website of the manufacturer "Sofiderm";
- expiration date of the product;
- Lot code;
- CE marking
- Short name of embodiment
- Volume.

4.5. The following information is applied to each blister pack:

- Embodiment
- Volume
- Composition
- Sodium hyaluronate concentration
- Storage condition
- Sign "sterile" "non toxic"
- name and address of the manufacturer "Sofiderm"
- name and address of the authorized representative of the manufacturer in the EU;
- symbols of clause 4.8.

4.6. At each Stickers for patient cards applied the following characters:

4 6. Explanation of symbols on the label

The prohibition on re-use

The expiration date of the product shelf life

Sterilization method - steam or dry heat

LOT

REF

Lot code

Manufacturer

Catalog number

Do not re-sterilize

4.7. Explanation of symbols on the needle marking

Needle marking	
	Manufacturer
LOT	LOT number
STERILE R	method of sterilization - Radiation
EXP	Expiration date
	The prohibition on re-use
CE 2460	CE marking

4.8 Each consumer (individual) package is the following characters:

4.8. Explanation of symbols on the label and blister (EN ISO 15223-1-2012)

The prohibition on re-use

The expiration date of the product shelf life

STERILE I

Sterilization method - steam or dry heat

STERILE R

method of sterilization - Radiation

LOT

Lot code

date of manufacture


Manufacturer

EC REP

Authorized Representative in the European Community

Refer to the instructions for use

temperature limits

Do not re-sterilize

Do not use if packaging is damaged

Avoid exposure to sunlight

Protect from moisture

CE marking and the technical expertise Record ID Number Authority.
Certification according to MDD (93/42 / EEC)

Example Consumer labeling:

Sterile intradermal implant of hyaluronic acid «Sofiderm»
An embodiment of the Derm 1.0 with needles
Syringe 1/1 **STERILE I**

LOT Party № 180724-1

The date of manufacture 06.2018; Use up to 06.2020
Registration certificate number _____ date _____

Name and address of the manufacturer:
 " Hangzhou Techderm Biological products Co.Ltd»
China, Hangzhou, Binjiang District, Jiangling Road 88, East Building 6, Floor 2


Needles: 2x27G 1/2"(0,4mm x 13mm) **STERILE R**
Name and address of the manufacturer of needles:
 Becton Dickinson France SAS
Becton Dickinson France SAS, rue Aristide Berges, 38801 Le Pont-de-Claix, Cedex, France

EC REP **Lendlink GmbH**
Dorfstrasse 2/4, 79312, Emmendingen, Germany


Authorised representative: "TECHDERM" LLC.
107370, Russia, Moscow, Rokossovsky B-d, 6, Building 1, Apartment 26

Keep in the packaging at a temperature + 2 ° C to + 25 ° C,
and humidity of 10% to 80%, do not freeze.
STERILE! Non-pyrogenic! Non-toxic!
The concentration of sodium hyaluronate 20mg / ml

Example labeling for patient card:

Example consumer labeling for blister:

5. Storage and shelf life

"Sofiderm" is kept in a clean, dry, dark place out of the reach of children in the manufacturer's packaging at a temperature of from 2 to 25 and humidity from 10% to 80%. Do not freeze. Protect from sunlight. Transportation of a medical product is performed by any type of covered transport in accordance with the requirements for transportation and the current regulations for each type of transport at a temperature of 2 ° C to 25 ° C. During storage and transportation it is forbidden to squeeze the package. Avoid direct sunlight, freezing, and condensation. "Sofiderm" shelf life - 2 years from the date of manufacture indicated on the package.

6. Procedure for disposal and destruction

Empty syringes, packaging and «Sofiderm» in violation of the integrity and / or hermetically sealed, with a change of color and / or texture, expired, not used after opening the package, be disposed of solid waste, in accordance with all federal, state / provincial and local laws and regulations. Used products are disposed of as potentially infectious class B medical waste.

7. Warranty

Manufacturer - " Hangzhou Techderm Biological products Co. Ltd» China, Hangzhou. Binjiang District, Jiangling Road 88, East Building 6, Floor 2

authorized representative of the manufacturer - "TECHDERM RUSSIA" LLC (107370 Moscow, B-d Marshall Rokossovskogo, 6, Bldg 1, app. 26), provided that Sofiderm corresponds to the declared properties during shelf life and conforms to the rules of storage, use and transportation.

Period of validity of 2 years from the production date printed on the package.

Monitoring activities after the sale is carried out according to SOP QA-24 "Control after the sale."


8. The maintenance requirements and repairs of medical products. Service maintenance

Product after application cannot be used again, so that does not require the service, maintenance and repair.

9. Complaints

Organization authorized to accept complaints from consumers:

LLC "TECHDERM RUSSIA", 107370, Moscow, B-d Marshall Rokossovskogo, 6, building 1, Apt. 26

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