
Instructions for Use
MAK-NV
Dated July 27, 2020
Consisting of 14 pages

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
Instructions for use
MEDICAL DEVICE

«The non-vascular access kit MAK-NV™»

Manufactured by
Merit Medical Systems, Inc., the USA

DEVICE DESCRIPTION:

The non-vascular access kit MAK-NV™ consists of a coaxial introducer sheath (a pair of coupled introducer sheath / dilator), a introducer needle with a stylet and a mini guidewire, used to provide percutaneous introduction of a non-vascular Seldinger catheter into a body cavity or lumen during interventional radiology procedures, including drainage.

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- Introducer needle with a trocar stylet, diameter: 21 G (0,8 mm), length: 15 cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1,33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.(if necessary);
- guidewire, diameter 0,038" (0,97 mm), length: 80 cm, 150 cm - 1 pc. (if necessary);
- Instructions for use.

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INDICATIONS FOR USE

The non-vascular access kit MAK-NVTM is intended for percutaneous placement of a 0,035" (0,89 mm) or 0,038" (0,97 mm) guidewire with the introduction of a non-vascular catheter by the Seldinger method into a cavity or lumen of the body.

CONTRAINDICATIONS

There are no contraindications.

Px only: **DEVICE RESTRICTED SALES**

For professional use only.

Attention! The Federal (USA) device authorizes this sale only to or on the order of a physician trained and / or experienced in the use of the data.

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PRECAUTIONARY MEASURES

- The kit should only be used for non-vascular procedures.
- Read the manufacturer's instructions before use.
- The contents of unopened and undamaged packaging are sterilized using ethylene oxide and pyrogen-free.
- Do not use if packaging is open, damaged, or broken.

- Use for only one patient. Reuse, handling or sterilization is prohibited. Reuse, handling or sterilization may impair the integrity of the device and / or cause the device to malfunction and, as a result, result in injury, illness or death of the patient. Reuse, handling or sterilization can also create a risk of contamination of the device and / or lead to infection of the patient or cross-infection, including but not limited to the transfer of infectious diseases from patient to patient. Contamination of the device can be harmful to the patient, resulting in illness or death.
- Observe generally accepted precautions when introducing and maintaining this device. Since there is a risk of contracting bloodborne pathogens, healthcare providers should always follow standard precautions when performing procedures with blood and body fluids. Sterility should always be observed.
- Do not use after the expiration date.
- Store in a cool, dry place.

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- edema
- blood clot formation
- infection
- pseudo-aneurysm formation
- hematoma
- perforation or rupture of a vessel or internal organ
- bleeding
- embolization of the conductor or catheter
- hemothorax
- inflammation, necrosis, or scarring
- hydrothorax
- pain in the introduction area
- extravasation

INSTRUCTIONS FOR USE:

1. Determine the injection site and prepare the injection site using a proper aseptic technique and, if required, anesthesia.
2. Enter the 21-caliber Introducer needle using standard technology.
3. Carefully slide forward end of the flexible conductor diameter of 0.018 inches (0.46 mm) through the needle.
Slide the guide forward as far as it is necessary. Make sure the placement.
WARNING. Do not move forward in the case of the conductor resistance.
4. Remove the needle retaining location conductor diameter 0.018 "(0.46 mm).
WARNING. Conclusion Extruding back or manipulation of the distal tip of the conductor through the needle tip may cause damage or embolization.
To avoid damage to the conductor during manipulation with it, remove Introducer needle and proceed to the next step.

5. Enter coaxial introducer with dilator, with a rigid cannula or not, through a wire diameter of 0.018 "(0.46 mm) and push forward to the desired position.
6. Remove the dilator / rigid cannula, leaving the sheath in place. Guidewire diameter of 0.018 "(0.46 mm) can be removed or left in place.
7. Enter the diameter of the guidewire 0.038 "(0.97 mm) through the introducer.
8. Remove the sheath, leaving the wire in place.

Warning about the inadmissibility Reapply

The product is intended for single use only. Do not allow the re-use, recycling and sterilization products. Repeated application, processing and sterilization can disrupt the structural integrity of the product and / or lead to failure, which in turn may create a risk of injury, illness or death of the patient outcome. Repeated application, processing and sterilization may also create a risk of contamination of the product and / or cause an infection of the patient or cross-infection, including, but not limited to, the transmission of infection (infection) from one patient to another. Product contamination can lead to injury, illness or death of the patient outcome.

**ANNEX TO THE INSTRUCTIONS FOR USE
MEDICAL DEVICES**

«The non-vascular access kit MAK-NV™»

**Manufactured by
Merit Medical Systems, Inc., the USA**

1. DEVICE DESCRIPTION

The non-vascular access kit MAK-NV™ consists of a coaxial introducer sheath (a pair of coupled introducer sheath / dilator), a introducer needle with a stylet and a mini guidewire with a diameter of 0,018" (0,46 mm) and / or a diagnostic guidewire with a diameter of 0,038" (0,97 mm). The non-vascular access kit MAK-NV™ is available with an introducer needle with a trocar or Chiba stylet, a mini guidewire, and / or a diagnostic guidewire.

The composition of the medical device is presented below.

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The coaxial introducer is equipped with a Luer type connector with internal thread, which allows the system to function as a single product. A conductor needle with a stylet provides support for a pair of introducer / dilator during installation. The distal end of the introducer sheath is equipped with a radiopaque marker, which facilitates the visualization of the tip during fluoroscopy

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3. PRODUCT SPECIFICATIONS

3.1 Guidewire needle with stylet, diameter: 21 G (0,8 mm), length: 15 cm, 20 cm

The outer diameter of the needle	21 G (0,800 mm)
needle length	15 cm 20 cm
Inner diameter of a needle	0.580 mm
Tube wall type	thin-walled
Needle color coding	Dark green
Working length of the stylet	170 mm (with 150 mm (15 cm) needle tube length) 220 mm (with a needle tube length of 200 mm (20 cm)) The stylet must be longer than the needle when the two components are assembled into a single structure.
The outer diameter of the stylet	0,500 mm
Tip shape stylet	stylet trocar stilet type Chiba (chiba)
Sharpening angle of the stylet trocar type	28°
Stylet apex angle, trocar type	35°
Stiletto sharpening shape trocar type	pyramidal
Sharpening angle of the stylet type Chiba	32°
Surface	When carrying out a visual check on the external surface of the working length of the device no foreign substance or surface defects. On needlepoint no protrusions growths, hooks etc.
The protruding portion of the needle tip	In the assembled structure cannula does not extend beyond the tip of the stylet, or its length does not exceed the length of the stylet more than 2 mm.
Explorer compatibility	Gidewire diameter of 0.018 inches (0.46 mm)
Radiopacity	The distal end of the Introducer needle is equipped with a radiopaque marker, 10mm in length, which facilitates the visualization of the tip during fluoroscopy.
The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)	Length of female conical element: 7,500 -10,500 mm Luer Taper: 6% Smaller outer lip diameter: 5,515-6,730 mm Larger outer ledge diameter: 7,730-7,830 mm

3.2 Coaxial sheath diameter 6F

External diameter	The length of the exposed portion	Length products
6F (2,00 mm)	6 mm	260 mm

- dilator diameter 4F,

External diameter	The length of the exposed portion	Length products
4F (1,33 mm)	6 mm	250 mm
Compatibility with the guidewire of 0.038 inches (0.97 mm)	Must pass through the dilator conductor 0.038 inches (0.97 mm)	
The geometric shape of the distal end of the dilator	The presence of tapering, polished the edges. Edge should be sharp. Without cracks, scratches.	
radio-opacity	It should be viewed with fluoroscopy	
The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)	Length of female conical element: 7,500 -10,500 mm Luer Taper: 6% Smaller outer lip diameter: 5,515-6.730 mm Larger outer ledge diameter: 7,730-7,830 mm	
The main characteristics of the connector (connector)) of the locking (plug) type Luer-Lock (according to EN ISO 80369-7: 2016)	Maximum internal thread diameter (root diameter): 8,000 mm Minimum internal thread diameter (thread crest diameter): 7,000 mm Luer Lock Taper: 6% Length of plug cone (e): 7,500 -10,500 mm	

- a metal cannula is rigid,

Characteristic	
The outer diameter of the cannula	0.810 mm
The length of the exposed portion	238 mm
The shape of the cannula tip	At the distal end of the cannula tip should be straight
Surface	On the outer surface should be free of surface defects and process.
The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)	Length of female conical element: 7,500 -10,500 mm Luer Taper: 6% Smaller outer lip diameter: 5,515-6.730 mm Larger outer ledge diameter: 7,730-7,830 mm
The main characteristics of the connector (connector)) of the locking	Maximum internal thread diameter (root diameter): 8,000 mm Minimum internal thread diameter (thread crest

(plug) type Luer-Lock (according to EN ISO 80369-7: 2016)	diameter): 7,000 mm Luer Lock Taper: 6% Length of plug cone (e): 7,500 -10,500 mm
--	---

3.3 guidewire diameter of 0.018 "(0.46 mm), 60 cm long

The diameter of the guidewire	0.018 " (0.46 mm)
The length of the guidewire	600 mm

3.4 guidewire diameter of 0.038 "(0.97 mm) and a length of 80 cm/ 150 cm, c J-shaped tip

The diameter of the guidewire	0.038" (0.97 mm)
The length of the guidewire	1500 mm 800 mm

Material Specification

Component	Description of materials
Needle wires 21G	
Luer-type connector	polycarbonate CALIBRE™ Megarad™ 2081
Needle	Stainless steel AISI 304
Stylet trocar / stylet type Chiba	Stainless steel alloys AISI 304 /AISI 316
Sleeve	polycarbonate CALIBRE™ Megarad™ 2081 Dye - pelletsCC10130299 TRANS GREEN HALF PC PolyOne™ Color: Green, classification code Pantone 354U
coaxial introducer	
Luer connector	Polyester polycarbonate blend Xenoy HX5600HP, Pigment - titanium oxide solid fraction Xenoy EXN X5600HP Color: bright white, code classification Pantone 11-0601 TCX Green ink INK-PG Pad Printing TYPE UV 2002 IWT
nut	The thermoplastic ABS resin (ABS) Dye - pellets 10200443 TRANS GREEN PolyOne™
dilator tube	polyethylene copolymer Alathon H5618, barium sulfate Pigment - titanium oxide solid fraction A32329C PolyOne™ Color: bright white, code classification Pantone 11-0601 TCX
The rigid cannula	Stainless steel alloys AISI 304 /AISI 316
guidewire diameter diameter 0.018 "(0.46 mm) diameter and 0.038" (0.97 mm)	
	Stainless steel alloys AISI 304 /AISI 316 Polycarbonate CALIBRE™ Megarad™ 2081 Polypropylene Bormed HD810MO
guidewire diameter diameter 0.018 "(0.46 mm)	
	Nitinol- platinum Inconel X Polycarbonate CALIBRE™ Megarad™ 2081 Polypropylene Bormed HD810MO

4. DESCRIPTION OF PACKAGING WITH INDICATION OF NUMBER OF PRODUCTS AND ACCESSORIES IN PACKING

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The non-vascular access kit MAK-NV™

1. The non-vascular access kit MAK-NV™, consisting of:

- A introducer needle with a trocar stylet, diameter: 21 G (0,8 mm), length: 15 cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1,33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018 " (0,46 mm), length 60 cm - 1 pc., (if necessary);
- guidewire, diameter 0,038 " (0,97 mm), length: 80 cm, 150 cm - 1 pc., (if necessary).

2. The non-vascular access kit MAK-NV™, consisting of:



- A introducer needle with Chiba stylet, diameter: 21 G (0.8 mm), length: 15cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1,33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018 " (0,46 mm), length 60 cm - 1 pc., (if necessary);
- guidewire, diameter 0,038" (0,97 mm), length: 80 cm, 150 cm- 1 pc., (if necessary).

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
5. INFORMATION ON LABELING OF PRODUCTS

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	Symbol	Decoding the symbol
1.		Authorized Representative in the European Community
2.		Manufacturer
3.		It sold only to doctors and hospitals
4.		A ban on re-use. Do not re-sterilize
5.		Number REF.
6.		use up

7.		lot number
8.		sterilized using ethylene oxide
9.		Attention! Refer to the instructions for use
10.		pyrogen-free
11.		Compatibility with the conductor
12.		Do not use if packaging is damaged
13.		Protect from moisture
14.		This side up
15.		Fragile, handle with care
16.		date of manufacture
17.		Conformity mark
18.		Thin-walled needle tube lettering
19.		The inscription "non-toxic"

6. Compliance with Relevant Regulations and Technical Standards

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The qualification activities that were performed for the product were based on standards that were in effect at the time of the testing. Standards applied include:
 USP United States Pharmacopoeia. ISO 13485, ISO 14971, EN 556, ISO 15223-1, ISO 11135, ISO 11607-1, ISO 11607-2, ISO 10993-1, EN 1617, ISO 594-1, ISO 594-2, EN ISO 80369-7 : 2017.

7. TRANSPORTATION, STORAGE AND HANDLING

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Usage environment:

Ambient temperature not exceeding 25 °C.

Temperature range for use - the body temperature of 32 °C to 42 °C

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9. Limited Warranty and Disclaimer

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device The non-vascular access kit MAK-NV™ is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial

number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

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Limited Liability Company "Merit Technologies" (LLC "Merit Technologies"),

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

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d

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Cover Letter
Non-Pyrogenicity symbol
Dated July 27, 2020
Consisting of 16 pages

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

Notary Public

My Commission Expires: June 27, 2021

The Russian Federal Service
for Surveillance in Healthcare (Roszdravnadzor)

Concerning the non-pyrogenicity
symbol on the product label

Cover letter


Merit Medical Systems, Inc. (further - MMSI) located at 1600 West Merit Parkway South Jordan, Utah, 84095, USA in response to the Notification No. 10-31233/20 dated June 10, 2020 claims that all risk class 2a and higher medical devices, manufactured by the Company, are mandatory tested for biological safety, including pyrogenicity.

MMSI asks The Russian Federal Service for Surveillance in Healthcare to attach the reports on the biological safety assessment tests for the medical device «The non-vascular access kit MAK-NV™» kit (Appendix No. 1) to the regulatory dossier (income No. № 31711 dated 14.05.2020)

The reports state that MMSI's product "The non-vascular access kit MAK-NV™" was tested for systemic toxicity according to ISO 10993-11 (including rabbit pyrogen test*). Based on the protocol criteria and the data obtained, the products meet the requirements for pyrogenicity and are considered as non-pyrogenic.

Merit Medical Systems, Inc., guided by the principles of confidence and transparency and privileges the safety of manufacturing products, informs potential users about the non-pyrogenicity of The non-vascular access kit MAK-NV™ device and marking symbol of pyrogenicity(~~XX~~) on the label of the product.

With the kind regards

Cory Marsh
Manager, Regulatory Operations
27 July 2020

* The pyrogenicity test was developed to determine the presence of chemical pyrogens in extracts of solid materials in order to limit the risks of a febrile reaction in a patient to an acceptable level after the product introduction. The test involves measuring the temperature rise of the New Zealand white rabbits after intravenous injection of the testing product's extract.

APPENDIX 1
PYROGENICITY TEST-REPORTS FOR
The non-vascular access kit MAK-NV™

Bacterial Endotoxins Test Final Report

Study Number: 1135287-S01

Test Article: MERIT MINI ACCESS KIT
Part N. MAK-NV-004/A
Lot. I1481596

Purchase Order: 408662

Study Received Date: 21 Dec 2018

Testing Facility: Nelson Laboratories, LLC
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.Test Procedure(s): Standard Test Protocol (STP) Number: STP0046 Rev 15
Customer Specification Sheet (CSS) Number: 201705958 Rev 2

Deviation(s): None

Summary: The Bacterial Endotoxins Test (BET), or *Limulus* Amebocyte Lysate (LAL) test, is an *in vitro* assay to detect and quantify bacterial endotoxin, a component of the cell wall of Gram negative bacteria. Standard controls and a positive product control (PPC) demonstrate a compliant assay. A PPC recovery within the 50%-200% range indicates that the test solution is free of interfering factors given the specific conditions of the test. If applicable, dilutions are calculated into the reported endotoxin level. All test method acceptance criteria were met.

The testing was conducted in accordance with the following regulatory documents: ANSI/AAMI ST72:2011/(R)2016, USP <161>, USP <85>, EP 2.6.14, and JP 4.01. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Kinetic Turbidimetric Results:

Sample ID	Extraction Volume	Bacterial Endotoxin	PPC Recovery
10 (Pooled)	59 mL/device	<0.00500 EU/mL or <0.297 EU/device	125%

Endotoxin Limit: The endotoxin limit provided by the sponsor is not more than 20 EU/device.

Maximum Extraction Volume: With an endotoxin limit of 20 EU/device, tested with an assay sensitivity of 0.005 EU/mL, the maximum amount of extraction liquid that can be used is 4000 mL/device.

Darcy M. Rawson electronically approved for

Study Director

Emily Spackman, B.S., RM(NRCM)

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com

28 Dec 2018 07:12 (+00:00)

Study Completion Date and Time

Page 1 of 2

Test Method Acceptance Criteria: The following conditions were met as part of a compliant assay:

- The absolute value of the correlation coefficient of the standard curve is ≥ 0.980
- The blank value is less than the endotoxin detection limit of the lysate reagent employed
- The coefficient of variation (CV) between replicates is $\leq 10\%$

Validation: The validity of the test requires demonstration that the test solution does not inhibit or enhance the assay. Validation of the kinetic BET is accomplished with the PPC. The PPC recovery must be between 50-200%. A PPC recovery of $<50\%$ suggests inhibition and a recovery of $>200\%$ suggests enhancement.

Preparation: The extraction was performed by immersing the test article in endotoxin free water and placing it on an orbital shaker in an incubator for not less than 60 minutes at $37-40^{\circ}\text{C}$.

Test Procedure: The assay was performed at a sensitivity of 0.005 EU/mL using Charles River reagents. For valid results, the combination of lysate and test solution was verified to be within 6-8 pH units.

TOXIKON FINAL REPORT: 03-3171-G6
RABBIT PYROGEN TEST (MATERIAL MEDIATED) - ISO

Author
Kathleen Krenzer O.D., Ph.D.

Final Report Date
July 16, 2003

MANAGEMENT OF THE STUDY

Performing Laboratory
Toxikon Corporation
15 Wiggins Avenue
Bedford, MA 01730

Sponsor
Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, UT 84094

REVIEWED BY:

SIGNATURE:

DATE:

[Signature]
7/21/03

TABLE OF CONTENTS

Title Page

Table of Contents

Study Summary

Quality Assurance Statement

Study Director Signature and Verification Dates

- 1.0 Purpose
- 2.0 References
- 3.0 Compliance
- 4.0 Identification of Test and Control Articles
- 5.0 Identification of Test System
- 6.0 Justification of Test System and Route of Administration
- 7.0 Experimental Design
- 8.0 Dosage
- 9.0 Evaluation Criteria
- 10.0 Results
- 11.0 Conclusion
- 12.0 Records
- 13.0 Confidentiality Agreement
- 14.0 Policy on Pain and Suffering in Animals
- 15.0 Animal Usage

Table I: Pyrogen Test Data

STUDY SUMMARY

The 0.9% USP Sodium Chloride for Injection (NaCl) extract of the test article, Mini Access Kit, was evaluated for its potential to produce a pyrogenic response when tested in New Zealand White rabbits. Based upon the criteria of the protocol, the test article meets the requirements of the Pyrogen Test and is therefore considered non-pyrogenic.

QUALITY ASSURANCE STATEMENT

This study was conducted in compliance with U.S. Food and Drug Administration regulations set forth in 21 CFR, Part 58.

The sections of the regulations not performed by or under the direction of Toxikon Corporation, exempt from this Good Laboratory Practice Statement, included characterization and stability of the test article and its mixture with carriers, 21 CFR, Parts 58.105 and 58.113.

The Quality Assurance Unit conducted inspections on the following dates. The findings were reported to the Study Director and to Toxikon's Management.

INSPECTIONS	DATE OF INSPECTION	DATE REPORTED STUDY DIRECTOR	DATE REPORTED MANAGEMENT
SCORING	07/03/03	07/03/03	07/03/03
RAW DATA	07/16/03	07/16/03	07/16/03
FINAL REPORT	07/16/03	07/16/03	07/16/03

John Christensen, B.S.
Quality Assurance Officer

07/16/03
Date

STUDY DIRECTOR SIGNATURE AND VERIFICATION DATES


This study meets the technical requirements of the protocol. The study also meets with the requirements of the Good Laboratory Practice Regulations, 21 CFR, Part 58, with the exemptions as stated in the Quality Assurance Statement.

Protocol Number: MMS/VIVO/004-03/000

Study Director: Kathleen Krenzer O.D., Ph.D.

Company: Toxikon Corporation

Signature:

Date:

07/16/03

Study Supervisor: William Entwistle, AS, LATG

VERIFICATION DATES:

The Toxikon Protocol Effective Date is the GLP Study Initiation Date. The study dates were as follows:

Protocol Effective Date:	04/22/03
Test Article Receipt:	06/25/03
Project Log Date:	06/25/03
Stock Training/Sham Test:	06/30/03 - 07/02/03
Extraction Dates:	06/30/03 - 07/03/03
Technical Initiation:	07/03/03
Technical Completion:	07/03/03

1.0 PURPOSE

The pyrogen test was designed to determine the presence of chemical pyrogens in extracts of solid materials in order to limit the risks of febrile reaction in a patient to an acceptable level following administration of the product concerned. The test involves measuring the rise in temperature of New Zealand White rabbits following the intravenous injection of a test article. It is designed for products that can be tolerated by the test rabbit in a dose not to exceed 10 mL per Kg, within a period of not more than 10 minutes.

2.0 REFERENCES

The test was based on the following references:

- 2.1 Biological Evaluation of Medical Devices - Part 11: Tests for Systemic Toxicity, ISO 10993-11, 1993.
- 2.2 Biological Evaluation of Medical Devices - Part 12: Sample Preparation and Reference Materials, ISO 10993-12, 2002.
- 2.3 General Requirements for the Competence of Testing and Calibration Laboratories, ISO 17025, 1999.

3.0 COMPLIANCE

The study conformed to all applicable laws and regulations. Specific regulatory requirements included the current FDA 21 CFR, Part 58 - Good Laboratory Practice for Nonclinical Laboratory Studies; AAALAC, International; "Guide for the Care and Use of Laboratory Animals," National Research Council, 1996. (NIH) (OLAW), "Public Health Service Policy on Humane Care and Use of Laboratory Animals," Health Research Extension Act of 1985 (Public Law 99-158 November 20, 1985), Reprinted 1996; USDA, Department of Agriculture, Animal and Plant Health Inspection Service, 9 CFR Ch.1 (1/1/95 edition), Subchapter A-Animal Welfare, ISO 10993-2, 1992, Biological Evaluation of Medical Devices-Part 2: Animal Welfare Requirements.

4.0 IDENTIFICATION OF TEST AND CONTROL ARTICLES

The following information was supplied by the Sponsor wherever applicable; it did not apply to confidential information. The Sponsor was responsible for all test article characterization data as specified in the GLP regulations. Test and control articles (exclusive of extracts) that are mixed with carriers require verification of concentration, homogeneity, and stability. Samples of test and control article extracts were not returned to the Sponsor for characterization and verification.

4.1 Test Article Name: Mini Access Kit

CAS/Code #: MAK401NV

Lot/Batch #: MQ04990

Physical State: Not Supplied by Sponsor (N/S)

Color: N/S

Density: N/S

Stability: N/S

Solubility: N/S

Storage Conditions: Room Temperature

Safety Precautions: Standard Toxikon Laboratory Safety Precautions

4.2 Negative Control Article (Toxikon Supplied)

Control Article Name: 0.9% USP Sodium Chloride for Injection (NaCl)

QC Inventory #: CSC-03-04-006-VIV

Physical State: Liquid

Color: Clear

Storage Conditions: Room Temperature

Safety Precautions: Standard Laboratory Safety Precautions

5.0 IDENTIFICATION OF TEST SYSTEM

5.1 Animals Used in the Test

Number and species: 4 New Zealand White rabbits (*Oryctolagus cuniculus*)

Sex: 2 males and 2 females

Weight/Age range: 2.09 - 2.16 kilograms / at least 10 weeks old (adult)
Weighed to the nearest 10 gm

Health status: healthy, not previously used in other experimental procedures

Animal purchase: Millbrook Breeding Labs, Amherst, MA

Animal selection: selected from a larger pool of animals and examined to ensure lack of adverse clinical signs

Animal identification: ear marker

Acclimatization: minimum 5 days under the same conditions as the actual study

5.2 Animal Care and Maintenance

Housing: individually housed

Cages: suspended stainless steel

Bedding: hardwood chips, P.W.I. Industries, St-Hyacinthe, Quebec, Canada
(non-contact)

Animal room temperature: $68 \pm 5^{\circ}\text{F}$

Animal room relative humidity: 30-70%

Air exchanges per hour: 10 to 15

Lights: 12-hour light/dark cycle, full spectrum fluorescent lights

The laboratory and animal rooms were maintained as limited-access facilities.

Animal rations: TEK 8630 Rabbit Diet, Harlan Teklad, Madison, WI, *ad libitum*

Water: tap water, *ad libitum*

There were no known contaminants present in the feed, water, or bedding expected to interfere with the test data.

6.0 JUSTIFICATION OF TEST SYSTEM AND ROUTE OF ADMINISTRATION

6.1 New Zealand White rabbits were used in this study because they historically have been used in pyrogen evaluation studies. The animal species, number and route are recommended by the ISO 10993-11 guidelines.

6.2 The test article was administered through a solution compatible with the test system. This was the optimal route of administration available in this test system.

7.0 EXPERIMENTAL DESIGN

7.1 Prior to using the rabbits in the pyrogen test, they were conditioned not more than 7 days before use in a Sham Test. This included all the steps as directed in the pyrogen test, except the injection.

7.2 Test and control articles were prepared as described in Section 8.0. All apparatus was prepared according to the requirements of ISO 10993-11.

7.3 The control article was 0.9% USP Sodium Chloride for Injection (NaCl), treated in the same manner as the test article.

7.4 Single-use disposable syringes and hypodermic needles were used for the administration of the test article. Glassware used in testing was rendered free from pyrogens by heating at 200°C for not less than 3 hours.

7.5 A Digital Thermometer with thermistor probes (accuracy of 0.1°C) was used to measure rectal temperatures. A probe was inserted into the rectum of the test rabbit to a depth of not less than 7.5 cm. After a period of temperature stabilization (approximately 1 minute), the body temperature was taken.

7.6 The test was performed in the Pyrogen Test Room, designated solely for the pyrogen testing, under environmental conditions similar to those under which rabbits are housed and free from disturbances, which may excite them. The rabbits were restrained for the course of the study and were not offered food or water during the three hour test period. Since rectal temperature measuring probes were inserted throughout the testing period, the rabbits were restrained with light-fitting neck stocks that allowed the animals to assume a natural resting position.

7.7 Not more than 30 minutes prior to injection of the test article, the baseline temperatures were determined. This served as a baseline for determining any increase in temperature resulting from the injection of the test article. None of the test animals exhibited a baseline temperature which varied by more than 1°C from each other or which was greater than 39.8°C.

7.8 Body temperatures were recorded at 0 hour and then at 30 minute intervals between 1 and 3 hours subsequent to injection.

7.9 At the end of the observation period, the animals were returned to the general colony.

8.0 DOSAGE

8.1 Preparation of Test Article

The test and control articles were prepared in NaCl at ratios specified by ISO 10993-12.

The samples were sterilized by the Sponsor (3 x EO). Per Sponsor's specifications, prior to extraction, the needle cover and wire hoop were removed and discarded. The test article was extracted at a ratio of 0.2 gm per 1.0 mL of 0.9% USP Sodium Chloride for Injection (NaCl) at 50±2°C for 72 hours.

Properly prepared test article was placed in an extraction bottle and NaCl was added. The extraction medium completely covered the test article.

The extracting medium (control article) was prepared for parallel treatments and comparisons. The control article was prepared in the same manner as the test article.

Each extract was agitated vigorously prior to administration.

8.2 The test or control article was injected into the marginal ear vein of each rabbit at a dose based on the body weight of each animal. Test article extracts were injected at a dose of 10 mL per Kg. Each injection was completed within 10 minutes after the start of administration.

9.0 EVALUATION CRITERIA

9.1 Temperature decreases are considered as zero rise.

9.2 If no rabbit shows an individual rise in temperature of 0.5°C or more above the baseline temperature, the test article meets the requirements for the absence of pyrogens.

9.3 If any rabbit shows an individual temperature rise of 0.5°C or more, a retest should be conducted using five other rabbits. If not more than three of the eight rabbits show individual rises in temperature of 0.5°C or more, and if the sum of the eight individual maximum temperature rises does not exceed 3.3°C, the test article meets the ISO requirements for the absence of pyrogens.

10.0 RESULTS

The temperature increases for the treated animals were 0.1, 0.0, and 0.0°C. The increases did not exceed the test limit for the maximum individual temperature rise. The temperature increase of the control animal was 0.0°C.

11.0 CONCLUSION

The 0.9% USP Sodium Chloride for Injection (NaCl) extract of the test article, Mini Access Kit, was evaluated for its potential to produce a pyrogenic response when tested in New Zealand White rabbits. Based upon the criteria of the protocol, the test article meets the requirements of the Pyrogen Test and is therefore considered non-pyrogenic.

12.0 RECORDS

- 12.1 Original raw data is archived at Toxikon Corporation.
- 12.2 A copy of the final report and any report amendments are archived at Toxikon Corporation.
- 12.3 The original final report and a copy of any protocol amendments or deviations were forwarded to the Sponsor.
- 12.4 All unused test article shall be disposed of by Toxikon, per Sponsor's specifications.
- 12.5 Final reports shall not be reproduced, except in full, without the written authorization/approval from Toxikon.

13.0 CONFIDENTIALITY AGREEMENT

Statements of confidentiality were as agreed upon prior to study initiation.

14.0 POLICY ON PAIN AND SUFFERING IN ANIMALS

No evidence of pain and suffering was observed in the animals.

15.0 ANIMAL USAGE

The Sponsor assures that, to the best of their knowledge, this study did not unnecessarily duplicate previous testing.

TABLE I
PYROGEN TEST DATA

Test Article: Mini Access Kit

Lot/Batch No: MQ04990

Test Parameters	Animals			
	30651 (Treated)	30652 (Treated)	30653 (Treated)	30654 (Control)
Weight (Kg)	2.11	2.09	2.14	2.16
Dose (mL)	21.1	20.9	21.4	21.6
Probe #	9	6	18	4
Baseline Temp (°C)	39.2	39.6	39.2	39.4
1.0 Hour Temp	39.3	39.3	38.9	39.3
1.5 Hour Temp	39.2	39.1	38.9	39.4
2.0 Hour Temp	39.1	39.3	39.0	39.4
2.5 Hour Temp	38.9	39.3	39.2	39.2
3.0 Hour Temp	38.9	39.5	39.0	39.1
Temp Increase*	0.1	0.0	0.0	0.0

* Temperature increases based upon maximum temperatures recorded at 30 minute intervals between 1 and 3 hour observations.

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
78
Technical Documentation
MAK-NV
Dated July 27, 2020
Consisting of 31 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 27th day of July 2020, before me, Jennifer Fordham, a notary public, personally appeared Cory Marsh, Manager, Regulatory Operations proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.


Notary Public

My Commission Expires: June 27, 2021

Cory Marsh
Manager, Regulatory Operations
Merit Medical Systems
27 July 2020

**Extract from technical documentation for
MEDICAL DEVICE**

«The non-vascular access kit MAK-NV™»

Manufactured by

Merit Medical Systems, Inc., the USA

Manufactured by:

Merit Medical Systems, Inc.

1600 West Merit Parkway South Jordan, Utah, 84095, USA

The place of production of the medical device

1) Merit Medical Systems, Inc., 1600 West Merit Parkway South Jordan, Utah, 84095, USA

2) Merit Maquiladora Mexico, S. DE R.L. DE C.V., Avenida Sor Juana Ines de la Cruz
19970 Interior B, Edificio 2, Parque Industrial Frontera Tijuana, Baja California, MEXICO
22630

The authorized representative of the manufacturer on the territory of the Russian Federation:

"Merit technologies"

Of the Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

phone: +7 (495) 221-89-02

1. DEVICE DESCRIPTION

The non-vascular access kit MAK-NV™ consists of a coaxial introducer sheath (a pair of coupled introducer sheath / dilator), a introducer needle with a stylet and a mini guidewire with a diameter of 0,018" (0,46 mm) and / or a diagnostic guidewire with a diameter of 0,038" (0,97 mm). The non-vascular access kit MAK-NV™ is available with an introducer needle with a trocar or Chiba stylet, a mini guidewire, and / or a diagnostic guidewire. The composition of the medical device is presented below.

The non-vascular access kit MAK-NV™:

1. The non-vascular access kit MAK-NV™, with a guide needle with a trocar stylet, diameter: 21 G (0.8 mm), consisting of:

- Introducer needle with a trocar stylet, diameter: 21 G (0,8 mm), length: 15 cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1,33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018 " (0,46 mm), length 60 cm - 1 pc.(if necessary);
- guidewire, diameter 0,038 " (0,97 mm), length: 80 cm, 150 cm - 1 pc. (if necessary);
- Instructions for use.

2. The non-vascular access kit MAK-NV™, with a guide needle with a Chiba stylet, diameter: 21 G (0.8 mm), consisting of:

- Introducer needle with Chiba stylet, diameter: 21 G (0.8 mm), length: 15cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1,33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018" (0,46 mm), length 60 cm - 1 pc.(if necessary);
- guidewire, diameter 0,038" (0,97 mm), length: 80 cm, 150 cm - 1 pc. (if necessary);
- Instructions for use.

The non-vascular access kit MAK-NV™ used to provide percutaneous introduction of a non-vascular Seldinger catheter into a body cavity or lumen during interventional radiology procedures, including drainage.

The coaxial introducer is equipped with a Luer type connector with internal thread, which allows the system to function as a single product. A conductor needle with a stylet provides support for a pair of introducer / dilator during installation. The distal end of the introducer sheath is equipped with a radiopaque marker, which facilitates the visualization of the tip during fluoroscopy.


Figure 1. The non-vascular access kit MAK-NV™ (with Chiba type stylet)


Figure 2. The non-vascular access kit MAK-NV™ (guidewire needle with trocar stylet)

2. Intended use / Indications for use, contraindications.

INDICATIONS FOR USE

The non-vascular access kit MAK-NV™ is intended for percutaneous placement of a 0,035" (0,89 mm) or 0,038" (0,97 mm) guidewire with the introduction of a non-vascular catheter by the Seldinger method into a cavity or lumen of the body.

CONTRAINDICATIONS AND WARNINGS:

There are no contraindications.

Px Only: DEVICE SALE RESTRICTION

For professional use only.

Attention! Federal law (USA) authorizes the sale of this device only to a doctor or by order of a doctor who has been trained to use this device and / or has experience using it.

PRECAUTIONARY MEASURES

- The kit should only be used for non-vascular procedures.
- Read the manufacturer's instructions before use.
- The contents of unopened and undamaged packaging are sterilized using ethylene oxide and pyrogen-free.
- Do not use if packaging is open, damaged, or broken.
- Use for only one patient. Reuse, handling or sterilization is prohibited. Reuse, handling or sterilization may impair the integrity of the device and / or cause the device to malfunction and, as a result, result in injury, illness or death of the patient. Reuse, handling or sterilization can also create a risk of contamination of the device and / or lead to infection of the patient or cross-infection, including but not limited to the transfer of infectious diseases from patient to patient. Contamination of the device can be harmful to the patient, resulting in illness or death.
- Observe generally accepted precautions when introducing and maintaining this device. Since there is a risk of contracting bloodborne pathogens, healthcare providers should always follow standard precautions when performing procedures with blood and body fluids. Sterility should always be observed.
- Do not use after the expiration date.
- Store in a cool, dry place.

WARNING

Removing, pulling, or manipulating the distal tip of the conductor through the tip of the needle can cause breakage or embolization.

The guidewire should not be advanced if resistance is felt.

ADVERSE REACTIONS

The non-vascular access kit MAK-NV™ should be used by risk-informed individuals who are qualified to perform these procedures. Possible complications include risks commonly associated with percutaneous diagnostic and / or interventional procedures. Other complications may include, but are not limited to:

- edema
- blood clot formation
- infection
- pseudo-aneurysm formation

- hematoma
- perforation or rupture of a vessel or internal organ
- bleeding
- embolization of the conductor or catheter
- hemothorax
- inflammation, necrosis, or scarring
- hydrothorax
- pain in the area of application
- extravasation

3. CLASSIFICATION:

The non-vascular access kit MAK-NV™ is designated as class IIa products in accordance with Annex II of the Medical Device Directive 93/42 / EEC on medical devices.

4. Sterilization Method

Merit Medical Systems, Inc., utilizes Ethylene Oxide (EO) sterilization for the products described in this Technical File. The EO sterilization cycle is validated according to European International Standard ISO 11135. Annual re-validation occurs assessing a worst-case reference load based on density and boxing configuration, which allows routine sterilization of a mixed-pallet load. Merit's products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals, no more than 10^{-6} .


Protocols and reports, validated sterilization process, are stored in the control system of documentation of "Merit Medical Systems, Inc.," And indicated in the summary table for the sterilization of "Merit Medical Systems, Inc."

This product does not provide a user re-sterilization.

5. PRODUCT SPECIFICATIONS

5.1 Guiding needle with stylet, diameter: 21 G (0.8 mm), length: 15 cm, 20 cm

Outside diameter of the needle	21 G (0.800mm+ 0.020mm)
Inner diameter of the needle	0.580 mm $\pm$ 0.020 mm
Tube wall type	thin-walled
Needle color coding	Dark green
Needle tube length	15.0 cm $\pm$ 0.5 cm 20.0 cm $\pm$ 0.5 cm
Needle tip sharpening angle	$32^{\circ} \pm 2^{\circ}$
Effective stylet length	170 mm $\pm$ 5 mm (with 150 mm (15 cm) needle tube length) 220 mm $\pm$ 5 mm (with a needle tube length of 200 mm (20 cm)) The stylet should be longer than the needle when the two components are assembled into one structure.
	0.500mm $\pm$ 0.020mm

Stiletto tip shape	stiletto trocar stiletto chiba type
Sharpening angle of the stylet trocar type	$28^{\circ} \pm 2^{\circ}$
Stylet apex angle trocar type	$35^{\circ} \pm 2^{\circ}$
Stylet shape, trocar type	pyramidal
Sharpening angle of the stylet type Chiba	$32^{\circ} \pm 2^{\circ}$
Tensile strength of the stylet-sleeve structure (fastening strength)	140 N (31.5 lb) minimum
product Weight	$2,45 \pm 0,5$ g
Tensile strength at break structure "needle-hub" (mounting strength)	No less than 350 N (78,7 lb)
Surface	When carrying out a visual check on the external surface of the working length of the device no foreign substance or surface defects. On needlepoint no protrusions growths, hooks etc.
The protruding portion of the needle tip	In the assembled structure cannula does not extend beyond the tip of the stylet, or its length does not exceed the length of the stylet more than 2 mm.
Allowable angle deviations needle without compromising the integrity of	$20^{\circ} \pm 1^{\circ}$
Explorer compatibility	Gidewire diameter of 0.018 inches (0.46 mm)
The hardness of the needle Rockwell	52-54 HRC
The radius of the working portion of blunting needle	not more than 0,03 mm
needle roughness	no more than Ra 0,4 microns
Radiopacity	The distal end of the Introducer needle is equipped with a radiopaque marker, $10\text{mm} \pm 1\text{mm}$ in length, which facilitates the visualization of the tip during fluoroscopy.
The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)	 Length of female conical element (E): 7,500 - 10,500 mm

	Luer Taper: 6% Profile angle of the outer protrusion on the load-bearing surface providing resistance to separation (degrees), Σ: 25 ° -30 °. Smaller outer rib diameter ($\emptyset J$): 5,515-6,730 mm Larger outer diameter ($\emptyset H$): 7.730-7.830 mm Radius of rounding or chamfer at the outer tip of the socket conical element: no more than 0.500 mm Inner diameter of the open end of the female cone at a distance of 0.75 mm from the tip of the female cone (according to EN ISO 80369-7: 2016): 4.198-4.298 mm Inner diameter of the narrow end of the female cone at a distance of 7.5 mm from the tip of the female cone, (according to EN ISO 80369-7: 2016): 3.793-3.893 mm
--	---


5.2 Coaxial introducer, diameter 6F

Outside diameter	Taper length	Effective length
6F $\pm$ 1F (2.000mm $\pm$ 0.300mm)	6 mm $\pm$ 1 mm	260 mm $\pm$ 10 mm

- dilator, 4F diameter,

Outside diameter	Taper length	Effective length
4F $\pm$ 1F (1.330mm $\pm$ 0.300mm)	6 mm $\pm$ 1 mm	250 mm $\pm$ 10 mm
Tensile strength	15 N (3.375 lb) minimum	
The inner diameter of the distal (narrowed) end of the dilator	For diameter 4F	0.041 " $\pm$ 0.002 " (1.04mm $\pm$ 0.05mm)
Guidewire compatibility 0.038 in. (0.97 mm)	The dilator should have a 0.038 inch (0.97 mm) guidewire	
Geometric shape of the distal end of the dilator	The presence of a tapered taper, ground the edges. The edge should not be sharp. No cracks scratches.	
Weight	1.60 g $\pm$ 0.5 g	
Radiopacity	should be visible on fluoroscopy	

The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)

Length of female conical element (E): 7,500 -10,500 mm

Luer Taper: 6%

Profile angle of the outer protrusion on the load-bearing surface providing resistance to separation (degrees), Σ : 25 ° -30 °.

Smaller outer rib diameter (ØJ): 5,515-6,730 mm


Larger outer diameter (ØH): 7.730-7.830 mm

Radius of rounding or chamfer at the outer tip of the socket conical element: no more than 0.500 mm

Inner diameter of the open end of the female cone at a distance of 0.75 mm from the tip of the female cone (according to EN ISO 80369-7: 2016): 4.198-4.298 mm

Inner diameter of the narrow end of the female cone at a distance of 7.5 mm from the tip of the female cone, (according to EN ISO 80369-7: 2016): 3.793-3.893 mm

The main characteristics of the connector (connector) of the locking (plug) type Luer-Lock (according to EN ISO 80369-7: 2016)


Diameter of the smallest cylinder covering the outer surfaces of the outer elements of the coupling (Øw): 8,800-11,500 mm

Maximum internal thread diameter (root diameter) (Øh): 8,000 ± 0.100 mm

	Minimum internal thread diameter (thread crest diameter) ($\emptyset j$): $7,000 \pm 0.200$ mm Luer Lock Taper: 6% Luer-Lock thread pitch (p): 2,500 mm Part of connector tip extending beyond the threaded sleeve (c): 2,100-2,573 mm Inner diameter of the tip of the conical plug ($\emptyset f$): 2,100-2,900 mm Length of plug cone (e): 7,500-10,500 mm Diameter of the narrow base of the outer cone of the plug cone: 3.970 - 4.072 mm Diameter of the wide base of the outer cone of the plug-in cone: 4,376-4,476 mm Radius of rounding or chamfer at the outer tip of the conical plug element: max. 0.500 mm
--	--

- rigid metal cannula,

Characteristic	
The outer diameter of the cannula	$0.810 \text{ mm} +0.0127 \text{ mm} / \text{mm} -0.0254$
The inner diameter of the cannula	$0.560 \text{ mm}, + 0.0127 \text{ mm} / \text{mm} -0.0254$
The length of the exposed portion (1 mm + / -2 mm)	238
The strength of attachment (lbs)	At least 10 h (2.25 lb)
Fracture strength	At least 10 h (2.25 lb)
Rockwell hardness	53-55 HRC
roughness Ra	$\leq 0,4 \text{ m}$
The allowable deviation angle, without compromising the integrity of Surface	$15^0 \pm 1^0$
Weight	$1.20 \text{ g} \pm 0.5 \text{ g}$
The main characteristics of the connector (connector) female lock type Luer (according to EN ISO 80369-7: 2016)	 Length of female conical element (E): 7,500 -10,500 mm Luer Taper: 6% Profile angle of the outer protrusion on the load-bearing

	surface providing resistance to separation (degrees), Σ: 25 ° -30 °. Smaller outer rib diameter ($\emptyset J$): 5,515-6,730 mm Larger outer diameter ($\emptyset H$): 7.730-7.830 mm Radius of rounding or chamfer at the outer tip of the socket conical element: no more than 0.500 mm Inner diameter of the open end of the female cone at a distance of 0.75 mm from the tip of the female cone (according to EN ISO 80369-7: 2016): 4.198-4.298 mm Inner diameter of the narrow end of the female cone at a distance of 7.5 mm from the tip of the female cone, (according to EN ISO 80369-7: 2016): 3.793-3.893 mm
The main characteristics of the connector (connector)) of the locking (plug) type Luer-Lock (according to EN ISO 80369-7: 2016)	Diameter of the smallest cylinder covering the outer surfaces of the outer elements of the coupling ($\emptyset w$): 8,800-11,500 mm Maximum internal thread diameter (root diameter) ($\emptyset h$): 8,000 $\pm$ 0.100 mm Minimum internal thread diameter (thread crest diameter) ($\emptyset j$): 7,000 $\pm$ 0.200 mm Luer Lock Taper: 6% Luer-Lock thread pitch (p): 2,500 mm Part of connector tip extending beyond the threaded sleeve (c): 2,100-2,573 mm Inner diameter of the tip of the conical plug ($\emptyset f$): 2,100-2,900 mm Length of plug cone (e): 7,500-10,500 mm Diameter of the narrow base of the outer cone of the plug cone: 3.970 - 4.072 mm Diameter of the wide base of the outer cone of the plug-in cone: 4,376-4,476 mm Radius of rounding or chamfer at the outer tip of the conical plug element: max. 0.500 mm

5.3 guidewire diameter of 0.018 "(0.46 mm), 60 cm long

The diameter of the guidewire	0.018 " (+0005"/ -0.0010") (0.460 mm +0.0127 mm / mm -0.0254)
The length of the guidewire	600 mm $\pm$ 25 mm
The diameter of the guard guidewire	4 mm $\pm$ 1 mm
The radius of the distal tip	3 mm $\pm$ 1 mm
Lubricity rod coating	Coefficient of friction should be not more than 0.386
Durability distal tip	Not less than 22.25 N (5.0 lbs)
Durability proximal tip	Not less than 12.23 N (2.75 lbs)
The allowable deviation angle, without compromising the integrity of	90° $\pm$ 1°
The weight of the guidewire in a protective casing	3.20 $\pm$ 0.5 g

5.4 guidewire diameter of 0.038 "(0.97 mm) and a length of 150 cm, c J-shaped tip

The diameter of the guidewire	0.038" (+0005"/ -0.0010") (0.970 mm, +0.0127 mm / mm -0.0254)
The length of the guidewire	1500 mm $\pm$ 125 mm 800 mm $\pm$ 55 mm
The diameter of the guard guidewire	7 mm $\pm$ 1 mm
The radius of the distal tip	3 mm $\pm$ 1 mm
Lubricity rod coating	Coefficient of friction should be not more than 0.386
Durability distal tip	Not less than 22.25 N (5.0 lbs)
Durability proximal tip	Not less than 12.23 N (2.75 lbs)
The allowable deviation angle, without compromising the integrity of	90° $\pm$ 1°
The weight of the guidewire in a protective casing, length 80 cm	13,0 $\pm$ 0.5 g
The weight of the guidewire in a protective casing, length 150 cm	16,5 $\pm$ 0,5 g

5.5 Requirements for the packaging of the terminal sterilized product (main characteristics).

Packing width	245 ± 5 mm
Length packaging	495 ± 5 mm
Adhesive width	10 ± 2 mm
Tensile strength of the adhesive layer	≥ 2 N
Packing material density	0.92 ± 0.05 g / cm ³
Longitudinal tensile strength	> 0.7 kN / cm ²
Tensile strength in the transverse direction	> 0.9 kN / cm ²
Bursting strength	> 650 kPa

Certification of packaging for products of "Merit Medical" conducted in accordance with the requirements of the following standards:

ISO 11607 "Packaging of medical devices to be terminally sterilized"

ASTM 4169 "Test performance packaging qualities"

ISO 2233 "Packaging - Complete, filled transport packages and unit loads - Conditioning for testing"

6. MATERIALS OF CONSTRUCTION

Component	Description of materials
Introducer needle with a trocar stylet, diameter: 21 G	
Luer-type connector	polycarbonate CALIBRE™ Megarad™ 2081 manufacturer Trinseo LLC, USA
Needle	Stainless steel AISI 304 (composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) manufacturer of Ulbrich Stainless Steels & Special Metals, Inc, USA
Stylet trocar / stylet type Chiba	Stainless steel alloys AISI 304 (composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 16.0 / 18.0%, Ni - 10.0 / 14.0%, Mo - 2.0-3.0%, Fe - 68.92-61.92%) manufacturer of Ulbrich Stainless Steels & Special Metals, Inc, USA
Sleeve	polycarbonate CALIBRE™ Megarad™ 2081 manufacturer Trinseo LLC, USA
	Dye - pellets CC10130299 TRANS GREEN HALF PC PolyOne™ Color: Green, classification code Pantone 354U manufacturer POLYONE CORPORATION, USA
coaxial introducer	
Luer connector	Polyester polycarbonate blend Xenoy HX5600HP,

Component	Description of materials
	manufacturer of SABIC Innovative Plastics BV Pigment - titanium oxide solid fraction Xenoy EXN X5600HP Color: bright white, code classification Pantone 11-0601 TCX manufacturer of SABIC Innovative Plastics BV Green ink INK-PG Pad Printing TYPE UV 2002 IWT Trans Tech, manufacturer of Trans Tech America, USA
nut	The thermoplastic ABS resin (ABS) Lustran® 248 NR, manufacturer Dow Plastic, USA Dye - pellets 10200443 TRANS GREEN PolyOne™, manufacturer PolyOne technologies, US Color: Green, Pantone 354U Classification code (for diameter 6,5F)
dilator tube	polyethylene copolymer Alathon H5618, manufacturer Trinseo LLC, USA. Pigment - titanium oxide solid fraction A32329C PolyOne™, manufacturer PolyOne technologies, US Color: bright white, code classification Pantone 11-0601 TCX Barium sulfate - 12%, the manufacturer of Sachtleben Chemie GmbH
The rigid cannula	Stainless steel alloys AISI 304 (composition: C -0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 16.0 / 18.0%, Ni - 10.0 / 14.0%, Mo - 2.0-3.0%, Fe - 68.92-61.92%), manufacturer of Ulbrich Stainless Steels & Special Metals, Inc, USA
guidewire diameter diameter 0.018 "(0.46 mm) diameter and 0.038" (0.97 mm)	
	Stainless steel alloys AISI 304 (composition: C -0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 16.0 / 18.0%, Ni - 10.0 / 14.0%, Mo - 2.0-3.0%, Fe - 68.92-61.92%), manufacturer of Ulbrich Stainless Steels & Special Metals, Inc, USA Polycarbonate CALIBRE™ Megarad™ 2081 manufacturer Trinseo LLC, USA Polypropylene Bormed HD810MO, manufacturer BOREALIS AG
guidewire diameter diameter 0.018 "(0.46 mm)	
	Nitinol- platinum Inconel X-750 manufacturerDow Corning Corporation Polycarbonate CALIBRE™ Megarad™ 2081 manufacturer Trinseo LLC, USA Polypropylene Bormed HD810MO, manufacturer

Component	Description of materials
	BOREALIS AG

** The manufacture of this medical product is used a chemical element barium sulfate. Barium sulfate serves as a component in the manufacture of these products, in this case it is not a dosage form or substance. Barium sulfate - one of the main X-ray contrast materials, which are added to the polymeric medical devices. Barium sulfate has no pharmacologic and metabolic effects on the body. Chemical formula - BaSO₄. Barium sulphate due to the low solubility is not toxic to the body substance and may therefore its use as a radiopaque agent. In applying these products to the destination BaSO₄ entering the human body is excluded. Barium sulfate is used as a part of MI MI for visualization during surgery.*

Kits are packaged and labeled as sterile single use products. sterility protection system defined in accordance with standard ISO 11607 for preventing penetration of microorganisms and provide an aseptic product presentation in the place of use. As part of the production process is carried out verification of the integrity of the packaging.

packaging component	Material
Package	High-density polyethylene (NDPE) DuPont™ Tyvek® 1073B, manufacturer «DuPont», USA
Cardboard box	Corrugated
shipping container	Corrugated

7. Compliance with Relevant Regulations and Technical Standards

The product is in compliance with the Council Directive 93/42/EEC, Medical Device Directive of the European Union. The applicable CE certificate is CE 541900.

The product is legally marketed according to the requirements of the United States Food and Drug Administration in compliance with the premarket notification 510(k) process. The 510(k) clearance letter is K141408.

The qualification activities that were performed for the product were based on standards that were in effect at the time of the testing. Standards applied include:

Standard	header
Common standards	
EU Regulation 2017/745	EU regulations on medical devices (Regulation of the European Parliament and of the Council on medical products 2017/745 on April 5, 2017)
ISO 14971	Medical devices. Application of risk management to medical devices.
Council Directive 93/42 / EEC	93/42 / EEC Council directive on medical device European Union
USP	U.S.P
ISO 13485	Medical devices - Quality management systems - Requirements for regulatory purposes
ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and research as part of the risk management process.
ISO 11135	Sterilization of health care products. Ethylene oxide. Requirements for the

	development, validation and routine control of a sterilization process of medical devices.
ISO 11607-1	Packaging of medical devices to be sterilized at the final stage of production - Part 1: Requirements for materials, sterile protection systems and packaging systems.
ISO 11607-2	Packaging medical devices to be sterilized in the final stages of production - Part 2: Requirements for validation processes forming, sealing and assembly.
ASTM D4169	Standard methods of testing the efficiency of transport of containers and systems.
ASTM F1980	Standard Guide for accelerated aging protect the sterility of medical devices systems.
BS EN 1041	Information supplied by the manufacturer of medical devices.
ISO 15223-1	Medical devices. Symbols used on labels for labeling of medical devices and in the information provided - Part 1: General requirements.
BS EN 15986	Symbols used in the labeling of medical devices. Labeling Requirements medical devices containing phthalates.
EN 556	Sterilization of medical devices. Requirements for medical devices subject to final sterilization category "sterile."
Specific standards relevant to the product	
EN 1617	Sterile drainage catheters and accessories for single use them.
ISO 594-1	Details of the connection with a taper of 6% (Luer) taper for syringes, needles and other medical equipment. Part 1: General requirements
ISO 594-2	Details of connecting with the taper 6% (Luer) taper for syringes, needles and certain other medical equipment - Part 2: latching connection.
EN ISO 80369-7:2017	Small bore connectors for liquids and gases in medical applications

8. Risk Analysis System

Potential risks are evaluated during the design / product development and manufacturing processes according to ISO 14971. The appropriate risk assessment tool (FMEA, Clinical Risk Assessment - CRA.) is applied in order to identify potential failure modes and potential effect on patients or users. Appropriate control and prevention mechanisms are defined and implemented. Risk Analysis reports are revisited and updated as necessary, and controlled under Merit's Document Control system. A cross-functional team, including persons with clinical experience, reviewed the Risk Analysis reports to ensure that all potential risks were addressed.

9. Biocompatible

The product is certified as meeting the requirements of ISO 10993-1 and the leadership of the US Office of Quality Control of Food and Drug Administration "Required training for biocompatibility and toxicological profiles for the evaluation of medical devices" (US Department of quality control of food and medicines, Memorandum Blue G95-1 book).

<u>PRODUCT</u>	<u>CLASSIFICATION</u>
----------------	-----------------------

PRODUCT	CLASSIFICATION
The non-vascular access kit MAK-NV™	An article externally attachable and contacting with the tissue for a limited time period (≤ 24 hours)

A brief report on the results of tests on biocompatibility

Tests Performed to Evaluate the Biological Safety of the MAK-NV per Patient Contacting Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Polyethylene Resin (PN 000178001), Marker Tip Resin (PN 000251001)* - All samples were prepared in accordance with ISO 10993-12. Sterilization information is unavailable.		
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution	Non-Cytotoxic	Toxikon GLP Report 07-4673-G3
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	Toxikon GLP Report 07-4673-G2
Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Toxikon GLP Report 07-4673-G5
Systemic Toxicity (Acute)/ ISO 10993-11/ Acute Systemic Injection	Negative	Toxikon GLP Report 07-4673-G4
Systemic Toxicity/ ISO 10993-11/ Material Mediated Rabbit Pyrogenicity	Non-Pyrogenic	Toxikon GLP Report 07-4673-G7
Tests Performed to Evaluate the Biological Safety of the MAK-NV per Patient Contacting Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Implantation/ ISO 10993-6/ 2 Week Implant Study	Negative	Toxikon GLP Report 07-4673-G8
Implantation/ ISO 10993-6/ 4 Week Implant Study	Negative	Toxikon GLP Report 07-4673-G9
Subchronic Toxicity/ ISO10993-11/ 14 Day Toxicity Study	Negative	Toxikon GLP Report 07-4673-G1
Genotoxicity/ ISO 10993-3/ Salmonella Typhimurium and Escherichia Coli Reverse Mutation Assay	Non-Genotoxic	Toxikon GLP Report 07-4673-G6
Hemocompatibility/ ISO 10993-4/ Indirect Contact Hemolysis	Non-Hemolytic	Toxikon GLP Report 07-4673-G14
Hemocompatibility/ ISO 10993-4/ Direct Contact Hemolysis	Non-Hemolytic	Toxikon GLP Report 07-4673-G15

Hemocompatibility/ ISO 10993-4/ In Vitro Hemocompatibility	Meets Acceptance Criteria	Toxikon GLP Report 07-4673-G12
Hemocompatibility/ ISO 10993-4/ Complement Activation	Meets Acceptance Criteria	Toxikon GLP Report 07-4673-G13
Hemocompatibility/ ISO10993-4/ Prothrombin Time	Meets Acceptance Criteria	Toxikon GLP Report 07-4673-G16
Hemocompatibility/ ISO10993-4/ Unactivated Partial Thromboplastin Time	Meets Acceptance Criteria	Toxikon GLP Report 07-4673-G17
Hemocompatibility/ ISO10993-4/ Lee-White Clotting Time	Meets Acceptance Criteria	Toxikon GLP Report 07-4673-G11
Chemical Characterization/ USP Physicochemical Tests for Plastics	Pass	Toxikon GLP Report 07-4673-G18
Tests Performed to Evaluate the Biological Safety of the MAK-NV per Patient Contacting Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
White Color Additive (PN 000020), Polyethylene Tubing (PN 702185002), Polyethylene Tubing (PN 702183002) - All samples were prepared in accordance with ISO 10993-12. Sterilization information is unavailable.		
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution	Non-Cytotoxic	Toxikon GLP Report 03-3171-G1
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	Toxikon GLP Report 03-3171-G2
Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Toxikon GLP Report 03-3171-G3
Systemic Toxicity (Acute)/ ISO 10993-11/ Acute Systemic Injection	Negative	Toxikon GLP Report 03-3171-G4
Systemic Toxicity/ ISO 10993-11/ Material Mediated Rabbit Pyrogenicity	Non-Pyrogenic	Toxikon GLP Report 03-3171-G6
Genotoxicity/ ISO 10993-3/ Salmonella Typhimurium and Escherichia Coli Reverse Mutation Assay	Non-Genotoxic	Nelson Laboratories Report 239315
Hemocompatibility/ ISO 10993-4/ Direct Contact Hemolysis	Non-Hemolytic	Toxikon GLP Report 03-3171-G5
Stainless Steel (PN 790030003) - All samples were prepared in accordance with ISO 10993-12 and were sterilized by ethylene oxide.		
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution	Non-Cytotoxic	Nelson Laboratories Report 490191

Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	Toxikon GLP Report 10-2857-G1
Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Toxikon GLP Report 10-2857-G2
Tests Performed to Evaluate the Biological Safety of the MAK-NV per Patient Contacting Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Systemic Toxicity (Acute)/ ISO 10993-11/ Acute Systemic Injection	Negative	Toxikon GLP Report 10-2857-G3
Systemic Toxicity/ ISO 10993-11/ Material Mediated Rabbit Pyrogenicity	Non-Pyrogenic	Toxikon GLP Report 03-3171-G6
Hemocompatibility/ ISO 10993-4/ Direct Contact Hemolysis	Non-Hemolytic	Toxikon GLP Report 10-2857-G5
Hemocompatibility/ ISO 10993-4/ In Vitro Hemocompatibility	Meets Acceptance Criteria	Toxikon GLP Report 10-2857-G6
Hemocompatibility/ ISO 10993-4/ Complement Activation	Meets Acceptance Criteria	Toxikon GLP Report 10-2857-G9
Hemocompatibility/ ISO 10993-4/ Prothrombin Time	Meets Acceptance Criteria	Toxikon GLP Report 10-2857-G7
Hemocompatibility/ ISO 10993-4/ Unactivated Partial Thromboplastin Time	Meets Acceptance Criteria	Toxikon GLP Report 10-2857-G8
Polyethylene Resin (PN 000186001) - All samples were prepared in accordance with ISO 10993-12. Sterilization information is unavailable.		
Cytotoxicity/ ISO 10993-5/ L929 MEM Elution	Non-Cytotoxic	Toxikon GLP Report 06-4959-G1
Sensitization/ ISO 10993-10/ Kligman Maximization Test	Non-Sensitizer	Toxikon GLP Report 06-4959-G2
Intracutaneous Reactivity/ ISO 10993-10/ Intracutaneous Injection Test	Negligible Irritant	Toxikon GLP Report 06-4959-G4
Systemic Toxicity (Acute)/ ISO 10993-11/ Acute Systemic Injection	Negative	Toxikon GLP Report 06-4959-G3

Tests Performed to Evaluate the Biological Safety of the MAK-NV per Patient Contacting Component Listed in the Materials Table. Acceptance criteria were met on all tests.		
Test Method	Results	Test Report
Hemocompatibility/ ISO 10993-4/ Indirect Contact Hemolysis	Non-Hemolytic	Toxikon GLP Report 06-4959-G5

Chemical Characterization/ USP Physicochemical Tests for Plastics	Pass	Nelson Laboratories Report 351460
*PN 000251001 was determined to be equivalent to PN 000232001 used in Merit's Prelude Sheath Introducers. The testing presented for this resin was conducted on the Prelude marker tip tubes. In November of 2011 a material change was made to the polypropylene material used in this resin. Subsequently the part number for the resin was changed to 000251002. Full biocompatibility testing was performed on this new polypropylene resin with no adverse responses identified. As the polypropylene material was only a subcomponent of the resin, the testing information presented for 000251001 applies to 000251002.		

Performance and functionality of the product have been tested and certified in the process of designing and developing products in accordance with the approved protocols and reports. The product meets the requirements specified in the design process. Requirements and acceptance criteria performance test articles were based on standard EN 1617 "Sterile catheters and drainage aids for single use", the respective items of auxiliary products for drainage catheters. Completed tests confirm that the products are suitable for their intended use. Below is a brief overview of the tests for the verification / validation of the design of medical products.

Verification / validation of the design

LINK	TITLE
The non-vascular access kit MAK-NV™	
08-0121	Preliminary and final examination tilting operation of the marker and installation certification for the introducer 6 F MAK-NV™
70091-2008-003	Verification report of the casting process of the plug hub for MAK-NV dilator™, Introducer cannula and rigid, using the equipment for injection molding Van Dorn
08-0123	operating certification protocol MAK-NV introducer™ after tilting
70091-2008-015	AE protocol for the study of accelerated aging introducers assembly sets for non-vascular mini-access (MAK-NV™)
07-0227	Certification nut dilator IAC art. № 101251
07-0262	Compliance HDPE molded luer cannulas Alathon H5618 standard ISO 594-2
07-0263	Compliance HDPE company "Dow" molded luer cannulas 10462N Standard ISO 594-2
70091-2008-002	Check the insertion molding process for dilator hub MAK-NV™ using equipment for injection molding Van Dorn
09-0096	AE tilting introducer 6 F MAK-NV™
09-0109	Validation of the new tube of tungsten carbide for finishing tilting introducer 6 F MAK-NV™
09-0225	Certification operation (AE) tilting marker introducer 6 F MAK-NV™
10-0197	Certification operation (AE) in the transfer process, tilting the marker introducer 6 F MAK-NV™ Structures of 87 South Jordan, in West Jordan

Conclusion: Based on the data evaluated in the certification protocols, products MAK-NV™ meet performance specifications.

10. Stability and durability

Merit determines shelf life based on accelerated aging studies performed according to the Arrhenius model for both the packaging and the product. The Arrhenius model $[Q_{10} = 2^{\frac{T_2 - T_1}{10}}]$ quantifies the doubling of reaction rate for every 10°C increase in temperature; providing the number of weeks of aging needed for a material under specified conditions. It is widely accepted as a standard in industry.

The Arrhenius model relates the variables of time and temperature to the deterioration of materials and therefore is useful in calculating accelerated thermal aging parameters. It is preferred for this purpose over the n-degree rule in that it has a firmer theoretical foundation and has been verified for many materials. (reference: Electric Power Research Institute EPRI NP-1558; Project 890-01; Final Report September 1980. "A Review of Equipment Aging Theory and Technology")

Aging studies are followed by real-time aging of product. Merit's internal procedure, STP0002, is derived from the above standard and is utilized in the aging studies of the majority of Merit's products.

Products mini access MAK-NV™ is labeled for a three (3) year shelf life.

11. DETAILS OF CLINICAL STUDIES

This device requires periodic clinical evaluation.

Products mini access MAK-NV™ used in practice since 2006, with a minimum of adverse effects in patients referred to in the Report on clinical evaluation. Assessment of the risks and complications associated with the use of these products was carried out by experienced doctors in accordance with the ISO 14155 series of standards.

These complaints and claims in the operation and clinical evaluation of this product indicates that the value of the product when used in a clinical setting outweighs the risk of harm to the patient if used for its intended purpose, that the performance requirements have been met, that the products show the declared characteristics and are one of the most modern devices.

12. **Manufacturing Process**

The block diagram showing the production process is shown in Appendix A. All products MAK-NV™ produced in a suitably controlled environment as described in the following standard operating procedures:

- GPS 999.031 to environmental action requirements in the clean room environment
- GPS 999.044 controlled environmental conditions
- Procedure QAP 999-950 environmental monitoring

Validation of processes and equipment used for the production of medical products, in accordance with established protocols, and the results are documented and approved. Records are maintained on validated / certified processes, equipment and personnel; they include the use operating instructions, data and test date of the process, the identification of persons who operate or control operation, and, if necessary, the equipment used. The block diagram shown in Appendix A, provides a visual representation of the production process.

Quality control

The criteria for the release of the finished product

Quality control is conducted at regular intervals during the manufacturing process in accordance with the parameters defined in "Merit" procedures and specifications in the technological device map. (See Appendix A. A block diagram (s) of the production process).

Requirements for the verification of systems MAK-NVTM include the following:

Visual inspection includes:

Pollutants and foreign material.

Missing, incorrectly assembled, incorrect or damaged components.

Casting defects, deformed parts, lack of shading, sinks, discoloration, burrs.

Inadequate seams on the package (ie. E. Cavity and excessive wrinkling, etc.).

The bends the tube, the cut (s), dent (s), protuberance (s) on the tube fading.

Damage to the tip, incomplete assembly, an air pocket, sink, hot spots, dirt, line leakage, etc. bulge. (Including marker tip).

Cracked, broken, deformed, or otherwise damaged components or parts.

Check the nut to make sure it fits snugly and is free to rotate.

Check the information printed on the hub in order to avoid the lack of, irregularity, promiscuity or non-printed information.

Check sizes include:

Check the critical dimensions of the molded components.

Tight seal components.

The length of the metal cannula.

The length of the stylet.

Measure the length of protrusion of the needle of the stylet. Do not use the product if this length is outside the specification according to the terms of reference to the drawings.

The distance from the edge tip marker bands to the catheter tip is within the specification

Functional test includes:

The tensile strength with a compound of the stylet sleeve and the tube with the sleeve.
packet burst test.

Leak Testing.

Components / assemblies, which do not meet the requirements specified in the procedures for the revision removed or destroyed.

Upon completion of the production of finished products sent for packaging and sterilization. The sterilized product was placed in a quarantine area after sterilization until verified records sterilization. The product is released for distribution, signed by specialists to ensure quality.

In accordance with the production process and Certificate No. FM 534441 on the conformity of the quality management system of the company "Merit Medical Systems, Inc.", USA (Merit Medical Systems, Inc.), in the field of design, development, production of medical devices to the applied ISO 13485: 2016 standard "Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes", some of the items included in the package are manufactured and packaged at the production site - Merit Maquiladora Mexico, S. DE RL DE CV, Avenida Sor Juana Ines de la Cruz 19970 Interior B, Edificio 2, Parque Industrial Frontera Tijuana, Baja California, MEXICO 22630. Manufacturing of other products and final assembly and packaging of finished products is carried out at the production site: Merit Medical Systems, Inc., 1600 West Merit Parkway South Jordan, Utah 84095 USA. The country of origin is indicated on the label.

13. DESCRIPTION OF PACKAGING

For products "Merit Medical Systems, Inc." Permissible standard warehousing, if the package is not given special storage conditions, in accordance with the risk management reports. The product is packaged and labeled as a sterile disposable article. Sterile barrier systems specified in the design are designed to prevent the ingress of microorganisms and provide an

aseptic preservation of products to end use. Checking the integrity of the packaging is carried out within the production process.

The non-vascular access kit MAK-NV™ placed in a package of the Tyvek® material, the printing process using a validated thermal welding, intended to cover and protect the product during sterilization time of loading, transportation, handling and storage. On the package label applied to the total, which shows the kit.

Needle with stylet wires and coaxial introducer packed in individual packing as a package of Tyvek® material sealable by means of a validated thermal welding process.

The non-vascular access kit MAK-NV™:

1. The non-vascular access kit MAK-NV™, with a guide needle with a trocar stylet, diameter: 21 G (0.8 mm), consisting of:

- Introducer needle with a trocar stylet, diameter: 21 G (0.8 mm), length: 15 cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1.33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018" (0.46 mm), length 60 cm - 1 pc.(if necessary);
- guidewire, diameter 0,038" (0.97 mm), length: 80 cm, 150 cm - 1 pc. (if necessary);
- Instructions for use.

2. The non-vascular access kit MAK-NV™, with a guide needle with a Chiba stylet, diameter: 21 G (0.8 mm), consisting of:

- Introducer needle with Chiba stylet, diameter: 21 G (0.8 mm), length: 15cm, 20 cm;
- Introducer coaxial in composition, diameter 6F (2.0 mm) (if necessary):
 - dilator, diameter 4F (1.33 mm) - 1 pc.,
 - hard metal cannula - 1 pc.,
- guidewire, diameter 0,018" (0.46 mm), length 60 cm - 1 pc.(if necessary);
- guidewire, diameter 0,038" (0.97 mm), length: 80 cm, 150 cm - 1 pc. (if necessary);
- Instructions for use.

Instructions for use can be embedded in packaging group sets. Sets for percutaneous extravascular mini-access collective package is placed in a shipping box.

Certification of product packaging MAK-NV™ of "Merit" It was conducted in accordance with the requirements of the following standards:

ISO 11607-1: packaging of medical devices to be sterilized at the final stage of production - Part 1: Requirements for materials, sterile protection systems and packaging systems.

ISO 11607-2: packaging of medical devices to be sterilized in the final stages of production - Part 2: Requirements for validation processes forming, sealing and assembly.

ASTM D4169: Standard method of testing the efficiency of transport of containers and systems.

ISO 2233 "Packaging - Complete, filled transport packages and unit loads - Conditioning for testing".


14. INFORMATION ON LABELING OF PRODUCTS

Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

The company "Merit Medical Systems, Inc." Uses the following standards for the maintenance and formation of marking symbols in accordance with the current requirements of the Medical Devices Directive:

- EN 1041: Information supplied by the manufacturer.
- ISO 15223-1: medical products. Symbols used on labels for labeling of medical devices and in the information provided - Part 1: General requirements.
- ISO 15223-2: medical products. The symbols used on medical devices, labeling and the information provided - Part 2: Design, selection and validation of characters.
- FDA guidance on the labeling of medical devices for patients; final guidance for industry and FDA reviews (19 April 2001)

Description of the marking symbol is represented in the table.

	Symbol	Decoding the symbol
1.		Authorized Representative in the European Community
2.		Manufacturer
3.		It sold only to doctors and hospitals
4.		A ban on re-use. Do not re-sterilize
5.		Number REF.
6.		use up
7.		lot number
8.		sterilized using ethylene oxide
9.		Attention! Refer to the instructions for use
10.		pyrogen-free
11.		Compatibility with the conductor
12.		Do not use if packaging is damaged
13.		Protect from moisture
14.		This side up
15.		Fragile, handle with care
16.		date of manufacture

17.		Conformity mark
18.		Thin-walled needle tube lettering
19.		The inscription "non-toxic"

Examples of marking are shown below:

Fig. 3 The non-vascular access kit MAK-NVTM needle with stylet wire, Chiba, diameter: 21 G (0,8 mm) length of 15 cm


Fig. 4 The non-vascular access kit MAK-NV™ needle wires with stylet, trocar, diameter: 21 G (0,8 mm) length of 15 cm

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Fig. 5 An example of the layout marking on a medical device, The non-vascular access kit MAK-NV™ needle wire with stylet, trocar, Diameter: 21 G (0,8 mm) with a length of 15 cm to refer to the territory of the Russian Federation

Fig. 5 EXAMPLE

transport packaging

15. TRANSPORTATION, STORAGE AND HANDLING

The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures. Rotate inventory so that device is used prior to the Use By date on the package label. Products released from the transport packaging should be kept at an ambient condition:

Storage conditions:

Temperature +5°C to +27 °C (environment)

Relative humidity 10% to 90 % (non-condensing)

Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Transport conditions:

Temperature -35°C to +55 °C (environment)
Relative humidity 10% to 90 % (non-condensing)
Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Usage environment:

Ambient temperature not exceeding 25 °C.
Temperature range for use - the body temperature of 32 °C to 42 °C
Relative humidity: 100% (blood contact)
Atmospheric pressure 86-106 kPa.

The shelf life of sterile devices - 36 months from the date of sterilization

16. DISPOSAL

To prevent the formation of health risks when disposal of medical devices and also formation of environment risks connected with disposal of medical devices, the personnel have to dispose MD per protocol adopted in particular organization.

Medical Device "The non-vascular access kit MAK-NV™" according to the Sanitary Rules and Regulations 2.1.7.2790-10 is related to Class B contaminated and potentially pathological wastes (epidemiologically hazardous wastes).

17. WARRANTY

Limited Warranty and Disclaimer

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device The non-vascular access kit MAK-NV™ is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

18. Complaint

The organization, which takes claims from consumers about the product quality on the Russian Federation territory:

Limited Liability Company "Merit Technologies"
(LLC "Merit Technologies"),

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6
Tel / fax: +7 (495) 221-89-02

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GPS994.033
Sub Molding
Molded
Subassembly or Finished good
Subassembly
Introducer or Dilator/Stiffner

Introducer
Hub Print/Curing
MPS780.031
Sub Molding
Tipping
MPS780.029
MPS780.033
MPS780.061
In process inspection
MPS780.033
MPS780.061
Fail
Scrap
Packaging Subs
MPS795.023
QA inspection
TBL243
Satisfactory inspection?

No
Generate NC
GPS994.003
Yes
Send materia to kanban
GPS994.033
Sub Assy
Dilator/Stiffner
Sub Molding
Nut Thumping
MPS780.032
Tipping
MPS780.029
MPS780.033
MPS780.061
In process inspection
MPS780.033
MPS780.061
Packaging Subs
MPS795.023
QA inspection
TBL243
Satisfactory inspection?

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MPS999.270
Sub Assy
Inner Seal
MPS999.270
MPS999.203
Apply Label
MPS999.270
MPS994.062
Additional Components
MPS999.270
Outer Seal
MPS999.270
MPS999.203
Boxing
MPS994.062
MPS999.271
QA inspection
TBL243/QAP994.012
Satisfactory inspection?

No
Generate NC
GPS994.003
Yes
Shipping
GPS994.013
Needs sterilization?
BowTie
Package BowTie in TJ
MPS999.818
(Only MAK / SMAK-BT)
Guidewire from GAL
Needle from SLC
HVA Valve from SLC (Only PAK)
Retort
QSP994.004
Pre Sterilizer
QSP994.004

Sterilization
QSP994.004
Post Sterilizer
QSP994.037
Locate
GPS994.009
End
Symbols Legend
QA inspection
Operation
Decision
Scrap/Disposition
Generate Docs

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