BAUSCH+LOMB

See better. Live better.

1400 North Goodman Street
Rochester, NY 14609
585.338.6000
www.bausch.com

Certified True Copy Letter

I, Isabelle B. Lefebvre, hereby swear (or affirm) that the attached Instruction for Use Of Medical Device for Akreos Single Use Insertion Device is true and correct.

Isabelle B. Lefebvre

Isabelle B. Lefebvre, MSc.RA
Vice President, Regulatory Affairs
Bausch & Lomb Incorporated

Date 01 July 2017

State of New Jersey

County of Somerset

Subscribed and sworn to (or affirmed) before me on this 1 day of June, 2017, by ISABELLE B. LEFEBVRE, proved to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Christine A. Hnatko

Signature of Notary Public

CHRISTINE A. HNATKO
Commission # 50001786
Notary Public, State of New Jersey
My Commission Expires
August 13, 2019

INSTRUCTION FOR USE OF MEDICAL DEVICE

Akreos Single Use Insertion Device

Design facility:

Bausch & Lomb Incorporated, USA

1400 North Goodman Street Rochester, NY 14609, USA

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Bausch & Lomb Incorporated, USA

1400 North Goodman Street Rochester, NY 14609, USA

Authorized Representative of Manufacturer:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

The legal entity addressed to whom the registration certificate can be issued

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Addresses of manufacturing facility of medical device:

1. 1400 North Goodman Street Rochester, NY 14609, USA
1. 21 Park Place Blvd. N. Clearwater, FL 33759, USA

DESCRIPTION AND THE DESIGN OF MEDICAL DEVICE

The Akreos Single Use Insertion device system consists of syringe shaped tube (includes lens loading deck and lens cap), plunger, and transition cell. The syringe tube and transition cell are made of polypropylene material. The plunger is made of polyphthalamide material. The Akreos Insertion Device is sterilized by ethylene oxide (EO).

The Akreos Single-Use Insertion Device is used for folding and delivering the Bausch + Lomb Akreos Advanced Optics Aspheric IOL into the eye during cataract surgery. The inserter consists of a syringe-shaped body a transition cell and is molded out of polypropylene resin. The primary components are the transition cell, the body, insert, and plunger. The Akreos Single-Use Insertion Device is a sterile, disposable, single-use device.


Fig. 1. Akreos Single Use Insertion Device. Photo.

MATERIALS USED IN MEDICAL DEVICE CONSTRUCTION AND ITS COMPONENTS

Akreos Single Use Insertion Device (and therefore each part of Akreos Single Use Insertion Device) is sterile and sterilized by Ethylene Oxide (SAL (Sterility Assurance Level) is 10^{-6}).

All the parts of Akreos Single Use Insertion Device have direct or indirect contact with human body (internal damaged tissues or/ and organs) by viscoelastic or balanced salt solution.

Tube/Body, Transition cell, Plunger: Basell PF531 homopolymer compounded with Pationic 901 additive

Plunger: Amodel A1240L NT Polyphthalamide

O-ring:

70 durometer silicone, red colored (the mark and composition of the material is highly confidential, this information is confirmed by the Letter of Confidentiality from the Manufacturer)

Spring:

Stainless Steel 302:

Fe	base
C	0.15 % Max
Cr	17 – 19 %
Ni	8 – 10 %
Si	1.0 % Max
Mo	1.0 % Max
Mn	2.0 % Max
P	0.045 % Max
S	0.03 % Max)

The blister package has only direct contact with the skin of medical personnel if the medical personnel works without gloves:

The lid: Tyvek ®

The tray: PETG Eastar Copolyester 6763

The conjunction of the lid and tray is thermos

The ink on the lid of the package: the mark and composition of the material is highly confidential, this information is confirmed by the Letter of Confidentiality from the Manufacturer

THE MAIN GUARANTEED CHARACTERISTICS OF MEDICAL PRODUCT

The main overall dimensions of the medical device (of its main parts) are introduced in the drawings and extracts of the drawings below, all the dimensions are introduced in millimeters:


Fig. 2. Body and transition deck
(extract from the Manufacturer drawing, the dimension are introduced in millimeters)


Fig. 3. Plunger

(extract from the Manufacturer drawing, the dimension are introduced in millimeters)


Fig. 4. Plunger

(extract from the Manufacturer drawing, the dimension are introduced in millimeters)

The overall mass of Akreos Single Use Insertion Device is 5,5 gram $\pm$ 1 gram.

The minimal force needed for the Plunger shifting during the medical device use is 20,0 – 40,0 N.

The period of validity of Akreos Single Use Insertion Device is 1 year (12 months) from date of manufacture (from the month of sterilization) if storage and transportations conditions are observed.

INDICATIONS

The Akreos™ Single Use Insertion Device is used for folding and delivering the Bausch & Lomb one-piece Akreos™ Advanced Optics Aspheric Lens into the eye. The device provides a small tubular pathway in which the Akreos Advanced Optics Aspheric Lens can be placed into the eye with one continuous forward motion.

The Akreos Single Use Insertion Device should only be used with the Bausch & Lomb Akreos Advanced Optics Aspheric Lens from 0 to 30 diopters in power.

CONTRAINDICATIONS

There are no contraindications for the Akreos Single Use Insertion Device.

ADVERSE EVENTS

1. IOL haptic and/or optic damage.
2. Anaphylactic Shock.
3. Microbial Infection.
4. Posterior Capsule Rupture.
5. Corneal Endothelium Damage.
6. Zonules, Torn, Weakened.
7. Elevated Intraocular Pressure.
8. Enlarged Incision Size/ Inadvertent Tissue Damage.
9. Surgical Procedure Impact, Delay (<30 min).

Adverse reaction reporting

Adverse reactions and/or potentially sight-threatening complications that maybe reasonably regarded as Akreos Single Use Insertion Device related and that were not previously expected in nature, severity or degree of incidence should be reported to your customer service affiliate.

INSTRUCTIONS FOR USE

1. Open the package using standard sterile procedures. Place contents onto the sterile field.
2. Fill the rectangular end of the transition cell with approximately 0.2 mL of a Bausch & Lomb viscoelastic, using the line as a guide (**Fig. 5**).


Fig. 5.

3. Select the syringe-shaped tube and plunger. Pull the plunger until the tip of the plunger is just visible in the lens loading deck. Place a drop of viscoelastic directly in the center of the loading deck. Open the vial and remove the lens holder. Position the lens holder with its jaws facing down (**Fig. 6a**). With a blunt, non-toothed forceps, remove the lens from the lens holder by pushing the lens down. Transfer the lens to the recessed area in the injector loading deck (**Fig. 6a**), ensuring that the lens is not flipped during the transfer. The lens is loaded in correct orientation if its haptic features are as shown in **Fig. 6b**.


Fig. 6a.


Fig. 6b.

4. Advance the plunger until the fork fully captures the lens edge. Place a drop of viscoelastic in the center of the lens. Verify that the lens is inside of the recessed area of the loading deck, prior to closing the lens cap.
5. Close the lens cap (**Fig. 7**).


Fig. 7.

6. Align the key on the lens cap with the slot in the transition cell (**Fig. 8**). Push together with a firm motion until the tube with the cap and transition cell snap together (**Fig. 9**). Viscoelastic material will flow and fill both the lens chamber and the transition cell. **Once assembled, pull on the transition cell to verify that it cannot be separated from the tube.**


Fig. 8.


Fig. 9.

7. Fill the distal end of the device with viscoelastic material or balanced salt solution to reduce the possibility of introducing air pockets into the eye during implantation.
8. The lens is now ready to be delivered. Turn the device over, so that the bevel is face down, and insert the beveled portion of the transition cell through the incision in the eye. When the tip is positioned, apply continuous uniform forward motion to the plunger. A slow smooth forward motion of the plunger is essential for proper transition of the lens into position. Reverse movement of the plunger may result in the loss of the control of the lens or damage to the lens. Continue forward motion until the lens is fully expressed from the tip.
9. The edge of the unfolded lens may remain engaged in the plunger fork until it is released into the eye. To disengage the lens, release pressure on the plunger end allowing it to retract into the tip.
10. Discard the device after use.

OPERATING CONDITIONS

The device must only be used at adequate facilities by a healthcare professional after training on this device.

Ambient Operating conditions	Temperature: 21 – 24 (21) °C Humidity: 30–60 % (30% to 60% Relative) The velocity of air motion: not more than 0,2 meters per second
------------------------------	--

STORAGE AND TRANSPORTATION CONDITIONS

Store and transport the Akreos Single Use Insertion Device sterile sealed package at the following conditions:

Ambient Storage/Transport conditions	Temperature: +15 – +30 °C Humidity: 30–70 % (30% to 70% Relative)
--------------------------------------	--

SAFETY REQUIREMENTS (INCLUDING SAFETY METHODS AND TOOLS)

Store and transport the Akreos Single Use Insertion Device sterile sealed package at the following conditions:

Ambient Storage/Transport conditions	Temperature: +15 – +30 °C Humidity: 30–70 % (30% to 70% Relative)
--------------------------------------	--

Contents are sterile unless package is opened or damaged. Do not resterilize. Do not reuse. The device is for Single use only. Discontinue lens delivery if excessive force is needed to advance the plunger and lens.

CAUTION: Federal Law (USA) restricts this device to sale by, or on the order of, a physician.

SURGICAL PROCEDURE

Proper surgical technique is the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based upon their medical training and experience.

WARRANTY

Bausch & Lomb Incorporated warrants that the Akreos Single Use Insertion Device, when delivered, will conform to the manufacturer's then current version of the published specifications for such Akreos Single Use Insertion Device in all material respects and shall be free from defects in material or workmanship for a period equal to the Akreos Single Use Insertion Device expiration date.

Bausch & Lomb Incorporated warrants that the Akreos Single Use Insertion Device, when delivered, will conform to the manufacturer's then current version of the published specifications for such Akreos Single Use Insertion Device in all material respects and shall be free from defects in material or workmanship for a period equal to the Akreos Single Use Insertion Device expiration date.

Bausch & Lomb Incorporated excludes all other warranties, whether expressed, implied or by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for a particular purpose. Bausch & Lomb Incorporated shall not be liable for any incidental, consequential or exemplary loss, damage or expense, directly or indirectly arising from the use of the Akreos Single Use Insertion Device even if Bausch & Lomb Incorporated has been advised of the possibility of such loss, damage or expense.

RETURN POLICY

Bausch + Lomb will, within the return period as specified on your invoice, from the date of invoice, accept return of this product for a full refund less any handling and shipping charges incurred by Bausch + Lomb. Customer must call their local Bausch + Lomb customer service representative to request a Return Good Authorization prior to expiration of the return period. It is the Customer's responsibility to properly pack all items being returned.

Please call +1-800-338-2020 for the authorized product return and detailed information obtaining.
The Authorized Representative of Bausch + Lomb in Russia:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Tel.: +7 495 510 2879, fax: +7 495 510 2879

Office.RU@valeant.com

ADVERSE REACTION REPORTING

Adverse reactions and/or potentially sight-threatening complications that maybe reasonably regarded as Akreos Single Use Insertion Device related and that were not previously expected in nature, severity or degree of incidence should be reported to your customer service affiliate.

The Authorized Representative of Bausch + Lomb in Russia:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Tel.: +7 495 510 2879, fax: +7 495 510 2879

Office.RU@valeant.com

SEVICE INFORMATION

Bausch & Lomb Incorporated warrants that the Akreos Single Use Insertion Device, when delivered, will conform to the manufacturer's then current version of the published specifications for such Akreos Single Use Insertion Device in all material respects and shall be free from defects in material or workmanship for a period equal to the Akreos Single Use Insertion Device expiration date.

REQUIREMENTS FOR THE RECYCLING (THE PROTECTION OF THE ENVIRONMENT IN THE APPLICATION OF MEDICAL PRODUCTS)

Akreos Single Use Insertion Device after its using is classified as a class B of waste (according Russian Legislation) and should be disposed according to SanPin 2.1.2790-10 ("Sanitary rules and norms #2.1.7.2790-10").

PACKAGING AND LABELING INFORMATION

On the package and label of Akreos Single Use Insertion Device the following symbols are used:


Fig. 8. The symbols used on the package and label.

Each Akreos Single Use Insertion Device in individual blister package is packaged in carton package.

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