Cover Page for Instructions for Use documentation, the following documentation is provided:

- IFU in English (DEV-01731)
- IFU in Russian translated from the English IFU (DEV-01768)

Handwritten signature

Handwritten signature of Søren Olesen

Signature

Handwritten date: 2021-07-08

Date

Søren Olesen,
Global Business Director

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Document Title **IFU english for Russian registration**
Document ID **DEV-01731**
Version **7.0**
Release Date **07-Jul-2021**
Lifecycle State **Approved**
Internal Legacy Number

Content Approval

The Electronic signature of Content Approval means that this document has been reviewed and approved for compliance with the end user's and owner's requirements.

Signed by **ELOSMEDTECH\lazfar.rehman**

Decision date **07-jul-2021 10:17**

QA Authorisation

The Electronic Signature of QA Authorisation means that this document has been reviewed and approved for compliance with business quality standards and SOPs.

Signed by **ELOSMEDTECH\lazfar.rehman**

Decision date **07-jul-2021 10:25**

Torque Wrench & Driver

Instruction for use

Medtech Pinol A/S
Engvej 33
3330 Gørlese, Denmark
www.medtech.com

Batch code

Catalogue number

Manufacturer

Date of manufacture

By prescription only

Consult instructions for use

Do not use if package is damaged

Medical device

Unique Device Identifier

Non sterile

Keep dry

Keep away from sunlight

Russian certification symbol according to GOST

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
Controlled Copy Location Printed by Nise Souto (ELOS MEDTECH/nise.souto) on 07-Jul-2021

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Electronic IFU can be accessed at <http://medtech.com/IFU/>
Document ID: DEV-01731 version 7.0 DRAFT (based on DEV-00275 version 7.0) Released yyyy/mm/dd

Medtech Pinol A/S
Engvej 33
3330 Gårlev, Denmark
www.medtech.com

Instruction for use - English

Torque Wrench EV & Driver (35, 45, 70 & 100Ncm)

1 Intended use

A dental torque wrench for placement and adjustment of dental implants, abutments, abutment screw and prosthetic screws during oral surgery and prosthetic procedures. By use of exchangeable screwdriver tips and bits the tool can be used for Astra Tech Implant System.

2 Product Description

The Torque Wrench (main body, head and Driver) is made of lightweight titanium and stainless steel. The Torque Wrench & Driver product design, finish and quality make the tool well suited for placement and adjustment of implants, abutments and screws during surgery and other clinical procedures. The Torque Wrench & Driver is easy to prepare for cleaning, disinfection and sterilization.

Exchangeable screwdrivers and the Driver make the tool compatible to various implant systems. The driver is available for a variety of connections. For specific product descriptions, please refer to individual product labels.

3 Use

The Torque Wrench & Driver is delivered non-sterile. Prior to use it must be cleaned with water and brush followed by sterilization according to section 5. The use of the Torque Wrench & Driver is described below.

1. After cleaning make sure the first line on the scale of the Torque Wrench arm aligns with the arrow (see figure 1.A).
2. Attach the head to the body of the Torque Wrench by pushing the components together and turning them in opposite direction until an audible click is achieved (see figure 2).
3. The arrow on the head of the Torque Wrench shows the direction in which the Torque Wrench & Driver is functioning, and the head is labeled "IN" (see figure 3.A) and "OUT" (see figure 3.B). "IN" denotes clockwise rotation; "OUT" denotes counterclockwise rotation. In order to change direction, flip the Torque Wrench to the opposite side (see figure 3) and push the Driver into the head of the Torque Wrench until an audible click is achieved (see figure 4).
4. Insert the tool into the Driver.
5. Place the tool in the implant/screw and turn the Torque Wrench arm in the direction of the arrow until the desired torque is achieved (see figure 6). The applied torque must be in accordance with the screw/implant manufacturers' instructions.
6. After each use clean and sterilize the Torque Wrench & Driver according to section 5.
7. Inspect instruments under normal lighting for the removal of visible soil and for damage or wear. Any instrument showing visible damage or wear must be disposed of according to section 10.

4 Warnings and precautions

Potential adverse occurrences may include but are not limited to the following list. Proper placement and use of this instrument may help to mitigate these risks.

• If overloaded or dropped

If overloaded, dropped or in other way mishandled, the product must no longer be used since correct functionality can no longer be guaranteed.

• Do not use hydrogen peroxide

This product must be sterilized prior to use and it must not be cleaned by using hydrogen peroxide.

• Scale alignment

The Torque Wrench & Driver does not require recalibration but before each use, make sure that the Torque Wrench & Driver is fully functional and that the first line on the scale aligns with the arrow (see figure 1.A).

• Ensure correct assembly

Before each use, make sure that all parts are correct assembled (no loose parts are allowed).

• Do not exceed max scale

The arm of the Torque Wrench must not go beyond the end of the scale (see figure 6.B), as this could result in permanent deformation of the Torque Wrench arm i.e. future inaccurate torque readings (see figure 1.B).

• No scale indication when used as ratchet wrench

If the Torque Wrench & Driver is used as a ratchet wrench it does NOT provide the user any indication of what torque that is used.

Notice: the torque applied must not exceed 150Ncm when the Torque Wrench & Driver is used as a ratchet wrench.

• Use screw/screwdriver/Driver/implant manufacturer torque recommendations

Please check the individual screw/screwdriver/ Driver/implant manufacturers' instruction for use or label regarding the maximum allowed torque. Do not under any circumstances apply higher torque than stated even if the Torque Wrench & Driver is used as a ratchet wrench (without using the Torque Wrench arm).

• Cleaning – do not use metal brush or steel wool

Do not use a metal brush or steel wool to clean the Torque Wrench or Driver as this will damage the instrument.

• Only for use by professionals

The Torque Wrench & Driver must only be used by professionals.

5 Cleaning and sterilization

The following cleaning and disinfection method are recommended. Further instruction can be found in "Cleaning and sterilization guideline" at <https://elosmedtech.com/IFU/>

1. Disassemble the Torque Wrench (see figure 5)
2. The separated parts are now ready to be cleaned according to the "Cleaning and sterilization guideline".
3. Sterilization according to the "Cleaning and sterilization guideline". Summary below:

Note: For US: Steam Sterilization Cycle 132°C (270°F)/4min.

Procedure	Dynamic-air-removal sterilization Cycle	
Exposure time	4 min.	3 min.
Temperature	132°C (270°F)	134° C (273°F)
Drying time	20 min.	

Caution: Accompanying packaging material CANNOT be autoclaved unless it is explicitly stated that the material is suitable for sterilization.

6 Further information

For additional information about the use of the Torque Wrench & Driver, please contact your local sales representative.

7 Validity

Upon publication of this instruction for use, all previous versions are superseded.

8 Storage and Handling

Recommended storage temperature range is 10-40°C (50-104°F) and the recommended relative humidity range is 0-80%. The shelf life is maximum 10 years.

9 Transportation conditions:

Products should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device's nature at a temperature range from -50 to +50°C, with a humidity up to 80 %.

10 Disposal

When disposed of The Torque Wrench & Driver may present an infection or microbial hazard and should therefore be disposed of as biological waste, according to national requirements.

11 Serious Incidents:

If any serious incidents occur in relation with this device, they should be reported to Elos Medtech Pinol A/S and to the competent authorities in the country where the incident occurred. A serious incident is a situation where the device directly or indirectly led or might have led to death, a serious deterioration in the health of a person or a serious public health threat. Contact Elos Medtech Pinol A/S at: complaintempi@medtech.com

Caution: U.S. Federal Law restricts this device to sale by or on the order of a licensed dentist or physician.

12 Warranty:

Elos warrants that Products are manufactured in accordance with ISO 13485 and specifications agreed in writing and that Products will be free from defects in workmanship for a period of 2 years for instruments from the date of Elos' delivery.

13 Authorized representative in Russia

Authorized representative of the manufacturer in the territory of the Russian Federation:

LLC "Dentsply Sirona", Russian Federation, 115432 Moscow, Andropov, Prospekt 18, Building 6 tel.: +7 (495) 725 1087

Medtech Pinol A/S
Engvej 33
3330 Gørlese, Denmark
www.medtech.com

Document Title IFU Russian for the Russian registration
Document ID DEV-01768
Version 6.0
Release Date 06-Jul-2021
Lifecycle State Approved
Internal Legacy Number

Content Approval

The Electronic signature of Content Approval means that this document has been reviewed and approved for compliance with the end user's and owner's requirements.

Signed by ELOSMEDTECH\lazfar.rehman

Decision date 06-jul-2021 18:29

QA Authorisation

The Electronic Signature of QA Authorisation means that this document has been reviewed and approved for compliance with business quality standards and SOPs.

Signed by ELOSMEDTECH\lazfar.rehman

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Controlled Copy Location Printed by Nise Souto (ELOS MEDTECH/nise_souto) on 07-Jul-2021

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